



DICOM Conformance Statement for
LISENDO 880

Company Name: FUJIFILM Corporation

Product Name: Diagnostic Ultrasound System LISENDO 880

Version: Version 6.3.0

Date: 2024.07.01

Internal document number: SFT-116833

(C) FUJIFILM Corporation 2024.

The equipment described in this document has been validated in accordance with the DICOM standard. FUJIFILM Corporation shall not be liable for errors contained herein or consequential damages in connection with the furnishing, performance or use of this document. FUJIFILM Corporation reserves the right to revise this document and to make changes to its contents at any time, without obligation to notify any person or entity of such revisions and changes.

1. CONFORMANCE STATEMENT OVERVIEW

This product LISEND0 880 implements the necessary DICOM services to download worklists from an information system, save acquired Ultrasound and Ultrasound Multi-frame images to a network storage device or storage media, print Ultrasound images to a networked hardcopy device and inform the information system about the work actually done.

Table1-1 provides an overview of the network services supported by LISEND0 880.

Table 1-1 NETWORK SERVICES

SOP Classes	User of Service (SCU)	Provider of Service (SCP)
Transfer		
Ultrasound Image Storage	Yes	No
Ultrasound Multi-frame Image Storage	Yes	No
Comprehensive SR ¹	Yes	No
Encapsulated PDF Storage ²	Yes	No
Secondary Capture Image Storage ³	Yes	No
Multi-frame True Color Secondary Capture Image Storage ³	Yes	No
Workflow Management		
Modality Worklist	Yes	No
Modality Performed Procedure Step	Yes	No
Storage Commitment Push Model	Yes	No
Verification	Yes	Yes
Print Management		
Basic Grayscale Print Management Meta	Yes	No
Basic Color Print Management Meta	Yes	No
Query/Retrieve ⁴		
Study Root Information Model FIND	Yes	No
Study Root Information Model MOVE	Yes	No
Ultrasound Image Storage	No	Yes
Ultrasound Multi-frame Image Storage	No	Yes
CT Image Storage	No	Yes
MR Image Storage	No	Yes
Positron Emission Tomography Image Storage	No	Yes
Encapsulated PDF Storage	No	Yes

Note: 1. SOP-LISEND0880-21 is necessary.

Note: 2. SOP-LISEND0880-70 is necessary to create Encapsulated PDF report.

Note: 3. SOP-LISEND0880-112, SOP-LISEND0880-113, or SOP-LISEND0880-114 is necessary to create Secondary Capture Image.

Note: 4. SOP-LISENDO880-59 is necessary.

Table 1-2 provides an overview of the Media Storage Application Profiles supported by LISENDO 880.

Table 1-2 Media Services

Media Storage Application Profile	Write Files (FSC or FSU)	Read Files (FSR) ¹
DVD-RAM & Compact Disk – Recordable		
Ultrasound Image Display	Yes	Yes
Ultrasound Image Spatial Calibration	Yes	Yes ²
Ultrasound Image Combined Calibration	No	No
Ultrasound Multi-frame Image Display	Yes	Yes
Ultrasound Multi-frame Image Spatial Calibration	Yes	Yes
Ultrasound Multi-frame Image Combined Calibration	No	No

Note: 1. Structured Reports, SC images, and SC-mf images cannot be imported.

Note: 2. Not support Ultrasound images from other products.

2. TABLE OF CONTENTS

1. CONFORMANCE STATEMENT OVERVIEW.....	2
2. TABLE OF CONTENTS	4
3. INTRODUCTION	7
3.1 REVISION HISTORY	7
3.2 AUDIENCE	7
3.3 REMARKS.....	7
3.4 DEFINITIONS, TERMS AND ABBREVIATIONS	8
3.5 REFERENCES	9
4. NETWORKING.....	10
4.1 IMPLEMENTATION MODEL	10
4.1.1 <i>Application Data Flow</i>	10
4.1.2 <i>Functional Definition of AEs</i>	11
4.1.2.1 Functional Definition of Workflow Application Entity	11
4.1.2.2 Functional Definition of Storage Application Entity	11
4.1.2.3 Functional Definition of Hardcopy Application Entity	12
4.1.2.4 Functional Definition of Query/Retrieve Application Entity	12
4.1.3 <i>Sequencing of Real-World Activities</i>	13
4.2 AE SPECIFICATIONS	14
4.2.1 <i>Workflow Application Entity Specification</i>	14
4.2.1.1 SOP Classes	14
4.2.1.2 Association Policies.....	14
4.2.1.2.1 General.....	14
4.2.1.2.2 Number of Associations	15
4.2.1.2.3 Asynchronous Nature.....	15
4.2.1.2.4 Implementation Identifying Information	15
4.2.1.3 Association Initiation Policy.....	15
4.2.1.3.1 Activity – Find Worklist.....	15
4.2.1.3.2 Activity – End Exam.....	20
4.2.1.4 Association Acceptance Policy	27
4.2.1.4.1 Activity – Receive Storage Commitment Response.....	27
4.2.2 <i>Storage Application Entity Specification</i>	30
4.2.2.1 SOP Classes	30
4.2.2.2 Association Policies.....	30
4.2.2.2.1 General.....	30
4.2.2.2.2 Number of Associations	30
4.2.2.2.3 Asynchronous Nature.....	30
4.2.2.2.4 Implementation Identifying Information	30
4.2.2.3 Association Initiation Policy.....	31
4.2.2.3.1 Activity – Send Images, Structured Reports, and Encapsulated PDF reports	31
4.2.3 <i>Hardcopy Application Entity Specification</i>	37
4.2.3.1 SOP Classes	37
4.2.3.2 Association Policies.....	37
4.2.3.2.1 General.....	37
4.2.3.2.2 Number of Associations	37
4.2.3.2.3 Asynchronous Nature.....	37
4.2.3.2.4 Implementation Identifying Information	37
4.2.3.3 Association Initiation Policy.....	38

4.2.3.3.1	Activity – Print Images	38
4.2.3.4	Association Acceptance Policy	46
4.2.4	<i>Query/Retrieve Application Entity Specification</i>	47
4.2.4.1	SOP Classes	47
4.2.4.2	Association Policies	47
4.2.4.2.1	General	47
4.2.4.2.2	Number of Associations	47
4.2.4.2.3	Asynchronous Nature	47
4.2.4.2.4	Implementation Identifying Information	48
4.2.4.3	Association Initiation Policy	48
4.2.4.3.1	Activity – Query Image Information and Retrieve Images	48
4.3	<i>NETWORK INTERFACES</i>	55
4.3.1	<i>Physical Network Interface</i>	55
4.3.2	<i>Additional Protocols</i>	55
4.3.3	<i>IPv4 and IPv6 Support</i>	55
4.4	<i>CONFIGURATION</i>	56
4.4.1	<i>AE Title/Presentation Address Mapping</i>	56
4.4.1.1	Local AE Titles	56
4.4.1.2	Remote AE Title/Presentation Address Mapping	56
4.4.1.2.1	Workflow	56
4.4.1.2.2	Storage	56
4.4.1.2.3	Hardcopy	57
4.4.1.2.4	Query/Retrieve	57
4.4.2	<i>Parameters</i>	57
5.	MEDIA INTERCHANGE	60
5.1	<i>IMPLEMENTATION MODEL</i>	60
5.1.1	<i>Application Data Flow</i>	60
5.1.2	<i>Functional Definition of AEs</i>	60
5.1.2.1	Functional Definition of Offline-Media Application Entity	60
5.1.3	<i>Sequencing of Real-World Activities</i>	60
5.1.4	<i>File Meta Information Options</i>	60
5.2	<i>AE SPECIFICATIONS</i>	61
5.2.1	<i>Offline-Media Application Entity Specification</i>	61
5.2.1.1	File Meta Information for the Application Entity	61
5.2.1.2	Real-World Activities	61
5.2.1.2.1	Activity – Export to Media	61
5.3	<i>AUGMENTED AND PRIVATE APPLICATION PROFILES</i>	62
5.4	<i>MEDIA CONFIGURATION</i>	62
6.	SUPPORT OF CHARACTER SETS	63
7.	SECURITY	64
7.1	<i>SECURITY PROFILES</i> ¹	64
8.	ANNEXES	65
8.1	<i>IOD CONTENTS</i>	65
8.1.1	<i>Created SOP Instances</i>	65
8.1.1.1	Ultrasound and Ultrasound Multi-frame Image IODs	65
8.1.1.2	Secondary Capture and Multi-frame True Color Secondary Capture Image IODs	66
8.1.1.3	Comprehensive SR IOD	67
8.1.1.4	Encapsulated PDF IOD	68
8.1.1.5	Common Modules	68
8.1.1.6	Ultrasound Modules	76
8.1.1.7	SR Document Modules	86
8.1.1.8	Encapsulated Document Modules	88

8.1.1.9	Secondary Capture Modules.....	89
8.1.2	<i>Used Fields in received IOD by application</i>	90
8.1.3	<i>Attribute mapping</i>	90
8.1.4	<i>Coerced/Modified Fields</i>	91
8.2	DATA DICTIONARY OF PRIVATE ATTRIBUTES	91
8.3	CODED TERMINOLOGY AND TEMPLATES	99
8.4	STANDARD EXTENDED / SPECIALIZED / PRIVATE SOP CLASSES.....	99
8.4.1	<i>Ultrasound and Ultrasound Multi-frame Image Storage SOP Class</i>	99
8.4.2	<i>Comprehensive SR Storage SOP Class</i>	99
8.5	PRIVATE TRANSFER SYNTAXES	99
8.6	STRUCTURED REPORTS	100
8.6.1	<i>Applications and Generated Templates</i>	100
8.6.1.1	<i>Image and Spatial Coordinates</i>	100
8.6.2	<i>Templates</i>	101
8.6.2.1	TID 5000 OB-GYN Ultrasound Procedure Report	101
8.6.2.2	TID 5100 Vascular Ultrasound Report.....	111
8.6.2.3	TID 5200 Echocardiography Procedure Report (Non-Advanced Cardiac)	129
8.6.2.4	TID 5200 Echocardiography Procedure Report (Advanced Cardiac)	208
8.6.3	<i>Context Groups (Non-Advanced Cardiac)</i>	286
8.6.4	<i>Private Code Definitions</i>	332
8.6.5	<i>Standard Extended and Private Template Definitions (Non-Advanced Cardiac)</i>	369

The company names, product names and other names described in this document are trademarks or registered trademarks of FUJIFILM Corporation or its subsidiaries.

ALOKA is a registered trademark or trademark of ALOKA Co., Ltd. in Japan and other countries.

DICOM is a registered trademark of the National Electrical Manufacturers Association in the United States for its standards publications relating to digital communications of medical information.

3. INTRODUCTION

3.1 REVISION HISTORY

Document Version	Date of Issue
1.0	2016.09.30
1.1	2016.10.07
1.2	2016.10.14
1.3	2017.02.17
2.0	2017.09.08
2.1	2018.01.10
2.2	2018.04.20
3.0	2018.10.12
3.2	2019.05.29
4.0	2019.12.17
4.1	2020.12.01
4.1.1	2021.04.01
4.2	2021.10.04
6.3	2024.07.01

3.2 AUDIENCE

This document is intended for hospital staff, health system integrators, software designers or implementers. It is assumed that the reader has a working understanding of DICOM.

3.3 REMARKS

DICOM, by itself, does not guarantee interoperability. However, the Conformance Statement facilitates a first-level validation for interoperability between different applications supporting the same DICOM functionality.

This Conformance Statement is not intended to replace validation with other DICOM equipment to ensure proper exchange of information intended.

The scope of this Conformance Statement is to facilitate communication with LISENDO 880 and other vendors' Medical equipment. The Conformance Statement should be read and understood in conjunction with the DICOM Standard [DICOM]. However, by itself it is not guaranteed to ensure the desired interoperability and a successful interconnectivity.

The user should be aware of the following important issues:

- The comparison of different conformance statements is the first step towards assessing interconnectivity between the equipments produced by different manufacturers.
- Test procedures should be defined to validate the desired level of connectivity.

- The DICOM standard will evolve to meet the users' future requirements. FUJIFILM Corporation reserves the right to make changes to its products or to discontinue its delivery.

3.4 DEFINITIONS, TERMS AND ABBREVIATIONS

Definitions, terms and abbreviations used in this document are defined within the different parts of the DICOM standard.

Abbreviations and terms are as follows:

AE	DICOM Application Entity
AET	Application Entity Title
ACSE	Association Control Service Element
CD-R	Compact Disk Recordable
CID	Context Identifier
CM	Code Meaning
CSD	Coding Scheme Designator
CSE	Customer Service Engineer
CV	Code Value
DVD	A trademark of the DVD Forum that is not an abbreviation
DVD-RAM	DVD Random Access Memory
FSC	File-Set Creator
FSU	File-Set Updater
FSR	File-Set Reader
GUI	Graphical User Interface
HDD	Hard Disk Drive
IHE	Integrating the Healthcare Enterprise
IOD	(DICOM) Information Object Definition
ISO	International Standard Organization
LCD	Liquid Crystal Display
MPPS	Modality Performed Procedure Step
MSPS	Modality Scheduled Procedure Step
NEMA	National Electrical Manufacturers Association
R	Required Key Attribute
O	Optional Key Attribute
PDU	DICOM Protocol Data Unit
PHI	Protected Health Information
SCU	DICOM Service Class User (DICOM client)
SCP	DICOM Service Class Provider (DICOM server)
SOP	DICOM Service-Object Pair

SR	Structured Reporting
TID	Template Identifier
U	Unique Key Attribute
UID	Unique Identifier
USB	Universal Serial Bus
VR	Value Representation

3.5 REFERENCES

[DICOM]Digital Imaging and Communications in Medicine (DICOM), NEMA PS 3.1-3.18, 2013

4. NETWORKING

4.1 IMPLEMENTATION MODEL

4.1.1 Application Data Flow




Figure 4.1-1 APPLICATION DATA FLOW DIAGRAM

- The Workflow Application Entity receives Worklist information from and sends MPPS information to a remote AE. It is associated with the local real-world activities "Find Worklist" and "End Exam". When the "Find Worklist" local real-world activity is requested by an operator the Workflow Application Entity queries a remote AE for worklist items and provides the set of worklist items matching the query

request. When "Acquire Images" local real-world activity is performed on the selected patient for the first time, the Workflow Application Entity automatically creates Modality Performed Procedure Step instance managed by a remote AE. When "End Exam" local real-world activity is requested by the operator, the MPPS Completed or Discontinued updates the MPPS instance managed by the remote AE. If the remote AE is configured as an archive device the Workflow AE will request Storage Commitment and if a commitment is successfully obtained the Workflow AE will record this information in the local database.

- The Storage Application Entity stores images(Ultrasound images, Ultrasound Multi-frame images, Secondary Capture images, Multi-frame True Color Secondary Capture images), Encapsulated PDF Reports, and Structured Reports to a remote AE. It is associated with the local real-world activity "Send Images", "Send PDF Reports", and "Send Structured Reports", respectively. An ultrasound modality displays image in real-time, and the operator acquires it by pressing the freeze button. A "Frozen" image may be sent to a remote AE or may be stored in the local HDD, DVD-RAM, CD-R buffer or USB Media for review and batch send. "Send Images" is performed upon user request for each study completed or for specific images selected. Sending Structured Reports is performed automatically at "End Exam" local real-world activity.
- The Hardcopy Application Entity prints images on a remote AE (Printer). It is associated with the local real-world activity "Print Images". "Print Images" creates a print-session within the print queue containing one or more virtual film sheets composed from images selected by the user.
- The Query/Retrieve Application Entity queries instances information of Ultrasound and/or Ultrasound Multi-frame images or CT Images or MR Images or PET Images or Documents, and retrieves the instances from a Remote AE. It is associated with local real-world activity "Query Instance Information" and "Retrieve Instance". When the "Query Instance Information" local real-world activity is requested by an operator, the Query/Retrieve Application Entity queries a remote AE for image information and provides the information list matching the query request. When "Retrieve Instance" local real-world activity is requested by an operator, the Query/Retrieve Application Entity requests a remote AE to transfer the instances.

4.1.2 Functional Definition of AEs

4.1.2.1 Functional Definition of Workflow Application Entity

The "Search Worklist" local activity is provided in the New Patient Registration GUI initiated by pressing the "NEW PATIENT" button on the console. The Patient ID, Patient Name, Accession Number, Scheduled Procedure Step Start Date, Scheduled Performing Physicians Name, and/or Requested Procedure ID may optionally be supplied before clicking the "Search" button in the GUI. Other default keys are the Modality (US), Scheduled Station AE Title (local AET), and Scheduled Procedure Step Start Date (Date of the day). When the "Search" button is clicked, the Workflow AE tries to open an association to a remote AE. If the Workflow AE establishes an association to a remote AE, it will transfer all worklist items via the open Association. During receiving the worklist response, the Workflow AE counts items and cancels the query processing, if the built-in limit of items (500) is reached. The results will be displayed in a list, which will be cleared with the next "Search Worklist" activity.

The Workflow AE automatically creates MPPS Instance when an image is sent to remote AE or stored in local drive for the first time in an examination. Further updates on the MPPS data can be performed from the "End Exam" user interface. The MPPS "Completed" or "Discontinued" states can only be set by the operator interaction. After a successful update of the MPPS Completed, the Workflow AE will issue a Storage Commitment request on images and Structured Reports already sent in the examination.

4.1.2.2 Functional Definition of Storage Application Entity

The Storage Application Entity can be requested in two modes. After the proper Worklist Item is selected or the patient identification is supplied by the operator, pressing the "STORE" button will directly send an image to the remote storage AE when it is configured to send to network. Or when it is configured to store the image in the local drive, the image is written in the drive for later reference. The images selected by the operator will be sent to the remote storage AE. A storage association will be initiated by clicking the "Send" button on the screen. or touching the "Send" button on the LCD touch panel.

A storage association for sending Structured Reports will be automatically initiated at the end of an examination if reports have been created in the examination. When measurements are performed under specific applications, Structured Reports will be created automatically at the end of the examination. Structured Reports can also be created manually by the operator.

The color of I-mark and R-mark overlaid on an image icon or the color of SR-mark overlaid on a report icon indicates the status of the image or report, respectively: **Green** - original, **Light Blue** - stored to media, **Orange** - sent to an Image Archive, and **Blue** - storage committed.

4.1.2.3 Functional Definition of Hardcopy Application Entity

The Hardcopy Application Entity locates in the "Print" button. By clicking the thumbnail the operator may select or deselect images to be printed by the remote AE. A print association will be initiated by clicking the "Print" button on the screen or touching the "Print(DICOM)" button on the LCD touch panel. After an association is established with the printer, its status is determined. If the printer is operating normally, the film sheets composed of selected images will be printed. If the printer is not operating normally, or the print returns a failure status during the association, the error is reported to the user.

4.1.2.4 Functional Definition of Query/Retrieve Application Entity

The Query/Retrieve Application Entity is provided in the Import GUI initialized by clicking the "Import" button on the "Find" Window. The Patient ID, Patient's Name, Accession Number, and/or Study Date may optionally be supplied before clicking the "Search" button. Other default key is Modality (US, CT, MR , PET or DOC). When the "Search" button is clicked, the Query/Retrieve AE will try to open an association to a remote AE. If the Query/Retrieve AE establishes an association to the remote AE, it will transfer all instance information items via the open Association. During receiving the response, the Query/Retrieve AE cancels the query processing, if the "Cancel" button is clicked by the operator. The results will be displayed in a separate list, which will be cleared with the next "Query Instance Information" activity.

When the operator selects instance information items in the list and clicks "Preview" or "Import" button, the Query/Retrieve AE will try to open an association to the remote AE for "Retrieve instance" activity. If the Query/Retrieve AE establishes an association to the remote AE, the Query/Retrieve AE opens a TCP/IP port and waits for another association requesting for the instance transfer from the remote AE. If the Remote AE establishes an association to the Query/Retrieve AE, it will transfer all instances via the Open Association from the remote AE. During receiving the response, the Query/Retrieve AE cancels the processing, if the "Cancel" button is clicking by the operator.

4.1.3 Sequencing of Real-World Activities



Figure 4.1-2 SEQUENCING CONSTRAINTS

Under normal scheduled workflow conditions the sequencing constraints illustrated in Figure 4.1-2 apply:

1. Query Worklist
2. Receive Worklist of Modality Scheduled Procedure Steps (MSPS)
3. Select Workitem (MSPS) from Worklist
4. Acquire Images

5. Create MPPS at the first Acquisition of Image, Encapsulated PDF, or Structured Report
6. Store acquired image instances
7. Create and store Structured Reports by "Send SR" user interface (optional step)
8. Create Encapsulated PDF reports (optional step)
9. Store Encapsulated PDF report instances (optional step)
10. Select Complete or Discontinue in "End Exam" user interface
11. Store unset image instances and store Structured Report instances (depends on conditions). Images are sent to remote storage AE over a separate association from Structured Report storage.
12. Complete acquisition and finalize MPPS
13. The Workflow AE will request Storage Commitment for the image and report instances, if the Image Manager is configured as an archive device.
14. Print acquired images (optional step)
15. Query Instance(Image or Document) Information (optional step)
16. Retrieve Instances (optional step)

Other workflow situations (e.g. unscheduled procedure steps) will have other sequencing constraints. Printing could equally take place after the acquired images have been stored. Printing could be omitted completely if no printer is connected or hardcopies are not required.

4.2 AE SPECIFICATIONS

4.2.1 Workflow Application Entity Specification

4.2.1.1 SOP Classes

This product provides Standard Conformance to the following SOP Classes:

Table 4.2-1 SOP CLASSES FOR AE WORKFLOW

SOP Class Name	SOP Class UID	SCU	SCP
Modality Worklist Information Model – FIND	1.2.840.10008.5.1.4.31	Yes	No
Modality Performed Procedure Step	1.2.840.10008.3.1.2.3.3	Yes	No
Storage Commitment Push Model	1.2.840.10008.1.20.1	Yes	No
Verification	1.2.840.10008.1.1	Yes	Yes

4.2.1.2 Association Policies

4.2.1.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed:

Table 4.2-2 DICOM APPLICATION CONTEXT FOR AE WORKFLOW

Application Context Name	1.2.840.10008.3.1.1.1
--------------------------	-----------------------

4.2.1.2.2 Number of Associations

This product initiates one Association at a time for a Workflow requests.

Table 4.2-3 NUMBER OF ASSOCIATIONS INITIATED FOR AE WORKFLOW

Maximum number of simultaneous Associations	1
---	---

4.2.1.2.3 Asynchronous Nature

This product does not support asynchronous communication (multiple outstanding transactions over a single Association).

Table 4.2-4 ASYNCHRONOUS NATURE AS A SCU FOR AE WORKFLOW

Maximum number of outstanding asynchronous transactions	1
---	---

4.2.1.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

Table 4.2-5 DICOM IMPLEMENTATION CLASS AND VERSION FOR AE WORKFLOW

Implementation Class UID	1.2.392.200039.117
Implementation Version Name	20180308 (subject to change without notice)

4.2.1.3 Association Initiation Policy**4.2.1.3.1 Activity – Find Worklist****4.2.1.3.1.1 Description and Sequencing of Activities**

An interactive query for Worklist is initiated by pressing the "Find" button in the "NEW PATIENT" Registration GUI. The built-in query keys are the Modality (US), Scheduled Procedure Step Start Date (actual date). The Scheduled Station AE Title may be excluded by the configuration. Additional "Patient-based" query keys, Patient's ID, Patient Name, Accession Number, and/or Requested Procedure ID, may be supplied in the dialog by the operator.

Upon initiation of the request, this product will build an Identifier for the C-FIND request, will initiate an Association to send the request and will wait for Worklist responses. After retrieval of all responses, this product will access the local database to update patient demographic data. To protect the system from overflow, this product will limit the number of processed worklist responses to 500. During receiving the worklist response items are counted and the query processing is canceled by issuing a C-FIND-CANCEL if the limit of items is reached. The results will be displayed in a separate list.

The retrieved Worklist items are stored locally during the day, which will be cleared with the next worklist update. If the list is a latest and additional examination is to be performed on a patient, or the equipment is disconnected for mobile examination, the stored worklist items may be referenced by pressing

"Worklist" button in the "NEW PATIENT" Registration GUI. The additional examination using the same MSPS generates a second series of images coping with the Append Case among the IHE use cases.




Figure 4.2-1 SEQUENCING OF ACTIVITY – WORKLIST QUERY

A possible sequence of interactions between the Workflow AE and a Departmental Scheduler (e.g. a device such as a RIS or HIS which supports the Modality Worklist SOP Class as an SCP) is illustrated in the Figure above:

1. The Workflow AE opens an association with the Department Scheduler.
2. The Workflow AE sends a C-FIND request to the Department Scheduler containing the Worklist Query attributes.
3. The Department Scheduler returns a C-FIND response containing the requested attributes of the first matching Worklist Item.
4. The Department Scheduler returns another C-FIND response containing the requested attributes of the second matching Worklist Item.
5. The Department Scheduler returns another C-FIND response with status Success indicating that no further matching Worklist Items exist. This example assumes that only 2 Worklist items match the Worklist Query.
6. The Workflow AE closes the association with the Department Scheduler.
7. The user selects a Worklist Item from the Worklist database and prepares to acquire images.

4.2.1.3.1.2 Proposed Presentation Contexts

This product will propose Presentation Contexts as shown in the following table:

Table 4.2-6 PROPOSED PRESENTATION CONTEXTS FOR ACTIVITY WORKLIST QUERY

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Ext. Neg.
Name	UID	Name List	UID List		
Modality Worklist Information Model – FIND	1.2.840.10008.5.1.4.31	Implicit VR Little Endian Explicit VR Little Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.1	SCU	None
Verification	1.2.840.10008.1.1	Implicit VR Little Endian Explicit VR Little Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.1	SCU	None

4.2.1.3.1.3 SOP Specific Conformance for Modality Worklist

The behavior of this product when encountering status codes in a Modality Worklist C-FIND response is summarized in the Table below. If any other SCP response status than "Success" or "Pending" is received by this product, a Worklist Error Message will appear on the user interface.

Table 4.2-7 MODALITY WORKLIST C-FIND RESPONSE STATUS HANDLING BEHAVIOR

Service Status	Further Meaning	Error Code	Behavior
Success	Matching is complete	0000	The SCP has completed the matches. Worklist items are available for display or further processing.
Refused	Out of Resources	A700	The Association is closed using A-RELEASE and the worklist query is failed.
Failed	Identifier does not match SOP Class	A900	The Association is closed using A-RELEASE and the worklist query is failed.
Failed	Unable to Process	C000 – CFFF	The Association is closed using A-RELEASE and the worklist query is failed.
Cancel	Matching terminated due to Cancel request	FE00	If the query was cancelled due to too may worklist items then the SCP has completed the matches. Worklist items are available for display or further processing. Otherwise, the Association is aborted using A-ABORT and the worklist query is failed.
Pending	Matches are continuing	FF00	The worklist item contained in the Identifier is collected for later display or further processing.
Pending	Matches are continuing – Warning that one or more Optional Keys were not supported	FF01	The worklist item contained in the Identifier is collected for later display or further processing.
*	*	Any other status code.	The Association is closed using A-RELEASE and the worklist is failed.

The behavior of this product during communication failure is summarized in the Table below.

Table 4.2-8 MODALITY WORKLIST COMMUNICATION FAILURE BEHAVIOR

Exception	Behavior
Timeout	The Association is aborted using A-ABORT and the worklist query is failed. A Worklist Error is reported to the user.
Association aborted by the SCP or network layers	The worklist query is failed. A Worklist Error is reported to the user.

Acquired images will always use the Study Instance UID specified for the Scheduled Procedure Step (if available). If an acquisition is unscheduled, a Study Instance UID will be generated locally.

The Table below provides a description of the Worklist Request Identifier of this product and specifies the attributes that are copied into the images. Unexpected attributes returned in a C-FIND response are ignored.

Requested return attributes not supported by the SCP are set to have no value. Non-matching responses returned by the SCP due to unsupported optional matching keys are ignored. No attempt is made to filter out possible duplicate entries.

Table 4.2-9 WORKLIST REQUEST IDENTIFIERS

Module Name Attribute Name	Tag	VR	M	R	Q	D	IOD
SOP Common Specific Character Set	(0008,0005)	CS		x			x

Scheduled Procedure Step							
Scheduled Procedure Step Sequence	(0040,0100)	SQ		x			
> Modality	(0008,0060)	CS	(S)	x		x	x
> Requested Contrast Agent	(0032,1070)	LO		x		x	
> Scheduled Station AET	(0040,0001)	AE	(S)	x		x	
> Scheduled Procedure Step Start Date	(0040,0002)	DA	R	x	x	x	x
> Scheduled Procedure Step Start Time	(0040,0003)	TM		x		x	x
> Scheduled Performing Physician's Name	(0040,0006)	PN		x	x	x	
> Scheduled Procedure Step Description	(0040,0007)	LO		x		x	x
> Scheduled Protocol Code Sequence	(0040,0008)	SQ		x		x	x
> Scheduled Procedure Step ID	(0040,0009)	SH		x		x	x
> Scheduled Station Name	(0040,0010)	SH		x		x	
> Scheduled Procedure Step Location	(0040,0011)	SH		x		x	
> Pre-Medication	(0040,0012)	LO		x		x	
Requested Procedure							
Referenced Study Sequence	(0008,1110)	SQ		x			x
Referenced Patient Sequence	(0008,1120)	SQ		x			
Study Instance UID	(0020,000D)	UI		x		x	x
Requested Procedure Description	(0032,1060)	LO		x		x	x
Requested Procedure Code Sequence	(0032,1064)	SQ		x		x	x
Requested Procedure ID	(0040,1001)	SH		x	x	x	x
Requested Procedure Priority	(0040,1003)	SH		x		x	
Patient Transport Arrangements	(0040,1004)	LO		x		x	

Imaging Service Request							
Accession Number	(0008,0050)	SH		x	x	x	x
Referring Physician's Name	(0008,0090)	PN		x		x	x
Requesting Physician	(0032,1032)	PN		x		x	x
Requesting Service	(0032,1033)	LO		x		x	
Admission ID	(0038,0010)	LO		x		x	
Patient Identification							
Patient Name	(0010,0010)	PN		x	x	x	x
Patient ID	(0010,0020)	LO		x	x	x	x
Patient Demographic							
Patient's Birth Date	(0010,0030)	DA		x		x	x
Patient's Sex	(0010,0040)	CS		x		x	x
Patient's Size	(0010,1020)	DS		x		x	x
Patient's Weight	(0010,1030)	DS		x		x	x
Ethnic Group	(0010,2160)	SH		x		x	
Occupation	(0010,2180)	SH		x		x	x
Patient Comments	(0010,4000)	LT		x		x	
Confidentiality Constraint on Patient Data Description	(0040,3001)	LO		x		x	
Patient Medical							
Medical Alerts	(0010,2000)	LO		x		x	x
Contrast Allergies	(0010,2110)	LO		x		x	x
Pregnancy Status	(0010,21C0)	US		x		x	
Last Menstrual Date	(0010,21D0)	DA		x		x	x
Special Needs	(0038,0050)	LO		x		x	
Patient State	(0038,0500)	LO		x		x	
Visit Status							
Current Patient Location	(0038,0300)	LO		x		x	

The above table should be read as follows:

Module Name: The name of the associated module for supported worklist attributes.

Attribute Name: Attributes supported to build a Worklist Request Identifier of this product.

Tag: DICOM tag for this attribute.

VR: DICOM VR for this attribute.

M: Matching keys for Worklist Query. An "S" will indicate that this product will supply an attribute value for Single Value Matching, an "R" will indicate Range Matching and a "*" will denote wildcard matching. An "(S)" will indicate that keys are configurable either Single Value Matching or Universal Matching.

R: Return keys. An "x" will indicate that this product will supply this attribute as Return Key with zero length for Universal Matching.

Q: Interactive Query Key. An "x" will indicate that this product will supply this attribute as matching key, if entered in the New Patient dialog. For example, the Patient Name can be entered thereby restricting Worklist responses to Procedure Steps scheduled for the patient.

D: Displayed keys. An "x" indicates that this worklist attribute is displayed to the user during a patient registration dialog. For example, Patient Name will be displayed when registering the patient prior to an examination.

IOD: An "x" indicates that this Worklist attribute is included into all Object Instances created during performance of the related Procedure Step.

The Modality (0008,0060) matching key is configurable to "US" or zero length (Universal matching), the Scheduled Station AET (0040,0001) key is configurable to the local AET (ex. LISENDO 880) or zero length, and the Scheduled Procedure Start Date (0040,0002) key is configurable to a date, date range, or zero length.

4.2.1.3.2 Activity – End Exam

4.2.1.3.2.1 Description and Sequencing of Activities

After the "NEW PATIENT" Registration, this product is awaiting local image storage or image transmission to remote Archive. The trigger to create a MPPS SOP Instance is derived from this event. An Association to the configured MPPS SCP system is established immediately and the related MPPS SOP Instance will be created.

The MPPS user interface is initiated by pressing the "END EXAM" button on the console. Only a manual update can be performed by the operator where is it possible to set the final state of the MPPS to "COMPLETED" or "DISCONTINUED". In the "Discontinued" case the user can also select the discontinuation reason from a pick list corresponding to Context Group 9300. An MPPS Instance that has been sent with a state of "COMPLETED" or "DISCONTINUED" can no longer be updated.

This product will support creation of "unscheduled cases" by allowing MPPS Instances to be communicated for locally registered Patients. This product only supports a 0-to-1 relationship between Scheduled and Performed Procedure Steps.

This product will initiate an Association to issue an:

- N-CREATE request according to the CREATE Modality Performed Procedure Step SOP Instance operation or a
- N-SET request to update the contents and state of the MPPS according to the SET Modality Performed Procedure Step Information operation.




Figure 4.2-2 SEQUENCING OF ACTIVITY – ACQUIRE IMAGES AND END EXAM

A possible sequence of interactions between the Workflow AE and a Department Scheduler (e.g. a device such as a RIS or HIS which supports the MPPS SOP Class as an SCP) is illustrated in Figure 4.2-2:

1. First image is acquired for the MSPS.
2. The Workflow AE opens an association with the Department Scheduler.
3. The Workflow AE sends an N-CREATE request to the Department Scheduler to create an MPPS instance with status of “IN PROGRESS” and create all necessary attributes. The Department Scheduler acknowledges the MPPS creation with an N-CREATE response (status success).
4. The Workflow AE closes the association with the Department Scheduler.
5. After all images are sent to Image Archive or stored in the local database, the operator requests "End Exam" and closes database the examination. At this time Structured Reports, if created, are sent to Image Archive automatically.
6. The Workflow AE opens a second association with the Department Scheduler.
7. The Workflow AE sends an N-SET request to the Department Scheduler to update the MPPS instance with status of “COMPLETED” and set all necessary attributes. The Department Scheduler acknowledges the MPPS update with an N-SET response (status success).

8. The Workflow AE closes the association with the Department Scheduler.
9. If the images/reports associated with the examination have been sent to the Image Manager, the Workflow AE opens independent association with the Image Manager.
10. The Workflow AE sends an N-ACTION request to the Image Manager to request the Images/reports are storage-committed. The Image Manager acknowledges the Storage Commitment Request with an N-ACTION response (status success).
11. The Workflow AE closes the association with the Image Manager.

4.2.1.3.2.2 Proposed Presentation Contexts

This product will propose Presentation Contexts as shown in the following table:

**Table 4.2-10
PROPOSED PRESENTATION CONTEXTS FOR REAL-WORLD ACTIVITY**

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Ext. Neg.
Name	UID	Name List	UID List		
Modality Performed Procedure Step	1.2.840.10008.3.1.2.3.3	Implicit VR Little Endian Explicit VR Little Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.1	SCU	None
Verification	1.2.840.10008.1.1	Implicit VR Little Endian Explicit VR Little Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.1	SCU	None

4.2.1.3.2.3 SOP Specific Conformance for MPPS

The behavior of this product when encountering status codes in an MPPS N-CREATE or N-SET response is summarized in Table 4.2-11. If any other SCP response status than "Success" is received by this product, a message "Network communication error" will appear on the user interface.

Table 4.2-11 MPPS N-CREATE / N-SET RESPONSE STATUS HANDLING BEHAVIOR

Service Status	Further Meaning	Error Code	Behavior
Success	Success	0000	The SCP has completed the operation successfully.
Failure	Processing Failure – Performed Procedure Step Object may no longer be updated	0110	The Association is closed using A-RELEASE. The error status is reported to the user.
Warning	Attribute Value Out of Range	0116H	The Association is closed using A-RELEASE. The error status is reported to the user.
*	*	Any other status code.	The Association is closed using A-RELEASE. The error status is reported to the user.

The behavior of this product during communication failure is summarized in the Table below:

Table 4.2-12 MPPS COMMUNICATION FAILURE BEHAVIOR

Exception	Behavior
Timeout	The Association is aborted using A-ABORT and the timeout status is reported to the user.
Association aborted by the SCP or network layers	The MPPS is failed and the status is reported to the user.

Table 4.2-13 provides a description of the MPPS N-CREATE and N-SET request identifiers sent by this product. Empty cells in the N-CREATE and N-SET columns indicate that the attribute is not sent. An “x” indicates that an appropriate value will be sent. A “Zero length” attribute will be sent with zero length.

Table 4.2-13 MPPS N-CREATE / N-SET REQUEST IDENTIFIER

Attribute Name	Tag	VR	N-CREATE	N-SET
Specific Character Set	(0008,0005)	CS	ISO_IR 100	-
Referenced Patient Sequence	(0008,1120)	SQ	Zero length	-
Patient's Name	(0010,0010)	PN	From Modality Worklist or user input. The user can modify values provided via Modality Worklist.	-
Patient ID	(0010,0020)	LO	From Modality Worklist or user input. The user can modify values provided via Modality Worklist.	-
Patient's Birth Date	(0010,0030)	DA	From Modality Worklist or user input. The user can modify values provided via Modality Worklist.	-
Patient's Sex	(0010,0040)	CS	From Modality Worklist or user input. The user can modify values provided via Modality Worklist.	-
Scheduled Step Attributes Sequence	(0040,0270)	SQ	If the procedure step creates an SOP Instance	-
> Accession Number	(0008,0050)	SH	From Modality Worklist or user input. The user can modify values provided via Modality Worklist.	-
>Referenced Study Sequence	(0008,1110)	SQ	From Modality Worklist	-
>>Referenced SOP Class UID	(0008,1150)	UI	From Modality Worklist	-
>>Referenced SOP Instance UID	(0008,1155)	UI	From Modality Worklist	-
>Study Instance UID	(0020,000D)	UI	From Modality Worklist or internally generated	-
>Requested Procedure Description	(0032,1060)	LO	From Modality Worklist	-

>Scheduled Procedure Step Description	(0040,0007)	LO	From Modality Worklist	-
>Scheduled Protocol Code Sequence	(0040,0008)	SQ	From Modality Worklist	-
>Scheduled Procedure Step ID	(0040,0009)	SH	From Modality Worklist	-
>Requested Procedure ID	(0040,1001)	SH	From Modality Worklist or user input. The user can modify values provided via Modality Worklist.	-
Procedure Code Sequence	(0008,1032)	SQ	From Modality Worklist, mapped from Requested Procedure Code Sequence (0032,1064)	-
Performed Station AE Title	(0040,0241)	AE	AE Title of the Equipment	-
Performed Station Name	(0040,0242)	SH	Station Name of the Equipment	-
Performed Location	(0040,0243)	SH	Zero length	-
Performed Procedure Step Start Date	(0040,0244)	DA	Actual start date	-
Performed Procedure Step Start Time	(0040,0245)	TM	Actual start time	-
Performed Procedure Step End Date	(0040,0250)	DA	Zero length	Actual end date
Performed Procedure Step End Time	(0040,0251)	TM	Zero length	Actual end time
Performed Procedure Step Status	(0040,0252)	CS	IN PROGRESS	COMPLETED or DISCONTINUED
Performed Procedure Step ID	(0040,0253)	SH	Automatically created.	-
Performed Procedure Step Description	(0040,0254)	LO	Input by the user at "Study Description" in New Patient Registration.	-
Performed Procedure Type Description	(0040,0255)	LO	Zero length	-
Performed Procedure Step Discontinuation Reason Code Sequence	(0040,0281)	SQ	Zero length	If Performed Procedure Step Status (0040,0252) is "DISCONTINUED" and a discontinuation reason is selected by the user, then a single item will be present containing a user-selected entry drawn from Context Group 9300.
Modality	(0008,0060)	CS	US	-
Study ID	(0020,0010)	SH	Copied from Requested Procedure ID (0040,1001) in	-

			MWL. The user can modify values provided via Modality Worklist.	
Performed Protocol Code Sequence	(0040,0260)	SQ	Zero length	Zero or more items
Performed Series Sequence	(0040,0340)	SQ	Zero length	One or more items
>Retrieve AE Title	(0008,0054)	AE		Zero length
>Series Description	(0008,103E)	LO		x
>Performing Physician's Name	(0008,1050)	PN		x
>Operator's Name	(0008,1070)	PN		x
>Referenced Image Sequence	(0008,1140)	SQ		Zero or more items
>>Referenced SOP Class UID	(0008,1150)	UI		x
>>Referenced SOP Instance UID	(0008,1155)	UI		x
>Protocol Name	(0018,1030)	LO		x
>Series Instance UID	(0020,000E)	UI		x
>Referenced Non-image Composite SOP Instance Sequence	(0040,0220)	SQ		Zero or more items
>>Referenced SOP Class UID	(0008,1150)	UI		x
>>Referenced SOP Instance UID	(0008,1155)	UI		x

4.2.1.3.2.4 SOP Specific Conformance for Storage Commitment SOP Class

4.2.1.3.2.4.1 Storage Commitment Operations (N-ACTION)

The Workflow AE will request storage commitment for instances of the Ultrasound Image, Ultrasound Multi-frame Image and Structured Report Storage Classes, if the remote AE is configured as an Archive Device and a presentation context for the Storage Commitment Push Model has been accepted.

The Workflow AE will consider Storage Commitment failed if no N-EVENT-REPORT is received for a Transaction UID within a configurable time period after receiving a successful N-ACTION response (duration of applicability for a Transaction UID).

The Workflow AE does not send the optional Storage Media File Set ID & UID Attributes or the Referenced Performed Procedure Step Sequence Attribute in the N-ACTION.

The behavior of Storage AE when encountering status codes in an N-ACTION response is summarized in the Table below:

**Table 4.2-14
STORAGE COMMITMENT N-ACTION RESPONSE STATUS HANDLING BEHAVIOR**

Service Status	Further Meaning	Error Code	Behavior
Success	Success	0000	The request for storage commitment is considered successfully sent. A timer is started which will expire if no N-EVENT-REPORT for the Transaction UID is received within a configurable timeout period.
*	*	Any other status code.	The Association is closed using A-RELEASE and the failure is reported to the user.

The behavior of Workflow AE during communication failure is summarized in the Table below:

Table 4.2-15 STORAGE COMMITMENT COMMUNICATION FAILURE BEHAVIOR

Exception	Behavior
Timeout	The Association is aborted using A-ABORT and the timeout error is reported to the user.
Association aborted by the SCP or network layers	The association failure is reported to the user.

4.2.1.3.2.4.2 Storage Commitment Notifications (N-EVENT-REPORT)

The Workflow AE is capable of receiving an N-EVENT-REPORT notification if it has successfully negotiated a Presentation Context for the Storage Commitment Push Model (i.e. only associations established with archive devices).

Upon receipt of an N-EVENT-REPORT the timer associated with the Transaction UID will be canceled.

The behavior of Workflow AE when receiving Event Types within the N-EVENT-REPORT is summarized in the Table below.

Table 4.2-16 STORAGE COMMITMENT N-EVENT-REPORT BEHAVIOUR

Event Type Name	Event Type ID	Behavior
Storage Commitment Request Successful	1	The Referenced SOP Instances under Referenced SOP Sequence (0008,1199) are marked within the database as committed. Successfully committed SOP Instances are candidates for deletion from the local database. The least recently accessed SOP Instances are deleted first.
Storage Commitment Request Complete – Failures Exist	2	The Referenced SOP Instances under Referenced SOP Sequence (0008,1199) are treated in the same way as in the success case (Event Type 1). The Referenced SOP Instances under Failed SOP Sequence (0008,1198) are not marked as committed within the database.

The reasons for returning specific status codes in an N-EVENT-REPORT response are summarized in the Table below:

Table 4.2-17 STORAGE COMMITMENT N-EVENT-REPORT RESPONSE STATUS REASONS

Service Status	Further Meaning	Error Code	Reasons
Success	Success	0000	The storage commitment result has been successfully received.
Failure	Unrecognized Operation	0211H	The Transaction UID in the N-EVENT-REPORT request is not recognized (was never issued within an N-ACTION request).
Failure	Resource Limitation	0213H	The Transaction UID in the N-EVENT-REPORT request has expired (no N-EVENT-REPORT was received within a configurable time limit).
Failure	No Such Event Type	0113H	An invalid Event Type ID was supplied in the N-EVENT-REPORT request.
Failure	Processing Failure	0110H	An internal error occurred during processing of the N-EVENT-REPORT. A short description of the error will be returned in Error Comment (0000,0902).
Failure	Invalid Argument Value	0115H	One or more SOP Instance UIDs with the Referenced SOP Sequence (0008,1199) or Failed SOP Sequence (0008,1198) was not included in the Storage Commitment Request associated with this Transaction UID. The unrecognized SOP Instance UIDs will be returned within the Event Information of the N-EVENT-REPORT response.

4.2.1.4 Association Acceptance Policy

This product accepts Associations to receive N-EVENT-REPORT notifications for the Storage Commitment Push Model SOP Class.

Table 4.2-18 NUMBER OF ASSOCIATIONS ACCEPTED FOR WORKFLOW AE

Maximum number of simultaneous Associations accepted	1
--	---

4.2.1.4.1 Activity – Receive Storage Commitment Response**4.2.1.4.1.1 Description and Sequencing of Activities**

The Workflow AE will accept associations in order to receive responses to a Storage Commitment N-EVENT-REPORT Request.




Figure 4.2-3 SEQUENCING OF ACTIVITY - RECEIVE STORAGE COMMITMENT RESPONSE

A possible sequence of interactions between the Workflow AE and an Image Manager (e.g. a storage or archive device supporting Storage Commitment SOP Classes as an SCP) is illustrated in the Figure above:

1. The Image Manager opens a new association with the Workflow AE.
2. The Image Manager sends an N-EVENT-REPORT request notifying the Workflow AE of the status of a previous Storage Commitment Request. The Workflow AE replies with an N-EVENT-REPORT response confirming receipt.
3. The Image Manager closes the association with the Workflow AE.

The Workflow AE may reject association attempts as shown in the Table below. The Result, Source and Reason/Diag columns represent the values returned in the appropriate fields of an ASSOCIATE-RJ PDU (see PS 3.8, Section 9.3.4). The contents of the Source column is abbreviated to save space and the meaning of the abbreviations are:

- a) 1 – DICOM UL service-user
- b) 2 – DICOM UL service-provider (ASCE related function)
- c) 3 – DICOM UL service-provider (Presentation related function)

Table 4.2-19 ASSOCIATION REJECTION REASONS

Result	Source	Reason/Diag	Explanation
2 – rejected-transient	c	2 – local-limit-exceeded	The maximum number of simultaneous associations has been reached. An association request with the same parameters may succeed at a later time.
2 – rejected-transient	c	1 – temporary-congestion	No associations can be accepted at this time because insufficient resources are available (e.g. memory, processes, or threads). An association request with the same parameters may succeed at a later time.
1 – rejected-permanent	a	2 – application-context-name-not-supported	The association request contained an unsupported Application Context Name. An association request with the same parameters will not succeed at a later time.
1 – rejected-permanent	a	7 – called-AE-title-not-recognized	The association request contained an unrecognized Called AE Title. An association request with the same parameters will not succeed at a later time unless configuration changes are made. This rejection reason normally occurs when the

			association initiator is incorrectly configured and attempts to address the association acceptor using the wrong AE Title.
1 – rejected-permanent	a	3 – calling-AE-title-not-recognized	The association request contained an unrecognized Calling AE Title. An association request with the same parameters will not succeed at a later time unless configuration changes are made. This rejection reason normally occurs when the association acceptor has not been configured to recognize the AE Title of the association initiator.
1 – rejected-permanent	b	1 – no-reason-given	The association request could not be parsed. An association request with the same format will not succeed at a later time.

4.2.1.4.1.2 Accepted Presentation Contexts

The Workflow AE will accept Presentation Contexts as shown in the Table below:

**Table 4.2-20 ACCEPTABLE PRESENTATION CONTEXTS
FOR ACTIVITY RECEIVE STORAGE COMMITMENT RESPONSE**

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Ext. Neg.
Name	UID	Name List	UID List		
Storage Commitment Push Model	1.2.840.10008.1.20.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
Verification	1.2.840.10008.1.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		

The Workflow AE will prefer to select the Explicit VR Little Endian Transfer Syntax if multiple transfer syntaxes are offered. The Workflow AE will only accept the SCU role (which must be proposed via SCP/SCU Role Selection Negotiation) within a Presentation Context for the Storage Commitment Push Model SOP Class.

4.2.1.4.1.3 SOP Specific Conformance for Storage Commitment SOP Class

4.2.1.4.1.3.1 Storage Commitment Notifications (N-EVENT-REPORT)

Upon receipt of an N-EVENT-REPORT the timer associated with the Transaction UID will be canceled.

The behavior of Workflow AE when receiving Event Types within the N-EVENT-REPORT is summarized in Table 4.2-16.

The reasons for returning specific status codes in an N-EVENT-REPORT response are summarized in Table 4.2-17.

4.2.1.4.1.4 SOP Specific Conformance for Verification SOP Class

The Workflow AE provides standard conformance to the Verification SOP Class as an SCP. If the C-ECHO request was successfully received, a 0000 (Success) status code will be returned in the C-ECHO response. Otherwise, the Association is aborted by A-ABORT by the Verification SCP.

4.2.2 Storage Application Entity Specification

4.2.2.1 SOP Classes

This product provides the following Conformance to the SOP Classes listed below:

Table 4.2-21 SOP CLASSES FOR AE STORAGE

SOP Class Name	SOP Class UID	Conformance	SCU	SCP
Ultrasound Image Storage	1.2.840.10008.5.1.4.1.1.6.1	Standard Extended	Yes	No
Ultrasound Multi-frame Image Storage	1.2.840.10008.5.1.4.1.1.3.1	Standard Extended	Yes	No
Comprehensive SR	1.2.840.10008.5.1.4.1.1.88.33	Standard	Yes	No
Encapsulated PDF Storage	1.2.840.10008.5.1.4.1.1.104.1	Standard	Yes	No
Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7	Standard Extended	Yes	No
Multi-frame True Color Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.4	Standard Extended	Yes	No
Verification	1.2.840.10008.1.1	Standard	Yes	No

4.2.2.2 Association Policies

4.2.2.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed:

Table 4.2-22 DICOM APPLICATION CONTEXT FOR AE STORAGE

Application Context Name	1.2.840.10008.3.1.1.1
--------------------------	-----------------------

4.2.2.2.2 Number of Associations

This product initiates one Association at a time for the destination to which a transfer request is being processed. It does not automatically request associations to multiple destinations.

Table 4.2-23 NUMBER OF ASSOCIATIONS INITIATED FOR AE STORAGE

Maximum number of simultaneous Associations initiated	1
---	---

Note: However, each AE for image storage, Structured Report storage, and Encapsulated PDF storage is independent. Therefore two associations for image storage and report storage may be initiated simultaneously.

4.2.2.2.3 Asynchronous Nature

This product does not support asynchronous communication (multiple outstanding transactions over a single Association).

Table 4.2-24 ASYNCHRONOUS NATURE AS A SCU FOR AE STORAGE

Maximum number of outstanding asynchronous transactions	1
---	---

4.2.2.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

Table 4.2-25 DICOM IMPLEMENTATION CLASS AND VERSION FOR AE STORAGE

Implementation Class UID	1.2.392.200039.117
Implementation Version Name	20180308 (subject to change without notice)

4.2.2.3 Association Initiation Policy

4.2.2.3.1 Activity – Send Images, Structured Reports, and Encapsulated PDF reports

4.2.2.3.1.1 Description and Sequencing of Activities

The Storage AE for sending images may be invoked in two ways. After the New Patient registration is completed selecting the worklist MSPS or by manual entry, ultrasound image is displayed in real-time. When the operator presses the "FREEZE" button on the console, the image is frozen and ready to send the image to remote storage AE or to store in a local drive. If the operator presses the "STORE" button in real-time, a multi-frame image may be acquired to send. The "STORE" button may be configured for either NET (DICOM), HDD, DVD-RAM, CD-R buffer or USB Media.

When it is configured for "NET (DICOM)" storage, pressing the "STORE" button will open an association to the remote Storage SCP. If the association is accepted by the SCP, local Storage AE will send an ultrasound Image Instance or a multi-frame image instance to the SCP, and if the SCP responds with a success status the association is closed normally. If the association is rejected or it is not responded within the configured time interval, or the Storage AE receives failure status, the association is aborted and the failure is reported to the user. The Storage AE will retain a copy of the Image Instance sent directly from the "STORE" button.

When the "STORE" button is configured for HDD storage, an Image Instance will be stored in the local HDD each time it is pressed. The Image Instances are reviewed on the "Tile" screen page or "Full screen" screen page. After the user selects one or more images, clicking the "Send" button on the screen or pressing the "Send" button on the LCD touch panel will open a storage association to the remote Archive. If the association is accepted, the Storage AE will send the all instances of the selected images to the remote Storage SCP within the single association. Each Instance successfully stored is indicated by an **orange** mark on the icon, and the storage-committed Instance is indicated by a **blue** mark. If an Image Instance is not responded by the successful status, or the Storage SCP does not respond within the configured interval, the association is aborted and the status is reported to the user. The Storage AE can be configured to retry automatically.

To create Encapsulated PDF reports, Advanced Cardiac feature needs to be configured. When the "Create DICOM PDF" button is pressed in the Advanced Reporting Window, Encapsulated PDF report will be created. The created PDF reports can be reviewed on the "Find" Window. After the user selects one or more PDF reports, clicking the "Send" button on the screen or pressing the "Send" button on the LCD touch panel will open a storage association to the remote Archive. If the association is accepted, the Storage AE will send the all instances of the selected PDF reports to the remote Storage SCP within the single association. Indication color scheme of PDF report on the icon is the same as one of image. If an Encapsulated PDF Instance is not responded by the successful status, or the Storage SCP does not respond within the configured interval, the association is aborted and the status is reported to the user. The Storage AE can be configured to retry automatically.

Images, Structured Reports, and PDF Reports are sent to remote storage AE over a separate association with each other.




Figure 4.2-4 SEQUENCING OF ACTIVITY – SEND IMAGES AND REPORTS

A possible sequence of interactions between the Storage AE and an Image Manager (e.g. a storage or archive device supporting the Storage SOP Classes as an SCP) is illustrated in Figure 4.2-4:

1. The Storage AE opens an association with the Image Manager
2. An acquired US image is transmitted to the Image Manager using a C-STORE request and the Image Manager replies with a C-STORE response (status success).
3. Another acquired US Multi-frame image is transmitted to the Image Manager using a C-STORE request and the Image Manager replies with a C-STORE response (status success).
4. The Storage AE closes the association with the Image Manager.
5. The Storage AE opens another association with the Image Manager.
6. An acquired SC image is transmitted to the Image Manager using a C-STORE request and the Image Manager replies with a C-STORE response (status success).
7. Another acquired SC Multi-frame image is transmitted to the Image Manager using a C-STORE request and the Image Manager replies with a C-STORE response (status success).
8. The Storage AE closes the association with the Image Manager.
9. The Storage AE opens another association with the Image Manager.
10. Structured Reports are transmitted to the Image Manager using C-STORE request and the Image Manager replies with a C-STORE response (status success).
11. The Storage AE closes the association with the Image Manager.
12. The Storage AE opens another association with the Image Manager.
13. Encapsulated PDF reports are transmitted to the Image Manager using C-STORE request and the Image Manager replies with a C-STORE response (status success).
14. The Storage AE closes the association with the Image Manager.

4.2.2.3.1.2 Proposed Presentation Contexts

This product is capable of proposing the Presentation Contexts shown in the following table:

Table 4.2-26 PROPOSED PRESENTATION CONTEXTS FOR ACTIVITY SEND IMAGES / REPORTS

Send Ultrasound Image Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Ext. Neg.
Name	UID	Name List	UID List		
Ultrasound Image Storage	1.2.840.10008.5.1.4.1.1.6.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		JPEG Baseline Compression	1.2.840.10008.1.2.4.50		
		RLE Compression	1.2.840.10008.1.2.5		

Ultrasound Multi-frame Image Storage	1.2.840.10008. 5.1.4.1.1.3.1	Implicit VR Little Endian Explicit VR Little Endian JPEG Baseline Compression ¹ RLE Compression	1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.5	SCU	None
Send Secondary Capture Image Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Ext. Neg.
Name	UID	Name List	UID List		
Secondary Capture Image Storage	1.2.840.10008. 5.1.4.1.1.7	Implicit VR Little Endian Explicit VR Little Endian JPEG Baseline Compression RLE Compression	1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.5	SCU	None
Multi-frame True Color Secondary Capture Image Storage	1.2.840.10008. 5.1.4.1.1.7.4	Implicit VR Little Endian Explicit VR Little Endian JPEG Baseline Compression RLE Compression	1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.5	SCU	None
Send Structured Report Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Ext. Neg.
Name	UID	Name List	UID List		
Comprehensive SR	1.2.840.10008. 5.1.4.1.1.88.33	Implicit VR Little Endian Explicit VR Little Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.1	SCU	None
Enhanced SR ²	1.2.840.10008. 5.1.4.1.1.88.22	Implicit VR Little Endian Explicit VR Little Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.1	SCU	None
Send Encapsulated PDF Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Ext. Neg.
Name	UID	Name List	UID List		
Encapsulated PDF Storage	1.2.840.10008. 5.1.4.1.1.104.1	Implicit VR Little Endian Explicit VR Little Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.1	SCU	None

Note: 1. except for images acquired for 3D Reconstruction.

Note: 2 only Non-Advanced Cardiac.

Presentation Contexts for Ultrasound Image Storage and Ultrasound Multi-frame Image Storage will always be proposed. An error message will be issued at sending the SOP Instance of the Presentation Context of which Abstract Syntax has been rejected by the remote AE.

Presentation Contexts for Comprehensive SR Storage and Enhanced SR Storage will always be proposed. When the remote AE accepts Comprehensive SR the reports will be sent as Comprehensive SR, or if the remote AE rejects Comprehensive SR the reports may be sent as Enhanced SR.

Presentation Contexts for Encapsulated PDF Storage will always be proposed. An error message will be issued at sending the SOP Instance of the Presentation Context of which Abstract Syntax has been rejected by the remote AE.

4.2.2.3.1.3 SOP Specific Conformance for Image and Report Storage SOP Classes

All Image and Structured Report Storage SOP Classes supported by the Storage AE exhibit the same behavior, except where stated, and are described together in this section.

If Ultrasound Multi-frame Image Storage SOP Instances are included in the Send Job and a corresponding Presentation Context is not accepted then the Association is aborted using A-ABORT and the send job is marked as failed. The job failure is reported to the user via the dialog window. The remaining Ultrasound (single frame) Image Storage SOP Instances should be selected and retried by the operator.

Table 4.2-27 STORAGE C-STORE RESPONSE STATUS HANDLING BEHAVIOR

Service Status	Further Meaning	Error Code	Behavior
Success	Success	0000	The SCP has successfully stored the SOP Instance.
Refused	Out of Resources	A700-A7FF	The Association is aborted using A-ABORT and the failure is reported to the user. This is a transient failure.
Error	Data Set does not match SOP Class	A900-A9FF	The Association is aborted using A-ABORT and the failure is reported to the user.
Error	Cannot Understand	C000-CFFF	The Association is aborted using A-ABORT and the failure is reported to the user.
Warning	Coercion of Data Elements	B000	Image transmission is considered successful therefore the warning is not reported to the user.
Warning	Data Set does not match SOP Class	B007	Image transmission is considered successful therefore the warning is not reported to the user.
Warning	Elements Discarded	B006	Image transmission is considered successful therefore the warning is not reported to the user.
Warning ¹	Attribute list error	0107	Object instance transmission is considered successful therefore the warning is not reported to the user.
Warning ¹	Attribute Value Out of Range	0116	Object instance transmission is considered successful therefore the warning is not reported to the user.
*	*	Any other status code.	The Association is aborted using A-ABORT and the failure is reported to the user.

Note: 1. Image storage excepted.

The behavior of Storage AE during communication failure is summarized in the Table below:

Table 4.2-28 STORAGE COMMUNICATION FAILURE BEHAVIOR

Exception	Behavior
Timeout	The Association is aborted using A-ABORT and the timeout error is reported to the user.
Association aborted by the SCP or network layers	The association failure is reported to the user.

A failed storage association may automatically be restarted if so configured.

4.2.3 Hardcopy Application Entity Specification

4.2.3.1 SOP Classes

This product provides Standard Conformance to the following SOP Classes:

Table 4.2-29 SOP CLASSES FOR AE HARDCOPY

SOP Class Name	SOP Class UID	SCU	SCP
Basic Grayscale Print Management Meta	1.2.840.10008.5.1.1.9	Yes	No
Basic Color Print Management Meta	1.2.840.10008.5.1.1.18	Yes	No
Verification	1.2.840.10008.1.1	Yes	No

4.2.3.2 Association Policies

4.2.3.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed:

Table 4.2-30 DICOM APPLICATION CONTEXT FOR AE HARDCOPY

Application Context Name	1.2.840.10008.3.1.1.1
--------------------------	-----------------------

4.2.3.2.2 Number of Associations

This product initiates one Association at a time for each configured hardcopy device.

Table 4.2-31 NUMBER OF ASSOCIATIONS INITIATED FOR AE HARDCOPY

Maximum number of simultaneous Associations	1
---	---

4.2.3.2.3 Asynchronous Nature

This product does not support asynchronous communication (multiple outstanding transactions over a single Association).

**Table 4.2-32
ASYNCHRONOUS NATURE AS A SCU FOR AE HARDCOPY**

Maximum number of outstanding asynchronous transactions	1
---	---

4.2.3.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

Table 4.2-33 DICOM IMPLEMENTATION CLASS AND VERSION FOR AE HARDCOPY

Implementation Class UID	1.2.392.200039.117
Implementation Version Name	20180308

4.2.3.3 Association Initiation Policy**4.2.3.3.1 Activity – Print Images****4.2.3.3.1.1 Description and Sequencing of Activities**

When a user selects images and requests to print them, the images are sent to the PRINTER_QUEUE folder. The virtual film sheets are composed according to the pre-defined film format. The film sheets are requests to be sent to a specific hardcopy device. The user can select the desired film format, number of copies, and other printing conditions in the Print Property GUI.




Figure 4.2-5
SEQUENCING OF ACTIVITY – FILM IMAGES

A typical sequence of DIMSE messages sent over an association between Hardcopy AE and a Printer is illustrated in Figure 4.2-5:

1. Hardcopy AE opens an association with the Printer
2. N-GET on the Printer SOP Class is used to obtain current printer status information. If the Printer reports a status of FAILURE, the print-job is switched to a failed state and the user informed.
3. N-CREATE on the Film Session SOP Class creates a Film Session.
4. N-CREATE on the Film Box SOP Class creates a Film Box linked to the Film Session. Pre-configured number of Image Boxes will be created as the result of this operation.

5. Each N-SET on Image Box SOP Class transfers the requested image to the printer at the successive position on the film sheet. The Hardcopy does not support the Presentation LUT SOP Class.
6. N-ACTION on the Film Box SOP Class instructs the printer to print the Film Box
7. The printer prints the requested number of film sheets. The sequence 4 through 6 may be repeated when the single film sheet is not enough to contain the requested images.
8. The Printer asynchronously reports its status via N-EVENT-REPORT notification (Printer SOP Class). The printer can send this message at any time. Hardcopy AE does not require the N-EVENT-REPORT to be sent. Hardcopy AE is capable of receiving an N-EVENT-REPORT notification at any time during an association. If the Printer reports a status of FAILURE, the Print Session is terminated and the user informed.
9. N-DELETE on the Film Session SOP Class deletes the complete Film Session SOP Instance hierarchy.
10. Hardcopy AE closes the association with the Printer.

If any Response from the remote Application contains a status other than Success or Warning, the Association is aborted and the related Film Session is terminated and the status is user informed.

4.2.3.3.1.2 Proposed Presentation Contexts

This product is capable of proposing the Presentation Contexts shown in the Table below:

Table 4.2-34 PROPOSED PRESENTATION CONTEXTS FOR ACTIVITY PRINT IMAGES

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Ext. Neg.
Name	UID	Name List	UID List		
Basic Grayscale Print Management Meta	1.2.840.10008.5.1.1.9	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
Basic Color Print Management Meta	1.2.840.10008.5.1.1.18	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
Verification	1.2.840.10008.1.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		

4.2.3.3.1.3 Common SOP Specific Conformance for all Print SOP Classes

The general behavior of Hardcopy AE during communication failure is summarized in the Table below. This behavior is common for all SOP Classes supported by Hardcopy AE.

**Table 4.2-35
HARDCOPY COMMUNICATION FAILURE BEHAVIOR**

Exception	Behavior
Timeout	The Association is aborted using A-ABORT and the print-job is terminated. The reason is reported to the user.
Association aborted by the SCP or network layers	The print-job is terminated and the print-job is terminated. The reason is reported to the user.

4.2.3.3.1.4 SOP Specific Conformance for the Printer SOP Class

Hardcopy AE supports the following DIMSE operations and notifications for the Printer SOP Class:

- N-GET
- N-EVENT-REPORT

Details of the supported attributes and status handling behavior are described in the following subsections.

4.2.3.3.1.4.1 Printer SOP Class Operations (N-GET)

Hardcopy AE uses the Printer SOP Class N-GET operation to obtain information about the current printer status. The attributes obtained via N-GET are listed in the Table below:

Table 4.2-36
PRINTER SOP CLASS N-GET REQUEST ATTRIBUTES

Attribute Name	Tag	VR	Value	Presence of Value	Source
Manufacturer	(0008,0070)	LO	Provided by Printer	ANAP	Printer
Manufacturer's Model Name	(0008,1090)	LO	Provided by Printer	ANAP	Printer
Device Serial Number	(0018,1000)	LO	Provided by Printer	ANAP	Printer
Software Versions	(0018,1020)	LO	Provided by Printer	ANAP	Printer
Printer Status	(2110,0010)	CS	Provided by Printer	ALWAYS	Printer
Printer Status Info	(2110,0020)	CS	Provided by Printer	ALWAYS	Printer
Printer Name	(2110,0030)	LO	Provided by Printer	ANAP	Printer

The Printer Status information is evaluated as follows:

1. If Printer status (2110,0010) is FAILURE, the Hardcopy AE is terminated and status is user informed.
2. If Printer status (2110,0010) is NORMAL or WARNING, the Hardcopy AE continues to print.

The behavior of Hardcopy AE when encountering status codes in a N-GET response is summarized in the Table below:

Table 4.2-37
PRINTER SOP CLASS N-GET RESPONSE STATUS HANDLING BEHAVIOR

Service Status	Further Meaning	Error Code	Behavior
Success	Success	0000	The request to get printer status information was success.
*	*	Any other status code.	The Association is aborted using A-ABORT and the print-job is terminated. The status is reported to the user.

4.2.3.3.1.4.2 Printer SOP Class Notifications (N-EVENT-REPORT)

Hardcopy AE is capable of receiving an N-EVENT-REPORT request at any time during an association.

The behavior of Hardcopy AE when receiving Event Types within the N-EVENT-REPORT is summarized in the Table below:

Table 4.2-38 PRINTER SOP CLASS N-EVENT-REPORT BEHAVIOUR

Event Type Name	Event Type ID	Behavior
Normal	1	The print session continues to be printed.
Warning	2	The print session continues to be printed. The Warning status is not reported to user.
Failure	3	The print session is terminated. The Failure is reported to user.
*	*	The print session is terminated. The Failure is reported to user.

The reasons for returning specific status codes in an N-EVENT-REPORT response are summarized in the Table below:

**Table 4.2-39
PRINTER SOP CLASS N-EVENT-REPORT RESPONSE STATUS REASONS**

Service Status	Further Meaning	Error Code	Reasons
Success	Success	0000	The notification event has been successfully received.
Failure	No Such Event Type	0113H	An invalid Event Type ID was supplied in the N-EVENT-REPORT request.
Failure	Processing Failure	0110H	An internal error occurred during processing of the N-EVENT-REPORT. A short description of the error will be returned in Error Comment (0000,0902).

4.2.3.3.1.5 SOP Specific Conformance for the Film Session SOP Class

Hardcopy AE supports the following DIMSE operations for the Film Session SOP Class:

- N-CREATE
- N-DELETE

Details of the supported attributes and status handling behavior are described in the following subsections.

4.2.3.3.1.5.1 Film Session SOP Class Operations (N-CREATE)

The attributes supplied in an N-CREATE Request are listed in the Table below. The values are typical and may be configured by the CSE:

Table 4.2-40 FILM SESSION SOP CLASS N-CREATE REQUEST ATTRIBUTES

Attribute Name	Tag	VR	Value	Presence of Value	Source
Number of Copies	(2000,0010)	IS	1 .. 9	ALWAYS	CONFIG
Print Priority	(2000,0020)	CS	HIGH, MED, or LOW	ALWAYS	CONFIG
Medium Type	(2000,0030)	CS	BLUE FILM, CLEAR FILM, or PAPER	ALWAYS	CONFIG
Film Destination	(2000,0040)	CS	MAGAZINE or PROCESSOR	ALWAYS	CONFIG
Memory Allocation	(2000,0060)	IS		ANAP	CONFIG

The behavior of Hardcopy AE when encountering status codes in an N-CREATE response is summarized in the Table below:

**Table 4.2-41
FILM SESSION SOP CLASS N-CREATE RESPONSE STATUS HANDLING BEHAVIOR**

Service Status	Further Meaning	Error Code	Behavior
Success	Success	0000	The SCP has completed the operation successfully.
Warning	Attribute Value Out of Range	0116H	The N-CREATE operation is considered successful and the status is not reported to the user.
Warning	Attribute List Error	0107H	The N-CREATE operation is considered successful and the status is not reported to the user.
*	*	Any other status code.	The Association is aborted using A-ABORT and the print session is terminated. The status is reported to the user.

4.2.3.3.1.5.2 Film Session SOP Class Operations (N-DELETE)

The behavior of Hardcopy AE when encountering status codes in an N-DELETE response is summarized in the Table below:

**Table 4.2-42
FILM SESSION SOP CLASS N-DELETE RESPONSE STATUS HANDLING BEHAVIOR**

Service Status	Further Meaning	Error Code	Behavior
Success	Success	0000	The SCP has completed the operation successfully.
*	*	Any other status code.	The Association is aborted using A-ABORT and the print session is terminated.

4.2.3.3.1.6 SOP Specific Conformance for the Film Box SOP Class

Hardcopy AE supports the following DIMSE operations for the Film Box SOP Class:

- N-CREATE
- N-ACTION

Details of the supported attributes and status handling behavior are described in the following subsections.

4.2.3.3.1.6.1 Film Box SOP Class Operations (N-CREATE)

The attributes supplied in an N-CREATE Request are listed in the Table below. The values are typical and may be configured by the CSE:

Table 4.2-43 FILM BOX SOP CLASS N-CREATE REQUEST ATTRIBUTES

Attribute Name	Tag	VR	Value	Presence of Value	Source
Image Display Format	(2010,0010)	ST	STANDARD\m,n	ALWAYS	CONFIG
Film Orientation	(2010,0040)	CS	PORTRAIT or LANDSCAPE	ALWAYS	CONFIG
Film Size ID	(2010,0050)	CS	14INX17IN,14INX14IN, 11INX14IN,11INX11IN, 85INX11IN, 8INX10IN	ALWAYS	CONFIG

Magnification Type	(2010,0060)	CS	REPLICATE, BILINEAR, CUBIC or NONE	ALWAYS	CONFIG
Smoothing Type	(2010,0080)	CS		ANAP	CONFIG
Border Density	(2010,0100)	CS	BLACK or WHITE	ALWAYS	CONFIG
Empty Image Density	(2010,0110)	CS	BLACK or WHITE	ALWAYS	CONFIG
Min Density	(2010,0120)	US	0 .. 329	ALWAYS	CONFIG
Max Density	(2010,0130)	US	1 .. 330	ALWAYS	CONFIG
Trim	(2010,0140)	CS	YES or NO	ANAP	CONFIG
Configuration Information	(2010,0150)	ST	Set if requested by Printer	ANAP	CONFIG
Referenced Film Session Sequence	(2010,0500)	SQ		ALWAYS	AUTO
>Referenced SOP Class UID	(0008,1150)	UI	1.2.840.10008.5.1.1.1	ALWAYS	AUTO
>Referenced SOP Instance UID	(0008,1155)	UI	From created Film Session SOP Instance	ALWAYS	AUTO

The behavior of Hardcopy AE when encountering status codes in an N-CREATE response is summarized in the Table below:

Table 4.44 FILM BOX SOP CLASS N-CREATE RESPONSE STATUS HANDLING BEHAVIOR

Service Status	Further Meaning	Error Code	Behavior
Success	Success	0000	The SCP has completed the operation successfully.
Warning	Requested Min Density or Max Density outside of printer's operating range	B605H	The N-CREATE operation is considered successful and the status is not reported to the user.
*	*	Any other status code.	The Association is aborted using A-ABORT and the print-job is terminated. The status is reported to the user.

4.2.3.3.1.6.2 Film Box SOP Class Operations (N-ACTION)

An N-ACTION Request is issued to instruct the Print SCP to print the contents of the Film Box. The Action Reply argument in an N-ACTION response is not evaluated.

The behavior of Hardcopy AE when encountering status codes in an N-ACTION response is summarized in the Table below:

**Table 4.2-45
FILM BOX SOP CLASS N-ACTION RESPONSE STATUS HANDLING BEHAVIOR**

Service Status	Further Meaning	Error Code	Behavior
Success	Success	0000	The SCP has completed the operation successfully. The film has been accepted for printing.
Warning	Film Box SOP Instance hierarchy does not contain Image Box SOP Instances (empty	B603H	The Association is aborted using A-ABORT and the print-job is

	page)		terminated. The status is reported to the user.
Warning	Image size is larger than Image Box size. The image has been demagnified.	B604H	The N-ACTION operation is considered successful and the status is not reported to the user.
Warning	Image size is larger than Image Box size. The image has been cropped to fit.	B609H	The N-ACTION operation is considered successful and the status is not reported to the user.
Warning	Image size or Combined Print Image Size is larger than Image Box size. The image or combined Print Image has been decimated to fit.	B60AH	The N-ACTION operation is considered successful and the status is not reported to the user.
Failure	Unable to create Print Job SOP Instance; print queue is full.	C602	The Association is aborted using A-ABORT and the status is reported to the user.
Failure	Image size is larger than Image Box size.	C603	The Association is aborted using A-ABORT and the status is reported to the user.
Failure	Combined Print Image Size is larger than Image Box size.	C613	The Association is aborted using A-ABORT and the status is reported to the user.
*	*	Any other status code.	The Association is aborted using A-ABORT and the status is reported to the user.

4.2.3.3.1.7 SOP Specific Conformance for the Image Box SOP Class

Hardcopy AE supports the following DIMSE operations for the Basic Grayscale and Basic Color Image Box SOP Classes:

— N-SET

Details of the supported attributes and status handling behavior are described in the following subsections.

4.2.3.3.1.7.1 Basic Grayscale Image Box SOP Class Operations (N-SET)

The attributes supplied in an N-SET Grayscale Image Box Request are listed in the Table below:

Table 4.2-46 BASIC GRAYSCALE IMAGE BOX SOP CLASS N-SET REQUEST ATTRIBUTES

Attribute Name	Tag	VR	Value	Presence of Value	Source
Image Position	(2020,0010)	US	1 to maximum image position allowed for the Image Display Format (2010, 0010)	ALWAYS	AUTO
Polarity	(2020,0020)	CS	NORMAL or REVERSE	ALWAYS	CONFIG
Basic Grayscale Image Sequence	(2020,0110)	SQ		ALWAYS	AUTO
>Samples Per Pixel	(0028,0002)	US	1	ALWAYS	AUTO

>Photometric Interpretation	(0028,0004)	CS	MONOCHROME2	ALWAYS	AUTO
>Rows	(0028,0010)	US	Copied from source image	ALWAYS	AUTO
>Columns	(0028,0011)	US	Copied from source image	ALWAYS	AUTO
>Pixel Aspect Ratio	(0028,0034)	IS	Copied from source image	ANAP	AUTO
>Bits Allocated	(0028,0100)	US	8	ALWAYS	AUTO
>Bits Stored	(0028,0101)	US	8	ALWAYS	AUTO
>High Bit	(0028,0102)	US	7	ALWAYS	AUTO
>Pixel Representation	(0028,0103)	US	0	ALWAYS	AUTO
>Pixel Data	(7FE0,0010)	OW	Pixels from source image	ALWAYS	AUTO

The behavior of Hardcopy AE when encountering status codes in an N-SET response is summarized in Table 4.2-48.

4.2.3.3.1.7.2 Basic Color Image Box SOP Class Operations (N-SET)

The attributes supplied in an N-SET Color Image Box Request are listed in the Table below:

Table 4.2-47 BASIC COLOR IMAGE BOX SOP CLASS N-SET REQUEST ATTRIBUTES

Attribute Name	Tag	VR	Value	Presence of Value	Source
Image Position	(2020,0010)	US	1 to maximum image position allowed for the Image Display Format (2010, 0010)	ALWAYS	AUTO
Polarity	(2020,0020)	CS	Printer may ignore the value	ALWAYS	AUTO
Preformatted Color Image Sequence	(2020,0111)	SQ		ALWAYS	AUTO
>Samples Per Pixel	(0028,0002)	US	3	ALWAYS	AUTO
>Photometric Interpretation	(0028,0004)	CS	RGB	ALWAYS	AUTO
>Planar Configuration	(0028,0006)	US	0 or 1 (default is 1)	ALWAYS	CONFIG ^[1]
>Rows	(0028,0010)	US	Copied from source image	ALWAYS	AUTO
>Columns	(0028,0011)	US	Copied from source image	ALWAYS	AUTO
>Pixel Aspect Ratio	(0028,0034)	IS	Copied from source image	ANAP	AUTO
>Bits Allocated	(0028,0100)	US	8	ALWAYS	AUTO
>Bits Stored	(0028,0101)	US	8	ALWAYS	AUTO
>High Bit	(0028,0102)	US	7	ALWAYS	AUTO
>Pixel Representation	(0028,0103)	US	0	ALWAYS	AUTO
>Pixel Data	(7FE0,0010)	OW	Pixels from source image	ALWAYS	AUTO

^[1] The Planar Configuration could be configured by the CSE referencing the DICOM Conformance Statement of the remote printer.

The behavior of Hardcopy AE when encountering status codes in an N-SET response is summarized in Table 4.2-48.

Table 4.2-48 IMAGE BOX SOP CLASSES N-SET RESPONSE STATUS HANDLING BEHAVIOR

Service Status	Further Meaning	Error Code	Behavior
Success	Success	0000	The SCP has completed the operation successfully. Image successfully stored in Image Box.
Warning	Image size is larger than Image Box size. The image has been demagnified.	B604	The N-SET operation is considered successful and the status is not reported.
Warning	Requested Min Density or Max Density outside of printer's operating range.	B605	The N-SET operation is considered successful and the status is not reported.
Warning	Image size is larger than Image Box size. The image has been cropped to fit.	B609	The N-SET operation is considered successful and the status is not reported.
Warning	Image size or Combined Print Image Size is larger than Image Box size. The image or combined Print Image has been decimated to fit.	B60A	The N-SET operation is considered successful and the status is not reported.
Failure	Image size is larger than Image Box size.	C603	The Association is aborted using A-ABORT and the status is reported to the user.
Failure	Insufficient memory in printer to store the image.	C605	The Association is aborted using A-ABORT and the status is reported to the user.
Failure	Combined Print Image Size is larger than Image Box size.	C613	The Association is aborted using A-ABORT and the status is reported to the user.
*	*	Any other status code.	The Association is aborted using A-ABORT and the status is reported to the user.

4.2.3.4 Association Acceptance Policy

The Hardcopy Application Entity does not accept Associations.

4.2.4 Query/Retrieve Application Entity Specification

4.2.4.1 SOP Classes

This product provides Standard Conformance to the following SOP Classes:

Table 4.2-49 SOP CLASSES FOR AE QUERY/RETRIEVE

SOP Class Name	SOP Class UID	SCU	SCP
Study Root Query/Retrieve Information Model – FIND	1.2.840.10008.5.1.4.1.2.2.1	Yes	No
Study Root Query/Retrieve Information Model – MOVE	1.2.840.10008.5.1.4.1.2.2.2	Yes	No
Ultrasound Image Storage	1.2.840.10008.5.1.4.1.1.6.1	No	Yes
Ultrasound Multi-frame Image Storage	1.2.840.10008.5.1.4.1.1.3.1	No	Yes
CT Image Storage	1.2.840.10008.5.1.4.1.1.2	No	Yes
MR Image Storage	1.2.840.10008.5.1.4.1.1.4	No	Yes
Positron Emission Tomography Image Storage	1.2.840.10008.5.1.4.1.1.128	No	Yes
Encapsulated PDF Storage	1.2.840.10008.5.1.4.1.1.104.1	No	Yes

4.2.4.2 Association Policies

4.2.4.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed:

Table 4.2-50 DICOM APPLICATION CONTEXT FOR AE QUERY/RETRIEVE

Application Context Name	1.2.840.10008.3.1.1.1
--------------------------	-----------------------

4.2.4.2.2 Number of Associations

This product initiates one Association at a time for a Workflow requests.

Table 4.2-51 NUMBER OF ASSOCIATIONS INITIATED FOR AE QUERY/RETRIEVE

Maximum number of simultaneous Associations	1
---	---

4.2.4.2.3 Asynchronous Nature

This product does not support asynchronous communication (multiple outstanding transactions over a single Association).

Table 4.2-52 ASYNCHRONOUS NATURE AS A SCU FOR AE QUERY/RETRIEVE

Maximum number of outstanding asynchronous transactions	1
---	---

4.2.4.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

Table 4.2-53 DICOM IMPLEMENTATION CLASS AND VERSION FOR AE QUERY/RETRIEVE

Implementation Class UID	1.2.276.0.7230010.3.0.3.5.4
Implementation Version Name	OFFIS_DCMTK_354 (subject to change without notice)

4.2.4.3 Association Initiation Policy

4.2.4.3.1 Activity – Query Image Information and Retrieve Images

4.2.4.3.1.1 Description and Sequencing of Activities

An interactive query for Image Information is initiated by clicking the "Search" button in the "Import" Registration GUI. In the window, Patient ID, Patient's Name, Accession Number, Study Date, and/or Modalities in Study (US, CT, MR, PET, or DOC) may be supplied by the operator as the query keys.

When the "Search" button is clicked by the operator, this product sends the C-FIND request several times with changing Query/Retrieve Level to "STUDY", "SERIES" or "IMAGE". After retrieval of all responses, this product displays the result in a separate list. During receiving the response, the query processing is canceled by issuing a C-FIND-CANCEL if the operator clicks the "Cancel" button.

When the operator selects the instance information items in the list and clicks "Preview" or "Import" button, this product sends the C-MOVE requests and waits for receiving C-Store requests from the image manager. After retrieval of all C-MOVE responses, this product displays the received instances and copies them to local HDD. During receiving the C-MOVE response, the processing is canceled by issuing a C-MOVE-CANCEL if the operator clicks the "Cancel" button.




Figure 4.2-5
SEQUENCING OF ACTIVITY – QUERY INSTANCE INFORMATION AND RETRIEVE INSTANCE

A possible sequence of interactions between the Query/Retrieve AE and an Image Manager is illustrated in Figure 4.2-5:

1. The operator requests “Query Instance Information”.
2. The Query/Retrieve AE (SCU) opens the association with the Image Manager (SCP).
3. The Query/Retrieve AE (SCU) sends a C-FIND request to the Image Manager (SCP) containing the Instance Information Query attributes.
4. The Image Manager (SCP) returns a C-FIND response containing the requested attributes of the first matching Instance Information Item.
5. The Image Manager (SCP) returns a C-FIND response containing the requested attributes of the second matching Instance Information Item.

6. The Image Manager (SCP) returns another C-FIND response with status Success indicating that no further matching Instance Information Items exist. This example assumes that only 2 Instance Information items match the Query.
7. The Query/Retrieve AE (SCU) closes the association with the Image Manager.
8. The operator requests "Retrieve Image".
9. The Query/Retrieve AE (SCU) opens the association with the Image Manager (SCP).
10. The Query/Retrieve AE (SCU) sends a C-MOVE request to the Image Manager (SCP).
11. The Image Manager (SCP) opens the association with the Query/Retrieve AE (SCP).
12. The Image Manager (SCP) transfers instance to the Query/Retrieve AE (SCP) using a C-STORE request.
13. The Query/Retrieve AE (SCP) replies with a C-STORE response (status success).
14. The Image Manager (SCP) sends C-MOVE response to the Query/Retrieve AE (SCU) with status Pending.
15. The Image Manager (SCP) closes the association with the Query/Retrieve AE (SCU).
16. The Image Manager (SCP) sends C-MOVE response to the Query/Retrieve AE (SCU) with status Success indicating that no further instances should be transferred.
17. The Query/Retrieve AE (SCU) closes the association with the Image Manager (SCP).

4.2.4.3.1.2 Proposed Presentation Contexts

This product will propose Presentation Contexts as shown in the following table:

Table 4.2-54
PROPOSED PRESENTATION CONTEXTS FOR ACTIVITY QUERY IMAGE INFORMATION

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Ext. Neg.
Name	UID	Name List	UID List		
Study Root Query/Retrieve Information Model - FIND	1.2.840.10008.5.1.4.1.2.2.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
Study Root Query/Retrieve Information Model - MOVE	1.2.840.10008.5.1.4.1.2.2.2	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

**Table 4.2-55
PROPOSED PRESENTATION CONTEXTS FOR ACTIVITY RETRIEVE IMAGES**

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Ext. Neg.
Name	UID	Name List	UID List		
Ultrasound Image Storage	1.2.840.10008.5.1.4.1.1.6.1	Implicit VR Little Endian Explicit VR Little Endian JPEG Baseline Compression RLE Compression	1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.5	SCP	None
Ultrasound Multi-frame Image Storage	1.2.840.10008.5.1.4.1.1.3.1	Implicit VR Little Endian Explicit VR Little Endian JPEG Baseline Compression ¹ RLE Compression	1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.5	SCP	None
CT Image Storage	1.2.840.10008.5.1.4.1.1.2	Implicit VR Little Endian Explicit VR Little Endian JPEG Lossless, Nonhierarchical, First- Order Prediction	1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.4.70	SCP	None
MR Image Storage	1.2.840.10008.5.1.4.1.1.4	Implicit VR Little Endian Explicit VR Little Endian JPEG Lossless, Nonhierarchical, First- Order Prediction	1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.4.70	SCP	None
Positron Emission Tomography Image Storage	1.2.840.10008.5.1.4.1.1.128	Implicit VR Little Endian Explicit VR Little Endian JPEG Lossless, Nonhierarchical, First- Order Prediction	1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.4.70	SCP	None
Encapsulated PDF Storage	1.2.840.10008.5.1.4.1.1.104.1	Implicit VR Little Endian Explicit VR Little Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.1	SCP	None

4.2.4.3.1.3 SOP Specific Conformance for Query/Retrieve SOP Classes

The behavior of Query/Retrieve AE when encountering status codes in C-FIND response is summarized in the table below.

**Table 4.2-56 QUERY INSTANCE INFORMATION C-FIND RESPONSE STATUS
HANDLING BEHAVIOR**

Service Status	Further Meaning	Error Code	Behavior
Success	Matching is complete	0000	The SCP has completed the matches. Instance Information items are available for display or further processing.
Cancel	Matching terminated due to	FE00	If the query was cancelled due to user's operation then the SCP has completed the matches. Instance Information

	Cancel request		items are available for display or further processing. Otherwise, the Association is aborted using A-ABORT and the Image Information query is failed.
Failed	Out of Resources	A700	The Association is aborted using A-ABORT and the Query Instance Information is failed.
Failed	Identifier does not match SOP Class	A900	The Association is aborted using A-ABORT and the Query Instance Information is failed.
Failed	Unable to Process	C000 – CFFF	The Association is aborted using A-ABORT and the Query Instance Information is failed.
Pending	Matches are continuing	FF00	The Instance Information item contained in the Identifier is collected for later display or further processing.
Pending	Matches are continuing – Warning that one or more Optional Keys were not supported	FF01	The Instance Information item contained in the Identifier is collected for later display or further processing.
Failed	*	Any other status code.	The Association is aborted using A-ABORT and the Query Instance Information is failed.

The behavior of Query/Retrieve AE when encountering status codes in C-MOVE response is summarized in the table below.

Table 4.2-57 RETRIEVE INSTANCE INFORMATION C-MOVE RESPONSE STATUS HANDLING BEHAVIOR

Service Status	Further Meaning	Error Code	Behavior
Success	Sub-operations Complete – No Failures	0000	The SCP has completed the sub-operations. All requested instances are available for display or further processing.
Warning	Sub-operations Complete – One or more Failures	B000	The Association is aborted using A-ABORT and Retrieve Instances is failed.
Cancel	Sub-operations terminated due to Cancel Indication	FE00	The Association is aborted using A-ABORT and Retrieve Instances is failed.
Failed	Unable to calculate number of matches	A701	The Association is aborted using A-ABORT and Retrieve Instances is failed.
Failed	Move Destination unknown	A801	The Association is aborted using A-ABORT and Retrieve Instances is failed.
Failed	Unable to perform sub-operations	A702	The Association is aborted using A-ABORT and Retrieve Instances is failed.
Failed	Unable to Process	C000 – CFFF	The Association is aborted using A-ABORT and Retrieve Instances is failed.
Pending	Sub-operations are continuing	FF00	It succeeded to receive an instances and the retrieve process is continuing.

Failed	*	Any other status code.	The Association is aborted using A-ABORT and Retrieve Instances is failed.
--------	---	------------------------	--

The behavior of this product during communication failure is summarized in the Table below.

Table 4.2-58 QUERY/RETRIEVE COMMUNICATION FAILURE BEHAVIOR

Exception	Behavior
Timeout	The Association is aborted using A-ABORT and the process is failed. An error is reported to the user.
Association aborted by the SCP or network layers	The process is failed. An error is reported to the user.

The table below provides a description of the Study Root Query/Retrieve Information Model. Unexpected attributes returned in a C-FIND response are ignored.

Table 4.2-59 STUDY ROOT QUERY/RETRIEVE INFORMATION MODEL

Attribute Name	Tag	VR	M	R	Q	D
Study Level Keys						
Specific Character Set	(0008,0005)	CS		x		
Study Date	(0008,0020)	DA	R	x	x	x
Study Time	(0008,0030)	TM		x	x	x
Accession Number	(0008,0050)	SH	*	x	x	x
Modalities In Study	(0008,0061)	CS		x		
Referring Physicians Name	(0008,0090)	PN		x		x
Study Description	(0008,1030)	LO		x		x
Patient's Name	(0010,0010)	PN	*	x	x	x
Patient ID	(0010,0020)	LO	*	x	x	x
Patient's Birth Date	(0010,0030)	DA		x		x
Patient's Sex	(0010,0040)	CS		x		x
Study ID	(0020,0010)	SH		x		x
Study Instance UID	(0020,000D)	UI		x	x	
Requesting Service	(0032,1033)	LO		x		
Series Level Keys						
Series Date	(0008,0021)	DA		x		
Series Time	(0008,0031)	TM		x		x
Modality	(0008,0060)	SH	S	x	x	x
Series Description	(0008,103E)	LO		x		x
Body Part Examined	(0018,0015)	CS		x		x
Series Instance UID	(0020,000E)	UI		x	x	
Series Number	(0020,0011)	IS		x		x

Image Level Keys						
SOP Instance UID	(0008,0018)	CS		x	x	
Instance Number	(0020,0013)	IS		x		x

The above table should be read as follows:

- Attribute Name: Attributes supported to build an Instance Information Request Identifier of this product.
- Tag: DICOM tag for this attribute.
- VR: DICOM VR for this attribute.
- M: Matching keys for Query Instance Information. An "S" will indicate that this product will supply an attribute value for Single Value Matching, an "R" will indicate Range Matching and a "*" will denote wildcard matching.
- R: Return keys. An "x" will indicate that this product will supply this attribute as Return Key with zero length for Universal Matching.
- Q: Interactive Query Key. An "x" will indicate that this product will supply this attribute as matching key, if entered in the Search item or selected item on a Import window.
- D: Displayed keys. An "x" indicates that this attribute is displayed to the user during a Import window.

4.3 NETWORK INTERFACES

4.3.1 Physical Network Interface

This product supports both wired and wireless network interface. The wired network interface is always used if wired and wireless network interfaces are connected at the same time. It is necessary to assign IP Address on each network interface.

Table 4.3-1
SUPPORTED PHYSICAL NETWORK INTERFACES

Ethernet 100baseTX or 10baseT is automatically detected
Wireless 802.11ac/a/n/g/b

4.3.2 Additional Protocols

This product may be configured to get the local configuration via the DHCP and to synchronize the system time with the NTP server. However it does not conforms to other System Management Profiles as DNS nor LDAP.

4.3.3 IPv4 and IPv6 Support

This product supports both IPv4 and IPv6 ¹ connections.

Note: 1. Except for Storage Application Entity for Encapsulated PDF Storage.

4.4 CONFIGURATION

4.4.1 AE Title/Presentation Address Mapping

4.4.1.1 Local AE Titles

All local applications share the same AE Titles and TCP/IP Address configured at the Local AET registry in the Network Setting and DICOM-Common GUI. The AE Title, Station Name, IP Address, and Port Number must be configured during installation, since no default values are provided and not automatically configured.

Table 4.4-1 AE TITLE CONFIGURATION TABLE

Application Entity	Default AE Title	Default TCP/IP Port
Storage	No Default	104
Workflow	same as Storage	104
Hardcopy	same as Storage	104
Query/Retrieve	same as Storage	2350

4.4.1.2 Remote AE Title/Presentation Address Mapping

The AE Titles, Station Names, IP Addresses, and Port numbers of remote applications are configured through the DICOM-Server/Worklist, DICOM-QR, DICOM-MPPS/Commitment, DICOM-SR, DICOM-Printer and DICOM-PDF GUI.

Some characters allowed for VR AE may not be used to represent the local and remote AE Titles as shown in the table below:

Table 4.4-2 FORBIDDEN CHARACTERS FOR AE TITLE

Forbidden characters for AE Title	; : " < > * \ ? ,
-----------------------------------	---------------------

4.4.1.2.1 Workflow

The Application Entity Title, Station Name, IP Address, and the Port Number of the remote Modality Worklist SCP are registered at the Worklist row of DICOM-Server/Worklist GUI. Optionally the Application Entity Titles, Station Names, IP Addresses, and the Port numbers of remote MPPS server and Storage Commitment (Image Manager) should be registered in the DICOM-MPPS/Commitment GUI before each service is used.

Only one Storage Commitment server can be activated although the destinations of image storage and SR storage may be registered separately.

4.4.1.2.2 Storage

The Application Entity Titles, Station Names, IP Addresses, and the Port Numbers of the remote Storage SCPs receiving images are registered at DICOM-Server/Worklist GUI. The remote Storage SCPs are able to register to five. Multiple storage destinations for single frame and/or multi-frame images can be activated by selecting the check box in the GUI.

The Application Entity Titles, Station Names, IP Addresses, and the Port Numbers of the remote Storage SCPs receiving Structured Reports are registered at DICOM-SR GUI. The destinations of Structured Report storage may not be same as the destinations of other storage.

The Application Entity Titles, Station Names, IP Addresses, and the Port Numbers of the remote Storage SCPs receiving PDF Reports are registered at DICOM-PDF GUI. The destinations of PDF Report storage may not be same as the destinations of other storage.

4.4.1.2.3 Hardcopy

The Application Entity Titles, Station Names, IP Addresses, and the Port Numbers of the remote Hardcopy SCPs are registered in the DICOM-Printer GUI. Although the Multiple Hardcopy destinations may be registered, only one destination must be activated by selected the radio button in the GUI. The same Application Entity may be registered to configure different Hardcopy settings for selections by the user's preference.

4.4.1.2.4 Query/Retrieve

The Application Entity Titles, Station Names, IP Addresses, and the Port Numbers of the remote Query/Retrieve SCPs receiving images are registered at DICOM-QR GUI. The remote Query/Retrieve SCPs are able to register to five. The local Port Number of receiving instances from the Image Manager(SCU) is also registered in the GUI

4.4.2 Parameters

A large number of parameters related to acquisition and general operation can be configured using the DICOM configuration user interface. The Table below shows those configuration parameters relevant to DICOM communication.

Table 4.4-3 CONFIGURATION PARAMETERS TABLE

Parameter	Configurable (Yes/No)	Default Value
General Parameters		
Max PDU Receive Size	No	16K bytes
Max PDU Send Size (larger PDUs will never be sent, even if the receiver supports a larger Max PDU Receive Size. If the receiver supports a smaller Max PDU Receive Size then the Max PDU Send Size will be reduced accordingly for the duration of the Association. Max PDU Receive Size information is exchanged during DICOM Association Negotiation in the Maximum Length Sub-Item of the A-ASSOCIATION-RQ and A-ASSOCIATE-AC)	No	16K Bytes
Time-out waiting for an acceptance or rejection response to an Association Request (Application Level Timeout)	No	INFINITE
Time-out waiting for a response to an Association release request (Application Level Timeout)	No	INFINITE
Time-out waiting for completion of a TCP/IP connect request (Low-level timeout)	No	20 s
Time-out waiting a Response to a DIMSE Request (Low-Level Timeout)	Yes	30 s
Time-out for waiting for data between TCP/IP-packets (Low Level Timeout)	Yes	30 s
Modality Worklist Parameters		
Modality Worklist SCU time-out waiting for the final response to a C-FIND-RQ	Yes	30 s
Maximum number of Worklist Items	No	500
Supported Transfer Syntaxes for Modality Worklist	Yes	Implicit VR Little Endian Explicit VR Little Endian

Parameter	Configurable (Yes/No)	Default Value
Query Worklist for specific Scheduled Station AE Title	Yes	No default
Query Worklist for specific Modality Value	No	US
MPPS Parameters		
MPPS SCU time-out waiting for a response to a N-CREATE-RQ	Yes	30 s
MPPS SCU time-out waiting for a response to a N-SET-RQ	Yes	30 s
Supported Transfer Syntaxes for MPPS	Yes	Implicit VR Little Endian Explicit VR Little Endian
Storage Commitment Parameters		
Timeout waiting for a Storage Commitment Notification (maximum duration of applicability for a Storage Commitment Transaction UID).	Yes	1 hours
Maximum number of simultaneously accepted Associations by the Storage AE	No	1
Delay association release after sending a Storage Commitment Request (wait for a Storage Commitment Notification over the same association).	Yes	5 s
Storage Parameters		
Storage SCU time-out waiting for a response to a C-STORE-RQ	Yes	30 s
Maximum number of simultaneously initiated Associations by the Storage AE	Yes	1
Supported Transfer Syntaxes (separately configurable for each Presentation Context)	Yes	Implicit VR Little Endian Explicit VR Little Endian RLE Lossless ¹ JPEG Baseline ¹
Print Parameters		
Print SCU time-out waiting for a response to a N-CREATE-RQ	Yes	30 s
Print SCU time-out waiting for a response to a N-SET-RQ	Yes	30 s
Print SCU time-out waiting for a response to a N-ACTION-RQ	Yes	30 s
Supported Transfer Syntaxes	Yes	Implicit VR Little Endian Explicit VR Little Endian
Query/Retrieve Parameters		
Time-out waiting for an acceptance or rejection response to an Association Request	Yes	4s

Parameter	Configurable (Yes/No)	Default Value
Query/Retrieve SCU time-out waiting for a response to a C-FIND-RQ	Yes	150s
Query/Retrieve SCU time-out waiting for a response to a C-MOVE-RQ	Yes	150s
Query/Retrieve SCU time-out waiting for a response to a C-FIND-CANCEL-RQ	Yes	1s
Query/Retrieve SCU time-out waiting for a response to a C-MOVE-CANCEL-RQ	Yes	1s
Time-out waiting for a response to an Association release request	Yes	30s

Note: 1. RLE or JPEG is not used for Structured Report storage.

5. MEDIA INTERCHANGE

5.1 IMPLEMENTATION MODEL

5.1.1 Application Data Flow




Figure 5.1-1 APPLICATION DATA FLOW DIAGRAM FOR MEDIA STORAGE

- The Offline-Media Application Entity exports Ultrasound and Ultrasound Multi-frame images, and Structured Reports to some kinds of media. It is associated with the local real-world activity "STORE" button or "Copy" button.

5.1.2 Functional Definition of AEs

5.1.2.1 Functional Definition of Offline-Media Application Entity

The Media Application Entity can be requested in two modes. After the patient identification is supplied by the operator, pressing the "STORE" button will directly pass an image to the Offline-Media Application Entity. Or when it is configured to store the image in the local drive, the image is buffered in the drive for later reference. By pressing the "Copy" button, the images selected by the operator will be sent to the Offline-Media Application Entity. A storage association will be initiated by clicking the "Copy" button on the screen or pressing the "Copy" button on the LCD touch panel.

"CD Burning" button will export the images to the CD-R media from CD-R Buffer.

5.1.3 Sequencing of Real-World Activities

At least one object instance must exist and be selected before the Offline-Media Application Entity can be invoked. The operator can insert DVD-RAM or a new CD-R media at any time before or after invocation of the Offline-Media Application Entity. The Offline-Media Application Entity will wait indefinitely for a media to be inserted before starting to write to the CD-R device. If no CD-R media is available, the export job can be canceled from the job queue.

"STORE" button will export the acquired image. "Copy" button will export the images selected by the operator. Storage device can be selected at menu setting(DVD/USB). "CD Burning" button will export the images to the CD-R media from CD-R Buffer.

5.1.4 File Meta Information Options

The implementation information written to the File Meta Header in each file is:

Table 5.1-1 DICOM IMPLEMENTATION CLASS AND VERSION FOR MEDIA STORAGE

Implementation Class UID	1.2.392.200039.117
Implementation Version Name	ALOKA20161222

	(subject to change without notice)
--	------------------------------------

5.2 AE SPECIFICATIONS

5.2.1 Offline-Media Application Entity Specification

The Offline-Media Application Entity provides standard conformance to the DICOM Interchange Option of the Media Storage Service Class. The Application Profiles and roles are listed below:

Table 5.2-1 APPLICATION PROFILES, ACTIVITIES AND ROLES FOR OFFLINE-MEDIA

Application Profiles Supported	Single Frame	Multi-Frame	Real World Activity	Role
Image Display	STD-US-ID-SF-xxxx	STD-US-ID-MF-xxxx	"Store (USB)"	FSC
Spatial Calibration	STD-US-SC-SF-xxxx	STD-US-SC-MF-xxxx	"Store(CD-R buffer)", "STORE(DVD)", "Copy(USB)", "Copy(CD-R buffer)", "Copy(DVD)", "Burning CD"	

5.2.1.1 File Meta Information for the Application Entity

The Source Application Entity Title included in the File Meta Header is same that of local Storage AET.

5.2.1.2 Real-World Activities

5.2.1.2.1 Activity – Export to Media

The Offline-Media Application Entity acts as an FSC using the interchange option when requested to export SOP Instances from the local database to a removable storage medium, DVD-RAMs, CD-Rs, or USBs.

Object Instances can be added to a storage medium and if it reaches full, the user will be prompted to replace to an empty medium. The DicomDIR will be updated each time an object instance is successfully written to it.

When the object instance is exported to CD-R, the user will be prompted to insert an empty CD-R for each export. The contents of the export job will be written together with a corresponding DicomDIR to a single-session CD-R. Writing in multi-session mode is not supported.

5.2.1.2.1.1 Media Storage Application Profiles

The supported Offline-Media Application Profiles are listed in Table 5.2-1.

5.2.1.2.1.1.1 Options

The Offline-Media Application Entity supports the SOP Classes and Transfer Syntaxes listed in the Table below:

Table 5.2-2 IODS, SOP CLASSES AND TRANSFER SYNTAXES FOR OFFLINE MEDIA

Information Object Definition	SOP Class UID	Transfer Syntax	Transfer Syntax UID
Media Storage Directory Storage	1.2.840.10008.1.3.10	Explicit VR Little Endian	1.2.840.10008.1.2.1
Ultrasound Image Storage	1.2.840.10008.5.1.4.1 .1.6.1	Explicit VR Little Endian Implicit VR Little Endian JPEG Baseline Compression RLE Lossless Compression	1.2.840.10008.1.2.1 1.2.840.10008.1.2.1 ¹ 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.5
Ultrasound Multi-frame Image Storage	1.2.840.10008.5.1.4.1 .1.3.1	JPEG Baseline Compression	1.2.840.10008.1.2.4.50
Ultrasound Multi-frame Image Storage for 3D	1.2.840.10008.5.1.4.1 .1.3.1	Explicit VR Little Endian Implicit VR Little Endian RLE Lossless Compression	1.2.840.10008.1.2.1 1.2.840.10008.1.2 1.2.840.10008.1.2.5
Secondary Capture Image Storage	1.2.840.10008.5.1.4.1 .1.7	Explicit VR Little Endian Implicit VR Little Endian JPEG Baseline Compression RLE Lossless Compression	1.2.840.10008.1.2.1 1.2.840.10008.1.2 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.5
Multi-frame True Color Secondary Capture Image Storage	1.2.840.10008.5.1.4.1 .1.7.4	JPEG Baseline Compression	1.2.840.10008.1.2.4.50
Comprehensive SR	1.2.840.10008.5.1.4.1 .1.88.33	Explicit VR Little Endian	1.2.840.10008.1.2.1
Encapsulated PDF Storage	1.2.840.10008.5.1.4.1 .1.104.1	Explicit VR Little Endian	1.2.840.10008.1.2.1

Note: 1. Cardiac 3D Volume does not support Implicit Little Endian for offline media.

5.3 AUGMENTED AND PRIVATE APPLICATION PROFILES

This product does not support any augmented or private Application Profiles.

5.4 MEDIA CONFIGURATION

The Application Entity Title for Media Services is same that is configured for Storage Service:

Table 5.4-1 AE TITLE CONFIGURATION TABLE

Application Entity	Default AE Title
Offline-Media	No default

6. SUPPORT OF CHARACTER SETS

DICOM applications of PDF reports and the Structured Reports with Advanced Cardiac feature of this product support the character set below:

ISO_IR 192 (Unicode in UTF-8)

DICOM application of Structured Report with Non-Advanced Cardiac feature of this product supports the character set below:

- ISO_IR 144 (ISO 8859-5:Cyrillic) for Russian language configuration.
- ISO_IR 100 (ISO 8859-1:1987 Latin Alphabet No.1 supplementary set) for other language configuration.

All other DICOM applications of this product support the

ISO_IR 100 (ISO 8859-1:1987 Latin Alphabet No. 1 supplementary set)

7. SECURITY

This product supports security measures described in section 7.1.

In addition to the section, this product supports below.

- a. Incoming TCP packets to ports other than configured by Storage AE and Query/Retrieve AE listed in Table 4.4-1 are blocked.

It is recommended that any communication with external hosts and services outside the locally secured environment use appropriate secure network channels (e.g. such as a Virtual Private Network (VPN)).

Other network security procedures such as automated intrusion detection may be appropriate in some environments. Additional security features may be established by the local security policy and are beyond the scope of this conformance statement.

7.1 SECURITY PROFILES ¹

This product conforms to the Non-Downgrading BCP 195 TLS Secure Transport Connection Profile and AES TLS Secure Transport Connection Profile. The protocol has features listed in the Table below:

Table 7.1-1 SECURITY PROFILES FEATURES

TLS version	Cipher Suites
TLS 1.2 ²	TLS_DHE_RSA_WITH_AES_128_GCM_SHA256 TLS_ECDHE_RSA_WITH_AES_128_GCM_SHA256 TLS_DHE_RSA_WITH_AES_256_GCM_SHA384 TLS_ECDHE_RSA_WITH_AES_256_GCM_SHA384 TLS_ECDHE_ECDSA_WITH_AES_256_GCM_SHA384 TLS_ECDHE_ECDSA_WITH_AES_128_GCM_SHA256
TLS 1.0 ³	TLS_RSA_WITH_AES_128_CBC_SHA TLS_RSA_WITH_3DES_EDE_CBC_SHA

When an integrity check fails, the connection will be dropped per the TLS protocol, causing both the sender and the receiver to issue an A-P-ABORT indication to the upper layers with an implementation-specific provider reason. This product does not define the reason.

Note: 1. Only products for USA.

Note: 2. Not supported Storage Commitment Application Entity, Query/Retrieve Application Entity nor Hardcopy Application Entity.

Note: 3. Not supported Hardcopy Application Entity.

8. ANNEXES

8.1 IOD CONTENTS

8.1.1 Created SOP Instances

Table 8.1-1 specifies the attributes of an Ultrasound and Ultrasound Multi-frame Image transmitted by the storage application of this product.

Table 8.1-2 specifies the attributes of an Secondary Capture and Multi-frame True Color Secondary Capture Image transmitted by the storage application of this product.

Table 8.1-3 specifies the attributes of a Structured Report transmitted by the storage application of this product.

Table 8.1-4 specifies the attributes of an Encapsulated PDF transmitted by the storage application of this product.

The following tables use a number of abbreviations. The abbreviations used in the “Presence of ...” column are:

ALWAYS	Always Present
ANAP	Attribute Not Always Present
VNAP	Value Not Always Present (attribute sent zero length if no value is present)
EMPTY	Attribute is sent without a value

The abbreviations used in the “Source” column:

MWL	the attribute value source Modality Worklist
USER	the attribute value source is from User input
AUTO	the attribute value is generated automatically
MPPS	the attribute value is the same as that use for Modality Performed Procedure Step
CONFIG	the attribute value source is a configurable parameter

Attributes in *Italic* are additions to the Standard Information Entity Modules.

The values of private attributes may be removed if the Teaching File option is set.

8.1.1.1 Ultrasound and Ultrasound Multi-frame Image IODs

Table 8.1-1 IOD OF ULTRASOUND AND ULTRASOUND MULTI-FRAME IMAGE SOP INSTANCES

IE	Module	Reference	Presence of Module	
			US	US-mf
Patient	Patient	Table 8.1-5	ALWAYS	ALWAYS
	Clinical Trial Subject	-	Not used	Not used

Study	General Study	Table 8.1-6	ALWAYS	ALWAYS
	Patient Study	Table 8.1-7	ALWAYS	ALWAYS
	Clinical Trial Study	-	Not used	Not used
Series	General Series	Table 8.1-8	ALWAYS	ALWAYS
	Clinical Trial Series	-	Not used	Not used
Frame of Reference	Frame of Reference	-	Not used	Not used
	Synchronization	-	Not used	Not used
Equipment	General Equipment	Table 8.1-9	ALWAYS	ALWAYS
Image	General Image	Table 8.1-10	ALWAYS	ALWAYS
	Contrast/bolus	Table 8.1-11	ANAP	ANAP
	Cine	Table 8.1-12	Not used	ALWAYS
	Frame Pointers	Table 8.1-13	Not used	ALWAYS
	Multi-Frame	Table 8.1-14	Not used	ALWAYS
	US Region Calibration	Table 8.1-19	ANAP	ANAP
	US Image	Table 8.1-20	ALWAYS	ALWAYS
	Palette Color Lookup Table	-	Not used	Not used
	Image Pixel	Table 8.1-15	ALWAYS	ALWAYS
	Overlay Plane	-	Not used	Not used
	VOI LUT	Table 8.1-16	Only if (0028,0004) equals to "MONOCHROME2"	Only if (0028,0004) equals to "MONOCHROME2"
	SOP Common	Table 8.1-17	ALWAYS	ALWAYS

8.1.1.2 Secondary Capture and Multi-frame True Color Secondary Capture Image IODs

Table 8.1-2 IOD OF SECONDARY CAPTURE AND MULTI-FRAME TRUE COLOR SECONDARY CAPTURE IMAGE SOP INSTANCES

IE	Module	Reference	Presence of Module	
			SC	SC-mf
Patient	Patient	Table 8.1-5	ALWAYS	ALWAYS
	Clinical Trial Subject	-	Not used	Not used
Study	General Study	Table 8.1-6	ALWAYS	ALWAYS
	Patient Study	Table 8.1-7	ALWAYS	ALWAYS
	Clinical Trial Study	-	Not used	Not used
Series	General Series	Table 8.1-8	ALWAYS	ALWAYS
	Clinical Trial Series	-	Not used	Not used
Frame of Reference	Frame of Reference	-	Not used	Not used
	Synchronization	-	Not used	Not used
Equipment	General Equipment	Table 8.1-9	ALWAYS	ALWAYS
	SC Equipment	Table 8.1-26	ALWAYS	ALWAYS

Image	General Image	Table 8.1-10	ALWAYS	ALWAYS
	General Reference	-	Not used	Not used
	Image Pixel	Table 8.1-15	ALWAYS	ALWAYS
	Cine	Table 8.1-12	Not used	ALWAYS
	Multi-frame	Table 8.1-14	Not used	ALWAYS
	Frame Pointers	Table 8.1-13	Not used	ALWAYS
	Device	-	Not used	Not used
	Multi-frame Functional Groups	-	Not used	Not used
	Multi-frame Dimension	-	Not used	Not used
	Specimen	-	Not used	Not used
	SC Image	Table 8.1-27	ALWAYS	ALWAYS
	SC Multi-frame Image	Table 8.1-28	Not used	ALWAYS
	SC Multi-frame Vector	-	Not used	Not used
	Overlay Plane	-	Not used	Not used
	Modality LUT	-	Not used	Not used
	VOI LUT	-	Not used	Not used
	ICC Profile	-	Not used	Not used
	SOP Common	Table 8.1-17	ALWAYS	ALWAYS
	Common Instance Reference	-	Not used	Not used
	Frame Extraction	-	Not used	Not used

8.1.1.3 Comprehensive SR IOD

Table 8.1-3 IOD OF CREATED COMPREHENSIVE SR SOP INSTANCES

IE	Module	Reference	Presence of Module
Patient	Patient	Table 8.1-5	ALWAYS
	Specimen Identification	-	Not used
	Clinical Trial Subject	-	Not used
Study	General Study	Table 8.1-6	ALWAYS
	Patient Study	Table 8.1-7	ALWAYS ¹
	Clinical Trial Study	-	Not used
Series	SR Document Series	Table 8.1-21	ALWAYS
	Clinical Trial Series	-	Not used
Frame of Reference	Frame of Reference	-	Not used
	Synchronization	-	Not used
Equipment	General Equipment	Table 8.1-9	ALWAYS
Document	SR Document General	Table 8.1-22	ALWAYS
	SR Document Content	Table 8.1-23	ALWAYS
	SOP Common	Table 8.1-17	ALWAYS

Note: 1. ANAP for Advanced Cardiac feature.

8.1.1.4 Encapsulated PDF IOD

Table 8.1-4 IOD OF CREATED ENCAPSULATED PDF SOP INSTANCES

IE	Module	Reference	Presence of Module
Patient	Patient	Table 8.1-5	ALWAYS
	Clinical Trial Subject	-	Not used
Study	General Study	Table 8.1-6	ALWAYS
	Patient Study	Table 8.1-7	ANAP
	Clinical Trial Study	-	Not used
Series	Encapsulated Document Series	Table 8.1-24	ALWAYS
	Clinical Trial Series	-	Not used
Equipment	General Equipment	Table 8.1-9	ALWAYS
	SC Equipment	Table 8.1-18	ALWAYS
Encapsulated Document	Encapsulated Document	Table 8.1-25	ALWAYS
	SOP Common	Table 8.1-17	ALWAYS

8.1.1.5 Common Modules

Table 8.1-5 PATIENT MODULE OF CREATED SOP INSTANCES

Tag	VR	Attribute Name	Value	Presence of Value	Source
(0010,0010)	PN	Patient's Name	From Modality Worklist or input by user. Values supplied via Modality Worklist will be entered as received. Values supplied via user input will contain all 5 components (some possibly empty). Maximum 64 characters. ¹	VNAP	MWL/USER
(0010,0020)	LO	Patient ID	From Modality Worklist or input by user. Maximum 64 characters. ²	ALWAYS	MWL/USER
(0010,0030)	DA	Patient's Birth Date	From Modality Worklist or input by user ¹	VNAP	MWL/USER
(0010,0040)	CS	Patient's Sex	From Modality Worklist or input by user ¹	VNAP	MWL/USER
(0012,0062)	CS	Patient Identity Removed	YES = Teaching File option NO = Otherwise	ANAP ³	AUTO
(0012,0063)	LO	De-identification Method	Set "Dummy\Removed\Zero Length" if the Teaching File option is set.	ANAP ³	AUTO

Note: 1. The value may be replaced with a zero length if the Teaching File option is set.

Note: 2. The value may be replaced with a non-zero length value that may be a dummy value and consistent with the VR if the Teaching File option is set.

Note: 3. Attribute present only in image.

Table 8.1-6 GENERAL STUDY MODULE OF CREATED SOP INSTANCES

Tag	VR	Attribute Name	Value	Presence of Value	Source
(0008,0020)	DA	Study Date	"yyyymmdd"	ALWAYS	AUTO
(0008,0030)	TM	Study Time	"hhmmss"	ALWAYS	AUTO
(0008,0050)	SH	Accession Number	From Modality Worklist or input by user ²	VNAP	MWL/USER
(0008,0090)	PN	Referring Physician's Name	From Modality Worklist or Input by user ²	VNAP	MWL/USER
(0008,1030)	LO	Study Description	From Modality Worklist or input by user at Comment text box in study list. Maximum 64 bytes. Configurable in "ID Option" to be copied from Requested Procedure Description (0032,1060) or from Scheduled Procedure Step Description (0040,0007). ³	ANAP	MWL/USER
(0008,1032)	SQ	Procedure Code Sequence	From Modality Worklist, mapped from Requested Procedure Code Sequence (0032,1064) ³	ANAP ¹	MWL
>Include "Code Sequence Macro"					
(0008,1110)	SQ	Referenced Study Sequence	From Modality Worklist	ANAP ⁶	MWL
>(0008,1150)	UI	Referenced SOP Class UID	From Modality Worklist	VNAP	MWL
>(0008,1155)	UI	Referenced SOP Instance UID	From Modality Worklist	VNAP	MWL
(0008,1060)	PN	Name of Physician(s) reading Study.	Entered as "Reporting Phys" in the New Patient Registration ³	ANAP ⁵	USER
(0020,000D)	UI	Study Instance UID	From Modality Worklist or generated by device	ALWAYS	MWL/AUTO
(0020,0010)	SH	Study ID	Copied from Requested Procedure ID (0040,1001) in Worklist or generated by device. User may modify the value. ⁴	ALWAYS	MWL/AUTO/USER
(0032,1060)	LO	<i>Requested Procedure Description</i>	<i>From Modality Worklist</i> ²	<i>VNAP</i> ¹	<i>MWL</i>

Note: 1. Attribute Not Present in Structured Report, PDF, SC Image, and SC-mf image.

Note: 2. The value may be replaced with a zero length if the Teaching File option is set for image.

Note: 3. The attribute may be removed if the Teaching File option is set for image.

Note: 4. The value may be replaced with a non-zero length value that may be a dummy value and consistent with the VR if the Teaching File option is set for image.

Note: 5. Attribute Not Present in PDF and analysis condition file with SC image.

Note: 6. Attribute Not Present in Structured Report, PDF, and analysis condition file with SC image.

Table 8.1-7 PATIENT STUDY MODULE OF CREATED SOP INSTANCES

Tag	VR	Attribute Name	Value	Presence of Value	Source
(0010,1010)	AS	Patient's Age	Calculated from DoB input on base of actual Date ³	ANAP ⁵	AUTO
(0010,1020)	DS	Patient's Size	From Modality Worklist or user input ³	ANAP ⁵	MWL/ USER
(0010,1030)	DS	Patient's Weight	From Modality Worklist or user input ³	ANAP ⁵	MWL/ USER
(0010,2000)	LO	Medical Alerts	From Modality Worklist ²	VNAP ^{1, 5}	MWL
(0010,2110)	LO	Contrast Allergies	From Modality Worklist ²	VNAP ^{1, 5}	MWL
(0010,2180)	SH	Occupation	From Modality Worklist or User input ³	ANAP	MWL/ USER
(0010,21D0)	DA	Last Menstrual Date	From Modality Worklist or User input ³	ANAP ^{4, 5}	MWL/ USER
(0038,0010)	LO	Admission ID	Empty	VNAP ⁶	AUTO
(0038,0050)	LO	Special Needs	From Modality Worklist ²	VNAP ^{1, 5}	MWL

Note: 1. Attribute Not Present in Structured Report and PDF.

Note: 2. The value may be replaced with a zero length if the Teaching File option is set for image.

Note: 3. The attribute may be removed if the Teaching File option is set for image.

Note: 4. Attribute Not Present in Structured Report, PDF, SC Image, and SC-mf image.

Note: 5. Attribute Not Present in analysis condition file with SC image.

Note: 6. Attribute presents only in analysis condition file with SC image.

Table 8.1-8 GENERAL SERIES MODULE OF CREATED SOP INSTANCES

Tag	VR	Attribute Name	Value	Presence of Value	Source
(0008,0021)	DA	Series Date	"yyyymmdd"	ALWAYS	AUTO
(0008,0031)	TM	Series Time	"hhmmss"	ALWAYS	AUTO
(0008,0060)	CS	Modality	"US"	ALWAYS	AUTO
(0008,103E)	LO	Series Description	From Modality Worklist or input by user: Maximum 64 bytes. Configurable in "ID Option" to be copied from Scheduled Procedure Step Description (0040,0007) or from Requested Procedure Description (0032,1060). ²	ANAP ⁵	MWL/ USER

(0008,1070)	PN	Operator's Name	Sonographer field in Study list. Maximum 64 characters.	ANAP ⁵	USER
(0008,1111)	SQ	Referenced Performed Procedure Step Sequence	Identifies the MPPS SOP Instance to which this image is related ²	ANAP ⁵	MPPS
>(0008,1150)	UI	Referenced SOP Class UID	MPPS SOP Class UID	ALWAYS	MPPS
>(0008,1155)	UI	Referenced SOP Instance UID	MPPS SOP Instance UID	ALWAYS	MPPS
(0018,0015)	CS	Body Part Examined	Set by user from a pick list	ANAP ⁴	USER
(0018,1030)	LO	Protocol Name	Application type selected by user or Stress Echo Protocol Name.	ANAP ⁴	AUTO
(0020,000E)	UI	Series Instance UID	Generated by device	ALWAYS	AUTO
(0020,0011)	IS	Series Number	Generated by device ³	ALWAYS	AUTO
(0020,0060)	CS	Laterality	Set by user from a pick list	ANAP ⁴	USER
(0032,1032)	PN	<i>Requesting Physician</i>	<i>From Modality Worklist ¹</i>	VNAP ⁴	MWL
(0040,0244)	DA	Performed Procedure Step Start Date	Same as Series Date (0008,0021) ²	ALWAYS	MPPS
(0040,0245)	TM	Performed Procedure Step Start Time	Same as Series Time (0008,0031) ²	ALWAYS	MPPS
(0040,0253)	SH	Performed Procedure Step ID	Sequence Number ²	ALWAYS	AUTO
(0040,0254)	LO	Performed Procedure Step Description	Input by user Same as MPPS. From user input as Study Description. Maximum 64 characters. ²	ANAP	USER
(0040,0260)	SQ	Performed Protocol Code Sequence	Derived from Scheduled Protocol Code Sequence or overridden by Ultrasound Protocol Types (CID 12001) when Stress Echo Mode is requested. May be modified by user. ²	ANAP	AUTO/ USER
>Include "Code Sequence Macro"					
(0040,0275)	SQ	Request Attributes Sequence	Zero or 1 item will be present ²	ANAP	AUTO
>(0032,1060)	LO	Requested Procedure Description	From Modality Worklist	ANAP	MWL
>(0040,0007)	LO	Scheduled Procedure Step Description	From Modality Worklist	ANAP	MWL
>(0040,0008)	SQ	Scheduled Protocol Code Sequence	From Modality Worklist	ANAP	MWL
>>Include "Code Sequence Macro"					
>(0040,0009)	SH	Scheduled Procedure Step ID	From Modality Worklist	ANAP	MWL

>(0040,1001)	SH	Requested Procedure ID	From Modality Worklist or Input by user.	ANAP	MWL/USER
--------------	----	------------------------	--	------	----------

- Note: 1. The value may be replaced with a zero length if the Teaching File option is set for image.
- Note: 2. The attribute may be removed if the Teaching File option is set for image.
- Note: 3. The value may be replaced with a non-zero length value that may be a dummy value and consistent with the VR if the Teaching File option is set for image.
- Note: 4. Attribute Not Present in SC and SC-mf image.
- Note: 5. Attribute Not Present in analysis condition file with SC image.

Table 8.1-9 GENERAL EQUIPMENT MODULE OF CREATED SOP INSTANCES

Tag	VR	Attribute Name	Value	Presence of Value	Source
(0008,0070)	LO	Manufacturer	"FUJIFILM Corporation"	ALWAYS	AUTO
(0008,0080)	LO	Institution Name	From Configuration ¹	VNAP	CONFIG
(0008,0081)	ST	Institution Address	From Configuration ^{1, 2}	ANAP	CONFIG
(0008,1010)	SH	Station Name	From Configuration ¹	ANAP	CONFIG
(0008,1040)	LO	Institutional Department Name	From Configuration ^{1, 2}	ANAP	CONFIG
(0008,1090)	LO	Manufacturer's Model Name	"LISENDO 880" ¹	ALWAYS	AUTO
(0018,1000)	LO	Device Serial Number	Built-in ³	ANAP	AUTO
(0018,1020)	LO	Software Version(s)	Built-in	ALWAYS	AUTO

- Note: 1. The value may be replaced with a zero length if the Teaching File option is set for image.
- Note: 2. Attribute Not Present in Structured Report, PDF, SC image, and SC-mf image.
- Note: 3. Attribute Not Present in SC and SC-mf image.

Table 8.1-10 GENERAL IMAGE MODULE OF CREATED SOP INSTANCES

Tag	VR	Attribute Name	Value	Presence of Value	Source
(0008,0008)	CS	Image Type	ORIGINAL or DERIVED for value 1, PRIMARY or SECONDARY for value 2. See Table 8.1-19 for value 3 and 4.	ANAP ¹	AUTO
(0008,0023)	DA	Content Date	"yyyymmdd"	ALWAYS	AUTO
(0008,0033)	TM	Content Time	"hhmmss"	ALWAYS	AUTO
(0020,0013)	IS	Instance Number	Generated by device ²	ALWAYS	AUTO
(0020,0020)	CS	Patient Orientation	From Pull Down Menu or Input by user.	ANAP	USER
(0020,4000)	LT	Image Comments	Input by user	ANAP ³	USER
(0028,2110)	CS	Lossy Image Compression	Generated by device	ANAP	AUTO

(0028,0301)	CS	Burned In Annotation	YES = Image may contain sufficient burned in annotation to identify the patient. NO = Image does not contain sufficient burned in annotation to identify the patient.	ANAP	AUTO
-------------	----	----------------------	--	------	------

- Note:
1. Attribute Not Present in analysis condition file with SC image.
 2. The value may be replaced with a non-zero length value that may be a dummy value and consistent with the VR if the Teaching File option is set for image.
 3. Attribute presents only in analysis condition file with SC image.

Table 8.1-11 CONTRAST/BOLUS MODULE OF CREATED SOP INSTANCES

Tag	VR	Attribute Name	Value	Presence of Value	Source
(0018,0010)	LO	Contrast/Bolus Agent	From MWL or Input by user ¹	ANAP	MWL/ USER
(0018,1040)	LO	Contrast/Bolus Route	May be input by user	ANAP	USER
(0018,1041)	DS	Contrast/Bolus Volume	May be input by user	ANAP	USER
(0018,1042)	TM	Contrast/Bolus Start Time	Generated by device	ANAP	AUTO
(0018,1044)	DS	Contrast/Bolus Total Dose	May be input by user	ANAP	USER

- Note:
1. The value may be replaced with a zero length if the Teaching File option is set for image.

Table 8.1-12 CINE MODULE OF CREATED SOP INSTANCES

Tag	VR	Attribute Name	Value	Presence of Value	Source
(0008,2142)	IS	Start Trim	1	ANAP ¹	AUTO
(0008,2143)	IS	Stop Trim	Last frame in Multi-frame IOD	ANAP ¹	AUTO
(0008,2144)	IS	Recommended Display Frame Rate	frames per second [fps]	ANAP ¹	AUTO
(0009,0010)	LO	Private Identification Code	ALOKA:1.2.392.200039.117.2	ANAP	AUTO
(0009,1004)	US		Private attribute without PHI	ANAP	AUTO
(0009,1006)	US		Private attribute without PHI	ANAP	AUTO
(0018,0040)	IS	Cine Rate	frames per second [fps]	ANAP ¹	AUTO
(0018,0072)	DS	Effective Duration	0	ANAP ¹	AUTO
(0018,1063)	DS	Frame Time	If Frame Increment Pointer (0028,0009) is Frame Time	ANAP	AUTO
(0018,1065)	DS	Frame Time Vector	If Frame Increment Pointer (0028,0009) is Frame Time Vector	ANAP	AUTO
(0018,1066)	DS	Frame Delay	0	ANAP ¹	AUTO

(0018,1242)	IS	Actual Frame Duration		ANAP ¹	AUTO
(0018,1244)	US	Preferred Playback Sequence	0=Looping	ANAP ¹	AUTO
(0028,6040)	US	R Wave Pointer		ANAP	AUTO

Note: 1. Attribute Not Present in SC-mf Image.

Table 8.1-13 FRAME POINTERS MODULE OF CREATED SOP INSTANCES

Tag	VR	Attribute Name	Value	Presence of Value	Source
(0028,6010)	US	Representative Frame Number	May be used in multi-frame image	ANAP	AUTO
(0028,6020)	US	Frame Number of Interest (FOI)		ANAP	AUTO
(0028,6022)	LO	Frame(s) of Interest Description		ANAP	AUTO
(0028,6023)	CS	Frame of Interest Type		ANAP	AUTO

Table 8.1-14 MULTI-FRAME MODULE OF CREATED SOP INSTANCES

Tag	VR	Attribute Name	Value	Presence of Value	Source
(0028,0008)	IS	Number of Frames		ALWAYS	AUTO
(0028,0009)	AT	Frame Increment Pointer	(0018,1063) or (0018,1065)	ALWAYS	AUTO

Table 8.1-15 IMAGE PIXEL MODULE OF CREATED SOP INSTANCES

Tag	VR	Attribute Name	Value	Presence of Value	Source
(0028,0002)	US	Samples per Pixel	1=PALETTE COLOR, MONOCHROME2, 3=RGB, YBR_FULL_422	ALWAYS	AUTO
(0028,0004)	CS	Photometric Interpretation	"PALETTE COLOR", "RGB", "YBR_FULL_422", "MONOCHROME2" (MONOCHROME2 is not applicable to media storage except for images acquired for 3D Reconstruction.)	ALWAYS	AUTO/ CONFIG
(0028,0006)	US	Planar Configuration	0=color-by-pixel for YBR_FULL_422, RGB 1=color-by-plane for RGB (0=color-by-pixel for RGB is not applicable to compression transfer syntax or media storage.)	ANAP	AUTO/ CONFIG

(0028,0010)	US	Rows	832, 768, 600, or 480 for US and US-mf image. (may vary if image is acquired for 3D Reconstruction and image is stored on Analysis screen page) 971 for SC and SC-mf image.	ALWAYS	AUTO
(0028,0011)	US	Columns	1152, 1024, 800 or 640 for US and US-mf image. (may vary if image is acquired for 3D Reconstruction and image is stored on Analysis screen page) 1648 for SC image. 1650 for SC-mf image.	ALWAYS	AUTO
(0028,0034)	IS	Pixel Aspect Ratio	Only if pixel aspect ratio is not 1:1.	ANAP	AUTO
(0028,0100)	US	Bits Allocated	8 for RGB, YBR_FULL_422, MONOCHROME2, 16 for PALETTE COLOR	ALWAYS	AUTO
(0028,0101)	US	Bits Stored	8 for RGB, YBR_FULL_422, MONOCHROME2, 16 for PALETTE COLOR	ALWAYS	AUTO
(0028,0102)	US	High Bit	7 for RGB, YBR_FULL_422, MONOCHROME2, 15 for PALETTE COLOR	ALWAYS	AUTO
(0028,0103)	US	Pixel Representation	0	ALWAYS	AUTO
(7FE0,0010)	OW /OB	Pixel Data	The Pixel Data contains burned-in annotation (Patient ID, Patient's Name, Scale Mark etc.)	ALWAYS	AUTO

Table 8.1-16 VOI LUT MODULE OF CREATED SOP INSTANCES

Tag	VR	Attribute Name	Value	Presence of Value	Source
(0028,3010)	SQ	VOI LUT Sequence	One item	ANAP	AUTO
>(0028,3002)	US	LUT Descriptor	<256,0,8>	ALWAYS	AUTO
>(0028,3006)	US	LUT Data	LUT	ALWAYS	AUTO
(0028,1050)	DS	Window Center	0-255	ANAP	USER

(0028,1051)	DS	Window Width	1-256	ANAP	USER
-------------	----	--------------	-------	------	------

Table 8.1-17 SOP COMMON MODULE OF CREATED SOP INSTANCES

Tag	VR	Attribute Name	Value	Presence of Value	Source
(0008,0005)	CS	Specific Character Set	See chapter 6	ANAP	MWL/AUTO
(0008,0012)	DA	Instance Creation Date	Instance Date	ANAP ¹	AUTO
(0008,0013)	TM	Instance Creation Time	Instance Time	ANAP ¹	AUTO
(0008,0014)	UI	Instance Creator UID	Built-in	ANAP ¹	AUTO
(0008,0016)	UI	SOP Class UID	1.2.840.10008.5.1.4.1.1.6.1 1.2.840.10008.5.1.4.1.1.3.1 1.2.840.10008.5.1.4.1.1.88.3 1.2.840.10008.5.1.4.1.1.104.1 1.2.840.10008.5.1.4.1.1.7 1.2.840.10008.5.1.4.1.1.7.4	ALWAYS	AUTO
(0008,0018)	UI	SOP Instance UID	Generated by device	ALWAYS	AUTO

Note: 1. Attribute Not Present in Structured Report with Non-Advanced Cardiac feature, SC image, and SC-mf image.

Table 8.1-18 SC EQUIPMENT OF CREATED SOP INSTANCES

Tag	VR	Attribute Name	Value	Presence of Value	Source
(0008,0064)	CS	Conversion Type	WSD	ALWAYS	AUTO

8.1.1.6 Ultrasound Modules

Table 8.1-19 US REGION CALIBRATION MODULE OF CREATED SOP INSTANCES

Tag	VR	Attribute Name	Value	Presence of Value	Source
(0018,6011)	SQ	Sequence of Ultrasound Regions	If Region Calibration is turned on	ANAP	CONFIG
>(0018,6012)	US	Region Spatial Format	0=none, 1=cross section, 2=M-mode, 3=spectral	ALWAYS	AUTO
>(0018,6014)	US	Region Data Type	0=none, 1=tissue, 2=color flow velocity, 3=PW, 4=CW, D=Gray Scale, E=Color Scale	ALWAYS	AUTO
>(0018,6016)	UL	Region Flag	1=transparent, 2=scale protected, 4=Doppler represented in frequency.	ALWAYS	AUTO
>(0018,6018)	UL	Region Location Min. x0	X and Y coordinates of upper	ALWAYS	AUTO

>(0018,601A)	UL	Region Location Min. y0	left corner of the region	ALWAYS	AUTO
>(0018,601C)	UL	Region Location Max. x1	X and Y coordinates of lower right corner of the region	ALWAYS	AUTO
>(0018,601E)	UL	Region Location Max. y1		ALWAYS	AUTO
>(0018,6020)	SL	Referenced Pixel X0	B: Transducer surface center, M: Transducer surface left, and D: Baseline left	ANAP	AUTO
>(0018,6022)	SL	Referenced Pixel Y0		ANAP	AUTO
>(0018,6024)	US	Physical Unit X Direction	Physical units for X and Y directions	ALWAYS	AUTO
>(0018,6026)	US	Physical Unit Y Direction		ALWAYS	AUTO
>(0018,6028)	FD	Referenced Pixel Physical Value X	Reference pixel physical values in X and Y direction physical units, respectively	ANAP	AUTO
>(0018,602A)	FD	Referenced Pixel Physical Value Y		ANAP	AUTO
>(0018,602C)	FD	Physical Delta X	Physical values in X and Y direction physical units, respectively	ALWAYS	AUTO
>(0018,602E)	FD	Physical Delta Y		ALWAYS	AUTO
>(0018,6030)	UL	Transducer Frequency	In [kHz]	ANAP	AUTO
>(0018,6032)	UL	Pulse Repetition Frequency	PRF in [kHz]	ANAP	AUTO
>(0018,6034)	FD	Doppler Correction Angle	Used if Doppler image [deg]	ANAP	AUTO
>(0018,6036)	FD	Steering Angle	Used if Color Flow image [deg]	ANAP	AUTO
>(0019,0010)	LO	Private Identification Code	ALOKA:1.2.392.200039.117.2	ANAP	AUTO
>(0019,1008)	FD		Private attribute without PHI	ANAP	AUTO
>(0019,100C)	CS		Private attribute without PHI	ANAP	AUTO
>(0019,100E)	DS		Private attribute without PHI	ANAP	AUTO
>(0019,1018)	SL		Private attribute without PHI	ANAP	AUTO
>(0019,101A)	SL		Private attribute without PHI	ANAP	AUTO
>(0019,1040)	SS		Private attribute without PHI	ANAP	AUTO
>(0019,1046)	US		Private attribute without PHI	ANAP	AUTO
>(0019,1050)	SL		Private attribute without PHI	ANAP	AUTO
>(0019,1052)	DS		Private attribute without PHI	ANAP	AUTO
>(0019,1054)	DS		Private attribute without PHI	ANAP	AUTO
>(0019,1056)	FD		Private attribute without PHI	ANAP	AUTO
>(0019,1060)	US		Private attribute without PHI	ANAP	AUTO
>(0019,1061)	UL		Private attribute without PHI	ANAP	AUTO

Table 8.1-20 US IMAGE MODULE OF CREATED SOP INSTANCES

Tag	VR	Attribute Name	Value	Presence of Value	Source
(0008,0008)	CS	Image Type	Value 3 is picked by the user from the Pull Down List, and Value 4 is set automatically.	VNAP	AUTO
(0008,2120)	SH	Stage Name	If Stress Echo protocol used	ANAP	AUTO
(0008,2122)	IS	Stage Number	If Stress Echo protocol used	ANAP	AUTO
(0008,2124)	IS	Number of Stages	If Stress Echo protocol used	ANAP	AUTO
(0008,2127)	SH	View Name	If Stress Echo protocol used	ANAP	AUTO
(0008,2128)	IS	View Number	If Stress Echo protocol used	ANAP	AUTO
(0008,2129)	IS	Number of Event Timers	If multi-frame and ECG Sync	ANAP	AUTO
(0008,212A)	IS	Number of Views in Stage	If Stress Echo protocol used	ANAP	AUTO
(0008,2130)	DS	Event Elapsed Time(s)		ANAP	AUTO
(0008,2132)	LO	Event Timer Name(s)		ANAP	AUTO
(0009,0010)	LO	Private Identification Code	ALOKA:1.2.392.200039.117.2	ANAP	AUTO
(0009,1000)	SH		Private attribute without PHI	ANAP	AUTO
(0009,1004)	US		Private attribute without PHI	ANAP	AUTO
(0009,1006)	US		Private attribute without PHI	ANAP	AUTO
(0009,100A)	SH		Private attribute without PHI	ANAP	AUTO
(0009,1012)	US		Private attribute without PHI	ANAP	AUTO
(0009,1014)	FD		Private attribute without PHI	ANAP	AUTO
(0009,1016)	DA		Private attribute without PHI	ANAP	AUTO
(0009,1018)	TM		Private attribute without PHI	ANAP	AUTO
(0009,101A)	LO		Private attribute without PHI	ANAP	AUTO
(0009,1020)	CS		Private attribute without PHI	ANAP	AUTO
(0009,1022)	CS		Private attribute without PHI	ANAP	AUTO
(0009,1024)	CS		Private attribute without PHI	ANAP	AUTO
(0009,1026)	IS		Private attribute without PHI	ANAP	AUTO
(0009,1028)	IS		Private attribute without PHI	ANAP	AUTO
(0009,102A)	DS		Private attribute without PHI	ANAP	AUTO
(0009,102C)	FD		Private attribute without PHI	ANAP	AUTO
(0009,102E)	FD		Private attribute without PHI	ANAP	AUTO
(0009,1030)	FD		Private attribute without PHI	ANAP	AUTO

(0009,1032)	DS		Private attribute without PHI	ANAP	AUTO
(0009,1034)	CS		Private attribute without PHI	ANAP	AUTO
(0009,1036)	FD		Private attribute without PHI	ANAP	AUTO
(0009,1038)	CS		Private attribute without PHI	ANAP	AUTO
(0009,103A)	FD		Private attribute without PHI	ANAP	AUTO
(0009,103C)	FD		Private attribute without PHI	ANAP	AUTO
(0009,103E)	FD		Private attribute without PHI	ANAP	AUTO
(0009,1040)	FD		Private attribute without PHI	ANAP	AUTO
(0009,1042)	FD		Private attribute without PHI	ANAP	AUTO
(0009,1044)	FD		Private attribute without PHI	ANAP	AUTO
(0009,1046)	US		Private attribute without PHI	ANAP	AUTO
(0009,1048)	FD		Private attribute without PHI	ANAP	AUTO
(0009,104A)	FD		Private attribute without PHI	ANAP	AUTO
(0009,104C)	FD		Private attribute without PHI	ANAP	AUTO
(0009,104E)	FD		Private attribute without PHI	ANAP	AUTO
(0009,1050)	US		Private attribute without PHI	ANAP	AUTO
(0009,1052)	US		Private attribute without PHI	ANAP	AUTO
(0009,1054)	LO		Private attribute without PHI	ANAP	AUTO
(0009,1056)	IS		Private attribute without PHI	ANAP	AUTO
(0009,1058)	CS		Private attribute without PHI	ANAP	AUTO
(0009,105A)	UL		Private attribute without PHI	ANAP	AUTO
(0009,105C)	CS		Private attribute without PHI	ANAP	AUTO
(0009,105E)	DA		Private attribute without PHI	ANAP	AUTO
(0009,1060)	TM		Private attribute without PHI	ANAP	AUTO
(0009,1062)	CS		Private attribute without PHI	ANAP	AUTO
(0018,1060)	DS	Trigger Time	If ECG is connected	ANAP	AUTO
(0018,1062)	IS	Nominal Interval	If ECG is connected	ANAP	AUTO
(0018,1088)	IS	Heart Rate	If ECG signal is input	ANAP	AUTO
(0018,5000)	SH	Output Power		ANAP	AUTO
(0018,5010)	LO	Transducer Data		ALWAYS	AUTO
(0018,5022)	DS	Mechanical Index		ANAP	AUTO
(0018,5024)	DS	Bone Thermal Index		ANAP	AUTO
(0018,5026)	DS	Cranial Thermal Index		ANAP	AUTO
(0018,5027)	DS	Soft Tissue Thermal Index		ANAP	AUTO
(0018,5050)	IS	Depth of Scan Field		ANAP	AUTO
(0018,6031)	CS	Transducer Type		ALWAYS	AUTO
(0028,0002)	US	Samples per Pixel	3= RGB, YBR_FULL_422 1= PALETTE COLOR, MONOCHROME2	ALWAYS	AUTO
(0028,0004)	CS	Photometric Interpretation	PALETTE COLOR, RGB,	ALWAYS	AUTO/

			YBR_FULL_422, MONOCHROME2 (MONOCHROME2 is not applicable to media storage except for images acquired for 3D Reconstruction.)		CONFIG
(0028,0014)	US	Ultrasound Color Data Present	Only if color image acquired	ANAP	AUTO
(0028,0100)	US	Bits Allocated	8= RGB, YBR_FULL_422, MONOCHROME2, 16= PALETTE COLOR	ALWAYS	AUTO
(0028,0101)	US	Bits Stored	8= RGB, YBR_FULL_422, MONOCHROME2, 16= PALETTE COLOR	ALWAYS	AUTO
(0028,0102)	US	High Bit	7= RGB, YBR_FULL_422, MONOCHROME2, 15= PALETTE COLOR	ALWAYS	AUTO
(0028,0103)	US	Pixel Representation	0000H	ALWAYS	AUTO
(0028,2110)	CS	Lossy Image Compression	If lossy image compression applied	ANAP	AUTO
(0028,6040)	US	R Wave Pointer	Only if multi-frame	ANAP	AUTO
(0040,000A)	SQ	Stage Code Sequence	If Stress Echo protocol used	ANAP	AUTO
(0054,0220)	SQ	View Code Sequence	If Stress Echo protocol used	ANAP	AUTO
(830F,0010)	LO	<i>Private Identification Code</i>	<i>ALOKA:1.2.392.200039.117.2</i>	<i>ANAP</i>	<i>AUTO</i>
(830F,1000)	US		<i>Private attribute without PHI</i>	<i>ANAP</i>	<i>AUTO</i>
(830F,10FF)	OB		<i>Private attribute without PHI</i>	<i>ANAP</i>	<i>AUTO</i>
(833D,0010)	LO	<i>Private Identification Code</i>	<i>ALOKA:1.2.392.200039.117.2</i>	<i>ANAP</i>	<i>AUTO</i>
(833D,0020)	LO	<i>Private Identification Code</i>	<i>ALOKA:1.2.392.200039.117.2</i>	<i>ANAP</i>	<i>AUTO</i>
(833D,1040)	LO		<i>Private attribute without PHI</i>	<i>ANAP</i>	<i>AUTO</i>
(833D,1042)	UL		<i>Private attribute without PHI</i>	<i>ANAP</i>	<i>AUTO</i>
(833D,104F)	OB		<i>Private attribute without PHI</i>	<i>ANAP</i>	<i>AUTO</i>
(833D,105F)	OB		<i>Private attribute without PHI</i>	<i>ANAP</i>	<i>AUTO</i>
(833D,106F)	OB		<i>Private attribute without PHI</i>	<i>ANAP</i>	<i>AUTO</i>
(833D,1070)	SQ		<i>Private attribute without PHI</i>	<i>ANAP</i>	<i>AUTO</i>
>(833D,107F)	OB		<i>Private attribute without PHI</i>	<i>ANAP</i>	<i>AUTO</i>
(833D,108F)	OB		<i>Private attribute without PHI</i>	<i>ANAP</i>	<i>AUTO</i>
(833D,109F)	OB		<i>Private attribute without PHI</i>	<i>ANAP</i>	<i>AUTO</i>

(833D,10AF)	OB		Private attribute without PHI	ANAP	AUTO
(833D,2000)	OB		Private attribute without PHI	ANAP	AUTO
(833D,2002)	OB		Private attribute without PHI	ANAP	AUTO
(833D,2004)	OB		Private attribute without PHI	ANAP	AUTO
(833D,2050)	SQ		Private attribute without PHI	ANAP	AUTO
>(833D,205E)	OB		Private attribute without PHI	ANAP	AUTO
>(833D,205F)	OB		Private attribute without PHI	ANAP	AUTO
(833D,2070)	OB		Private attribute without PHI	ANAP	AUTO
(833D,2072)	OB		Private attribute without PHI	ANAP	AUTO
(833D,2074)	OB		Private attribute without PHI	ANAP	AUTO
(83FF,0010)	LO	Private Identification Code	ALOKA:1.2.392.200039.117.2	ANAP	AUTO
(83FF,0020)	LO	Private Identification Code	ALOKA:1.2.392.200039.117.2	ANAP	AUTO
(83FF,0021)	LO	Private Identification Code	ALOKA:1.2.392.200039.117.2	ANAP	AUTO
(83FF,1000)	US		Private attribute without PHI	ANAP	AUTO
(83FF,1002)	UL		Private attribute without PHI	ANAP	AUTO
(83FF,1004)	US		Private attribute without PHI	ANAP	AUTO
(83FF,1006)	UL		Private attribute without PHI	ANAP	AUTO
(83FF,1008)	UL		Private attribute without PHI	ANAP	AUTO
(83FF,100A)	UL		Private attribute without PHI	ANAP	AUTO
(83FF,1050)	SH		Private attribute without PHI	ANAP	AUTO
(83FF,1060)	SQ		Private attribute without PHI	ANAP	AUTO
>(83FF,1061)	US		Private attribute without PHI	ANAP	AUTO
>(83FF,1062)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,2000)	DS		Private attribute without PHI	ANAP	AUTO
>(83FF,2002)	DS		Private attribute without PHI	ANAP	AUTO
>(83FF,2004)	UL		Private attribute without PHI	ANAP	AUTO
>(83FF,2006)	UL		Private attribute without PHI	ANAP	AUTO
>(83FF,2008)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,200A)	US		Private attribute without PHI	ANAP	AUTO
>(83FF,200C)	US		Private attribute without PHI	ANAP	AUTO
>(83FF,200E)	US		Private attribute without PHI	ANAP	AUTO
>(83FF,2010)	US		Private attribute without PHI	ANAP	AUTO
>(83FF,2012)	US		Private attribute without PHI	ANAP	AUTO
>(83FF,2014)	US		Private attribute without PHI	ANAP	AUTO
>(83FF,2016)	DS		Private attribute without PHI	ANAP	AUTO
>(83FF,2018)	US		Private attribute without PHI	ANAP	AUTO

>(83FF,201A)	FD		Private attribute without PHI	ANAP	AUTO
>(83FF,201C)	FD		Private attribute without PHI	ANAP	AUTO
>(83FF,2020)	FD		Private attribute without PHI	ANAP	AUTO
>(83FF,2022)	FD		Private attribute without PHI	ANAP	AUTO
>(83FF,2024)	FD		Private attribute without PHI	ANAP	AUTO
>(83FF,2026)	FD		Private attribute without PHI	ANAP	AUTO
>(83FF,2028)	FD		Private attribute without PHI	ANAP	AUTO
>(83FF,202A)	FD		Private attribute without PHI	ANAP	AUTO
>(83FF,202C)	FD		Private attribute without PHI	ANAP	AUTO
>(83FF,202E)	FD		Private attribute without PHI	ANAP	AUTO
>(83FF,2030)	FD		Private attribute without PHI	ANAP	AUTO
>(83FF,2032)	FD		Private attribute without PHI	ANAP	AUTO
>(83FF,2034)	FD		Private attribute without PHI	ANAP	AUTO
>(83FF,2036)	FD		Private attribute without PHI	ANAP	AUTO
>(83FF,2038)	FL		Private attribute without PHI	ANAP	AUTO
>(83FF,203A)	FL		Private attribute without PHI	ANAP	AUTO
>(83FF,203C)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,203E)	SS		Private attribute without PHI	ANAP	AUTO
>(83FF,2040)	FL		Private attribute without PHI	ANAP	AUTO
>(83FF,2042)	US		Private attribute without PHI	ANAP	AUTO
>(83FF,2044)	US		Private attribute without PHI	ANAP	AUTO
>(83FF,2046)	US		Private attribute without PHI	ANAP	AUTO
>(83FF,2048)	US		Private attribute without PHI	ANAP	AUTO
>(83FF,204A)	SS		Private attribute without PHI	ANAP	AUTO
>(83FF,204C)	SS		Private attribute without PHI	ANAP	AUTO
>(83FF,204E)	FD		Private attribute without PHI	ANAP	AUTO
>(83FF,2050)	FD		Private attribute without PHI	ANAP	AUTO
>(83FF,2052)	US		Private attribute without PHI	ANAP	AUTO
>(83FF,2054)	US		Private attribute without PHI	ANAP	AUTO
>(83FF,2056)	US		Private attribute without PHI	ANAP	AUTO
>(83FF,2058)	US		Private attribute without PHI	ANAP	AUTO
>(83FF,205A)	FL		Private attribute without PHI	ANAP	AUTO
>(83FF,205C)	US		Private attribute without PHI	ANAP	AUTO
>(83FF,205E)	FL		Private attribute without PHI	ANAP	AUTO
>(83FF,2060)	US		Private attribute without PHI	ANAP	AUTO
>(83FF,2062)	US		Private attribute without PHI	ANAP	AUTO
>(83FF,2064)	US		Private attribute without PHI	ANAP	AUTO
>(83FF,2066)	US		Private attribute without PHI	ANAP	AUTO

>(83FF,2068)	US		Private attribute without PHI	ANAP	AUTO
>(83FF,206A)	DS		Private attribute without PHI	ANAP	AUTO
>(83FF,206C)	DS		Private attribute without PHI	ANAP	AUTO
>(83FF,206E)	US		Private attribute without PHI	ANAP	AUTO
>(83FF,2070)	DS		Private attribute without PHI	ANAP	AUTO
>(83FF,2072)	DS		Private attribute without PHI	ANAP	AUTO
>(83FF,2074)	US		Private attribute without PHI	ANAP	AUTO
>(83FF,2076)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,2078)	FL		Private attribute without PHI	ANAP	AUTO
>(83FF,207A)	FL		Private attribute without PHI	ANAP	AUTO
>(83FF,207C)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,207E)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,2080)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,2082)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,2084)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,2086)	SL		Private attribute without PHI	ANAP	AUTO
>(83FF,2088)	SL		Private attribute without PHI	ANAP	AUTO
>(83FF,208A)	FL		Private attribute without PHI	ANAP	AUTO
>(83FF,208C)	FL		Private attribute without PHI	ANAP	AUTO
>(83FF,208E)	FL		Private attribute without PHI	ANAP	AUTO
>(83FF,2090)	FL		Private attribute without PHI	ANAP	AUTO
>(83FF,2092)	US		Private attribute without PHI	ANAP	AUTO
>(83FF,2094)	SS		Private attribute without PHI	ANAP	AUTO
>(83FF,2096)	SL		Private attribute without PHI	ANAP	AUTO
>(83FF,2098)	FL		Private attribute without PHI	ANAP	AUTO
>(83FF,209A)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,209C)	SL		Private attribute without PHI	ANAP	AUTO
>(83FF,209E)	UL		Private attribute without PHI	ANAP	AUTO
>(83FF,20A0)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,20A2)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,20A4)	US		Private attribute without PHI	ANAP	AUTO
>(83FF,20A6)	FD		Private attribute without PHI	ANAP	AUTO
>(83FF,20A8)	FD		Private attribute without PHI	ANAP	AUTO
>(83FF,20AA)	FD		Private attribute without PHI	ANAP	AUTO
>(83FF,20AC)	US		Private attribute without PHI	ANAP	AUTO
>(83FF,20AE)	US		Private attribute without PHI	ANAP	AUTO
>(83FF,20B0)	UL		Private attribute without PHI	ANAP	AUTO
>(83FF,20B2)	FL		Private attribute without PHI	ANAP	AUTO

>(83FF,20B4)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,20B6)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,20B8)	US		Private attribute without PHI	ANAP	AUTO
>(83FF,20B9)	UL		Private attribute without PHI	ANAP	AUTO
>(83FF,20BA)	FD		Private attribute without PHI	ANAP	AUTO
>(83FF,20BC)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,20BD)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,20BE)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,20BF)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,20C0)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,20C1)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,20C2)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,20C3)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,20C4)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,2100)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,2101)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,2102)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,2103)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,2104)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,2105)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,2106)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,2107)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,2108)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,2109)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,210A)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,210B)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,210C)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,210D)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,210E)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,2110)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,2112)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,2114)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,2116)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,2118)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,211A)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,211C)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,211E)	OB		Private attribute without PHI	ANAP	AUTO
>(83FF,2120)	OB		Private attribute without PHI	ANAP	AUTO

>(83FF,216E)	OB		Private attribute without PHI	ANAP	AUTO
(83FF,1070)	OB		Private attribute without PHI	ANAP	AUTO
(83FF,1072)	UL		Private attribute without PHI	ANAP	AUTO
(83FF,1074)	FD		Private attribute without PHI	ANAP	AUTO
(83FF,1076)	US		Private attribute without PHI	ANAP	AUTO
(83FF,1078)	OB		Private attribute without PHI	ANAP	AUTO
(83FF,10A0)	SH		Private attribute without PHI	ANAP	AUTO
(83FF,10B1)	SQ		Private attribute without PHI	ANAP	AUTO
(83FF,10B2)	SQ		Private attribute without PHI	ANAP	AUTO
(83FF,10B3)	SQ		Private attribute without PHI	ANAP	AUTO
(83FF,10B4)	SQ		Private attribute without PHI	ANAP	AUTO
(83FF,10C0)	SQ		Private attribute without PHI	ANAP	AUTO
(83FF,10C1)	SQ		Private attribute without PHI	ANAP	AUTO
>(83FF,10FF)	OB		Private attribute without PHI	ANAP	AUTO

8.1.1.7 SR Document Modules

Table 8.1-21 SR DOCUMENT SERIES MODULE OF CREATED SOP INSTANCES

Tag	VR	Attribute Name	Value	Presence of Value	Source
(0008,0021)	DA	Series Date	"yyyymmdd"	ANAP ¹	AUTO
(0008,0031)	TM	Series Time	"hhmmss"	ANAP ¹	AUTO
(0008,0060)	CS	Modality	"SR"	ALWAYS	AUTO
(0008,1111)	SQ	Referenced Performed Procedure Step Sequence	Identifies the MPPS SOP Instance to which this SOP instance is related	VNAP	MPPS
>(0008,1150)	UI	Referenced SOP Class UID	MPPS SOP Class UID	ALWAYS	MPPS
>(0008,1155)	UI	Referenced SOP Instance UID	MPPS SOP Instance UID	ALWAYS	MPPS
(0020,000E)	UI	Series Instance UID	Generated by device	ALWAYS	AUTO
(0020,0011)	IS	Series Number	Generated by device	ALWAYS	AUTO

Note: 1. Attribute presents only for Advanced Cardiac feature.

Table 8.1-22 SR DOCUMENT GENERAL MODULE OF CREATED SOP INSTANCES

Tag	VR	Attribute Name	Value	Presence of Value	Source
(0008,0023)	DA	Content Date	"yyyymmdd"	ALWAYS	AUTO
(0008,0033)	TM	Content Time	"hhmmss"	ALWAYS	AUTO
(0020,0013)	IS	Instance Number	Generated by device	ALWAYS	AUTO
(0040,A370)	SQ	Referenced Request Sequence	From Modality Worklist, or absent if unscheduled case	ANAP	AUTO

>(0008,0050)	SH	Accession Number	From Modality Worklist	VNAP	MWL
>(0008,1110)	SQ	Referenced Study Sequence	From Modality Worklist	VNAP	MWL
>>(0008,1150)	UI	Referenced SOP Class UID	From Modality Worklist	ALWAYS	MWL
>>(0008,1155)	UI	Referenced SOP Instance UID	From Modality Worklist	ALWAYS	MWL
>(0020,000D)	UI	Study Instance UID	From Modality Worklist	ALWAYS	MWL
>(0032,1060)	LO	Requested Procedure Description	From Modality Worklist	VNAP	MWL
>(0032,1064)	SQ	Requested Procedure Code Sequence	From Modality Worklist	VNAP ¹	MWL
>(0040,1001)	SH	Requested Procedure ID	From Modality Worklist	VNAP	MWL
>(0040,2016)	LO	Placer Order Number / Imaging Service Request	Zero length	EMPTY	AUTO
>(0040,2017)	LO	Filler Order Number / Imaging Service Request	Zero length	EMPTY	AUTO
(0040,A372)	SQ	Performed Procedure Code Sequence	From Modality Worklist, mapped from Requested Procedure Code Sequence (0032,1064)	VNAP	MWL
(0040,A491)	CS	Completion Flag	PARTIAL	ALWAYS	AUTO
(0040,A493)	CS	Verification Flag	UNVERIFIED	ALWAYS	AUTO

Note: 1. Attribute Not Present for Advanced Cardiac feature.

Table 8.1-23 SR DOCUMENT CONTENT MODULE OF CREATED SOP INSTANCES

Tag	VR	Attribute Name	Value	Presence of Value	Source
(0040,A040)	CS	Value Type	CONTAINER	ALWAYS	AUTO
(0040,A043)	SQ	Concept Name Code Sequence	Document Title	ALWAYS	AUTO
>(0008,0100)	SH	Code Value	125000, 125100 or 125200	ALWAYS	AUTO
>(0008,0102)	SH	Coding Scheme Designator	DCM	ALWAYS	AUTO
>(0008,0104)	LO	Code Meaning	"OB-GYN Ultrasound Procedure Report", "Vascular Ultrasound Procedure Report" or "Adult Echocardiography Procedure Report"	ALWAYS	AUTO
(0040,A050)	CS	Continuity of Content	SEPARATE	ALWAYS	AUTO
(0040,A504)	SQ	Content Template Sequence		ALWAYS	AUTO
>(0008,0105)	CS	Mapping Resource	DCMR	ALWAYS	AUTO
>(0040,DB00)	CS	Template Identifier	5000, 5100 or 5200	ALWAYS	AUTO

(0040,A032)	DT	Observation Date Time	"yyyymmddhhmmss" Same as Study Date, Study Time	ALWAYS	AUTO
(0040,A730)	SQ	Content Sequence		ALWAYS	AUTO
>(0040,A010)	CS	Relationship Type	See TID 5000 OB-GYN Ultrasound Procedure Report, TID 5100 Vascular Ultrasound Report and TID 5200 Echocardiography Procedure Report	ALWAYS	AUTO
> ...		Document Relationship Macro	See TID 5000 OB-GYN Ultrasound Procedure Report, TID 5100 Vascular Ultrasound Report and TID 5200 Echocardiography Procedure Report	ALWAYS	AUTO
> ...		Document Content Macro	See TID 5000 OB-GYN Ultrasound Procedure Report, TID 5100 Vascular Ultrasound Report and TID 5200 Echocardiography Procedure Report	ALWAYS	AUTO
>(0040,DB73)	UL	Referenced Content Item Identifier	Not used	-	-

8.1.1.8 Encapsulated Document Modules

Table 8.1-24 ENCAPSULATED DOCUMENT SERIES OF CREATED SOP INSTANCES

Tag	VR	Attribute Name	Value	Presence of Value	Source
(0008,0060)	CS	Modality	"DOC"	ALWAYS	AUTO
(0008,1111)	SQ	Referenced Performed Procedure Step Sequence	Identifies the MPPS SOP Instance to which this document is related	VNAP	MPPS
>(0008,1150)	UI	Referenced SOP Class UID	MPPS SOP Class UID	ALWAYS	MPPS
>(0008,1155)	UI	Referenced SOP Instance UID	MPPS SOP Instance UID	ALWAYS	MPPS
(0020,000E)	UI	Series Instance UID	Generated by device	ALWAYS	AUTO
(0020,0011)	IS	Series Number	Generated by device	ALWAYS	AUTO
(0040,0275)	SQ	Request Attributes Sequence	Zero or 1 item will be present	VNAP	AUTO
>(0040,0007)	LO	Scheduled Procedure Step Description	From Modality Worklist	VNAP	MWL
>(0040,0009)	SH	Scheduled Procedure Step ID	From Modality Worklist	VNAP	MWL
>(0040,1001)	SH	Requested Procedure ID	From Modality Worklist or Input by user	VNAP	MWL/ USER

Table 8.1-25 ENCAPSULATED DOCUMENT OF CREATED SOP INSTANCES

Tag	VR	Attribute Name	Value	Presence of Value	Source
(0008,0023)	DA	Content Date	"yyyymmdd"	ALWAYS	AUTO
(0008,0033)	TM	Content Time	"hhmmss"	ALWAYS	AUTO
(0008,002A)	DT	Acquisition DateTime	Zero length	EMPTY	AUTO
(0020,0013)	IS	Instance Number	Generated by device	ALWAYS	AUTO
(0028,0301)	CS	Burned In Annotation	"YES"	ALWAYS	AUTO
(0040,A043)	SQ	Concept Name Code Sequence	Zero length	EMPTY	AUTO
(0042,0010)	ST	Document Title	"Report"	VNAP	AUTO
(0042,0011)	OB	Encapsulated Document	Encapsulated document stream encoded by MIME Type of Encapsulated Document (0040,0012).	ALWAYS	AUTO
(0042,0012)	LO	MIME Type of Encapsulated Document	"application/pdf"	ALWAYS	AUTO

8.1.1.9 Secondary Capture Modules**Table 8.1-26 SC EQUIPMENT MODULE OF CREATED SOP INSTANCES**

Tag	VR	Attribute Name	Value	Presence of Value	Source
(0008,0064)	CS	Conversion Type	"WSD"	ALWAYS	AUTO
(0008,0060)	CS	Modality	"OT"	ALWAYS	AUTO

Table 8.1-27 SC IMAGE MODULE OF CREATED SOP INSTANCES

Tag	VR	Attribute Name	Value	Presence of Value	Source
(0009,0010)	LO	Private Identification Code	ALOKA:1.2.392.200039.11 7.2	ANAP	AUTO
(0009,1000)	SH		Private attribute without PHI	ANAP	AUTO
(0009,1012)	US		Private attribute without PHI	ANAP	AUTO
(0018,1012)	DA	Date of Secondary Capture	"yyyymmdd"	ALWAYS	AUTO
(0018,1014)	TM	Time of Secondary Capture	"hhmmss"	ALWAYS	AUTO
(7FDF,0050)	LO	Private Identification Code	"TomTec"	ANAP	AUTO
(7FDF,5050)	OB		Private attribute without PHI	ANAP	AUTO
(7FDF,5051)	OB		Private attribute without PHI	ANAP	AUTO

Table 8.1-28 SC MULTI-FRAME IMAGE MODULE OF CREATED SOP INSTANCES

Tag	VR	Attribute Name	Value	Presence of Value	Source
(0028,0301)	CS	Burned In Annotation	"NO"	ALWAYS	AUTO
(0028,0009)	AT	Frame Increment Pointer	(0018,1063)	ANAP	AUTO

8.1.2 Used Fields in received IOD by application

The Query/Retrieve application of this product may receive SOP Instances described in section 4.2.4.1. To "Import" the received SOP Instance, value of Patient ID (0010, 0020), Study ID (0020,0010), Series Number (0020,0011), and Instance Number (0020,0013) are necessary in it. The patient identification in the local database is done by Patient ID (0010,0020).

The usage of attributes received via Modality Worklist is described in section 4.2.1.3.1.3.

8.1.3 Attribute mapping

The relationships between attributes received via Modality Worklist, stored in acquired images and communicated via MPPS are summarized in Table 8.1-29.

Table 8.1-29 ATTRIBUTE MAPPING BETWEEN MODALITY WORKLIST, IMAGE AND MPPS

Modality Worklist	Image IOD	MPPS IOD
Patient Name	Patient Name	Patient Name
Patient ID	Patient ID	Patient ID
Patient's Birth Date	Patient's Birth Date	Patient's Birth Date
Patient's Sex	Patient's Sex	Patient's Sex
Patient's Size	Patient's Size	----
Patient's Weight	Patient's Weight	----
Referring Physician's Name	Referring Physician's Name	
----	----	Scheduled Step Attributes Sequence
Study Instance UID	Study Instance UID	>Study Instance UID
Referenced Study Sequence	Referenced Study Sequence	>Referenced Study Sequence
Accession Number	Accession Number	>Accession Number
----	Request Attributes Sequence	----
Requesting Physician	Requesting Physician	----
Requested Procedure ID	>Requested Procedure ID	>Requested Procedure ID
Requested Procedure Description	>Requested Procedure Description	>Requested Procedure Description
Scheduled Procedure Step ID	>Scheduled Procedure Step ID	>Scheduled Procedure Step ID
Scheduled Procedure Step Description	>Scheduled Procedure Step Description	>Scheduled Procedure Step Description
Scheduled Protocol Code Sequence	>Scheduled Protocol Code Sequence	>Scheduled Protocol Code Sequence
----	Performed Protocol Code Sequence	Performed Protocol Code Sequence

Modality Worklist	Image IOD	MPPS IOD
Requested Procedure ID	Study ID	Study ID
----	Performed Procedure Step ID	Performed Procedure Step ID
----	Performed Procedure Step Start Date	Performed Procedure Step Start Date
----	Performed Procedure Step Start Time	Performed Procedure Step Start Time
----	Performed Procedure Step Description	Performed Procedure Step Description
----	----	Performed Series Sequence
----	Name of Physician's Reading Study	>Performing Physician's Name
Requested Procedure Code Sequence	Procedure Code Sequence	Procedure Code Sequence
----	Referenced Performed Procedure Step Sequence	----
----	>Referenced SOP Class UID	SOP Class UID
----	>Referenced SOP Instance UID	SOP Instance UID
----	Protocol Name	Protocol Name
Medical Alerts	Medical Alerts	----
Contrast Allergies	Contrast Allergies	----
Occupation	Occupation	----

8.1.4 Coerced/Modified Fields

The Modality Worklist AE will truncate attribute values received in the response to a Modality Worklist Query if the value length is longer than the maximum length permitted by the attribute's VR.

8.2 DATA DICTIONARY OF PRIVATE ATTRIBUTES

The Private Attributes added to create SOP Instances are listed in the Table below. This product reserves blocks of private attributes in groups 0009, 0019, 830F, 833D and 83FF. Further details on usage of these private attributes are contained in Section 8.1. Groups 830F, 833D, 83FF are used in DICOM+Raw and 3D Reconstruction files.

This product does not support reading SOP instance that the group of Private Attributes has been relocated

Table 8.2-1 DATA DICTIONARY OF PRIVATE ATTRIBUTES

Tag	VR	VM	Attribute Name
(0009,0010)	LO	1	
(0009,1000)	SH	1	
(0009,1004)	US	1-n	
(0009,1006)	US	1-n	
(0009,100A)	SH	1	
(0009,1012)	US	1	

Tag	VR	VM	Attribute Name
(0009,1014)	FD	1	
(0009,1016)	DA	1	
(0009,1018)	TM	1	
(0009,101A)	LO	1	
(0009,1020)	CS	1	
(0009,1022)	CS	1	
(0009,1024)	CS	1	
(0009,1026)	IS	1	
(0009,1028)	IS	1	
(0009,102A)	DS	1	
(0009,102C)	FD	6-6n	
(0009,102E)	FD	3-3n	
(0009,1030)	FD	1	
(0009,1032)	DS	1	
(0009,1034)	CS	1	
(0009,1036)	FD	1	
(0009,1038)	CS	1	
(0009,103A)	FD	1	
(0009,103C)	FD	1	
(0009,103E)	FD	1	
(0009,1040)	FD	1	
(0009,1042)	FD	1	
(0009,1044)	FD	1	
(0009,1046)	US	1	
(0009,1048)	FD	1	
(0009,104A)	FD	1	
(0009,104C)	FD	1	
(0009,104E)	FD	1	
(0009,1050)	US	1	
(0009,1052)	US	1	
(0009,1054)	LO	1	
(0009,1056)	IS	1-n	
(0009,1058)	CS	1	
(0009,105A)	UL	1	
(0009,105C)	CS	1	
(0009,105E)	DA	1	
(0009,1060)	TM	1	
(0009,1062)	US	1	

Tag	VR	VM	Attribute Name
(0019,0010)	LO	1	
(0019,1008)	FD	1	
(0019,100C)	CS	1	
(0019,100E)	DS	1	
(0019,1018)	SL	1	
(0019,101A)	SL	1	
(0019,1040)	SS	1	
(0019,1046)	US	1	
(0019,1050)	SL	1	
(0019,1052)	DS	1	
(0019,1054)	DS	1	
(0019,1056)	FD	1	
(0019,1060)	US	1	
(0019,1061)	UL	1	
(7FDF,0050)	LO	1	
(7FDF,5050)	OB	1	
(7FDF,5051)	OB	1	
(830F,0010)	LO	1	
(830F,1000)	US	1	
(830F,10FF)	OB	1	
(833D,0010)	LO	1	
(833D,0020)	OB	1	
(833D,1040)	LO	1	
(833D,1042)	UL	1	
(833D,104F)	OB	1	
(833D,105F)	OB	1	
(833D,106F)	OB	1	
(833D,1070)	SQ	1	
>(833D,107F)	OB	1	
(833D,108F)	OB	1	
(833D,109F)	OB	1	
(833D,10AF)	OB	1	
(833D,2000)	OB	1	
(833D,2002)	OB	1	
(833D,2004)	OB	1	
(833D,2050)	SQ	1	
>(833D,205E)	OB	1	
>(833D,205F)	OB	1	

Tag	VR	VM	Attribute Name
(833D,2070)	OB	1	
(833D,2072)	OB	1	
(833D,2074)	OB	1	
(83FF,0010)	LO	1	
(83FF,0020)	LO	1	
(83FF,0021)	LO	1	
(83FF,1000)	US	1	
(83FF,1002)	UL	1	
(83FF,1004)	US	1	
(83FF,1006)	UL	1	
(83FF,1008)	UL	4	
(83FF,100A)	UL	1	
(83FF,1050)	SH	1	
(83FF,1060)	SQ	1	
(83FF,1061)	US	1	
(83FF,1062)	OB	1	
(83FF,1070)	OB	1	
(83FF,1072)	UL	1	
(83FF,1074)	FD	1	
(83FF,1076)	US	1	
(83FF,1078)	OB	1	
(83FF,10A0)	SH	1	
(83FF,10B1)	SQ	1	
(83FF,10B2)	SQ	1	
(83FF,10B3)	SQ	1	
(83FF,10B4)	SQ	1	
(83FF,10C0)	SQ	1	
(83FF,10C1)	SQ	1	
(83FF,10FF)	OB	1	
(83FF,2000)	DS	1	
(83FF,2002)	DS	1	
(83FF,2004)	UL	1	
(83FF,2006)	UL	1	
(83FF,2008)	OB	1	
(83FF,200A)	US	1	
(83FF,200C)	US	1	
(83FF,200E)	US	1	
(83FF,2010)	US	1	

Tag	VR	VM	Attribute Name
(83FF,2012)	US	1	
(83FF,2014)	US	1	
(83FF,2016)	DS	1	
(83FF,2018)	US	1	
(83FF,201A)	FD	1	
(83FF,201C)	FD	1	
(83FF,2020)	FD	1	
(83FF,2022)	FD	1	
(83FF,2024)	FD	1	
(83FF,2026)	FD	1	
(83FF,2028)	FD	7	
(83FF,202A)	FD	7	
(83FF,202C)	FD	7	
(83FF,202E)	FD	7	
(83FF,2030)	FD	7	
(83FF,2032)	FD	7	
(83FF,2034)	FD	7	
(83FF,2036)	FD	7	
(83FF,2038)	FL	1	
(83FF,203A)	FL	256	
(83FF,203C)	OB	1	
(83FF,203E)	SS	1	
(83FF,2040)	FL	1	
(83FF,2042)	US	1	
(83FF,2044)	US	1	
(83FF,2046)	US	1	
(83FF,2048)	US	1	
(83FF,204A)	SS	1	
(83FF,204C)	SS	1	
(83FF,204E)	FD	1	
(83FF,2050)	FD	1	
(83FF,2052)	US	4	
(83FF,2054)	US	4	
(83FF,2056)	US	4	
(83FF,2058)	US	4	
(83FF,205A)	FL	1	
(83FF,205C)	US	1	
(83FF,205E)	FL	1-n	

Tag	VR	VM	Attribute Name
(83FF,2060)	US	1-n	
(83FF,2062)	US	1-n	
(83FF,2064)	US	1-n	
(83FF,2066)	US	1	
(83FF,2068)	US	1-n	
(83FF,206A)	DS	1	
(83FF,206C)	DS	1	
(83FF,206E)	US	1	
(83FF,2070)	DS	1	
(83FF,2072)	DS	1	
(83FF,2074)	US	1	
(83FF,2076)	OB	1	
(83FF,2078)	FL	16	
(83FF,207A)	FL	16	
(83FF,207C)	OB	1	
(83FF,207E)	OB	1	
(83FF,2080)	OB	1	
(83FF,2082)	OB	1	
(83FF,2084)	OB	1	
(83FF,2086)	SL	1	
(83FF,2088)	SL	1	
(83FF,208A)	FL	1	
(83FF,208C)	FL	1	
(83FF,208E)	FL	7	
(83FF,2090)	FL	7	
(83FF,2092)	US	9	
(83FF,2094)	SS	9	
(83FF,2096)	SL	1	
(83FF,2098)	FL	1	
(83FF,209A)	OB	1	
(83FF,209C)	SL	1	
(83FF,209E)	UL	3	
(83FF,20A0)	OB	1	
(83FF,20A2)	OB	1	
(83FF,20A4)	US	1	
(83FF,20A6)	FD	1	
(83FF,20A8)	FD	1	
(83FF,20AA)	FD	1	

Tag	VR	VM	Attribute Name
(83FF,20AC)	US	1	
(83FF,20AE)	US	2	
(83FF,20B0)	UL	4	
(83FF,20B2)	FL	1	
(83FF,20B4)	OB	1	
(83FF,20B6)	OB	1	
(83FF,20B8)	US	1	
(83FF,20B9)	UL	1	
(83FF,20BA)	FD	16	
(83FF,20BC)	OB	1	
(83FF,20BD)	OB	3	
(83FF,20BE)	OB	1	
(83FF,20BF)	OB	1	
(83FF,20C0)	OB	1	
(83FF,20C1)	OB	1	
(83FF,20C2)	OB	1	
(83FF,20C3)	OB	1	
(83FF,20C4)	OB	1	
(83FF,2100)	OB	1	
(83FF,2101)	OB	1	
(83FF,2102)	OB	1	
(83FF,2103)	OB	1	
(83FF,2104)	OB	1	
(83FF,2105)	OB	1	
(83FF,2106)	OB	1	
(83FF,2107)	OB	1	
(83FF,2108)	OB	1	
(83FF,2109)	OB	1	
(83FF,210A)	OB	1	
(83FF,210B)	OB	1	
(83FF,210C)	OB	1	
(83FF,210D)	OB	1	
(83FF,210E)	OB	1	
(83FF,2110)	OB	1	
(83FF,2112)	OB	1	
(83FF,2114)	OB	1	
(83FF,2116)	OB	1	
(83FF,2118)	OB	1	

Tag	VR	VM	Attribute Name
(83FF,211A)	OB	1	
(83FF,211C)	OB	1	
(83FF,211E)	OB	1	
(83FF,2120)	OB	1	
(83FF,2122)	OB	1	
(83FF,2124)	OB	1	
(83FF,2126)	OB	1	
(83FF,2128)	OB	1	
(83FF,212A)	OB	1	
(83FF,212C)	OB	1	
(83FF,212E)	OB	1	
(83FF,2130)	OB	1	
(83FF,2132)	OB	1	
(83FF,2134)	OB	1	
(83FF,2136)	OB	1	
(83FF,2138)	OB	1	
(83FF,213A)	OB	1	
(83FF,213C)	OB	1	
(83FF,213E)	OB	1	
(83FF,2140)	OB	1	
(83FF,2142)	OB	1	
(83FF,2144)	OB	1	
(83FF,2146)	OB	1	
(83FF,2148)	OB	1	
(83FF,214A)	OB	1	
(83FF,214C)	OB	1	
(83FF,214E)	OB	1	
(83FF,2150)	OB	1	
(83FF,2152)	OB	1	
(83FF,2154)	OB	1	
(83FF,2156)	OB	1	
(83FF,2158)	OB	1	
(83FF,215A)	OB	1	
(83FF,215C)	OB	1	
(83FF,215E)	OB	1	
(83FF,2160)	OB	1	
(83FF,2162)	OB	1	
(83FF,2164)	OB	1	

Tag	VR	VM	Attribute Name
(83FF,2166)	OB	1	
(83FF,2168)	OB	1	
(83FF,216A)	OB	1	
(83FF,216C)	OB	1	
(83FF,216E)	OB	1	

8.3 CODED TERMINOLOGY AND TEMPLATES

The contents of Performed Procedure Step Discontinuation Reason Code Sequence (0040,0281) for a discontinued MPPS will be filled with a code selected by the user from a fixed list corresponding to Context Group 9300.

The Structured Reports use the Standard Templates and Context Groups supplied by DCMR (DICOM Content Mapping Resource), LOINC, and SRT. As a private definition, "99IMPROOF" for Advanced Cardiac feature and "99ALOKA" for other features are used. See section 8.6 for details.

8.4 STANDARD EXTENDED / SPECIALIZED / PRIVATE SOP CLASSES

8.4.1 Ultrasound and Ultrasound Multi-frame Image Storage SOP Class

The Ultrasound and Ultrasound Multi-frame Image Storage SOP Classes are extended to create a Standard Extended SOP Classes by addition of standard and private attributes to the created SOP Instances as documented in section 8.1.

8.4.2 Comprehensive SR Storage SOP Class

The Comprehensive SR Storage SOP Class which is generated as "Advanced Cardiac Report" is extended to create a Standard Extended SOP Class by addition of standard and private attributes to the created SOP Instances as documented in section 8.1.

8.5 PRIVATE TRANSFER SYNTAXES

No Private Transfer Syntaxes are supported.

8.6 STRUCTURED REPORTS

The equipment supports Standard Templates TID 5000 OB-GYN, TID 5100 Vascular and TID 5200 Echocardiography. TIDs are described in DICOM Part 16.

There is no attempt to provide all tags that will be sent in the SR object. For those specifics, please refer to DICOM PS3.3, PS3.16 and PS3.17 for complete discussion on Enhanced and/or Comprehensive SR Storage SOP Classes support and the Templates described here.

DICOM PS 3.17 of the DICOM Standard includes tree diagrams showing graphic examples of the structure of each template.

8.6.1 Applications and Generated Templates

Application	Template ID	Template Name
OB	5000	OB-GYN Ultrasound Procedure Report
GYN	5000	OB-GYN Ultrasound Procedure Report
Vascular	5100	Vascular Ultrasound Report
Abdomen	5100	Vascular Ultrasound Report
Cardio	5200	Echocardiography Procedure Report
other applications	SR is not generated.	-

OB application generates an OB-GYN Ultrasound Procedure Report for OB measurements, GYN application generates an OB-GYN Ultrasound Procedure Report for GYN measurements. Vascular application generates a Vascular Ultrasound Report for peripheral vascular measurements, Abdomen application generates a Vascular Ultrasound Report for abdominal vascular measurements. Cardio application generates an Echocardiography Procedure Report. No SR is generated for any other applications.

Note that DICOM Standard provides TID 5200 for an "Adult" Echocardiography Procedure Report, however Cardio application will generate TID 5200 Adult Echocardiography Procedure Report regardless of patient's age or the Preset Application like "Pediatric Heart".

This product is able to generate two types of TID 5200 Structured Report below:

1. Non-Advanced Cardiac feature:
Some user-defined measurement items or user-defined tables/equations are not included in SR. All private Codes Values use the Coding Scheme Designator "99ALOKA".
2. Advanced Cardiac feature:
No user-defined measurement item is included in SR. All private Code Values use the Coding Scheme Designator "99IMPROOF".

Some user-defined measurement items or user-defined tables/equations are not included in SR.

8.6.1.1 Image and Spatial Coordinates

It is configurable to output the referenced single images and spatial coordinates in distance-typed measurement results except for Advanced Cardiac feature. The referenced image and spatial coordinates are confirmed by pressing the "STORE" button on the image whose the measurements have been done. TID 320 row 3 and row 4 are used to store spatial coordinates and image respectively. The

Value Multiplicity for TID 320 is up to 10 on each result. The value of "POLYLINE" is used for Graphic Type (0070,0023).

8.6.2 Templates

8.6.2.1 TID 5000 OB-GYN Ultrasound Procedure Report

TID 5000 row 3 - TID 1001 Observation Context

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
Patient Name	DCM	121029	Subject Name	TID 1001 row 3 - TID 1006 row 2 - TID 1007 row 2

TID 5000 row 4 - TID 5001 OB-GYN Patient Characteristics

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
Height	LN	8302-2	Patient Height	TID 5001 row 3
Weight	LN	29463-7	Patient Weight	TID 5001 row 4
GRAV	LN	11996-6	Gravida	TID 5001 row 5
PARA	LN	11977-6	Para	TID 5001 row 6
AB	LN	11612-9	Aborta	TID 5001 row 7
ECTO	LN	33065-4	Ectopic Pregnancies	TID 5001 row 8

TID 5000 row 7 - TID 5002 OB-GYN Procedure Summary Section

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
LMP	LN	11955-2	LMP	TID 5002 row 2
BBT	LN	11976-8	Ovulation date	TID 5002 row 2
LMP-EDC	LN	11779-6	EDD from LMP	TID 5002 row 2
BBT-EDC	LN	11780-4	EDD from ovulation date	TID 5002 row 2
Composite US-EDC ¹	LN	11781-2	EDD from average ultrasound age	TID 5002 row 2
-	LN	11878-6	Number of Fetuses	TID 5002 row 3
Comments	DCM	121106	Comment	TID 5002 row 4
LMP-GA	LN	11885-1	Gestational Age by LMP	TID 5002 row 6 - TID 5003 row 5

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
Composite US-GA	LN	11888-5	Composite Ultrasound Age	TID 5002 row 6 - TID 5003 row 5

Note: 1. The earliest Composite US-EDC is included when more than one fetus are observed.

OB measurements :

TID 5000 row 7 - TID 5002 row 6 - TID 5003 OB-GYN Fetus Summary
 TID 5000 row 8 - TID 5004 Fetal Biometry Ratio Section
 TID 5000 row 9 - TID 5005 Fetal Biometry Section
 TID 5000 row 10 - TID 5006 Fetal Long Bones Section
 TID 5000 row 11 - TID 5007 Fetal Cranium Section
 TID 5000 row 12 - TID 5009 Fetal Biophysical Profile Section
 TID 5000 row 13 - TID 5011 Early Gestation Section
 TID 5000 row 14 - TID 5010 Amniotic Sac Section
 TID 5000 row 15 - TID 5015 Pelvis and Uterus Section
 TID 5000 row 21 - TID 5025 OB-GYN Fetal Vascular Ultrasound Measurement Group
 TID 5000 row 24 - TID 5026 OB-GYN Pelvic Vascular Ultrasound Measurement Group

TID 5000 row 25(Extension) - TID ALOKA-5001 OB-GYN Original Measurement Section
 TID 5000 row 26(Extension) - TID ALOKA-5000 OB Original Fetal Measurement Section
 TID 5000 row 27(Extension) - TID ALOKA-5000 OB Original Fetal Measurement Section
 TID 5000 row 28(Extension) - TID ALOKA-5002 OB-GYN User Defined Equation Section

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
FHR	LN	11948-7	Fetal Heart Rate	TID 5003 row 5
PreHR(Amnio)	99ALOKA	A12019-001	Fetal Heart Rate before Biopsy	TID 5003 row 5
PstHR(Amnio)	99ALOKA	A12019-002	Fetal Heart Rate after Biopsy	TID 5003 row 5
FW (up to 5 items)	LN	11727-5	Estimated Weight	TID 5003 row 5
EFW Ratio (up to 5 items)	99ALOKA	A12019-004	Estimated Fetal Weight Ratio	TID 5003 row 5
- (FW %ile rank by Doubilet)	LN	11767-1	EFW percentile rank	TID 5003 row 5
Auto FHR	99ALOKA	A12019-003	Auto Fetal Heart Rate	TID 5003 row 5
CI (BPD/OFD)	LN	11823-2	Cephalic Index	TID 5004 row 3

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
CI (BPD _o /OFD _o)	99ALOKA	A12004-001	Cephalic Index(BPD _o /OFD _o)	TID 5004 row 3
FL/AC	LN	11871-1	FL/AC	TID 5004 row 3
FL/BPD	LN	11872-9	FL/BPD	TID 5004 row 3
FL/HC	LN	11873-7	FL/HC	TID 5004 row 3
HC/AC	LN	11947-9	HC/AC	TID 5004 row 3
LVW/HW	99ALOKA	A12004-002	LVW/HW	TID 5004 row 3
AC	LN	11979-2	Abdominal Circumference	TID 5005 row 3
AC2	99ALOKA	A12005-011	Abdominal Circumference for BMUS	TID 5005 row 3
AD	99ALOKA	A12005-006	Abdominal Diameter	TID 5005 row 3
AF Pocket	99ALOKA	A12005-009	Amniotic Fluid Volume	TID 5005 row 3
AFV	99ALOKA	A12005-009	Amniotic Fluid Volume	TID 5005 row 3
MVP	99ALOKA	A12005-012	Maximum Vertical Pocket	TID 5005 row 3
APD	LN	11818-2	Anterior-Posterior Abdominal Diameter	TID 5005 row 3
APTD	LN	11818-2	Anterior-Posterior Abdominal Diameter	TID 5005 row 3
AXT	LN	33191-8	APAD * TAD	TID 5005 row 3
BD	99ALOKA	A12005-001	Binocular Distance	TID 5005 row 3
BPD	LN	11820-8	Biparietal Diameter	TID 5005 row 3
BPD _o	99ALOKA	A12005-002	Biparietal Diameter outer-to-outer	TID 5005 row 3
CD	LN	11863-8	Trans Cerebellar Diameter	TID 5007 row 3
CRL	LN	11957-8	Crown Rump Length	TID 5011 row 3
EES	99ALOKA	A12009-001	Early Embryonic Size	TID 5011 row 3
FIB	LN	11964-4	Fibula length	TID 5006 row 3
FL	LN	11963-6	Femur Length	TID 5006 row 3
FTA	99ALOKA	A12005-003	Fetal Trunk Cross Sectional Area	TID 5005 row 3
GS	LN	11850-5	Gestational Sac Diameter	TID 5011 row 3
HC	LN	11984-2	Head Circumference	TID 5005 row 3
HC2	99ALOKA	A12005-008	Head Circumference for Merz, Hansmann	TID 5005 row 3
HC3	99ALOKA	A12005-010	Head Circumference for BMUS	TID 5005 row 3
HL	LN	11966-9	Humerus length	TID 5006 row 3

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
HW	LN	12170-7	Width of Hemisphere	TID 5007 row 3
IOD	LN	33070-4	Inner Orbital Diameter	TID 5007 row 3
LV	99ALOKA	A12005-004	Length of Vertebrae	TID 5005 row 3
LVW	LN	12171-5	Lateral Ventricular width	TID 5007 row 3
mGS, D1	99ALOKA	A12009-002	Gestational Sac Diameter 1	TID 5011 row 3
mGS, D2	99ALOKA	A12009-003	Gestational Sac Diameter 2	TID 5011 row 3
mGS, D3	99ALOKA	A12009-004	Gestational Sac Diameter 3	TID 5011 row 3
mGS, mGS	99ALOKA	A12009-005	Mean Gestational Sac Diameter	TID 5011 row 3
NBL	99ALOKA	A12006-001	Nasal Bone Length	TID 5006 row 3
NT	LN	33069-6	Nuchal Translucency	TID 5011 row 3
OFD	LN	11851-3	Occipital-Frontal Diameter	TID 5005 row 3
OFDo	99ALOKA	A12005-005	Occipital-Frontal Diameter outer-to-outer	TID 5005 row 3
OOD	LN	11629-3	Outer Orbital Diameter	TID 5007 row 3
RAD	LN	11967-7	Radius length	TID 5006 row 3
TAD	LN	11862-0	Tranverse Abdominal Diameter	TID 5005 row 3
TC	LN	11988-3	Thoracic Circumference	TID 5005 row 3
TIB	LN	11968-5	Tibia length	TID 5006 row 3
TL	99ALOKA	A12005-007	Thoracic Length	TID 5005 row 3
TTD	LN	11862-0	Tranverse Abdominal Diameter	TID 5005 row 3
ULNA	LN	11969-3	Ulna length	TID 5006 row 3
BPP, Breathing	LN	11632-7	Fetal Breathing	TID 5009 row 4
BPP, Movement	LN	11631-9	Gross Body Movement	TID 5009 row 3
BPP, Tone	LN	11635-0	Fetal Tone	TID 5009 row 5
BPP, Fluid	LN	11630-1	Amniotic Fluid Volume	TID 5009 row 7
BPP, Non-Stress Test	LN	11635-5	Fetal Heart Reactivity	TID 5009 row 6
BPP, Total Score	LN	11634-3	Biophysical Profile Sum Score	TID 5009 row 8
AFI, AFI ¹	LN	11627-7	Amniotic Fluid Index	TID 5010 row 3
AFI, Q1	LN	11624-4	First Quadrant Diameter	TID 5010 row 4

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
AFI, Q2	LN	11626-9	Second Quadrant Diameter	TID 5010 row 4
AFI, Q3	LN	11625-1	Third Quadrant Diameter	TID 5010 row 4
AFI, Q4	LN	11623-6	Fourth Quadrant Diameter	TID 5010 row 4
Cervix	LN	11961-0	Cervix Length	TID 5015 row 3
MCA	SRT	T-45600	Middle Cerebral Artery	TID 5025 row 1
MCA, PI	LN	12008-9	Pulsatility Index	TID 5025 row 4
MCA, RI	LN	12023-8	Resistivity Index	TID 5025 row 4
MCA, S/D	LN	12144-2	Systolic to Diastolic Velocity Ratio	TID 5025 row 4
MCA, PSV	LN	11726-7	Peak Systolic Velocity	TID 5025 row 4
MCA, EDV	LN	11653-3	End Diastolic Velocity	TID 5025 row 4
MCA, MnV	LN	11692-1	Time averaged peak velocity	TID 5025 row 4
UmA ¹	SRT	T-F1810	Umbilical Artery	TID 5026 row 1, TID ALOKA-5000 row3 - TID ALOKA-300 row1
UmA, PI	LN	12008-9	Pulsatility Index	TID 5026 row 4, TID ALOKA-300 row3
UmA, RI	LN	12023-8	Resistivity Index	TID 5026 row 4, TID ALOKA-300 row3
UmA, S/D	LN	12144-2	Systolic to Diastolic Velocity Ratio	TID 5026 row 4, TID ALOKA-300 row3
UmA, PSV	LN	11726-7	Peak Systolic Velocity	TID 5026 row 4, TID ALOKA-300 row3
UmA, EDV	LN	11653-3	End Diastolic Velocity	TID 5026 row 4, TID ALOKA-300 row3
UmA, MnV	LN	11692-1	Time averaged peak velocity	TID 5026 row 4, TID ALOKA-300 row3
Rt.UtA	SRT	T-46820	Uterine Artery	TID ALOKA-5001 row2 -TID ALOKA-300 row1
Rt.UtA, PI	LN	12008-9	Pulsatility Index	TID ALOKA-300 row3
Rt.UtA, RI	LN	12023-8	Resistivity Index	TID ALOKA-300 row3
Rt.UtA, S/D	LN	12144-2	Systolic to Diastolic Velocity Ratio	TID ALOKA-300 row3

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
Rt.UtA, PSV	LN	11726-7	Peak Systolic Velocity	TID ALOKA-300 row3
Rt.UtA, EDV	LN	11653-3	End Diastolic Velocity	TID ALOKA-300 row3
Rt.UtA, MnV	LN	11692-1	Time averaged peak velocity	TID ALOKA-300 row3
Lt.UtA	SRT	T-46820	Uterine Artery	TID ALOKA-5001 row2 -TID ALOKA-300 row1
Lt.UtA, PI	LN	12008-9	Pulsatility Index	TID ALOKA-300 row3
Lt.UtA, RI	LN	12023-8	Resistivity Index	TID ALOKA-300 row3
Lt.UtA, S/D	LN	12144-2	Systolic to Diastolic Velocity Ratio	TID ALOKA-300 row3
Lt.UtA, PSV	LN	11726-7	Peak Systolic Velocity	TID ALOKA-300 row3
Lt.UtA, EDV	LN	11653-3	End Diastolic Velocity	TID ALOKA-300 row3
Lt.UtA, MnV	LN	11692-1	Time averaged peak velocity	TID ALOKA-300 row3
D-Ao	SRT	T-D0765	Descending Aorta	TID 5025 row 1
D-Ao, PI	LN	12008-9	Pulsatility Index	TID 5025 row 4
D-Ao, RI	LN	12023-8	Resistivity Index	TID 5025 row 4
D-Ao, S/D	LN	12144-2	Systolic to Diastolic Velocity Ratio	TID 5025 row 4
D-Ao, PSV	LN	11726-7	Peak Systolic Velocity	TID 5025 row 4
D-Ao, EDV	LN	11653-3	End Diastolic Velocity	TID 5025 row 4
D-Ao, MnV	LN	11692-1	Time averaged peak velocity	TID 5025 row 4
OB Dop1	99ALOKA	ALOKA_001-001	OB User Definition Doppler1	TID ALOKA-5000 row3 -TID ALOKA-300 row1
OB Dop1, PI	LN	12008-9	Pulsatility Index	TID ALOKA-300 row3
OB Dop1, RI	LN	12023-8	Resistivity Index	TID ALOKA-300 row3
OB Dop1, S/D	LN	12144-2	Systolic to Diastolic Velocity Ratio	TID ALOKA-300 row3
OB Dop1, PSV	LN	11726-7	Peak Systolic Velocity	TID ALOKA-300 row3
OB Dop1, EDV	LN	11653-3	End Diastolic Velocity	TID ALOKA-300 row3
OB Dop1, MnV	LN	11692-1	Time averaged peak velocity	TID ALOKA-300 row3
OB Dop2	99ALOKA	ALOKA_001-002	OB User Definition Doppler2	TID ALOKA-5000 row3 -TID ALOKA-300 row1
OB Dop2, PI	LN	12008-9	Pulsatility Index	TID ALOKA-300 row3
OB Dop2, RI	LN	12023-8	Resistivity Index	TID ALOKA-300 row3

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
OB Dop2, S/D	LN	12144-2	Systolic to Diastolic Velocity Ratio	TID ALOKA-300 row3
OB Dop2, PSV	LN	11726-7	Peak Systolic Velocity	TID ALOKA-300 row3
OB Dop2, EDV	LN	11653-3	End Diastolic Velocity	TID ALOKA-300 row3
OB Dop2, MnV	LN	11692-1	Time averaged peak velocity	TID ALOKA-300 row3
OB Dop3	99ALOKA	ALOKA_001-003	OB User Definition Doppler3	TID ALOKA-5000 row3 -TID ALOKA-300 row1
OB Dop3, PI	LN	12008-9	Pulsatility Index	TID ALOKA-300 row3
OB Dop3, RI	LN	12023-8	Resistivity Index	TID ALOKA-300 row3
OB Dop3, S/D	LN	12144-2	Systolic to Diastolic Velocity Ratio	TID ALOKA-300 row3
OB Dop3, PSV	LN	11726-7	Peak Systolic Velocity	TID ALOKA-300 row3
OB Dop3, EDV	LN	11653-3	End Diastolic Velocity	TID ALOKA-300 row3
OB Dop3, MnV	LN	11692-1	Time averaged peak velocity	TID ALOKA-300 row3
DV(Ductus Venosus)	99ALOKA	A12141-001	Ductus Venosus	TID 5025 row 1
DV(Ductus Venosus), PI	LN	12008-9	Pulsatility Index	TID 5025 row 4
DV(Ductus Venosus), a/S	99ALOKA	A12121-002	atrial contraction to Systolic peak Velocity Ratio	TID 5025 row 4
DV(Ductus Venosus), S	LN	11726-7	Peak Systolic Velocity	TID 5025 row 4
DV(Ductus Venosus), a	99ALOKA	A12120-001	atrial contraction	TID 5025 row 4
DV(Ductus Venosus), MnV	LN	11692-1	Time averaged peak velocity	TID 5025 row 4
DV(Ductus Venosus), D	99ALOKA	A12120-002	Peak Diastolic Velocity	TID 5025 row 4
DV(Ductus Venosus), S/a	99ALOKA	A12121-003	Systolic peak to atrial contraction Velocity Ratio	TID 5025 row 4
DV(Ductus Venosus), PVIV	99ALOKA	A12121-004	Peak velocity index for veins	TID 5025 row 4
User-defined Equation Result Parameter	99ALOKA	ALOKA_002-001	Calculated result.	TID ALOKA-5002 row2 -TID ALOKA-301 row4 -TID ALOKA-302 row1

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
User-defined Equation Parameter	99ALOKA	ALOKA_002-002	Equation parameter	TID ALOKA-5002 row2 -TID ALOKA-301 row4 -TID ALOKA-302 row1

Note: 1. UmA is included in both TID 5026(not TID 5025) and TID ALOKA-5000.

When UmA is measured on multiple fetuses, the precedence to be included in TID 5026 is Fetus a, b, then c. TID ALOKA-5000 includes UmA measured on each fetus.

TID 1008 Subject Context, Fetus

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
Fetus a/b/c	LN	11951-1	Fetus ID	TID 1008 row 4

Note: Fetus ID is included in Structured Report as "A", "B", "C". When singleton fetus, Fetus ID is always "A".

GYN measurements :

TID 5000 row 15 - TID 5015 Pelvis and Uterus Section

TID 5000 row 16 - TID 5012 Ovaries Section

TID 5000 row 17 - TID 5013 Follicles Section

TID 5000 row 18 - TID 5013 Follicles Section

TID 5000 row 29(Extension) - TID ALOKA-301 User Defined Equation Section

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
Uterus, Length	LN	11842-2	Uterus Length	TID 5015 row 2 - TID 5016 row 3
Uterus, AP	LN	11859-6	Uterus Height	TID 5015 row 2 - TID 5016 row 5
Uterus, Width	LN	11865-3	Uterus Width	TID 5015 row 2 - TID 5016 row 4
Uterus, Volume	LN	33192-6	Uterus Volume	TID 5015 row 2 - TID 5016 row 2
Cervix, Length	LN	11961-0	Cervix Length	TID 5015 row 3
Cervix, AP	99ALOKA	A12011-002	Cervix Antero-Posterior Diameter	TID 5015 row 3

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
Cervix, Width	99ALOKA	A12011-001	Cervix Width	TID 5015 row 3
Endom-T	LN	12145-9	Endometrium Thickness	TID 5015 row 3
Pre Bldrvol, Length	99ALOKA	A12011-003	Pre Void Bladder Length	TID 5015 row 3
Pre Bldrvol, AP	99ALOKA	A12011-004	Pre Void Bladder Antero-Posterior Diameter	TID 5015 row 3
Pre Bldrvol, Width	99ALOKA	A12011-005	Pre Void Bladder Width	TID 5015 row 3
Pre Bldrvol, Volume	99ALOKA	A12011-006	Pre Void Bladder Volume	TID 5015 row 3
Pst Bldrvol, Length	99ALOKA	A12011-007	Post Void Bladder Length	TID 5015 row 3
Pst Bldrvol, AP	99ALOKA	A12011-008	Post Void Bladder Antero-Posterior Diameter	TID 5015 row 3
Pst Bldrvol, Width	99ALOKA	A12011-009	Post Void Bladder Width	TID 5015 row 3
Pst Bldrvol, Volume	99ALOKA	A12011-010	Post Void Bladder Volume	TID 5015 row 3
Void Volume	99ALOKA	A12011-011	Bladder Void Volume	TID 5015 row 3
Left Ovary, Length	LN	11840-6	Left Ovary Length	TID 5012 row 3 - TID 5016 row 3
Left Ovary, AP	LN	11857-0	Left Ovary Height	TID 5012 row 3 - TID 5016 row 5
Left Ovary, Width	LN	11829-9	Left Ovary Width	TID 5012 row 3 - TID 5016 row 4
Left Ovary, Volume	LN	12164-0	Left Ovary Volume	TID 5012 row 3 - TID 5016 row 2
Right Ovary, Length	LN	11841-4	Right Ovary Length	TID 5012 row 4 - TID 5016 row 3
Right Ovary, AP	LN	11858-8	Right Ovary Height	TID 5012 row 4 - TID 5016 row 5
Right Ovary, Width	LN	11830-7	Right Ovary Width	TID 5012 row 4 - TID 5016 row 4
Right Ovary, Volume	LN	12165-7	Right Ovary Volume	TID 5012 row 4 - TID 5016 row 2
Follicles (up to 10 items)	-	-	-	-
Left Follicles, D1	LN	11793-7	Follicle Diameter	TID 5013 row 5 - TID 5014 row 4
Left Follicles, D2	LN	11793-7	Follicle Diameter	TID 5013 row 5 - TID 5014 row 4
Left Follicles, Average	LN	11793-7	Follicle Diameter	TID 5013 row 5 - TID 5014 row 4

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
Left Follicles, Volume	SRT	G-D705	Volume	TID 5013 row 5 - TID 5014 row 3
Right Follicles, D1	LN	11793-7	Follicle Diameter	TID 5013 row 5 - TID 5014 row 4
Right Follicles, D2	LN	11793-7	Follicle Diameter	TID 5013 row 5 - TID 5014 row 4
Right Follicles, Average	LN	11793-7	Follicle Diameter	TID 5013 row 5 - TID 5014 row 4
Right Follicles, Volume	SRT	G-D705	Volume	TID 5013 row 5 - TID 5014 row 3
Follicles Volume (up to 10 items)	-	-	-	-
Left Follicles, D1	LN	11793-7	Follicle Diameter	TID 5013 row 5 - TID 5014 row 4
Left Follicles, D2	LN	11793-7	Follicle Diameter	TID 5013 row 5 - TID 5014 row 4
Left Follicles, D3	LN	11793-7	Follicle Diameter	TID 5013 row 5 - TID 5014 row 4
Left Follicles, Average	LN	11793-7	Follicle Diameter	TID 5013 row 5 - TID 5014 row 4
Left Follicles, Volume	SRT	G-D705	Volume	TID 5013 row 5 - TID 5014 row 3
Right Follicles, D1	LN	11793-7	Follicle Diameter	TID 5013 row 5 - TID 5014 row 4
Right Follicles, D2	LN	11793-7	Follicle Diameter	TID 5013 row 5 - TID 5014 row 4
Right Follicles, D3	LN	11793-7	Follicle Diameter	TID 5013 row 5 - TID 5014 row 4
Right Follicles, Average	LN	11793-7	Follicle Diameter	TID 5013 row 5 - TID 5014 row 4
Right Follicles, Volume	SRT	G-D705	Volume	TID 5013 row 5 - TID 5014 row 3
Rt.Uterine Artery	SRT	T-46820	Uterine Artery	TID ALOKA_5001 row2 - TID ALOKA_300 row1
Rt.Uterine Artery, PI	LN	12008-9	Pulsatility Index	TID ALOKA_300 row3
Rt.Uterine Artery, RI	LN	12023-8	Resistivity Index	TID ALOKA_300 row3

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
Rt.Uterine Artery, S/D	LN	12144-2	Systolic to Diastolic Velocity Ratio	TID ALOKA _300 row3
Rt.Uterine Artery, PSV	LN	11726-7	Peak Systolic Velocity	TID ALOKA _300 row3
Rt.Uterine Artery, EDV	LN	11653-3	End Diastolic Velocity	TID ALOKA _300 row3
Rt.Uterine Artery, MnV	LN	11692-1	Time averaged peak velocity	TID ALOKA _300 row3
User-defined Equation Result Parameter	99ALOKA	ALOKA_002-001	Calculated result.	TID ALOKA-301 row4 -TID ALOKA-302 row1
User-defined Equation Parameter	99ALOKA	ALOKA_002-002	Equation parameter	TID ALOKA-301 row4 -TID ALOKA-302 row1

8.6.2.2 TID 5100 Vascular Ultrasound Report

TID 5100 row 4 - TID 1001 Observation Context

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
Patient Name	DCM	121029	Subject Name	TID 1001 row 3 - TID 1006 row 2 - TID 1007 row 2

TID 5100 row 5 - TID 5101 Vascular Patient Characteristics

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
Age	DCM	121033	Subject Age	TID 5101 row 2
Sex	DCM	121032	Subject Sex	TID 5101 row 3
Blood Pressure (Systolic)	SRT	F-008EC	Systolic Blood Pressure	TID 5101 row 5
Blood Pressure (Diastolic)	SRT	F-008ED	Diastolic Blood Pressure	TID 5101 row 6

TID 5100 row 8 - TID 5102 Vascular Procedure Summary Section

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
Comments	DCM	121106	Comment	TID 5102 row 2

Vascular measurements :

TID 5100 row 9 - TID 5103 Vascular Ultrasound Section (Blood Vessel of Head, Left)
 TID 5100 row 10 - TID 5103 Vascular Ultrasound Section (Blood Vessel of Head, Right)
 TID 5100 row 11 - TID 5103 Vascular Ultrasound Section (Blood Vessel of Head, Unilateral)
 TID 5100 row 12 - TID 5103 Vascular Ultrasound Section (Artery of neck, Left)
 TID 5100 row 13 - TID 5103 Vascular Ultrasound Section (Artery of neck, Right)
 TID 5100 row 14 - TID 5103 Vascular Ultrasound Section (Artery of Lower Extremity, Left)
 TID 5100 row 15 - TID 5103 Vascular Ultrasound Section (Artery of Lower Extremity, Right)
 TID 5100 row 16 - TID 5103 Vascular Ultrasound Section (Vein of Lower Extremity, Left)
 TID 5100 row 17 - TID 5103 Vascular Ultrasound Section (Vein of Lower Extremity, Right)
 TID 5100 row 18 - TID 5103 Vascular Ultrasound Section (Artery Of Upper Extremity, Left)
 TID 5100 row 19 - TID 5103 Vascular Ultrasound Section (Artery Of Upper Extremity, Right)
 TID 5100 row 20 - TID 5103 Vascular Ultrasound Section (Vein Of Upper Extremity, Left)
 TID 5100 row 21 - TID 5103 Vascular Ultrasound Section (Vein Of Upper Extremity, Right)
 TID 5100 row 31(Extension) - TID ALOKA-301 User Defined Equation Section

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
(Modifier) ¹				
Left	SRT	G-A101	Left	TID 5103 row 3
Right	SRT	G-A100	Right	TID 5103 row 3
-	SRT	G-A103	Unilateral	TID 5103 row 3
- (prox)	SRT	G-A118	Proximal	TID 5103 row 4 - TID 5104 row 2
- (mid)	SRT	G-A188	Mid-longitudinal	TID 5103 row 4 - TID 5104 row 2
- (distal)	SRT	G-A119	Distal	TID 5103 row 4 - TID 5104 row 2
(Measurements) ²				
PSV pV	LN	11726-7	Peak Systolic Velocity	TID 5103 row 4 - TID 5104 row 4
EDV	LN	11653-3	End Diastolic Velocity	TID 5103 row 4 - TID 5104 row 4
MnV	LN	11692-1	Time averaged peak velocity	TID 5103 row 4 - TID 5104 row 4

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
PI	LN	12008-9	Pulsatility Index	TID 5103 row 4 - TID 5104 row 4
RI	LN	12023-8	Resistivity Index	TID 5103 row 4 - TID 5104 row 4
S/D	LN	12144-2	Systolic to Diastolic Velocity Ratio	TID 5103 row 4 - TID 5104 row 4
(Vessels)				
(row 9, 10)				
ACA	SRT	T-45540	Anterior Cerebral Artery	TID 5103 row 4 - TID 5104 row 1
MCA	SRT	T-45600	Middle Cerebral Artery	TID 5103 row 4 - TID 5104 row 1
PCA	SRT	T-45900	Posterior Cerebral Artery	TID 5103 row 4 - TID 5104 row 1
PCoA	SRT	T-45320	Posterior Communicating Artery	TID 5103 row 4 - TID 5104 row 1
TICA	SRT	R-102BD	Terminal internal carotid artery	TID 5103 row 4 - TID 5104 row 1
(row 11)				
ACoA ³	SRT	T-45530	Anterior Communicating Artery	TID 5103 row 4 - TID 5104 row 1
BA	SRT	T-45800	Basilar Artery	TID 5103 row 4 - TID 5104 row 1
(row 12, 13)				
BIFUR	SRT	T-45160	Carotid Bifurcation	TID 5103 row 4 - TID 5104 row 1
CCA prox CCA mid CCA distal	SRT	T-45100	Common Carotid Artery	TID 5103 row 4 - TID 5104 row 1
ECA	SRT	T-45200	External Carotid Artery	TID 5103 row 4 - TID 5104 row 1
ICA ICA prox ICA mid ICA distal	SRT	T-45300	Internal Carotid Artery	TID 5103 row 4 - TID 5104 row 1
VA ⁴ VERT	SRT	T-45700	Vertebral Artery	TID 5103 row 4 - TID 5104 row 1
(row 14, 15)				

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
ATA	SRT	T-47700	Anterior Tibial Artery	TID 5103 row 4 - TID 5104 row 1
CFA	SRT	T-47400	Common Femoral Artery	TID 5103 row 4 - TID 5104 row 1
CIA	SRT	T-46710	Common Iliac Artery	TID 5103 row 4 - TID 5104 row 1
DFA (Deep Femoral Artery)	SRT	T-47440	Profunda Femoris Artery	TID 5103 row 4 - TID 5104 row 1
DPA	SRT	T-47741	Dorsalis Pedis Artery	TID 5103 row 4 - TID 5104 row 1
EIA	SRT	T-46910	External Iliac Artery	TID 5103 row 4 - TID 5104 row 1
IIA	SRT	T-46740	Internal Iliac Artery	TID 5103 row 4 - TID 5104 row 1
PerA	SRT	T-47630	Peroneal Artery	TID 5103 row 4 - TID 5104 row 1
PopA	SRT	T-47500	Popliteal Artery	TID 5103 row 4 - TID 5104 row 1
PTA	SRT	T-47600	Posterior Tibial Artery	TID 5103 row 4 - TID 5104 row 1
SFA	SRT	T-47403	Superficial Femoral Artery	TID 5103 row 4 - TID 5104 row 1
Lwr Art.1	99ALOKA	A12109-001	User Definition Artery1	TID 5103 row 4 - TID 5104 row 1
Lwr Art.2	99ALOKA	A12109-002	User Definition Artery2	TID 5103 row 4 - TID 5104 row 1
Lwr Art.3	99ALOKA	A12109-003	User Definition Artery3	TID 5103 row 4 - TID 5104 row 1
Lwr Art.4	99ALOKA	A12109-004	User Definition Artery4	TID 5103 row 4 - TID 5104 row 1
Lwr Art.5	99ALOKA	A12109-005	User Definition Artery5	TID 5103 row 4 - TID 5104 row 1
Lwr Art.6	99ALOKA	A12109-006	User Definition Artery6	TID 5103 row 4 - TID 5104 row 1
Lwr Art.7	99ALOKA	A12109-007	User Definition Artery7	TID 5103 row 4 - TID 5104 row 1
Lwr Art.8	99ALOKA	A12109-008	User Definition Artery8	TID 5103 row 4 - TID 5104 row 1
(row 16, 17)				
ATV	SRT	T-49630	Anterior Tibial Vein	TID 5103 row 4 - TID 5104 row 1

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
CFV	SRT	G-035B	Common Femoral Vein	TID 5103 row 4 - TID 5104 row 1
CIV	SRT	T-48920	Common Iliac Vein	TID 5103 row 4 - TID 5104 row 1
DFV (Deep Femoral Vein)	SRT	T-49660	Profunda Femoris Vein	TID 5103 row 4 - TID 5104 row 1
EIV	SRT	T-48930	External Iliac Vein	TID 5103 row 4 - TID 5104 row 1
GSV	SRT	T-49530	Great Saphenous Vein	TID 5103 row 4 - TID 5104 row 1
IIV	SRT	T-48940	Internal iliac vein	TID 5103 row 4 - TID 5104 row 1
LSV	SRT	T-49550	Lesser Saphenous Vein	TID 5103 row 4 - TID 5104 row 1
PerV	SRT	T-49650	Peroneal Vein	TID 5103 row 4 - TID 5104 row 1
PopV	SRT	T-49640	Popliteal Vein	TID 5103 row 4 - TID 5104 row 1
PTV	SRT	T-49620	Posterior Tibial Vein	TID 5103 row 4 - TID 5104 row 1
SFV	SRT	G-035A	Superficial Femoral Vein	TID 5103 row 4 - TID 5104 row 1
Lwr Vein.1	99ALOKA	A12110-001	User Definition Vein1	TID 5103 row 4 - TID 5104 row 1
Lwr Vein.2	99ALOKA	A12110-002	User Definition Vein2	TID 5103 row 4 - TID 5104 row 1
Lwr Vein.3	99ALOKA	A12110-003	User Definition Vein3	TID 5103 row 4 - TID 5104 row 1
Lwr Vein.4	99ALOKA	A12110-004	User Definition Vein4	TID 5103 row 4 - TID 5104 row 1
Lwr Vein.5	99ALOKA	A12110-005	User Definition Vein5	TID 5103 row 4 - TID 5104 row 1
Lwr Vein.6	99ALOKA	A12110-006	User Definition Vein6	TID 5103 row 4 - TID 5104 row 1
Lwr Vein.7	99ALOKA	A12110-007	User Definition Vein7	TID 5103 row 4 - TID 5104 row 1
Lwr Vein.8	99ALOKA	A12110-008	User Definition Vein8	TID 5103 row 4 - TID 5104 row 1
(row 18, 19)				
AA	SRT	T-47100	Axillary Artery	TID 5103 row 4 - TID 5104 row 1

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
BA	SRT	T-47160	Brachial Artery	TID 5103 row 4 - TID 5104 row 1
BasA	99ALOKA	A12107-002	Basilic Artery	TID 5103 row 4 - TID 5104 row 1
DBA	99ALOKA	A12107-001	Deep Brachial Artery	TID 5103 row 4 - TID 5104 row 1
RA	SRT	T-47300	Radial Artery	TID 5103 row 4 - TID 5104 row 1
ScA	SRT	T-46100	Subclavian Artery	TID 5103 row 4 - TID 5104 row 1
SPA	SRT	T-47240	Superficial Palmar Arch	TID 5103 row 4 - TID 5104 row 1
UA	SRT	T-47200	Ulnar Artery	TID 5103 row 4 - TID 5104 row 1
Upr Art.1	99ALOKA	A12107-003	User Definition Artery1	TID 5103 row 4 - TID 5104 row 1
Upr Art.2	99ALOKA	A12107-004	User Definition Artery2	TID 5103 row 4 - TID 5104 row 1
Upr Art.3	99ALOKA	A12107-005	User Definition Artery3	TID 5103 row 4 - TID 5104 row 1
Upr Art.4	99ALOKA	A12107-006	User Definition Artery4	TID 5103 row 4 - TID 5104 row 1
Upr Art.5	99ALOKA	A12107-007	User Definition Artery5	TID 5103 row 4 - TID 5104 row 1
Upr Art.6	99ALOKA	A12107-008	User Definition Artery6	TID 5103 row 4 - TID 5104 row 1
Upr Art.7	99ALOKA	A12107-009	User Definition Artery7	TID 5103 row 4 - TID 5104 row 1
Upr Art.8	99ALOKA	A12107-010	User Definition Artery8	TID 5103 row 4 - TID 5104 row 1
(row 20, 21)				
AV	SRT	T-49110	Axillary vein	TID 5103 row 4 - TID 5104 row 1
BasV	SRT	T-48052	Basilic vein	TID 5103 row 4 - TID 5104 row 1
BV	SRT	T-49350	Brachial vein	TID 5103 row 4 - TID 5104 row 1
CV	SRT	T-49240	Cephalic vein	TID 5103 row 4 - TID 5104 row 1
DBV	99ALOKA	A12108-001	Deep Brachial vein	TID 5103 row 4 - TID 5104 row 1

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
IJV	SRT	T-48170	Internal Jugular vein	TID 5103 row 4 - TID 5104 row 1
RV	SRT	T-49340	Radial vein	TID 5103 row 4 - TID 5104 row 1
ScV	SRT	T-48330	Subclavian vein	TID 5103 row 4 - TID 5104 row 1
UV	SRT	T-49330	Ulnar vein	TID 5103 row 4 - TID 5104 row 1
Upr Vein.1	99ALOKA	A12108-002	User Definition Vein1	TID 5103 row 4 - TID 5104 row 1
Upr Vein.2	99ALOKA	A12108-003	User Definition Vein2	TID 5103 row 4 - TID 5104 row 1
Upr Vein.3	99ALOKA	A12108-004	User Definition Vein3	TID 5103 row 4 - TID 5104 row 1
Upr Vein.4	99ALOKA	A12108-005	User Definition Vein4	TID 5103 row 4 - TID 5104 row 1
Upr Vein.5	99ALOKA	A12108-006	User Definition Vein5	TID 5103 row 4 - TID 5104 row 1
Upr Vein.6	99ALOKA	A12108-007	User Definition Vein6	TID 5103 row 4 - TID 5104 row 1
Upr Vein.7	99ALOKA	A12108-008	User Definition Vein7	TID 5103 row 4 - TID 5104 row 1
Upr Vein.8	99ALOKA	A12108-009	User Definition Vein8	TID 5103 row 4 - TID 5104 row 1
User-defined Equation Result Parameter	99ALOKA	ALOKA_002-001	Calculated result.	TID ALOKA-301 row4 -TID ALOKA-302 row1
User-defined Equation Parameter	99ALOKA	ALOKA_002-002	Equation parameter	TID ALOKA-301 row4 -TID ALOKA-302 row1

- Note:
1. Prox / mid / distal is included in Structured Report only for the measurements of CCA prox, CCA mid, CCA distal, ICA prox, ICA mid and ICA distal. Otherwise not included.
 2. PSV, EDV, MnV, PI, RI and S/D are included in Structured Report for the measurements of artery, pV for the measurements of vein.
 3. Although DICOM Standard defines ACoA (Anterior Communicating Artery) in CID 12105, it has no laterality and is included in row 11 "Unilateral" section of TID 5100.
 4. When both VA and VERT are measured, only VERT is included in Structured Report, not VA.

Abdomen measurements :

TID 5100 row 22 - TID 5103 Vascular Ultrasound Section (Vascular Structure Of Kidney, Left)
TID 5100 row 23 - TID 5103 Vascular Ultrasound Section (Vascular Structure Of Kidney, Right)
TID 5100 row 24 - TID 5103 Vascular Ultrasound Section (Artery of Abdomen, Left)
TID 5100 row 25 - TID 5103 Vascular Ultrasound Section (Artery of Abdomen, Right)

TID 5100 row 26 - TID 5103 Vascular Ultrasound Section (Artery of Abdomen, Unilateral)
TID 5100 row 29 - TID 5103 Vascular Ultrasound Section (Vein of Abdomen, Unilateral)
TID 5100 row 31(Extension) - TID ALOKA-301 User Defined Equation Section
TID 5100 row 32(Extension) - TID ALOKA-5100 Vascular Original Measurements Section

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
(Modifier) ¹				
Left	SRT	G-A101	Left	TID 5103 row 3
Right	SRT	G-A100	Right	TID 5103 row 3
-	SRT	G-A103	Unilateral	TID 5103 row 3
- (prox)	SRT	G-A118	Proximal	TID 5103 row 4 - TID 5104 row 2
- (mid)	SRT	G-A188	Mid-longitudinal	TID 5103 row 4 - TID 5104 row 2
- (distal)	SRT	G-A119	Distal	TID 5103 row 4 - TID 5104 row 2
Pre Prandial	99ALOKA	A-001	Pre-prandial	TID 5103 row 4 - TID 5104 row 6
Post Prandial	SRT	G-A491	Post-prandial	TID 5103 row 4 - TID 5104 row 6
Upper	SRT	G-A116	Superior	TID 5104 Row 3
Mid	SRT	G-A109	Medial	TID 5104 Row 3
Lower	SRT	G-A115	Inferior	TID 5104 Row 3
(Measurements) ²				
PSV pV	LN	11726-7	Peak Systolic Velocity	TID 5103 row 4 - TID 5104 row 4
EDV	LN	11653-3	End Diastolic Velocity	TID 5103 row 4 - TID 5104 row 4
MnV	LN	11692-1	Time averaged peak velocity	TID 5103 row 4 - TID 5104 row 4
PI	LN	12008-9	Pulsatility Index	TID 5103 row 4 - TID 5104 row 4
RI	LN	12023-8	Resistivity Index	TID 5103 row 4 - TID 5104 row 4
S/D	LN	12144-2	Systolic to Diastolic Velocity Ratio	TID 5103 row 4 - TID 5104 row 4
AccT	LN	20168-1	Acceleration Time	TID 5103 row 4 - TID 5104 row 4
ACC	LN	20167-3	Acceleration Index	TID 5103 row 4 - TID 5104 row 4
FlowT	99ALOKA	A12122-001	Flow Time	TID 5103 row 4 - TID 5104 row 4

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
AccT/FlowT	99ALOKA	A12121-001	Acceleration Time to Flow Time Ratio	TID 5103 row 4 - TID 5104 row 4
RAR	LN	33869-9	Renal Artery/Aorta velocity ratio	TID 5103 row 5
(Vessels)				
(row 22, 23)				
Renal-A	SRT	T-46600	Renal Artery	TID 5103 row 4 - TID 5104 row 1
Renal-A hilum	SRT	G-035C	Hilar Artery	TID 5103 row 4 - TID 5104 row 1
Seg-A	SRT	T-46659	Segmental Artery	TID 5103 row 4 - TID 5104 row 1
(row 24, 25)				
CIA ³	SRT	T-46710	Common Iliac Artery	TID 5103 row 4 - TID 5104 row 1
(row 26)				
A-Ao	SRT	T-42000	Aorta	TID 5103 row 4 - TID 5104 row 1
CA	SRT	T-46400	Celiac Axis	TID 5103 row 4 - TID 5104 row 1
CHA	SRT	T-46421	Common Hepatic Artery	TID 5103 row 4 - TID 5104 row 1
SA	SRT	T-46460	Splenic Artery	TID 5103 row 4 - TID 5104 row 1
SMA Prandial SMA	SRT	T-46510	Superior Mesenteric Artery	TID 5103 row 4 - TID 5104 row 1
IMA	SRT	T-46520	Inferior Mesenteric Artery	TID 5103 row 4 - TID 5104 row 1
HA, Left	SRT	T-46427	Left Branch of Hepatic Artery	TID 5103 row 4 - TID 5104 row 1
HA, Right	SRT	T-46423	Right Branch of Hepatic Artery	TID 5103 row 4 - TID 5104 row 1
Artery1	99ALOKA	A12112-001	User Definition Artery1	TID 5103 row 4 - TID 5104 row 1
Artery2	99ALOKA	A12112-002	User Definition Artery2	TID 5103 row 4 - TID 5104 row 1
Artery3	99ALOKA	A12112-003	User Definition Artery3	TID 5103 row 4 - TID 5104 row 1
(row 29)				

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
Main PV	SRT	T-48810	Portal Vein	TID 5103 row 4 - TID 5104 row 1
Lt.PV	SRT	T-4881F	Left Main Branch of Portal Vein	TID 5103 row 4 - TID 5104 row 1
Rt.PV	SRT	T-4882A	Right Main Branch of Portal Vein	TID 5103 row 4 - TID 5104 row 1
Prox Shunt Mid Shunt Distal Shunt	SRT	G-036C	Transjugular Intrahepatic Portosystemic Shunt	TID 5103 row 4 - TID 5104 row 1
(row 31)				
User-defined Equation Result Parameter	99ALOKA	ALOKA_002-001	Calculated result.	TID ALOKA-301 row4 -TID ALOKA-302 row1
User-defined Equation Parameter	99ALOKA	ALOKA_002-002	Equation parameter	TID ALOKA-301 row4 -TID ALOKA-302 row1
(row 32)				
%Stenosis Diam.	99ALOKA	ALOKA_005-001	Vessel Diameter Stenosis Ratio	TID ALOKA-5100 row 2 - TID ALOKA-300 row 3
%Stenosis Area	99ALOKA	ALOKA_005-002	Vessel Area Stenosis Ratio	TID ALOKA-5100 row 2 - TID ALOKA-300 row 3
Gallbladder	SRT	T-63000	Gall Bladder	TID ALOKA-5100 row 2 - TID ALOKA-300 row 3
GB Wall-T	99ALOKA	ALOKA_005-003	Gall Bladder Wall	TID ALOKA-5100 row 2 - TID ALOKA-300 row 3
Right Liver Left Liver	SRT	T-62000	Liver	TID ALOKA-5100 row 2 - TID ALOKA-300 row 3
CDB	99ALOKA	ALOKA_005-004	Common Bile Duct	TID ALOKA-5100 row 2 - TID ALOKA-300 row 3

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
Pancreas	SRT	T-65000	Pancreas	TID ALOKA-5100 row 2 - TID ALOKA-300 row 3
P-Duct	SRT	T-65010	Pancreatic duct	TID ALOKA-5100 row 2 - TID ALOKA-300 row 3
Right Renal Volume Left Renal Volume	99ALOKA	ALOKA_005-005	Renal Volume	TID ALOKA-5100 row 2 - TID ALOKA-300 row 3
Spleen	SRT	T-C3000	Spleen	TID ALOKA-5100 row 2 - TID ALOKA-300 row 3
SOL1	99ALOKA	ALOKA_005-006	Space Occupying Lesion1	TID ALOKA-5100 row 2 - TID ALOKA-300 row 3
SOL2	99ALOKA	ALOKA_005-007	Space Occupying Lesion2	TID ALOKA-5100 row 2 - TID ALOKA-300 row 3
Aorta Diam	99ALOKA	ALOKA_005-008	Aorta Diameter	TID ALOKA-5100 row 2 - TID ALOKA-300 row 3
PV Diam	99ALOKA	ALOKA_005-009	Portal Vein Diameter	TID ALOKA-5100 row 2 - TID ALOKA-300 row 3

- Note:
1. Prox / mid / distal is included in Structured Report only for the measurements of Prox Shunt, Mid Shunt and Distal Shunt. Otherwise not included.
Pre Prandial / Post Prandial is included only for the measurement of Prandial SMA. Otherwise not included.
 2. PSV, EDV, MnV, PI, RI, S/D, AccT, ACC, FlowT and AccT/FlowT are included in Structured Report for the measurements of Renal-A, Renal-A hilum and Seg-A (TID 5100 row 22, 23).
RAR is included for the measurements of Renal-A (TID 5100 row 22, 23).
PSV, EDV, MnV, PI, RI, S/D, AccT and ACC are included for the measurements of Artery of Abdomen (TID 5100 row 24, 25, 26).
pV is included for the measurements of Vein of Abdomen (TID 5100 row 29).
 3. Although DICOM Standard defines CIA (Common Iliac Artery) in CID 12112 (unilateral), it has laterality and is included in row 24 "Left" and row 25 "Right" sections of TID 5100.

max-IMT measurements :

TID 5100 row 12 - TID 5103 Vascular Ultrasound Section (Artery of neck, Left)

TID 5100 row 13 - TID 5103 Vascular Ultrasound Section (Artery of neck, Right)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
(Modifier)				
Left	SRT	G-A101	Left	TID 5103 row 3
Right	SRT	G-A100	Right	TID 5103 row 3
(Measurements)				
1 max-IMT	99ALOKA	A12122-002	Intima-media thickness 1	TID 5103 row 4 - TID 5104 row 4
2 max-IMT	99ALOKA	A12122-003	Intima-media thickness 2	TID 5103 row 4 - TID 5104 row 4
3 max-IMT	99ALOKA	A12122-004	Intima-media thickness 3	TID 5103 row 4 - TID 5104 row 4
4 max-IMT	99ALOKA	A12122-005	Intima-media thickness 4	TID 5103 row 4 - TID 5104 row 4
5 max-IMT	99ALOKA	A12122-006	Intima-media thickness 5	TID 5103 row 4 - TID 5104 row 4
6 max-IMT	99ALOKA	A12122-007	Intima-media thickness 6	TID 5103 row 4 - TID 5104 row 4
7 max-IMT	99ALOKA	A12122-008	Intima-media thickness 7	TID 5103 row 4 - TID 5104 row 4
8 max-IMT	99ALOKA	A12122-009	Intima-media thickness 8	TID 5103 row 4 - TID 5104 row 4
9 max-IMT	99ALOKA	A12122-010	Intima-media thickness 9	TID 5103 row 4 - TID 5104 row 4
10 max-IMT	99ALOKA	A12122-011	Intima-media thickness 10	TID 5103 row 4 - TID 5104 row 4
max-IMT	99ALOKA	A12122-012	Maximum Intima-media thickness	TID 5103 row 4 - TID 5104 row 4
(Vessels)				

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
(row 12, 13)				
max-IMT	99ALOKA	A12104-001	Maximum Intima-media thickness	TID 5103 row 4 - TID 5104 row 1
Cmax-IMT	SRT	T-45100	Common Carotid Artery	TID 5103 row 4 - TID 5104 row 1
Bmax-IMT	SRT	T-45160	Carotid Bifurcation	TID 5103 row 4 - TID 5104 row 1
Imax-IMT	SRT	T-45300	Internal Carotid Artery	TID 5103 row 4 - TID 5104 row 1

mean-IMT measurements :

TID 5100 row 12 - TID 5103 Vascular Ultrasound Section (Artery of neck, Left)

TID 5100 row 13 - TID 5103 Vascular Ultrasound Section (Artery of neck, Right)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
(Modifier)				
Left	SRT	G-A101	Left	TID 5103 row 3
Right	SRT	G-A100	Right	TID 5103 row 3
(Measurements)				
1 IMT	99ALOKA	A12122-002	Intima-media thickness 1	TID 5103 row 4 - TID 5104 row 4
2 IMT	99ALOKA	A12122-003	Intima-media thickness 2	TID 5103 row 4 - TID 5104 row 4
3 IMT	99ALOKA	A12122-004	Intima-media thickness 3	TID 5103 row 4 - TID 5104 row 4
mean-IMT	99ALOKA	A12122-013	Mean Intima-media thickness	TID 5103 row 4 - TID 5104 row 4
(Vessels)				

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
(row 12, 13)				
mean-IMT	99ALOKA	A12104-002	Mean Intima-media thickness	TID 5103 row 4 - TID 5104 row 1

IMT measurements :

TID 5100 row 12 - TID 5103 Vascular Ultrasound Section (Artery of neck, Left)

TID 5100 row 13 - TID 5103 Vascular Ultrasound Section (Artery of neck, Right)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
(Modifier)				
Left	SRT	G-A101	Left	TID 5103 row 3
Right	SRT	G-A100	Right	TID 5103 row 3
- (Farwall)	99ALOKA	A12117-001	far wall	TID 5103 row 4 - TID 5104 row 3
- (Nearwall)	99ALOKA	A12117-002	near wall	TID 5103 row 4 - TID 5104 row 3
(Measurements)				
D-max	99ALOKA	A12122-014	Maximum Intima-media thickness at diastole	TID 5103 row 4 - TID 5104 row 4
D-min	99ALOKA	A12122-015	Minimum Intima-media thickness at diastole	TID 5103 row 4 - TID 5104 row 4
D-average	99ALOKA	A12122-016	Average Intima-media thickness at diastole	TID 5103 row 4 - TID 5104 row 4
D-SD	99ALOKA	A12122-018	Standard Deviation of Intima-media thickness at diastole	TID 5103 row 4 - TID 5104 row 4
D-points	99ALOKA	A12122-019	Measured number of Intima-media thickness points at diastole	TID 5103 row 4 - TID 5104 row 4

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
D-width	99ALOKA	A12122-020	Width of ROI at diastole	TID 5103 row 4 - TID 5104 row 4
S-max	99ALOKA	A12122-021	Maximum Intima-media thickness at systole	TID 5103 row 4 - TID 5104 row 4
S-min	99ALOKA	A12122-022	Minimum Intima-media thickness at systole	TID 5103 row 4 - TID 5104 row 4
S-average	99ALOKA	A12122-023	Average Intima-media thickness at systole	TID 5103 row 4 - TID 5104 row 4
S-SD	99ALOKA	A12122-025	Standard Deviation of Intima-media thickness at systole	TID 5103 row 4 - TID 5104 row 4
S-points	99ALOKA	A12122-026	Measured number of Intima-media thickness points at systole	TID 5103 row 4 - TID 5104 row 4
S-width	99ALOKA	A12122-027	Width of ROI at systole	TID 5103 row 4 - TID 5104 row 4
(Vessels)				
(row 12, 13)				
CCA-IMT (Auto)	99ALOKA	A12104-003	Intima-media thickness of CCA by auto-detection	TID 5103 row 4 - TID 5104 row 1
ICA-IMT (Auto)	99ALOKA	A12104-004	Intima-media thickness of ICA by auto-detection	TID 5103 row 4 - TID 5104 row 1
ECA-IMT (Auto)	99ALOKA	A12104-005	Intima-media thickness of ECA by auto-detection	TID 5103 row 4 - TID 5104 row 1
BIF-IMT (Auto)	99ALOKA	A12104-006	Intima-media thickness of BIFUR by auto-detection	TID 5103 row 4 - TID 5104 row 1

Cmean-IMT measurements :

TID 5100 row 12 - TID 5103 Vascular Ultrasound Section (Artery of neck, Left)

TID 5100 row 13 - TID 5103 Vascular Ultrasound Section (Artery of neck, Right)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
(Modifier)				
Left	SRT	G-A101	Left	TID 5103 row 3
Right	SRT	G-A100	Right	TID 5103 row 3
- (Farwall)	99ALOKA	A12117-001	far wall	TID 5103 row 4 - TID 5104 row 3
- (Nearwall)	99ALOKA	A12117-002	near wall	TID 5103 row 4 - TID 5104 row 3
(Measurements)				
D-mean-IMT	99ALOKA	A12122-031	Mean Intima-media thickness at diastole	TID 5103 row 4 - TID 5104 row 4
D-1 IMT	99ALOKA	A12122-028	Intima-media thickness 1 at diastole	TID 5103 row 4 - TID 5104 row 4
D-2 IMT	99ALOKA	A12122-029	Intima-media thickness 2 at diastole	TID 5103 row 4 - TID 5104 row 4
D-3 IMT	99ALOKA	A12122-030	Intima-media thickness 3 at diastole	TID 5103 row 4 - TID 5104 row 4
D-points	99ALOKA	A12122-019	Measured number of Intima-media thickness points at diastole	TID 5103 row 4 - TID 5104 row 4
D-width	99ALOKA	A12122-020	Width of ROI at diastole	TID 5103 row 4 - TID 5104 row 4
S-mean-IMT	99ALOKA	A12122-035	Mean Intima-media thickness at systole	TID 5103 row 4 - TID 5104 row 4
S-1 IMT	99ALOKA	A12122-032	Intima-media thickness 1 at systole	TID 5103 row 4 - TID 5104 row 4
S-2 IMT	99ALOKA	A12122-033	Intima-media thickness 2 at systole	TID 5103 row 4 - TID 5104 row 4
S-3 IMT	99ALOKA	A12122-034	Intima-media thickness 3 at systole	TID 5103 row 4 - TID 5104 row 4
S-points	99ALOKA	A12122-026	Measured number of Intima-media thickness points at systole	TID 5103 row 4 - TID 5104 row 4

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
S-width	99ALOKA	A12122-027	Width of ROI at systole	TID 5103 row 4 - TID 5104 row 4
(Vessels)				
(row 12, 13)				
CmeanIMT (Auto)	99ALOKA	A12104-007	Mean Carotid Intima-media thickness by auto-detection	TID 5103 row 4 - TID 5104 row 1

IMT-C10 measurements :

TID 5100 row 12 - TID 5103 Vascular Ultrasound Section (Artery of neck, Left)

TID 5100 row 13 - TID 5103 Vascular Ultrasound Section (Artery of neck, Right)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
(Modifier)				
Left	SRT	G-A101	Left	TID 5103 row 3
Right	SRT	G-A100	Right	TID 5103 row 3
(Measurements)				
C10	99ALOKA	A12122-036	IMT-C10	TID 5103 row 4 - TID 5104 row 4
(Vessels)				
(row 12, 13)				
IMT-C10	99ALOKA	A12104-008	IMT-C10	TID 5103 row 4 - TID 5104 row 1

Auto IMT-C10 measurements :

TID 5100 row 12 - TID 5103 Vascular Ultrasound Section (Artery of neck, Left)

TID 5100 row 13 - TID 5103 Vascular Ultrasound Section (Artery of neck, Right)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
(Modifier)				
Left	SRT	G-A101	Left	TID 5103 row 3
Right	SRT	G-A100	Right	TID 5103 row 3
(Measurements)				
D-C10	99ALOKA	A12122-017	IMT-C10 at diastole	TID 5103 row 4 - TID 5104 row 4
D-max	99ALOKA	A12122-014	Maximum Intima-media thickness at diastole	TID 5103 row 4 - TID 5104 row 4
D-min	99ALOKA	A12122-015	Minimum Intima-media thickness at diastole	TID 5103 row 4 - TID 5104 row 4
D-avg	99ALOKA	A12122-016	Average Intima-media thickness at diastole	TID 5103 row 4 - TID 5104 row 4
D-SD	99ALOKA	A12122-018	Standard Deviation of Intima-media thickness at diastole	TID 5103 row 4 - TID 5104 row 4
D-points	99ALOKA	A12122-019	Measured number of Intima-media thickness points at diastole	TID 5103 row 4 - TID 5104 row 4
D-width	99ALOKA	A12122-020	Width of ROI at diastole	TID 5103 row 4 - TID 5104 row 4
S-C10	99ALOKA	A12122-024	IMT-C10 at systole	TID 5103 row 4 - TID 5104 row 4
S-max	99ALOKA	A12122-021	Maximum Intima-media thickness at systole	TID 5103 row 4 - TID 5104 row 4
S-min	99ALOKA	A12122-022	Minimum Intima-media thickness at systole	TID 5103 row 4 - TID 5104 row 4
S-avg	99ALOKA	A12122-023	Average Intima-media thickness at systole	TID 5103 row 4 - TID 5104 row 4

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
S-SD	99ALOKA	A12122-025	Standard Deviation of Intima-media thickness at systole	TID 5103 row 4 - TID 5104 row 4
S-points	99ALOKA	A12122-026	Measured number of Intima-media thickness points at systole	TID 5103 row 4 - TID 5104 row 4
S-width	99ALOKA	A12122-027	Width of ROI at systole	TID 5103 row 4 - TID 5104 row 4
(Vessels)				
(row 12, 13)				
Auto IMT-C10	99ALOKA	A12104-009	IMT-C10 by auto-detection	TID 5103 row 4 - TID 5104 row 1

8.6.2.3 TID 5200 Echocardiography Procedure Report (Non-Advanced Cardiac)

TID 5200 row 3 - TID 1001 Observation Context

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
Patient Name	DCM	121029	Subject Name	TID 1001 row 3 - TID 1006 row 2 - TID 1007 row 2

TID 5200 row 4 - TID 5201 Echocardiography Patient Characteristics

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
Age	DCM	121033	Subject Age	TID 5201 row 2
Sex	DCM	121032	Subject Sex	TID 5201 row 3
Blood Pressure (Systolic)	SRT	F-008EC	Systolic Blood Pressure	TID 5201 row 5
Blood Pressure (Diastolic)	SRT	F-008ED	Diastolic Blood Pressure	TID 5201 row 6

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Mapping
BSA	LN	8277-6	Body Surface Area	TID 5201 row 7
BSA Equation, DuBois	DCM	122241	$BSA = 0.007184 * WT^{0.425} * HT^{0.725}$	TID 5201 row 8
BSA Equation, Boyd	99ALOKA	A3663-001	$BSA = 0.0003207 * WT^{(0.7285 - 0.0188 \log(WT))} * HT^{0.3}$	TID 5201 row 8
BSA Equation, Shintani	99ALOKA	A3663-002	$BSA = 0.007358 * HT^{0.725} * WT^{0.425}$	TID 5201 row 8

Cardio measurements :

TID 5200 row 9 - TID 5202 Echo Section (Echocardiography Left Ventricle)
 TID 5200 row 10 - TID 5202 Echo Section (Echocardiography Right Ventricle)
 TID 5200 row 11 - TID 5202 Echo Section (Echocardiography Left Atrium)
 TID 5200 row 12 - TID 5202 Echo Section (Echocardiography Right Atrium)
 TID 5200 row 13 - TID 5202 Echo Section (Echocardiography Aortic Valve)
 TID 5200 row 14 - TID 5202 Echo Section (Echocardiography Mitral Valve)
 TID 5200 row 15 - TID 5202 Echo Section (Echocardiography Pulmonic Valve)
 TID 5200 row 16 - TID 5202 Echo Section (Echocardiography Tricuspid Valve)
 TID 5200 row 17 - TID 5202 Echo Section (Echocardiography Aorta)
 TID 5200 row 18 - TID 5202 Echo Section (Echocardiography Pulmonary Artery)
 TID 5200 row 19 - TID 5202 Echo Section (Echocardiography Vena Cavae)
 TID 5200 row 20 - TID 5202 Echo Section (Echocardiography Pulmonary Veins)
 TID 5200 row 22 - TID 5202 Echo Section (Echocardiography Cardiac Shunt)
 TID 5200 row 25(Extension) - TID ALOKA-301 User Defined Equation Section

The rightmost column "row" denotes the row number of TID 5200. In the "Modifier" column, "Image Mode(Group)" means TID 5202 row 4, "Image Mode" means TID 5203 row 5.

Label	Concept Name	Modifier	row
(Area-Length, BP-Ellipse, Bullet, Simpson)			
LVLd	(18077-8, LN, "Left Ventricle diastolic major axis")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode")	9
LVLs	(18076-0, LN, "Left Ventricle systolic major axis")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode")	9

Label	Concept Name	Modifier	row
(Simpson(Disc))			
LVL2d	(18077-8, LN, "Left Ventricle diastolic major axis")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19B, SRT, "Apical two chamber")	9
LVL2s	(18076-0, LN, "Left Ventricle systolic major axis")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19B, SRT, "Apical two chamber")	9
LVL4d	(18077-8, LN, "Left Ventricle diastolic major axis")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19C, SRT, "Apical four chamber")	9
LVL4s	(18076-0, LN, "Left Ventricle systolic major axis")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19C, SRT, "Apical four chamber")	9
(LV Volume (Pombo, Teichholz, Gibson))			
LVIDd	(29436-3, LN, "Left Ventricle Internal End Diastolic Dimension")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode")	9
LVIDs	(29438-9, LN, "Left Ventricle Internal Systolic Dimension")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode")	9
(LV Function (Pombo, Teichholz, Gibson))			
LVIDd ²	(29436-3, LN, "Left Ventricle Internal End Diastolic Dimension")	Image Mode(Group)=(G-0394, SRT, "M mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-0394, SRT, "M mode")	9
LVIDs ²	(29438-9, LN, "Left Ventricle Internal Systolic Dimension")	Image Mode(Group)=(G-0394, SRT, "M mode"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-0394, SRT, "M mode")	9
(BP-Ellipse)			

Label	Concept Name	Modifier	row
LVSLMVd	(A12201-001, 99ALOKA, "Left Ventricular Short Axis Length at Mitral Valve")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-039A, SRT, "Parasternal short axis at the Mitral Valve level")	9
LVSLMVd	(A12201-001, 99ALOKA, "Left Ventricular Short Axis Length at Mitral Valve")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-039A, SRT, "Parasternal short axis at the Mitral Valve level")	9
(Area-Length, BP-Ellipse)			
LVLAd	(G-0375, SRT, "Left Ventricular Diastolic Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125220, DCM, "Planimetry"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode")	9
LVLAs	(G-0374, SRT, "Left Ventricular Systolic Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125220, DCM, "Planimetry"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode")	9
(Simpson(Disc))			
LVLA2d	(G-0375, SRT, "Left Ventricular Diastolic Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125220, DCM, "Planimetry"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19B, SRT, "Apical two chamber")	9
LVLA2s	(G-0374, SRT, "Left Ventricular Systolic Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125220, DCM, "Planimetry"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19B, SRT, "Apical two chamber")	9
LVLA4d	(G-0375, SRT, "Left Ventricular Diastolic Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125220, DCM, "Planimetry"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19C, SRT, "Apical four chamber")	9

Label	Concept Name	Modifier	row
LVLA4s	(G-0374, SRT, "Left Ventricular Systolic Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125220, DCM, "Planimetry"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19C, SRT, "Apical four chamber")	9
(BP-Ellipse, Simpson)			
LVSAMVd	(G-0375, SRT, "Left Ventricular Diastolic Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125220, DCM, "Planimetry"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-039A, SRT, "Parasternal short axis at the Mitral Valve level")	9
LVSAMVs	(G-0374, SRT, "Left Ventricular Systolic Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125220, DCM, "Planimetry"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-039A, SRT, "Parasternal short axis at the Mitral Valve level")	9
(Bullet, Simpson)			
LVSAPMd	(G-0375, SRT, "Left Ventricular Diastolic Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125220, DCM, "Planimetry"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-039B, SRT, "Parasternal short axis at the Papillary Muscle level")	9
LVSAPMs	(G-0374, SRT, "Left Ventricular Systolic Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125220, DCM, "Planimetry"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-039B, SRT, "Parasternal short axis at the Papillary Muscle level")	9
(LV Mass (AL))			
Aepi	(G-0379, SRT, "Left Ventricle Epicardial Diastolic Area, psax pap view")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125220, DCM, "Planimetry"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-039B, SRT, "Parasternal short axis at the Papillary Muscle level")	9

Label	Concept Name	Modifier	row
Aend	(G-0375, SRT, "Left Ventricular Diastolic Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125220, DCM, "Planimetry"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-039B, SRT, "Parasternal short axis at the Papillary Muscle level")	9
thick	(A12201-002, 99ALOKA, "Mean Wall Thickness")	Image Mode(Group)=(G-03A2, SRT, "2D mode")	9
LVM(AL)	(18087-7, LN, "Left Ventricle Mass")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12232-001, 99ALOKA, "Left Ventricle Mass by Area Length")	9
LVM(AL)/BSA	(A12203-001, 99ALOKA, "Left Ventricular Mass Index")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12232-001, 99ALOKA, "Left Ventricle Mass by Area Length")	9
(B - LV Mass (Devereux))			
LVM(Devereux) ³	(18087-7, LN, "Left Ventricle Mass")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12232-002, 99ALOKA, "Left Ventricle Mass by Penn")	9
LVM(Devereux)/B SA ³	(A12203-001, 99ALOKA, "Left Ventricular Mass Index")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12232-002, 99ALOKA, "Left Ventricle Mass by Penn")	9
(M - LV Mass (Devereux))			
LVM(Devereux) ³	(18087-7, LN, "Left Ventricle Mass")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(A12232-002, 99ALOKA, "Left Ventricle Mass by Penn")	9
LVM(Devereux)/B SA ³	(A12203-001, 99ALOKA, "Left Ventricular Mass Index")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(A12232-002, 99ALOKA, "Left Ventricle Mass by Penn")	9
(B - LV Mass (ASE))			
LVM(ASE)	(18087-7, LN, "Left Ventricle Mass")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125221, DCM, "Left Ventricle Mass by M-mode")	9
LVM(ASE)/BSA ³	(A12203-001, 99ALOKA, "Left Ventricular Mass Index")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125221, DCM, "Left Ventricle Mass by M-mode")	9
(M - LV Mass (ASE))			

Label	Concept Name	Modifier	row
LVM(ASE)	(18087-7, LN, "Left Ventricle Mass")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(125221, DCM, "Left Ventricle Mass by M-mode"), Image Mode=(G-0394, SRT, "M mode")	9
LVM(ASE)/BSA ³	(A12203-001, 99ALOKA, "Left Ventricular Mass Index")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(125221, DCM, "Left Ventricle Mass by M-mode"), Image Mode=(G-0394, SRT, "M mode")	9
(Simpson(Disc))			
EDV	(18026-5, LN, "Left Ventricular End Diastolic Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125207, DCM, "Method of Disks, Biplane"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	9
ESV	(18148-7, LN, "Left Ventricular End Systolic Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125207, DCM, "Method of Disks, Biplane"), Cardiac Cycle Point=(109070, DCM, "End Systole")	9
SV	(F-32120, SRT, "Stroke Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125207, DCM, "Method of Disks, Biplane")	9
SVI	(F-00078, SRT, "Stroke Index")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125207, DCM, "Method of Disks, Biplane")	9
CO	(F-32100, SRT, "Cardiac Output")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125207, DCM, "Method of Disks, Biplane")	9
COI	(F-32110, SRT, "Cardiac Index")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125207, DCM, "Method of Disks, Biplane")	9
EF	(18043-0, LN, "Left Ventricular Ejection Fraction")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125207, DCM, "Method of Disks, Biplane")	9
GLS	(A12203-013, 99ALOKA, "Global Longitudinal Strain")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125207, DCM, "Method of Disks, Biplane")	9
EDV(ap4C)	(18026-5, LN, "Left Ventricular End Diastolic Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125208, DCM, "Method of Disks, Single Plane"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image View=(G-A19C, SRT, "Apical four chamber")	9
ESV(ap4C)	(18148-7, LN, "Left Ventricular End Systolic Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125208, DCM, "Method of Disks, Single Plane"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image View=(G-A19C, SRT, "Apical four chamber")	9

Label	Concept Name	Modifier	row
EDV(ap2C)	(18026-5, LN, "Left Ventricular End Diastolic Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-001, 99ALOKA, "Method of Disks, Single Plane with Apical two chamber"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image View=(G-A19B, SRT, "Apical two chamber")	9
ESV(ap2C)	(18148-7, LN, "Left Ventricular End Systolic Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-001, 99ALOKA, "Method of Disks, Single Plane with Apical two chamber"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image View=(G-A19B, SRT, "Apical two chamber")	9
SV(ap4C)	(F-32120, SRT, "Stroke Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125208, DCM, "Method of Disks, Single Plane"), Image View=(G-A19C, SRT, "Apical four chamber")	9
SVI(ap4C)	(F-00078, SRT, "Stroke Index")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125208, DCM, "Method of Disks, Single Plane"), Image View=(G-A19C, SRT, "Apical four chamber")	9
CO(ap4C)	(F-32100, SRT, "Cardiac Output")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125208, DCM, "Method of Disks, Single Plane"), Image View=(G-A19C, SRT, "Apical four chamber")	9
COI(ap4C)	(F-32110, SRT, "Cardiac Index")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125208, DCM, "Method of Disks, Single Plane"), Image View=(G-A19C, SRT, "Apical four chamber")	9
EF(ap4C)	(18043-0, LN, "Left Ventricular Ejection Fraction")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125208, DCM, "Method of Disks, Single Plane"), Image View=(G-A19C, SRT, "Apical four chamber")	9
areaEF4	(G-0376, SRT, "Left Ventricular Fractional Area Change")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125208, DCM, "Method of Disks, Single Plane"), Image View=(G-A19C, SRT, "Apical four chamber")	9
SV(ap2C)	(F-32120, SRT, "Stroke Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-001, 99ALOKA, "Method of Disks, Single Plane with Apical two chamber"), Image View=(G-A19B, SRT, "Apical two chamber")	9

Label	Concept Name	Modifier	row
SVI(ap2C)	(F-00078, SRT, "Stroke Index")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-001, 99ALOKA, "Method of Disks, Single Plane with Apical two chamber"), Image View=(G-A19B, SRT, "Apical two chamber")	9
CO(ap2C)	(F-32100, SRT, "Cardiac Output")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-001, 99ALOKA, "Method of Disks, Single Plane with Apical two chamber"), Image View=(G-A19B, SRT, "Apical two chamber")	9
COI(ap2C)	(F-32110, SRT, "Cardiac Index")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-001, 99ALOKA, "Method of Disks, Single Plane with Apical two chamber"), Image View=(G-A19B, SRT, "Apical two chamber")	9
EF(ap2C)	(18043-0, LN, "Left Ventricular Ejection Fraction")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-001, 99ALOKA, "Method of Disks, Single Plane with Apical two chamber"), Image View=(G-A19B, SRT, "Apical two chamber")	9
areaEF2	(G-0376, SRT, "Left Ventricular Fractional Area Change")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-001, 99ALOKA, "Method of Disks, Single Plane with Apical two chamber"), Image View=(G-A19B, SRT, "Apical two chamber")	9
%difD	(A12203-003, 99ALOKA, "Long Axis(at End Diastole or End Systole) Length % Difference")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	9
%difS	(A12203-003, 99ALOKA, "Long Axis(at End Diastole or End Systole) Length % Difference")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(109070, DCM, "End Systole")	9
GLS(ap4C)	(A12203-013, 99ALOKA, "Global Longitudinal Strain")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125208, DCM, "Method of Disks, Single Plane"), Image View=(G-A19C, SRT, "Apical four chamber")	9
GLS(ap2C)	(A12203-013, 99ALOKA, "Global Longitudinal Strain")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-001, 99ALOKA, "Method of Disks, Single Plane with Apical two chamber"), Image View=(G-A19B, SRT, "Apical two chamber")	9
(LV Volume (Teichholz))			
EDV ²	(18026-5, LN, "Left Ventricular End Diastolic Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125209, DCM, "Teichholz"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	9

Label	Concept Name	Modifier	row
ESV ²	(18148-7, LN, "Left Ventricular End Systolic Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125209, DCM, "Teichholz"), Cardiac Cycle Point=(109070, DCM, "End Systole")	9
SV ²	(F-32120, SRT, "Stroke Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125209, DCM, "Teichholz")	9
SVI ²	(F-00078, SRT, "Stroke Index")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125209, DCM, "Teichholz")	9
CO ²	(F-32100, SRT, "Cardiac Output")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125209, DCM, "Teichholz")	9
COI ²	(F-32110, SRT, "Cardiac Index")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125209, DCM, "Teichholz")	9
EF ²	(18043-0, LN, "Left Ventricular Ejection Fraction")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125209, DCM, "Teichholz")	9
(LV Volume (Pombo, Teichholz, Gibson))			
FS ²	(18051-3, LN, "Left Ventricular Fractional Shortening")	Image Mode(Group)=(G-03A2, SRT, "2D mode"),	9
mFS ²	(A12203-004, 99ALOKA, "Midwall Fractional Shortening")	Image Mode(Group)=(G-03A2, SRT, "2D mode"),	9
(LV Function (Teichholz))			
EDV ²	(18026-5, LN, "Left Ventricular End Diastolic Volume")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(125209, DCM, "Teichholz"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	9
ESV ²	(18148-7, LN, "Left Ventricular End Systolic Volume")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(125209, DCM, "Teichholz"), Cardiac Cycle Point=(109070, DCM, "End Systole")	9
SV ²	(F-32120, SRT, "Stroke Volume")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(125209, DCM, "Teichholz")	9
SVI ²	(F-00078, SRT, "Stroke Index")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(125209, DCM, "Teichholz")	9
CO ²	(F-32100, SRT, "Cardiac Output")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(125209, DCM, "Teichholz")	9

Label	Concept Name	Modifier	row
COI ²	(F-32110, SRT, "Cardiac Index")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(125209, DCM, "Teichholz")	9
EF ²	(18043-0, LN, "Left Ventricular Ejection Fraction")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(125209, DCM, "Teichholz")	9
(LV Function (Pombo, Teichholz, Gibson))			
FS ²	(18051-3, LN, "Left Ventricular Fractional Shortening")	Image Mode(Group)=(G-0394, SRT, "M mode")	9
mFS ²	(A12203-004, 99ALOKA, "Midwall Fractional Shortening")	Image Mode(Group)=(G-0394, SRT, "M mode")	9
(LV Volume (Pombo))			
EDV ²	(18026-5, LN, "Left Ventricular End Diastolic Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125206, DCM, "Cube Method"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	9
ESV ²	(18148-7, LN, "Left Ventricular End Systolic Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125206, DCM, "Cube Method"), Cardiac Cycle Point=(109070, DCM, "End Systole")	9
SV ²	(F-32120, SRT, "Stroke Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125206, DCM, "Cube Method")	9
SVI ²	(F-00078, SRT, "Stroke Index")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125206, DCM, "Cube Method")	9
CO ²	(F-32100, SRT, "Cardiac Output")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125206, DCM, "Cube Method")	9
COI ²	(F-32110, SRT, "Cardiac Index")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125206, DCM, "Cube Method")	9
EF ²	(18043-0, LN, "Left Ventricular Ejection Fraction")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125206, DCM, "Cube Method")	9
(LV Function (Pombo))			
EDV ²	(18026-5, LN, "Left Ventricular End Diastolic Volume")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(125206, DCM, "Cube Method"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	9

Label	Concept Name	Modifier	row
ESV ²	(18148-7, LN, "Left Ventricular End Systolic Volume")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(125206, DCM, "Cube Method"), Cardiac Cycle Point=(109070, DCM, "End Systole")	9
SV ²	(F-32120, SRT, "Stroke Volume")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(125206, DCM, "Cube Method")	9
SVI ²	(F-00078, SRT, "Stroke Index")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(125206, DCM, "Cube Method")	9
CO ²	(F-32100, SRT, "Cardiac Output")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(125206, DCM, "Cube Method")	9
COI ²	(F-32110, SRT, "Cardiac Index")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(125206, DCM, "Cube Method")	9
EF ²	(18043-0, LN, "Left Ventricular Ejection Fraction")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(125206, DCM, "Cube Method")	9
(LV Volume (Gibson))			
EDV ²	(18026-5, LN, "Left Ventricular End Diastolic Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-002, 99ALOKA, "Gibson"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	9
ESV ²	(18148-7, LN, "Left Ventricular End Systolic Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-002, 99ALOKA, "Gibson"), Cardiac Cycle Point=(109070, DCM, "End Systole")	9
SV ²	(F-32120, SRT, "Stroke Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-002, 99ALOKA, "Gibson")	9
SVI ²	(F-00078, SRT, "Stroke Index")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-002, 99ALOKA, "Gibson")	9
CO ²	(F-32100, SRT, "Cardiac Output")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-002, 99ALOKA, "Gibson")	9
COI ²	(F-32110, SRT, "Cardiac Index")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-002, 99ALOKA, "Gibson")	9
EF ²	(18043-0, LN, "Left Ventricular Ejection Fraction")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-002, 99ALOKA, "Gibson")	9
(LV Function (Gibson))			

Label	Concept Name	Modifier	row
EDV ²	(18026-5, LN, "Left Ventricular End Diastolic Volume")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(A12228-002, 99ALOKA, "Gibson"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	9
ESV ²	(18148-7, LN, "Left Ventricular End Systolic Volume")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(A12228-002, 99ALOKA, "Gibson"), Cardiac Cycle Point=(109070, DCM, "End Systole")	9
SV ²	(F-32120, SRT, "Stroke Volume")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(A12228-002, 99ALOKA, "Gibson")	9
SVI ²	(F-00078, SRT, "Stroke Index")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(A12228-002, 99ALOKA, "Gibson")	9
CO ²	(F-32100, SRT, "Cardiac Output")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(A12228-002, 99ALOKA, "Gibson")	9
COI ²	(F-32110, SRT, "Cardiac Index")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(A12228-002, 99ALOKA, "Gibson")	9
EF ²	(18043-0, LN, "Left Ventricular Ejection Fraction")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(A12228-002, 99ALOKA, "Gibson")	9
(Area-Length)			
EDV	(18026-5, LN, "Left Ventricular End Diastolic Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125205, DCM, "Area-Length Single Plane"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	9
ESV	(18148-7, LN, "Left Ventricular End Systolic Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125205, DCM, "Area-Length Single Plane"), Cardiac Cycle Point=(109070, DCM, "End Systole")	9
SV	(F-32120, SRT, "Stroke Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125205, DCM, "Area-Length Single Plane")	9
SVI	(F-00078, SRT, "Stroke Index")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125205, DCM, "Area-Length Single Plane")	9
CO	(F-32100, SRT, "Cardiac Output")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125205, DCM, "Area-Length Single Plane")	9
COI	(F-32110, SRT, "Cardiac Index")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125205, DCM, "Area-Length Single Plane")	9
EF	(18043-0, LN, "Left Ventricular Ejection Fraction")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125205, DCM, "Area-Length Single Plane")	9

Label	Concept Name	Modifier	row
areaEF	(G-0376, SRT, "Left Ventricular Fractional Area Change")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125205, DCM, "Area-Length Single Plane")	9
(BP-Ellipse)			
EDV	(18026-5, LN, "Left Ventricular End Diastolic Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125204, DCM, "Area-Length Biplane"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	9
ESV	(18148-7, LN, "Left Ventricular End Systolic Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125204, DCM, "Area-Length Biplane"), Cardiac Cycle Point=(109070, DCM, "End Systole")	9
SV	(F-32120, SRT, "Stroke Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125204, DCM, "Area-Length Biplane")	9
SVI	(F-00078, SRT, "Stroke Index")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125204, DCM, "Area-Length Biplane")	9
CO	(F-32100, SRT, "Cardiac Output")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125204, DCM, "Area-Length Biplane")	9
COI	(F-32110, SRT, "Cardiac Index")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125204, DCM, "Area-Length Biplane")	9
EF	(18043-0, LN, "Left Ventricular Ejection Fraction")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125204, DCM, "Area-Length Biplane")	9
areaEFIx	(A12203-005, 99ALOKA, "Area Ejection Fraction at Long Axis View")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125204, DCM, "Area-Length Biplane")	9
areaEFsx	(A12203-006, 99ALOKA, "Area Ejection Fraction at Short Axis View")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125204, DCM, "Area-Length Biplane")	9
(Simpson)			
EDV	(18026-5, LN, "Left Ventricular End Diastolic Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-003, 99ALOKA, "Modified Simpson's"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	9
ESV	(18148-7, LN, "Left Ventricular End Systolic Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-003, 99ALOKA, "Modified Simpson's"), Cardiac Cycle Point=(109070, DCM, "End Systole")	9
SV	(F-32120, SRT, "Stroke Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-003, 99ALOKA, "Modified Simpson's")	9

Label	Concept Name	Modifier	row
SVI	(F-00078, SRT, "Stroke Index")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-003, 99ALOKA, "Modified Simpson's")	9
CO	(F-32100, SRT, "Cardiac Output")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-003, 99ALOKA, "Modified Simpson's")	9
COI	(F-32110, SRT, "Cardiac Index")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-003, 99ALOKA, "Modified Simpson's")	9
EF	(18043-0, LN, "Left Ventricular Ejection Fraction")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-003, 99ALOKA, "Modified Simpson's")	9
areaEFmv	(G-0376, SRT, "Left Ventricular Fractional Area Change")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-003, 99ALOKA, "Modified Simpson's"), Image View=(G-039A, SRT, "Parasternal short axis at the Mitral Valve level")	9
areaEFpm	(G-0376, SRT, "Left Ventricular Fractional Area Change")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-003, 99ALOKA, "Modified Simpson's"), Image View=(G-039B, SRT, "Parasternal short axis at the Papillary Muscle level")	9
(Bullet)			
EDV	(18026-5, LN, "Left Ventricular End Diastolic Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-004, 99ALOKA, "Bullet"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	9
ESV	(18148-7, LN, "Left Ventricular End Systolic Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-004, 99ALOKA, "Bullet"), Cardiac Cycle Point=(109070, DCM, "End Systole")	9
SV	(F-32120, SRT, "Stroke Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-004, 99ALOKA, "Bullet")	9
SVI	(F-00078, SRT, "Stroke Index")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-004, 99ALOKA, "Bullet")	9
CO	(F-32100, SRT, "Cardiac Output")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-004, 99ALOKA, "Bullet")	9
COI	(F-32110, SRT, "Cardiac Index")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-004, 99ALOKA, "Bullet")	9
EF	(18043-0, LN, "Left Ventricular Ejection Fraction")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-004, 99ALOKA, "Bullet")	9

Label	Concept Name	Modifier	row
areaEF	(G-0376, SRT, "Left Ventricular Fractional Area Change")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-004, 99ALOKA, "Bullet")	9
(LV Volume, LV Function)			
HR (LV Volume (Simpson(Disc))) ¹	(8867-4, LN, "Heart rate")		9
HR (LV Volume) ^{1 2}	(8867-4, LN, "Heart rate")	Image Mode(Group)=(G-03A2, SRT, "2D mode")	9
HR (LV Function) ^{1 2}	(8867-4, LN, "Heart rate")	Image Mode(Group)=(G-0394, SRT, "M mode")	9
(LV Function (Pombo, Teichholz, Gibson))			
ET	(18041-4, LN, "Aortic Valve Ejection Time")	Image Mode(Group)=(G-0394, SRT, "M mode")	13
MVCF	(A12203-009, 99ALOKA, "Mean Velocity of Circumferential Fiber Shortening")	Image Mode(Group)=(G-0394, SRT, "M mode")	9
(MVA)			
MVA	(G-038E, SRT, "Cardiovascular Orifice Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125220, DCM, "Planimetry"), Flow Direction=(R-42047, SRT, "Antegrade Flow"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-039A, SRT, "Parasternal short axis at the Mitral Valve level")	14
(AVA)			
AVA	(G-038E, SRT, "Cardiovascular Orifice Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125220, DCM, "Planimetry"), Flow Direction=(R-42047, SRT, "Antegrade Flow"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-0398, SRT, "Parasternal short axis at the aortic valve level")	13
(RVD, LV Volume (Pombo, Teichholz, Gibson))			
RVDd ²	(20304-2, LN, "Right Ventricular Internal Diastolic Dimension")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	10
RVDs ²	(20305-9, LN, "Right Ventricular Internal Systolic Dimension")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(109070, DCM, "End Systole")	10

Label	Concept Name	Modifier	row
(LV Function (Pombo, Teichholz, Gibson))			
RVDd ²	(20304-2, LN, "Right Ventricular Internal Diastolic Dimension")	Image Mode(Group)=(G-0394, SRT, "M mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	10
RVDs ²	(20305-9, LN, "Right Ventricular Internal Systolic Dimension")	Image Mode(Group)=(G-0394, SRT, "M mode"), Cardiac Cycle Point=(109070, DCM, "End Systole")	10
(Ratio, LV Volume (Pombo, Teichholz, Gibson))			
IVSd ²	(18154-5, LN, "Interventricular Septum Diastolic Thickness")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	9
IVSs ²	(18158-6, LN, "Interventricular Septum Systolic Thickness")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(109070, DCM, "End Systole")	9
%IVSTF ²	(18054-7, LN, "Interventricular Septum % Thickening")	Image Mode(Group)=(G-03A2, SRT, "2D mode")	9
LVPWd ²	(18152-9, LN, "Left Ventricle Posterior Wall Diastolic Thickness")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	9
LVPWs ²	(18156-0, LN, "Left Ventricle Posterior Wall Systolic Thickness")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(109070, DCM, "End Systole")	9
%PWTF ²	(18053-9, LN, "Left Ventricle Posterior Wall % Thickening")	Image Mode(Group)=(G-03A2, SRT, "2D mode")	9
IVS/LVPW ²	(18155-2, LN, "Interventricular Septum to Posterior Wall Thickness Ratio")	Image Mode(Group)=(G-03A2, SRT, "2D mode")	9
(LV Function (Pombo, Teichholz, Gibson))			
IVSd ²	(18154-5, LN, "Interventricular Septum Diastolic Thickness")	Image Mode(Group)=(G-0394, SRT, "M mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	9
IVSs ²	(18158-6, LN, "Interventricular Septum Systolic Thickness")	Image Mode(Group)=(G-0394, SRT, "M mode"), Cardiac Cycle Point=(109070, DCM, "End Systole")	9
%IVSTF ²	(18054-7, LN, "Interventricular Septum % Thickening")	Image Mode(Group)=(G-0394, SRT, "M mode")	9

Label	Concept Name	Modifier	row
LVPWd ²	(18152-9, LN, "Left Ventricle Posterior Wall Diastolic Thickness")	Image Mode(Group)=(G-0394, SRT, "M mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	9
LVPWs ²	(18156-0, LN, "Left Ventricle Posterior Wall Systolic Thickness")	Image Mode(Group)=(G-0394, SRT, "M mode"), Cardiac Cycle Point=(109070, DCM, "End Systole")	9
%PWTF ²	(18053-9, LN, "Left Ventricle Posterior Wall % Thickening")	Image Mode(Group)=(G-0394, SRT, "M mode")	9
IVS/LVPW ²	(18155-2, LN, "Interventricular Septum to Posterior Wall Thickness Ratio")	Image Mode(Group)=(G-0394, SRT, "M mode")	9
(LVOT Flow)			
LVOT	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
CSA(LVOT)	(G-038E, SRT, "Cardiovascular Orifice Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
(RVOT Flow)			
RVOT	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract"), Image Mode=(G-03A2, SRT, "2D mode")	10
CSA(RVOT)	(G-038E, SRT, "Cardiovascular Orifice Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	10
(RVD)			
RVAWd	(18153-7, LN, "Right Ventricular Anterior Wall Diastolic Thickness")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	10
RVAWs	(18157-8, LN, "Right Ventricular Anterior Wall Systolic Thickness")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(109070, DCM, "End Systole")	10
(B - IVC)			
Insp ²	(18006-7, LN, "Inferior Vena Cava Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Respiratory Cycle Point=(F-20010, SRT, "During Inspiration")	19
Exp ²	(18006-7, LN, "Inferior Vena Cava Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Respiratory Cycle Point=(F-20020, SRT, "During Expiration")	19
%Collapse ²	(18050-5, LN, "Inferior Vena Cava % Collapse")	Image Mode(Group)=(G-03A2, SRT, "2D mode")	19
(M - IVC)			
Insp ²	(18006-7, LN, "Inferior Vena Cava Diameter")	Image Mode(Group)=(G-0394, SRT, "M mode"), Respiratory Cycle Point=(F-20010, SRT, "During Inspiration")	19

Label	Concept Name	Modifier	row
Exp ²	(18006-7, LN, "Inferior Vena Cava Diameter")	Image Mode(Group)=(G-0394, SRT, "M mode"), Respiratory Cycle Point=(F-20020, SRT, "During Expiration")	19
%Collapse ²	(18050-5, LN, "Inferior Vena Cava % Collapse")	Image Mode(Group)=(G-0394, SRT, "M mode")	19
(B - LA/AO)			
LADd ²	(A12205-001, 99ALOKA, "Left Atrium Antero-posterior Diastolic Dimension")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	11
LADs ²	(29469-4, LN, "Left Atrium Antero-posterior Systolic Dimension")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(109070, DCM, "End Systole")	11
AODd ²	(18015-8, LN, "Aortic Root Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	17
AODs ²	(18015-8, LN, "Aortic Root Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(109070, DCM, "End Systole")	17
LADs/AODd ²	(17985-3, LN, "Left Atrium to Aortic Root Ratio")	Image Mode(Group)=(G-03A2, SRT, "2D mode")	11
(M - LA/AO)			
LADd ²	(A12205-001, 99ALOKA, "Left Atrium Antero-posterior Diastolic Dimension")	Image Mode(Group)=(G-0394, SRT, "M mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	11
LADs ²	(29469-4, LN, "Left Atrium Antero-posterior Systolic Dimension")	Image Mode(Group)=(G-0394, SRT, "M mode"), Cardiac Cycle Point=(109070, DCM, "End Systole")	11
AODd ²	(18015-8, LN, "Aortic Root Diameter")	Image Mode(Group)=(G-0394, SRT, "M mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	17
AODs ²	(18015-8, LN, "Aortic Root Diameter")	Image Mode(Group)=(G-0394, SRT, "M mode"), Cardiac Cycle Point=(109070, DCM, "End Systole")	17
AVDs ²	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-0394, SRT, "M mode"), Finding Site=(T-35410, SRT, "Aortic Valve Ring"), Flow Direction=(R-42047, SRT, "Antegrade Flow"), Cardiac Cycle Point=(F-32020, SRT, "Systole")	13
LADs/AODd	(17985-3, LN, "Left Atrium to Aortic Root Ratio")	Image Mode(Group)=(G-0394, SRT, "M mode")	11
(Mitral V)			
C-Eamp	(A12207-001, 99ALOKA, "Mitral Valve Dimension of C point to E point by M-Mode")	Image Mode(Group)=(G-0394, SRT, "M mode"), Image Mode=(G-0394, SRT, "M mode")	14

Label	Concept Name	Modifier	row
C-Aamp	(A12207-002, 99ALOKA, "Mitral Valve Dimension of C point to A point by M-Mode")	Image Mode(Group)=(G-0394, SRT, "M mode"), Image Mode=(G-0394, SRT, "M mode")	14
EPSS	(18036-4, LN, "Mitral Valve EPSS, E wave")	Image Mode(Group)=(G-0394, SRT, "M mode"), Image Mode=(G-0394, SRT, "M mode")	14
E-Fslop	(18040-6, LN, "Mitral Valve E-F Slope by M-Mode")	Image Mode(Group)=(G-0394, SRT, "M mode"), Image Mode=(G-0394, SRT, "M mode")	14
A/E	(A12207-003, 99ALOKA, "Mitral Valve C-A Dimension to C-E Dimension Ratio by M-Mode")	Image Mode(Group)=(G-0394, SRT, "M mode")	14
E/A	(A12207-004, 99ALOKA, "Mitral Valve C-E Dimension to C-A Dimension Ratio by M-Mode")	Image Mode(Group)=(G-0394, SRT, "M mode")	14
(Tricuspid V)			
C-Eamp	(A12208-001, 99ALOKA, "Tricuspid Valve Dimension of C point to E point by M-Mode")	Image Mode(Group)=(G-0394, SRT, "M mode"), Image Mode=(G-0394, SRT, "M mode")	16
C-Aamp	(A12208-002, 99ALOKA, "Tricuspid Valve Dimension of C point to A point by M-Mode")	Image Mode(Group)=(G-0394, SRT, "M mode"), Image Mode=(G-0394, SRT, "M mode")	16
D-Eamp	(A12208-003, 99ALOKA, "Tricuspid Valve Dimension of D point to E point by M-Mode")	Image Mode(Group)=(G-0394, SRT, "M mode"), Image Mode=(G-0394, SRT, "M mode")	16
E-Fslop	(A12208-004, 99ALOKA, "Tricuspid Valve Velocity from E point to F point by M-Mode")	Image Mode(Group)=(G-0394, SRT, "M mode"), Image Mode=(G-0394, SRT, "M mode")	16
D-Eslop	(A12208-005, 99ALOKA, "Tricuspid Valve Velocity from D point to E point by M-Mode")	Image Mode(Group)=(G-0394, SRT, "M mode"), Image Mode=(G-0394, SRT, "M mode")	16
A/E	(A12208-006, 99ALOKA, "Tricuspid Valve C-A Dimension to C-E Dimension Ratio by M-Mode")	Image Mode(Group)=(G-0394, SRT, "M mode")	16
E/A	(A12208-007, 99ALOKA, "Tricuspid Valve C-E Dimension to C-A Dimension Ratio by M-Mode")	Image Mode(Group)=(G-0394, SRT, "M mode")	16
(Pulmonary V)			

Label	Concept Name	Modifier	row
A wave amp	(A12209-001, 99ALOKA, "Pulmonic Valve Dimension of F point to A point by M-Mode")	Image Mode(Group)=(G-0394, SRT, "M mode"), Image Mode=(G-0394, SRT, "M mode")	15
B-Camp	(A12209-002, 99ALOKA, "Pulmonic Valve Dimension of B point to C point by M-Mode")	Image Mode(Group)=(G-0394, SRT, "M mode"), Image Mode=(G-0394, SRT, "M mode")	15
E-Fslop	(A12209-003, 99ALOKA, "Pulmonic Valve Velocity from E point to F point by M-Mode")	Image Mode(Group)=(G-0394, SRT, "M mode"), Image Mode=(G-0394, SRT, "M mode")	15
B-Cslop	(A12209-004, 99ALOKA, "Pulmonic Valve Velocity from B point to C point by M-Mode")	Image Mode(Group)=(G-0394, SRT, "M mode"), Image Mode=(G-0394, SRT, "M mode")	15
(LVOT Flow)			
pV	(11726-7, LN, "Peak Velocity")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
MnV	(20352-1, LN, "Mean Velocity")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
MPG	(20256-4, LN, "Mean Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
VTI	(20354-7, LN, "Velocity Time Integral")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
AccT	(20168-1, LN, "Acceleration Time")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
PEP	(A12222-001, 99ALOKA, "Pre-Ejection Period")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
ET	(18041-4, LN, "Aortic Valve Ejection Time")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	13
HR	(8867-4, LN, "Heart rate")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
AccT/ET	(G-0382, SRT, "Ratio of Aortic Valve Acceleration Time to Ejection Time")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	13
PG	(20247-3, LN, "Peak Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
SV(LVOT)	(F-32120, SRT, "Stroke Volume")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
CO(LVOT)	(F-32100, SRT, "Cardiac Output")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
PEP/ET	(A12222-002, 99ALOKA, "PEP/ET")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9

Label	Concept Name	Modifier	row
SVI(LVOT)	(F-00078, SRT, "Stroke Index")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
COI(LVOT)	(F-32110, SRT, "Cardiac Index")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
(RVOT Flow)			
pV	(11726-7, LN, "Peak Velocity")	Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	10
MnV	(20352-1, LN, "Mean Velocity")	Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	10
MPG	(20256-4, LN, "Mean Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	10
VTI	(20354-7, LN, "Velocity Time Integral")	Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	10
AccT	(20168-1, LN, "Acceleration Time")	Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	10
PEP	(A12222-001, 99ALOKA, "Pre-Ejection Period")	Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	10
ET	(18042-2, LN, "Pulmonic Valve Ejection Time")	Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	15
HR	(8867-4, LN, "Heart rate")	Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	10
AccT/ET	(G-0388, SRT, "Ratio of Pulmonic Valve Acceleration Time to Ejection Time")	Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	15
PG	(20247-3, LN, "Peak Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	10
SV(RVOT)	(F-32120, SRT, "Stroke Volume")	Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	10
CO(RVOT)	(F-32100, SRT, "Cardiac Output")	Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	10
PEP/ET	(A12222-002, 99ALOKA, "PEP/ET")	Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	10
Qp/Qs	(29462-9, LN, "Pulmonary-to-Systemic Shunt Flow Ratio")		22
SVI(RVOT)	(F-00078, SRT, "Stroke Index")	Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	10
COI(RVOT)	(F-32110, SRT, "Cardiac Index")	Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	10
(Trans M Flow)			
eV	(18037-2, LN, "Mitral Valve E-Wave Peak Velocity")		14

Label	Concept Name	Modifier	row
aV	(17978-8, LN, "Mitral Valve A-Wave Peak Velocity")		14
MnV	(20352-1, LN, "Mean Velocity")	Image Mode(Group)=(R-409E4, SRT, "Doppler Pulsed"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	14
MPG	(20256-4, LN, "Mean Gradient")	Image Mode(Group)=(R-409E4, SRT, "Doppler Pulsed"), Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	14
AccT	(20168-1, LN, "Acceleration Time")		14
P1/2T	(20280-4, LN, "Pressure Half-Time")	Image Mode(Group)=(R-409E4, SRT, "Doppler Pulsed"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	14
IRT	(18071-1, LN, "Left Ventricular Isovolumic Relaxation Time")		9
DecT	(G-0384, SRT, "Mitral Valve E-Wave Deceleration Time")		14
Edur	(A12207-005, 99ALOKA, "Mitral Valve E-Wave Duration")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	14
Adur	(G-0385, SRT, "Mitral Valve A-Wave Duration")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	14
VTI	(20354-7, LN, "Velocity Time Integral")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	14
LVDF	(A12203-010, 99ALOKA, "Left Ventricle Diastole Filling Time")		9
RR	(122182, DCM, "R-R interval")		9
MVA(P1/2T)	(G-038E, SRT, "Cardiovascular Orifice Area")	Image Mode(Group)=(R-409E4, SRT, "Doppler Pulsed"), Measurement Method=(125210, DCM, "Area by Pressure Half-Time"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	14
E/A	(18038-0, LN, "Mitral Valve E to A Ratio")		14
A/E	(A12207-006, 99ALOKA, "Mitral Valve A to E Ratio")		14
EPG	(A12207-007, 99ALOKA, "Mitral Valve E-wave Peak Pressure Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli")	14
APG	(A12207-008, 99ALOKA, "Mitral Valve A-wave Peak Pressure Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli")	14

Label	Concept Name	Modifier	row
E/Em	(A12203-011, 99ALOKA, "Ratio of MV E-Wave Peak Vel. to Early Diastolic Myocardium Vel.")		9
PVAdur-Adur	(A12207-009, 99ALOKA, "Subtraction of A-wave Duration from PVA-wave Duration")		14
LVDF/RR	(A12203-012, 99ALOKA, "Ratio of Left Ventricle Diastole Filling Time to R-R interval")		9
(AS Flow)			
pV	(11726-7, LN, "Peak Velocity")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
MnV	(20352-1, LN, "Mean Velocity")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
MPG	(20256-4, LN, "Mean Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
VTI	(20354-7, LN, "Velocity Time Integral")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
AVA	(G-038E, SRT, "Cardiovascular Orifice Area")	Measurement Method=(125215, DCM, "Continuity Equation by Velocity Time Integral"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
PG	(20247-3, LN, "Peak Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
AVA/BSA	(79959-3, LN, "Aortic valve area (Continuity VTI) / BSA")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
(PS Flow)			
pV	(11726-7, LN, "Peak Velocity")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	15
MnV	(20352-1, LN, "Mean Velocity")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	15
MPG	(20256-4, LN, "Mean Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	15
PG	(20247-3, LN, "Peak Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	15
(MS Flow)			
pV	(11726-7, LN, "Peak Velocity")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	14

Label	Concept Name	Modifier	row
MnV	(20352-1, LN, "Mean Velocity")	Image Mode(Group)=(R-409E3, SRT, "Doppler Continuous Wave"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	14
MPG	(20256-4, LN, "Mean Gradient")	Image Mode(Group)=(R-409E3, SRT, "Doppler Continuous Wave"), Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	14
FlowT	(A12222-003, 99ALOKA, "Flow Time")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	14
P1/2T	(20280-4, LN, "Pressure Half-Time")	Image Mode(Group)=(R-409E3, SRT, "Doppler Continuous Wave"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	14
PG	(20247-3, LN, "Peak Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	14
MVA(P1/2T)	(G-038E, SRT, "Cardiovascular Orifice Area")	Image Mode(Group)=(R-409E3, SRT, "Doppler Continuous Wave"), Measurement Method=(125210, DCM, "Area by Pressure Half-Time"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	14
MVA(VTI)	(G-038E, SRT, "Cardiovascular Orifice Area")	Measurement Method=(125215, DCM, "Continuity Equation by Velocity Time Integral"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	14
VTI	(20354-7, LN, "Velocity Time Integral")	Image Mode(Group)=(R-409E3, SRT, "Doppler Continuous Wave"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	14
(TS Flow)			
pV	(11726-7, LN, "Peak Velocity")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
MnV	(20352-1, LN, "Mean Velocity")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
MPG	(20256-4, LN, "Mean Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
FlowT	(A12222-003, 99ALOKA, "Flow Time")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
P1/2T	(20280-4, LN, "Pressure Half-Time")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
PG	(20247-3, LN, "Peak Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
(AR Flow)			

Label	Concept Name	Modifier	row
pV	(11726-7, LN, "Peak Velocity")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	13
MnV	(20352-1, LN, "Mean Velocity")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	13
MPG	(20256-4, LN, "Mean Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	13
P1/2T	(20280-4, LN, "Pressure Half-Time")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	13
PG	(20247-3, LN, "Peak Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	13
(MR Flow)			
pV	(11726-7, LN, "Peak Velocity")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	14
MnV	(20352-1, LN, "Mean Velocity")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	14
MPG	(20256-4, LN, "Mean Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	14
FlowT	(A12222-003, 99ALOKA, "Flow Time")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	14
PG	(20247-3, LN, "Peak Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	14
dP/dt	(18035-6, LN, "Mitral Regurgitation dP/dt derived from Mitral Reg. velocity")		14
Reg. Volume	(A12207-014, 99ALOKA, "Mitral valve Regurgitant volume")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	14
Reg. Frac	(G-0390, SRT, "Regurgitant Fraction")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	14
EROA	(G-038E, SRT, "Cardiovascular Orifice Area")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	14
(PR Flow)			
pV	(11726-7, LN, "Peak Velocity")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	15
MnV	(20352-1, LN, "Mean Velocity")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	15
MPG	(20256-4, LN, "Mean Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	15

Label	Concept Name	Modifier	row
PG	(20247-3, LN, "Peak Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	15
(PV Flow)			
PVS	(29450-4, LN, "Pulmonary Vein Systolic Peak Velocity")		20
PVD	(29451-2, LN, "Pulmonary Vein Diastolic Peak Velocity")		20
PVA	(29453-8, LN, "Pulmonary Vein Atrial Contraction Reversal Peak Velocity")		20
S-VTI	(G-038C, SRT, "Pulmonary Vein S-Wave Velocity Time Integral")		20
D-VTI	(G-038D, SRT, "Pulmonary Vein D-Wave Velocity Time Integral")		20
PVAdur	(G-038B, SRT, "Pulmonary Vein A-Wave Duration")		20
DecT	(A12214-001, 99ALOKA, "Deceleration Time of D-Wave Flow")		20
S/D	(29452-0, LN, "Pulmonary Vein Systolic to Diastolic Ratio")		20
SF	(A12214-002, 99ALOKA, "Systolic Fraction")		20
(AR Vol. PISA)			
PISA Radius	(A12222-004, 99ALOKA, "Radius of Flow Convergence")		13
AR Alias V (Vr)	(A12222-005, 99ALOKA, "Aliasing Velocity")		13
Angle (PISA)	(A12222-006, 99ALOKA, "Proximal Isovelocity Surface Area Angle")		13
VTI(AR)	(A12211-001, 99ALOKA, "Velocity Time Integral of Aortic Regurgitant Flow")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	13
AR Vol	(33878-0, LN, "Volume Flow")	Measurement Method=(125216, DCM, "Proximal Isovelocity Surface Area"), Finding Site=(T-35410, SRT, "Aortic Valve Ring"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	13

Label	Concept Name	Modifier	row
AR EROA	(G-038E, SRT, "Cardiovascular Orifice Area")	Measurement Method=(125216, DCM, "Proximal Isovelocity Surface Area"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	13
AR Flow Rt	(34141-2, LN, "Peak Instantaneous Flow Rate")	Measurement Method=(125216, DCM, "Proximal Isovelocity Surface Area"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	13
AR PISA	(A12211-003, 99ALOKA, "Aortic Regurgitant Proximal Isovelocity Surface Area")	Measurement Method=(125216, DCM, "Proximal Isovelocity Surface Area"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	13
AR RF	(G-0390, SRT, "Regurgitant Fraction")	Measurement Method=(125216, DCM, "Proximal Isovelocity Surface Area"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	13
(MR Vol. PISA)			
PISA Radius	(A12222-004, 99ALOKA, "Radius of Flow Convergence")		14
MR Alias V (Vr)	(A12222-005, 99ALOKA, "Aliasing Velocity")		14
Angle (PISA)	(A12222-006, 99ALOKA, "Proximal Isovelocity Surface Area Angle")		14
VTI(MR)	(A12207-010, 99ALOKA, "Velocity Time Integral of Mitral Regurgitant Flow")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	14
MV Diam	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Finding Site=(T-35313, SRT, "Mitral Annulus"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	14
VTI(MVannu)	(20354-7, LN, "Velocity Time Integral")	Finding Site=(T-35313, SRT, "Mitral Annulus"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	14
MV DiamA4C	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Finding Site=(T-35313, SRT, "Mitral Annulus"), Flow Direction=(R-42047, SRT, "Antegrade Flow"), Image View=(G-A19C, SRT, "Apical four chamber")	14
MV DiamA2C	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Finding Site=(T-35313, SRT, "Mitral Annulus"), Flow Direction=(R-42047, SRT, "Antegrade Flow"), Image View=(G-A19B, SRT, "Apical two chamber")	14
MR Vol	(33878-0, LN, "Volume Flow")	Measurement Method=(125216, DCM, "Proximal Isovelocity Surface Area"), Finding Site=(T-35313, SRT, "Mitral Annulus"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	14

Label	Concept Name	Modifier	row
MR EROA	(G-038E, SRT, "Cardiovascular Orifice Area")	Measurement Method=(125216, DCM, "Proximal Isovelocity Surface Area"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	14
MR Flow Rt	(34141-2, LN, "Peak Instantaneous Flow Rate")	Measurement Method=(125216, DCM, "Proximal Isovelocity Surface Area"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	14
MR PISA	(A12207-012, 99ALOKA, "Mitral Regurgitant Proximal Isovelocity Surface Area")	Measurement Method=(125216, DCM, "Proximal Isovelocity Surface Area"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	14
MR RF	(G-0390, SRT, "Regurgitant Fraction")	Measurement Method=(125216, DCM, "Proximal Isovelocity Surface Area"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	14
SV(MV)	(A12207-013, 99ALOKA, "Flow Volume of Mitral Valve Annulus in Flow")	Measurement Method=(125219, DCM, "Doppler Volume Flow")	14
(TR Vol. PISA)			
PISA Radius	(A12222-004, 99ALOKA, "Radius of Flow Convergence")		16
TR Alias V (Vr)	(A12222-005, 99ALOKA, "Aliasing Velocity")		16
Angle (PISA)	(A12222-006, 99ALOKA, "Proximal Isovelocity Surface Area Angle")		16
VTI(TR)	(A12208-008, 99ALOKA, "Velocity Time Integral of Tricuspid Regurgitant Flow")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	16
TV Diam	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Finding Site=(T-35111, SRT, "Tricuspid Annulus"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
VTI(TVannu)	(20354-7, LN, "Velocity Time Integral")	Finding Site=(T-35111, SRT, "Tricuspid Annulus"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
TR Vol	(33878-0, LN, "Volume Flow")	Measurement Method=(125216, DCM, "Proximal Isovelocity Surface Area"), Finding Site=(T-35111, SRT, "Tricuspid Annulus"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	16
TR EROA	(G-038E, SRT, "Cardiovascular Orifice Area")	Measurement Method=(125216, DCM, "Proximal Isovelocity Surface Area"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	16

Label	Concept Name	Modifier	row
TR Flow Rt	(34141-2, LN, "Peak Instantaneous Flow Rate")	Measurement Method=(125216, DCM, "Proximal Isovelocity Surface Area"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	16
TR PISA	(A12208-010, 99ALOKA, "Tricuspid Regurgitant Proximal Isovelocity Surface Area")	Measurement Method=(125216, DCM, "Proximal Isovelocity Surface Area"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	16
TR RF	(G-0390, SRT, "Regurgitant Fraction")	Measurement Method=(125216, DCM, "Proximal Isovelocity Surface Area"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	16
SV(TV)	(A12208-011, 99ALOKA, "Flow Volume of Tricuspid Valve Annulus in Flow")	Measurement Method=(125219, DCM, "Doppler Volume Flow")	16
(PR Vol. PISA)			
PISA Radius	(A12222-004, 99ALOKA, "Radius of Flow Convergence")		15
PR Alias V (Vr)	(A12222-005, 99ALOKA, "Aliasing Velocity")		15
Angle (PISA)	(A12222-006, 99ALOKA, "Proximal Isovelocity Surface Area Angle")		15
VTI(PR)	(A12209-005, 99ALOKA, "Velocity Time Integral of Pulmonic Regurgitant Flow")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	15
PR Vol	(33878-0, LN, "Volume Flow")	Measurement Method=(125216, DCM, "Proximal Isovelocity Surface Area"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	15
PR EROA	(G-038E, SRT, "Cardiovascular Orifice Area")	Measurement Method=(125216, DCM, "Proximal Isovelocity Surface Area"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	15
PR Flow Rt	(34141-2, LN, "Peak Instantaneous Flow Rate")	Measurement Method=(125216, DCM, "Proximal Isovelocity Surface Area"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	15
PR PISA	(A12209-007, 99ALOKA, "Pulmonic Regurgitant Proximal Isovelocity Surface Area")	Measurement Method=(125216, DCM, "Proximal Isovelocity Surface Area"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	15
PR RF	(G-0390, SRT, "Regurgitant Fraction")	Measurement Method=(125216, DCM, "Proximal Isovelocity Surface Area"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	15

Label	Concept Name	Modifier	row
(LA Volume Simpson(Disc))			
LAL4s	(A12205-002, 99ALOKA, "Left Atrium systolic major axis")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-005, 99ALOKA, "Method of Disks, Biplane of LA"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19C, SRT, "Apical four chamber")	11
LALA4s	(17977-0, LN, "Left Atrium Systolic Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-005, 99ALOKA, "Method of Disks, Biplane of LA"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19C, SRT, "Apical four chamber")	11
LAL2s	(A12205-002, 99ALOKA, "Left Atrium systolic major axis")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-005, 99ALOKA, "Method of Disks, Biplane of LA"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19B, SRT, "Apical two chamber")	11
LALA2s	(A12205-003, 99ALOKA, "Left Atrium Systolic Area by Apical two chamber")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-005, 99ALOKA, "Method of Disks, Biplane of LA"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19B, SRT, "Apical two chamber")	11
LAS4s	(A12205-007, 99ALOKA, "Left Atrium systolic minor axis")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-005, 99ALOKA, "Method of Disks, Biplane of LA"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19C, SRT, "Apical four chamber")	11
LAS2s	(A12205-007, 99ALOKA, "Left Atrium systolic minor axis")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-005, 99ALOKA, "Method of Disks, Biplane of LA"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19B, SRT, "Apical two chamber")	11
LA Volume(Bi-plane)	(A12205-004, 99ALOKA, "Left Atrial Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-005, 99ALOKA, "Method of Disks, Biplane of LA"), Cardiac Cycle Point=(109070, DCM, "End Systole")	11

Label	Concept Name	Modifier	row
LA Volume/BSA(Bi-plane)	(A12205-005, 99ALOKA, "Left Atrial Volume divided by Body Surface Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-005, 99ALOKA, "Method of Disks, Biplane of LA"), Cardiac Cycle Point=(109070, DCM, "End Systole")	11
LA Volume(ap4C)	(A12205-004, 99ALOKA, "Left Atrial Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-006, 99ALOKA, "Method of Disks, Single Plane with Apical four chamber of LA"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image View=(G-A19C, SRT, "Apical four chamber")	11
LA Volume/BSA(ap4C)	(A12205-005, 99ALOKA, "Left Atrial Volume divided by Body Surface Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-006, 99ALOKA, "Method of Disks, Single Plane with Apical four chamber of LA"), Cardiac Cycle Point=(109070, DCM, "End Systole")	11
LA Volume(ap2C)	(A12205-004, 99ALOKA, "Left Atrial Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-007, 99ALOKA, "Method of Disks, Single Plane with Apical two chamber of LA"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image View=(G-A19B, SRT, "Apical two chamber")	11
LA Volume/BSA(ap2C)	(A12205-005, 99ALOKA, "Left Atrial Volume divided by Body Surface Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-007, 99ALOKA, "Method of Disks, Single Plane with Apical two chamber of LA"), Cardiac Cycle Point=(109070, DCM, "End Systole")	11
%difS	(A12205-006, 99ALOKA, "Long Axis at End Systole Length % Difference of Left Atrium")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-005, 99ALOKA, "Method of Disks, Biplane of LA"), Cardiac Cycle Point=(109070, DCM, "End Systole")	11
(RA Volume Simpson(Disc))			
RAL4s	(A12206-001, 99ALOKA, "Right Atrium systolic major axis")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-008, 99ALOKA, "Method of Disks, Biplane of RA"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19C, SRT, "Apical four chamber")	12

Label	Concept Name	Modifier	row
RALA4s	(17988-7, LN, "Right Atrium Systolic Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-008, 99ALOKA, "Method of Disks, Biplane of RA"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19C, SRT, "Apical four chamber")	12
RAL2s	(A12206-001, 99ALOKA, "Right Atrium systolic major axis")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-008, 99ALOKA, "Method of Disks, Biplane of RA"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19B, SRT, "Apical two chamber")	12
RALA2s	(A12206-002, 99ALOKA, "Right Atrium Systolic Area by Apical two chamber")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-008, 99ALOKA, "Method of Disks, Biplane of RA"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19B, SRT, "Apical two chamber")	12
RAS4s	(A12206-006, 99ALOKA, "Right Atrium systolic minor axis")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-008, 99ALOKA, "Method of Disks, Biplane of RA"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19C, SRT, "Apical four chamber")	12
RAS2s	(A12206-006, 99ALOKA, "Right Atrium systolic minor axis")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-008, 99ALOKA, "Method of Disks, Biplane of RA"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19B, SRT, "Apical two chamber")	12
RA Volume(Bi-plane)	(A12206-003, 99ALOKA, "Right Atrial Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-008, 99ALOKA, "Method of Disks, Biplane of RA"), Cardiac Cycle Point=(109070, DCM, "End Systole")	12
RA Volume/BSA(Bi-plane)	(A12206-004, 99ALOKA, "Right Atrial Volume divided by Body Surface Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-008, 99ALOKA, "Method of Disks, Biplane of RA"), Cardiac Cycle Point=(109070, DCM, "End Systole")	12

Label	Concept Name	Modifier	row
RA Volume(ap4C)	(A12206-003, 99ALOKA, "Right Atrial Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-009, 99ALOKA, "Method of Disks, Single Plane with Apical four chamber of RA"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image View=(G-A19C, SRT, "Apical four chamber")	12
RA Volume/BSA(ap4C)	(A12206-004, 99ALOKA, "Right Atrial Volume divided by Body Surface Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-009, 99ALOKA, "Method of Disks, Single Plane with Apical four chamber of RA"), Cardiac Cycle Point=(109070, DCM, "End Systole")	12
RA Volume(ap2C)	(A12206-003, 99ALOKA, "Right Atrial Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-010, 99ALOKA, "Method of Disks, Single Plane with Apical two chamber of RA"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image View=(G-A19B, SRT, "Apical two chamber")	12
RA Volume/BSA(ap2C)	(A12206-004, 99ALOKA, "Right Atrial Volume divided by Body Surface Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-010, 99ALOKA, "Method of Disks, Single Plane with Apical two chamber of RA"), Cardiac Cycle Point=(109070, DCM, "End Systole")	12
%difS	(A12206-005, 99ALOKA, "Long Axis at End Systole Length % Difference of Right Atrium")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-008, 99ALOKA, "Method of Disks, Biplane of RA"), Cardiac Cycle Point=(109070, DCM, "End Systole")	12
(LA Volume Area-Length)			
LAL4s	(A12205-002, 99ALOKA, "Left Atrium systolic major axis")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-011, 99ALOKA, "Area-Length Biplane of Left Atrium"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19C, SRT, "Apical four chamber")	11
LALA4s	(17977-0, LN, "Left Atrium Systolic Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-011, 99ALOKA, "Area-Length Biplane of Left Atrium"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19C, SRT, "Apical four chamber")	11

Label	Concept Name	Modifier	row
LAL2s	(A12205-002, 99ALOKA, "Left Atrium systolic major axis")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-011, 99ALOKA, "Area-Length Biplane of Left Atrium"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19B, SRT, "Apical two chamber")	11
LALA2s	(A12205-003, 99ALOKA, "Left Atrium Systolic Area by Apical two chamber")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-011, 99ALOKA, "Area-Length Biplane of Left Atrium"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19B, SRT, "Apical two chamber")	11
LAS4s	(A12205-007, 99ALOKA, "Left Atrium systolic minor axis")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-011, 99ALOKA, "Area-Length Biplane of Left Atrium"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19C, SRT, "Apical four chamber")	11
LAS2s	(A12205-007, 99ALOKA, "Left Atrium systolic minor axis")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-011, 99ALOKA, "Area-Length Biplane of Left Atrium"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19B, SRT, "Apical two chamber")	11
LA Volume	(A12205-004, 99ALOKA, "Left Atrial Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-011, 99ALOKA, "Area-Length Biplane of Left Atrium"), Cardiac Cycle Point=(109070, DCM, "End Systole")	11
LA Volume/BSA	(A12205-005, 99ALOKA, "Left Atrial Volume divided by Body Surface Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-011, 99ALOKA, "Area-Length Biplane of Left Atrium"), Cardiac Cycle Point=(109070, DCM, "End Systole")	11
%difS	(A12205-006, 99ALOKA, "Long Axis at End Systole Length % Difference of Left Atrium")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-011, 99ALOKA, "Area-Length Biplane of Left Atrium"), Cardiac Cycle Point=(109070, DCM, "End Systole")	11
(RA Volume Area-Length)			

Label	Concept Name	Modifier	row
RAL4s	(A12206-001, 99ALOKA, "Right Atrium systolic major axis")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-012, 99ALOKA, "Area-Length Biplane of Right Atrium"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19C, SRT, "Apical four chamber")	12
RALA4s	(17988-7, LN, "Right Atrium Systolic Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-012, 99ALOKA, "Area-Length Biplane of Right Atrium"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19C, SRT, "Apical four chamber")	12
RAL2s	(A12206-001, 99ALOKA, "Right Atrium systolic major axis")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-012, 99ALOKA, "Area-Length Biplane of Right Atrium"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19B, SRT, "Apical two chamber")	12
RALA2s	(A12206-002, 99ALOKA, "Right Atrium Systolic Area by Apical two chamber")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-012, 99ALOKA, "Area-Length Biplane of Right Atrium"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19B, SRT, "Apical two chamber")	12
RAS4s	(A12206-006, 99ALOKA, "Right Atrium systolic minor axis")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-012, 99ALOKA, "Area-Length Biplane of Right Atrium"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19C, SRT, "Apical four chamber")	12
RAS2s	(A12206-006, 99ALOKA, "Right Atrium systolic minor axis")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-012, 99ALOKA, "Area-Length Biplane of Right Atrium"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19B, SRT, "Apical two chamber")	12
RA Volume	(A12206-003, 99ALOKA, "Right Atrial Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-012, 99ALOKA, "Area-Length Biplane of Right Atrium"), Cardiac Cycle Point=(109070, DCM, "End Systole")	12

Label	Concept Name	Modifier	row
RA Volume/BSA	(A12206-004, 99ALOKA, "Right Atrial Volume divided by Body Surface Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-012, 99ALOKA, "Area-Length Biplane of Right Atrium"), Cardiac Cycle Point=(109070, DCM, "End Systole")	12
%difS	(A12206-005, 99ALOKA, "Long Axis at End Systole Length % Difference of Right Atrium")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-012, 99ALOKA, "Area-Length Biplane of Right Atrium"), Cardiac Cycle Point=(109070, DCM, "End Systole")	12
(TAPSE)			
TAPSE	(A12208-012, 99ALOKA, "Tricuspid Annular Plane Systolic Excursion")	Image Mode(Group)=(G-0394, SRT, "M mode")	16
(FAC)			
RV Area ED	(A12204-001, 99ALOKA, "Right Ventricular Diastolic Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125220, DCM, "Planimetry"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	10
RV Area ES	(A12204-002, 99ALOKA, "Right Ventricular Systolic Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125220, DCM, "Planimetry"), Cardiac Cycle Point=(R-FAB5B, SRT, "End Systole")	10
RV FAC	(A12204-003, 99ALOKA, "Right Ventricular Fractional Area Change")	Image Mode(Group)=(G-03A2, SRT, "2D mode")	10
(2DTT(LV))			
Peak Global	(ALOKA_004-001, 99ALOKA, "Peak GLS")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-013, 99ALOKA, "2D Speckle Tracking, Single Plane"), Image View=(G-A19C, SRT, "Apical four chamber")	9
Peak Global	(ALOKA_004-001, 99ALOKA, "Peak GLS")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-013, 99ALOKA, "2D Speckle Tracking, Single Plane"), Image View=(G-A19B, SRT, "Apical two chamber")	9
Peak Global	(ALOKA_004-001, 99ALOKA, "Peak GLS")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-013, 99ALOKA, "2D Speckle Tracking, Single Plane"), Image View=(G-0395, SRT, "Apical long axis")	9
Peak Global	(ALOKA_004-002, 99ALOKA, "Peak GLS of Overlay")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-014, 99ALOKA, "2D Speckle Tracking, Multi Plane")	9

Label	Concept Name	Modifier	row
Time from R	(ALOKA_004-006, 99ALOKA, "Time from R to Peak GLS phase")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-013, 99ALOKA, "2D Speckle Tracking, Single Plane"), Image View=(G-A19C, SRT, "Apical four chamber")	9
Time from R	(ALOKA_004-006, 99ALOKA, "Time from R to Peak GLS phase")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-013, 99ALOKA, "2D Speckle Tracking, Single Plane"), Image View=(G-A19B, SRT, "Apical two chamber")	9
Time from R	(ALOKA_004-006, 99ALOKA, "Time from R to Peak GLS phase")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-013, 99ALOKA, "2D Speckle Tracking, Single Plane"), Image View=(G-0395, SRT, "Apical long axis")	9
Time from R	(ALOKA_004-007, 99ALOKA, "Time from R to Peak GLS phase of Overlay")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-014, 99ALOKA, "2D Speckle Tracking, Multi Plane")	9
-	(ALOKA_004-003, 99ALOKA, "Averaged Peak GLS")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-014, 99ALOKA, "2D Speckle Tracking, Multi Plane")	9
L. Strain(EN)[%]	(ALOKA_004-004, 99ALOKA, "GLS-Sel.")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-013, 99ALOKA, "2D Speckle Tracking, Single Plane"), Image View=(G-A19C, SRT, "Apical four chamber")	9
L. Strain(EN)[%]	(ALOKA_004-004, 99ALOKA, "GLS-Sel.")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-013, 99ALOKA, "2D Speckle Tracking, Single Plane"), Image View=(G-A19B, SRT, "Apical two chamber")	9
L. Strain(EN)[%]	(ALOKA_004-004, 99ALOKA, "GLS-Sel.")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-013, 99ALOKA, "2D Speckle Tracking, Single Plane"), Image View=(G-0395, SRT, "Apical long axis")	9
L. Strain(EN)[%]	(ALOKA_004-005, 99ALOKA, "GLS-Sel.of Overlay")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-014, 99ALOKA, "2D Speckle Tracking, Multi Plane")	9
Time from R	(ALOKA_004-008, 99ALOKA, "Time from R to GLS-Sel. phase")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-013, 99ALOKA, "2D Speckle Tracking, Single Plane"), Image View=(G-A19C, SRT, "Apical four chamber")	9

Label	Concept Name	Modifier	row
Time from R	(ALOKA_004-008, 99ALOKA, "Time from R to GLS-Sel. phase")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-013, 99ALOKA, "2D Speckle Tracking, Single Plane"), Image View=(G-A19B, SRT, "Apical two chamber")	9
Time from R	(ALOKA_004-008, 99ALOKA, "Time from R to GLS-Sel. phase")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-013, 99ALOKA, "2D Speckle Tracking, Single Plane"), Image View=(G-0395, SRT, "Apical long axis")	9
Time from R	(ALOKA_004-009, 99ALOKA, "Time from R to GLS-Sel. phase of Overlay")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-014, 99ALOKA, "2D Speckle Tracking, Multi Plane")	9
(LV Mass T-E)			
LVM	(18087-7, LN, "Left Ventricle Mass")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125222, DCM, "Left Ventricle Mass by Truncated Ellipse")	9
LVM/BSA	(A12203-001, 99ALOKA, "Left Ventricular Mass Index")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125222, DCM, "Left Ventricle Mass by Truncated Ellipse")	9
LVAd PM Epi	(G-0379, SRT, "Left Ventricle Epicardial Diastolic Area, psax pap view")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125222, DCM, "Left Ventricle Mass by Truncated Ellipse"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-039B, SRT, "Parasternal short axis at the Papillary Muscle level")	9
LVAd PM Endo	(A12240-001, 99ALOKA, "Left Ventricular Endocardial Diastolic Area, psax pap view")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125222, DCM, "Left Ventricle Mass by Truncated Ellipse"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-039B, SRT, "Parasternal short axis at the Papillary Muscle level")	9
a	(399063007, SCT, "Left Ventricle Semi-major Axis Diastolic Dimension")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125222, DCM, "Left Ventricle Mass by Truncated Ellipse"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode")	9

Label	Concept Name	Modifier	row
d	(399309003, SCT, "Left Ventricle Truncated Semi-major Axis Diastolic Dimension")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125222, DCM, "Left Ventricle Mass by Truncated Ellipse"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode")	9
t	(A12201-002, 99ALOKA, "Mean Wall Thickness")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125222, DCM, "Left Ventricle Mass by Truncated Ellipse")	9
b	(A12201-004, 99ALOKA, "Short Axis Radius Calculated from Short Axis Cavity Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125222, DCM, "Left Ventricle Mass by Truncated Ellipse")	9
(LV Mass Pedi)			
LVM	(18087-7, LN, "Left Ventricle Mass")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12227-003, 99ALOKA, "Left Ventricle Mass for Pediatric")	9
LVM/BSA	(A12203-001, 99ALOKA, "Left Ventricular Mass Index")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12227-003, 99ALOKA, "Left Ventricle Mass for Pediatric")	9
EDL	(A12201-003, 99ALOKA, "Left Ventricular End Diastolic Length")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12227-003, 99ALOKA, "Left Ventricle Mass for Pediatric"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode")	9
EDA	(G-0375, SRT, "Left Ventricular Diastolic Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12227-003, 99ALOKA, "Left Ventricle Mass for Pediatric"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode")	9
EDLepi	(A12201-003, 99ALOKA, "Left Ventricular End Diastolic Length")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12227-003, 99ALOKA, "Left Ventricle Mass for Pediatric"), Finding Site=(A12236-006, 99ALOKA, "Epicardium"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode")	9
EDAepi	(G-0379, SRT, "Left Ventricle Epicardial Diastolic Area, psax pap view")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12227-003, 99ALOKA, "Left Ventricle Mass for Pediatric"), Finding Site=(A12236-006, 99ALOKA, "Epicardium"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode")	9

Label	Concept Name	Modifier	row
EDV	(18026-5, LN, "Left Ventricular End Diastolic Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	9
EDVepi	(18026-5, LN, "Left Ventricular End Diastolic Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12227-003, 99ALOKA, "Left Ventricle Mass for Pediatric"), Finding Site=(A12236-006, 99ALOKA, "Epicardium"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	9
(LV Mass Pedi Simpson)			
LVM	(18087-7, LN, "Left Ventricle Mass")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12227-004, 99ALOKA, "Left Ventricle Mass for Pediatric by Modified Simpson")	9
LVM/BSA	(A12203-001, 99ALOKA, "Left Ventricular Mass Index")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12227-004, 99ALOKA, "Left Ventricle Mass for Pediatric by Modified Simpson")	9
EDV(A4C)	(18026-5, LN, "Left Ventricular End Diastolic Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125208, DCM, "Method of Disks, Single Plane"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image View=(G-A19C, SRT, "Apical four chamber")	9
EDV(A2C)	(18026-5, LN, "Left Ventricular End Diastolic Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-001, 99ALOKA, "Method of Disks, Single Plane with Apical two chamber"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image View=(G-A19B, SRT, "Apical two chamber")	9
EDV	(18026-5, LN, "Left Ventricular End Diastolic Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125207, DCM, "Method of Disks, Biplane"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	9
EDVepi(A4C)	(18026-5, LN, "Left Ventricular End Diastolic Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12227-004, 99ALOKA, "Left Ventricle Mass for Pediatric by Modified Simpson"), Finding Site=(A12236-006, 99ALOKA, "Epicardium"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image View=(G-A19C, SRT, "Apical four chamber")	9

Label	Concept Name	Modifier	row
EDVepi(A2C)	(18026-5, LN, "Left Ventricular End Diastolic Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12227-004, 99ALOKA, "Left Ventricle Mass for Pediatric by Modified Simpson"), Finding Site=(A12236-006, 99ALOKA, "Epicardium"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image View=(G-A19B, SRT, "Apical two chamber")	9
EDVepi	(18026-5, LN, "Left Ventricular End Diastolic Volume")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12227-004, 99ALOKA, "Left Ventricle Mass for Pediatric by Modified Simpson"), Finding Site=(A12236-006, 99ALOKA, "Epicardium"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	9
%difD	(A12203-003, 99ALOKA, "Long Axis(at End Diastole or End Systole) Length % Difference")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	9
%difDepi	(A12203-003, 99ALOKA, "Long Axis(at End Diastole or End Systole) Length % Difference")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12227-004, 99ALOKA, "Left Ventricle Mass for Pediatric by Modified Simpson"), Finding Site=(A12236-006, 99ALOKA, "Epicardium")	9
(RV(2D))			
RV WT	(A12204-007, 99ALOKA, "Right Ventricular Wall Thickness")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode")	10
RVDd base	(A12204-004, 99ALOKA, "End-diastolic Basal Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode")	10
RVDd mid	(A12204-005, 99ALOKA, "End-diastolic Mid-cavity Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode")	10
RVLd	(A12204-006, 99ALOKA, "Right Ventricular End-diastolic Length")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode")	10
RVOT	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract"), Image Mode=(G-03A2, SRT, "2D mode")	10
RVOT VTI	(20354-7, LN, "Velocity Time Integral")	Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	10

Label	Concept Name	Modifier	row
RVAWd	(18153-7, LN, "Right Ventricular Anterior Wall Diastolic Thickness")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	10
RVAWs	(18157-8, LN, "Right Ventricular Anterior Wall Systolic Thickness")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(109070, DCM, "End Systole")	10
RVDd lax	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(443083000, SCT, "Parasternal long axis view of right ventricular outflow tract")	10
RVDd prox	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(A12236-001, 99ALOKA, "Proximal Right Ventricle Outflow Tract"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(399195005, SCT, "Right Ventricular Outflow Tract View")	10
RVDd dist	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(A12236-002, 99ALOKA, "Distal Right Ventricle Outflow Tract"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(399195005, SCT, "Right Ventricular Outflow Tract View")	10
RVDs	(20305-9, LN, "Right Ventricular Internal Systolic Dimension")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(109070, DCM, "End Systole")	10
RV EDA	(A12204-001, 99ALOKA, "Right Ventricular Diastolic Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125220, DCM, "Planimetry"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	10
RV ESA	(A12204-002, 99ALOKA, "Right Ventricular Systolic Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125220, DCM, "Planimetry"), Cardiac Cycle Point=(R-FAB5B, SRT, "End Systole")	10
(AV(2D))			
LVOT	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
AV Ann diam	(A12211-004, 99ALOKA, "Aortic Annulus Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Flow Direction=(R-42047, SRT, "Antegrade Flow"), Image Mode=(G-03A2, SRT, "2D mode")	13

Label	Concept Name	Modifier	row
AVA(Trace)	(G-038E, SRT, "Cardiovascular Orifice Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125220, DCM, "Planimetry"), Flow Direction=(R-42047, SRT, "Antegrade Flow"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-0398, SRT, "Parasternal short axis at the aortic valve level")	13
AR VCW	(A12211-009, 99ALOKA, "Aortic Regurgitation Vena Contracta Width")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(G-03A2, SRT, "2D mode")	13
Jet A	(A12222-056, 99ALOKA, "Regurgitation Jet Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(G-03A2, SRT, "2D mode")	13
Jet H	(A12222-057, 99ALOKA, "Regurgitation Jet Height")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(G-03A2, SRT, "2D mode")	13
Jet W	(A12222-058, 99ALOKA, "Regurgitation Jet Width")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(G-03A2, SRT, "2D mode")	13
Jet W/LVOT	(A12211-008, 99ALOKA, "Aortic Regurgitant Jet Width/LVOT Width")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	13
Jet A/LVOT(CSA)	(A12211-007, 99ALOKA, "Aortic Regurgitant Jet Area/LVOT Area")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	13
(MV Annulus)			
Ann diam comm	(A12222-033, 99ALOKA, "Commissural Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(65197004, SCT, "Mitral Annulus"), Image Mode=(G-03A2, SRT, "2D mode")	14
Ann diam a-p	(A12222-022, 99ALOKA, "Anterior-Posterior Annulus Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(65197004, SCT, "Mitral Annulus"), Image Mode=(G-03A2, SRT, "2D mode")	14
Ann Area	(A12207-019, 99ALOKA, "Mitral Annulus Area")	Measurement Method=(125220, DCM, "Planimetry"), Finding Site=(65197004, SCT, "Mitral Annulus")	14
MVA(Trace)	(G-038E, SRT, "Cardiovascular Orifice Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125220, DCM, "Planimetry"), Flow Direction=(R-42047, SRT, "Antegrade Flow"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-039A, SRT, "Parasternal short axis at the Mitral Valve level")	14

Label	Concept Name	Modifier	row
MV Ann VTI	(20354-7, LN, "Velocity Time Integral")	Image Mode(Group)=(R-409E4, SRT, "Doppler Pulsed"), Measurement Method=(125215, DCM, "Continuity Equation by Velocity Time Integral"), Finding Site=(65197004, SCT, "Mitral Annulus"), Image Mode=(R-409E4, SRT, "Doppler Pulsed")	14
Ann Vmax e'	(399133000, SCT, "Left Ventricular Peak Early Diastolic Tissue Velocity")	Image Mode(Group)=(R-409E4, SRT, "Doppler Pulsed"), Image Mode=(R-409E4, SRT, "Doppler Pulsed")	9
Ann Vmax a'	(399007006, SCT, "LV Peak Diastolic Tissue Velocity During Atrial Systole")	Image Mode(Group)=(R-409E4, SRT, "Doppler Pulsed"), Cardiac Cycle Point=(R-FAB5B, SRT, "End Systole"), Image Mode=(R-409E4, SRT, "Doppler Pulsed")	9
MV Ann E/e'	(399140004, SCT, "Ratio of MV Peak Velocity to LV Peak Tissue Velocity E-Wave")		9
Ann PLAX	(A12207-020, 99ALOKA, "Mitral Annulus Diameter PLAX")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(65197004, SCT, "Mitral Annulus"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-0396, SRT, "Parasternal long axis")	14
(TV Annulus)			
Ann diam comm	(A12222-033, 99ALOKA, "Commissural Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(113259005, SCT, "Tricuspid Annulus"), Image Mode=(G-03A2, SRT, "2D mode")	16
Ann diam a-p	(A12222-022, 99ALOKA, "Anterior-Posterior Annulus Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(113259005, SCT, "Tricuspid Annulus"), Image Mode=(G-03A2, SRT, "2D mode")	16
Ann Area	(A12208-022, 99ALOKA, "Tricuspid Annulus Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125220, DCM, "Planimetry"), Finding Site=(113259005, SCT, "Tricuspid Annulus")	16
Ann Vmax e'	(A12208-017, 99ALOKA, "Right Ventricular Peak Early Diastolic Tissue Velocity")	Image Mode(Group)=(R-409E4, SRT, "Doppler Pulsed"), Finding Site=(113259005, SCT, "Tricuspid Annulus"), Image Mode=(R-409E4, SRT, "Doppler Pulsed")	16
Ann Vmax a'	(A12208-020, 99ALOKA, "RV Peak Diastolic Tissue Velocity During Atrial Systole")	Image Mode(Group)=(R-409E4, SRT, "Doppler Pulsed"), Finding Site=(113259005, SCT, "Tricuspid Annulus"), Image Mode=(R-409E4, SRT, "Doppler Pulsed")	16
TV Ann E/e'	(59111-5, LN, "E Velocity to Annulus E Velocity Ratio")	Finding Site=(113259005, SCT, "Tricuspid Annulus")	16

Label	Concept Name	Modifier	row
Ann PLAX	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(113259005, SCT, "Tricuspid Annulus"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-0396, SRT, "Parasternal long axis")	16
Ann PSAX	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(113259005, SCT, "Tricuspid Annulus"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-0397, SRT, "Parasternal short axis")	16
(PV Annulus)			
Ann diam a-p	(A12222-022, 99ALOKA, "Anterior-Posterior Annulus Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Image Mode=(G-03A2, SRT, "2D mode")	15
Ann diam trans	(A12209-009, 99ALOKA, "Transverse Annulus Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Image Mode=(G-03A2, SRT, "2D mode")	15
Ann Area	(A12209-016, 99ALOKA, "Purmonary Annulus Area")	Measurement Method=(125220, DCM, "Planimetry")	15
Ann PLAX	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-0396, SRT, "Parasternal long axis")	15
Ann PSAX	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-0397, SRT, "Parasternal short axis")	15
(TAVI/TAVR)			
Ann dmin	(A12211-006, 99ALOKA, "Aortic Annulus Minimum Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(A12236-016, 99ALOKA, "Aortic Valve Annulus"), Image Mode=(G-03A2, SRT, "2D mode")	13
Ann dmax	(A12211-005, 99ALOKA, "Aortic Annulus Maximum Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(A12236-016, 99ALOKA, "Aortic Valve Annulus"), Image Mode=(G-03A2, SRT, "2D mode")	13
Ann Area	(A12211-010, 99ALOKA, "Aortic Valve Annulus Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(A12236-016, 99ALOKA, "Aortic Valve Annulus"), Image Mode=(G-03A2, SRT, "2D mode")	13
Ao Ring	(A12211-004, 99ALOKA, "Aortic Annulus Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(T-35410, SRT, "Aortic Valve Ring"), Image Mode=(G-03A2, SRT, "2D mode")	13
Ao SV	(A12212-009, 99ALOKA, "Aortic Sinuses of Valsalva Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Image Mode=(G-03A2, SRT, "2D mode")	17

Label	Concept Name	Modifier	row
Ao STJ	(A12212-008, 99ALOKA, "Aortic Sinotubular Junction Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Image Mode=(G-03A2, SRT, "2D mode")	17
Ann-Ost Lt.	(A12222-020, 99ALOKA, "Annulus-Ostium Left diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Image Mode=(G-03A2, SRT, "2D mode")	17
Ann-Ost Rt.	(A12222-021, 99ALOKA, "Annulus-Ostium Right diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Image Mode=(G-03A2, SRT, "2D mode")	17
(MV(M))			
C-Eamp	(A12207-001, 99ALOKA, "Mitral Valve Dimension of C point to E point by M-Mode")	Image Mode(Group)=(G-0394, SRT, "M mode"), Image Mode=(G-0394, SRT, "M mode")	14
C-Aamp	(A12207-002, 99ALOKA, "Mitral Valve Dimension of C point to A point by M-Mode")	Image Mode(Group)=(G-0394, SRT, "M mode"), Image Mode=(G-0394, SRT, "M mode")	14
D-Eamp	(A12207-017, 99ALOKA, "D-E amplitude")	Image Mode(Group)=(G-0394, SRT, "M mode"), Image Mode=(G-0394, SRT, "M mode")	14
D-Eslop (F)	(59127-1, LN, "D-E Slope")	Image Mode(Group)=(G-0394, SRT, "M mode"), Image Mode=(G-0394, SRT, "M mode")	14
E-Fslop	(18040-6, LN, "Mitral Valve E-F Slope by M-Mode")	Image Mode(Group)=(G-0394, SRT, "M mode"), Image Mode=(G-0394, SRT, "M mode")	14
EPSS	(18036-4, LN, "Mitral Valve EPSS, E wave")	Image Mode(Group)=(G-0394, SRT, "M mode"), Image Mode=(G-0394, SRT, "M mode")	14
MAPSE	(A12207-018, 99ALOKA, "Mitral Annular Plane Systolic Excursion")	Image Mode(Group)=(G-0394, SRT, "M mode"), Finding Site=(65197004, SCT, "Mitral Annulus"), Image Mode=(G-0394, SRT, "M mode")	14
A/E	(A12207-003, 99ALOKA, "Mitral Valve C-A Dimension to C-E Dimension Ratio by M-Mode")	Image Mode(Group)=(G-0394, SRT, "M mode")	14
E/A	(A12207-004, 99ALOKA, "Mitral Valve C-E Dimension to C-A Dimension Ratio by M-Mode")	Image Mode(Group)=(G-0394, SRT, "M mode")	14
(RV(M))			
RV WT	(A12204-007, 99ALOKA, "Right Ventricular Wall Thickness")	Image Mode(Group)=(G-0394, SRT, "M mode"), Image Mode=(G-0394, SRT, "M mode")	10
TAPSE	(A12208-012, 99ALOKA, "Tricuspid Annular Plane Systolic Excursion")	Image Mode(Group)=(G-0394, SRT, "M mode")	16

Label	Concept Name	Modifier	row
RVDd	(20304-2, LN, "Right Ventricular Internal Diastolic Dimension")	Image Mode(Group)=(G-0394, SRT, "M mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	10
RVOT	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-0394, SRT, "M mode"), Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract"), Image Mode=(G-0394, SRT, "M mode")	10
(LVOT)			
VTI	(20354-7, LN, "Velocity Time Integral")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
pV	(11726-7, LN, "Peak Velocity")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
MnV	(20352-1, LN, "Mean Velocity")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
PG	(20247-3, LN, "Peak Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
MPG	(20256-4, LN, "Mean Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
LVOT	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
IVRT	(18071-1, LN, "Left Ventricular Isovolumic Relaxation Time")		9
LVET	(18041-4, LN, "Aortic Valve Ejection Time")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	13
PEP	(A12222-001, 99ALOKA, "Pre-Ejection Period")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
AccS	(20167-3, LN, "Acceleration Slope")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
AccT	(20168-1, LN, "Acceleration Time")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
PEP/ET	(A12222-002, 99ALOKA, "PEP/ET")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
AccT/ET	(G-0382, SRT, "Ratio of Aortic Valve Acceleration Time to Ejection Time")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	13
dPEP	(A12222-038, 99ALOKA, "Interventricular Dyssynchrony - Interventricular Delay")	Image Mode(Group)=(R-409E4, SRT, "Doppler Pulsed"), Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract"), Image Mode=(R-409E4, SRT, "Doppler Pulsed")	9
LVOT/AV(pV)	(A12203-014, 99ALOKA, "pV Ratio of LVOT to AV")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
LVOT/AV(VTI)	(A12203-015, 99ALOKA, "VTI Ratio of LVOT to AV")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9

Label	Concept Name	Modifier	row
SV	(F-32120, SRT, "Stroke Volume")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
CO	(F-32100, SRT, "Cardiac Output")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
SVi	(F-00078, SRT, "Stroke Index")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
COi	(F-32110, SRT, "Cardiac Index")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
AVA (VTI)	(G-038E, SRT, "Cardiovascular Orifice Area")	Measurement Method=(125215, DCM, "Continuity Equation by Velocity Time Integral"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
AVAi (VTI)	(79959-3, LN, "Aortic valve area (Continuity VTI) / BSA")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
AVA (pV)	(G-038E, SRT, "Cardiovascular Orifice Area")	Measurement Method=(125214, DCM, "Continuity Equation by Peak Velocity")	13
AVAi (pV)	(A12222-032, 99ALOKA, "Cardiovascular Orifice Area Index")	Measurement Method=(125214, DCM, "Continuity Equation by Peak Velocity")	13
HR	(8867-4, LN, "Heart rate")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
Qp/Qs	(29462-9, LN, "Pulmonary-to-Systemic Shunt Flow Ratio")		22
Qp-Qs	(A12217-001, 99ALOKA, "Pulmonary minus Systemic Shunt Flow - Doppler Volume Flow")		22
(AS Flow)			
pV	(11726-7, LN, "Peak Velocity")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
PG	(20247-3, LN, "Peak Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
MnV	(20352-1, LN, "Mean Velocity")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
MPG	(20256-4, LN, "Mean Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
VTI	(20354-7, LN, "Velocity Time Integral")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
LVOT	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
CSA	(G-038E, SRT, "Cardiovascular Orifice Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9

Label	Concept Name	Modifier	row
VTI (LVOT)	(20354-7, LN, "Velocity Time Integral")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
pV (LVOT)	(11726-7, LN, "Peak Velocity")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
PG (LVOT)	(20247-3, LN, "Peak Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
MnV (LVOT)	(20352-1, LN, "Mean Velocity")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
MPG (LVOT)	(20256-4, LN, "Mean Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
AVA (VTI)	(G-038E, SRT, "Cardiovascular Orifice Area")	Measurement Method=(125215, DCM, "Continuity Equation by Velocity Time Integral"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
AVAi (VTI)	(79959-3, LN, "Aortic valve area (Continuity VTI) / BSA")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
AVA (pV)	(G-038E, SRT, "Cardiovascular Orifice Area")	Measurement Method=(125214, DCM, "Continuity Equation by Peak Velocity")	13
AVAi (pV)	(A12222-032, 99ALOKA, "Cardiovascular Orifice Area Index")	Measurement Method=(125214, DCM, "Continuity Equation by Peak Velocity")	13
VR	(A12211-011, 99ALOKA, "Velocity Ratio of LVOT to AV")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
(AV)			
VTI	(20354-7, LN, "Velocity Time Integral")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
pV	(11726-7, LN, "Peak Velocity")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
MnV	(20352-1, LN, "Mean Velocity")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
PG	(20247-3, LN, "Peak Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
MPG	(20256-4, LN, "Mean Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
AccSlope	(20167-3, LN, "Acceleration Slope")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
AccT	(20168-1, LN, "Acceleration Time")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
LVET	(18041-4, LN, "Aortic Valve Ejection Time")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	13

Label	Concept Name	Modifier	row
AccT/ET	(G-0382, SRT, "Ratio of Aortic Valve Acceleration Time to Ejection Time")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
DecSlope	(20216-8, LN, "Deceleration Slope")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
DecT	(20217-6, LN, "Deceleration Time")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
LVOT	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
LVOT VTI	(20354-7, LN, "Velocity Time Integral")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
LVOT SV	(F-32120, SRT, "Stroke Volume")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
LVOT SVi	(F-00078, SRT, "Stroke Index")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
AVA(VTI)	(G-038E, SRT, "Cardiovascular Orifice Area")	Measurement Method=(125215, DCM, "Continuity Equation by Velocity Time Integral"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
AVAi (VTI)	(79959-3, LN, "Aortic valve area (Continuity VTI) / BSA")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	13
AVA (pV)	(G-038E, SRT, "Cardiovascular Orifice Area")	Measurement Method=(125214, DCM, "Continuity Equation by Peak Velocity")	13
AVAi (pV)	(A12222-032, 99ALOKA, "Cardiovascular Orifice Area Index")	Measurement Method=(125214, DCM, "Continuity Equation by Peak Velocity")	13
HR	(8867-4, LN, "Heart rate")	Image Mode(Group)=(R-409E4, SRT, "Doppler Pulsed"), Image Mode=(R-409E4, SRT, "Doppler Pulsed")	13
(AR)			
pV	(11726-7, LN, "Peak Velocity")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	13
PG	(20247-3, LN, "Peak Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	13
Ved	(11653-3, LN, "End Diastolic Velocity")	Image Mode(Group)=(R-409E3, SRT, "Doppler Continuous Wave"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(R-409E3, SRT, "Doppler Continuous Wave")	13

Label	Concept Name	Modifier	row
PGed	(A12222-034, 99ALOKA, "Diastolic Pressure Gradient")	Image Mode(Group)=(R-409E3, SRT, "Doppler Continuous Wave"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(R-409E3, SRT, "Doppler Continuous Wave")	13
DecSlope	(20216-8, LN, "Deceleration Slope")	Image Mode(Group)=(R-409E3, SRT, "Doppler Continuous Wave"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(R-409E3, SRT, "Doppler Continuous Wave")	13
DecT	(20217-6, LN, "Deceleration Time")	Image Mode(Group)=(R-409E3, SRT, "Doppler Continuous Wave"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(R-409E3, SRT, "Doppler Continuous Wave")	13
P1/2T	(20280-4, LN, "Pressure Half-Time")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	13
VTI(AR)	(A12211-001, 99ALOKA, "Velocity Time Integral of Aortic Regurgitant Flow")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	13
MnV	(20352-1, LN, "Mean Velocity")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	13
MPG	(20256-4, LN, "Mean Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	13
AccSlope	(20167-3, LN, "Acceleration Slope")	Image Mode(Group)=(R-409E3, SRT, "Doppler Continuous Wave"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(R-409E3, SRT, "Doppler Continuous Wave")	13
AccT	(20168-1, LN, "Acceleration Time")	Image Mode(Group)=(R-409E3, SRT, "Doppler Continuous Wave"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(R-409E3, SRT, "Doppler Continuous Wave")	13
LVOT	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
LVOT VTI	(20354-7, LN, "Velocity Time Integral")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
LVOT SV	(F-32120, SRT, "Stroke Volume")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
LVOT SVi	(F-00078, SRT, "Stroke Index")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
MV SV	(A12207-013, 99ALOKA, "Flow Volume of Mitral Valve Annulus in Flow")	Measurement Method=(125219, DCM, "Doppler Volume Flow")	14

Label	Concept Name	Modifier	row
Reg.Vol	(33878-0, LN, "Volume Flow")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	13
Reg.Frac	(G-0390, SRT, "Regurgitant Fraction")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	13
Jet A	(A12222-056, 99ALOKA, "Regurgitation Jet Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(G-03A2, SRT, "2D mode")	13
JetH	(A12222-057, 99ALOKA, "Regurgitation Jet Height")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(G-03A2, SRT, "2D mode")	13
JetW	(A12222-058, 99ALOKA, "Regurgitation Jet Width")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(G-03A2, SRT, "2D mode")	13
(MV)			
eV	(18037-2, LN, "Mitral Valve E-Wave Peak Velocity")		14
aV	(17978-8, LN, "Mitral Valve A-Wave Peak Velocity")		14
E/A	(18038-0, LN, "Mitral Valve E to A Ratio")		14
A/E	(A12207-006, 99ALOKA, "Mitral Valve A to E Ratio")		14
E/e'	(A12203-011, 99ALOKA, "Ratio of MV E-Wave Peak Vel. to Early Diastolic Myocardium Vel.")		9
TE-e'	(A12207-024, 99ALOKA, "Time E to e'")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(G-0391, SRT, "Medial Mitral Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
DecS	(20216-8, LN, "Deceleration Slope")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	14
DecT	(G-0384, SRT, "Mitral Valve E-Wave Deceleration Time")		14
Edur	(A12207-005, 99ALOKA, "Mitral Valve E-Wave Duration")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	14
Adur	(G-0385, SRT, "Mitral Valve A-Wave Duration")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	14
IVRT	(18071-1, LN, "Left Ventricular Isovolumic Relaxation Time")		9

Label	Concept Name	Modifier	row
MV C-Odur	(399104001, SCT, "Mitral Valve Closure to Opening Time")		14
LVET	(18041-4, LN, "Aortic Valve Ejection Time")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	13
VTI	(20354-7, LN, "Velocity Time Integral")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	14
pV	(11726-7, LN, "Peak Velocity")	Image Mode(Group)=(R-409E4, SRT, "Doppler Pulsed"), Flow Direction=(R-42047, SRT, "Antegrade Flow"), Image Mode=(R-409E4, SRT, "Doppler Pulsed")	14
MnV	(20352-1, LN, "Mean Velocity")	Image Mode(Group)=(R-409E4, SRT, "Doppler Pulsed"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	14
PG	(20247-3, LN, "Peak Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli")	14
MPG	(20256-4, LN, "Mean Gradient")	Image Mode(Group)=(R-409E4, SRT, "Doppler Pulsed"), Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	14
P1/2T	(20280-4, LN, "Pressure Half-Time")	Image Mode(Group)=(R-409E4, SRT, "Doppler Pulsed"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	14
MVA (P1/2T)	(G-038E, SRT, "Cardiovascular Orifice Area")	Image Mode(Group)=(R-409E4, SRT, "Doppler Pulsed"), Measurement Method=(125210, DCM, "Area by Pressure Half-Time"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	14
MVA (VTI)	(G-038E, SRT, "Cardiovascular Orifice Area")	Image Mode(Group)=(R-409E4, SRT, "Doppler Pulsed"), Measurement Method=(125215, DCM, "Continuity Equation by Velocity Time Integral")	14
Diam.A4C	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Finding Site=(T-35313, SRT, "Mitral Annulus"), Flow Direction=(R-42047, SRT, "Antegrade Flow"), Image View=(G-A19C, SRT, "Apical four chamber")	14
Diam.A2C	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Finding Site=(T-35313, SRT, "Mitral Annulus"), Flow Direction=(R-42047, SRT, "Antegrade Flow"), Image View=(G-A19B, SRT, "Apical two chamber")	14
SV (MV)	(A12207-026, 99ALOKA, "Flow Volume of Mitral Valve Annulus in Flow (1)")	Measurement Method=(125219, DCM, "Doppler Volume Flow")	14
CO (MV)	(F-32100, SRT, "Cardiac Output")	Measurement Method=(125207, DCM, "Method of Disks, Biplane")	14
SVi (MV)	(F-00078, SRT, "Stroke Index")	Measurement Method=(125207, DCM, "Method of Disks, Biplane")	14

Label	Concept Name	Modifier	row
COi (MV)	(F-32110, SRT, "Cardiac Index")	Measurement Method=(125207, DCM, "Method of Disks, Biplane")	14
Acc Slope	(20167-3, LN, "Acceleration Slope")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	14
AccT	(20168-1, LN, "Acceleration Time")		14
PVAdur-Adur	(A12207-009, 99ALOKA, "Subtraction of A-wave Duration from PVA-wave Duration")		14
HR	(8867-4, LN, "Heart rate")	Image Mode(Group)=(R-409E4, SRT, "Doppler Pulsed"), Image Mode=(R-409E4, SRT, "Doppler Pulsed")	14
(MV(Valsalva))			
eV(Val)	(18037-2, LN, "Mitral Valve E-Wave Peak Velocity")	Measurement Method=(A12227-001, 99ALOKA, "Valsalva")	14
aV(Val)	(17978-8, LN, "Mitral Valve A-Wave Peak Velocity")	Measurement Method=(A12227-001, 99ALOKA, "Valsalva")	14
E/A(Val)	(18038-0, LN, "Mitral Valve E to A Ratio")	Measurement Method=(A12227-001, 99ALOKA, "Valsalva")	14
(MR Flow)			
VTI(MR)	(A12207-010, 99ALOKA, "Velocity Time Integral of Mitral Regurgitant Flow")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	14
pV	(11726-7, LN, "Peak Velocity")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	14
MnV	(20352-1, LN, "Mean Velocity")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	14
PG	(20247-3, LN, "Peak Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	14
MPG	(20256-4, LN, "Mean Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	14
FlowT	(A12222-003, 99ALOKA, "Flow Time")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	14
dP/dt	(18035-6, LN, "Mitral Regurgitation dP/dt derived from Mitral Reg. velocity")		14
dt	(A12222-059, 99ALOKA, "Time")	Image Mode(Group)=(R-409E3, SRT, "Doppler Continuous Wave"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(R-409E3, SRT, "Doppler Continuous Wave")	14

Label	Concept Name	Modifier	row
PG1	(A12222-010, 99ALOKA, "Peak Gradient 1")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	14
V1	(A12222-024, 99ALOKA, "Blood Flow Velocity 1")	Image Mode(Group)=(R-409E3, SRT, "Doppler Continuous Wave"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(R-409E3, SRT, "Doppler Continuous Wave")	14
PG2	(A12222-011, 99ALOKA, "Peak Gradient 2")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	14
V2	(A12222-027, 99ALOKA, "Blood Flow Velocity 2")	Image Mode(Group)=(R-409E3, SRT, "Doppler Continuous Wave"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(R-409E3, SRT, "Doppler Continuous Wave")	14
MR VCV	(A12222-007, 99ALOKA, "Vena Contracta Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode")	14
Jet A	(A12222-056, 99ALOKA, "Regurgitation Jet Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(G-03A2, SRT, "2D mode")	14
Jet H	(A12222-057, 99ALOKA, "Regurgitation Jet Height")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(G-03A2, SRT, "2D mode")	14
Jet W	(A12222-058, 99ALOKA, "Regurgitation Jet Width")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(G-03A2, SRT, "2D mode")	14
VTI(MV annu)	(20354-7, LN, "Velocity Time Integral")	Finding Site=(T-35313, SRT, "Mitral Annulus"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	14
Diam.A4C	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Finding Site=(T-35313, SRT, "Mitral Annulus"), Flow Direction=(R-42047, SRT, "Antegrade Flow"), Image View=(G-A19C, SRT, "Apical four chamber")	14
Diam.A2C	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Finding Site=(T-35313, SRT, "Mitral Annulus"), Flow Direction=(R-42047, SRT, "Antegrade Flow"), Image View=(G-A19B, SRT, "Apical two chamber")	14
MV SV	(A12207-013, 99ALOKA, "Flow Volume of Mitral Valve Annulus in Flow")	Measurement Method=(125219, DCM, "Doppler Volume Flow")	14
MV Leaf Sep	(59109-9, LN, "Leaflet Separation")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(G-03A2, SRT, "2D mode")	14

Label	Concept Name	Modifier	row
Reg.Vol	(A12207-014, 99ALOKA, "Mitral valve Regurgitant volume")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	14
Reg.Frac	(G-0390, SRT, "Regurgitant Fraction")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	14
EROA	(G-038E, SRT, "Cardiovascular Orifice Area")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	14
(PVA)			
VTI	(20354-7, LN, "Velocity Time Integral")		15
pV	(11726-7, LN, "Peak Velocity")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	15
PG	(20247-3, LN, "Peak Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	15
MnV	(20352-1, LN, "Mean Velocity")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	15
MPG	(20256-4, LN, "Mean Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	15
RVOT VTI	(20354-7, LN, "Velocity Time Integral")	Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	10
RVOT pV	(11726-7, LN, "Peak Velocity")	Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	10
RVOT PG	(20247-3, LN, "Peak Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	10
RVOT MnV	(20352-1, LN, "Mean Velocity")	Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	10
RVOT MPG	(20256-4, LN, "Mean Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	10
RVOT	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract"), Image Mode=(G-03A2, SRT, "2D mode")	10
PVA(VTI)	(G-038E, SRT, "Cardiovascular Orifice Area")	Measurement Method=(125215, DCM, "Continuity Equation by Velocity Time Integral")	15
PVAi(VTI)	(A12222-032, 99ALOKA, "Cardiovascular Orifice Area Index")	Measurement Method=(125215, DCM, "Continuity Equation by Velocity Time Integral")	15
PVA(pV)	(G-038E, SRT, "Cardiovascular Orifice Area")	Measurement Method=(125214, DCM, "Continuity Equation by Peak Velocity")	15

Label	Concept Name	Modifier	row
PVAi(pV)	(A12222-032, 99ALOKA, "Cardiovascular Orifice Area Index")	Measurement Method=(125214, DCM, "Continuity Equation by Peak Velocity")	15
CSA	(G-038E, SRT, "Cardiovascular Orifice Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	10
VR	(A12209-017, 99ALOKA, "Velocity Ratio of RVOT velocity to Pulmonic Valve velocity")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	15
(PV)			
pV	(11726-7, LN, "Peak Velocity")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	15
VTI	(20354-7, LN, "Velocity Time Integral")		15
MnV	(20352-1, LN, "Mean Velocity")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	15
PG	(20247-3, LN, "Peak Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	15
MPG	(20256-4, LN, "Mean Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	15
ET	(18042-2, LN, "Pulmonic Valve Ejection Time")	Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	15
PEP	(A12222-001, 99ALOKA, "Pre-Ejection Period")	Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	10
PEP/ET	(A12222-002, 99ALOKA, "PEP/ET")	Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	10
RV AccSlope	(20167-3, LN, "Acceleration Slope")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	15
RV AccT	(20168-1, LN, "Acceleration Time")		15
RV AccT/ET	(G-0388, SRT, "Ratio of Pulmonic Valve Acceleration Time to Ejection Time")		15
PA AccSlope	(20167-3, LN, "Acceleration Slope")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	18
PA AccT	(20168-1, LN, "Acceleration Time")		18
PA P1/2T	(20280-4, LN, "Pressure Half-Time")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	18
Q to PV close	(20295-2, LN, "Time from Q wave to Pulmonic Valve Closes")		15

Label	Concept Name	Modifier	row
PV diam	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Image Mode=(G-03A2, SRT, "2D mode")	15
RVOT	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract"), Image Mode=(G-03A2, SRT, "2D mode")	10
dPEP	(A12222-038, 99ALOKA, "Interventricular Dyssynchrony - Interventricular Delay")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	15
(PR)			
pV	(11726-7, LN, "Peak Velocity")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	15
PG	(20247-3, LN, "Peak Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	15
MnV	(20352-1, LN, "Mean Velocity")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	15
MPG	(20256-4, LN, "Mean Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	15
VTI(PR)	(A12209-005, 99ALOKA, "Velocity Time Integral of Pulmonic Regurgitant Flow")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	15
PRend pV	(11653-3, LN, "End Diastolic Velocity")	Image Mode(Group)=(R-409E3, SRT, "Doppler Continuous Wave"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(R-409E3, SRT, "Doppler Continuous Wave")	15
PRend PG	(A12222-034, 99ALOKA, "Diastolic Pressure Gradient")	Image Mode(Group)=(R-409E3, SRT, "Doppler Continuous Wave"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	15
RAP	(18070-3, LN, "Right Atrium Systolic Pressure")		12
PADP	(A12209-013, 99ALOKA, "Pulmonary Artery Pressure")	Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	15
sPAP	(A12209-013, 99ALOKA, "Pulmonary Artery Pressure")	Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Cardiac Cycle Point=(R-FAB5B, SRT, "End Systole")	15
mPAP	(A12209-013, 99ALOKA, "Pulmonary Artery Pressure")	Measurement Method=(A12227-002, 99ALOKA, "Mean Pulmonary Arterial Pressure"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	15

Label	Concept Name	Modifier	row
DecSlope	(20216-8, LN, "Deceleration Slope")	Image Mode(Group)=(R-409E3, SRT, "Doppler Continuous Wave"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(R-409E3, SRT, "Doppler Continuous Wave")	15
DecT	(20217-6, LN, "Deceleration Time")	Image Mode(Group)=(R-409E3, SRT, "Doppler Continuous Wave"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(R-409E3, SRT, "Doppler Continuous Wave")	15
P1/2T	(20280-4, LN, "Pressure Half-Time")	Image Mode(Group)=(R-409E3, SRT, "Doppler Continuous Wave"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(R-409E3, SRT, "Doppler Continuous Wave")	15
PR dur	(A12209-014, 99ALOKA, "Pulmonary Valve Regurgitation Duration")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	15
dur diastole	(A12209-010, 99ALOKA, "Duration of Diastole")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	15
PR index	(A12209-011, 99ALOKA, "Ratio of PR dur to PR Index")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	15
PR VCW	(A12209-015, 99ALOKA, "Pulmonary Valve Regurgitation Vena Contracta Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(G-03A2, SRT, "2D mode")	15
Jet A	(A12222-056, 99ALOKA, "Regurgitation Jet Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(G-03A2, SRT, "2D mode")	15
Jet H	(A12222-057, 99ALOKA, "Regurgitation Jet Height")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(G-03A2, SRT, "2D mode")	15
Jet W	(A12222-058, 99ALOKA, "Regurgitation Jet Width")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(G-03A2, SRT, "2D mode")	15
PR VCW/Ann diam	(A12209-012, 99ALOKA, "Ratio of PR VCW to Ann diam a-p")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	15
RVOT	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract"), Image Mode=(G-03A2, SRT, "2D mode")	10
RVOT VTI	(20354-7, LN, "Velocity Time Integral")	Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	10
RVOT SV	(F-32120, SRT, "Stroke Volume")	Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	10

Label	Concept Name	Modifier	row
LVOT	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
LVOT VTI	(20354-7, LN, "Velocity Time Integral")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
LVOT SV	(F-32120, SRT, "Stroke Volume")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9
Reg.Vol	(33878-0, LN, "Volume Flow")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	15
Reg.Frac	(G-0390, SRT, "Regurgitant Fraction")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	15
(TV)			
eV	(59080-2, LN, "E-Wave Peak Velocity")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
EPG	(A12208-015, 99ALOKA, "E-Wave Peak Gradient")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
aV	(59081-0, LN, "A-Wave Peak Velocity")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
APG	(A12222-023, 99ALOKA, "A-Wave Peak Gradient")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
E/A	(59104-0, LN, "Peak E wave/Peak A wave by US")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
A/E	(A12208-013, 99ALOKA, "Peak A wave/Peak E wave by US")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
E/e'(l)	(59111-5, LN, "E Velocity to Annulus E Velocity Ratio")	Finding Site=(A12236-005, 99ALOKA, "Lateral Tricuspid Annulus"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
DecSlope	(20216-8, LN, "Deceleration Slope")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
DecT	(20217-6, LN, "Deceleration Time")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
Edur	(A12208-014, 99ALOKA, "E-Wave Duration")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
Adur	(59105-7, LN, "A-Wave Duration")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
RVET	(18042-2, LN, "Pulmonic Valve Ejection Time")	Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	15
VTI	(20354-7, LN, "Velocity Time Integral")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
pV	(11726-7, LN, "Peak Velocity")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
MnV	(20352-1, LN, "Mean Velocity")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	16

Label	Concept Name	Modifier	row
PG	(20247-3, LN, "Peak Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
MPG	(20256-4, LN, "Mean Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
P1/2T	(20280-4, LN, "Pressure Half-Time")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
TVA(P1/2T)	(G-038E, SRT, "Cardiovascular Orifice Area")	Measurement Method=(125210, DCM, "Area by Pressure Half-Time"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
Q to TV open	(20296-0, LN, "Time from Q wave to Tricuspid Valve Opens")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
AccSlope	(20167-3, LN, "Acceleration Slope")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
AccT	(20168-1, LN, "Acceleration Time")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
IVRT	(A12222-042, 99ALOKA, "Isovolumic Relaxation Time")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
FlowT	(A12222-003, 99ALOKA, "Flow Time")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
(TV Valsalva)			
eV(Val)	(59080-2, LN, "E-Wave Peak Velocity")	Measurement Method=(A12227-001, 99ALOKA, "Valsalva")	16
aV(Val)	(59081-0, LN, "A-Wave Peak Velocity")	Measurement Method=(A12227-001, 99ALOKA, "Valsalva")	16
E/A(Val)	(59104-0, LN, "Peak E wave/Peak A wave by US")	Measurement Method=(A12227-001, 99ALOKA, "Valsalva")	16
(TVA(VTI))			
TV VTI	(20354-7, LN, "Velocity Time Integral")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
RVOT VTI	(20354-7, LN, "Velocity Time Integral")	Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	10
RVOT diam	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract"), Image Mode=(G-03A2, SRT, "2D mode")	10
TVA(VTI)	(G-038E, SRT, "Cardiovascular Orifice Area")	Measurement Method=(125215, DCM, "Continuity Equation by Velocity Time Integral")	16
TVAi(VTI)	(A12222-032, 99ALOKA, "Cardiovascular Orifice Area Index")	Measurement Method=(125215, DCM, "Continuity Equation by Velocity Time Integral")	16
(TVA(pV))			

Label	Concept Name	Modifier	row
TV pV	(11726-7, LN, "Peak Velocity")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
RVOT pV	(11726-7, LN, "Peak Velocity")	Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	10
RVOT diam	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract"), Image Mode=(G-03A2, SRT, "2D mode")	10
TVA (pV)	(G-038E, SRT, "Cardiovascular Orifice Area")	Measurement Method=(125214, DCM, "Continuity Equation by Peak Velocity")	16
TVAi(pV)	(A12222-032, 99ALOKA, "Cardiovascular Orifice Area Index")	Measurement Method=(125214, DCM, "Continuity Equation by Peak Velocity")	16
(TR Flow)			
pV	(11726-7, LN, "Peak Velocity")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	16
PG	(20247-3, LN, "Peak Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	16
MnV	(20352-1, LN, "Mean Velocity")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	16
MPG	(20256-4, LN, "Mean Gradient")	Measurement Method=(125218, DCM, "Simplified Bernoulli"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	16
VTI	(A12208-008, 99ALOKA, "Velocity Time Integral of Tricuspid Regurgitant Flow")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	16
dP/dt	(18034-9, LN, "Tricuspid Regurgitation dP/dt")		16
PG1	(A12222-012, 99ALOKA, "Peak Pressure Gradient 1")	Image Mode(Group)=(R-409E3, SRT, "Doppler Continuous Wave"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	16
V1	(A12222-030, 99ALOKA, "Blood Flow Velocity at Peak Pressure Gradient 1")	Image Mode(Group)=(R-409E3, SRT, "Doppler Continuous Wave"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(R-409E3, SRT, "Doppler Continuous Wave")	16
PG2	(A12222-013, 99ALOKA, "Peak Pressure Gradient 2")	Image Mode(Group)=(R-409E3, SRT, "Doppler Continuous Wave"), Flow Direction=(R-42E61, SRT, "Retrograde Flow")	16
V2	(A12222-031, 99ALOKA, "Blood Flow Velocity at Peak Pressure Gradient 2")	Image Mode(Group)=(R-409E3, SRT, "Doppler Continuous Wave"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(R-409E3, SRT, "Doppler Continuous Wave")	16

Label	Concept Name	Modifier	row
dt	(A12208-021, 99ALOKA, "Time Difference between v1 and v2")	Image Mode(Group)=(R-409E3, SRT, "Doppler Continuous Wave"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(R-409E3, SRT, "Doppler Continuous Wave")	16
FlowT	(A12222-003, 99ALOKA, "Flow Time")	Flow Direction=(R-42E61, SRT, "Retrograde Flow")	16
RVSP	(G-0380, SRT, "Right Ventricular Peak Systolic Pressure")	Cardiac Cycle Point=(F-32020, SRT, "Systole")	10
RAP	(18070-3, LN, "Right Atrium Systolic Pressure")		12
TR VCW	(A12222-007, 99ALOKA, "Vena Contracta Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode")	16
Jet A	(A12222-056, 99ALOKA, "Regurgitation Jet Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(G-03A2, SRT, "2D mode")	16
Jet H	(A12222-057, 99ALOKA, "Regurgitation Jet Height")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(G-03A2, SRT, "2D mode")	16
Jet W	(A12222-058, 99ALOKA, "Regurgitation Jet Width")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Flow Direction=(R-42E61, SRT, "Retrograde Flow"), Image Mode=(G-03A2, SRT, "2D mode")	16
(M Annulus(m))			
e'(m)	(78185-6, LN, "MV med ann Peak e' wave US.TDI")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
a'(m)	(A12222-045, 99ALOKA, "Peak Diastolic Tissue Velocity During Atrial Systole")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
e'/a'(m)	(A12222-035, 99ALOKA, "Ea to Aa Ratio")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
E/e'(m)	(A12203-011, 99ALOKA, "Ratio of MV E-Wave Peak Vel. to Early Diastolic Myocardium Vel.")		9
s'(m)	(A12222-048, 99ALOKA, "Peak Systolic Tissue Velocity 1")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
ACC	(20167-3, LN, "Acceleration Slope")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging")	14

Label	Concept Name	Modifier	row
AccT	(20168-1, LN, "Acceleration Time")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
DEC	(20216-8, LN, "Deceleration Slope")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(G-0391, SRT, "Medial Mitral Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
DecT	(20217-6, LN, "Deceleration Time")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(G-0391, SRT, "Medial Mitral Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
IVRT	(A12222-042, 99ALOKA, "Isovolumic Relaxation Time")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
IVCT	(A12222-039, 99ALOKA, "Isovolumic Contraction Time")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
TE-e'(m)	(A12207-024, 99ALOKA, "Time E to e'")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(G-0391, SRT, "Medial Mitral Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
IVC Vel(m)	(A12207-022, 99ALOKA, "Peak MV Annular Velocity in Isovolumic Contraction")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(G-0391, SRT, "Medial Mitral Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
IVC TPV(m)	(A12207-025, 99ALOKA, "Time to peak velocity in isovolumic contractuion")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(G-0391, SRT, "Medial Mitral Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
IVA(m)	(A12207-015, 99ALOKA, "Isovolumic Acceleration")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(G-0391, SRT, "Medial Mitral Annulus")	14
(M Annulus(l))			
e'(l)	(G-037A, SRT, "Left Ventricular Peak Early Diastolic Tissue Velocity")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(G-0392, SRT, "Lateral Mitral Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9

Label	Concept Name	Modifier	row
a'(l)	(G-037C, SRT, "LV Peak Diastolic Tissue Velocity During Atrial Systole")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(G-0392, SRT, "Lateral Mitral Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
e'/a'(l)	(A12207-021, 99ALOKA, "Mitral Valve Ea to Aa Ratio")	Finding Site=(G-0392, SRT, "Lateral Mitral Annulus")	14
E/e'(l)	(G-037B, SRT, "Ratio of MV Peak Velocity to LV Peak Tissue Velocity E-Wave")	Finding Site=(G-0392, SRT, "Lateral Mitral Annulus")	9
av e'	(A12207-016, 99ALOKA, "Average of e'(l) and e'(s)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
E/av e'	(A12207-023, 99ALOKA, "Ratio of MV Peak Velocity to LV av e'")		14
s'(l)	(G-037D, SRT, "Left Ventricular Peak Systolic Tissue Velocity")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(G-0392, SRT, "Lateral Mitral Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
ACC	(20167-3, LN, "Acceleration Slope")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(G-0392, SRT, "Lateral Mitral Annulus")	14
AccT	(20168-1, LN, "Acceleration Time")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(G-0392, SRT, "Lateral Mitral Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
DEC	(20216-8, LN, "Deceleration Slope")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(G-0392, SRT, "Lateral Mitral Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
DecT	(20217-6, LN, "Deceleration Time")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(G-0392, SRT, "Lateral Mitral Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
IVRT	(A12222-042, 99ALOKA, "Isovolumic Relaxation Time")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(G-0392, SRT, "Lateral Mitral Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
IVCT	(A12222-039, 99ALOKA, "Isovolumic Contraction Time")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(G-0392, SRT, "Lateral Mitral Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14

Label	Concept Name	Modifier	row
TE-e'(l)	(A12207-024, 99ALOKA, "Time E to e'")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(G-0392, SRT, "Lateral Mitral Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
IVC Vel(l)	(A12207-022, 99ALOKA, "Peak MV Annular Velocity in Isovolumic Contraction")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(G-0392, SRT, "Lateral Mitral Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
IVC TPV(l)	(A12207-025, 99ALOKA, "Time to peak velocity in isovolumic contractuion")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(G-0392, SRT, "Lateral Mitral Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
IVA(l)	(A12207-015, 99ALOKA, "Isovolumic Acceleration")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(G-0392, SRT, "Lateral Mitral Annulus")	14
(M Annulus(s))			
e'(s)	(G-037A, SRT, "Left Ventricular Peak Early Diastolic Tissue Velocity")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(A12236-003, 99ALOKA, "Septal Mitral Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
a'(s)	(G-037C, SRT, "LV Peak Diastolic Tissue Velocity During Atrial Systole")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(A12236-003, 99ALOKA, "Septal Mitral Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
e'/a'(s)	(A12207-021, 99ALOKA, "Mitral Valve Ea to Aa Ratio")	Finding Site=(A12236-003, 99ALOKA, "Septal Mitral Annulus")	14
E/e'(s)	(G-037B, SRT, "Ratio of MV Peak Velocity to LV Peak Tissue Velocity E-Wave")	Finding Site=(A12236-003, 99ALOKA, "Septal Mitral Annulus")	9
av e'	(A12207-016, 99ALOKA, "Average of e'(l) and e'(s)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
E/av e'	(A12207-023, 99ALOKA, "Ratio of MV Peak Velocity to LV av e'")		14
s'(s)	(G-037D, SRT, "Left Ventricular Peak Systolic Tissue Velocity")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(A12236-003, 99ALOKA, "Septal Mitral Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
ACC	(20167-3, LN, "Acceleration Slope")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(A12236-003, 99ALOKA, "Septal Mitral Annulus")	14

Label	Concept Name	Modifier	row
AccT	(20168-1, LN, "Acceleration Time")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(A12236-003, 99ALOKA, "Septal Mitral Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
DEC	(20216-8, LN, "Deceleration Slope")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(A12236-003, 99ALOKA, "Septal Mitral Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
DecT	(20217-6, LN, "Deceleration Time")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(A12236-003, 99ALOKA, "Septal Mitral Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
IVRT	(A12222-042, 99ALOKA, "Isovolumic Relaxation Time")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(A12236-003, 99ALOKA, "Septal Mitral Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
IVCT	(A12222-039, 99ALOKA, "Isovolumic Contraction Time")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(A12236-003, 99ALOKA, "Septal Mitral Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
TE-e'(s)	(A12207-024, 99ALOKA, "Time E to e'")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(A12236-003, 99ALOKA, "Septal Mitral Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
IVC Vel(s)	(A12207-022, 99ALOKA, "Peak MV Annular Velocity in Isovolumic Contraction")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(A12236-003, 99ALOKA, "Septal Mitral Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
IVC TPV(s)	(A12207-025, 99ALOKA, "Time to peak velocity in isovolumic contractuion")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(A12236-003, 99ALOKA, "Septal Mitral Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
IVA(s)	(A12207-015, 99ALOKA, "Isovolumic Acceleration")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(A12236-003, 99ALOKA, "Septal Mitral Annulus")	14
(T Annulus(s))			
TV e'(s)	(A12208-017, 99ALOKA, "Right Ventricular Peak Early Diastolic Tissue Velocity")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(A12236-004, 99ALOKA, "Septal Tricuspid Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	16

Label	Concept Name	Modifier	row
TV a'(s)	(A12208-020, 99ALOKA, "RV Peak Diastolic Tissue Velocity During Atrial Systole")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(A12236-004, 99ALOKA, "Septal Tricuspid Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	16
TV s'(s)	(A12208-018, 99ALOKA, "Right Ventricular Peak Systolic Tissue Velocity")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(A12236-004, 99ALOKA, "Septal Tricuspid Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	16
TV E/e'(s)	(A12208-016, 99ALOKA, "Ratio of TV Peak Velocity to RV Peak Tissue Velocity E-Wave")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(A12236-004, 99ALOKA, "Septal Tricuspid Annulus")	16
(T Annulus(l))			
TV e'(l)	(A12208-017, 99ALOKA, "Right Ventricular Peak Early Diastolic Tissue Velocity")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(A12236-005, 99ALOKA, "Lateral Tricuspid Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	16
TV a'(l)	(A12208-020, 99ALOKA, "RV Peak Diastolic Tissue Velocity During Atrial Systole")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(A12236-005, 99ALOKA, "Lateral Tricuspid Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	16
TV s'(l)	(A12208-018, 99ALOKA, "Right Ventricular Peak Systolic Tissue Velocity")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Finding Site=(A12236-005, 99ALOKA, "Lateral Tricuspid Annulus"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	16
TV E/e'(l)	(59111-5, LN, "E Velocity to Annulus E Velocity Ratio")	Finding Site=(A12236-005, 99ALOKA, "Lateral Tricuspid Annulus"), Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
(LV MPI)			
ET	(18041-4, LN, "Aortic Valve Ejection Time")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	13
MV C-Odur	(399104001, SCT, "Mitral Valve Closure to Opening Time")		14
IVCT	(399051002, SCT, "Left Ventricular Isovolumic Contraction Time")	Image Mode(Group)=(R-409E4, SRT, "Doppler Pulsed"), Image Mode=(R-409E4, SRT, "Doppler Pulsed")	9
IVRT	(18071-1, LN, "Left Ventricular Isovolumic Relaxation Time")		9
LV MPI	(399266005, SCT, "Left Ventricular Index of Myocardial Performance")	Finding Site=(T-32650, SRT, "Left Ventricle Outflow Tract")	9

Label	Concept Name	Modifier	row
(RV MPI)			
RVET	(18042-2, LN, "Pulmonic Valve Ejection Time")	Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	15
TV C-Odur	(A12208-023, 99ALOKA, "Tricuspid Valve Closure to Opening Time")	Image Mode(Group)=(R-409E4, SRT, "Doppler Pulsed"), Image Mode=(R-409E4, SRT, "Doppler Pulsed")	16
IVCT	(A12222-039, 99ALOKA, "Isovolumic Contraction Time")	Image Mode(Group)=(R-409E4, SRT, "Doppler Pulsed"), Image Mode=(R-409E4, SRT, "Doppler Pulsed")	10
IVRT	(A12222-042, 99ALOKA, "Isovolumic Relaxation Time")	Flow Direction=(R-42047, SRT, "Antegrade Flow")	16
RV MPI ⁴	(399154007, SCT, "Right Ventricular Index of Myocardial Performance")	Finding Site=(T-32550, SNM3, "Right Ventricle Outflow Tract")	10
(RV MPI(Tissue))			
ET	(18042-2, LN, "Pulmonic Valve Ejection Time")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	10
TV C-Odur	(A12208-023, 99ALOKA, "Tricuspid Valve Closure to Opening Time")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	10
IVCT	(A12222-039, 99ALOKA, "Isovolumic Contraction Time")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	10
IVRT	(A12222-042, 99ALOKA, "Isovolumic Relaxation Time")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	10
RV MPI ⁴	(399154007, SCT, "Right Ventricular Index of Myocardial Performance")		10
(Ascending Ao Doppler)			
pV	(11726-7, LN, "Peak Velocity")	Finding Site=(A12236-007, 99ALOKA, "Ascending Aorta")	17
VTI	(20354-7, LN, "Velocity Time Integral")	Finding Site=(A12236-007, 99ALOKA, "Ascending Aorta")	17
MnV	(20352-1, LN, "Mean Velocity")	Finding Site=(A12236-007, 99ALOKA, "Ascending Aorta")	17
PG	(20247-3, LN, "Peak Gradient")	Finding Site=(A12236-007, 99ALOKA, "Ascending Aorta")	17
MPG	(20256-4, LN, "Mean Gradient")	Finding Site=(A12236-007, 99ALOKA, "Ascending Aorta")	17

Label	Concept Name	Modifier	row
(Descending Ao Doppler)			
pV	(11726-7, LN, "Peak Velocity")	Finding Site=(A12236-008, 99ALOKA, "Descending Aorta")	17
VTI	(20354-7, LN, "Velocity Time Integral")	Finding Site=(A12236-008, 99ALOKA, "Descending Aorta")	17
MnV	(20352-1, LN, "Mean Velocity")	Finding Site=(A12236-008, 99ALOKA, "Descending Aorta")	17
PG	(20247-3, LN, "Peak Gradient")	Finding Site=(A12236-008, 99ALOKA, "Descending Aorta")	17
MPG	(20256-4, LN, "Mean Gradient")	Finding Site=(A12236-008, 99ALOKA, "Descending Aorta")	17
(Transverse Ao Doppler)			
pV	(11726-7, LN, "Peak Velocity")	Finding Site=(A12236-009, 99ALOKA, "Transverse Aorta")	17
VTI	(20354-7, LN, "Velocity Time Integral")	Finding Site=(A12236-009, 99ALOKA, "Transverse Aorta")	17
MnV	(20352-1, LN, "Mean Velocity")	Finding Site=(A12236-009, 99ALOKA, "Transverse Aorta")	17
PG	(20247-3, LN, "Peak Gradient")	Finding Site=(A12236-009, 99ALOKA, "Transverse Aorta")	17
MPG	(20256-4, LN, "Mean Gradient")	Finding Site=(A12236-009, 99ALOKA, "Transverse Aorta")	17
(Abdominal Ao Doppler)			
pV	(11726-7, LN, "Peak Velocity")	Finding Site=(A12236-010, 99ALOKA, "Abdominal Aorta")	17
VTI	(20354-7, LN, "Velocity Time Integral")	Finding Site=(A12236-010, 99ALOKA, "Abdominal Aorta")	17
MnV	(20352-1, LN, "Mean Velocity")	Finding Site=(A12236-010, 99ALOKA, "Abdominal Aorta")	17
PG	(20247-3, LN, "Peak Gradient")	Finding Site=(A12236-010, 99ALOKA, "Abdominal Aorta")	17
MPG	(20256-4, LN, "Mean Gradient")	Finding Site=(A12236-010, 99ALOKA, "Abdominal Aorta")	17
(ASD)			
pV	(11726-7, LN, "Peak Velocity")	Finding Site=(D4-31220, SRT, "Atrial Septal Defect")	23
PG	(20247-3, LN, "Peak Gradient")	Finding Site=(D4-31220, SRT, "Atrial Septal Defect")	23
MnV	(20352-1, LN, "Mean Velocity")	Finding Site=(D4-31220, SRT, "Atrial Septal Defect")	23
MPG	(20256-4, LN, "Mean Gradient")	Finding Site=(D4-31220, SRT, "Atrial Septal Defect")	23

Label	Concept Name	Modifier	row
VTI	(20354-7, LN, "Velocity Time Integral")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(D4-31220, SRT, "Atrial Septal Defect"), Image Mode=(G-03A2, SRT, "2D mode")	23
ASD	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(D4-31220, SRT, "Atrial Septal Defect"), Image Mode=(G-03A2, SRT, "2D mode")	23
ASD lax	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(D4-31220, SRT, "Atrial Septal Defect"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-0396, SRT, "Parasternal long axis")	23
ASD sax	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(D4-31220, SRT, "Atrial Septal Defect"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-0397, SRT, "Parasternal short axis")	23
ASD Ao Rim	(A12218-001, 99ALOKA, "Atrial Septal Defect Rim diameter(Ao Position)")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(D4-31220, SRT, "Atrial Septal Defect"), Image Mode=(G-03A2, SRT, "2D mode")	23
ASD SVC Rim	(A12218-002, 99ALOKA, "Atrial Septal Defect Rim diameter(SVC Position)")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(D4-31220, SRT, "Atrial Septal Defect"), Image Mode=(G-03A2, SRT, "2D mode")	23
ASD IVC Rim	(A12218-003, 99ALOKA, "Atrial Septal Defect Rim diameter(IVC Position)")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(D4-31220, SRT, "Atrial Septal Defect"), Image Mode=(G-03A2, SRT, "2D mode")	23
(VSD)			
pV	(11726-7, LN, "Peak Velocity")	Finding Site=(D4-31150, SRT, "Ventricular Septal Defect")	23
PG	(20247-3, LN, "Peak Gradient")	Finding Site=(D4-31150, SRT, "Ventricular Septal Defect")	23
MnV	(20352-1, LN, "Mean Velocity")	Finding Site=(D4-31150, SRT, "Ventricular Septal Defect")	23
MPG	(20256-4, LN, "Mean Gradient")	Finding Site=(D4-31150, SRT, "Ventricular Septal Defect")	23
VTI	(20354-7, LN, "Velocity Time Integral")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(D4-31150, SRT, "Ventricular Septal Defect"), Image Mode=(G-03A2, SRT, "2D mode")	23
VSD	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(D4-31150, SRT, "Ventricular Septal Defect"), Image Mode=(G-03A2, SRT, "2D mode")	23

Label	Concept Name	Modifier	row
VSD lax	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(D4-31150, SRT, "Ventricular Septal Defect"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-0396, SRT, "Parasternal long axis")	23
VSD sax	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(D4-31150, SRT, "Ventricular Septal Defect"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-0397, SRT, "Parasternal short axis")	23
VSD EROA	(G-038E, SRT, "Cardiovascular Orifice Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(D4-31150, SRT, "Ventricular Septal Defect"), Image Mode=(G-03A2, SRT, "2D mode")	23
(PFO)			
PFO	(G-038F, SRT, "Cardiovascular Orifice Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(A12236-011, 99ALOKA, "Patent Foramen Ovale"), Image Mode=(G-03A2, SRT, "2D mode")	23
PFO pV	(11726-7, LN, "Peak Velocity")	Finding Site=(A12236-011, 99ALOKA, "Patent Foramen Ovale")	23
(PDA)			
PDA	(81827009, SCT, "Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(A12236-012, 99ALOKA, "Patent Ductus Arteriosus"), Image Mode=(G-03A2, SRT, "2D mode")	23
dias pV	(11726-7, LN, "Peak Velocity")	Finding Site=(A12236-012, 99ALOKA, "Patent Ductus Arteriosus"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	23
dias MnV	(20352-1, LN, "Mean Velocity")	Finding Site=(A12236-012, 99ALOKA, "Patent Ductus Arteriosus"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	23
dias VTI	(20354-7, LN, "Velocity Time Integral")	Finding Site=(A12236-012, 99ALOKA, "Patent Ductus Arteriosus"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole")	23
sys pV	(11726-7, LN, "Peak Velocity")	Finding Site=(A12236-012, 99ALOKA, "Patent Ductus Arteriosus"), Cardiac Cycle Point=(109070, DCM, "End Systole")	23
sys MnV	(20352-1, LN, "Mean Velocity")	Finding Site=(A12236-012, 99ALOKA, "Patent Ductus Arteriosus"), Cardiac Cycle Point=(109070, DCM, "End Systole")	23
sys VTI	(20354-7, LN, "Velocity Time Integral")	Finding Site=(A12236-012, 99ALOKA, "Patent Ductus Arteriosus"), Cardiac Cycle Point=(109070, DCM, "End Systole")	23
(PA)			

Label	Concept Name	Modifier	row
MPA	(18020-8, LN, "Main Pulmonary Artery Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(A12236-013, 99ALOKA, "Main Pulmonary Artery"), Image Mode=(G-03A2, SRT, "2D mode")	18
LPA	(18019-0, LN, "Left Pulmonary Artery Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(A12236-014, 99ALOKA, "Left Pulmonary Artery"), Image Mode=(G-03A2, SRT, "2D mode")	18
RPA	(18021-6, LN, "Right Pulmonary Artery Diameter")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Finding Site=(A12236-015, 99ALOKA, "Right Pulmonary Artery"), Image Mode=(G-03A2, SRT, "2D mode")	18
RAP	(18070-3, LN, "Right Atrium Systolic Pressure")		12
MPA pV	(11726-7, LN, "Peak Velocity")	Finding Site=(A12236-013, 99ALOKA, "Main Pulmonary Artery")	18
MPA PG	(20247-3, LN, "Peak Gradient")	Finding Site=(A12236-013, 99ALOKA, "Main Pulmonary Artery")	18
MPA MnV	(20352-1, LN, "Mean Velocity")	Finding Site=(A12236-013, 99ALOKA, "Main Pulmonary Artery")	18
MPA MPG	(20256-4, LN, "Mean Gradient")	Finding Site=(A12236-013, 99ALOKA, "Main Pulmonary Artery")	18
MPA VTI	(20354-7, LN, "Velocity Time Integral")	Finding Site=(A12236-013, 99ALOKA, "Main Pulmonary Artery")	18
LPA pV	(11726-7, LN, "Peak Velocity")	Finding Site=(A12236-014, 99ALOKA, "Left Pulmonary Artery")	18
LPA PG	(20247-3, LN, "Peak Gradient")	Finding Site=(A12236-014, 99ALOKA, "Left Pulmonary Artery")	18
LPA MnV	(20352-1, LN, "Mean Velocity")	Finding Site=(A12236-014, 99ALOKA, "Left Pulmonary Artery")	18
LPA MPG	(20256-4, LN, "Mean Gradient")	Finding Site=(A12236-014, 99ALOKA, "Left Pulmonary Artery")	18
LPA VTI	(20354-7, LN, "Velocity Time Integral")	Finding Site=(A12236-014, 99ALOKA, "Left Pulmonary Artery")	18
RPA pV	(11726-7, LN, "Peak Velocity")	Finding Site=(A12236-015, 99ALOKA, "Right Pulmonary Artery")	18
RPA PG	(20247-3, LN, "Peak Gradient")	Finding Site=(A12236-015, 99ALOKA, "Right Pulmonary Artery")	18
RPA MnV	(20352-1, LN, "Mean Velocity")	Finding Site=(A12236-015, 99ALOKA, "Right Pulmonary Artery")	18
RPA MPG	(20256-4, LN, "Mean Gradient")	Finding Site=(A12236-015, 99ALOKA, "Right Pulmonary Artery")	18
RPA VTI	(20354-7, LN, "Velocity Time Integral")	Finding Site=(A12236-015, 99ALOKA, "Right Pulmonary Artery")	18
(TDI PW MA)			

Label	Concept Name	Modifier	row
Sm1	(A12222-048, 99ALOKA, "Peak Systolic Tissue Velocity 1")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
Sm2	(A12222-051, 99ALOKA, "Peak Systolic Tissue Velocity 2")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
sMnV	(20352-1, LN, "Mean Velocity")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
Em	(78185-6, LN, "MV med ann Peak e' wave US.TDI")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
Am	(A12222-045, 99ALOKA, "Peak Diastolic Tissue Velocity During Atrial Systole")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
dMnV	(20352-1, LN, "Mean Velocity")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
RIVRT	(A12222-042, 99ALOKA, "Isovolumic Relaxation Time")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
RIVCT	(A12222-039, 99ALOKA, "Isovolumic Contraction Time")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
time	(A12222-059, 99ALOKA, "Time")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
v1	(A12222-024, 99ALOKA, "Blood Flow Velocity 1")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
v2	(A12222-027, 99ALOKA, "Blood Flow Velocity 2")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
AccT	(20168-1, LN, "Acceleration Time")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	14

Label	Concept Name	Modifier	row
Em/Am	(A12222-035, 99ALOKA, "Ea to Aa Ratio")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
E/Em	(A12203-011, 99ALOKA, "Ratio of MV E-Wave Peak Vel. to Early Diastolic Myocardium Vel.")		9
ACC	(20167-3, LN, "Acceleration Slope")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging")	14
(TDI PW 1)			
Sm1	(A12222-049, 99ALOKA, "Peak Systolic Tissue Velocity 1(1)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
Sm2	(A12222-052, 99ALOKA, "Peak Systolic Tissue Velocity 2(1)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
sMnV	(A12222-008, 99ALOKA, "Mean Velocity (1)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
Em	(A12222-018, 99ALOKA, "Ann Peak e' wave (1)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
Am	(A12222-046, 99ALOKA, "Peak Diastolic Tissue Velocity During Atrial Systole (1)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
dMnV	(A12222-008, 99ALOKA, "Mean Velocity (1)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
RIVRT	(A12222-043, 99ALOKA, "Isovolumic Relaxation Time (1)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
RIVCT	(A12222-040, 99ALOKA, "Isovolumic Contraction Time (1)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
time	(A12222-060, 99ALOKA, "Time (1)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9

Label	Concept Name	Modifier	row
v1	(A12222-025, 99ALOKA, "Blood Flow Velocity 1(1)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
v2	(A12222-028, 99ALOKA, "Blood Flow Velocity 2(1)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
AccT	(A12222-016, 99ALOKA, "Acceleration Time (1)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
Em/Am	(A12222-036, 99ALOKA, "Ea to Aa Ratio (1)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
E/Em	(A12222-054, 99ALOKA, "Ratio of MV E-Wave Peak V to Early Diastolic Myocardium V (1)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
ACC	(A12222-014, 99ALOKA, "Acceleration Slope (1)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
(TDI PW 2)			
Sm1	(A12222-050, 99ALOKA, "Peak Systolic Tissue Velocity 1(2)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
Sm2	(A12222-053, 99ALOKA, "Peak Systolic Tissue Velocity 2(2)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
sMnV	(A12222-009, 99ALOKA, "Mean Velocity (2)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
Em	(A12222-019, 99ALOKA, "Ann Peak e' wave (2)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
Am	(A12222-047, 99ALOKA, "Peak Diastolic Tissue Velocity During Atrial Systole (2)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
dMnV	(A12222-009, 99ALOKA, "Mean Velocity (2)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9

Label	Concept Name	Modifier	row
RIVRT	(A12222-044, 99ALOKA, "Isovolumic Relaxation Time (2)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
RIVCT	(A12222-041, 99ALOKA, "Isovolumic Contraction Time (2)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
time	(A12222-061, 99ALOKA, "Time (2)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
v1	(A12222-026, 99ALOKA, "Blood Flow Velocity 1(2)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
v2	(A12222-029, 99ALOKA, "Blood Flow Velocity 2(2)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
AccT	(A12222-017, 99ALOKA, "Acceleration Time (2)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging"), Image Mode=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
Em/Am	(A12222-037, 99ALOKA, "Ea to Aa Ratio (2)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
E/Em	(A12222-055, 99ALOKA, "Ratio of MV E-Wave Peak V to Early Diastolic Myocardium V (2)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
ACC	(A12222-015, 99ALOKA, "Acceleration Slope (2)")	Image Mode(Group)=(P5-B0128, SRT, "Tissue Doppler Imaging")	9
(GLS3P)			
GLS3	(A12203-013, 99ALOKA, "Global Longitudinal Strain")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125208, DCM, "Method of Disks, Single Plane"), Image View=(A12226-001, 99ALOKA, "Apical three chamber")	9
GLS(3P)	(A12203-013, 99ALOKA, "Global Longitudinal Strain")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(A12228-015, 99ALOKA, "Method of Disks, 3 planes")	9
ML3d	(A12203-019, 99ALOKA, "Myocardial Length")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(A12226-001, 99ALOKA, "Apical three chamber")	9

Label	Concept Name	Modifier	row
ML3s	(A12203-019, 99ALOKA, "Myocardial Length")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(A12226-001, 99ALOKA, "Apical three chamber")	9
ML4d	(A12203-019, 99ALOKA, "Myocardial Length")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19C, SRT, "Apical four chamber")	9
ML4s	(A12203-019, 99ALOKA, "Myocardial Length")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19C, SRT, "Apical four chamber")	9
ML2d	(A12203-019, 99ALOKA, "Myocardial Length")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19B, SRT, "Apical two chamber")	9
ML2s	(A12203-019, 99ALOKA, "Myocardial Length")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(G-A19B, SRT, "Apical two chamber")	9
LVL3d	(18077-8, LN, "Left Ventricle diastolic major axis")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(A12226-001, 99ALOKA, "Apical three chamber")	9
LVL3d	(G-0375, SRT, "Left Ventricular Diastolic Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125220, DCM, "Planimetry"), Cardiac Cycle Point=(F-32011, SRT, "End Diastole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(A12226-001, 99ALOKA, "Apical three chamber")	9
LVL3s	(18076-0, LN, "Left Ventricle systolic major axis")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(A12226-001, 99ALOKA, "Apical three chamber")	9

Label	Concept Name	Modifier	row
LVL3s	(G-0374, SRT, "Left Ventricular Systolic Area")	Image Mode(Group)=(G-03A2, SRT, "2D mode"), Measurement Method=(125220, DCM, "Planimetry"), Cardiac Cycle Point=(109070, DCM, "End Systole"), Image Mode=(G-03A2, SRT, "2D mode"), Image View=(A12226-001, 99ALOKA, "Apical three chamber")	9
(LV Function (Pombo, Teichholz, Gibson))			
LVMi2.7(Devereux)	(A12203-018, 99ALOKA, "Left Ventricle Mass adjusted by Height(2.7)")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(A12232-002, 99ALOKA, "Left Ventricle Mass by Penn")	9
LVMi2.16(Devereux)	(A12203-017, 99ALOKA, "Left Ventricle Mass adjusted by Height(2.16)")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(A12232-002, 99ALOKA, "Left Ventricle Mass by Penn")	9
LVMi1.7(Devereux)	(A12203-016, 99ALOKA, "Left Ventricle Mass adjusted by Height(1.7)")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(A12232-002, 99ALOKA, "Left Ventricle Mass by Penn")	9
LVMi2.7(ASE)	(A12203-018, 99ALOKA, "Left Ventricle Mass adjusted by Height(2.7)")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(125221, DCM, "Left Ventricle Mass by M-mode")	9
LVMi2.16(ASE)	(A12203-017, 99ALOKA, "Left Ventricle Mass adjusted by Height(2.16)")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(125221, DCM, "Left Ventricle Mass by M-mode")	9
LVMi1.7(ASE)	(A12203-016, 99ALOKA, "Left Ventricle Mass adjusted by Height(1.7)")	Image Mode(Group)=(G-0394, SRT, "M mode"), Measurement Method=(125221, DCM, "Left Ventricle Mass by M-mode")	9

- Note:
1. When these HRs are measured at a time, only one HR is included in Structured Report. The precedence is LV Volume(Simpson(Disc)), LV Volume, then LV Function.
 2. It is configurable not to include Image Mode(Group) and/or Image Mode modifiers. In this case, only a value of either B mode or M mode is included by priority.
 3. It is configurable to include a value of either B mode or M mode by priority.
 4. If RV MPI of RV MPI measurement and RV MPI of RV MPI(Tissue) measurement are measured, two (399154007, SCT, "Right Ventricular Index of Myocardial Performance") with the same modifiers are included in Structured Report. The former one is for RV MPI measurement and the latter one is for RV MPI(Tissue) measurement.

8.6.2.4 TID 5200 Echocardiography Procedure Report (Advanced Cardiac)

The following tables present information used in Structured Reports for this template. The tables are sorted by the label value, which corresponds to the cardiac measurement label displayed on the Worksheet window.

The tables include the actual Coding Scheme Designator, Code Value and Code Meaning used for that label.

Body Surface Area Equations

Mosteller is used as Body Surface Area to calculate indexed value.

[Mosteller, R.D. 'Simplified Calculation of Body SurfaceArea.' N Engl J Med 1987 Oct 22;317(17):1098].

Left Ventricle

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
Teichholz (2D)	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
RVAWd	LN	18153-7	Right Ventricular Anterior Wall Diastolic Thickness	SRT	T-32500	Right Ventricle	SRT	G-03A2	2D mode (Image Mode)
RVIDd	LN	20304-2	Right Ventricular Internal Diastolic Dimension	SRT	T-32500	Right Ventricle	SRT	G-03A2	2D mode (Image Mode)
IVSd	LN	18154-5	Interventricular Septum Diastolic Thickness	SRT	T-32600	Left Ventricle	SRT	G-03A2	2D mode (Image Mode)
LVIDd	LN	29436-3	Left Ventricle Internal End Diastolic Dimension	SRT	T-32600	Left Ventricle	SRT	G-03A2	2D mode (Image Mode)
LVPWd	LN	18152-9	Left Ventricle Posterior Wall Diastolic Thickness	SRT	T-32600	Left Ventricle	SRT	G-03A2	2D mode (Image Mode)
IVSs	LN	18158-6	Interventricular Septum Systolic Thickness	SRT	T-32600	Left Ventricle	SRT	G-03A2	2D mode (Image Mode)

LVIDs	LN	29438-9	Left Ventricle Internal Systolic Dimension	SRT	T-32600	Left Ventricle	SRT	G-03A2	2D mode (Image Mode)
LVPWs	LN	18156-0	Left Ventricle Posterior Wall Systolic Thickness	SRT	T-32600	Left Ventricle	SRT	G-03A2	2D mode (Image Mode)
HR	LN	8867-4	Heart rate	SRT	T-32600	Left Ventricle	DCM	125209	Teichholz (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
EDV	LN	18026-5	Left Ventricular End Diastolic Volume	SRT	T-32600	Left Ventricle	DCM	125209	Teichholz (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
ESV	LN	18148-7	Left Ventricular End Systolic Volume	SRT	T-32600	Left Ventricle	DCM	125209	Teichholz (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
LV EDVI	LN	18026-5	Left Ventricular End Diastolic Volume	SRT	T-32600	Left Ventricle	DCM	125209	Teichholz (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							LN	8277-6	Body Surface Area (Index)
LV ESVI	LN	18148-7	Left Ventricular End Systolic Volume	SRT	T-32600	Left Ventricle	DCM	125209	Teichholz (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							LN	8277-6	Body Surface Area (Index)

LV EDVI (ht)	LN	18026-5	Left Ventricular End Diastolic Volume	SRT	T-32600	Left Ventricle	DCM	125209	Teichholz (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							LN	8302-2	Patient Height (Index)
LV ESVI (ht)	LN	18148-7	Left Ventricular End Systolic Volume	SRT	T-32600	Left Ventricle	DCM	125209	Teichholz (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							LN	8302-2	Patient Height (Index)
SV	SRT	F-32120	Stroke Volume	SRT	T-32600	Left Ventricle	DCM	125209	Teichholz (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
CO	SRT	F-32100	Cardiac Output	SRT	T-32600	Left Ventricle	DCM	125209	Teichholz (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
EF	LN	18043-0	Left Ventricular Ejection Fraction	SRT	T-32600	Left Ventricle	DCM	125209	Teichholz (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
FS	LN	18051-3	Left Ventricular Fractional Shortening	SRT	T-32600	Left Ventricle	SRT	G-03A2	2D mode (Image Mode)
SI	SRT	F-00078	Stroke Index	SRT	T-32600	Left Ventricle	DCM	125209	Teichholz (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)

CI	SRT	F-32110	Cardiac Index	SRT	T-32600	Left Ventricle	DCM	125209	Teichholz (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
LV MASS-c	LN	18087-7	Left Ventricle Mass	SRT	T-32600	Left Ventricle	SRT	G-03A2	2D mode (Image Mode)
							DCM	125221	Left Ventricle Mass by M-mode (Measurement Method)
LV MASSI	LN	18087-7	Left Ventricle Mass	SRT	T-32600	Left Ventricle	SRT	G-03A2	2D mode (Image Mode)
							LN	8277-6	Body Surface Area (Index)
LV MASS Penn	LN	18087-7	Left Ventricle Mass	SRT	T-32600	Left Ventricle	SRT	G-03A2	2D mode (Image Mode)
							99IM PRO OF	MASSPENN	Left Ventricle Mass by Penn (Measurement Method)
LV RWT	99IM PRO OF	RWT	Relative Wall Thickness	SRT	T-32600	Left Ventricle	SRT	G-03A2	2D mode (Image Mode)
LVIDd sax PM ant-post	LN	29436-3	Left Ventricle Internal End Diastolic Dimension	SRT	T-32600	Left Ventricle	SRT	G-03A2	2D mode (Image Mode)
							DCM	122675	Anterior-Posterior (Finding Site)
							SRT	G-039B	Parasternal short axis at the Papillary Muscle level (Image View)
LVIDd sax PM sep-lat	LN	29436-3	Left Ventricle Internal End Diastolic Dimension	SRT	T-32600	Left Ventricle	SRT	G-03A2	2D mode (Image Mode)
							99IM PRO OF	SEPLAT	Septal-Lateral

							SRT	G-039B	Parasternal short axis at the Papillary Muscle level (Image View)
LVEId	99IM PRO OF	LVEI	LV Eccentricity Index	SRT	T-32600	Left Ventricle	SRT	F-32011	End Diastole (Cardiac Cycle Point)
LVIDs sax PM ant-post	LN	29438-9	Left Ventricle Internal Systolic Dimension	SRT	T-32600	Left Ventricle	SRT	G-03A2	2D mode (Image Mode)
							DCM	122675	Anterior-Posterior (Finding Site)
							SRT	G-039B	Parasternal short axis at the Papillary Muscle level (Image View)
LVIDs sax PM sep-lat	LN	29438-9	Left Ventricle Internal Systolic Dimension	SRT	T-32600	Left Ventricle	SRT	G-03A2	2D mode (Image Mode)
							99IM PRO OF	SEPLAT	Septal-Lateral
							SRT	G-039B	Parasternal short axis at the Papillary Muscle level (Image View)
LVEIs	99IM PRO OF	LVEI	LV Eccentricity Index	SRT	T-32600	Left Ventricle	SRT	R-FAB5B	End Systole (Cardiac Cycle Point)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
RVAWd	LN	18153-7	Right Ventricular Anterior Wall Diastolic Thickness	SRT	T-32500	Right Ventricle	SRT	G-0394	MMode (Image Mode)

RVIDd	LN	20304-2	Right Ventricular Internal Diastolic Dimension	SRT	T-32500	Right Ventricle	DCM	125209	Teichholz (Measurement Method)
							SRT	G-0394	MMode (Image Mode)
IVSd	LN	18154-5	Interventricular Septum Diastolic Thickness	SRT	T-32600	Left Ventricle	SRT	G-0394	MMode (Image Mode)
LVIDd	LN	29436-3	Left Ventricle Internal End Diastolic Dimension	SRT	T-32600	Left Ventricle	SRT	G-0394	MMode (Image Mode)
LVPWd	LN	18152-9	Left Ventricle Posterior Wall Diastolic Thickness	SRT	T-32600	Left Ventricle	SRT	G-0394	MMode (Image Mode)
IVSs	LN	18158-6	Interventricular Septum Systolic Thickness	SRT	T-32600	Left Ventricle	SRT	G-0394	MMode (Image Mode)
LVIDs	LN	29438-9	Left Ventricle Internal Systolic Dimension	SRT	T-32600	Left Ventricle	SRT	G-0394	MMode (Image Mode)
LVPWs	LN	18156-0	Left Ventricle Posterior Wall Systolic Thickness	SRT	T-32600	Left Ventricle	SRT	G-0394	MMode (Image Mode)
SPWMD	99IMPROOF	SPWMD	Septal-to-Posterior-Wall Motion Delay	SRT	T-32600	Left Ventricle	SRT	G-0394	MMode (Image Mode)
HR	LN	8867-4	Heart rate	SRT	T-32600	Left Ventricle	DCM	125209	Teichholz (Measurement Method)
							SRT	G-0394	MMode (Image Mode)

EDV	LN	18026-5	Left Ventricular End Diastolic Volume	SRT	T-32600	Left Ventricle	DCM	125209	Teichholz (Measurement Method)
							SRT	G-0394	MMode (Image Mode)
ESV	LN	18148-7	Left Ventricular End Systolic Volume	SRT	T-32600	Left Ventricle	DCM	125209	Teichholz (Measurement Method)
							SRT	G-0394	MMode (Image Mode)
EDVI	LN	18026-5	Left Ventricular End Diastolic Volume	SRT	T-32600	Left Ventricle	DCM	125209	Teichholz (Measurement Method)
							SRT	G-0394	MMode (Image Mode)
							LN	8277-6	Body Surface Area (Index)
ESVI	LN	18148-7	Left Ventricular End Systolic Volume	SRT	T-32600	Left Ventricle	DCM	125209	Teichholz (Measurement Method)
							SRT	G-0394	MMode (Image Mode)
							LN	8277-6	Body Surface Area (Index)
EDVI (ht)	LN	18026-5	Left Ventricular End Diastolic Volume	SRT	T-32600	Left Ventricle	DCM	125209	Teichholz (Measurement Method)
							SRT	G-0394	MMode (Image Mode)
							LN	8302-2	Patient Height (Index)
ESVI (ht)	LN	18148-7	Left Ventricular End Systolic Volume	SRT	T-32600	Left Ventricle	DCM	125209	Teichholz (Measurement Method)
							SRT	G-0394	MMode (Image Mode)

							LN	8302-2	Patient Height (Index)
SV	SRT	F-32120	Stroke Volume	SRT	T-32600	Left Ventricle	DCM	125209	Teichholz (Measurement Method)
							SRT	G-0394	MMode (Image Mode)
CO	SRT	F-32100	Cardiac Output	SRT	T-32600	Left Ventricle	DCM	125209	Teichholz (Measurement Method)
							SRT	G-0394	MMode (Image Mode)
EF	LN	18043-0	Left Ventricular Ejection Fraction	SRT	T-32600	Left Ventricle	DCM	125209	Teichholz (Measurement Method)
							SRT	G-0394	MMode (Image Mode)
FS	LN	18051-3	Left Ventricular Fractional Shortening	SRT	T-32600	Left Ventricle	SRT	G-0394	MMode (Image Mode)
SI	SRT	F-00078	Stroke Index	SRT	T-32600	Left Ventricle	DCM	125209	Teichholz (Measurement Method)
							SRT	G-0394	MMode (Image Mode)
CI	SRT	F-32110	Cardiac Index	SRT	T-32600	Left Ventricle	DCM	125209	Teichholz (Measurement Method)
							SRT	G-0394	MMode (Image Mode)
mVcf	99IM PRO OF	CA_LV_V CF	Velocity of Circumferential Fiber Shortening	SRT	T-32600	Left Ventricle	SRT	G-0394	MMode (Image Mode)
LV MASS-c	LN	18087-7	Left Ventricle Mass	SRT	T-32600	Left Ventricle	DCM	125221	Left Ventricle Mass by M-mode (Measurement Method)

LV MASSI	LN	18087-7	Left Ventricle Mass	SRT	T-32600	Left Ventricle	SRT	G-0394	MMode (Image Mode)
							LN	8277-6	Body Surface Area (Index)
LV MASS Penn	LN	18087-7	Left Ventricle Mass	SRT	T-32600	Left Ventricle	SRT	G-0394	MMode (Image Mode)
							99IM PRO OF	MASSPEN	Left Ventricle Mass by Penn (Measurement Method)
LV RWT	99IM PRO OF	RWT	Relative Wall Thickness	SRT	T-32600	Left Ventricle	SRT	G-0394	MMode (Image Mode)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
Simpson	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
EDV	LN	18026-5	Left Ventricular End Diastolic Volume	SRT	T-32600	Left Ventricle	DCM	125207	Method of Disks, Biplane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
ESV	LN	18148-7	Left Ventricular End Systolic Volume	SRT	T-32600	Left Ventricle	DCM	125207	Method of Disks, Biplane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
EDVI	LN	18026-5	Left Ventricular End Diastolic Volume	SRT	T-32600	Left Ventricle	DCM	125207	Method of Disks, Biplane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)

							LN	8277-6	Body Surface Area (Index)
ESVI	LN	18148-7	Left Ventricular End Systolic Volume	SRT	T-32600	Left Ventricle	DCM	125207	Method of Disks, Biplane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							LN	8277-6	Body Surface Area (Index)
EDVI (ht)	LN	18026-5	Left Ventricular End Diastolic Volume	SRT	T-32600	Left Ventricle	DCM	125207	Method of Disks, Biplane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							LN	8302-2	Patient Height (Index)
ESVI (ht)	LN	18148-7	Left Ventricular End Systolic Volume	SRT	T-32600	Left Ventricle	DCM	125207	Method of Disks, Biplane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							LN	8302-2	Patient Height (Index)
SV	SRT	F-32120	Stroke Volume	SRT	T-32600	Left Ventricle	DCM	125207	Method of Disks, Biplane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)

CO	SRT	F-32100	Cardiac Output	SRT	T-32600	Left Ventricle	DCM	125207	Method of Disks, Biplane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
EF	LN	18043-0	Left Ventricular Ejection Fraction	SRT	T-32600	Left Ventricle	DCM	125207	Method of Disks, Biplane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
SI	SRT	F-00078	Stroke Index	SRT	T-32600	Left Ventricle	DCM	125207	Method of Disks, Biplane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
CI	SRT	F-32110	Cardiac Index	SRT	T-32600	Left Ventricle	DCM	125207	Method of Disks, Biplane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
HR	LN	8867-4	Heart rate	SRT	T-32600	Left Ventricle	DCM	125207	Method of Disks, Biplane (Measurement Method)
EDV (A4C)	LN	18026-5	Left Ventricular End Diastolic Volume	SRT	T-32600	Left Ventricle	DCM	125208	Method of Disks, Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19C	Apical four chamber (Image View)

ESV (A4C)	LN	18148-7	Left Ventricular End Systolic Volume	SRT	T-32600	Left Ventricle	DCM	125208	Method of Disks, Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19C	Apical four chamber (Image View)
EDVI (A4C)	LN	18026-5	Left Ventricular End Diastolic Volume	SRT	T-32600	Left Ventricle	DCM	125208	Method of Disks, Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19C	Apical four chamber (Image View)
							LN	8277-6	Body Surface Area (Index)
ESVI (A4C)	LN	18148-7	Left Ventricular End Systolic Volume	SRT	T-32600	Left Ventricle	DCM	125208	Method of Disks, Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19C	Apical four chamber (Image View)
							LN	8277-6	Body Surface Area (Index)
EDVI (ht) (A4C)	LN	18026-5	Left Ventricular End Diastolic Volume	SRT	T-32600	Left Ventricle	DCM	125208	Method of Disks, Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)

							SRT	G-A19C	Apical four chamber (Image View)
							LN	8302-2	Patient Height (Index)
ESVI (ht) (A4C)	LN	18148-7	Left Ventricular End Systolic Volume	SRT	T-32600	Left Ventricle	DCM	125208	Method of Disks, Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19C	Apical four chamber (Image View)
							LN	8302-2	Patient Height (Index)
SV (A4C)	SRT	F-32120	Stroke Volume	SRT	T-32600	Left Ventricle	DCM	125208	Method of Disks, Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19C	Apical four chamber (Image View)
CO (A4C)	SRT	F-32100	Cardiac Output	SRT	T-32600	Left Ventricle	DCM	125208	Method of Disks, Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19C	Apical four chamber (Image View)
EF (A4C)	LN	18043-0	Left Ventricular Ejection Fraction	SRT	T-32600	Left Ventricle	DCM	125208	Method of Disks, Single Plane (Measurement Method)

							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19C	Apical four chamber (Image View)
SI (A4C)	SRT	F-00078	Stroke Index	SRT	T-32600	Left Ventricle	DCM	125208	Method of Disks, Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19C	Apical four chamber (Image View)
CI (A4C)	SRT	F-32110	Cardiac Index	SRT	T-32600	Left Ventricle	DCM	125208	Method of Disks, Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19C	Apical four chamber (Image View)
EDV (A2C)	LN	18026-5	Left Ventricular End Diastolic Volume	SRT	T-32600	Left Ventricle	DCM	125208	Method of Disks, Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19B	Apical two chamber (Image View)
ESV (A2C)	LN	18148-7	Left Ventricular End Systolic Volume	SRT	T-32600	Left Ventricle	DCM	125208	Method of Disks, Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)

							SRT	G-A19B	Apical two chamber (Image View)
EDVI (A2C)	LN	18026-5	Left Ventricular End Diastolic Volume	SRT	T-32600	Left Ventricle	DCM	125208	Method of Disks, Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19B	Apical two chamber (Image View)
							LN	8277-6	Body Surface Area (Index)
ESVI (A2C)	LN	18148-7	Left Ventricular End Systolic Volume	SRT	T-32600	Left Ventricle	DCM	125208	Method of Disks, Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19B	Apical two chamber (Image View)
							LN	8277-6	Body Surface Area (Index)
EDVI (ht) (A2C)	LN	18026-5	Left Ventricular End Diastolic Volume	SRT	T-32600	Left Ventricle	DCM	125208	Method of Disks, Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19B	Apical two chamber (Image View)
							LN	8302-2	Patient Height (Index)

ESVI (ht) (A2C)	LN	18148-7	Left Ventricular End Systolic Volume	SRT	T-32600	Left Ventricle	DCM	125208	Method of Disks, Single Plane (Measureme nt Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19B	Apical two chamber (Image View)
							LN	8302-2	Patient Height (Index)
SV (A2C)	SRT	F-32120	Stroke Volume	SRT	T-32600	Left Ventricle	DCM	125208	Method of Disks, Single Plane (Measureme nt Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19B	Apical two chamber (Image View)
CO (A2C)	SRT	F-32100	Cardiac Output	SRT	T-32600	Left Ventricle	DCM	125208	Method of Disks, Single Plane (Measureme nt Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19B	Apical two chamber (Image View)
EF (A2C)	LN	18043-0	Left Ventricular Ejection Fraction	SRT	T-32600	Left Ventricle	DCM	125208	Method of Disks, Single Plane (Measureme nt Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19B	Apical two chamber (Image View)

SI (A2C)	SRT	F-00078	Stroke Index	SRT	T-32600	Left Ventricle	DCM	125208	Method of Disks, Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19B	Apical two chamber (Image View)
CI (A2C)	SRT	F-32110	Cardiac Index	SRT	T-32600	Left Ventricle	DCM	125208	Method of Disks, Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19B	Apical two chamber (Image View)
difD	99IM PRO OF	CA_LV_D IF	Major axis length difference	SRT	T-32600	Left Ventricle	DCM	125207	Method of Disks, Biplane (Measurement Method)
							SRT	F-32011	End Diastole (Cardiac Cycle Point)
							SRT	G-03A2	2D mode (Image Mode)
difS	99IM PRO OF	CA_LV_D IF	Major axis length difference	SRT	T-32600	Left Ventricle	DCM	125207	Method of Disks, Biplane (Measurement Method)
							SRT	R-FAB5B	End Systole (Cardiac Cycle Point)
							SRT	G-03A2	2D mode (Image Mode)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
Area Length	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
LVLD (A4C)	LN	18077-8	Left Ventricle Diastolic Major axis	SRT	T-32600	Left Ventricle	DCM	125205	Area-Length Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19C	Apical four chamber (Image View)
LVAd (A4C)	SRT	G-0375	Left Ventricular Diastolic Area	SRT	T-32600	Left Ventricle	DCM	125205	Area-Length Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19C	Apical four chamber (Image View)
LVLs (A4C)	LN	18076-0	Left Ventricle Systolic Major Axis	SRT	T-32600	Left Ventricle	DCM	125205	Area-Length Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19C	Apical four chamber (Image View)
LVAs (A4C)	SRT	G-0374	Left Ventricular Systolic Area	SRT	T-32600	Left Ventricle	DCM	125205	Area-Length Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19C	Apical four chamber (Image View)

EDV (A4C)	LN	18026-5	Left Ventricular End Diastolic Volume	SRT	T-32600	Left Ventricle	DCM	125205	Area-Length Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19C	Apical four chamber (Image View)
ESV (A4C)	LN	18148-7	Left Ventricular End Systolic Volume	SRT	T-32600	Left Ventricle	DCM	125205	Area-Length Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19C	Apical four chamber (Image View)
SV (A4C)	SRT	F-32120	Stroke Volume	SRT	T-32600	Left Ventricle	DCM	125205	Area-Length Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19C	Apical four chamber (Image View)
CO (A4C)	SRT	F-32100	Cardiac Output	SRT	T-32600	Left Ventricle	DCM	125205	Area-Length Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19C	Apical four chamber (Image View)
EF (A4C)	LN	18043-0	Left Ventricular Ejection Fraction	SRT	T-32600	Left Ventricle	DCM	125205	Area-Length Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)

							SRT	G-A19C	Apical four chamber (Image View)
SI (A4C)	SRT	F-00078	Stroke Index	SRT	T-32600	Left Ventricle	DCM	125205	Area-Length Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19C	Apical four chamber (Image View)
CI (A4C)	SRT	F-32110	Cardiac Index	SRT	T-32600	Left Ventricle	DCM	125205	Area-Length Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19C	Apical four chamber (Image View)
LVLd (A2C)	LN	18077-8	Left Ventricle Diastolic Major axis	SRT	T-32600	Left Ventricle	DCM	125205	Area-Length Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19B	Apical two chamber (Image View)
LVAd (A2C)	SRT	G-0375	Left Ventricular Diastolic Area	SRT	T-32600	Left Ventricle	DCM	125205	Area-Length Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19B	Apical two chamber (Image View)
LVLs (A2C)	LN	18076-0	Left Ventricle Systolic Major Axis	SRT	T-32600	Left Ventricle	DCM	125205	Area-Length Single Plane (Measurement Method)

							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19B	Apical two chamber (Image View)
LVAs (A2C)	SRT	G-0374	Left Ventricular Systolic Area	SRT	T-32600	Left Ventricle	DCM	125205	Area-Length Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19B	Apical two chamber (Image View)
EDV (A2C)	LN	18026-5	Left Ventricular End Diastolic Volume	SRT	T-32600	Left Ventricle	DCM	125205	Area-Length Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19B	Apical two chamber (Image View)
ESV (A2C)	LN	18148-7	Left Ventricular End Systolic Volume	SRT	T-32600	Left Ventricle	DCM	125205	Area-Length Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19B	Apical two chamber (Image View)
SV (A2C)	SRT	F-32120	Stroke Volume	SRT	T-32600	Left Ventricle	DCM	125205	Area-Length Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19B	Apical two chamber (Image View)

CO (A2C)	SRT	F-32100	Cardiac Output	SRT	T-32600	Left Ventricle	DCM	125205	Area-Length Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19B	Apical two chamber (Image View)
EF (A2C)	LN	18043-0	Left Ventricular Ejection Fraction	SRT	T-32600	Left Ventricle	DCM	125205	Area-Length Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19B	Apical two chamber (Image View)
SI (A2C)	SRT	F-00078	Stroke Index	SRT	T-32600	Left Ventricle	DCM	125205	Area-Length Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19B	Apical two chamber (Image View)
CI (A2C)	SRT	F-32110	Cardiac Index	SRT	T-32600	Left Ventricle	DCM	125205	Area-Length Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19B	Apical two chamber (Image View)
HR	LN	8867-4	Heart rate	SRT	T-32600	Left Ventricle	DCM	125205	Area-Length Single Plane (Measurement Method)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
Bullet	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
LVLd	LN	18077-8	Left Ventricle Diastolic Major Axis	SRT	T-32600	Left Ventricle	DCM	125228	Bullet Method (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
LVLS	LN	18076-0	Left Ventricle Systolic Major Axis	SRT	T-32600	Left Ventricle	DCM	125228	Bullet Method (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
LVAd sax MV	SRT	G-0375	Left Ventricular Diastolic Area	SRT	T-32600	Left Ventricle	DCM	125228	Bullet Method (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-0399	Parasternal short axis at the level of the mitral chords (Image View)
LVAs sax MV	SRT	G-0374	Left Ventricular Systolic Area	SRT	T-32600	Left Ventricle	DCM	125228	Bullet Method (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-0399	Parasternal short axis at the level of the mitral chords (Image View)

EDV	LN	18026-5	Left Ventricular End Diastolic Volume	SRT	T-32600	Left Ventricle	DCM	125228	Bullet Method (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
ESV	LN	18148-7	Left Ventricular End Systolic Volume	SRT	T-32600	Left Ventricle	DCM	125228	Bullet Method (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
SV	SRT	F-32120	Stroke Volume	SRT	T-32600	Left Ventricle	DCM	125228	Bullet Method (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
CO	SRT	F-32100	Cardiac Output	SRT	T-32600	Left Ventricle	DCM	125228	Bullet Method (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
EF	LN	18043-0	Left Ventricular Ejection Fraction	SRT	T-32600	Left Ventricle	DCM	125228	Bullet Method (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
SI	SRT	F-00078	Stroke Index	SRT	T-32600	Left Ventricle	DCM	125228	Bullet Method (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
CI	SRT	F-32110	Cardiac Index	SRT	T-32600	Left Ventricle	DCM	125228	Bullet Method (Measurement Method)

							SRT	G-03A2	2D mode (Image Mode)
HR	LN	8867-4	Heart rate	SRT	T-32600	Left Ventricle	DCM	125228	Bullet Method (Measureme nt Method)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
LVAd sax PM Epi	SRT	G-0379	Left Ventricle Epicardial Diastolic Area, psax pap view	SRT	T-32600	Left Ventricle	DCM	125222	Left Ventricle Mass by Truncated Ellipse (Measureme nt Method)
							SRT	G-03A2	2D mode (Image Mode)
LVAd sax PM Endo	SRT	G-0375	Left Ventricular Diastolic Area	SRT	T-32600	Left Ventricle	DCM	125222	Left Ventricle Mass by Truncated Ellipse (Measureme nt Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-039B	Parasternal short axis at the Papillary Muscle level (Image View)
a	SRT	G-0377	Left Ventricle Semi-major Axis Diastolic Dimension	SRT	T-32600	Left Ventricle	DCM	125222	Left Ventricle Mass by Truncated Ellipse (Measureme nt Method)
							SRT	G-03A2	2D mode (Image Mode)

d	SRT	G-0378	Left Ventricle Truncated Semi-major Axis Diastolic Dimension	SRT	T-32600	Left Ventricle	DCM	125222	Left Ventricle Mass by Truncated Ellipse (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
t	99IM PRO OF	THI_Myo	Myocardial Thickness	SRT	T-32600	Left Ventricle	DCM	125222	Left Ventricle Mass by Truncated Ellipse (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
b	99IM PRO OF	RADsax	Short axis radius calculated from short axis cavity area	SRT	T-32600	Left Ventricle	DCM	125222	Left Ventricle Mass by Truncated Ellipse (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
LV MASS	LN	18087-7	Left Ventricle Mass	SRT	T-32600	Left Ventricle	DCM	125222	Left Ventricle Mass by Truncated Ellipse (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
LV MASSI	LN	18087-7	Left Ventricle Mass	SRT	T-32600	Left Ventricle	DCM	125222	Left Ventricle Mass by Truncated Ellipse (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							LN	8277-6	Body Surface Area (Index)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
LV Mass A-L	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
LVAd sax PM Epi	SRT	G-0379	Left Ventricle Epicardial Diastolic Area, psax pap view	SRT	T-32600	Left Ventricle	DCM	125270	Left Ventricle Mass by Area Length (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
LVAd sax PM Endo	SRT	G-0375	Left Ventricular Diastolic Area	SRT	T-32600	Left Ventricle	DCM	125270	Left Ventricle Mass by Area Length (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-039B	Parasternal short axis at the Papillary Muscle level (Image View)
LVLD	LN	18077-8	Left Ventricle Diastolic Major axis	SRT	T-32600	Left Ventricle	DCM	125270	Left Ventricle Mass by Area Length (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
t	99IM PRO OF	THI_Myo	Myocardial Thickness	SRT	T-32600	Left Ventricle	DCM	125270	Left Ventricle Mass by Area Length (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
LV MASS	LN	18087-7	Left Ventricle Mass	SRT	T-32600	Left Ventricle	DCM	125270	Left Ventricle Mass by Area Length (Measurement Method)

							SRT	G-03A2	2D mode (Image Mode)
LV MASSI	LN	18087-7	Left Ventricle Mass	SRT	T-32600	Left Ventricle	DCM	125270	Left Ventricle Mass by Area Length (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							LN	8277-6	Body Surface Area (Index)

AV - LA

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
2D	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
RV diam	LN	20304-2	Right Ventricular Internal Diastolic Dimension	SRT	T-32500	Right Ventricle	SRT	G-03A2	2D mode (Image Mode)
							99IM PRO OF	AV_LA	AV/LA level
Ao diam systole	LN	18015-8	Aortic Root Diameter	SRT	T-42000	Aorta	SRT	F-32020	Systole (Cardiac Cycle Point)
							SRT	G-03A2	2D mode (Image Mode)
Ao diam diastole	LN	18015-8	Aortic Root Diameter	SRT	T-42000	Aorta	SRT	F-32010	Diastole (Cardiac Cycle Point)
							SRT	G-03A2	2D mode (Image Mode)
LA diam systole	SRT	M-02550	Diameter	SRT	T-32300	Left Atrium	SRT	F-32020	Systole (Cardiac Cycle Point)
							SRT	G-03A2	2D mode (Image Mode)

LA diam diastole	SRT	M-02550	Diameter	SRT	T-32300	Left Atrium	SRT	F-32010	Diastole (Cardiac Cycle Point)
							SRT	G-03A2	2D mode (Image Mode)
AO/LA systole	99IM PRO OF	AOs2LAs	Aortic Root to Left Atrium Ratio	SRT	T-32300	Left Atrium	SRT	G-03A2	2D mode (Image Mode)
LA/AO systole	LN	17985-3	Left Atrium to Aortic Root Ratio	SRT	T-32300	Left Atrium	SRT	G-03A2	2D mode (Image Mode)
LAs/AOd	99IM Proof	LAs2AOd	Left Atrium Systole to Aortic Root Diastole Diameter Ratio	SRT	T-32300	Left Atrium	SRT	G-03A2	2D mode (Image Mode)
LA diam systole Index	SRT	M-02550	Diameter	SRT	T-32300	Left Atrium	SRT	F-32020	Systole (Cardiac Cycle Point)
							SRT	G-03A2	2D mode (Image Mode)
							LN	8277-6	Body Surface Area (Index)
ACS	LN	17996-0	Aortic Valve Cusp Separation	SRT	T-35400	Aortic Valve	SRT	G-03A2	2D mode (Image Mode)
LA VMax App.	LN	29486-8	Left Atrial Appendage Peak Velocity	SRT	T-32300	Left Atrium			
Ao STJ diam	99IM PRO OF	CA_DIAM_SINO	Aortic Sinotubular Junction Diameter	SRT	T-42000	Aorta	SRT	G-03A2	2D mode (Image Mode)
Ao Arch diam	LN	18011-7	Aortic Arch Diameter	SRT	T-42000	Aorta	SRT	G-03A2	2D mode (Image Mode)
Ao SV diam	99IM PRO OF	CA_DIAM_SINUS	Aortic Sinuses of Valsalva Diameter	SRT	T-42000	Aorta	SRT	G-03A2	2D mode (Image Mode)

Ao Ring diam	99IM PRO OF	CA_DIAM _ANN	Aortic Annulus Diameter	SRT	T-42000	Aorta	SRT	G-03A2	2D mode (Image Mode)
Ao Asc diam	LN	18012-5	Ascending Aortic Diameter	SRT	T-42000	Aorta	SRT	G-03A2	2D mode (Image Mode)
Ao Desc diam	LN	18013-3	Descending Aortic Diameter	SRT	T-42000	Aorta	SRT	G-03A2	2D mode (Image Mode)
Ao Abd diam	SRT	M-02550	Diameter	SRT	T-42500	Abdominal Aorta	SRT	G-03A2	2D mode (Image Mode)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
MMode	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
RV diam	LN	20304-2	Right Ventricular Internal Diastolic Dimension	SRT	T-32500	Right Ventricle	SRT	G-0394	MMode (Image Mode)
							99IM PRO OF	AV_LA	AV/LA level
Ao diam systole	LN	18015-8	Aortic Root Diameter	SRT	T-42000	Aorta	SRT	F-32020	Systole (Cardiac Cycle Point)
							SRT	G-0394	MMode (Image Mode)
Ao diam diastole	LN	18015-8	Aortic Root Diameter	SRT	T-42000	Aorta	SRT	F-32010	Diastole (Cardiac Cycle Point)
							SRT	G-0394	MMode (Image Mode)
LA diam systole	SRT	M-02550	Diameter	SRT	T-32300	Left Atrium	SRT	F-32020	Systole (Cardiac Cycle Point)
							SRT	G-0394	MMode (Image Mode)
LA diam diastole	SRT	M-02550	Diameter	SRT	T-32300	Left Atrium	SRT	F-32010	Diastole (Cardiac Cycle Point)

							SRT	G-0394	MMode (Image Mode)
AO/LA systole	99IM PRO OF	AOs2LAs	Aortic Root to Left Atrium Ratio	SRT	T-32300	Left Atrium	SRT	G-0394	MMode (Image Mode)
LA/AO systole	LN	17985-3	Left Atrium to Aortic Root Ratio	SRT	T-32300	Left Atrium	SRT	G-0394	MMode (Image Mode)
LAs/AOd	99IM Proof	LAs2AOd	Left Atrium Systole to Aortic Root Diastole Diameter Ratio	SRT	T-32300	Left Atrium	SRT	G-0394	MMode (Image Mode)
LA diam systole Index	SRT	M-02550	Diameter	SRT	T-32300	Left Atrium	SRT	F-32020	Systole (Cardiac Cycle Point)
							SRT	G-0394	MMode (Image Mode)
							LN	8277-6	Body Surface Area (Index)
ACS	LN	17996-0	Aortic Valve Cusp Separation	SRT	T-35400	Aortic Valve	SRT	G-0394	MMode (Image Mode)
LVET	LN	18041-4	Aortic Valve Ejection Time	SRT	T-35400	Aortic Valve	SRT	G-0394	MMode (Image Mode)
LVPEP	LN	59085-1	Pre-Ejection Period	SRT	T-35400	Aortic Valve	SRT	G-0394	MMode (Image Mode)
LVPEP / LVET	LN	59088-5	Pre-Ejection Period/Ejection Time Ratio	SRT	T-35400	Aortic Valve	SRT	G-0394	MMode (Image Mode)

Mitral Valve

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
MV	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM

E pt	LN	18037-2	Mitral Valve E-Wave Peak Velocity	SRT	T-35300	Mitral Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
A pt	LN	17978-8	Mitral Valve A-Wave Peak Velocity	SRT	T-35300	Mitral Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
Dec. Slope	LN	20216-8	Deceleration Slope	SRT	T-35300	Mitral Valve			
E/A	LN	18038-0	Mitral Valve E to A Ratio	SRT	T-35300	Mitral Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
A/E	99IM PRO OF	CA_MV_A 2E	Mitral Valve A to E Ratio	SRT	T-35300	Mitral Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
Dec. Time	SRT	G-0384	Mitral Valve E-Wave Deceleration Time	SRT	T-35300	Mitral Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
E duration	99IM PRO OF	CA_MV_D URe	Mitral Valve E-Wave Duration	SRT	T-35300	Mitral Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
A duration	SRT	G-0385	Mitral Valve A-Wave Duration	SRT	T-35300	Mitral Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
IVRT	LN	18071-1	Left Ventricular Isovolumic Relaxation Time	SRT	T-32600	Left Ventricle			
MV C-O dur	SRT	G-0387	Mitral Valve Closure to Opening Time	SRT	T-35300	Mitral Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
LVET	LN	18041-4	Aortic Valve Ejection Time	SRT	T-35400	Aortic Valve			
MV VTI	LN	20354-7	Velocity Time Integral	SRT	T-35300	Mitral Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
Vmax	LN	11726-7	Peak Velocity	SRT	T-35300	Mitral Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
Vmean	LN	20352-1	Time averaged mean velocity	SRT	T-35300	Mitral Valve	SRT	R-42047	Antegrade Flow (Flow Direction)

PGmax	LN	20247-3	Peak Gradient	SRT	T-35300	Mitral Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
PGmean	LN	20256-4	Mean Gradient	SRT	T-35300	Mitral Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
PHT	LN	20280-4	Pressure Half-Time	SRT	T-35300	Mitral Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
MVA (PHT)	SRT	G-038E	Cardiovascular Orifice Area	SRT	T-35300	Mitral Valve	DCM	125210	Area by Pressure Half-Time (Measurement Method)
							SRT	R-42047	Antegrade Flow (Flow Direction)
MVA (VTI)	SRT	G-038E	Cardiovascular Orifice Area	SRT	T-35300	Mitral Valve	DCM	125215	Continuity Equation by Velocity Time Integral (Measurement Method)
							SRT	R-42047	Antegrade Flow (Flow Direction)
MV diam	SRT	G-038F	Cardiovascular Orifice Diameter	SRT	T-35300	Mitral Valve			
LVOT VTI	LN	20354-7	Velocity Time Integral	SRT	T-32650	Left Ventricle Outflow Tract	SRT	R-42047	Antegrade Flow (Flow Direction)
LVOT diam	SRT	G-038F	Cardiovascular Orifice Diameter	SRT	T-32650	Left Ventricle Outflow Tract			
HR	LN	8867-4	Heart rate	SRT	T-35300	Mitral Valve			
SV	SRT	F-32120	Stroke Volume	SRT	T-35300	Mitral Valve	DCM	125219	Doppler Volume Flow (Measurement Method)
CO	SRT	F-32100	Cardiac Output	SRT	T-35300	Mitral Valve	DCM	125219	Doppler Volume Flow (Measurement Method)

SI	SRT	F-00078	Stroke Index	SRT	T-35300	Mitral Valve	DCM	125219	Doppler Volume Flow (Measurement Method)
CI	SRT	F-32110	Cardiac Index	SRT	T-35300	Mitral Valve	DCM	125219	Doppler Volume Flow (Measurement Method)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
LV MPI	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
LVET	LN	18041-4	Aortic Valve Ejection Time	SRT	T-35400	Aortic Valve			
MV C-O dur	SRT	G-0387	Mitral Valve Closure to Opening Time	SRT	T-35300	Mitral Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
LV MPI	SRT	G-037F	Left Ventricular Index of Myocardial Performance	SRT	T-32600	Left Ventricle			

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
Annulus	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
Ann diam comm	99IM PRO OF	CA_MV_Dia_CC	Mitral Valve Commissural Diameter	SRT	T-35313	Mitral Annulus	SRT	F-32020	Systole (Cardiac Cycle Point)
							SRT	G-03A2	2D mode (Image Mode)
Ann diam ant-post	99IM PRO OF	CA_MV_Dia_AP	Mitral Valve Anterior-Posterior Annulus Diameter	SRT	T-35313	Mitral Annulus	SRT	F-32020	Systole (Cardiac Cycle Point)
							SRT	G-03A2	2D mode (Image Mode)
Ann. Area	SRT	G-A166	Area	SRT	T-35313	Mitral Annulus	SRT	F-32020	Systole (Cardiac Cycle Point)

							SRT	G-03A2	2D mode (Image Mode)
MVA (Trace)	SRT	G-038E	Cardiovascular Orifice Area	SRT	T-35300	Mitral Valve	DCM	125220	Planimetry (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
E'(m)	SRT	G-037A	Left Ventricular Peak Early Diastolic Tissue Velocity	SRT	T-32600	Left Ventricle	SRT	P5-B0128	Tissue Doppler Imaging (Measurement Method)
							SRT	G-0391	Medial Mitral Annulus (Finding Site)
A'(m)	SRT	G-037C	LV Peak Diastolic Tissue Velocity During Atrial Systole	SRT	T-32600	Left Ventricle	SRT	P5-B0128	Tissue Doppler Imaging (Measurement Method)
							SRT	G-0391	Medial Mitral Annulus (Finding Site)
E'/A'(m)	99IM PRO OF	CA_MV_Ea2Aa	Mitral Valve Ea to Aa Ratio	SRT	T-32600	Left Ventricle	SRT	P5-B0128	Tissue Doppler Imaging (Measurement Method)
							SRT	G-0391	Medial Mitral Annulus (Finding Site)
E/E'(m)	SRT	G-037B	Ratio of MV Peak Velocity to LV Peak Tissue Velocity E-Wave	SRT	T-32600	Left Ventricle	SRT	P5-B0128	Tissue Doppler Imaging (Measurement Method)
							SRT	G-0391	Medial Mitral Annulus (Finding Site)
MV Ann VTI	LN	20354-7	Velocity Time Integral	SRT	T-35313	Mitral Annulus	SRT	R-42047	Antegrade Flow (Flow Direction)

Ann Vmax E'	SRT	G-037A	Left Ventricular Peak Early Diastolic Tissue Velocity	SRT	T-35313	Mitral Annulus			
Ann Vmax A'	SRT	G-037C	LV Peak Diastolic Tissue Velocity During Atrial Systole	SRT	T-35313	Mitral Annulus			
MV Ann E/E'	SRT	G-037B	Ratio of MV Peak Velocity to LV Peak Tissue Velocity E-Wave	SRT	T-35313	Mitral Annulus			

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
MMode	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
CE amp	99IM PRO OF	CA_MV_E WaveAmpl	Amplitude E- Wave	SRT	T-35300	Mitral Valve	SRT	G-0394	MMode (Image Mode)
CA amp	99IM PRO OF	CA_MV_A WaveAmpl	Amplitude A-Wave	SRT	T-35300	Mitral Valve	SRT	G-0394	MMode (Image Mode)
DE slope (Flow Propagation)	LN	59127-1	D-E Slope	SRT	T-35300	Mitral Valve	SRT	G-0394	MMode (Image Mode)
DE amp	99IM PRO OF	CA_MV_D EWaveAmpl	Amplitude D-E Wave	SRT	T-35300	Mitral Valve	SRT	G-0394	MMode (Image Mode)
EPSS	LN	18036-4	Mitral Valve EPSS, E wave	SRT	T-35300	Mitral Valve	SRT	G-0394	MMode (Image Mode)
EF slope	LN	18040-6	Mitral Valve E-F Slope by M-Mode	SRT	T-35300	Mitral Valve	SRT	G-0394	MMode (Image Mode)
CA/CE	99IM PRO OF	MV_AMPL_A2AMPL_E	A-Wave Amplitude to E-Wave Amplitude	SRT	T-35300	Mitral Valve	SRT	G-0394	MMode (Image Mode)

CE/CA	99IM PRO OF	MV_AMPL _E2AMPL _A	E-Wave Amplitude to A-Wave Amplitude	SRT	T-35300	Mitral Valve	SRT	G-0394	MMode (Image Mode)
MAPSE	99IM PRO OF	MAPSE	Mitral Annular Plane Systolic Excursion	SRT	T-32600	Left Ventricle	SRT	G-0394	MMode (Image Mode)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
MR	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
MR VTI	LN	20354-7	Velocity Time Integral	SRT	T-35300	Mitral Valve	SRT	G-0367	Regurgitant Flow (Flow Direction)
MR Vmax	LN	11726-7	Peak Velocity	SRT	T-35300	Mitral Valve	SRT	G-0367	Regurgitant Flow (Flow Direction)
PISA Radius	LN	59102-4	Flow Radius	SRT	T-35300	Mitral Valve	DCM	125216	Proximal Isovelocity Surface Area (Measureme nt Method)
							SRT	G-0367	Regurgitant Flow (Flow Direction)
Aliasing Vel	LN	59130-5	Alias Velocity	SRT	T-35300	Mitral Valve	DCM	125216	Proximal Isovelocity Surface Area (Measureme nt Method)
							SRT	G-0367	Regurgitant Flow (Flow Direction)
MR Vmean	LN	20352-1	Time averaged mean velocity	SRT	T-35300	Mitral Valve	SRT	G-0367	Regurgitant Flow (Flow Direction)
MR PGmax	LN	20247-3	Peak Gradient	SRT	T-35300	Mitral Valve	SRT	G-0367	Regurgitant Flow (Flow Direction)
MR PGmean	LN	20256-4	Mean Gradient	SRT	T-35300	Mitral Valve	SRT	G-0367	Regurgitant Flow (Flow Direction)

MR dP/dt	LN	18035-6	Mitral Regurgitation dP/dt derived from Mitral Reg. velocity	SRT	T-35300	Mitral Valve	SRT	G-0367	Regurgitant Flow (Flow Direction)
MR dt	LN	20217-6	Deceleration Time	SRT	T-35300	Mitral Valve	SRT	G-0367	Regurgitant Flow (Flow Direction)
MR VC diam	99IM PRO OF	CA_DIAM_VC	Vena Contracta Diameter	SRT	T-35300	Mitral Valve	SRT	G-03A2	2D mode (Image Mode)
Jet Area	LN	59123-0	Jet Area	SRT	T-35300	Mitral Valve	DCM	125220	Planimetry (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
Jet Height	99IM PRO OF	HEIGHT_JET	Jet Height	SRT	T-35300	Mitral Valve	SRT	G-03A2	2D mode (Image Mode)
MV VTI	LN	20354-7	Velocity Time Integral	SRT	T-35300	Mitral Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
MV diam	SRT	G-038F	Cardiovascular Orifice Diameter	SRT	T-35300	Mitral Valve			
LVOT VTI	LN	20354-7	Velocity Time Integral	SRT	T-32650	Left Ventricle Outflow Tract	SRT	R-42047	Antegrade Flow (Flow Direction)
LVOT diam	SRT	G-038F	Cardiovascular Orifice Diameter	SRT	T-32650	Left Ventricle Outflow Tract			
SV	SRT	F-32120	Stroke Volume	SRT	T-35300	Mitral Valve	DCM	125219	Doppler Volume Flow (Measurement Method)
LVOT SV	SRT	F-32120	Stroke Volume	SRT	T-32650	Left Ventricle Outflow Tract	DCM	125219	Doppler Volume Flow (Measurement Method)
MR. Reg Volume	99IM PRO OF	CA_REG_VOL	Regurgitation Volume	SRT	T-35300	Mitral Valve	SRT	G-0367	Regurgitant Flow (Flow Direction)

MR. Reg Fraction	SRT	G-0390	Regurgitant Fraction	SRT	T-35300	Mitral Valve	SRT	G-0367	Regurgitant Flow (Flow Direction)
PISA Flow Rate	LN	34141-2	Peak Instantaneous Flow Rate	SRT	T-35300	Mitral Valve	DCM	125216	Proximal Isovelocity Surface Area (Measurement Method)
							SRT	G-0367	Regurgitant Flow (Flow Direction)
PISA ERO Area	SRT	G-038E	Cardiovascular Orifice Area	SRT	T-35300	Mitral Valve	DCM	125216	Proximal Isovelocity Surface Area (Measurement Method)
							SRT	G-0367	Regurgitant Flow (Flow Direction)
PISA. Reg Volume	LN	33878-0	Volume Flow	SRT	T-35300	Mitral Valve	DCM	125216	Proximal Isovelocity Surface Area (Measurement Method)
							SRT	G-0367	Regurgitant Flow (Flow Direction)
PISA. Reg Fraction	SRT	G-0390	Regurgitant Fraction	SRT	T-35300	Mitral Valve	DCM	125216	Proximal Isovelocity Surface Area (Measurement Method)
							SRT	G-0367	Regurgitant Flow (Flow Direction)

Aortic Valve

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
AV VTI	LN	20354-7	Velocity Time Integral	SRT	T-35400	Aortic Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
HR	LN	8867-4	Heart rate	SRT	T-35400	Aortic Valve			

AV Vmax	LN	11726-7	Peak Velocity	SRT	T-35400	Aortic Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
AV Vmean	LN	20352-1	Time averaged mean velocity	SRT	T-35400	Aortic Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
AV PGmax	LN	20247-3	Peak Gradient	SRT	T-35400	Aortic Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
AV PGmean	LN	20256-4	Mean Gradient	SRT	T-35400	Aortic Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
AR Vmax	LN	11726-7	Peak Velocity	SRT	T-35400	Aortic Valve	SRT	G-0367	Regurgitant Flow (Flow Direction)
AR PGmax	LN	20247-3	Peak Gradient	SRT	T-35400	Aortic Valve	SRT	G-0367	Regurgitant Flow (Flow Direction)
AR Ved	LN	11653-3	End Diastolic Velocity	SRT	T-35400	Aortic Valve	SRT	G-0367	Regurgitant Flow (Flow Direction)
AR PGed	99IM PRO OF	CA_AV_P	Aortic Valve Pressure	SRT	T-35400	Aortic Valve	SRT	F-32011	End Diastole (Cardiac Cycle Point)
							SRT	G-0367	Regurgitant Flow (Flow Direction)
AR Dec. Time	LN	20217-6	Deceleration Time	SRT	T-35400	Aortic Valve	DCM	125210	Area by Pressure Half-Time (Measurement Method)
							SRT	G-0367	Regurgitant Flow (Flow Direction)
AR Dec. Slope	LN	20216-8	Deceleration Slope	SRT	T-35400	Aortic Valve	DCM	125210	Area by Pressure Half-Time (Measurement Method)
							SRT	G-0367	Regurgitant Flow (Flow Direction)

AR PHT	LN	20280-4	Pressure Half-Time	SRT	T-35400	Aortic Valve	DCM	125210	Area by Pressure Half-Time (Measurement Method)
							SRT	G-0367	Regurgitant Flow (Flow Direction)
AR VTI	LN	20354-7	Velocity Time Integral	SRT	T-35400	Aortic Valve	SRT	G-0367	Regurgitant Flow (Flow Direction)
AR Vmean	LN	20352-1	Time averaged mean velocity	SRT	T-35400	Aortic Valve	SRT	G-0367	Regurgitant Flow (Flow Direction)
AR PGmean	LN	20256-4	Mean Gradient	SRT	T-35400	Aortic Valve	SRT	G-0367	Regurgitant Flow (Flow Direction)
PISA Radius	LN	59102-4	Flow Radius	SRT	T-35400	Aortic Valve	DCM	125216	Proximal Isovelocity Surface Area (Measurement Method)
							SRT	G-0367	Regurgitant Flow (Flow Direction)
PISA Aliasing Vel	LN	59130-5	Alias Velocity	SRT	T-35400	Aortic Valve	DCM	125216	Proximal Isovelocity Surface Area (Measurement Method)
							SRT	G-0367	Regurgitant Flow (Flow Direction)
PISA Flow Rate	LN	34141-2	Peak Instantaneous Flow Rate	SRT	T-35400	Aortic Valve	DCM	125216	Proximal Isovelocity Surface Area (Measurement Method)
							SRT	G-0367	Regurgitant Flow (Flow Direction)
PISA ERO Area	SRT	G-038E	Cardiovascular Orifice Area	SRT	T-35400	Aortic Valve	DCM	125216	Proximal Isovelocity Surface Area (Measurement Method)

							SRT	G-0367	Regurgitant Flow (Flow Direction)
PISA Reg. Volume	LN	33878-0	Volume Flow	SRT	T-35400	Aortic Valve	DCM	125216	Proximal Isovelocity Surface Area (Measurement Method)
							SRT	G-0367	Regurgitant Flow (Flow Direction)
PISA Reg. Fraction	SRT	G-0390	Regurgitant Fraction	SRT	T-35400	Aortic Valve	DCM	125216	Proximal Isovelocity Surface Area (Measurement Method)
							SRT	G-0367	Regurgitant Flow (Flow Direction)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
LVOT VTI	LN	20354-7	Velocity Time Integral	SRT	T-32650	Left Ventricle Outflow Tract	SRT	R-42047	Antegrade Flow (Flow Direction)
LVOT Vmax	LN	11726-7	Peak Velocity	SRT	T-32650	Left Ventricle Outflow Tract	SRT	R-42047	Antegrade Flow (Flow Direction)
LVOT Vmean	LN	20352-1	Time averaged mean velocity	SRT	T-32650	Left Ventricle Outflow Tract	SRT	R-42047	Antegrade Flow (Flow Direction)
LVOT PGmax	LN	20247-3	Peak Gradient	SRT	T-32650	Left Ventricle Outflow Tract	SRT	R-42047	Antegrade Flow (Flow Direction)
LVOT PGmean	LN	20256-4	Mean Gradient	SRT	T-32650	Left Ventricle Outflow Tract	SRT	R-42047	Antegrade Flow (Flow Direction)

LVOT diam	SRT	G-038F	Cardiovascular Orifice Diameter	SRT	T-32650	Left Ventricle Outflow Tract			
IVRT	LN	18071-1	Left Ventricular Isovolumic Relaxation Time	SRT	T-32600	Left Ventricle			
LVET	LN	18041-4	Aortic Valve Ejection Time	SRT	T-35400	Aortic Valve			
LVPEP	LN	59085-1	Pre-Ejection Period	SRT	T-35400	Aortic Valve			
LV Acc Time	LN	20168-1	Acceleration Time	SRT	T-35400	Aortic Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
LVPEP/LVET	LN	59088-5	Pre-Ejection Period/Ejection Time Ratio	SRT	T-35400	Aortic Valve			
dPEP	99IM PRO OF	dPEP	Interventricular Delay	99IM PRO OF	IVD	Interventricular Dyssynchrony			
LVOT Vmax/AV Vmax	99IM PRO OF	LVOTTOAoV_VP	Maximum velocity in the LVOT to Maximum velocity in AV	SRT	T-35400	Aortic Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
LVOT VTI/AV VTI	99IM PRO OF	LVOTTOAoV_VTI	VTI in the LVOT to VTI in AV	SRT	T-35400	Aortic Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
LV Acc T/LVET	SRT	G-0382	Ratio of Aortic Valve Acceleration Time to Ejection Time	SRT	T-35400	Aortic Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
LVOT SV	SRT	F-32120	Stroke Volume	SRT	T-32650	Left Ventricle Outflow Tract	DCM	125219	Doppler Volume Flow (Measurement Method)

LVOT CO	SRT	F-32100	Cardiac Output	SRT	T-32650	Left Ventricle Outflow Tract	DCM	125219	Doppler Volume Flow (Measurement Method)
LVOT SI	SRT	F-00078	Stroke Index	SRT	T-32650	Left Ventricle Outflow Tract	DCM	125219	Doppler Volume Flow (Measurement Method)
LVOT CI	SRT	F-32110	Cardiac Index	SRT	T-32650	Left Ventricle Outflow Tract	DCM	125219	Doppler Volume Flow (Measurement Method)
HR	LN	8867-4	Heart rate	SRT	T-35400	Aortic Valve			

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
2D	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
LVOT diam	SRT	G-038F	Cardiovascular Orifice Diameter	SRT	T-32650	Left Ventricle Outflow Tract			
AVA (Trace)	SRT	G-038E	Cardiovascular Orifice Area	SRT	T-35400	Aortic Valve	DCM	125220	Planimetry (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
AR VC diam	99IM PRO OF	CA_DIAM_VC	Vena Contracta Diameter	SRT	T-35400	Aortic Valve	SRT	G-03A2	2D mode (Image Mode)
Jet Area	LN	59123-0	Jet Area	SRT	T-35400	Aortic Valve	DCM	125220	Planimetry (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
Jet Height	99IM PRO OF	HEIGHT_JET	Jet Height	SRT	T-35400	Aortic Valve	SRT	G-03A2	2D mode (Image Mode)
Jet Width	99IM PRO OF	CA_AR_WIDTH_JET	Jet Width	SRT	T-35400	Aortic Valve	SRT	G-03A2	2D mode (Image Mode)

Jet Width / LVOT diam	99IM PRO OF	CA_AR_WJET2LV OTD	Jet Width to LVOT Diam	SRT	T-35400	Aortic Valve	SRT	G-03A2	2D mode (Image Mode)
-----------------------	-------------	-------------------	------------------------	-----	---------	--------------	-----	--------	----------------------

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
AVA (VTI)	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
AV VTI	LN	20354-7	Velocity Time Integral	SRT	T-35400	Aortic Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
LVOT VTI	LN	20354-7	Velocity Time Integral	SRT	T-32650	Left Ventricle Outflow Tract	SRT	R-42047	Antegrade Flow (Flow Direction)
LVOT diam	SRT	G-038F	Cardiovascular Orifice Diameter	SRT	T-32650	Left Ventricle Outflow Tract			
AVA (VTI)	SRT	G-038E	Cardiovascular Orifice Area	SRT	T-35400	Aortic Valve	DCM	125215	Continuity Equation by Velocity Time Integral (Measurement Method)
							SRT	R-42047	Antegrade Flow (Flow Direction)
AVAI (VTI)	SRT	G-038E	Cardiovascular Orifice Area	SRT	T-35400	Aortic Valve	DCM	125215	Continuity Equation by Velocity Time Integral (Measurement Method)
							SRT	R-42047	Antegrade Flow (Flow Direction)
							LN	8277-6	Body Surface Area (Index)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
------------------------------	-----------------------------	--	--	--------------------------------------	--	--	---	--	--

AVA (Vmax)	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
AV Vmax	LN	11726-7	Peak Velocity	SRT	T-35400	Aortic Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
LVOT Vmax	LN	11726-7	Peak Velocity	SRT	T-32650	Left Ventricle Outflow Tract	SRT	R-42047	Antegrade Flow (Flow Direction)
LVOT diam	SRT	G-038F	Cardiovascular Orifice Diameter	SRT	T-32650	Left Ventricle Outflow Tract			
AVA (Vmax)	SRT	G-038E	Cardiovascular Orifice Area	SRT	T-35400	Aortic Valve	DCM	125214	Continuity Equation by Peak Velocity (Measurement Method)
							SRT	R-42047	Antegrade Flow (Flow Direction)
AVAI (Vmax)	SRT	G-038E	Cardiovascular Orifice Area	SRT	T-35400	Aortic Valve	DCM	125214	Continuity Equation by Peak Velocity (Measurement Method)
							SRT	R-42047	Antegrade Flow (Flow Direction)
							LN	8277-6	Body Surface Area (Index)

Tricuspid Valve

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
TV	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
E pt	LN	18031-5	Tricuspid Valve E-Wave Peak Velocity	SRT	T-35100	Tricuspid valve	SRT	R-42047	Antegrade Flow (Flow Direction)

A pt	LN	18030-7	Tricuspid Valve A-Wave Peak Velocity	SRT	T-35100	Tricuspid valve	SRT	R-42047	Antegrade Flow (Flow Direction)
Dec. Slope	LN	20216-8	Deceleration Slope	SRT	T-35100	Tricuspid valve	SRT	R-42047	Antegrade Flow (Flow Direction)
Dec. Time	LN	20217-6	Deceleration Time	SRT	T-35100	Tricuspid valve	SRT	F-32010	Diastole (Cardiac Cycle Point)
PHT	LN	20280-4	Pressure Half-Time	SRT	T-35100	Tricuspid valve	SRT	R-42047	Antegrade Flow (Flow Direction)
TVA (PHT)	SRT	G-038E	Cardiovascular Orifice Area	SRT	T-35100	Tricuspid valve	DCM	125210	Area by Pressure Half-Time (Measurement Method)
							SRT	R-42047	Antegrade Flow (Flow Direction)
E/A	LN	18039-8	Tricuspid Valve E to A Ratio	SRT	T-35100	Tricuspid valve	SRT	R-42047	Antegrade Flow (Flow Direction)
A/E	99IM PRO OF	CA_TV_A 2E	Tricuspid Valve A to E Ratio	SRT	T-35100	Tricuspid valve	SRT	R-42047	Antegrade Flow (Flow Direction)
TV VTI	LN	20354-7	Velocity Time Integral	SRT	T-35100	Tricuspid valve	SRT	R-42047	Antegrade Flow (Flow Direction)
Vmax	LN	11726-7	Peak Velocity	SRT	T-35100	Tricuspid valve	SRT	R-42047	Antegrade Flow (Flow Direction)
Vmean	LN	20352-1	Time averaged mean velocity	SRT	T-35100	Tricuspid valve	SRT	R-42047	Antegrade Flow (Flow Direction)
PGmax	LN	20247-3	Peak Gradient	SRT	T-35100	Tricuspid valve	SRT	R-42047	Antegrade Flow (Flow Direction)
PGmean	LN	20256-4	Mean Gradient	SRT	T-35100	Tricuspid valve	SRT	R-42047	Antegrade Flow (Flow Direction)
Q_to_TVopen	LN	20296-0	Time from Q wave to Tricuspid Valve Opens	SRT	T-35100	Tricuspid valve	SRT	R-42047	Antegrade Flow (Flow Direction)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
RV MPI	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
RVET	LN	18042-2	Pulmonic Valve Ejection Time	SRT	T-35200	Pulmonary Valve			
TV C-O dur	SRT	G-0389	Tricuspid Valve Closure to Opening Time	SRT	T-35100	Tricuspid valve	SRT	R-42047	Antegrade Flow (Flow Direction)
RV MPI	SRT	G-0381	Right Ventricular Index of Myocardial Performance	SRT	T-32500	Right Ventricle			

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
Annulus	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
Ann diam comm	99IM PRO OF	TV_Dia_C C	Tricuspid Valve Commissural Diameter	SRT	T-35111	Tricuspid Annulus	SRT	F-32020	Systole (Cardiac Cycle Point)
							SRT	G-03A2	2D mode (Image Mode)
Ann diam ant-post	99IM PRO OF	TV_Dia_A P	Tricuspid Valve Anterior-Posterior Annulus Diameter	SRT	T-35111	Tricuspid Annulus	SRT	F-32020	Systole (Cardiac Cycle Point)
							SRT	G-03A2	2D mode (Image Mode)
E' (m)	99IM PRO OF	RV_Ea	Right Ventricular Peak Early Diastolic Tissue	SRT	T-32500	Right Ventricle	SRT	P5-B0128	Tissue Doppler Imaging (Measurement Method)

			Velocity				99IM PRO OF	MedTriAn n	Medial Tricuspid Annulus (Finding Site)
A' (m)	99IM PRO OF	RV_Aa	RV Peak Diastolic Tissue Velocity During Atrial Systole	SRT	T-32500	Right Ventricle	SRT	P5- B0128	Tissue Doppler Imaging (Measureme nt Method)
							99IM PRO OF	MedTriAn n	Medial Tricuspid Annulus (Finding Site)
S' (m)	99IM PRO OF	RV_Sa	Right Ventricular Peak Systolic Tissue Velocity	SRT	T-32500	Right Ventricle	SRT	P5- B0128	Tissue Doppler Imaging (Measureme nt Method)
							99IM PRO OF	MedTriAn n	Medial Tricuspid Annulus (Finding Site)
E' (l)	99IM PRO OF	RV_Ea	Right Ventricular Peak Early Diastolic Tissue Velocity	SRT	T-32500	Right Ventricle	SRT	P5- B0128	Tissue Doppler Imaging (Measureme nt Method)
							99IM PRO OF	LatTriAnn	Lateral Tricuspid Annulus (Finding Site)
A' (l)	99IM PRO OF	RV_Aa	RV Peak Diastolic Tissue Velocity During Atrial Systole	SRT	T-32500	Right Ventricle	SRT	P5- B0128	Tissue Doppler Imaging (Measureme nt Method)
							99IM PRO OF	LatTriAnn	Lateral Tricuspid Annulus (Finding Site)
S' (l)	99IM PRO OF	RV_Sa	Right Ventricular Peak Systolic Tissue	SRT	T-32500	Right Ventricle	SRT	P5- B0128	Tissue Doppler Imaging (Measureme nt Method)

			Velocity				99IM PRO OF	LatTriAnn	Lateral Tricuspid Annulus (Finding Site)
--	--	--	----------	--	--	--	-------------------	-----------	--

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
TR	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
TR VTI	LN	20354-7	Velocity Time Integral	SRT	T-35100	Tricuspid valve	SRT	G-0367	Regurgitant Flow (Flow Direction)
TR Vmax	LN	11726-7	Peak Velocity	SRT	T-35100	Tricuspid valve	SRT	G-0367	Regurgitant Flow (Flow Direction)
PISA Radius	LN	59102-4	Flow Radius	SRT	T-35100	Tricuspid valve	DCM	125216	Proximal Isovelocity Surface Area (Measureme nt Method)
							SRT	G-0367	Regurgitant Flow (Flow Direction)
Aliasing Vel	LN	59130-5	Alias Velocity	SRT	T-35100	Tricuspid valve	DCM	125216	Proximal Isovelocity Surface Area (Measureme nt Method)
							SRT	G-0367	Regurgitant Flow (Flow Direction)
TR Vmean	LN	20352-1	Time averaged mean velocity	SRT	T-35100	Tricuspid valve	SRT	G-0367	Regurgitant Flow (Flow Direction)
TR PGmax	LN	20247-3	Peak Gradient	SRT	T-35100	Tricuspid valve	SRT	G-0367	Regurgitant Flow (Flow Direction)
TR PGmean	LN	20256-4	Mean Gradient	SRT	T-35100	Tricuspid valve	SRT	G-0367	Regurgitant Flow (Flow Direction)
dP/dt	LN	18034-9	Tricuspid Regurgitati on dP/dt	SRT	T-35100	Tricuspid valve	SRT	G-0367	Regurgitant Flow (Flow Direction)

dt	LN	20217-6	Deceleration Time	SRT	T-35100	Tricuspid valve	SRT	G-0367	Regurgitant Flow (Flow Direction)
RAP	LN	18070-3	Right Atrium Systolic Pressure	SRT	T-32200	Right Atrium			
RVSP	SRT	G-0380	Right Ventricular Peak Systolic Pressure	SRT	T-32500	Right Ventricle			
TR VC diam	99IM PRO OF	CA_DIAM_VC	Vena Contracta Diameter	SRT	T-35100	Tricuspid valve	SRT	G-03A2	2D mode (Image Mode)
Jet Area	LN	59123-0	Jet Area	SRT	T-35100	Tricuspid valve	DCM	125220	Planimetry (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
Jet Height	99IM PRO OF	HEIGHT_JET	Jet Height	SRT	T-35100	Tricuspid valve	SRT	G-03A2	2D mode (Image Mode)
PISA Flow Rate	LN	34141-2	Peak Instantaneous Flow Rate	SRT	T-35100	Tricuspid valve	DCM	125216	Proximal Isovelocity Surface Area (Measurement Method)
							SRT	G-0367	Regurgitant Flow (Flow Direction)
PISA ERO Area	SRT	G-038E	Cardiovascular Orifice Area	SRT	T-35100	Tricuspid valve	DCM	125216	Proximal Isovelocity Surface Area (Measurement Method)
							SRT	G-0367	Regurgitant Flow (Flow Direction)
PISA Reg. Volume	LN	33878-0	Volume Flow	SRT	T-35100	Tricuspid valve	DCM	125216	Proximal Isovelocity Surface Area (Measurement Method)
							SRT	G-0367	Regurgitant Flow (Flow Direction)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
TVA (VTI)	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
TV VTI	LN	20354-7	Velocity Time Integral	SRT	T-35100	Tricuspid valve	SRT	R-42047	Antegrade Flow (Flow Direction)
RVOT VTI	LN	20354-7	Velocity Time Integral	SRT	T-32550	Right Ventricle Outflow Tract	SRT	R-42047	Antegrade Flow (Flow Direction)
RVOT diam	SRT	G-038F	Cardiovascular Orifice Diameter	SRT	T-32550	Right Ventricle Outflow Tract			
TVA (VTI)	SRT	G-038E	Cardiovascular Orifice Area	SRT	T-35100	Tricuspid valve	DCM	125215	Continuity Equation by Velocity Time Integral (Measurement Method)
							SRT	R-42047	Antegrade Flow (Flow Direction)
TVAI (VTI)	SRT	G-038E	Cardiovascular Orifice Area	SRT	T-35100	Tricuspid valve	DCM	125214	Continuity Equation by Peak Velocity (Measurement Method)
							SRT	R-42047	Antegrade Flow (Flow Direction)
							LN	8277-6	Body Surface Area (Index)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
TVA (Vmax)	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
Vmax	LN	11726-7	Peak Velocity	SRT	T-35100	Tricuspid valve	SRT	R-42047	Antegrade Flow (Flow Direction)

RVOT Vmax	LN	11726-7	Peak Velocity	SRT	T-32550	Right Ventricle Outflow Tract	SRT	R-42047	Antegrade Flow (Flow Direction)
RVOT diam	SRT	G-038F	Cardiovascular Orifice Diameter	SRT	T-32550	Right Ventricle Outflow Tract			
TVA (Vmax)	SRT	G-038E	Cardiovascular Orifice Area	SRT	T-35100	Tricuspid valve	DCM	125214	Continuity Equation by Peak Velocity (Measurement Method)
							SRT	R-42047	Antegrade Flow (Flow Direction)
TVAI (Vmax)	SRT	G-038E	Cardiovascular Orifice Area	SRT	T-35100	Tricuspid valve	DCM	125214	Continuity Equation by Peak Velocity (Measurement Method)
							SRT	R-42047	Antegrade Flow (Flow Direction)
							LN	8277-6	Body Surface Area (Index)

Pulmonic Valve

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
PV	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
PV VTI	LN	20354-7	Velocity Time Integral	SRT	T-35200	Pulmonary Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
PV Vmax	LN	11726-7	Peak Velocity	SRT	T-35200	Pulmonary Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
PV PGmax	LN	20247-3	Peak Gradient	SRT	T-35200	Pulmonary Valve	SRT	R-42047	Antegrade Flow (Flow Direction)

PV Vmean	LN	20352-1	Time averaged mean velocity	SRT	T-35200	Pulmonary Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
PV PGmean	LN	20256-4	Mean Gradient	SRT	T-35200	Pulmonary Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
RVET	LN	18042-2	Pulmonic Valve Ejection Time	SRT	T-35200	Pulmonary Valve			
RVPEP	LN	59085-1	Pre-Ejection Period	SRT	T-32500	Right Ventricle			
RV Acc Time	LN	20168-1	Acceleration Time	SRT	T-35200	Pulmonary Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
PA Acc Time	LN	20168-1	Acceleration Time	SRT	T-44000	Pulmonary Artery	SRT	R-42047	Antegrade Flow (Flow Direction)
PA PHT	LN	20280-4	Pressure Half-Time	SRT	T-44000	Pulmonary Artery	SRT	R-42047	Antegrade Flow (Flow Direction)
PA Acc Slope	LN	20167-3	Acceleration Slope	SRT	T-44000	Pulmonary Artery	SRT	R-42047	Antegrade Flow (Flow Direction)
Q_to_PVclose	LN	20295-2	Time from Q wave to Pulmonic Valve Closes	SRT	T-35200	Pulmonary Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
RVOT diam	SRT	G-038F	Cardiovascular Orifice Diameter	SRT	T-32550	Right Ventricle Outflow Tract			
RVPEP/RVET	LN	59088-5	Pre-Ejection Period/Ejection Time Ratio	SRT	T-35200	Pulmonary Valve			
RV Acc T/ET	SRT	G-0388	Ratio of Pulmonic Valve Acceleration Time to Ejection Time	SRT	T-35200	Pulmonary Valve	SRT	R-42047	Antegrade Flow (Flow Direction)

dPEP	99IM PRO OF	dPEP	Interventric ular Delay	99IM PRO OF	IVD	Interventric ular Dyssynchr ony			
HR	LN	8867-4	Heart rate	SRT	T-35200	Pulmonary Valve			

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
RVOT	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
RVOT VTI	LN	20354-7	Velocity Time Integral	SRT	T-32550	Right Ventricle Outflow Tract	SRT	R-42047	Antegrade Flow (Flow Direction)
RVOT Vmax	LN	11726-7	Peak Velocity	SRT	T-32550	Right Ventricle Outflow Tract	SRT	R-42047	Antegrade Flow (Flow Direction)
RVOT PGmax	LN	20247-3	Peak Gradient	SRT	T-32550	Right Ventricle Outflow Tract	SRT	R-42047	Antegrade Flow (Flow Direction)
RVOT Vmean	LN	20352-1	Time averaged mean velocity	SRT	T-32550	Right Ventricle Outflow Tract	SRT	R-42047	Antegrade Flow (Flow Direction)
RVOT PGmean	LN	20256-4	Mean Gradient	SRT	T-32550	Right Ventricle Outflow Tract	SRT	R-42047	Antegrade Flow (Flow Direction)
RVOT diam	SRT	G-038F	Cardiovasc ular Orifice Diameter	SRT	T-32550	Right Ventricle Outflow Tract			
HR	LN	8867-4	Heart rate	SRT	T-35200	Pulmonary Valve			
RVOT SV	SRT	F-32120	Stroke Volume	SRT	T-32550	Right Ventricle Outflow Tract	DCM	125219	Doppler Volume Flow (Measureme nt Method)
RVOT CO	SRT	F-32100	Cardiac Output	SRT	T-32550	Right Ventricle Outflow Tract	DCM	125219	Doppler Volume Flow (Measureme nt Method)

RVOT SI	SRT	F-00078	Stroke Index	SRT	T-32550	Right Ventricle Outflow Tract	DCM	125219	Doppler Volume Flow (Measurement Method)
RVOT CI	SRT	F-32110	Cardiac Index	SRT	T-32550	Right Ventricle Outflow Tract	DCM	125219	Doppler Volume Flow (Measurement Method)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
PV Vmax	LN	11726-7	Peak Velocity	SRT	T-35200	Pulmonary Valve	SRT	R-42047	Antegrade Flow (Flow Direction)
RVOT Vmax	LN	11726-7	Peak Velocity	SRT	T-32550	Right Ventricle Outflow Tract	SRT	R-42047	Antegrade Flow (Flow Direction)
RVOT diam	SRT	G-038F	Cardiovascular Orifice Diameter	SRT	T-32550	Right Ventricle Outflow Tract			
PVA (Vmax)	SRT	G-038E	Cardiovascular Orifice Area	SRT	T-35200	Pulmonary Valve	DCM	125214	Continuity Equation by Peak Velocity (Measurement Method)
							SRT	R-42047	Antegrade Flow (Flow Direction)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
PR VTI	LN	20354-7	Velocity Time Integral	SRT	T-35200	Pulmonary Valve	SRT	G-0367	Regurgitant Flow (Flow Direction)
PR Vmax	LN	11726-7	Peak Velocity	SRT	T-35200	Pulmonary Valve	SRT	G-0367	Regurgitant Flow (Flow Direction)

PR Vmean	LN	20352-1	Time averaged mean velocity	SRT	T-35200	Pulmonary Valve	SRT	G-0367	Regurgitant Flow (Flow Direction)
PR PGmax	LN	20247-3	Peak Gradient	SRT	T-35200	Pulmonary Valve	SRT	G-0367	Regurgitant Flow (Flow Direction)
PR PGmean	LN	20256-4	Mean Gradient	SRT	T-35200	Pulmonary Valve	SRT	G-0367	Regurgitant Flow (Flow Direction)
PR Ved	LN	11653-3	End Diastolic Velocity	SRT	T-35200	Pulmonary Valve	SRT	G-0367	Regurgitant Flow (Flow Direction)
RAP	LN	18070-3	Right Atrium Systolic Pressure	SRT	T-32200	Right Atrium			
dPAP	SRT	F-0212C	Pulmonary Artery Pressure	SRT	T-44000	Pulmonary Artery	SRT	F-32011	End Diastole (Cardiac Cycle Point)
PR VC diam	99IM PRO OF	CA_DIAM_VC	Vena Contracta Diameter	SRT	T-35200	Pulmonary Valve	SRT	G-03A2	2D mode (Image Mode)
Jet Area	LN	59123-0	Jet Area	SRT	T-35200	Pulmonary Valve	DCM	125220	Planimetry (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
Jet Height	99IM PRO OF	HEIGHT_JET	Jet Height	SRT	T-35200	Pulmonary Valve	SRT	G-03A2	2D mode (Image Mode)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
PISA (PR)	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
Radius	LN	59102-4	Flow Radius	SRT	T-35200	Pulmonary Valve	DCM	125216	Proximal Isovelocity Surface Area (Measurement Method)
							SRT	G-0367	Regurgitant Flow (Flow Direction)

PR VTI	LN	20354-7	Velocity Time Integral	SRT	T-35200	Pulmonary Valve	DCM	125216	Proximal Isovelocity Surface Area (Measurement Method)
							SRT	G-0367	Regurgitant Flow (Flow Direction)
PR Vmax	LN	11726-7	Peak Velocity	SRT	T-35200	Pulmonary Valve	DCM	125216	Proximal Isovelocity Surface Area (Measurement Method)
							SRT	G-0367	Regurgitant Flow (Flow Direction)
Aliasing Vel	LN	59130-5	Alias Velocity	SRT	T-35200	Pulmonary Valve	DCM	125216	Proximal Isovelocity Surface Area (Measurement Method)
							SRT	G-0367	Regurgitant Flow (Flow Direction)
PISA Flow Rate	LN	34141-2	Peak Instantaneous Flow Rate	SRT	T-35200	Pulmonary Valve	DCM	125216	Proximal Isovelocity Surface Area (Measurement Method)
							SRT	G-0367	Regurgitant Flow (Flow Direction)
PISA ERO Area	SRT	G-038E	Cardiovascular Orifice Area	SRT	T-35200	Pulmonary Valve	DCM	125216	Proximal Isovelocity Surface Area (Measurement Method)
							SRT	G-0367	Regurgitant Flow (Flow Direction)
PISA Reg. Volume	LN	33878-0	Volume Flow	SRT	T-35200	Pulmonary Valve	DCM	125216	Proximal Isovelocity Surface Area (Measurement Method)
							SRT	G-0367	Regurgitant Flow (Flow Direction)

PISA Reg. Fraction	SRT	G-0390	Regurgitant Fraction	SRT	T-35200	Pulmonary Valve	DCM	125216	Proximal Isovelocity Surface Area (Measurement Method)
							SRT	G-0367	Regurgitant Flow (Flow Direction)

Pulmonary Artery

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
PA	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
PA diam	LN	18020-8	Main Pulmonary Artery Diameter	SRT	T-44000	Pulmonary Artery	SRT	F-32020	Systole (Cardiac Cycle Point)
							SRT	G-03A2	2D mode (Image Mode)
PA Left diam	LN	18019-0	Left Pulmonary Artery Diameter	SRT	T-44000	Pulmonary Artery	SRT	F-32020	Systole (Cardiac Cycle Point)
							SRT	G-03A2	2D mode (Image Mode)
PA Right diam	LN	18021-6	Right Pulmonary Artery Diameter	SRT	T-44000	Pulmonary Artery	SRT	F-32020	Systole (Cardiac Cycle Point)
							SRT	G-03A2	2D mode (Image Mode)
RAP	LN	18070-3	Right Atrium Systolic Pressure	SRT	T-32200	Right Atrium			
TR Vmax	LN	11726-7	Peak Velocity	SRT	T-35100	Tricuspid valve	SRT	G-0367	Regurgitant Flow (Flow Direction)
PR Ved	LN	11653-3	End Diastolic Velocity	SRT	T-35200	Pulmonary Valve	SRT	G-0367	Regurgitant Flow (Flow Direction)

RVSP	SRT	G-0380	Right Ventricular Peak Systolic Pressure	SRT	T-32500	Right Ventricle			
sPAP	SRT	F-0212C	Pulmonary Artery Pressure	SRT	T-44000	Pulmonary Artery	SRT	F-32020	Systole (Cardiac Cycle Point)
dPAP	SRT	F-0212C	Pulmonary Artery Pressure	SRT	T-44000	Pulmonary Artery	SRT	F-32011	End Diastole (Cardiac Cycle Point)
mPAP	99IM PRO OF	MPAP	Mean Pulmonary Artery Pressure	SRT	T-44000	Pulmonary Artery			

Pulmonary Vein

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
PVE	CS D	CV	CM	CSD	CV	CM	CSD	CV	CM
S1 Vel	LN	29450-4	Pulmonary Vein Systolic Peak Velocity	SRT	T-48500	Pulmonary Vein	99IM PRO OF	Peak1	early systolic Peak S1 (Finding Site)
							SRT	R-42047	Antegrade Flow (Flow Direction)
D Vel	LN	29451-2	Pulmonary Vein Diastolic Peak Velocity	SRT	T-48500	Pulmonary Vein	SRT	R-42047	Antegrade Flow (Flow Direction)
A Vel	LN	29453-8	Pulmonary Vein Atrial Contractio n Reversal Peak Velocity	SRT	T-48500	Pulmonary Vein	SRT	F-32010	Diastole (Cardiac Cycle Point)
							SRT	R-42047	Antegrade Flow (Flow Direction)
A Dur	SRT	G-038B	Pulmonary Vein A- Wave Duration	SRT	T-48500	Pulmonary Vein	SRT	R-42047	Antegrade Flow (Flow Direction)

S VTI	SRT	G-038C	Pulmonary Vein S-Wave Velocity Time Integral	SRT	T-48500	Pulmonary Vein	SRT	R-42047	Antegrade Flow (Flow Direction)
D VTI	SRT	G-038D	Pulmonary Vein D-Wave Velocity Time Integral	SRT	T-48500	Pulmonary Vein	SRT	R-42047	Antegrade Flow (Flow Direction)
A VTI	99IM PRO OF	CA_VTI_A Wave_PVE	Pulmonary Vein A-Wave Velocity Time Integral	SRT	T-48500	Pulmonary Vein	SRT	R-42047	Antegrade Flow (Flow Direction)
S2 Vel	LN	29450-4	Pulmonary Vein Systolic Peak Velocity	SRT	T-48500	Pulmonary Vein	99IM PRO OF	Peak2	late systolic Peak S2 (Finding Site)
							SRT	R-42047	Antegrade Flow (Flow Direction)
S2/D	LN	29452-0	Pulmonary Vein Systolic to Diastolic Ratio	SRT	T-48500	Pulmonary Vein	SRT	R-42047	Antegrade Flow (Flow Direction)
D Dec. Time	LN	20217-6	Deceleration Time	SRT	T-48500	Pulmonary Vein	SRT	F-32010	Diastole (Cardiac Cycle Point)
							SRT	R-42047	Antegrade Flow (Flow Direction)
Sys Fraction	99IM PRO OF	CA_PVE_SF	Systolic Fraction	SRT	T-48500	Pulmonary Vein			
Diam Diastole	SRT	M-02550	Diameter	SRT	T-48500	Pulmonary Vein	SRT	F-32010	Diastole (Cardiac Cycle Point)
							SRT	G-03A2	2D mode (Image Mode)
Diam Systole	SRT	M-02550	Diameter	SRT	T-48500	Pulmonary Vein	SRT	F-32020	Systole (Cardiac Cycle Point)

							SRT	G-03A2	2D mode (Image Mode)
--	--	--	--	--	--	--	-----	--------	----------------------------

Right Ventricle

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
RV	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
RV WT	99IM PRO OF	RVWT	Right Ventricular Wall Thickness	SRT	T-32500	Right Ventricle	SRT	F-32011	End Diastole (Cardiac Cycle Point)
							SRT	G-03A2	2D mode (Image Mode)
RVDd base (RVD1)	LN	20304-2	Right Ventricular Internal Diastolic Dimension	SRT	T-32500	Right Ventricle	SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A123	Basal (Finding Site)
RVDd mid (RVD2)	LN	20304-2	Right Ventricular Internal Diastolic Dimension	SRT	T-32500	Right Ventricle	SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A109	Medial (Finding Site)
RVLd (RVD3)	99IM PRO OF	CA_RV_L_ BaseApex	Right Ventricular base-apex Length	SRT	T-32500	Right Ventricle	SRT	F-32011	End Diastole (Cardiac Cycle Point)
							SRT	G-03A2	2D mode (Image Mode)
RVOT diam	SRT	G-038F	Cardiovas cular Orifice Diameter	SRT	T-32550	Right Ventricle Outflow Tract	SRT	F-32011	End Diastole (Cardiac Cycle Point)
							SRT	G-03A2	2D mode (Image Mode)
RVOT VTI	LN	20354-7	Velocity Time Integral	SRT	T-32550	Right Ventricle Outflow Tract	SRT	R-42047	Antegrade Flow (Flow Direction)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
MMode	CS D	CV	CM	CSD	CV	CM	CSD	CV	CM
RV WT	99IM PRO OF	RVWT	Right Ventricular Wall Thickness	SRT	T-32500	Right Ventricle	SRT	F-32011	End Diastole (Cardiac Cycle Point)
							SRT	G-0394	MMode (Image Mode)
TAPSE	99IM PRO OF	TAPSE	Tricuspid Annular Plane Systolic Excursion	SRT	T-32500	Right Ventricle	SRT	G-0394	MMode (Image Mode)
RVIDd	LN	20304-2	Right Ventricular Internal Diastolic Dimension	SRT	T-32500	Right Ventricle	SRT	G-0394	MMode (Image Mode)
RVOT diam	SRT	G-038F	Cardiovascular Orifice Diameter	SRT	T-32550	Right Ventricle Outflow Tract	SRT	F-32011	End Diastole (Cardiac Cycle Point)
							SRT	G-0394	MMode (Image Mode)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
FAC	CS D	CV	CM	CSD	CV	CM	CSD	CV	CM
RV Area ED	99IM PRO OF	AREA_DIA_RV	Right Ventricular Diastolic Area	SRT	T-32500	Right Ventricle	SRT	G-03A2	2D mode (Image Mode)
RV Area ES	99IM PRO OF	AREA_SY S_RV	Right Ventricular Systolic Area	SRT	T-32500	Right Ventricle	SRT	G-03A2	2D mode (Image Mode)
RV FAC	99IM PRO OF	RVFAC	Right Ventricular Fractional Area Change	SRT	T-32500	Right Ventricle	SRT	G-03A2	2D mode (Image Mode)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
RV MPI	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
RVET	LN	18042-2	Pulmonic Valve Ejection Time	SRT	T-35200	Pulmonary Valve			
TV C-O dur	SRT	G-0389	Tricuspid Valve Closure to Opening Time	SRT	T-35100	Tricuspid valve	SRT	R-42047	Antegrade Flow (Flow Direction)
RV MPI	SRT	G-0381	Right Ventricular Index of Myocardial Performance	SRT	T-32500	Right Ventricle			

Venous Return

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
Vena Cava 2D	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
VC Inf diam	LN	18006-7	Inferior Vena Cava Diameter	SRT	T-48710	Inferior Vena Cava	SRT	G-03A2	2D mode (Image Mode)
VC Inf inspir diam	LN	18006-7	Inferior Vena Cava Diameter	SRT	T-48710	Inferior Vena Cava	SRT	F-20010	During Inspiration (Respiratory Cycle Point)
							SRT	G-03A2	2D mode (Image Mode)
VC Inf expir diam	LN	18006-7	Inferior Vena Cava Diameter	SRT	T-48710	Inferior Vena Cava	SRT	F-20020	During Expiration (Respiratory Cycle Point)
							SRT	G-03A2	2D mode (Image Mode)

VC Sup diam	99IM PRO OF	DIAM_VC_ Sup	Superior Vena Cava Diameter	SRT	T-48610	Superior Vena Cava	SRT	G-03A2	2D mode (Image Mode)
VC Sup inspir diam	99IM PRO OF	DIAM_VC_ Sup	Superior Vena Cava Diameter	SRT	T-48610	Superior Vena Cava	SRT	F-20010	During Inspiration(R espiratory Cycle Point)
							SRT	G-03A2	2D mode (Image Mode)
VC Sup expir diam	99IM PRO OF	DIAM_VC_ Sup	Superior Vena Cava Diameter	SRT	T-48610	Superior Vena Cava	SRT	F-20020	During Expiration(R espiratory Cycle Point)
							SRT	G-03A2	2D mode (Image Mode)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
	CS D	CV	CM	CSD	CV	CM	CSD	CV	CM
VC Inf inspir diam	LN	18006-7	Inferior Vena Cava Diameter	SRT	T-48710	Inferior Vena Cava	SRT	F-20010	During Inspiration(R espiratory Cycle Point)
							SRT	G-0394	MMode (Image Mode)
VC Inf expir diam	LN	18006-7	Inferior Vena Cava Diameter	SRT	T-48710	Inferior Vena Cava	SRT	F-20020	During Expiration(R espiratory Cycle Point)
							SRT	G-0394	MMode (Image Mode)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
	CS D	CV	CM	CSD	CV	CM	CSD	CV	CM
Hepatic Vein diam	SRT	M-02550	Diameter	SRT	T-48720	Hepatic Vein	SRT	G-03A2	2D mode (Image Mode)

Hepatic Vein Vmax	LN	11726-7	Peak Velocity	SRT	T-48720	Hepatic Vein			
Hepatic Vein PGmax	LN	20247-3	Peak Gradient	SRT	T-48720	Hepatic Vein			

Left Atrium

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
LA Vol (A-L)	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
LAA (A4C)	LN	17977-0	Left Atrium Systolic Area	SRT	T-32300	Left Atrium	SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19C	Apical four chamber (Image View)
LAA (A2C)	LN	17977-0	Left Atrium Systolic Area	SRT	T-32300	Left Atrium	SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19B	Apical two chamber (Image View)
LA L	SRT	G-A22A	Length	SRT	T-32300	Left Atrium	SRT	F-32020	Systole (Cardiac Cycle Point)
							SRT	G-03A2	2D mode (Image Mode)
LA Vol (A-L)	DCM	122408	Left Atrial End Systolic Volume	SRT	T-32300	Left Atrium	DCM	125204	Area-Length Biplane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
LA Voll (A-L)	DCM	122408	Left Atrial End Systolic Volume	SRT	T-32300	Left Atrium	DCM	125204	Area-Length Biplane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)

							LN	8277-6	Body Surface Area (Index)
--	--	--	--	--	--	--	----	--------	---------------------------

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
LA Vol (Simpson)	CS D	CV	CM	CSD	CV	CM	CSD	CV	CM
LA Vol (A4C)	DCM	122408	Left Atrial End Systolic Volume	SRT	T-32300	Left Atrium	DCM	125208	Method of Disks, Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19C	Apical four chamber (Image View)
LA Vol (A2C)	DCM	122408	Left Atrial End Systolic Volume	SRT	T-32300	Left Atrium	DCM	125208	Method of Disks, Single Plane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19B	Apical two chamber (Image View)
LA Vol (BP)	DCM	122408	Left Atrial End Systolic Volume	SRT	T-32300	Left Atrium	DCM	125207	Method of Disks, Biplane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)
LA Vol (BP)	DCM	122408	Left Atrial End Systolic Volume	SRT	T-32300	Left Atrium	DCM	125207	Method of Disks, Biplane (Measurement Method)
							SRT	G-03A2	2D mode (Image Mode)

							LN	8277-6	Body Surface Area (Index)
LAA (A4C)	LN	17977-0	Left Atrium Systolic Area	SRT	T-32300	Left Atrium	SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19C	Apical four chamber (Image View)
							99IM PRO OF	LA_Vol_SIMP	Measurement Package LA Vol(SIMP) (Measurement Method)
LAA (A2C)	LN	17977-0	Left Atrium Systolic Area	SRT	T-32300	Left Atrium	SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19B	Apical two chamber (Image View)
							99IM PRO OF	LA_Vol_SIMP	Measurement Package LA Vol(SIMP) (Measurement Method)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
LA Vol (3 Axis)	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
LA diam ant-post	LN	29469-4	Left Atrium Antero-posterior Systolic Dimension	SRT	T-32300	Left Atrium	SRT	G-03A2	2D mode (Image Mode)
LA diam sup-inf	99IM PRO OF	DIAM_SUP INF	Left Atrium Superior-inferior Systolic Dimension	SRT	T-32300	Left Atrium	SRT	G-03A2	2D mode (Image Mode)

LA diam mid-lat	99IM PRO OF	DIAM_MID LAT	Left Atrium Medial- lateral Systolic Dimensio n	SRT	T-32300	Left Atrium	SRT	G-03A2	2D mode (Image Mode)
LA Vol	SRT	G-0383	Left Atrium Systolic Volume	SRT	T-32300	Left Atrium	DCM	125211	Biplane Ellipse (Measureme nt Method)
							SRT	G-03A2	2D mode (Image Mode)

Right Atrium

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
	CS D	CV	CM	CSD	CV	CM	CSD	CV	CM
RAA (A4C)	L N	17988 -7	Right Atrium Systoli c Area	S RT	T- 32200	Right Atrium	S RT	G- 03A2	2D mode (Image Mode)
							SRT	G-A19C	Apical four chamber (Image View)
RA diam lax sys	99IM PRO OF	CA_RA_M AJAXES	Right Atrium Systolic Major Axis	SRT	T-32200	Right Atrium	SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19C	Apical four chamber (Image View)
RA diam sax sys	99IM PRO OF	CA_RA_MI NAXES	Right Atrium Systolic Minor Axis	SRT	T-32200	Right Atrium	SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19C	Apical four chamber (Image View)
RA Vol (A-L)	99IM PRO OF	CA_RA_V OLes	Right Atrial End Systolic Volume	SRT	T-32200	Right Atrium	DCM	125205	Area-Length Single Plane (Measureme nt Method)
							SRT	G-03A2	2D mode (Image Mode)

RA VolI (A-L)	99IM PRO OF	CA_RA_VO LEs	Right Atrial End Systolic Volume	SRT	T-32200	Right Atrium	DCM	125205	Area-Length Single Plane (Measureme nt Method)
							SRT	G-03A2	2D mode (Image Mode)
							LN	8277-6	Body Surface Area (Index)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
RA dimension	CS D	CV	CM	CSD	CV	CM	CSD	CV	CM
RA diam lax	SRT	G-A193	Major Axis	SRT	T-32200	Right Atrium	SRT	G-03A2	2D mode (Image Mode)
RA diam sax	SRT	G-A194	Minor Axis	SRT	T-32200	Right Atrium	SRT	G-03A2	2D mode (Image Mode)
RA Area ES	LN	17988-7	Right Atrium Systolic Area	SRT	T-32200	Right Atrium	SRT	G-03A2	2D mode (Image Mode)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
RA Vol (Simpson)	CS D	CV	CM	CSD	CV	CM	CSD	CV	CM
RA Vol (A4C)	99IM PRO OF	CA_RA_VO LEs	Right Atrial End Systolic Volume	SRT	T-32200	Right Atrium	DCM	125208	Method of Disks, Single Plane (Measureme nt Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19C	Apical four chamber (Image View)

RA Voll (A4C)	99IM PRO OF	CA_RA_VO Les	Right Atrial End Systolic Volume	SRT	T-32200	Right Atrium	DCM	125208	Method of Disks, Single Plane (Measureme nt Method)
							SRT	G-03A2	2D mode (Image Mode)
							SRT	G-A19C	Apical four chamber (Image View)
							LN	8277-6	Body Surface Area (Index)

Tissue Doppler

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
MV Annulus (m)	CS D	CV	CM	CSD	CV	CM	CSD	CV	CM
E'(m)	SRT	G-037A	Left Ventricula r Peak Early Diastolic Tissue Velocity	SRT	T-32600	Left Ventricle	SRT	P5-B0128	Tissue Doppler Imaging (Measureme nt Method)
							SRT	G-0391	Medial Mitral Annulus (Finding Site)
A'(m)	SRT	G-037C	LV Peak Diastolic Tissue Velocity During Atrial Systole	SRT	T-32600	Left Ventricle	SRT	P5-B0128	Tissue Doppler Imaging (Measureme nt Method)
							SRT	G-0391	Medial Mitral Annulus (Finding Site)
E'/A'(m)	99IM PRO OF	CA_MV_Ea 2Aa	Mitral Valve Ea to Aa Ratio	SRT	T-32600	Left Ventricle	SRT	P5-B0128	Tissue Doppler Imaging (Measureme nt Method)

							SRT	G-0391	Medial Mitral Annulus (Finding Site)
E/E'(m)	SRT	G-037B	Ratio of MV Peak Velocity to LV Peak Tissue Velocity E-Wave	SRT	T-32600	Left Ventricle	SRT	P5-B0128	Tissue Doppler Imaging (Measurement Method)
							SRT	G-0391	Medial Mitral Annulus (Finding Site)
AT E'(m)	99IM PRO OF	CA_MV_ATa	Mitral Valve Acceleration time of Ea	SRT	T-32600	Left Ventricle	SRT	P5-B0128	Tissue Doppler Imaging (Measurement Method)
							SRT	G-0391	Medial Mitral Annulus (Finding Site)
AR E'(m)	99IM PRO OF	CA_MV_A Ra	Mitral Valve Acceleration rate of Ea	SRT	T-32600	Left Ventricle	SRT	P5-B0128	Tissue Doppler Imaging (Measurement Method)
							SRT	G-0391	Medial Mitral Annulus (Finding Site)
DT E'(m)	99IM PRO OF	CA_MV_D Ta	Mitral Valve Deceleration time of Ea	SRT	T-32600	Left Ventricle	SRT	P5-B0128	Tissue Doppler Imaging (Measurement Method)
							SRT	G-0391	Medial Mitral Annulus (Finding Site)
DR E'(m)	99IM PRO OF	CA_MV_D Ra	Mitral Valve Deceleration rate of Ea	SRT	T-32600	Left Ventricle	SRT	P5-B0128	Tissue Doppler Imaging (Measurement Method)
							SRT	G-0391	Medial Mitral Annulus (Finding Site)

S'(m)	SRT	G-037D	Left Ventricular Peak Systolic Tissue Velocity	SRT	T-32600	Left Ventricle	SRT	P5-B0128	Tissue Doppler Imaging (Measurement Method)
							SRT	G-0391	Medial Mitral Annulus (Finding Site)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
MV Annulus (I)	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
E'(I)	SRT	G-037A	Left Ventricular Peak Early Diastolic Tissue Velocity	SRT	T-32600	Left Ventricle	SRT	P5-B0128	Tissue Doppler Imaging (Measurement Method)
							SRT	G-0392	Lateral Mitral Annulus (Finding Site)
A'(I)	SRT	G-037C	LV Peak Diastolic Tissue Velocity During Atrial Systole	SRT	T-32600	Left Ventricle	SRT	P5-B0128	Tissue Doppler Imaging (Measurement Method)
							SRT	G-0392	Lateral Mitral Annulus (Finding Site)
E'/A'(I)	99IM PROOF	CA_MV_Ea2Aa	Mitral Valve Ea to Aa Ratio	SRT	T-32600	Left Ventricle	SRT	P5-B0128	Tissue Doppler Imaging (Measurement Method)
							SRT	G-0392	Lateral Mitral Annulus (Finding Site)
E/E'(I)	SRT	G-037B	Ratio of MV Peak Velocity to LV Peak Tissue	SRT	T-32600	Left Ventricle	SRT	P5-B0128	Tissue Doppler Imaging (Measurement Method)

			Velocity E-Wave				SRT	G-0392	Lateral Mitral Annulus (Finding Site)
AT E'(I)	99IM PRO OF	CA_MV_A Ta	Mitral Valve Accelerati on time of Ea	SRT	T-32600	Left Ventricle	SRT	P5- B0128	Tissue Doppler Imaging (Measureme nt Method)
							SRT	G-0392	Lateral Mitral Annulus (Finding Site)
AR E'(I)	99IM PRO OF	CA_MV_A Ra	Mitral Valve Accelerati on rate of Ea	SRT	T-32600	Left Ventricle	SRT	P5- B0128	Tissue Doppler Imaging (Measureme nt Method)
							SRT	G-0392	Lateral Mitral Annulus (Finding Site)
DT E'(I)	99IM PRO OF	CA_MV_D Ta	Mitral Valve Deccelera tion time of Ea	SRT	T-32600	Left Ventricle	SRT	P5- B0128	Tissue Doppler Imaging (Measureme nt Method)
							SRT	G-0392	Lateral Mitral Annulus (Finding Site)
DR E'(I)	99IM PRO OF	CA_MV_D Ra	Mitral Valve Deccelera tion rate of Ea	SRT	T-32600	Left Ventricle	SRT	P5- B0128	Tissue Doppler Imaging (Measureme nt Method)
							SRT	G-0392	Lateral Mitral Annulus (Finding Site)
S'(I)	SRT	G-037D	Left Ventricula r Peak Systolic Tissue Velocity	SRT	T-32600	Left Ventricle	SRT	P5- B0128	Tissue Doppler Imaging (Measureme nt Method)
							SRT	G-0392	Lateral Mitral Annulus (Finding Site)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
MV Annulus (s)	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
E'(s)	SRT	G-037A	Left Ventricular Peak Early Diastolic Tissue Velocity	SRT	T-32600	Left Ventricle	SRT	P5-B0128	Tissue Doppler Imaging (Measurement Method)
							99IM Proof	SepMitAnn	Septal Mitral Annulus (Finding Site)
A'(s)	SRT	G-037C	LV Peak Diastolic Tissue Velocity During Atrial Systole	SRT	T-32600	Left Ventricle	SRT	P5-B0128	Tissue Doppler Imaging (Measurement Method)
							99IM Proof	SepMitAnn	Septal Mitral Annulus (Finding Site)
S'(s)	SRT	G-037D	Left Ventricular Peak Systolic Tissue Velocity	SRT	T-32600	Left Ventricle	SRT	P5-B0128	Tissue Doppler Imaging (Measurement Method)
							99IM Proof	SepMitAnn	Septal Mitral Annulus (Finding Site)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
TV Annulus (m)	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
E'(m)	99IM PROOF	RV_Ea	Right Ventricular Peak Early Diastolic Tissue Velocity	SRT	T-32500	Right Ventricle	SRT	P5-B0128	Tissue Doppler Imaging (Measurement Method)
							99IM PROOF	MedTriAnn	Medial Tricuspid Annulus (Finding Site)

A'(m)	99IM PRO OF	RV_Aa	RV Peak Diastolic Tissue Velocity During Atrial Systole	SRT	T-32500	Right Ventricle	SRT	P5- B0128	Tissue Doppler Imaging (Measureme nt Method)
							99IM PRO OF	MedTriA nn	Medial Tricuspid Annulus (Finding Site)
S'(m)	99IM PRO OF	RV_Sa	Right Ventricular Peak Systolic Tissue Velocity	SRT	T-32500	Right Ventricle	SRT	P5- B0128	Tissue Doppler Imaging (Measureme nt Method)
							99IM PRO OF	MedTriA nn	Medial Tricuspid Annulus (Finding Site)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
TV Annulus (l)	CS D	CV	CM	CSD	CV	CM	CSD	CV	CM
E'(l)	99IM PRO OF	RV_Ea	Right Ventricular Peak Early Diastolic Tissue Velocity	SRT	T-32500	Right Ventricle	SRT	P5- B0128	Tissue Doppler Imaging (Measureme nt Method)
							99IM PRO OF	LatTriAn n	Lateral Tricuspid Annulus (Finding Site)
A'(l)	99IM PRO OF	RV_Aa	RV Peak Diastolic Tissue Velocity During Atrial Systole	SRT	T-32500	Right Ventricle	SRT	P5- B0128	Tissue Doppler Imaging (Measureme nt Method)
							99IM PRO OF	LatTriAn n	Lateral Tricuspid Annulus (Finding Site)
S'(l)	99IM PRO OF	RV_Sa	Right Ventricular Peak Systolic Tissue	SRT	T-32500	Right Ventricle	SRT	P5- B0128	Tissue Doppler Imaging (Measureme nt Method)

			Velocity				99IM PRO OF	LatTriAn n	Lateral Tricuspid Annulus (Finding Site)
--	--	--	----------	--	--	--	-------------------	---------------	---

Congenital

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
ASD	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
ASD Vmax	LN	11726-7	Peak Velocity	SRT	D4- 31220	Atrial Septal Defect			

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
VSD	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
VSD Vmax	LN	11726-7	Peak Velocity	SRT	D4- 31150	Ventricular Septal Defect	SRT	R-42047	Antegrade Flow (Flow Direction)
VSD PGmax	LN	20247-3	Peak Gradient	SRT	D4- 31150	Ventricular Septal Defect	SRT	R-42047	Antegrade Flow (Flow Direction)

Label displayed in Worksheet	TID (5203) Echo Measurement			TID (5202) Echo Section Finding Site			Concept or Acquisition Context Modifier		
Qp/Qs	CSD	CV	CM	CSD	CV	CM	CSD	CV	CM
LVOT VTI	LN	20354-7	Velocity Time Integral	SRT	T-32650	Left Ventricle Outflow Tract	SRT	R-42047	Antegrade Flow (Flow Direction)
LVOT diam	SRT	G-038F	Cardiovascular Orifice Diameter	SRT	T-32650	Left Ventricle Outflow Tract			
HR	LN	8867-4	Heart rate	SRT	T-35400	Aortic Valve			
LVOT SV	SRT	F-32120	Stroke Volume	SRT	T-32650	Left Ventricle Outflow Tract	DCM	125219	Doppler Volume Flow (Measurement Method)

LVOT CO	SRT	F-32100	Cardiac Output	SRT	T-32650	Left Ventricle Outflow Tract	DCM	125219	Doppler Volume Flow (Measurement Method)
RVOT VTI	LN	20354-7	Velocity Time Integral	SRT	T-32550	Right Ventricle Outflow Tract	SRT	R-42047	Antegrade Flow (Flow Direction)
RVOT diam	SRT	G-038F	Cardiovascular Orifice Diameter	SRT	T-32550	Right Ventricle Outflow Tract			
HR	LN	8867-4	Heart rate	SRT	T-35200	Pulmonary Valve			
RVOT SV	SRT	F-32120	Stroke Volume	SRT	T-32550	Right Ventricle Outflow Tract	DCM	125219	Doppler Volume Flow (Measurement Method)
RVOT CO	SRT	F-32100	Cardiac Output	SRT	T-32550	Right Ventricle Outflow Tract	DCM	125219	Doppler Volume Flow (Measurement Method)
Qp/Qs	LN	29462-9	Pulmonary-to-Systemic Shunt Flow Ratio	SRT	P5-30031	Cardiac Shunt Study	DCM	125219	Doppler Volume Flow (Measurement Method)
Qp-Qs	99IM PRO OF	CA_CS_P F-SF	Pulmonary minus Systemic Shunt Flow	SRT	P5-30031	Cardiac Shunt Study	DCM	125219	Doppler Volume Flow (Measurement Method)

8.6.3 Context Groups (Non-Advanced Cardiac)

The Context Groups used in Structured Reports are shown below.

Terms in *Italic* are extensions to the Standard Context Groups. In the Context Group extension, the next attributes are always included besides CV/CSD/CM. (Attributes are not included for a term which is not an extension.)

Attributes used in Context Group Extension

Tag	VR	Attribute Name	Value
(0008,0105)	CS	Mapping Resource	"DCMR"
(0008,0106)	DT	Context Group Version	Described in PS 3.16
(0008,0107)	DT	Context Group Local Version	"YYYYMMDD" Described in each of Context Groups below
(0008,010B)	CS	Context Group Extension Flag	"Y"

Tag	VR	Attribute Name	Value
(0008,010D)	UI	Context Group Extension Creator UID	1.2.392.200039.106
(0008,010F)	CS	Context Identifier	Identifies Context Group to which a term is added

CID 42 Numeric Value Qualifier

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	DCM	114007	Measurement not attempted
-	DCM	114009	Value out of range

CID 223 Normal Range Values

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	SRT	R-0038B	Normal Range Upper Limit
-	SRT	R-10041	Normal Range Lower Limit

CID 227 Sample Statistical Descriptors

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	DCM	121416	Z-Score of measurement

CID 228 Equation or Table

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	DCM	121421	Equation Citation
-	DCM	121422	Table of Values Citation

CID 244 Laterality

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
Right	SRT	G-A100	Right

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
Left	SRT	G-A101	Left

CID 271 Observation Subject Class

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	DCM	121025	Patient

CID 3627 Measurement Type

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	SRT	R-00317	Mean
-	SRT	R-41D2D	Calculated
-	SRT	G-A437	Maximum

CID 3663 Body Surface Area Equations

(Local Version: 20060807)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
BSA Equation, DuBois	DCM	122241	$BSA = 0.007184 * WT^{0.425} * HT^{0.725}$
<i>BSA Equation, Boyd</i>	99ALOKA	A3663-001	$BSA = 0.0003207 * WT^{(0.7285 - 0.0188 \log(WT))} * HT^{0.3}$
<i>BSA Equation, Shintani</i>	99ALOKA	A3663-002	$BSA = 0.007358 * HT^{0.725} * WT^{0.425}$

CID 7455 Sex

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
Male	DCM	M	Male
Female	DCM	F	Female
Other	DCM	121103	Undetermined sex

CID 7456 Units of Measure for Age

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
- (Y)	UCUM	a	year
- (M)	UCUM	mo	month
- (W)	UCUM	wk	week
- (D)	UCUM	d	day

CID 12003 OB-GYN Dates

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
LMP-EDC	LN	11779-6	EDD from LMP
Composite US-EDC	LN	11781-2	EDD from average ultrasound age
BBT-EDC	LN	11780-4	EDD from ovulation date
LMP	LN	11955-2	LMP
BBT	LN	11976-8	Ovulation date

CID 12004 Fetal Biometry Ratios

(Local Version: 20051202)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
HC/AC	LN	11947-9	HC/AC
FL/AC	LN	11871-1	FL/AC
FL/BPD	LN	11872-9	FL/BPD
CI (BPD/OFD)	LN	11823-2	Cephalic Index
FL/HC	LN	11873-7	FL/HC
CI (BPD _o /OFD _o)	99ALOKA	A12004-001	Cephalic Index(BPD _o /OFD _o)
LVW/HW	99ALOKA	A12004-002	LVW/HW

CID 12005 Fetal Biometry Measurements

(Local Version: 20171130 , 20051202 , 20150515)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
AC	LN	11979-2	Abdominal Circumference
APD	LN	11818-2	Anterior-Posterior Abdominal Diameter

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
APTD	LN	11818-2	Anterior-Posterior Abdominal Diameter
BPD	LN	11820-8	Biparietal Diameter
HC	LN	11984-2	Head Circumference
OFD	LN	11851-3	Occipital-Frontal Diameter
TC	LN	11988-3	Thoracic Circumference
TAD	LN	11862-0	Tranverse Abdominal Diameter
TTD	LN	11862-0	Tranverse Abdominal Diameter
AXT	LN	33191-8	APAD * TAD
<i>BD</i>	<i>99ALOKA</i>	<i>A12005-001</i>	<i>Binocular Distance</i>
<i>BPD_o</i>	<i>99ALOKA</i>	<i>A12005-002</i>	<i>Biparietal Diameter outer-to-outer</i>
<i>FTA</i>	<i>99ALOKA</i>	<i>A12005-003</i>	<i>Fetal Trunk Cross Sectional Area</i>
<i>LV</i>	<i>99ALOKA</i>	<i>A12005-004</i>	<i>Length of Vertebrae</i>
<i>OFD_o</i>	<i>99ALOKA</i>	<i>A12005-005</i>	<i>Occipital-Frontal Diameter outer-to-outer</i>
<i>AD</i>	<i>99ALOKA</i>	<i>A12005-006</i>	<i>Abdominal Diameter</i>
<i>TL</i>	<i>99ALOKA</i>	<i>A12005-007</i>	<i>Thoracic Length</i>
<i>HC2</i>	<i>99ALOKA</i>	<i>A12005-008</i>	<i>Head Circumference for Merz, Hansmann</i>
<i>AF Pocket</i>	<i>99ALOKA</i>	<i>A12005-009</i>	<i>Amniotic Fluid Volume</i>
<i>AFV</i>	<i>99ALOKA</i>	<i>A12005-009</i>	<i>Amniotic Fluid Volume</i>
<i>HC3</i>	<i>99ALOKA</i>	<i>A12005-010</i>	<i>Head Circumference for BMUS</i>
<i>AC2</i>	<i>99ALOKA</i>	<i>A12005-011</i>	<i>Abdominal Circumference for BMUS</i>
<i>MVP</i>	<i>99ALOKA</i>	<i>A12005-012</i>	<i>Maximum Vertical Pocket</i>

CID 12006 Fetal Long Bones Biometry Measurements

(Local Version: 20051202)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
HL	LN	11966-9	Humerus length
RAD	LN	11967-7	Radius length
ULNA	LN	11969-3	Ulna length
TIB	LN	11968-5	Tibia length
FIB	LN	11964-4	Fibula length
FL	LN	11963-6	Femur Length
<i>NBL</i>	<i>99ALOKA</i>	<i>A12006-001</i>	<i>Nasal Bone Length</i>

CID 12007 Fetal Cranium

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
LVW	LN	12171-5	Lateral Ventricular width
IOD	LN	33070-4	Inner Orbital Diameter
OOD	LN	11629-3	Outer Orbital Diameter
CD	LN	11863-8	Trans Cerebellar Diameter
HW	LN	12170-7	Width of Hemisphere

CID 12008 OB-GYN Amniotic Sac

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
AFI, Q1	LN	11624-4	First Quadrant Diameter
AFI, Q2	LN	11626-9	Second Quadrant Diameter
AFI, Q3	LN	11625-1	Third Quadrant Diameter
AFI, Q4	LN	11623-6	Fourth Quadrant Diameter

CID 12009 Early Gestation Biometry Measurements

(Local Version: 20051202)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
CRL	LN	11957-8	Crown Rump Length
GS	LN	11850-5	Gestational Sac Diameter
NT	LN	33069-6	Nuchal Translucency
EES	99ALOKA	A12009-001	Early Embryonic Size
mGS, D1	99ALOKA	A12009-002	Gestational Sac Diameter 1
mGS, D2	99ALOKA	A12009-003	Gestational Sac Diameter 2
mGS, D3	99ALOKA	A12009-004	Gestational Sac Diameter 3
mGS, mGS	99ALOKA	A12009-005	Mean Gestational Sac Diameter

CID 12011 Ultrasound Pelvis and Uterus

(Local Version: 20051202)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
Cervix	LN	11961-0	Cervix Length
Cervix, Length	LN	11961-0	Cervix Length
Endom-T	LN	12145-9	Endometrium Thickness

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
<i>Cervix, Width</i>	99ALOKA	A12011-001	<i>Cervix Width</i>
<i>Cervix, AP</i>	99ALOKA	A12011-002	<i>Cervix Antero-Posterior Diameter</i>
<i>Pre Bldrvol, Length</i>	99ALOKA	A12011-003	<i>Pre Void Bladder Length</i>
<i>Pre Bldrvol, AP</i>	99ALOKA	A12011-004	<i>Pre Void Bladder Antero-Posterior Diameter</i>
<i>Pre Bldrvol, Width</i>	99ALOKA	A12011-005	<i>Pre Void Bladder Width</i>
<i>Pre Bldrvol, Volume</i>	99ALOKA	A12011-006	<i>Pre Void Bladder Volume</i>
<i>Pst Bldrvol, Length</i>	99ALOKA	A12011-007	<i>Post Void Bladder Length</i>
<i>Pst Bldrvol, AP</i>	99ALOKA	A12011-008	<i>Post Void Bladder Antero-Posterior Diameter</i>
<i>Pst Bldrvol, Width</i>	99ALOKA	A12011-009	<i>Post Void Bladder Width</i>
<i>Pst Bldrvol, Volume</i>	99ALOKA	A12011-010	<i>Post Void Bladder Volume</i>
<i>Void Volume</i>	99ALOKA	A12011-011	<i>Bladder Void Volume</i>

CID 12013 Gestational Age Equations and Tables

(Local Version: 20051202, 20200302)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
AC (Had-90%)	LN	11892-7	AC, Hadlock 1984
AC (Hadlock84)	LN	11892-7	AC, Hadlock 1984
AC (Hansmann)	LN	33073-8	AC, Hansmann1985
AC (Shinozuka)	LN	33076-1	AC, Shinozuka 1996
AXT (Shinozuka)	LN	33078-7	AxT, Shinozuka 1996
BPD (Hadlock84)	LN	11902-4	BPD, Hadlock 1984
BPD (Hansmann)	LN	11903-2	BPD, Hansmann 1985
BPD (Kurtz)	LN	11906-5	BPD, Kurtz 1980
BPD (Rempen)	LN	33083-7	BPD, Rempen 1991
BPD (Shinozuka)	LN	33084-5	BPD, Shinozuka 1996
BPD (Tokyo U)	LN	33085-2	BPD, Tokyo 1986
CRL (Daya)	LN	33091-0	CRL, Daya 1993
CRL (Hadlock)	LN	11910-7	CRL, Hadlock 1992
CRL (Hansmann)	LN	11911-5	CRL, Hansmann 1985

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
CRL (Nelson)	LN	11913-1	CRL, Nelson 1981
CRL (Rempen)	LN	33094-4	CRL, Rempen 1991
CRL (Robinson)	LN	11914-9	CRL, Robinson 1975
FL (Had-90%)	LN	11920-6	FL, Hadlock 1984
FL (Hadlock84)	LN	11920-6	FL, Hadlock 1984
FL (Hansmann)	LN	11921-4	FL, Hansmann 1985
FL (Hohler)	LN	11922-2	FL, Hohler 1982
FL (Jeanty)	LN	11923-0	FL, Jeanty 1984
FL (Shinozuka)	LN	33102-5	FL, Shinozuka 1996
FL (Tokyo U)	LN	33103-3	FL, Tokyo 1986
mGS (Hellman)	LN	11928-9	GS, Hellman 1969
mGS (Rempen)	LN	11929-7	GS, Rempen 1991
HC (Had-90%)	LN	11932-1	HC, Hadlock 1984
HC (Hadlock84)	LN	11932-1	HC, Hadlock 1984
HC2 (Hansmann)	LN	33112-4	HC, Hansmann 1985
HL (Jeanty)	LN	11936-2	Humerus, Jeanty 1984
LV (Tokyo U)	LN	33118-1	Length of Vertebra, Tokyo 1986
OFD (Hansmann)	LN	33544-8	OFD, Hansmann 1985
RAD (Jea-95%)	LN	33126-4	Radius, Jeanty 1983
TIB (Jeanty)	LN	11941-2	Tibia, Jeanty 1984
TC (Chitkara U)	LN	33131-4	ThC, Chitkara 1987
CD (Goldstein)	LN	33133-0	TCD, Goldstein 1987
CD (Hill)	LN	33134-8	TCD, Hill 1990
ULNA (Jeanty)	LN	11944-6	Ulna, Jeanty 1984
AC (Campbell)	99ALOKA	A12013-001	AC, Campbell
AC (Chitty)	99ALOKA	A12013-002	AC, Chitty 1994
AC (Hadlock)	99ALOKA	A12013-003	AC, Hadlock 1982
AC (Merz)	99ALOKA	A12013-005	AC, Merz 1996
APTD (Merz)	99ALOKA	A12013-006	APAD, Merz 1996
AXT (Tokyo U)	99ALOKA	A12013-007	AxT, Tokyo
BD (Jeanty)	99ALOKA	A12013-008	BD, Jeanty 1984
BPD (Campbell)	99ALOKA	A12013-009	BPD, Campbell
BPD (Chitty)	99ALOKA	A12013-010	BPD-oi, Chitty 1994
BPD (Hadlock)	99ALOKA	A12013-011	BPD, Hadlock 1982
BPD (Merz)	99ALOKA	A12013-012	BPD, Merz 1996

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
<i>BPD (Sabbagha)</i>	99ALOKA	A12013-013	<i>BPD, Sabbagha 1976</i>
<i>CRL (JSUM'03)</i>	99ALOKA	A12013-016	<i>CRL, JSUM 2003</i>
<i>CRL (Tokyo U)</i>	99ALOKA	A12013-017	<i>CRL, Tokyo</i>
<i>EES (Goldstein)</i>	99ALOKA	A12013-018	<i>EES, Goldstein 1994</i>
<i>FIB (Merz)</i>	99ALOKA	A12013-019	<i>FIB, Merz 1996</i>
<i>FL (Campbell)</i>	99ALOKA	A12013-020	<i>FL, Campbell</i>
<i>FL (Chitty)</i>	99ALOKA	A12013-021	<i>FL, Chitty 1994</i>
<i>FL (Hadlock)</i>	99ALOKA	A12013-022	<i>FL, Hadlock 1982</i>
<i>FL (Jea-95%)</i>	99ALOKA	A12013-025	<i>FL, Jeanty 95% 1983</i>
<i>FL (Merz)</i>	99ALOKA	A12013-026	<i>FL, Merz 1996</i>
<i>FL (Warda)</i>	99ALOKA	A12013-027	<i>FL, Warda 1985</i>
<i>GS (Tokyo U)</i>	99ALOKA	A12013-028	<i>GS, Tokyo</i>
<i>HC (Campbell)</i>	99ALOKA	A12013-029	<i>HC, Campbell</i>
<i>HC (Chitty)</i>	99ALOKA	A12013-030	<i>HC, Chitty 1994</i>
<i>HC (Hadlock)</i>	99ALOKA	A12013-031	<i>HC, Hadlock 1982</i>
<i>HC2 (Merz)</i>	99ALOKA	A12013-033	<i>HC, Merz 1996</i>
<i>HL (Hansmann)</i>	99ALOKA	A12013-034	<i>Humerus, Hansmann 1985</i>
<i>HL (Jea-95%)</i>	99ALOKA	A12013-035	<i>Humerus, Jeanty 95% 1983</i>
<i>HL (Merz)</i>	99ALOKA	A12013-036	<i>Humerus, Merz 1996</i>
<i>NBL (Sonek)</i>	99ALOKA	A12013-037	<i>NBL, Sonek 2003</i>
<i>OFD (Merz)</i>	99ALOKA	A12013-038	<i>OFD, Merz 1996</i>
<i>RAD (Merz)</i>	99ALOKA	A12013-039	<i>Radius, Merz 1996</i>
<i>TIB (Jea-95%)</i>	99ALOKA	A12013-040	<i>Tibia, Jeanty 95% 1983</i>
<i>TIB (Merz)</i>	99ALOKA	A12013-041	<i>Tibia, Merz 1996</i>
<i>TL (Chitkara U)</i>	99ALOKA	A12013-042	<i>TL, Chitkara 1987</i>
<i>TTD (Hansmann)</i>	99ALOKA	A12013-043	<i>TAD, Hansmann 1985</i>
<i>TTD (Merz)</i>	99ALOKA	A12013-044	<i>TAD, Merz 1996</i>
<i>ULNA (Jea-95%)</i>	99ALOKA	A12013-045	<i>Ulna, Jeanty 95% 1983</i>
<i>ULNA (Merz)</i>	99ALOKA	A12013-046	<i>Ulna, Merz 1996</i>
<i>AC (JSUM'03)</i>	99ALOKA	A12013-047	<i>AC, JSUM 2003</i>
<i>BPD (JSUM'03)</i>	99ALOKA	A12013-048	<i>BPD, JSUM 2003</i>
<i>BPD (Osaka U)</i>	99ALOKA	A12013-049	<i>BPD, Osaka</i>
<i>CRL (Osaka U)</i>	99ALOKA	A12013-050	<i>CRL, Osaka</i>
<i>FL (JSUM'03)</i>	99ALOKA	A12013-051	<i>FL, JSUM 2003</i>
<i>FL (O'Brien)</i>	99ALOKA	A12013-052	<i>FL, O'Brien 1981</i>
<i>FL (Osaka U)</i>	99ALOKA	A12013-053	<i>FL, Osaka</i>

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
<i>FTA (Osaka U)</i>	99ALOKA	A12013-054	<i>FTA, Osaka</i>
<i>HL (Osaka U)</i>	99ALOKA	A12013-055	<i>Humerus, Osaka</i>
<i>CRL (BMUS 2009)</i>	99ALOKA	A12013-056	<i>CRL, BMUS 2009</i>
<i>HC3 (BMUS 2009)</i>	99ALOKA	A12013-057	<i>HC, BMUS 2009</i>
<i>FL (BMUS 2009)</i>	99ALOKA	A12013-058	<i>FL, BMUS 2009</i>
<i>HC3 (BMUS-95% 2009)</i>	99ALOKA	A12013-059	<i>HC, BMUS-95% 2009</i>
<i>AC2 (BMUS-95% 2009)</i>	99ALOKA	A12013-060	<i>AC, BMUS-95% 2009</i>
<i>FL (BMUS-95% 2009)</i>	99ALOKA	A12013-061	<i>FL, BMUS-95% 2009</i>
<i>BPD (CFEF-97% 2006)</i>	99ALOKA	A12013-062	<i>BPD, CFEF-97% 2006</i>
<i>HC3 (CFEF-97% 2006)</i>	99ALOKA	A12013-063	<i>HC, CFEF-97% 2006</i>
<i>AC2 (CFEF-97% 2006)</i>	99ALOKA	A12013-064	<i>AC, CFEF-97% 2006</i>
<i>FL (CFEF-97% 2006)</i>	99ALOKA	A12013-066	<i>FL, CFEF-97% 2006</i>
<i>BPD (CFEF-90% 2006)</i>	99ALOKA	A12013-067	<i>BPD, CFEF-90% 2006</i>
<i>HC3 (CFEF-90% 2006)</i>	99ALOKA	A12013-068	<i>HC, CFEF-90% 2006</i>
<i>AC2 (CFEF-90% 2006)</i>	99ALOKA	A12013-069	<i>AC, CFEF-90% 2006</i>
<i>FL (CFEF-90% 2006)</i>	99ALOKA	A12013-071	<i>FL, CFEF-90% 2006</i>
<i>GA Table/Equ. (User1)</i>	99ALOKA	A12013-072	<i>GA, User defined Table or Equation1</i>
<i>GA Table/Equ. (User2)</i>	99ALOKA	A12013-073	<i>GA, User defined Table or Equation2</i>
<i>GA Table/Equ. (User3)</i>	99ALOKA	A12013-074	<i>GA, User defined Table or Equation3</i>
<i>GA Table/Equ. (User4)</i>	99ALOKA	A12013-075	<i>GA, User defined Table or Equation4</i>
<i>GA Table/Equ. (User5)</i>	99ALOKA	A12013-076	<i>GA, User defined Table or Equation5</i>
<i>GA Table/Equ. (User6)</i>	99ALOKA	A12013-077	<i>GA, User defined Table or Equation6</i>
<i>GA Table/Equ. (User7)</i>	99ALOKA	A12013-078	<i>GA, User defined Table or Equation7</i>

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
GA Table/Equ. (User8)	99ALOKA	A12013-079	GA, User defined Table or Equation8
GA Table/Equ. (User9)	99ALOKA	A12013-080	GA, User defined Table or Equation9
GA Table/Equ. (User10)	99ALOKA	A12013-081	GA, User defined Table or Equation10
GA Table/Equ. (User11)	99ALOKA	A12013-082	GA, User defined Table or Equation11
GA Table/Equ. (User12)	99ALOKA	A12013-083	GA, User defined Table or Equation12
GA Table/Equ. (User13)	99ALOKA	A12013-084	GA, User defined Table or Equation13
GA Table/Equ. (User14)	99ALOKA	A12013-085	GA, User defined Table or Equation14
GA Table/Equ. (User15)	99ALOKA	A12013-086	GA, User defined Table or Equation15
GA Table/Equ. (User16)	99ALOKA	A12013-087	GA, User defined Table or Equation16
GA Table/Equ. (User17)	99ALOKA	A12013-088	GA, User defined Table or Equation17
GA Table/Equ. (User18)	99ALOKA	A12013-089	GA, User defined Table or Equation18
GA Table/Equ. (User19)	99ALOKA	A12013-090	GA, User defined Table or Equation19
GA Table/Equ. (User20)	99ALOKA	A12013-091	GA, User defined Table or Equation20

CID 12014 OB Fetal Body Weight Equations and Tables

(Local Version: 20051202, 20200302)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
FW Equation (Campbell)	LN	11756-4	EFW by AC, Campbell 1975
FW Equation (Hadlock5)	LN	11732-5	EFW by AC, BPD, FL, HC, Hadlock 1985
FW Equation (Hadlock1)	LN	11751-5	EFW by AC, FL, Hadlock 1985
FW Equation (Hadlock2)	LN	11746-5	EFW by AC, FL, HC, Hadlock 1985
FW Equation (Hadlock4)	LN	11754-9	EFW by AC, HC Hadlock 1984

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
FW Equation (Hansmann)	LN	33139-7	EFW by BPD, TTD, Hansmann 1986
FW Equation (Shepard)	LN	11739-0	EFW by AC and BPD, Shepard 1982
FW Equation (Shinozuka)	LN	33142-1	EFW2 by Shinozuka 1996
FW Equation (Tokyo U)	LN	33144-7	EFW by BPD, APAD, TAD, FL, Tokyo 1987
<i>FW Equation (Hadlock3)</i>	99ALOKA	<i>A12014-001</i>	<i>EFW by BPD, AC, FL, Hadlock</i>
<i>FW Equation (JSUM'03)</i>	99ALOKA	<i>A12014-002</i>	<i>EFW by BPD, AC, FL, JSUM 2003</i>
<i>FW Equation (Osaka U)</i>	99ALOKA	<i>A12014-003</i>	<i>EFW by BPD, FTA, FL, Osaka</i>
<i>FW Equation (Warsof)</i>	99ALOKA	<i>A12014-004</i>	<i>EFW by BPD, AC, Warsof 1977</i>
<i>FW Equation (User1)</i>	99ALOKA	<i>A12014-005</i>	<i>EFW by User defined Equation1</i>
<i>FW Equation (User2)</i>	99ALOKA	<i>A12014-006</i>	<i>EFW by User defined Equation2</i>
<i>FW Equation (User3)</i>	99ALOKA	<i>A12014-007</i>	<i>EFW by User defined Equation3</i>
<i>FW Equation (User4)</i>	99ALOKA	<i>A12014-008</i>	<i>EFW by User defined Equation4</i>
<i>FW Equation (User5)</i>	99ALOKA	<i>A12014-009</i>	<i>EFW by User defined Equation5</i>

CID 12015 Fetal Growth Equations and Tables

(Local Version: 20051202, 20200302)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
CI (BPDo/OFD0) (Hadlock)	LN	33158-7	Cephalic Index by GA, Hadlock 1981
HC/AC (Campbell)	LN	33182-7	HC/AC by GA, Campbell 1977
<i>FL/AC (Hadlock)</i>	99ALOKA	<i>A12015-001</i>	<i>FL/AC by GA, Hadlock 1983</i>
<i>FL/HC (Hadlock)</i>	99ALOKA	<i>A12015-002</i>	<i>FL/HC by GA, Hadlock 1984</i>
<i>FL/BPD (Hohler)</i>	99ALOKA	<i>A12015-003</i>	<i>FL/BPD by GA, Hohler 1981</i>
<i>AFI (Jeng)</i>	99ALOKA	<i>A12015-004</i>	<i>AFI by GA, Jeng et al.</i>
<i>AFI (Moore)</i>	99ALOKA	<i>A12015-005</i>	<i>AFI by GA, Moore et al.</i>
<i>AFI (Phelan)</i>	99ALOKA	<i>A12015-006</i>	<i>AFI by GA, Phelan et al.</i>

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
AC (JSUM'03)	99ALOKA	A12015-007	AC by GA, JSUM 2003
BPD (JSUM'03)	99ALOKA	A12015-008	BPD by GA, JSUM 2003
BPD (Osaka U)	99ALOKA	A12015-009	BPD by GA, Osaka
CRL (Osaka U)	99ALOKA	A12015-010	CRL by GA, Osaka
FL (JSUM'03)	99ALOKA	A12015-011	FL by GA, JSUM 2003
FL (O'Brien)	99ALOKA	A12015-012	FL by GA, O'Brien 1981
FL (Osaka U)	99ALOKA	A12015-013	FL by GA, Osaka
FTA (Osaka U)	99ALOKA	A12015-014	FTA by GA, Osaka
HL (Osaka U)	99ALOKA	A12015-015	Humerus by GA, Osaka
BPD (CFEF-97% 2006)	99ALOKA	A12015-016	BPD by GA, CFEF-97% 2006
HC3 (CFEF-97% 2006)	99ALOKA	A12015-017	HC by GA, CFEF-97% 2006
AC2 (CFEF-97% 2006)	99ALOKA	A12015-018	AC by GA, CFEF-97% 2006
FL (CFEF-97% 2006)	99ALOKA	A12015-020	FL by GA, CFEF-97% 2006
BPD (CFEF-90% 2006)	99ALOKA	A12015-021	BPD by GA, CFEF-90% 2006
HC3 (CFEF-90% 2006)	99ALOKA	A12015-022	HC by GA, CFEF-90% 2006
AC2 (CFEF-90% 2006)	99ALOKA	A12015-023	AC by GA, CFEF-90% 2006
FL (CFEF-90% 2006)	99ALOKA	A12015-025	FL by GA, CFEF-90% 2006
BPD (Had 84_90%)	LN	33198-3	BPD by GA, Hadlock 1984
HC (Had 84_90%)	LN	33173-6	HC by GA, Hadlock 1984
AC (Had 84_90%)	LN	33146-2	AC by GA, Hadlock 1984
FL (Had 84_90%)	LN	33166-0	FL by GA, Hadlock 1984
BPD (Had 84_97%)	LN	33166-0	BPD by GA, Hadlock 1984
HC(Had 84_97%)	LN	33173-6	HC by GA, Hadlock 1984
AC (Had 84_97%)	LN	33146-2	AC by GA, Hadlock 1984
FL (Had 84_97%)	LN	33166-0	FL by GA, Hadlock 1984

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
<i>GA Table/Equ. (User1)</i>	99ALOKA	A12015-026	<i>GA by GA, User defined Table or Equation1</i>
<i>GA Table/Equ. (User2)</i>	99ALOKA	A12015-027	<i>GA by GA, User defined Table or Equation2</i>
<i>GA Table/Equ. (User3)</i>	99ALOKA	A12015-028	<i>GA by GA, User defined Table or Equation3</i>
<i>GA Table/Equ. (User4)</i>	99ALOKA	A12015-029	<i>GA by GA, User defined Table or Equation4</i>
<i>GA Table/Equ. (User5)</i>	99ALOKA	A12015-030	<i>GA by GA, User defined Table or Equation5</i>
<i>GA Table/Equ. (User6)</i>	99ALOKA	A12015-031	<i>GA by GA, User defined Table or Equation6</i>
<i>GA Table/Equ. (User7)</i>	99ALOKA	A12015-032	<i>GA by GA, User defined Table or Equation7</i>
<i>GA Table/Equ. (User8)</i>	99ALOKA	A12015-033	<i>GA by GA, User defined Table or Equation8</i>
<i>GA Table/Equ. (User9)</i>	99ALOKA	A12015-034	<i>GA by GA, User defined Table or Equation9</i>
<i>GA Table/Equ. (User10)</i>	99ALOKA	A12015-035	<i>GA by GA, User defined Table or Equation10</i>
<i>GA Table/Equ. (User11)</i>	99ALOKA	A12015-036	<i>GA by GA, User defined Table or Equation11</i>
<i>GA Table/Equ. (User12)</i>	99ALOKA	A12015-037	<i>GA by GA, User defined Table or Equation12</i>
<i>GA Table/Equ. (User13)</i>	99ALOKA	A12015-038	<i>GA by GA, User defined Table or Equation13</i>
<i>GA Table/Equ. (User14)</i>	99ALOKA	A12015-039	<i>GA by GA, User defined Table or Equation14</i>
<i>GA Table/Equ. (User15)</i>	99ALOKA	A12015-040	<i>GA by GA, User defined Table or Equation15</i>
<i>GA Table/Equ. (User16)</i>	99ALOKA	A12015-041	<i>GA by GA, User defined Table or Equation16</i>
<i>GA Table/Equ. (User17)</i>	99ALOKA	A12015-042	<i>GA by GA, User defined Table or Equation17</i>
<i>GA Table/Equ. (User18)</i>	99ALOKA	A12015-043	<i>GA by GA, User defined Table or Equation18</i>
<i>GA Table/Equ. (User19)</i>	99ALOKA	A12015-044	<i>GA by GA, User defined Table or Equation19</i>
<i>GA Table/Equ. (User20)</i>	99ALOKA	A12015-045	<i>GA by GA, User defined Table or Equation20</i>
<i>Fetus Ratio Table/Equation (User1)</i>	99ALOKA	A12015-046	<i>Fetus Ratio, User defined Table or Equation1</i>

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
<i>Fetus Ratio Table/Equation (User2)</i>	99ALOKA	A12015-047	<i>Fetus Ratio, User defined Table or Equation2</i>
<i>Fetus Ratio Table/Equation (User3)</i>	99ALOKA	A12015-048	<i>Fetus Ratio, User defined Table or Equation3</i>
<i>Fetus Ratio Table/Equation (User4)</i>	99ALOKA	A12015-049	<i>Fetus Ratio, User defined Table or Equation4</i>
<i>Fetus Ratio Table/Equation (User5)</i>	99ALOKA	A12015-050	<i>Fetus Ratio, User defined Table or Equation5</i>
<i>Fetus Ratio Table/Equation (User6)</i>	99ALOKA	A12015-051	<i>Fetus Ratio, User defined Table or Equation6</i>
<i>AFI Table/Equation (User1)</i>	99ALOKA	A12015-052	<i>AFI by User defined Table or Equation1</i>
<i>AFI Table/Equation (User2)</i>	99ALOKA	A12015-053	<i>AFI by User defined Table or Equation2</i>

CID 12016 Estimated Fetal Weight Percentile Equations and Tables (Local Version: 20051202, 2160617, 20200302)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
FW Growth (Hadlock)	LN	33183-5	FWP by GA, Hadlock 1991
FW Growth (Brenner)	LN	33189-2	FWP by GA, Brenner 1976
<i>FW Growth (Doubilet)</i>	99ALOKA	A12016-001	<i>FW, Doubilet 1997</i>
<i>FW Growth (Yarkoni(Twins))</i>	99ALOKA	A12016-005	<i>Twins FW, Yarkoni 1987</i>
CFEF_97%	99ALOKA	A12016-006	FW, CFEF-97% 2014
CFEF_90%	99ALOKA	A12016-007	FW, CFEF-90% 2014
<i>FW Growth Table/Equation (User1)</i>	99ALOKA	A12016-008	<i>FW, User defined Table or Equation1</i>
<i>FW Growth Table/Equation (User2)</i>	99ALOKA	A12016-009	<i>FW, User defined Table or Equation2</i>

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
<i>FW Growth Table/Equation (User3)</i>	99ALOKA	A12016-010	<i>FW, User defined Table or Equation3</i>
<i>FW Growth Table/Equation (User4)</i>	99ALOKA	A12016-011	<i>FW, User defined Table or Equation4</i>
<i>FW Growth Table/Equation (User5)</i>	99ALOKA	A12016-012	<i>FW, User defined Table or Equation5</i>

CID 12017 Growth Distribution Rank

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	DCM	125013	Growth Z-score

CID 12018 OB-GYN Summary

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	LN	11878-6	Number of Fetuses

CID 12019 OB-GYN Fetus Summary

(Local Version: 20151202, 20150515, 20151201)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
Composite US-GA	LN	11888-5	Composite Ultrasound Age
LMP-GA	LN	11885-1	Gestational Age by LMP
FW	LN	11727-5	Estimated Weight
- (FW %ile rank by Doubilet)	LN	11767-1	EFW percentile rank
FHR	LN	11948-7	Fetal Heart Rate
<i>PreHR(Amnio)</i>	99ALOKA	A12019-001	<i>Fetal Heart Rate before Biopsy</i>
<i>PstHR(Amnio)</i>	99ALOKA	A12019-002	<i>Fetal Heart Rate after Biopsy</i>
<i>Auto FHR</i>	99ALOKA	A12019-003	<i>Auto Fetal Heart Rate</i>
EFW Ratio (up to 5 items)	99ALOKA	A12019-004	<i>Estimated Fetal Weight Ratio</i>

CID 12101 Vascular Summary

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
Comments	DCM	121106	Comment

CID 12104 Extracranial Arteries

(Local Version: 20200424)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
BIFUR	SRT	T-45160	Carotid Bifurcation
CCA prox CCA mid CCA distal	SRT	T-45100	Common Carotid Artery
ECA	SRT	T-45200	External Carotid Artery
ICA ICA prox ICA mid ICA distal	SRT	T-45300	Internal Carotid Artery
VA VERT	SRT	T-45700	Vertebral Artery
<i>max-IMT</i>	99ALOKA	A12104-001	<i>Maximum Intima-media thickness</i>
<i>mean-IMT</i>	99ALOKA	A12104-002	<i>Mean Intima-media thickness</i>
<i>CCA-IMT (Auto)</i>	99ALOKA	A12104-003	<i>Intima-media thickness of CCA by auto-detection</i>
<i>ICA-IMT (Auto)</i>	99ALOKA	A12104-004	<i>Intima-media thickness of ICA by auto-detection</i>
<i>ECA-IMT (Auto)</i>	99ALOKA	A12104-005	<i>Intima-media thickness of ECA by auto-detection</i>
<i>BIF-IMT (Auto)</i>	99ALOKA	A12104-006	<i>Intima-media thickness of BIFUR by auto-detection</i>
<i>CmeanIMt (Auto)</i>	99ALOKA	A12104-007	<i>Mean Carotid Intima-media thickness by auto-detection</i>
<i>IMT-C10</i>	99ALOKA	A12104-008	<i>IMT-C10</i>
<i>Auto IMT-C10</i>	99ALOKA	A12104-009	<i>IMT-C10 by auto-detection</i>

CID 12105 Intracranial Cerebral Vessels

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
ACA	SRT	T-45540	Anterior Cerebral Artery

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
TICA	SRT	R-102BD	Terminal internal carotid artery
MCA	SRT	T-45600	Middle Cerebral Artery
PCA	SRT	T-45900	Posterior Cerebral Artery
PCoA	SRT	T-45320	Posterior Communicating Artery

CID 12106 Intracranial Cerebral Vessels (unilateral)

(Local Version: 20060807)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
BA	SRT	T-45800	Basilar Artery
ACoA	SRT	T-45530	Anterior Communicating Artery

CID 12107 Upper Extremity Arteries

(Local Version: 20060807)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
AA	SRT	T-47100	Axillary Artery
BA	SRT	T-47160	Brachial Artery
RA	SRT	T-47300	Radial Artery
ScA	SRT	T-46100	Subclavian Artery
SPA	SRT	T-47240	Superficial Palmar Arch
UA	SRT	T-47200	Ulnar Artery
DBA	99ALOKA	A12107-001	Deep Brachial Artery
BasA	99ALOKA	A12107-002	Basilic Artery
Upr Art.1	99ALOKA	A12107-003	User Definition Artery1
Upr Art.2	99ALOKA	A12107-004	User Definition Artery2
Upr Art.3	99ALOKA	A12107-005	User Definition Artery3
Upr Art.4	99ALOKA	A12107-006	User Definition Artery4
Upr Art.5	99ALOKA	A12107-007	User Definition Artery5
Upr Art.6	99ALOKA	A12107-008	User Definition Artery6
Upr Art.7	99ALOKA	A12107-009	User Definition Artery7
Upr Art.8	99ALOKA	A12107-010	User Definition Artery8

CID 12108 Upper Extremity Veins

(Local Version: 20060807)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
AV	SRT	T-49110	Axillary vein
BasV	SRT	T-48052	Basilic vein
BV	SRT	T-49350	Brachial vein
CV	SRT	T-49240	Cephalic vein
IJV	SRT	T-48170	Internal Jugular vein
RV	SRT	T-49340	Radial vein
ScV	SRT	T-48330	Subclavian vein
UV	SRT	T-49330	Ulnar vein
DBV	99ALOKA	A12108-001	Deep Brachial vein
Upr Vein.1	99ALOKA	A12108-002	User Definition Vein1
Upr Vein.2	99ALOKA	A12108-003	User Definition Vein2
Upr Vein.3	99ALOKA	A12108-004	User Definition Vein3
Upr Vein.4	99ALOKA	A12108-005	User Definition Vein4
Upr Vein.5	99ALOKA	A12108-006	User Definition Vein5
Upr Vein.6	99ALOKA	A12108-007	User Definition Vein6
Upr Vein..7	99ALOKA	A12108-008	User Definition Vein7
Upr Vein.8	99ALOKA	A12108-009	User Definition Vein8

CID 12109 Lower Extremity Arteries

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
CIA	SRT	T-46710	Common Iliac Artery
ATA	SRT	T-47700	Anterior Tibial Artery
CFA	SRT	T-47400	Common Femoral Artery
DPA	SRT	T-47741	Dorsalis Pedis Artery
EIA	SRT	T-46910	External Iliac Artery
IIA	SRT	T-46740	Internal Iliac Artery
PerA	SRT	T-47630	Peroneal Artery
PopA	SRT	T-47500	Popliteal Artery
PTA	SRT	T-47600	Posterior Tibial Artery
DFA (Deep Femoral Artery)	SRT	T-47440	Profunda Femoris Artery
SFA	SRT	T-47403	Superficial Femoral Artery
Lwr Art.1	99ALOKA	A12109-001	User Definition Artery1
Lwr Art.2	99ALOKA	A12109-002	User Definition Artery2

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
<i>Lwr Art.3</i>	99ALOKA	A12109-003	<i>User Definition Artery3</i>
<i>Lwr Art.4</i>	99ALOKA	A12109-004	<i>User Definition Artery4</i>
<i>Lwr Art.5</i>	99ALOKA	A12109-005	<i>User Definition Artery5</i>
<i>Lwr Art.6</i>	99ALOKA	A12109-006	<i>User Definition Artery6</i>
<i>Lwr Art.7</i>	99ALOKA	A12109-007	<i>User Definition Artery7</i>
<i>Lwr Art.8</i>	99ALOKA	A12109-008	<i>User Definition Artery8</i>

CID 12110 Lower Extremity Veins

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
ATV	SRT	T-49630	Anterior Tibial Vein
CFV	SRT	G-035B	Common Femoral Vein
CIV	SRT	T-48920	Common Iliac Vein
EIV	SRT	T-48930	External Iliac Vein
GSV	SRT	T-49530	Great Saphenous Vein
LSV	SRT	T-49550	Lesser Saphenous Vein
PerV	SRT	T-49650	Peroneal Vein
PopV	SRT	T-49640	Popliteal Vein
PTV	SRT	T-49620	Posterior Tibial Vein
DFV (Deep Femoral Vein)	SRT	T-49660	Profunda Femoris Vein
SFV	SRT	G-035A	Superficial Femoral Vein
IIV	SRT	T-48940	Internal iliac vein
<i>Lwr Vein.1</i>	99ALOKA	A12110-001	<i>User Definition Vein1</i>
<i>Lwr Vein.2</i>	99ALOKA	A12110-002	<i>User Definition Vein2</i>
<i>Lwr Vein.3</i>	99ALOKA	A12110-003	<i>User Definition Vein3</i>
<i>Lwr Vein.4</i>	99ALOKA	A12110-004	<i>User Definition Vein4</i>
<i>Lwr Vein.5</i>	99ALOKA	A12110-005	<i>User Definition Vein5</i>
<i>Lwr Vein.6</i>	99ALOKA	A12110-006	<i>User Definition Vein6</i>
<i>Lwr Vein.7</i>	99ALOKA	A12110-007	<i>User Definition Vein7</i>
<i>Lwr Vein.8</i>	99ALOKA	A12110-008	<i>User Definition Vein8</i>

CID 12111 Abdominal Arteries (lateral)

(Local Version: 20070424)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
CIA	SRT	T-46710	Common Iliac Artery

CID 12112 Abdominal Arteries (unilateral)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
A-Ao	SRT	T-42000	Aorta
CA	SRT	T-46400	Celiac Axis
CHA	SRT	T-46421	Common Hepatic Artery
IMA	SRT	T-46520	Inferior Mesenteric Artery
HA, Right	SRT	T-46423	Right Branch of Hepatic Artery
HA, Left	SRT	T-46427	Left Branch of Hepatic Artery
SA	SRT	T-46460	Splenic Artery
SMA Prandial SMA	SRT	T-46510	Superior Mesenteric Artery
Artery1	99ALOKA	A12112-001	User Definition Artery1
Artery2	99ALOKA	A12112-002	User Definition Artery2
Artery3	99ALOKA	A12112-003	User Definition Artery3

CID 12114 Abdominal Veins (unilateral)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
Main PV	SRT	T-48810	Portal Vein
Lt.PV	SRT	T-4881F	Left Main Branch of Portal Vein
Rt.PV	SRT	T-4882A	Right Main Branch of Portal Vein
Prox Shunt Mid Shunt Distal Shunt	SRT	G-036C	Transjugular Intrahepatic Portosystemic Shunt

CID 12115 Renal Vessels

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
Renal-A	SRT	T-46600	Renal Artery

CID 12116 Vessel Segment Modifiers

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
- (distal)	SRT	G-A119	Distal
- (mid)	SRT	G-A188	Mid-longitudinal
- (prox)	SRT	G-A118	Proximal

CID 12117 Vessel Branch Modifiers

(Local Version: 20200424)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
- (<i>Farwall</i>)	99ALOKA	A12117-001	<i>far wall</i>
- (<i>Nearwall</i>)	99ALOKA	A12117-002	<i>near wall</i>

CID 12120 Blood Velocity Measurements

(Local Version: 20150515, 20151201)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
EDV	LN	11653-3	End Diastolic Velocity
PSV pV	LN	11726-7	Peak Systolic Velocity
MnV	LN	11692-1	Time averaged peak velocity
DV(<i>Ductus Venosus</i>), a	99ALOKA	A12120-001	<i>atrial contraction</i>

CID 12121 Vascular Indices and Ratios

(Local Version: 20070424, 20150515, 20151201)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
ACC	LN	20167-3	Acceleration Index
PI	LN	12008-9	Pulsatility Index
RI	LN	12023-8	Resistivity Index
S/D	LN	12144-2	Systolic to Diastolic Velocity Ratio
AccT/FlowT	99ALOKA	A12121-001	<i>Acceleration Time to Flow Time Ratio</i>
DV(<i>Ductus Venosus</i>), a/S	99ALOKA	A12121-002	<i>atrial contraction to Systolic peak Velocity Ratio</i>
DV(<i>Ductus Venosus</i>), S/a	99ALOKA	A12121-003	<i>Systolic peak to atrial contraction Velocity Ratio</i>

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
AccT	LN	20168-1	Acceleration Time
FlowT	99ALOKA	A12122-001	Flow Time
1 max-IMT 1 IMT	99ALOKA	A12122-002	Intima-media thickness 1
2 max-IMT 2 IMT	99ALOKA	A12122-003	Intima-media thickness 2
3 max-IMT 3 IMT	99ALOKA	A12122-004	Intima-media thickness 3
4 max-IMT	99ALOKA	A12122-005	Intima-media thickness 4
5 max-IMT	99ALOKA	A12122-006	Intima-media thickness 5
6 max-IMT	99ALOKA	A12122-007	Intima-media thickness 6
7 max-IMT	99ALOKA	A12122-008	Intima-media thickness 7
8 max-IMT	99ALOKA	A12122-009	Intima-media thickness 8
9 max-IMT	99ALOKA	A12122-010	Intima-media thickness 9
10 max-IMT	99ALOKA	A12122-011	Intima-media thickness 10
max-IMT	99ALOKA	A12122-012	Maximum Intima-media thickness
mean-IMT	99ALOKA	A12122-013	Mean Intima-media thickness
D-max	99ALOKA	A12122-014	Maximum Intima-media thickness at diastole
D-min	99ALOKA	A12122-015	Minimum Intima-media thickness at diastole
D-average D-avg	99ALOKA	A12122-016	Average Intima-media thickness at diastole
D-C10	99ALOKA	A12122-017	IMT-C10 at diastole
D-SD	99ALOKA	A12122-018	Standard Deviation of Intima-media thickness at diastole
D-points	99ALOKA	A12122-019	Measured number of Intima-media thickness points at diastole
D-width	99ALOKA	A12122-020	Width of ROI at diastole
S-max	99ALOKA	A12122-021	Maximum Intima-media thickness at systole
S-min	99ALOKA	A12122-022	Minimum Intima-media thickness at systole
S-average S-avg	99ALOKA	A12122-023	Average Intima-media thickness at systole
S-C10	99ALOKA	A12122-024	IMT-C10 at systole
S-SD	99ALOKA	A12122-025	Standard Deviation of Intima-media thickness at systole
S-points	99ALOKA	A12122-026	Measured number of Intima-media thickness points at systole

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
<i>S-width</i>	99ALOKA	A12122-027	<i>Width of ROI at systole</i>
<i>D-1 IMT</i>	99ALOKA	A12122-028	<i>Intima-media thickness 1 at diastole</i>
<i>D-2 IMT</i>	99ALOKA	A12122-029	<i>Intima-media thickness 2 at diastole</i>
<i>D-3 IMT</i>	99ALOKA	A12122-030	<i>Intima-media thickness 3 at diastole</i>
<i>D-mean-IMT</i>	99ALOKA	A12122-031	<i>Mean Intima-media thickness at diastole</i>
<i>S-1 IMT</i>	99ALOKA	A12122-032	<i>Intima-media thickness 1 at systole</i>
<i>S-2 IMT</i>	99ALOKA	A12122-033	<i>Intima-media thickness 2 at systole</i>
<i>S-3 IMT</i>	99ALOKA	A12122-034	<i>Intima-media thickness 3 at systole</i>
<i>S-mean-IMT</i>	99ALOKA	A12122-035	<i>Mean Intima-media thickness at systole</i>
<i>C10</i>	99ALOKA	A12122-036	<i>IMT-C10</i>

CID 12124 Renal Ratios

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
RAR	LN	33869-9	Renal Artery/Aorta velocity ratio

CID 12140 Pelvic Vasculature Anatomical Location

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
UmA	SRT	T-F1810	Umbilical Artery
UtA	SRT	T-46820	Uterine Artery

CID 12141 Fetal Vasculature Anatomical Location

(Local Version: 20150515)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
MCA	SRT	T-45600	Middle Cerebral Artery
D-Ao	SRT	T-D0765	Descending Aorta
<i>DV(Ductus Venosus)</i>	99ALOKA	A12141-001	<i>Ductus Venosus</i>

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
LVIDd	LN	29436-3	Left Ventricle Internal End Diastolic Dimension
LVIDs	LN	29438-9	Left Ventricle Internal Systolic Dimension
FS	LN	18051-3	Left Ventricular Fractional Shortening
IVSd	LN	18154-5	Interventricular Septum Diastolic Thickness
IVS/LVPW	LN	18155-2	Interventricular Septum to Posterior Wall Thickness Ratio
%IVSTF	LN	18054-7	Interventricular Septum % Thickening
IVSs	LN	18158-6	Interventricular Septum Systolic Thickness
%PWTF	LN	18053-9	Left Ventricle Posterior Wall % Thickening
LVLd LVL4d LVL2d	LN	18077-8	Left Ventricle diastolic major axis
LVLs LVL4s LVL2s	LN	18076-0	Left Ventricle systolic major axis
LVPWs	LN	18156-0	Left Ventricle Posterior Wall Systolic Thickness
LVPWd	LN	18152-9	Left Ventricle Posterior Wall Diastolic Thickness
LVSLMVd LVSLMV _s	99ALOKA	A12201-001	Left Ventricular Short Axis Length at Mitral Valve
thick	99ALOKA	A12201-002	Mean Wall Thickness
EDL EDL _{epi}	99ALOKA	A12201-003	Left Ventricular End Diastolic Length
b	99ALOKA	A12201-004	Short Axis Radius Calculated from Short Axis Cavity Area

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
EDV	LN	18026-5	Left Ventricular End Diastolic Volume
ESV	LN	18148-7	Left Ventricular End Systolic Volume
EF	LN	18043-0	Left Ventricular Ejection Fraction

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
LVM	LN	18087-7	Left Ventricle Mass
IRT	LN	18071-1	Left Ventricular Isovolumic Relaxation Time
LVM/BSA	99ALOKA	A12203-001	Left Ventricular Mass Index
%difD %difS	99ALOKA	A12203-003	Long Axis(at End Diastole or End Systole) Length % Difference
mFS	99ALOKA	A12203-004	Midwall Fractional Shortening
areaEF _{lx}	99ALOKA	A12203-005	Area Ejection Fraction at Long Axis View
areaEF _{sx}	99ALOKA	A12203-006	Area Ejection Fraction at Short Axis View
MVCF	99ALOKA	A12203-009	Mean Velocity of Circumferential Fiber Shortening
LVDF _T	99ALOKA	A12203-010	Left Ventricle Diastole Filling Time
E/Em	99ALOKA	A12203-011	Ratio of MV E-Wave Peak Vel. to Early Diastolic Myocardium Vel.
LVDF _T /RR	99ALOKA	A12203-012	Ratio of Left Ventricle Diastole Filling Time to R-R interval
GLS GLS(ap4C) GLS(ap2C) GLS3 GLS(3P)	99ALOKA	A12203-013	Global Longitudinal Strain
LVOT/AV(pV)	99ALOKA	A12203-014	pV Ratio of LVOT to AV
LVOT/AV(VTI)	99ALOKA	A12203-015	VTI Ratio of LVOT to AV
LVMi1.7(Devereux) LVMi1.7(ASE)	99ALOKA	A12203-016	Left Ventricle Mass adjusted by Height(1.7)
LVMi2.16(Devereux) LVMi2.16(ASE)	99ALOKA	A12203-017	Left Ventricle Mass adjusted by Height(2.16)
LVMi2.7(Devereux) LVMi2.7(ASE)	99ALOKA	A12203-018	Left Ventricle Mass adjusted by Height(2.7)
ML3d ML3s ML4d ML4s ML2d ML2s	99ALOKA	A12203-019	Myocardial Length

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
RVDd	LN	20304-2	Right Ventricular Internal Diastolic Dimension
RVDs	LN	20305-9	Right Ventricular Internal Systolic Dimension
RVSP	SRT	G-0380	Right Ventricular Peak Systolic Pressure
RVAWd	LN	18153-7	Right Ventricular Anterior Wall Diastolic Thickness
RVAWs	LN	18157-8	Right Ventricular Anterior Wall Systolic Thickness
<i>RV Area ED</i>	99ALOKA	A12204-001	<i>Right Ventricular Diastolic Area</i>
<i>RV Area ES</i>	99ALOKA	A12204-002	<i>Right Ventricular Systolic Area</i>
<i>RV FAC</i>	99ALOKA	A12204-003	<i>Right Ventricular Fractional Area Change</i>
<i>RVDd base</i>	99ALOKA	A12204-004	<i>End-diastolic Basal Diameter</i>
<i>RVDd mid</i>	99ALOKA	A12204-005	<i>End-diastolic Mid-cavity Diameter</i>
<i>RVLd</i>	99ALOKA	A12204-006	<i>Right Ventricular End-diastolic Length</i>
<i>RV WT</i>	99ALOKA	A12204-007	<i>Right Ventricular Wall Thickness</i>

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
LADs	LN	29469-4	Left Atrium Antero-posterior Systolic Dimension
LADs/AODd	LN	17985-3	Left Atrium to Aortic Root Ratio
LALA4s	LN	17977-0	Left Atrium Systolic Area
<i>LADd</i>	99ALOKA	A12205-001	<i>Left Atrium Antero-posterior Diastolic Dimension</i>
<i>LAL4s</i> <i>LAL2s</i>	99ALOKA	A12205-002	<i>Left Atrium systolic major axis</i>
<i>LALA2s</i>	99ALOKA	A12205-003	<i>Left Atrium Systolic Area by Apical two chamber</i>
<i>LA Volume</i>	99ALOKA	A12205-004	<i>Left Atrial Volume</i>
<i>LA Volume/BSA</i>	99ALOKA	A12205-005	<i>Left Atrial Volume divided by Body Surface Area</i>
<i>%difS</i>	99ALOKA	A12205-006	<i>Long Axis at End Systole Length % Difference of Left Atrium</i>
<i>LAS4s</i> <i>LAS2s</i>	99ALOKA	A12205-007	<i>Left Atrium systolic minor axis</i>

CID 12206 Echocardiography Right Atrium

(Local Version: 20060807)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
RAP	LN	18070-3	Right Atrium Systolic Pressure
RALA4s	LN	17988-7	Right Atrium Systolic Area
RAL4s RAL2s	99ALOKA	A12206-001	Right Atrium systolic major axis
RALA2s	99ALOKA	A12206-002	Right Atrium Systolic Area by Apical two chamber
RA Volume	99ALOKA	A12206-003	Right Atrial Volume
RA Volume/BSA	99ALOKA	A12206-004	Right Atrial Volume divided by Body Surface Area
%difS	99ALOKA	A12206-005	Long Axis at End Systole Length % Difference of Right Atrium
RAS4s RAS2s	99ALOKA	A12206-006	Right Atrium systolic minor axis

CID 12207 Echocardiography Mitral Valve

(Local Version: 20060807, 20200424)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
aV	LN	17978-8	Mitral Valve A-Wave Peak Velocity
eV	LN	18037-2	Mitral Valve E-Wave Peak Velocity
E/A (Doppler)	LN	18038-0	Mitral Valve E to A Ratio
DecT	SRT	G-0384	Mitral Valve E-Wave Deceleration Time
E-Fslop	LN	18040-6	Mitral Valve E-F Slope by M-Mode
EPSS	LN	18036-4	Mitral Valve EPSS, E wave
Adur	SRT	G-0385	Mitral Valve A-Wave Duration
dP/dt	LN	18035-6	Mitral Regurgitation dP/dt derived from Mitral Reg. velocity
C-Eamp	99ALOKA	A12207-001	Mitral Valve Dimension of C point to E point by M-Mode
C-Aamp	99ALOKA	A12207-002	Mitral Valve Dimension of C point to A point by M-Mode
A/E	99ALOKA	A12207-003	Mitral Valve C-A Dimension to C-E Dimension Ratio by M-Mode
E/A	99ALOKA	A12207-004	Mitral Valve C-E Dimension to C-A Dimension Ratio by M-Mode
Edur	99ALOKA	A12207-005	Mitral Valve E-Wave Duration
A/E (Doppler)	99ALOKA	A12207-006	Mitral Valve A to E Ratio
EPG	99ALOKA	A12207-007	Mitral Valve E-wave Peak Pressure Gradient

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
APG	99ALOKA	A12207-008	Mitral Valve A-wave Peak Pressure Gradient
PVAdur-Adur	99ALOKA	A12207-009	Subtraction of A-wave Duration from PVA-wave Duration
VTI(MR)	99ALOKA	A12207-010	Velocity Time Integral of Mitral Regurgitant Flow
MR PISA	99ALOKA	A12207-012	Mitral Regurgitant Proximal Isovelocity Surface Area
SV(MV)	99ALOKA	A12207-013	Flow Volume of Mitral Valve Annulus in Flow
Reg. Volume	99ALOKA	A12207-014	Mitral valve Regurgitant volume
Em	LN	78185-6	MV med ann Peak e' wave US.TDI
IVA(l) IVA(s)	99ALOKA	A12207-015	Isovolumic Acceleration
av e'	99ALOKA	A12207-016	Average of e'(l) and e'(s)
D-Eamp	99ALOKA	A12207-017	D-E amplitude
MAPSE	99ALOKA	A12207-018	Mitral Annular Plane Systolic Excursion
Ann Area	99ALOKA	A12207-019	Mitral Annulus Area
Ann PLAX	99ALOKA	A12207-020	Mitral Annulus Diameter PLAX
e'/a'(l) e'/a'(s)	99ALOKA	A12207-021	Mitral Valve Ea to Aa Ratio
IVC Vel(m) IVC Vel(l) IVC Vel(s)	99ALOKA	A12207-022	Peak MV Annular Velocity in Isovolumic Contraction
E/av e'	99ALOKA	A12207-023	Ratio of MV Peak Velocity to LV av e'
TE-e' TE-e'(m) TE-e'(l) TE-e'(s)	99ALOKA	A12207-024	Time E to e'
IVC TPV(m) IVC TPV(l) IVC TPV(s)	99ALOKA	A12207-025	Time to peak velocity in isovolumic contractuion
SV (MV)	99ALOKA	A12207-026	Flow Volume of Mitral Valve Annulus in Flow (1)

CID 12208 Echocardiography Tricuspid Valve (Local Version: 20060807, 20140307, 20200424)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
dP/dt	LN	18034-9	Tricuspid Regurgitation dP/dt

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
<i>C-Eamp</i>	99ALOKA	A12208-001	<i>Tricuspid Valve Dimension of C point to E point by M-Mode</i>
<i>C-Aamp</i>	99ALOKA	A12208-002	<i>Tricuspid Valve Dimension of C point to A point by M-Mode</i>
<i>D-Eamp</i>	99ALOKA	A12208-003	<i>Tricuspid Valve Dimension of D point to E point by M-Mode</i>
<i>E-Fslop</i>	99ALOKA	A12208-004	<i>Tricuspid Valve Velocity from E point to F point by M-Mode</i>
<i>D-Eslop</i>	99ALOKA	A12208-005	<i>Tricuspid Valve Velocity from D point to E point by M-Mode</i>
<i>A/E</i>	99ALOKA	A12208-006	<i>Tricuspid Valve C-A Dimension to C-E Dimension Ratio by M-Mode</i>
<i>E/A</i>	99ALOKA	A12208-007	<i>Tricuspid Valve C-E Dimension to C-A Dimension Ratio by M-Mode</i>
<i>VTI(TR)</i>	99ALOKA	A12208-008	<i>Velocity Time Integral of Tricuspid Regurgitant Flow</i>
<i>TR PISA</i>	99ALOKA	A12208-010	<i>Tricuspid Regurgitant Proximal Isovelocity Surface Area</i>
<i>SV(TV)</i>	99ALOKA	A12208-011	<i>Flow Volume of Tricuspid Valve Annulus in Flow</i>
<i>TAPSE</i>	99ALOKA	A12208-012	<i>Tricuspid Annular Plane Systolic Excursion</i>
<i>A/E</i>	99ALOKA	A12208-013	<i>Peak A wave/Peak E wave by US</i>
<i>Edur</i>	99ALOKA	A12208-014	<i>E-Wave Duration</i>
<i>EPG</i>	99ALOKA	A12208-015	<i>E-Wave Peak Gradient</i>
<i>TV E/e'(s)</i>	99ALOKA	A12208-016	<i>Ratio of TV Peak Velocity to RV Peak Tissue Velocity E-Wave</i>
<i>Ann Vmax e'</i> <i>TV e'(s)</i> <i>TV e'(l)</i>	99ALOKA	A12208-017	<i>Right Ventricular Peak Early Diastolic Tissue Velocity</i>
<i>TV s'(s)</i> <i>TV s'(l)</i>	99ALOKA	A12208-018	<i>Right Ventricular Peak Systolic Tissue Velocity</i>
<i>Ann Vmax a'</i> <i>TV a'(s)</i> <i>TV a'(l)</i>	99ALOKA	A12208-020	<i>RV Peak Diastolic Tissue Velocity During Atrial Systole</i>
<i>dt</i>	99ALOKA	A12208-021	<i>Time Difference between v1 and v2</i>
<i>Ann Area</i>	99ALOKA	A12208-022	<i>Tricuspid Annulus Area</i>
<i>TV C-Odur</i>	99ALOKA	A12208-023	<i>Tricuspid Valve Closure to Opening Time</i>

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
ET	LN	18042-2	Pulmonic Valve Ejection Time
AccT/ET	SRT	G-0388	Ratio of Pulmonic Valve Acceleration Time to Ejection Time
A wave amp	99ALOKA	A12209-001	Pulmonic Valve Dimension of F point to A point by M-Mode
B-Camp	99ALOKA	A12209-002	Pulmonic Valve Dimension of B point to C point by M-Mode
E-Fslop	99ALOKA	A12209-003	Pulmonic Valve Velocity from E point to F point by M-Mode
B-Cslop	99ALOKA	A12209-004	Pulmonic Valve Velocity from B point to C point by M-Mode
VTI(PR)	99ALOKA	A12209-005	Velocity Time Integral of Pulmonic Regurgitant Flow
PR PISA	99ALOKA	A12209-007	Pulmonic Regurgitant Proximal Isovelocity Surface Area
Ann diam trans	99ALOKA	A12209-009	Transverse Annulus Diameter
dur diastole	99ALOKA	A12209-010	Duration of Diastole
PR index	99ALOKA	A12209-011	Ratio of PR dur to PR Index
PR VCW/Ann diam	99ALOKA	A12209-012	Ratio of PR VCW to Ann diam a-p
PADP sPAP mPAP	99ALOKA	A12209-013	Pulmonary Artery Pressure
PR dur	99ALOKA	A12209-014	Pulmonary Valve Regurgitation Duration
PR VCW	99ALOKA	A12209-015	Pulmonary Valve Regurgitation Vena Contracta Diameter
Ann Area	99ALOKA	A12209-016	Pulmonary Annulus Area
VR	99ALOKA	A12209-017	Velocity Ratio of RVOT velocity to Pulmonic Valve velocity

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
MPA MPG LPA MPG RPA MPG	LN	20256-4	Mean Gradient

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
MPA PG LPA PG RPA PG	LN	20247-3	Peak Gradient
MPA MnV LPA MnV RPA MnV	LN	20352-1	Mean Velocity
MPA pV LPA pV RPA pV	LN	11726-7	Peak Velocity
MPA VTI LPA VTI RPA VTI	LN	20354-7	Velocity Time Integral

CID 12211 Echocardiography Aortic Valve

(Local Version: 20060807, 20200424)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
ET	LN	18041-4	Aortic Valve Ejection Time
AccT/ET	SRT	G-0382	Ratio of Aortic Valve Acceleration Time to Ejection Time
VTI(AR)	99ALOKA	A12211-001	Velocity Time Integral of Aortic Regurgitant Flow
AR PISA	99ALOKA	A12211-003	Aortic Regurgitant Proximal Isovelocity Surface Area
AV Ann diam Ao Ring	99ALOKA	A12211-004	Aortic Annulus Diameter
Ann dmax	99ALOKA	A12211-005	Aortic Annulus Maximum Diameter
Ann dmin	99ALOKA	A12211-006	Aortic Annulus Minimum Diameter
Jet A/LVOT(CSA)	99ALOKA	A12211-007	Aortic Regurgitant Jet Area/LVOT Area
Jet W/LVOT	99ALOKA	A12211-008	Aortic Regurgitant Jet Width/LVOT Width
AR VCW	99ALOKA	A12211-009	Aortic Regurgitation Vena Contracta Width
Ann Area	99ALOKA	A12211-010	Aortic Valve Annulus Area
VR	99ALOKA	A12211-011	Velocity Ratio of LVOT to AV

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
AODd AODs	LN	18015-8	Aortic Root Diameter
MPG	LN	20256-4	Mean Gradient
PG	LN	20247-3	Peak Gradient
MnV	LN	20352-1	Mean Velocity
pV	LN	11726-7	Peak Velocity
VTI	LN	20354-7	Velocity Time Integral
Sinus of V	99ALOKA	A12212-001	Sinus of Valsalva
Sino tube	99ALOKA	A12212-002	Sinotubular junction
Aortic Arch	99ALOKA	A12212-003	Aortic Arch
Asc Aorta	99ALOKA	A12212-004	Ascending Aorta
Desc Aorta	99ALOKA	A12212-005	Descending Aorta
Abdom Aorta	99ALOKA	A12212-006	Abdominal Aorta
Hepatic Vein	99ALOKA	A12212-007	Hepatic Vein
Ao STJ	99ALOKA	A12212-008	Aortic Sinotubular Junction Diameter
Ao SV	99ALOKA	A12212-009	Aortic Sinuses of Valsalva Diameter

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
PVS	LN	29450-4	Pulmonary Vein Systolic Peak Velocity
PVD	LN	29451-2	Pulmonary Vein Diastolic Peak Velocity
S/D	LN	29452-0	Pulmonary Vein Systolic to Diastolic Ratio
PVA	LN	29453-8	Pulmonary Vein Atrial Contraction Reversal Peak Velocity
PVAdur	SRT	G-038B	Pulmonary Vein A-Wave Duration
D-VTI	SRT	G-038D	Pulmonary Vein D-Wave Velocity Time Integral
S-VTI	SRT	G-038C	Pulmonary Vein S-Wave Velocity Time Integral
DecT	99ALOKA	A12214-001	Deceleration Time of D-Wave Flow
SF	99ALOKA	A12214-002	Systolic Fraction

CID 12215 Echocardiography Vena Cavae

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
Insp Exp	LN	18006-7	Inferior Vena Cava Diameter
%Collapse	LN	18050-5	Inferior Vena Cava % Collapse

CID 12217 Echocardiography Cardiac Shunt

(Local Version: 20200424)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
Qp/Qs	LN	29462-9	Pulmonary-to-Systemic Shunt Flow Ratio
Qp-Qs	99ALOKA	A12217-001	<i>Pulmonary minus Systemic Shunt Flow - Doppler Volume Flow</i>

CID 12218 Echocardiography Congenital

(Local Version: 20200424)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
ASD Ao Rim	99ALOKA	A12218-001	<i>Atrial Septal Defect Rim diameter(Ao Position)</i>
ASD SVC Rim	99ALOKA	A12218-002	<i>Atrial Septal Defect Rim diameter(SVC Position)</i>
ASD IVC Rim	99ALOKA	A12218-003	<i>Atrial Septal Defect Rim diameter(IVC Position)</i>

CID 12220 Echocardiography Common Measurements

(Local Version: 20060807)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
HR	LN	8867-4	Heart rate
RR	DCM	122182	<i>R-R interval</i>

CID 12221 Flow Direction

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	SRT	R-42047	Antegrade Flow
-	SRT	R-42E61	Retrograde Flow

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
MR Vol TR Vol PR Vol AR Vol	LN	33878-0	Volume Flow
MR Flow Rt TR Flow Rt PR Flow Rt AR Flow Rt	LN	34141-2	Peak Instantaneous Flow Rate
CSA(LVOT) CSA(RVOT) MVA MVA(P1/2T) AVA MR EROA TR EROA PR EROA AR EROA	SRT	G-038E	Cardiovascular Orifice Area
LVOT RVOT MV Diam TV Diam AVDs	SRT	G-038F	Cardiovascular Orifice Diameter
MR RF TR RF PR RF AR RF	SRT	G-0390	Regurgitant Fraction
pV	LN	11726-7	Peak Velocity
MnV	LN	20352-1	Mean Velocity
PG	LN	20247-3	Peak Gradient
MPG	LN	20256-4	Mean Gradient
VTI VTI(MVannu) VTI(TVannu)	LN	20354-7	Velocity Time Integral
P1/2T	LN	20280-4	Pressure Half-Time
AccT	LN	20168-1	Acceleration Time

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
PEP	99ALOKA	A12222-001	Pre-Ejection Period
PEP/ET	99ALOKA	A12222-002	PEP/ET
FlowT	99ALOKA	A12222-003	Flow Time
PISA Radius	99ALOKA	A12222-004	Radius of Flow Convergence
MR Alias V (Vr) TR Alias V (Vr) PR Alias V (Vr) AR Alias V (Vr)	99ALOKA	A12222-005	Aliasing Velocity
Angle (PISA)	99ALOKA	A12222-006	Proximal Isovelocity Surface Area Angle
AR VC Diam MR VC Diam TR VC Diam PR VC Diam	99ALOKA	A12222-007	Vena Contracta Diameter
sMnV dMnV	99ALOKA	A12222-008	Mean Velocity (1)
sMnV dMnV	99ALOKA	A12222-009	Mean Velocity (2)
PG1	99ALOKA	A12222-010	Peak Gradient 1
PG2	99ALOKA	A12222-011	Peak Gradient 2
PG1	99ALOKA	A12222-012	Peak Pressure Gradient 1
PG2	99ALOKA	A12222-013	Peak Pressure Gradient 2
ACC	99ALOKA	A12222-014	Acceleration Slope (1)
ACC	99ALOKA	A12222-015	Acceleration Slope (2)
AccT	99ALOKA	A12222-016	Acceleration Time (1)
AccT	99ALOKA	A12222-017	Acceleration Time (2)
Em	99ALOKA	A12222-018	Ann Peak e' wave (1)
Em	99ALOKA	A12222-019	Ann Peak e' wave (2)
Ann-Ost Lt.	99ALOKA	A12222-020	Annulus-Ostium Left diameter
Ann-Ost Rt.	99ALOKA	A12222-021	Annulus-Ostium Right diameter
Ann diam a-p	99ALOKA	A12222-022	Anterior-Posterior Annulus Diameter
APG	99ALOKA	A12222-023	A-Wave Peak Gradient
V1 v1	99ALOKA	A12222-024	Blood Flow Velocity 1
v1	99ALOKA	A12222-025	Blood Flow Velocity 1(1)
v1	99ALOKA	A12222-026	Blood Flow Velocity 1(2)
V2 v2	99ALOKA	A12222-027	Blood Flow Velocity 2

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
v2	99ALOKA	A12222-028	Blood Flow Velocity 2(1)
v2	99ALOKA	A12222-029	Blood Flow Velocity 2(2)
V1	99ALOKA	A12222-030	Blood Flow Velocity at Peak Pressure Gradient 1
V2	99ALOKA	A12222-031	Blood Flow Velocity at Peak Pressure Gradient 2
AVAi(pV) PVAi(VTI) PVAi(pV) TVAi(VTI) TVAi(pV)	99ALOKA	A12222-032	Cardiovascular Orifice Area Index
Ann diam comm Ann diam comm	99ALOKA	A12222-033	Commissural Diameter
PGed PRend PG	99ALOKA	A12222-034	Diastolic Pressure Gradient
Em/Am	99ALOKA	A12222-035	Ea to Aa Ratio
Em/Am	99ALOKA	A12222-036	Ea to Aa Ratio (1)
Em/Am	99ALOKA	A12222-037	Ea to Aa Ratio (2)
dPEP	99ALOKA	A12222-038	Interventricular Dyssynchrony - Interventricular Delay
IVCT RIVCT	99ALOKA	A12222-039	Isovolumic Contraction Time
RIVCT	99ALOKA	A12222-040	Isovolumic Contraction Time (1)
RIVCT	99ALOKA	A12222-041	Isovolumic Contraction Time (2)
IVRT RIVRT	99ALOKA	A12222-042	Isovolumic Relaxation Time
RIVRT	99ALOKA	A12222-043	Isovolumic Relaxation Time (1)
RIVRT	99ALOKA	A12222-044	Isovolumic Relaxation Time (2)
Am	99ALOKA	A12222-045	Peak Diastolic Tissue Velocity During Atrial Systole
Am	99ALOKA	A12222-046	Peak Diastolic Tissue Velocity During Atrial Systole (1)
Am	99ALOKA	A12222-047	Peak Diastolic Tissue Velocity During Atrial Systole (2)
Sm1	99ALOKA	A12222-048	Peak Systolic Tissue Velocity 1
Sm1	99ALOKA	A12222-049	Peak Systolic Tissue Velocity 1(1)
Sm1	99ALOKA	A12222-050	Peak Systolic Tissue Velocity 1(2)
Sm2	99ALOKA	A12222-051	Peak Systolic Tissue Velocity 2
Sm2	99ALOKA	A12222-052	Peak Systolic Tissue Velocity 2(1)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
<i>Sm2</i>	99ALOKA	A12222-053	<i>Peak Systolic Tissue Velocity 2(2)</i>
<i>E/Em</i>	99ALOKA	A12222-054	<i>Ratio of MV E-Wave Peak V to Early Diastolic Myocardium V (1)</i>
<i>E/Em</i>	99ALOKA	A12222-055	<i>Ratio of MV E-Wave Peak V to Early Diastolic Myocardium V (2)</i>
<i>Jet A</i>	99ALOKA	A12222-056	<i>Regurgitation Jet Area</i>
<i>Jet H</i>	99ALOKA	A12222-057	<i>Regurgitation Jet Height</i>
<i>Jet W</i>	99ALOKA	A12222-058	<i>Regurgitation Jet Width</i>
<i>dt time</i>	99ALOKA	A12222-059	<i>Time</i>
<i>time</i>	99ALOKA	A12222-060	<i>Time (1)</i>
<i>time</i>	99ALOKA	A12222-061	<i>Time (2)</i>

CID 12223 Echocardiography Stroke Volume Origin

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	SNM3 ¹	T-32650	Left Ventricle Outflow Tract
-	SNM3	T-32550	Right Ventricle Outflow Tract

Note: 1. Actually "SRT" described in CID 12243 is included in Structured Report, not "SNM3".

CID 12224 Ultrasound Image Modes

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	SRT	G-03A2	2D mode
-	SRT	G-0394	M mode
-	SRT	R-409E4	Doppler Pulsed
-	SRT	R-409E3	Doppler Continuous Wave

CID 12226 Echocardiography Image View

(Local Version: 20200424)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	SRT	G-A19B	Apical two chamber
-	SRT	G-A19C	Apical four chamber
-	SRT	G-0398	Parasternal short axis at the aortic valve level

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	SRT	G-039A	Parasternal short axis at the Mitral Valve level
-	SRT	G-039B	Parasternal short axis at the Papillary Muscle level
-	99ALOKA	A12226-001	Apical three chamber

CID 12227 Echocardiography Measurement Method

(Local Version: 20200424)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	99ALOKA	A12227-001	Valsalva
-	99ALOKA	A12227-002	Mean Pulmonary Arterial Pressure
-	99ALOKA	A12227-003	Left Ventricle Mass for Pediatric
-	99ALOKA	A12227-004	Left Ventricle Mass for Pediatric by Modified Simpson

CID 12228 Volume Methods

(Local Version: 20190123, 20200424)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	DCM	125204	Area-Length Biplane
-	DCM	125205	Area-Length Single Plane
-	DCM	125206	Cube Method
-	DCM	125207	Method of Disks, Biplane
-	DCM	125208	Method of Disks, Single Plane
-	DCM	125209	Teichholz
-	99ALOKA	A12228-001	Method of Disks, Single Plane with Apical two chamber
-	99ALOKA	A12228-002	Gibson
-	99ALOKA	A12228-003	Modified Simpson's
-	99ALOKA	A12228-004	Bullet
-	99ALOKA	A12228-005	Method of Disks, Biplane of LA
-	99ALOKA	A12228-006	Method of Disks, Single Plane with Apical four chamber of LA
-	99ALOKA	A12228-007	Method of Disks, Single Plane with Apical two chamber of LA
-	99ALOKA	A12228-008	Method of Disks, Biplane of RA

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	99ALOKA	A12228-009	<i>Method of Disks, Single Plane with Apical four chamber of RA</i>
-	99ALOKA	A12228-010	<i>Method of Disks, Single Plane with Apical two chamber of RA</i>
-	99ALOKA	A12228-011	<i>Area-Length Biplane of Left Atrium</i>
-	99ALOKA	A12228-012	<i>Area-Length Biplane of Right Atrium</i>
-	99ALOKA	A12228-013	<i>2D Speckle Tracking, Single Plane</i>
-	99ALOKA	A12228-014	<i>2D Speckle Tracking, Multi Plane</i>
-	99ALOKA	A12228-015	<i>Method of Disks, 3 planes</i>

CID 12229 Area Methods

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	DCM	125210	Area by Pressure Half-Time
-	DCM	125215	Continuity Equation by Velocity Time Integral
-	DCM	125216	Proximal Isovelocity Surface Area
-	DCM	125220	Planimetry
-	99ALOKA	A12228-013	<i>2D Speckle Tracking, Single Plane</i>
-	99ALOKA	A12228-014	<i>2D Speckle Tracking, Multi Plane</i>

CID 12230 Gradient Methods

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	DCM	125218	Simplified Bernoulli

CID 12231 Volume Flow Methods

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	DCM	125219	Doppler Volume Flow
-	DCM	125216	Proximal Isovelocity Surface Area

CID 12232 Myocardium Mass Methods

(Local Version: 20060807)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	DCM	125221	Left Ventricle Mass by M-mode
-	99ALOKA	A12232-001	Left Ventricle Mass by Area Length
-	99ALOKA	A12232-002	Left Ventricle Mass by Penn

CID 12233 Cardiac Phase

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	SRT	F-32020	Systole
-	SRT	F-32011	End Diastole
-	DCM	109070	End Systole
-	SRT	R-FAB5B	End Systole

CID 12234 Respiration State

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	SRT	F-20010	During Inspiration
-	SRT	F-20020	During Expiration

CID 12235 Mitral Valve Anatomic Sites

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	SRT	T-35313	Mitral Annulus

CID 12236 Echo Anatomic Sites

(Local Version: 20200424)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	99ALOKA	A12236-001	Proximal Right Ventricle Outflow Tract
-	99ALOKA	A12236-002	Distal Right Ventricle Outflow Tract
-	99ALOKA	A12236-003	Septal Mitral Annulus
-	99ALOKA	A12236-004	Septal Tricuspid Annulus
-	99ALOKA	A12236-005	Lateral Tricuspid Annulus

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	99ALOKA	A12236-006	<i>Epicardium</i>
-	99ALOKA	A12236-007	<i>Ascending Aorta</i>
-	99ALOKA	A12236-008	<i>Descending Aorta</i>
-	99ALOKA	A12236-009	<i>Transverse Aorta</i>
-	99ALOKA	A12236-010	<i>Abdominal Aorta</i>
-	99ALOKA	A12236-011	<i>Patent Foramen Ovale</i>
-	99ALOKA	A12236-012	<i>Patent Ductus Arteriosus</i>
-	99ALOKA	A12236-013	<i>Main Pulmonary Artery</i>
-	99ALOKA	A12236-014	<i>Left Pulmonary Artery</i>
-	99ALOKA	A12236-015	<i>Right Pulmonary Artery</i>
-	99ALOKA	A12236-016	<i>Aortic Valve Annulus</i>

CID 12239 Cardiac Output Properties

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
SV SV(LVOT) SV(RVOT)	SRT	F-32120	Stroke Volume
CO CO(LVOT) CO(RVOT)	SRT	F-32100	Cardiac Output
COI COI(LVOT) COI(RVOT)	SRT	F-32110	Cardiac Index
SVI SVI(LVOT) SVI(RVOT)	SRT	F-00078	Stroke Index

CID 12240 Left Ventricle Area

(Local Version: 20200424)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
LVLAs LVLA4s LVLA2s LVSAMVs LVSAPMs	SRT	G-0374	Left Ventricular Systolic Area
LVLAd LVLA4d LVLA2d LVSAMVd LVSAPMd Aend	SRT	G-0375	Left Ventricular Diastolic Area
areaEF areaEF4 areaEF2 areaEFmv areaEFpm	SRT	G-0376	Left Ventricular Fractional Area Change
Aepi	SRT	G-0379	Left Ventricle Epicardial Diastolic Area, psax pap view
<i>LVAd PM Endo</i>	<i>99ALOKA</i>	<i>A12240-001</i>	<i>Left Ventricular Endocardial Diastolic Area, psax pap view</i>

CID 12241 Tricuspid Valve Finding Sites

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	SRT	T-35111	Tricuspid Annulus

CID 12242 Aortic Valve Finding Sites

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	SRT	T-35410	Aortic Valve Ring

CID 12243 Left Ventricle Finding Sites

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	SRT	T-32650	Left Ventricle Outflow Tract

CID ALOKA_001 OB Original Vascular Measurement Group

This Context Group is private definition.

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
OB Dop1	99ALOKA	ALOKA_001-001	OB User Definition Doppler1
OB Dop2	99ALOKA	ALOKA_001-002	OB User Definition Doppler2
OB Dop3	99ALOKA	ALOKA_001-003	OB User Definition Doppler3

CID ALOKA_002 User Defined Equation Group

This Context Group is private definition.

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
User defined equation result parameter	99ALOKA	ALOKA_002-001	Calculated result
User defined equation parameter	99ALOKA	ALOKA_002-002	Equation parameter

CID ALOKA_003 Common Measurement Items

This Context Group is private definition.

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
Length(Gallbladder) Length(Renal Volume)	SRT	G-D7FE	Length
Width(Gallbladder) Width(Renal Volume)	SRT	G-A220	Width

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
b(Spleen Index) Wall Thickness(GB Wall-T)	LN	59089-3	Thickness
Aorta Diam PV Diam	SRT	M-02550	Diameter
a(Spleen Index)	SRT	G-A185	Long axis
A-P(Renal Volume)	SRT	R-10206	Antero-posterior
Area(SOL1,2)	SRT	G-A166	Area
Volume(SOL1,2) Volume(Renal Volume)	SRT	G-D705	Volume
Vessel(%STEN O Diam)	SRT	G-0364	Vessel Lumen Diameter
Vessel(%STEN O Area)	SRT	G-0366	Vessel Lumen Cross-Sectional Area
<i>Residual(%STE NO Diam)</i>	99ALOKA	ALOKA_003-001	<i>Residual Vessel Lumen Diameter</i>
<i>%Diam Steno (%STENO Diam)</i>	99ALOKA	ALOKA_003-002	<i>Stenosis Ratio calculated by Vessel and Residual Diameter</i>
<i>Residual(%STE NO Area)</i>	99ALOKA	ALOKA_003-003	<i>Residual Vessel Lumen Cross-Sectional Area</i>
<i>%Area Steno (%STENO Area)</i>	99ALOKA	ALOKA_003-004	<i>Stenosis Ratio calculated by Vessel and Residual Area</i>
<i>Diam1(SOL1,2)</i>	99ALOKA	ALOKA_003-005	<i>Diameter1</i>
<i>Diam2(SOL1,2)</i>	99ALOKA	ALOKA_003-006	<i>Diameter2</i>
<i>Diam3(SOL1,2)</i>	99ALOKA	ALOKA_003-007	<i>Diameter3</i>
<i>L1(Liver)</i>	99ALOKA	ALOKA_003-008	<i>Length1</i>
<i>L2(Liver)</i>	99ALOKA	ALOKA_003-009	<i>Length2</i>

CID ALOKA_004 Echocardiography Analysis Measurement Items

(Local Version: 20190123)

This Context Group is private definition.

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
Peak Global	99ALOKA	ALOKA_004-001	Peak GLS
Peak Global	99ALOKA	ALOKA_004-002	Peak GLS of Overlay

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
-	99ALOKA	ALOKA_004-003	Averaged Peak GLS
L. Strain(EN)[%]	99ALOKA	ALOKA_004-004	GLS-Sel.
L. Strain(EN)[%]	99ALOKA	ALOKA_004-005	GLS-Sel. of Overlay
Time from R	99ALOKA	ALOKA_004-006	Time from R to Peak GLS phase
Time from R	99ALOKA	ALOKA_004-007	Time from R to Peak GLS phase of Overlay
Time from R	99ALOKA	ALOKA_004-008	Time from R to GLS-Sel. phase
Time from R	99ALOKA	ALOKA_004-009	Time from R to GLS-Sel. phase of Overlay

CID ALOKA_005 Abdomen Original Measurements Group (Local Version: 20190301)

This Context Group is private definition.

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
Gallbladder	SRT	T-63000	Gall Bladder
Liver	SRT	T-62000	Liver
Pancreas	SRT	T-65000	Pancreas
P-Duct	SRT	T-65010	Pancreatic duct
Spleen	SRT	T-C3000	Spleen
%STENO Diam	99ALOKA	ALOKA_005-001	Vessel Diameter Stenosis Ratio
%STENO Area	99ALOKA	ALOKA_005-002	Vessel Area Stenosis Ratio
GB Wall-T	99ALOKA	ALOKA_005-003	Gall Bladder Wall
CBD	99ALOKA	ALOKA_005-004	Common Bile Duct
Renal Volume	99ALOKA	ALOKA_005-005	Renal Volume
SOL1	99ALOKA	ALOKA_005-006	Space Occupying Lesion1
SOL2	99ALOKA	ALOKA_005-007	Space Occupying Lesion2
Aorta Diam	99ALOKA	ALOKA_005-008	Aorta Diameter
PV Diam	99ALOKA	ALOKA_005-009	Portal Vein Diameter

CID ALOKA_006 Abdomen Original Measurement Items (Local Version: 20190301)

This Context Group is private definition.

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
CBD	99ALOKA	ALOKA_006-001	Common Bile Duct
Head	99ALOKA	ALOKA_006-002	Pancreas Head

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
Body	99ALOKA	ALOKA_006-003	Pancreas Body
Tail	99ALOKA	ALOKA_006-004	Pancreas Tail
P-Duct	99ALOKA	ALOKA_006-005	Pancreas Duct
SI(AXB)	99ALOKA	ALOKA_006-006	Spleen Index

8.6.4 Private Code Definitions

This section specifies the meanings of private codes used in Structured Reports.

Private Code Definitions (Non-Advanced Cardiac, Coding Scheme Designator: "99ALOKA", Coding Scheme Version: not specified)

Code Value	Code Meaning	Definition
A-001	Pre-prandial	Pre-prandial
A-002	OB Original Fetal Measurement Group	Concept Name of TID ALOKA-5000 row1.
A-003	OB-GYN Original Measurement Group	Concept Name of TID ALOKA-5001 row1.
A-004	User Defined Equation	Concept Name of TID ALOKA-301 row1.
A-005	User Defined Equation Group Number	Concept Name of TID ALOKA-301 row2.
A-006	User Defined Equation Group Name	Concept Name of TID ALOKA-301 row3.
A-007	User Defined Parameter Number	Concept Name of TID ALOKA-302 row2.
A-008	User Defined Parameter Name	Concept Name of TID ALOKA-302 row3.
A-009	User Defined Equation	Concept Name of TID ALOKA-302 row4.
A-010	OB-GYN User Defined Measurement Group	Concept Name of TID ALOKA-5002 row1
A-011	Original Measurement Group of Abodomen	Concept Name of TID ALOKA-5100 row2
A3663-001	BSA = $0.0003207 \times WT^{(0.7285 - 0.0188 \log(WT))} \times HT^{0.3}$	BSA = $0.0003207 \times WT^{(0.7285 - 0.0188 \log(WT))} \times HT^{0.3}$ Body Surface Area computed from patient height(HT) and weight(WT). The formula is derived by Boyd : $BSA = 0.0003207 \times WT[g]^{(0.7285 - 0.0188 \log(WT[g]))} \times HT[cm]^{0.3}$ Reference : Boyd E. The growth of the surface area of the human body. Minneapolis:university of Minnesota Press, 1935.

Code Value	Code Meaning	Definition
A3663-002	BSA = 0.007358*HT^0.725*WT^0.425	BSA = 0.007358*HT^0.725*WT^0.425 Body Surface Area computed from patient height(HT) and weight(WT). The formula is derived by Shintani : BSA = 0.007358*HT[cm]^0.725*WT[kg]^0.425 Reference : 臨床検査法概要 29版, 金井泉 著, 金原出版
A12004-001	Cephalic Index(BPDo/OFDo)	Cephalic Index = BPDo/OFDo
A12004-002	LVW/HW	LVW/HW
A12005-001	Binocular Distance	Binocular Distance
A12005-002	Biparietal Diameter outer-to-outer	Biparietal Diameter outer-to-outer
A12005-003	Fetal Trunk Cross Sectional Area	Fetal Trunk Cross Sectional Area
A12005-004	Length of Vertebrae	Length of Vertebrae
A12005-005	Occipital-Frontal Diameter outer-to-outer	Occipital-Frontal Diameter outer-to-outer
A12005-006	Abdominal Diameter	Abdominal Diameter
A12005-007	Thoracic Length	Thoracic Length
A12005-008	Head Circumference for Merz, Hansmann	Head Circumference for Merz, Hansmann
A12005-009	Amniotic Fluid Volume	Amniotic Fluid Volume
A12005-010	Head Circumference for BMUS	Head Circumference for BMUS
A12005-011	Abdominal Circumference for BMUS	Abdominal Circumference for BMUS
A12005-012	Maximum Vertical Pocket	"Ultrasonography in Obstetrics and Gynecology 5th Edition" Peter Callen November 28, 2007. pp.760-767
A12006-001	Nasal Bone Length	Nasal Bone Length
A12009-001	Early Embryonic Size	Early Embryonic Size
A12009-002	Gestational Sac Diameter 1	Gestational Sac Diameter 1
A12009-003	Gestational Sac Diameter 2	Gestational Sac Diameter 2
A12009-004	Gestational Sac Diameter 3	Gestational Sac Diameter 3
A12009-005	Mean Gestational Sac Diameter	Mean Gestational Sac Diameter
A12011-001	Cervix Width	Cervix Width
A12011-002	Cervix Antero-Posterior Diameter	Cervix Antero-Posterior Diameter
A12011-003	Pre Void Bladder Length	Pre Void Bladder Length

Code Value	Code Meaning	Definition
A12011-004	Pre Void Bladder Antero-Posterior Diameter	Pre Void Bladder Antero-Posterior Diameter
A12011-005	Pre Void Bladder Width	Pre Void Bladder Width
A12011-006	Pre Void Bladder Volume	Pre Void Bladder Volume
A12011-007	Post Void Bladder Length	Post Void Bladder Length
A12011-008	Post Void Bladder Antero-Posterior Diameter	Post Void Bladder Antero-Posterior Diameter
A12011-009	Post Void Bladder Width	Post Void Bladder Width
A12011-010	Post Void Bladder Volume	Post Void Bladder Volume
A12011-011	Bladder Void Volume	Bladder Void Volume
A12013-001	AC, Campbell	Materials provided: Professor Campbell's Group at Harris Birthright Centre, King's College Hospital
A12013-002	AC, Chitty 1994	Charts of fetal size : 3. Abdominal measurements Lyn S Chitty. British Journal of Obstetrics and Gynaecology February 1994, Vol.101, pp.125-131 <Table 4>
A12013-003	AC, Hadlock 1982	Fetal Abdominal Circumference as a Predictor of Menstrual Age. Hadlock FP, Deter RL, Harrist RB, Park SK. AJR 139:367-370, August 1982
A12013-005	AC, Merz 1996	Das normale fetale Wachstumsprofil - ein einheitliches Modell zur Berechnung von Normkurven für die gängigen Kopf-und Abdomenparameter sowie die großen Extremitätenknochen. Ultraschall in Med. 17(1996),153 - 162 Table Data : 95 percentile data form «Growth format»
A12013-006	APAD, Merz 1996	Das normale fetale Wachstumsprofil - ein einheitliches Modell zur Berechnung von Normkurven für die gängigen Kopf-und Abdomenparameter sowie die großen Extremitätenknochen. Ultraschall in Med. 17(1996),153 - 162 Table Data : 95 percentile data form «Growth format»
A12013-007	AxT, Tokyo	胎児生理の総合的解析による新しい周産期管理へのアプローチ. 東京大学 岡井 崇他. 日本産婦人科学会雑誌 第38巻 第8号 別冊
A12013-008	BD, Jeanty 1984	Estimation of Gestational Age from Measurements of Fetal Long Bones. Jeanty P, Rodesch F, Delbeke D, Dumont JE. Journal of Ultrasound in Medicine 3:75-79, February 1984
A12013-009	BPD, Campbell	Materials provided: Professor Campbell's Group at Harris Birthright Centre, King's College Hospital
A12013-010	BPD-oi, Chitty 1994	Charts of fetal size : 2. Head measurements. Lyn S Chitty. British Journal of Obstetrics and Gynaecology February 1994, Vol.101, pp.35-43 <Table 4,7>
A12013-011	BPD, Hadlock 1982	Fetal Biparietal Diameter : A Critical Re-evaluation of the Relation to Menstrual Age by means of Real-time Ultrasound. Hadlock FP, Deter RL, Harrist RB, Park SK : Journal of Ultrasound in Medicine 1:97, 97-104

Code Value	Code Meaning	Definition
A12013-012	BPD, Merz 1996	Das normale fetale Wachstumsprofil - ein einheitliches Modell zur Berechnung von Normkurven für die gängigen Kopf-und Abdomenparameter sowie die großen Extremitätenknochen. Ultraschall in Med. 17(1996),153 - 162 Table Data : 95 percentile data form «Growth format»
A12013-013	BPD, Sabbagha 1976	Sonar Biparietal Diameter : I. Analysis of Percentile Growth Differences in Two Normal Populations Using Same Methodology. Sabbagha RE, Barton FB, Barton BA. American Journal of Obstetrics and Gynecology 126:479-484, October 1976
A12013-016	CRL, JSUM 2003	超音波胎児計測の標準化と日本人の基準値の公示について. J Med Ultrasonics Vol.30 No.3 2003
A12013-017	CRL, Tokyo	胎児生理の総合的解析による新しい周産期管理へのアプローチ. 東京大学 岡井 崇他. 日本産婦人科学会雑誌 第38巻 第8号 別冊
A12013-018	EES, Goldstein 1994	Endovaginal Ultrasonographic Measurement of Early Embryonic Size as a Means of Assessing Gestational Age. Steven R. Goldstein, MD, Robert Wolfson, MD, PhD. J. Ultrasound Med. 13:27-31, 1994. <Figure 3>
A12013-019	FIB, Merz 1996	Das normale fetale Wachstumsprofil - ein einheitliches Modell zur Berechnung von Normkurven für die gängigen Kopf-und Abdomenparameter sowie die großen Extremitätenknochen. Ultraschall in Med. 17(1996),153 - 162 Table Data : 95 percentile data form «Growth format»
A12013-020	FL, Campbell	Materials provided: Professor Campbell's Group at Harris Birthright Centre, King's College Hospital
A12013-021	FL, Chitty 1994	Charts of fetal size : 4. Femur length. Lyn S Chitty. British Journal of Obstetrics and Gynaecology February 1994, Vol.101, pp.132-135 <Table 2>
A12013-022	FL, Hadlock 1982	Fetal Femur Length as a Predictor of Menstrual Age : Sonographically Measured. Hadlock FP, Deter RL, Harrist RB, Park SK. AJR 138:875-878, May 1982
A12013-025	FL, Jeanty 95% 1983	『Fetal limb biometry』 Radiology 1983 ; 147 : 602. Table Data : 95 percentile data form «Growth format»
A12013-026	FL, Merz 1996	Das normale fetale Wachstumsprofil - ein einheitliches Modell zur Berechnung von Normkurven für die gängigen Kopf-und Abdomenparameter sowie die großen Extremitätenknochen. Ultraschall in Med. 17(1996),153 - 162 Table Data : 95 percentile data form «Growth format»
A12013-027	FL, Warda 1985	Fetal Femur Length : A Critical Reevaluation of the Relationship to Menstrual Age. Warda AH, Deter RL, Rossavik IK, Carpenter RJ, Hadlock FP. American Journal of Obstetrics and Gynecology 66(1):69-75, July 1985

Code Value	Code Meaning	Definition
A12013-028	GS, Tokyo	胎児生理の総合的解析による新しい周産期管理へのアプローチ. 東京大学 岡井 崇他. 日本産婦人科学会雑誌 第38巻 第8号 別冊
A12013-029	HC, Campbell	Materials provided: Professor Campbell's Group at Harris Birthright Centre, King's College Hospital
A12013-030	HC, Chitty 1994	Charts of fetal size : 2. Head measurements. Lyn S Chitty. British Journal of Obstetrics and Gynaecology February 1994, Vol.101, pp.35-43 <Table 4,7>
A12013-031	HC, Hadlock 1982	Fetal Head Circumference : Relation to Menstrual Age. Hadlock FP, Deter RL, Harrist RB, Park SK. AJR 138:649-653, April 1982
A12013-033	HC, Merz 1996	Das normale fetale Wachstumsprofil - ein einheitliches Modell zur Berechnung von Normkurven für die gängigen Kopf-und Abdomenparameter sowie die großen Extremitätenknochen. Ultraschall in Med. 17(1996),153 - 162 Table Data : 95 percentile data form «Growth format»
A12013-034	Humerus, Hansmann 1985	Ultrasound Diagnosis in Obstetrics and Gynecology. Hansmann M., Hackeloer B.J. and Staudach A. Springer-Verlag, Berlin, Heidelberg, New York, Tokyo 1985
A12013-035	Humerus, Jeanty 95% 1983	『Fetal limb biometry』 Radiology 1983 ; 147 : 602. Table Data : 95 percentile data form «Growth format»
A12013-036	Humerus, Merz 1996	Das normale fetale Wachstumsprofil - ein einheitliches Modell zur Berechnung von Normkurven für die gängigen Kopf-und Abdomenparameter sowie die großen Extremitätenknochen. Ultraschall in Med. 17(1996),153 - 162 Table Data : 95 percentile data form «Growth format»
A12013-037	NBL, Sonek 2003	Nasal bone length throughout gestation : normal ranges based on 3537 fetal ultrasound measurements. J. D. SONEK. Ultrasound Obstet Gynecol 2003 ; 21 ; 152-155
A12013-038	OFD, Merz 1996	Das normale fetale Wachstumsprofil - ein einheitliches Modell zur Berechnung von Normkurven für die gängigen Kopf-und Abdomenparameter sowie die großen Extremitätenknochen. Ultraschall in Med. 17(1996),153 - 162 Table Data : 95 percentile data form «Growth format»
A12013-039	Radius, Merz 1996	Das normale fetale Wachstumsprofil - ein einheitliches Modell zur Berechnung von Normkurven für die gängigen Kopf-und Abdomenparameter sowie die großen Extremitätenknochen. Ultraschall in Med. 17(1996),153 - 162 Table Data : 95 percentile data form «Growth format»
A12013-040	Tibia, Jeanty 95% 1983	『Fetal limb biometry』 Radiology 1983 ; 147 : 602. Table Data : 95 percentile data form «Growth format»

Code Value	Code Meaning	Definition
A12013-041	Tibia, Merz 1996	Das normale fetale Wachstumsprofil - ein einheitliches Modell zur Berechnung von Normkurven für die gängigen Kopf-und Abdomenparameter sowie die großen Extremitätenknochen. Ultraschall in Med. 17(1996),153 - 162 Table Data : 95 percentile data form 《Growth format》
A12013-042	TL, Chitkara 1987	Prenatal sonographic assessment of the fetal thorax: Normal values. Usha Chitkara, M.D., Joanne Rosenberg, R.D.M.S., Frank A. Chervenak, M.D., Gertrud S. Berkowitz, Ph.D., Rebecca Levine, M.A., Richard M. Fagerstrom, Ph.D., Barbara Walker, R.D.M.S., and Richard L. Berkowitz, M.D. American Journal of Obstetrics and Gynecology, Volume 156, Number 5, May 1987, pp.1069-1074. <Table 2>
A12013-043	TAD, Hansmann 1985	Ultrasound Diagnosis in Obstetrics and Gynecology. Hansmann M., Hackeloer B.J. and Staudach A. Springer-Verlag, Berlin, Heidelberg, New York, Tokyo 1985
A12013-044	TAD, Merz 1996	Das normale fetale Wachstumsprofil - ein einheitliches Modell zur Berechnung von Normkurven für die gängigen Kopf-und Abdomenparameter sowie die großen Extremitätenknochen. Ultraschall in Med. 17(1996),153 - 162 Table Data : 95 percentile data form 《Growth format》
A12013-045	Ulna, Jeanty 95% 1983	『Fetal limb biometry』 Radiology 1983 ; 147 : 602. Table Data : 95 percentile data form 《Growth format》
A12013-046	Ulna, Merz 1996	Das normale fetale Wachstumsprofil - ein einheitliches Modell zur Berechnung von Normkurven für die gängigen Kopf-und Abdomenparameter sowie die großen Extremitätenknochen. Ultraschall in Med. 17(1996),153 - 162 Table Data : 95 percentile data form 《Growth format》
A12013-047	AC, JSUM 2003	超音波胎児計測の標準化と日本人の基準値の公示について. J Med Ultrasonics Vol.30 No.3 2003
A12013-048	BPD, JSUM 2003	超音波胎児計測の標準化と日本人の基準値の公示について. J Med Ultrasonics Vol.30 No.3 2003
A12013-049	BPD, Osaka	超音波胎児体格計測によるIUGR診断. 青木 嶺夫. Perinatal Care Vol.9 No.5,(407-422)
A12013-050	CRL, Osaka	超音波胎児体格計測によるIUGR診断. 青木 嶺夫. Perinatal Care Vol.9 No.5,(407-422)
A12013-051	FL, JSUM 2003	超音波胎児計測の標準化と日本人の基準値の公示について. J Med Ultrasonics Vol.30 No.3 2003
A12013-052	FL, O'Brien 1981	Assessment of Gestational Age in the Second Trimester by Real-Time Ultrasound Measurement of the Femur Length. O'Brien GD, Queenan JT, Campbell S (American Journal of Obstetrics & Gynecology 139:540-545, Mar. 1981) Table Data : 《Growth format》

Code Value	Code Meaning	Definition
A12013-053	FL, Osaka	超音波胎児体格計測によるIUGR診断. 青木 嶺夫. Perinatal Care Vol.9 No.5,(407-422)
A12013-054	FTA, Osaka	超音波胎児体格計測によるIUGR診断. 青木 嶺夫. Perinatal Care Vol.9 No.5,(407-422)
A12013-055	Humerus, Osaka	超音波胎児体格計測によるIUGR診断. 青木 嶺夫. Perinatal Care Vol.9 No.5,(407-422)
A12013-056	CRL, BMUS 2009	[Fetal Size and Dating: Charts Recommended for Clinical Obstetric Practice] Pam Lougha, Lyn Chitty, Tony Evans & Trish Chudleigh Ultrasound 2009;17(3):160-166 British Medical Ultrasound Society:2009
A12013-057	HC, BMUS 2009	[Fetal Size and Dating: Charts Recommended for Clinical Obstetric Practice] Pam Lougha, Lyn Chitty, Tony Evans & Trish Chudleigh Ultrasound 2009;17(3):160-166 British Medical Ultrasound Society:2009
A12013-058	FL, BMUS 2009	[Fetal Size and Dating: Charts Recommended for Clinical Obstetric Practice] Pam Lougha, Lyn Chitty, Tony Evans & Trish Chudleigh Ultrasound 2009;17(3):160-166 British Medical Ultrasound Society:2009
A12013-059	HC, BMUS-95% 2009	[Fetal Size and Dating: Charts Recommended for Clinical Obstetric Practice] Pam Lougha, Lyn Chitty, Tony Evans & Trish Chudleigh Ultrasound 2009;17(3):160-166 British Medical Ultrasound Society:2009
A12013-060	AC, BMUS-95% 2009	[Fetal Size and Dating: Charts Recommended for Clinical Obstetric Practice] Pam Lougha, Lyn Chitty, Tony Evans & Trish Chudleigh Ultrasound 2009;17(3):160-166 British Medical Ultrasound Society:2009
A12013-061	FL, BMUS-95% 2009	[Fetal Size and Dating: Charts Recommended for Clinical Obstetric Practice] Pam Lougha, Lyn Chitty, Tony Evans & Trish Chudleigh Ultrasound 2009;17(3):160-166 British Medical Ultrasound Society:2009
A12013-062	BPD, CFEF-97% 2006	[French fetal biometry: reference equations and comparison with other charts.] Salomon LJ, Duyme M, Crequat J, Brodaty G, Talmant C, Fries N, Althuser M Ultrasound Obstet Gynecol. 2006 Aug;28(2):193-8

Code Value	Code Meaning	Definition
A12013-063	HC, CFEF-97% 2006	[French fetal biometry: reference equations and comparison with other charts.] Salomon LJ, Duyme M, Crequat J, Brodaty G, Talmant C, Fries N, Althuser M Ultrasound Obstet Gynecol. 2006 Aug;28(2):193-8
A12013-064	AC, CFEF-97% 2006	[French fetal biometry: reference equations and comparison with other charts.] Salomon LJ, Duyme M, Crequat J, Brodaty G, Talmant C, Fries N, Althuser M Ultrasound Obstet Gynecol. 2006 Aug;28(2):193-8
A12013-066	FL, CFEF-97% 2006	[French fetal biometry: reference equations and comparison with other charts.] Salomon LJ, Duyme M, Crequat J, Brodaty G, Talmant C, Fries N, Althuser M Ultrasound Obstet Gynecol. 2006 Aug;28(2):193-8
A12013-067	BPD, CFEF-90% 2006	[French fetal biometry: reference equations and comparison with other charts.] Salomon LJ, Duyme M, Crequat J, Brodaty G, Talmant C, Fries N, Althuser M Ultrasound Obstet Gynecol. 2006 Aug;28(2):193-8
A12013-068	HC, CFEF-90% 2006	[French fetal biometry: reference equations and comparison with other charts.] Salomon LJ, Duyme M, Crequat J, Brodaty G, Talmant C, Fries N, Althuser M Ultrasound Obstet Gynecol. 2006 Aug;28(2):193-8
A12013-069	AC, CFEF-90% 2006	[French fetal biometry: reference equations and comparison with other charts.] Salomon LJ, Duyme M, Crequat J, Brodaty G, Talmant C, Fries N, Althuser M Ultrasound Obstet Gynecol. 2006 Aug;28(2):193-8
A12013-071	FL, CFEF-90% 2006	[French fetal biometry: reference equations and comparison with other charts.] Salomon LJ, Duyme M, Crequat J, Brodaty G, Talmant C, Fries N, Althuser M Ultrasound Obstet Gynecol. 2006 Aug;28(2):193-8
A12013-072	User defined Table or Equation1	User defined Table or Equation1
A12013-073	User defined Table or Equation2	User defined Table or Equation2
A12013-074	User defined Table or Equation3	User defined Table or Equation3
A12013-075	User defined Table or Equation4	User defined Table or Equation4

Code Value	Code Meaning	Definition
A12013-076	User defined Table or Equation5	User defined Table or Equation5
A12013-077	User defined Table or Equation6	User defined Table or Equation6
A12013-078	User defined Table or Equation7	User defined Table or Equation7
A12013-079	User defined Table or Equation8	User defined Table or Equation8
A12013-080	User defined Table or Equation9	User defined Table or Equation9
A12013-081	User defined Table or Equation10	User defined Table or Equation10
A12013-082	User defined Table or Equation11	User defined Table or Equation11
A12013-083	User defined Table or Equation12	User defined Table or Equation12
A12013-084	User defined Table or Equation13	User defined Table or Equation13
A12013-085	User defined Table or Equation14	User defined Table or Equation14
A12013-086	User defined Table or Equation15	User defined Table or Equation15
A12013-087	User defined Table or Equation16	User defined Table or Equation16
A12013-088	User defined Table or Equation17	User defined Table or Equation17
A12013-089	User defined Table or Equation18	User defined Table or Equation18
A12013-090	User defined Table or Equation19	User defined Table or Equation19
A12013-091	User defined Table or Equation20	User defined Table or Equation20
A12014-001	EFW by BPD, AC, FL, Hadlock	Estimation of fetal weight with the use of head,body,and femur measurement - A prospective study. Frank P. Hadlock, R.B.Harrist, Ralph S.Sharman, Russel L Deter, and Seung K.Park. Am J Obstet Gynecol :Volume151 Number3: 333-337, February1,1985. Sonographic Estimation of Fetal weight. Frank P. Hadlock, R.B.Harrist, Robert J.Carpenter, Russel L Deter, Seung K.Park. Radiology Volume150 Number2:535-540
A12014-002	EFW by BPD, AC, FL, JSUM 2003	超音波胎児計測の標準化と日本人の基準値の公示について. J Med Ultrasonics Vol.30 No.3 2003
A12014-003	EFW by BPD, FTA, FL, Osaka	超音波胎児体格計測によるIUGR診断. 青木 嶺夫. Perinatal Care Vol.9 No.5,(407-422)

Code Value	Code Meaning	Definition
A12014-004	EFW by BPD, AC, Warsof 1977	The estimation of fetal weight by computer-assisted analysis. Steven L.Warsof, Parviz Gohari, Richard L. Berkowitz, John C.Hobbins. Am J Obstet Gynecol :Volume 128 Number 8 : 881-892, August 15, 1977
A12014-005	EFW by User defined Equation1	EFW by User defined Equation1
A12014-006	EFW by User defined Equation2	EFW by User defined Equation2
A12014-007	EFW by User defined Equation3	EFW by User defined Equation3
A12014-008	EFW by User defined Equation4	EFW by User defined Equation4
A12014-009	EFW by User defined Equation5	EFW by User defined Equation5
A12015-001	FL/AC by GA, Hadlock 1983	A Date-Independent Predictor of Intrauterine Growth Retardation : Femur Length/Abdominal Circumference Ratio. Hadlock FP, Deter RL, Harrist RB, Roecker E, Park SK. American Journal of Roentgenology 141:979-984, November 1983
A12015-002	FL/HC by GA, Hadlock 1984	The Femur Length/Head Circumference Relation in Obstetric Sonography. Frank P. Hadlock, MD, Ronald B. Harrist, PhD, Yogesh Shah, MD, Seung K. Park, MD. Journal of Ultrasound in Medicine, Volume 3, October 1984, pp.439-442. <Table 1>
A12015-003	FL/BPD by GA, Hohler 1981	Comparison of Ultrasound Femur Length and Biparietal Diameter in Late Pregnancy. Hohler CW, Quetel TA. American Journal of Obstetrics and Gynecology 141:759-762, December 1981
A12015-004	AFI by GA, Jeng et al.	Amniotic Fluid Index Measurement with the Four-Quadrant Technique During Pregnancy. Cherng-Jye Jeng, M.D., Tian-Jii Jou, M.D., Kuo-Gon Wang, M.D., Yuh-Cheng Yang, M.D., Yi-Nan Lee, M.D., Chung-Chi Lan, M.D. The Journal of Reproductive Medicine, Volume 35, Number 7, July 1990, pp.674-677. <Table 1>
A12015-005	AFI by GA, Moore et al.	The amniotic fluid index in normal human pregnancy. Thomas R. Moore, MD, and Jonathan E. Cayle, MD. American Journal of Obstetrics and Gynecology, Volume 162, Number 5, May 1990, pp.1168-1173. <Table 6>
A12015-006	AFI by GA, Phelan et al.	Amniotic Fluid Volume Assessment with the Four-Quadrant Technique at 36-42 Weeks' Gestation. Jeffrey P. Phelan, M.D., Carl Vernon Smith, M.D., Paula Broussard, R.N., Mary Small, M.D. The Journal of Reproductive Medicine, Volume 32, Number 7, July 1987, pp.540-542. <Table 1>
A12015-007	AC by GA, JSUM 2003	超音波胎児計測の標準化と日本人の基準値の公示について. J Med Ultrasonics Vol.30 No.3 2003

Code Value	Code Meaning	Definition
A12015-008	BPD by GA, JSUM 2003	超音波胎児計測の標準化と日本人の基準値の公示について. J Med Ultrasonics Vol.30 No.3 2003
A12015-009	BPD by GA, Osaka	超音波胎児体格計測によるIUGR診断. 青木 嶺夫. Perinatal Care Vol.9 No.5,(407-422)
A12015-010	CRL by GA, Osaka	超音波胎児体格計測によるIUGR診断. 青木 嶺夫. Perinatal Care Vol.9 No.5,(407-422)
A12015-011	FL by GA, JSUM 2003	超音波胎児計測の標準化と日本人の基準値の公示について. J Med Ultrasonics Vol.30 No.3 2003
A12015-012	FL by GA, O'Brien 1981	Assessment of Gestational Age in the Second Trimester by Real-Time Ultrasound Measurement of the Femur Length. O'Brien GD, Queenan JT, Campbell S (American Journal of Obstetrics & Gynecology 139:540-545, Mar. 1981) Table Data : 《Growth format》
A12015-013	FL by GA, Osaka	超音波胎児体格計測によるIUGR診断. 青木 嶺夫. Perinatal Care Vol.9 No.5,(407-422)
A12015-014	FTA by GA, Osaka	超音波胎児体格計測によるIUGR診断. 青木 嶺夫. Perinatal Care Vol.9 No.5,(407-422)
A12015-015	Humerus by GA, Osaka	超音波胎児体格計測によるIUGR診断. 青木 嶺夫. Perinatal Care Vol.9 No.5,(407-422)
A12015-016	BPD by GA, CFEF-97% 2006	[French fetal biometry: reference equations and comparison with other charts.] Salomon LJ, Duyme M, Crequat J, Brodaty G, Talmant C, Fries N, Althuser M Ultrasound Obstet Gynecol. 2006 Aug;28(2):193-8
A12015-017	HC by GA, CFEF-97% 2006	[French fetal biometry: reference equations and comparison with other charts.] Salomon LJ, Duyme M, Crequat J, Brodaty G, Talmant C, Fries N, Althuser M Ultrasound Obstet Gynecol. 2006 Aug;28(2):193-8
A12015-018	AC by GA, CFEF-97% 2006	[French fetal biometry: reference equations and comparison with other charts.] Salomon LJ, Duyme M, Crequat J, Brodaty G, Talmant C, Fries N, Althuser M Ultrasound Obstet Gynecol. 2006 Aug;28(2):193-8
A12015-020	FL by GA, CFEF-97% 2006	[French fetal biometry: reference equations and comparison with other charts.] Salomon LJ, Duyme M, Crequat J, Brodaty G, Talmant C, Fries N, Althuser M Ultrasound Obstet Gynecol. 2006 Aug;28(2):193-8
A12015-021	BPD by GA, CFEF-90% 2006	[French fetal biometry: reference equations and comparison with other charts.] Salomon LJ, Duyme M, Crequat J, Brodaty G, Talmant C, Fries N, Althuser M Ultrasound Obstet Gynecol. 2006 Aug;28(2):193-8

Code Value	Code Meaning	Definition
A12015-022	HC by GA, CFEF-90% 2006	[French fetal biometry: reference equations and comparison with other charts.] Salomon LJ, Duyme M, Crequat J, Brodaty G, Talmant C, Fries N, Althuser M Ultrasound Obstet Gynecol. 2006 Aug;28(2):193-8
A12015-023	AC by GA, CFEF-90% 2006	[French fetal biometry: reference equations and comparison with other charts.] Salomon LJ, Duyme M, Crequat J, Brodaty G, Talmant C, Fries N, Althuser M Ultrasound Obstet Gynecol. 2006 Aug;28(2):193-8
A12015-025	FL by GA, CFEF-90% 2006	[French fetal biometry: reference equations and comparison with other charts.] Salomon LJ, Duyme M, Crequat J, Brodaty G, Talmant C, Fries N, Althuser M Ultrasound Obstet Gynecol. 2006 Aug;28(2):193-8
A12015-026	GA by GA, User defined Table or Equation1	GA by GA, User defined Table or Equation1
A12015-027	GA by GA, User defined Table or Equation2	GA by GA, User defined Table or Equation2
A12015-028	GA by GA, User defined Table or Equation3	GA by GA, User defined Table or Equation3
A12015-029	GA by GA, User defined Table or Equation4	GA by GA, User defined Table or Equation4
A12015-030	GA by GA, User defined Table or Equation5	GA by GA, User defined Table or Equation5
A12015-031	GA by GA, User defined Table or Equation6	GA by GA, User defined Table or Equation6
A12015-032	GA by GA, User defined Table or Equation7	GA by GA, User defined Table or Equation7
A12015-033	GA by GA, User defined Table or Equation8	GA by GA, User defined Table or Equation8
A12015-034	GA by GA, User defined Table or Equation9	GA by GA, User defined Table or Equation9
A12015-035	GA by GA, User defined Table or Equation10	GA by GA, User defined Table or Equation10
A12015-036	GA by GA, User defined Table or Equation11	GA by GA, User defined Table or Equation11
A12015-037	GA by GA, User defined Table or Equation12	GA by GA, User defined Table or Equation12
A12015-038	GA by GA, User defined Table or Equation13	GA by GA, User defined Table or Equation13
A12015-039	GA by GA, User defined Table or Equation14	GA by GA, User defined Table or Equation14
A12015-040	GA by GA, User defined Table or Equation15	GA by GA, User defined Table or Equation15

Code Value	Code Meaning	Definition
A12015-041	GA by GA, User defined Table or Equation16	GA by GA, User defined Table or Equation16
A12015-042	GA by GA, User defined Table or Equation17	GA by GA, User defined Table or Equation17
A12015-043	GA by GA, User defined Table or Equation18	GA by GA, User defined Table or Equation18
A12015-044	GA by GA, User defined Table or Equation19	GA by GA, User defined Table or Equation19
A12015-045	GA by GA, User defined Table or Equation20	GA by GA, User defined Table or Equation20
A12015-046	Fetus Ratio, User defined Table or Equation1	Fetus Ratio, User defined Table or Equation1
A12015-047	Fetus Ratio, User defined Table or Equation2	Fetus Ratio, User defined Table or Equation2
A12015-048	Fetus Ratio, User defined Table or Equation3	Fetus Ratio, User defined Table or Equation3
A12015-049	Fetus Ratio, User defined Table or Equation4	Fetus Ratio, User defined Table or Equation4
A12015-050	Fetus Ratio, User defined Table or Equation5	Fetus Ratio, User defined Table or Equation5
A12015-051	Fetus Ratio, User defined Table or Equation6	Fetus Ratio, User defined Table or Equation6
A12015-052	AFI by User defined Table or Equation1	AFI by User defined Table or Equation1
A12015-053	AFI by User defined Table or Equation2	AFI by User defined Table or Equation2
A12016-001	FW, Doubilet 1997	Improved Birth Weight Table for Neonates Developed from Gestations Dated by Early Ultrasonography: Doubilet PM et al; J Ultrasound Med 16;241-249, 1997
A12016-005	Twins FW, Yarkoni 1987	Estimated Fetal Weight in the Evaluation of Growth in Twin Gestations: A Prospective Longitudinal Study. Shaul Yarkoni, MD, E. Albert Reece, MD, Theodore Holford, PhD, Theresa Z. O'Connor, MPH, And John C. Hobbins, MD. Obstetrics & Gynecology, Volume 69, Number 4, April 1987, pp.636-639.
A12016-006	FW, CFEF-97% 2014	Courbe d'estimation de poids fœtal 2014 par le Collège français d'échographie fœtale (CFEF) M. Massoud, M. Duyme, M. Fontanges, D. Combourieu, Collège français d'échographie fœtale (CFEF) Journal de Gynécologie Obstétrique et Biologie de la Reproduction - accepté le 15 janvier 2015.

Code Value	Code Meaning	Definition
A12016-007	FW, CFEF-90% 2014	Courbe d'estimation de poids fœtal 2014 par le Collège français d'échographie fœtale (CFEF) M. Massoud, M. Duyme, M. Fontanges, D. Combourieu, Collège français d'échographie fœtale (CFEF) Journal de Gynécologie Obstétrique et Biologie de la Reproduction - accepté le 15 janvier 2015.
A12016-008	FW, User defined Table or Equation1	FW, User defined Table or Equation1
A12016-009	FW, User defined Table or Equation2	FW, User defined Table or Equation2
A12016-010	FW, User defined Table or Equation3	FW, User defined Table or Equation3
A12016-011	FW, User defined Table or Equation4	FW, User defined Table or Equation4
A12016-012	FW, User defined Table or Equation5	FW, User defined Table or Equation5
A12019-001	Fetal Heart Rate before Biopsy	Fetal Heart Rate before Biopsy
A12019-002	Fetal Heart Rate after Biopsy	Fetal Heart Rate after Biopsy
A12019-003	Auto Fetal Heart Rate	Automated FHR Measurement
A12019-004	Estimated Fetal Weight Ratio	EFW Ratio = (Estimated Fetal Weight / Weightest Fetal Weight) * 100 [%]
A12104-001	Maximum Intima-media thickness	Maximum Intima-media thickness
A12104-002	Mean Intima-media thickness	Mean Intima-media thickness
A12104-003	Intima-media thickness of CCA by auto-detection	Intima-media thickness of Common Carotid Artery by auto-detection
A12104-004	Intima-media thickness of ICA by auto-detection	Intima-media thickness of Internal Carotid Artery by auto-detection
A12104-005	Intima-media thickness of ECA by auto-detection	Intima-media thickness of External Carotid Artery by auto-detection
A12104-006	Intima-media thickness of BIFUR by auto-detection	Intima-media thickness of Carotid Bifurcation by auto-detection
A12104-007	Mean Carotid Intima-media thickness by auto-detection	Mean Carotid Intima-media thickness by auto-detection
A12104-008	IMT-C10	IMT-C10 “超音波による頸動脈病変の標準的評価法 2017” JSUM0515 GUIDELINE
A12104-009	IMT-C10 by auto-detection	IMT-C10 by auto-detection
A12107-001	Deep Brachial Artery	Deep Brachial Artery
A12107-002	Basilic Artery	Basilic Artery

Code Value	Code Meaning	Definition
A12107-003	User Definition Artery1	Upr Art.1 (User Definition Artery)
A12107-004	User Definition Artery2	Upr Art.2 (User Definition Artery)
A12107-005	User Definition Artery3	Upr Art.3 (User Definition Artery)
A12107-006	User Definition Artery4	Upr Art.4 (User Definition Artery)
A12107-007	User Definition Artery5	Upr Art.5 (User Definition Artery)
A12107-008	User Definition Artery6	Upr Art.6 (User Definition Artery)
A12107-009	User Definition Artery7	Upr Art.7 (User Definition Artery)
A12107-010	User Definition Artery8	Upr Art.8 (User Definition Artery)
A12108-001	Deep Brachial vein	Deep Brachial Vein
A12108-002	User Definition Vein1	Upr Vein.1 (User Definition Vein)
A12108-003	User Definition Vein2	Upr Vein.2 (User Definition Vein)
A12108-004	User Definition Vein3	Upr Vein.3 (User Definition Vein)
A12108-005	User Definition Vein4	Upr Vein.4 (User Definition Vein)
A12108-006	User Definition Vein5	Upr Vein.5 (User Definition Vein)
A12108-007	User Definition Vein6	Upr Vein.6 (User Definition Vein)
A12108-008	User Definition Vein7	Upr Vein.7 (User Definition Vein)
A12108-009	User Definition Vein8	Upr Vein.8 (User Definition Vein)
A12109-001	User Definition Artery1	Lwr Art.1 (User Definition Artery)
A12109-002	User Definition Artery2	Lwr Art.2 (User Definition Artery)
A12109-003	User Definition Artery3	Lwr Art.3 (User Definition Artery)
A12109-004	User Definition Artery4	Lwr Art.4 (User Definition Artery)
A12109-005	User Definition Artery5	Lwr Art.5 (User Definition Artery)
A12109-006	User Definition Artery6	Lwr Art.6 (User Definition Artery)
A12109-007	User Definition Artery7	Lwr Art.7 (User Definition Artery)
A12109-008	User Definition Artery8	Lwr Art.8 (User Definition Artery)
A12110-001	User Definition Vein1	Lwr Vein.1 (User Definition Vein)
A12110-002	User Definition Vein2	Lwr Vein.2 (User Definition Vein)
A12110-003	User Definition Vein3	Lwr Vein.3 (User Definition Vein)
A12110-004	User Definition Vein4	Lwr Vein.4 (User Definition Vein)
A12110-005	User Definition Vein5	Lwr Vein.5 (User Definition Vein)
A12110-006	User Definition Vein6	Lwr Vein.6 (User Definition Vein)
A12110-007	User Definition Vein7	Lwr Vein.7 (User Definition Vein)
A12110-008	User Definition Vein8	Lwr Vein.8 (User Definition Vein)
A12112-001	User Definition Artery1	Artery1 (User Definition Artery (Abdom))
A12112-002	User Definition Artery2	Artery2 (User Definition Artery (Abdom))

Code Value	Code Meaning	Definition
A12112-003	User Definition Artery3	Artery3 (User Definition Artery (Abdom))
A12117-001	far wall	far wall “超音波による頸動脈病変の標準的評価法 2017” JSUM0515 GUIDELINE
A12117-002	near wall	near wall “超音波による頸動脈病変の標準的評価法 2017” JSUM0515 GUIDELINE
A12120-001	atrial contraction	日本超音波医学会用語・診断基準委員会 “胎児静脈血流波形基準値(2013) 胎児静脈管(ductus venosus:DV) PI, a/Sおよび胎児下大静脈(inferior vena cava:IVC) PLI(preload index), a/Sの計測法と週数毎の基準値” In Europe: K. Hecher, S. Campbell, R. Snijders and K. Nicolaides Ultrasound Obstet. Gynecol. 4 (1994) 381-390
A12120-002	Peak Diastolic Velocity	[Reference ranges for fetal venous and atrioventricular blood flow parameters] K. Hecher, S. Campbell, R. Snijders and K. Nicolaides Ultrasound Obstet. Gynecol. 4 (1994) 381-390
A12121-001	Acceleration Time to Flow Time Ratio	Acceleration Time to Flow Time Ratio
A12121-002	atrial contraction to Systolic peak Velocity Ratio	日本超音波医学会用語・診断基準委員会 “胎児静脈血流波形基準値(2013) 胎児静脈管(ductus venosus:DV) PI, a/Sおよび胎児下大静脈(inferior vena cava:IVC) PLI(preload index), a/Sの計測法と週数毎の基準値” In Europe: “Reference ranges for fetal venous and atrioventricular blood flow parameters” K. Hecher, S. Campbell, R. Snijders and K. Nicolaides Ultrasound Obstet. Gynecol. 4 (1994) 381-390

Code Value	Code Meaning	Definition
A12121-003	Systolic peak to atrial contraction Velocity Ratio	日本超音波医学会用語・診断基準委員会 “胎児静脈血流波形基準値(2013) 胎児静脈管(ductus venosus:DV) PI, a/Sおよび胎児下大静脈(inferior vena cava:IVC) PLI(preload index), a/Sの計測法と週数毎の基準値” In Europe: “Reference ranges for fetal venous and atrioventricular blood flow parameters” K. Hecher, S. Campbell, R. Snijders and K. Nicolaides Ultrasound Obstet. Gynecol. 4 (1994) 381-390
A12121-004	Peak velocity index for veins	[Reference ranges for fetal venous and atrioventricular blood flow parameters] K. Hecher, S. Campbell, R. Snijders and K. Nicolaides Ultrasound Obstet. Gynecol. 4 (1994) 381-390
A12122-001	Flow Time	Flow Time
A12122-002	Intima-media thickness 1	Intima-media thickness 1
A12122-003	Intima-media thickness 2	Intima-media thickness 2
A12122-004	Intima-media thickness 3	Intima-media thickness 3
A12122-005	Intima-media thickness 4	Intima-media thickness 4
A12122-006	Intima-media thickness 5	Intima-media thickness 5
A12122-007	Intima-media thickness 6	Intima-media thickness 6
A12122-008	Intima-media thickness 7	Intima-media thickness 7
A12122-009	Intima-media thickness 8	Intima-media thickness 8
A12122-010	Intima-media thickness 9	Intima-media thickness 9
A12122-011	Intima-media thickness 10	Intima-media thickness 10
A12122-012	Maximum Intima-media thickness	Mean Intima-media thickness
A12122-013	Mean Intima-media thickness	Mean Intima-media thickness
A12122-014	Maximum Intima-media thickness at diastole	Maximum Intima-media thickness at diastole
A12122-015	Minimum Intima-media thickness at diastole	Minimum Intima-media thickness at diastole
A12122-016	Average Intima-media thickness at diastole	Average Intima-media thickness at diastole
A12122-017	IMT-C10 at diastole	IMT-C10 at diastole
A12122-018	Standard Deviation of Intima-media thickness at diastole	Standard Deviation of Intima-media thickness at diastole

Code Value	Code Meaning	Definition
A12122-019	Measured number of Intima-media thickness points at diastole	Measured number of Intima-media thickness points at diastole
A12122-020	Width of ROI at diastole	Width of ROI at diastole
A12122-021	Maximum Intima-media thickness at systole	Maximum Intima-media thickness at systole
A12122-022	Minimum Intima-media thickness at systole	Minimum Intima-media thickness at systole
A12122-023	Average Intima-media thickness at systole	Average Intima-media thickness at systole
A12122-024	IMT-C10 at systole	IMT-C10 at systole
A12122-025	Standard Deviation of Intima-media thickness at systole	Standard Deviation of Intima-media thickness at systole
A12122-026	Measured number of Intima-media thickness points at systole	Measured number of Intima-media thickness points at systole
A12122-027	Width of ROI at systole	Width of ROI at systole
A12122-028	Intima-media thickness 1 at diastole	Intima-media thickness 1 at diastole
A12122-029	Intima-media thickness 2 at diastole	Intima-media thickness 2 at diastole
A12122-030	Intima-media thickness 3 at diastole	Intima-media thickness 3 at diastole
A12122-031	Mean Intima-media thickness at diastole	Mean Intima-media thickness at diastole
A12122-032	Intima-media thickness 1 at systole	Intima-media thickness 1 at systole
A12122-033	Intima-media thickness 2 at systole	Intima-media thickness 2 at systole
A12122-034	Intima-media thickness 3 at systole	Intima-media thickness 3 at systole
A12122-035	Mean Intima-media thickness at systole	Mean Intima-media thickness at systole
A12122-036	IMT-C10	IMT-C10 “超音波による頸動脈病変の標準的評価法 2017” JSUM0515 GUIDELINE

Code Value	Code Meaning	Definition
A12141-001	Ductus Venosus	日本超音波医学会用語・診断基準委員会 “胎児静脈血流波形基準値(2013) 胎児静脈管(ductus venosus:DV) PI, a/Sおよび胎児下大静脈(inferior vena cava:IVC) PLI(preload index), a/Sの計測法と週数毎の基準値” In Europe: “Reference ranges for fetal venous and atrioventricular blood flow parameters” K. Hecher, S. Campbell, R. Snijders and K. Nicolaides Ultrasound Obstet. Gynecol. 4 (1994) 381-390
A12201-001	Left Ventricular Short Axis Length at Mitral Valve	Left Ventricular Short Axis Length at Mitral Valve
A12201-002	Mean Wall Thickness	Mean Wall Thickness
A12201-003	Left Ventricular End Diastolic Length	Left Ventricular End Diastolic Length
A12201-004	Short Axis Radius Calculated from Short Axis Cavity Area	Short Axis Radius Calculated from Short Axis Cavity Area
A12203-001	Left Ventricular Mass Index	Left Ventricular Mass Index (LVMI). $LVMI = LVM / BSA$ LVMI unit is g/m^2, mass unit in gram and BSA unit in m^2.
A12203-003	Long Axis(at End Diastole or End Systole) Length % Difference	Long Axis(at End Diastole or End Systole) Length Percentage Difference. Length Percentage Difference unit is %.
A12203-004	Midwall Fractional Shortening	Midwall Fractional Shortening (mFS)
A12203-005	Area Ejection Fraction at Long Axis View	Area Ejection Fraction at Long Axis View
A12203-006	Area Ejection Fraction at Short Axis View	Area Ejection Fraction at Short Axis View
A12203-009	Mean Velocity of Circumferential Fiber Shortening	Mean Velocity of Circumferential Fiber Shortening
A12203-010	Left Ventricle Diastole Filling Time	Left Ventricle Diastole Filling Time
A12203-011	Ratio of MV E-Wave Peak Vel. to Early Diastolic Myocardium Vel.	Ratio of Mitral Valve E-Wave Peak Velocity to Early Diastolic Myocardium Velocity
A12203-012	Ratio of Left Ventricle Diastole Filling Time to R-R interval	Ratio of Left Ventricle Diastole Filling Time to R-R interval

Code Value	Code Meaning	Definition
A12203-013	Global Longitudinal Strain	Global Longitudinal Strain Bi-Plane GLS = { [(A4C ES Circum-length - A4C ED Circum-length) / A4C ED Circum-length] + [(A2C ES Circum-length - A2C ED Circum-length) / A2C ED Circum-length] } * 50 4ch GLS = { (A4C ES Circum-length - A4C ED Circum-length) / A4C ED Circum-length } * 100 2ch GLS = { (A2C ES Circum-length - A2C ED Circum-length) / A2C ED Circum-length } * 100
A12203-014	pV Ratio of LVOT to AV	pV Ratio of LVOT to AV
A12203-015	VTI Ratio of LVOT to AV	VTI Ratio of LVOT to AV
A12203-016	Left Ventricle Mass adjusted by Height(1.7)	Left Ventricle Mass adjusted by Height(1.7)
A12203-017	Left Ventricle Mass adjusted by Height(2.16)	Left Ventricle Mass adjusted by Height(2.16)
A12203-018	Left Ventricle Mass adjusted by Height(2.7)	Left Ventricle Mass adjusted by Height(2.7)
A12203-019	Myocardial Length	Myocardial Length
A12204-001	Right Ventricular Diastolic Area	Right Ventricular Diastolic Area
A12204-002	Right Ventricular Systolic Area	Right Ventricular Systolic Area
A12204-003	Right Ventricular Fractional Area Change	RV FAC = 100 * (RV Area ED - RV Area ES) / RV Area ED. Unit is %.
A12204-004	End-diastolic Basal Diameter	End-diastolic Basal Diameter
A12204-005	End-diastolic Mid-cavity Diameter	End-diastolic Mid-cavity Diameter
A12204-006	Right Ventricular End-diastolic Length	Right Ventricular End-diastolic Length
A12204-007	Right Ventricular Wall Thickness	Right Ventricular Wall Thickness
A12205-001	Left Atrium Antero-posterior Diastolic Dimension	Left Atrium Antero-posterior Diastolic Dimension
A12205-002	Left Atrium systolic major axis	Left Atrium systolic major axis
A12205-003	Left Atrium Systolic Area by Apical two chamber	Left Atrium Systolic Area by Apical two chamber
A12205-004	Left Atrial Volume	Left Atrial Volume
A12205-005	Left Atrial Volume divided by Body Surface Area	Left Atrial Volume divided by Body Surface Area
A12205-006	Long Axis at End Systole Length % Difference of Left Atrium	Long Axis at End Systole Length Percentage Difference of Left Atrium

Code Value	Code Meaning	Definition
A12205-007	Left Atrium systolic minor axis	Left Atrium systolic minor axis
A12206-001	Right Atrium systolic major axis	Right Atrium systolic major axis
A12206-002	Right Atrium Systolic Area by Apical two chamber	Right Atrium Systolic Area by Apical two chamber
A12206-003	Right Atrial Volume	Right Atrial Volume
A12206-004	Right Atrial Volume divided by Body Surface Area	Right Atrial Volume divided by Body Surface Area
A12206-005	Long Axis at End Systole Length % Difference of Right Atrium	Long Axis at End Systole Length Percentage Difference of Right Atrium
A12206-006	Right Atrium systolic minor axis	Right Atrium systolic minor axis
A12207-001	Mitral Valve Dimension of C point to E point by M-Mode	Mitral Valve Dimension of C point to E point by M-Mode
A12207-002	Mitral Valve Dimension of C point to A point by M-Mode	Mitral Valve Dimension of C point to A point by M-Mode
A12207-003	Mitral Valve C-A Dimension to C-E Dimension Ratio by M-Mode	Mitral Valve C-A Dimension to C-E Dimension Ratio by M-Mode
A12207-004	Mitral Valve C-E Dimension to C-A Dimension Ratio by M-Mode	Mitral Valve C-E Dimension to C-A Dimension Ratio by M-Mode
A12207-005	Mitral Valve E-Wave Duration	Mitral Valve E-Wave Duration
A12207-006	Mitral Valve A to E Ratio	Mitral Valve A to E Ratio
A12207-007	Mitral Valve E-wave Peak Pressure Gradient	Mitral Valve E-wave Peak Pressure Gradient
A12207-008	Mitral Valve A-wave Peak Pressure Gradient	Mitral Valve A-wave Peak Pressure Gradient
A12207-009	Subtraction of A-wave Duration from PVA-wave Duration	Subtraction of Mitral A-wave Duration from Pulmonary Vein Atrial Reversal Duration
A12207-010	Velocity Time Integral of Mitral Regurgitant Flow	Velocity Time Integral of Mitral Regurgitant Flow
A12207-012	Mitral Regurgitant Proximal Isovelocity Surface Area	Mitral Regurgitant Proximal Isovelocity Surface Area
A12207-013	Flow Volume of Mitral Valve Annulus in Flow	Flow Volume of Mitral Valve Annulus in Flow
A12207-014	Mitral valve Regurgitant volume	Mitral valve Regurgitant volume
A12207-015	Isovolumic Acceleration	Isovolumic Acceleration
A12207-016	Average of e'(l) and e'(s)	Average of e'(l) and e'(s)
A12207-017	D-E amplitude	D-E amplitude

Code Value	Code Meaning	Definition
A12207-018	Mitral Annular Plane Systolic Excursion	Mitral Annular Plane Systolic Excursion
A12207-019	Mitral Annulus Area	Mitral Annulus Area
A12207-020	Mitral Annulus Diameter PLAX	Mitral Annulus Diameter PLAX
A12207-021	Mitral Valve Ea to Aa Ratio	Mitral Valve Ea to Aa Ratio
A12207-022	Peak MV Annular Velocity in Isovolumic Contraction	Peak MV Annular Velocity in Isovolumic Contraction
A12207-023	Ratio of MV Peak Velocity to LV av e'	Ratio of MV Peak Velocity to LV av e'
A12207-024	Time E to e'	Time E to e'
A12207-025	Time to peak velocity in isovolumic contractuion	Time to peak velocity in isovolumic contractuion
A12207-026	Flow Volume of Mitral Valve Annulus in Flow (1)	Flow Volume of Mitral Valve Annulus in Flow (1)
A12208-001	Tricuspid Valve Dimension of C point to E point by M-Mode	Tricuspid Valve Dimension of C point to E point by M-Mode
A12208-002	Tricuspid Valve Dimension of C point to A point by M-Mode	Tricuspid Valve Dimension of C point to A point by M-Mode
A12208-003	Tricuspid Valve Dimension of D point to E point by M-Mode	Tricuspid Valve Dimension of D point to E point by M-Mode
A12208-004	Tricuspid Valve Velocity from E point to F point by M-Mode	Tricuspid Valve Velocity from E point to F point by M-Mode
A12208-005	Tricuspid Valve Velocity from D point to E point by M-Mode	Tricuspid Valve Velocity from D point to E point by M-Mode
A12208-006	Tricuspid Valve C-A Dimension to C-E Dimension Ratio by M-Mode	Tricuspid Valve C-A Dimension to C-E Dimension Ratio by M-Mode
A12208-007	Tricuspid Valve C-E Dimension to C-A Dimension Ratio by M-Mode	Tricuspid Valve C-E Dimension to C-A Dimension Ratio by M-Mode
A12208-008	Velocity Time Integral of Tricuspid Regurgitant Flow	Velocity Time Integral of Tricuspid Regurgitant Flow
A12208-010	Tricuspid Regurgitant Proximal Isovelocity Surface Area	Tricuspid Regurgitant Proximal Isovelocity Surface Area
A12208-011	Flow Volume of Tricuspid Valve Annulus in Flow	Flow Volume of Tricuspid Valve Annulus in Flow
A12208-012	Tricuspid Annular Plane Systolic Excursion	Tricuspid Annular Plane Systolic Excursion

Code Value	Code Meaning	Definition
A12208-013	Peak A wave/Peak E wave by US	Peak A wave/Peak E wave by US
A12208-014	E-Wave Duration	E-Wave Duration
A12208-015	E-Wave Peak Gradient	E-Wave Peak Gradient
A12208-016	Ratio of TV Peak Velocity to RV Peak Tissue Velocity E-Wave	Ratio of TV Peak Velocity to RV Peak Tissue Velocity E-Wave
A12208-017	Right Ventricular Peak Early Diastolic Tissue Velocity	Right Ventricular Peak Early Diastolic Tissue Velocity
A12208-018	Right Ventricular Peak Systolic Tissue Velocity	Right Ventricular Peak Systolic Tissue Velocity
A12208-020	RV Peak Diastolic Tissue Velocity During Atrial Systole	RV Peak Diastolic Tissue Velocity During Atrial Systole
A12208-021	Time Difference between v1 and v2	Time Difference between v1 and v2
A12208-022	Tricuspid Annulus Area	Tricuspid Annulus Area
A12208-023	Tricuspid Valve Closure to Opening Time	Tricuspid Valve Closure to Opening Time
A12209-001	Pulmonic Valve Dimension of F point to A point by M-Mode	Pulmonic Valve Dimension of F point to A point by M-Mode
A12209-002	Pulmonic Valve Dimension of B point to C point by M-Mode	Pulmonic Valve Dimension of B point to C point by M-Mode
A12209-003	Pulmonic Valve Velocity from E point to F point by M-Mode	Pulmonic Valve Velocity from E point to F point by M-Mode
A12209-004	Pulmonic Valve Velocity from B point to C point by M-Mode	Pulmonic Valve Velocity from B point to C point by M-Mode
A12209-005	Velocity Time Integral of Pulmonic Regurgitant Flow	Velocity Time Integral of Pulmonic Regurgitant Flow
A12209-007	Pulmonic Regurgitant Proximal Isovelocity Surface Area	Pulmonic Regurgitant Proximal Isovelocity Surface Area
A12209-009	Transverse Annulus Diameter	Transverse Annulus Diameter
A12209-010	Duration of Diastole	Duration of Diastole
A12209-011	Ratio of PR dur to PR Index	Ratio of PR dur to PR Index
A12209-012	Ratio of PR VCW to Ann diam a-p	Ratio of PR VCW to Ann diam a-p
A12209-013	Pulmonary Artery Pressure	Pulmonary Artery Pressure
A12209-014	Pulmonary Valve Regurgitation Duration	Pulmonary Valve Regurgitation Duration

Code Value	Code Meaning	Definition
A12209-015	Pulmonary Valve Regurgitation Vena Contracta Diameter	Pulmonary Valve Regurgitation Vena Contracta Diameter
A12209-016	Purmonary Annulus Area	Purmonary Annulus Area
A12209-017	Velocity Ratio of RVOT velocity to Pulmonic Valve velocity	Velocity Ratio of RVOT velocity to Pulmonic Valve velocity
A12211-001	Velocity Time Integral of Aortic Regurgitant Flow	Velocity Time Integral of Aortic Regurgitant Flow
A12211-003	Aortic Regurgitant Proximal Isovelocity Surface Area	Aortic Regurgitant Proximal Isovelocity Surface Area
A12211-004	Aortic Annulus Diameter	Aortic Annulus Diameter
A12211-005	Aortic Annulus Maximun Diamter	Aortic Annulus Maximun Diamter
A12211-006	Aortic Annulus Minimum Diamter	Aortic Annulus Minimum Diamter
A12211-007	Aortic Regurgitant Jet Area/LVOT Area	Aortic Regurgitant Jet Area/LVOT Area
A12211-008	Aortic Regurgitant Jet Width/LVOT Width	Aortic Regurgitant Jet Width/LVOT Width
A12211-009	Aortic Regurgitation Vena Contracta Width	Aortic Regurgitation Vena Contracta Width
A12211-010	Aortic Valve Annulus Area	Aortic Valve Annulus Area
A12211-011	Velocity Ratio of LVOT to AV	Velocity Ratio of LVOT to AV
A12212-001	Sinus of Valsalva	Sinus of Valsalva
A12212-002	Sinotubular junction	Sinotubular junction
A12212-003	Aortic Arch	Aortic Arch
A12212-004	Ascending Aorta	Ascending Aorta
A12212-005	Descending Aorta	Descending Aorta
A12212-006	Abdominal Aorta	Abdominal Aorta
A12212-007	Hepatic Vein	Hepatic Vein
A12212-008	Aortic Sinotubular Junction Diameter	Aortic Sinotubular Junction Diameter
A12212-009	Aortic Sinuses of Valsalva Diameter	Aortic Sinuses of Valsalva Diameter
A12214-001	Deceleration Time of D-Wave Flow	Deceleration Time of D-Wave Flow
A12214-002	Systolic Fraction	Systolic Fraction = $[S-VTI / (S-VTI + D-VTI)] * 100$. Unit is %.
A12217-001	Pulmonary minus Systemic Shunt Flow - Doppler Volume Flow	Pulmonary minus Systemic Shunt Flow - Doppler Volume Flow

Code Value	Code Meaning	Definition
A12218-001	Atrial Septal Defect Rim diameter(Ao Position)	Atrial Septal Defect Rim diameter(Ao Position)
A12218-002	Atrial Septal Defect Rim diameter(SVC Position)	Atrial Septal Defect Rim diameter(SVC Position)
A12218-003	Atrial Septal Defect Rim diameter(IVC Position)	Atrial Septal Defect Rim diameter(IVC Position)
A12222-001	Pre-Ejection Period	Pre-Ejection Period
A12222-002	PEP/ET	Ratio of Pre-Ejection Period to Ejection Time
A12222-003	Flow Time	Flow Time
A12222-004	Radius of Flow Convergence	Radius of Flow Convergence
A12222-005	Aliasing Velocity	Aliasing Velocity
A12222-006	Proximal Isovelocity Surface Area Angle	Proximal Isovelocity Surface Area Angle
A12222-007	Vena Contracta Diameter	Vena Contracta Diameter
A12222-008	Mean Velocity (1)	Mean Velocity (1)
A12222-009	Mean Velocity (2)	Mean Velocity (2)
A12222-010	Peak Gradient 1	Peak Gradient 1
A12222-011	Peak Gradient 2	Peak Gradient 2
A12222-012	Peak Pressure Gradient 1	Peak Pressure Gradient 1
A12222-013	Peak Pressure Gradient 2	Peak Pressure Gradient 2
A12222-014	Acceleration Slope (1)	Acceleration Slope (1)
A12222-015	Acceleration Slope (2)	Acceleration Slope (2)
A12222-016	Acceleration Time (1)	Acceleration Time (1)
A12222-017	Acceleration Time (2)	Acceleration Time (2)
A12222-018	Ann Peak e' wave (1)	Ann Peak e' wave (1)
A12222-019	Ann Peak e' wave (2)	Ann Peak e' wave (2)
A12222-020	Annulus-Ostium Left diameter	Annulus-Ostium Left diameter
A12222-021	Annulus-Ostium Right diameter	Annulus-Ostium Right diameter
A12222-022	Anterior-Posterior Annulus Diameter	Anterior-Posterior Annulus Diameter
A12222-023	A-Wave Peak Gradient	A-Wave Peak Gradient
A12222-024	Blood Flow Velocity 1	Blood Flow Velocity 1
A12222-025	Blood Flow Velocity 1(1)	Blood Flow Velocity 1(1)
A12222-026	Blood Flow Velocity 1(2)	Blood Flow Velocity 1(2)
A12222-027	Blood Flow Velocity 2	Blood Flow Velocity 2
A12222-028	Blood Flow Velocity 2(1)	Blood Flow Velocity 2(1)
A12222-029	Blood Flow Velocity 2(2)	Blood Flow Velocity 2(2)

Code Value	Code Meaning	Definition
A12222-030	Blood Flow Velocity at Peak Pressure Gradient 1	Blood Flow Velocity at Peak Pressure Gradient 1
A12222-031	Blood Flow Velocity at Peak Pressure Gradient 2	Blood Flow Velocity at Peak Pressure Gradient 2
A12222-032	Cardiovascular Orifice Area Index	Cardiovascular Orifice Area Index
A12222-033	Commissural Diameter	Commissural Diameter
A12222-034	Diastolic Pressure Gradient	Diastolic Pressure Gradient
A12222-035	Ea to Aa Ratio	Ea to Aa Ratio
A12222-036	Ea to Aa Ratio (1)	Ea to Aa Ratio (1)
A12222-037	Ea to Aa Ratio (2)	Ea to Aa Ratio (2)
A12222-038	Interventricular Dyssynchrony - Interventricular Delay	Interventricular Dyssynchrony - Interventricular Delay
A12222-039	Isovolumic Contraction Time	Isovolumic Contraction Time
A12222-040	Isovolumic Contraction Time (1)	Isovolumic Contraction Time (1)
A12222-041	Isovolumic Contraction Time (2)	Isovolumic Contraction Time (2)
A12222-042	Isovolumic Relaxation Time	Isovolumic Relaxation Time
A12222-043	Isovolumic Relaxation Time (1)	Isovolumic Relaxation Time (1)
A12222-044	Isovolumic Relaxation Time (2)	Isovolumic Relaxation Time (2)
A12222-045	Peak Diastolic Tissue Velocity During Atrial Systole	Peak Diastolic Tissue Velocity During Atrial Systole
A12222-046	Peak Diastolic Tissue Velocity During Atrial Systole (1)	Peak Diastolic Tissue Velocity During Atrial Systole (1)
A12222-047	Peak Diastolic Tissue Velocity During Atrial Systole (2)	Peak Diastolic Tissue Velocity During Atrial Systole (2)
A12222-048	Peak Systolic Tissue Velocity 1	Peak Systolic Tissue Velocity (1)
A12222-049	Peak Systolic Tissue Velocity 1(1)	Peak Systolic Tissue Velocity 1(1)
A12222-050	Peak Systolic Tissue Velocity 1(2)	Peak Systolic Tissue Velocity 1(2)
A12222-051	Peak Systolic Tissue Velocity 2	Peak Systolic Tissue Velocity (2)
A12222-052	Peak Systolic Tissue Velocity 2(1)	Peak Systolic Tissue Velocity 2(1)

Code Value	Code Meaning	Definition
A12222-053	Peak Systolic Tissue Velocity 2(2)	Peak Systolic Tissue Velocity 2(2)
A12222-054	Ratio of MV E-Wave Peak V to Early Diastolic Myocardium V (1)	Ratio of MV E-Wave Peak V to Early Diastolic Myocardium V (1)
A12222-055	Ratio of MV E-Wave Peak V to Early Diastolic Myocardium V (2)	Ratio of MV E-Wave Peak V to Early Diastolic Myocardium V (2)
A12222-056	Regurgitation Jet Area	Regurgitation Jet Area
A12222-057	Regurgitation Jet Height	Regurgitation Jet Height
A12222-058	Regurgitation Jet Width	Regurgitation Jet Width
A12222-059	Time	Time
A12222-060	Time (1)	Time (1)
A12222-061	Time (2)	Time (2)
A12226-001	Apical three chamber	Apical three chamber
A12227-001	Valsalva	Valsalva
A12227-002	Mean Pulmonary Arterial Pressure	Mean Pulmonary Arterial Pressure
A12227-003	Left Ventricle Mass for Pediatric	Left Ventricle Mass for Pediatric
A12227-004	Left Ventricle Mass for Pediatric by Modified Simpson	Left Ventricle Mass for Pediatric by Modified Simpson
A12228-001	Method of Disks, Single Plane with Apical two chamber	Method of Disks, Single Plane with Apical two chamber
A12228-002	Gibson	EDV and ESV are calculated as follows. $EDV = \pi / 6 * LVIDd^2 * (0.98 * LVIDd + 5.90)$ $ESV = \pi / 6 * LVIDs^2 * (1.14 * LVIDs + 4.18)$ Volume unit is milliliter and length in cm. Reference : Gibson, D. G. Measurement of left ventricular volumes in man by echocardiography – comparison with biplane angiographs. Br.Heart J. 1971;33:614.

Code Value	Code Meaning	Definition
A12228-003	Modified Simpson's	Volume = $(LVL / 9) * (4 * LVSAMV + 2 * LVSAPM + (LVSAMV * LVSAPM)^{(1 / 2)})$ Volume unit is milliliter, length in cm and area in cm². References : Folland, ED, et al. Assessment of Left Ventricular Ejection Fraction and Volumes by Real-Time, Two-Dimensional Echocardiography. Circulation,1979;60:760-766. A. F. Parisi, MD et al. Approaches to Determination of Left Ventricular Volume and Ejection Fraction by Real-Time Two-Dimensional Echocardiography. Clin.Cardiol.2,257-263(1979).
A12228-004	Bullet	Volume = $(5 * LVSAPM * LVL) / 6$ Volume unit is milliliter, length in cm and area in cm².
A12228-005	Method of Disks, Biplane of LA	Method of Disks, Biplane of Left Atrium
A12228-006	Method of Disks, Single Plane with Apical four chamber of LA	Method of Disks, Single Plane with Apical four chamber of Left Atrium
A12228-007	Method of Disks, Single Plane with Apical two chamber of LA	Method of Disks, Single Plane with Apical two chamber of Left Atrium
A12228-008	Method of Disks, Biplane of RA	Method of Disks, Biplane of Right Atrium
A12228-009	Method of Disks, Single Plane with Apical four chamber of RA	Method of Disks, Single Plane with Apical four chamber of Right Atrium
A12228-010	Method of Disks, Single Plane with Apical two chamber of RA	Method of Disks, Single Plane with Apical two chamber of Right Atrium
A12228-011	Area-Length Biplane of Left Atrium	Area-Length Biplane of Left Atrium
A12228-012	Area-Length Biplane of Right Atrium	Area-Length Biplane of Right Atrium
A12228-013	2D Speckle Tracking, Single Plane	2D Speckle Tracking, Single Plane display Speckle Tracking Definitions for a Common Standard for 2D Speckle Tracking Echocardiography: Consensus Document of the EACVI/ASE/Industry Task Force to Standardize Deformation Imaging(J Am Soc Echocardiogr 2015;28:183-93.)(European Heart Journal – Cardiovascular Imaging (2015) 16, 1–11)

Code Value	Code Meaning	Definition
A12228-014	2D Speckle Tracking, Multi Plane	2D Speckle Tracking, Multi Plane display Speckle Tracking Definitions for a Common Standard for 2D Speckle Tracking Echocardiography: Consensus Document of the EACVI/ASE/Industry Task Force to Standardize Deformation Imaging(J Am Soc Echocardiogr 2015;28:183-93.)(European Heart Journal – Cardiovascular Imaging (2015) 16, 1–11)
A12228-015	Method of Disks, 3 planes	Mean value calculated by results from 3 chamber view, 4 chamber view, and 2 chamber view.
A12232-001	Left Ventricle Mass by Area Length	$\text{Mass} = 1.05 * ((5 * A_{epi} * (LVLd + \text{thick}) / 6) - (5 * A_{end} * LVLd / 6))$ $\text{thick} = (A_{epi} / \pi)^{(1/2)} - (A_{end} / \pi)^{(1/2)}$ References : Nelson B. Schiller, MD, et al. Recommendations for Quantitation of the Left Ventricle by Two-Dimensional Echocardiography. American Society of Echocardiography Committee on Standards, Subcommittee on Quantitation of Two-Dimensional Echocardiograms. Journal of the American Society of Echocardiography Vol.2, No.5 September-October 1989.358-367. Nelson B. Schiller, MD, et al. Two-Dimensional Echocardiographic Determination of Left Ventricular Volume, Systolic Function, and Mass: Summary and Discussion of the 1989 Recommendations of the American Society of Echocardiography. Circulation Vol.84, No.3 1991;84 [Suppl I]: I -280- I -287.
A12232-002	Left Ventricle Mass by Penn	$\text{Mass} = 1.04 * ((IVSd + LVIDd + LVPWd)^3 - LVIDd^3) - 13.6$ Mass unit is grams and length in cm. References : Richard B. Devereux. Detection of Left Ventricular Hypertrophy by M-Mode Echocardiography. Anatomic Validation, Standardization, and Comparison to Other Methods. Hypertension 9 [Suppl II]; II - 19 - 26, 1987. Donald C. Wallerson and Richard B. Devereux. Reproducibility of Echocardiographic Left Ventricular Measurements. Hypertension 9 [Suppl II]; II - 6 - 18, 1987. American Society of Echocardiography Committee on Standards, Subcommittee on Quantitation of Two-Dimensional Echocardiograms. Recommendations for Quantitation of the Left Ventricle by Two-Dimensional Echocardiography. Journal of the American Society of Echocardiography Volume 2 Number 5 September-October 1989.

Code Value	Code Meaning	Definition
A12236-001	Proximal Right Ventricle Outflow Tract	Proximal Right Ventricle Outflow Tract
A12236-002	Distal Right Ventricle Outflow Tract	Distal Right Ventricle Outflow Tract
A12236-003	Septal Mitral Annulus	Septal Mitral Annulus
A12236-004	Septal Tricuspid Annulus	Septal Tricuspid Annulus
A12236-005	Lateral Tricuspid Annulus	Lateral Tricuspid Annulus
A12236-006	Epicardium	Epicardium
A12236-007	Ascending Aorta	Ascending Aorta
A12236-008	Descending Aorta	Descending Aorta
A12236-009	Transverse Aorta	Transverse Aorta
A12236-010	Abdominal Aorta	Abdominal Aorta
A12236-011	Patent Foramen Ovale	Patent Foramen Ovale
A12236-012	Patent Ductus Arteriosus	Patent Ductus Arteriosus
A12236-013	Main Pulmonary Artery	Main Pulmonary Artery
A12236-014	Left Pulmonary Artery	Left Pulmonary Artery
A12236-015	Right Pulmonary Artery	Right Pulmonary Artery
A12236-016	Aortic Valve Annulus	Aortic Valve Annulus
A12240-001	Left Ventricular Endocardial Diastolic Area, psax pap view	Left Ventricular Endocardial Diastolic Area, psax pap view
ALOKA_001-001	OB User Definition Doppler1	User Definition Doppler Measurement for OB Application.
ALOKA_001-002	OB User Definition Doppler2	User Definition Doppler Measurement for OB Application.
ALOKA_001-003	OB User Definition Doppler3	User Definition Doppler Measurement for OB Application.
ALOKA_001-005	Ductus Venosus	Ductus Venosus
ALOKA_002-001	Calculated result	Result Parameter of User Definition Equation Measurement for "User's calculation function".
ALOKA_002-002	Equation parameter	Parameter of User Definition Equation Measurement for "User's calculation function".
ALOKA_003-001	Residual Vessel Lumen Diameter	Residual Diameter
ALOKA_003-002	Stenosis Ratio calculated by Vessel and Residual Diameter	$\%STENO = 100 * (A-B) / A$ A : Vessel Lumen diameter B : Residual diameter
ALOKA_003-003	Residual Vessel Lumen Cross-Sectional Area	Residual Area

Code Value	Code Meaning	Definition
ALOKA_003-004	Stenosis Ratio calculated by Vessel and Residual Area	$\%STENO = 100 * (A-B) / A$ A : Vessel Lumen Area B : Residual Area
ALOKA_003-005	Diameter1	Diameter
ALOKA_003-006	Diameter2	Diameter
ALOKA_003-007	Diameter3	Diameter
ALOKA_003-008	Length1	Length
ALOKA_003-009	Length2	Length
ALOKA_004-001	Peak GLS	Peak of Global Longitudinal Strain $GLS(\%) = (MLs - MLd) / MLd;$
		References : Lang, Roberto M., et al. "Guidelines and Standards, Recommendations for cardiac chamber quantification by echocardiography in adults: an update from the American Society of Echocardiography and the European Association of Cardiovascular Imaging." European Heart Journal-Cardiovascular Imaging 16.3 (2015): 233-271. 2.3. Global Longitudinal Strain (GLS). Lagrangian strain is defined as the change in length of an object within a certain direction relative to its baseline length: $Strain(\%) = (Lt - L0) / L0;$ where Lt is the length at time t, and L0 is the initial length at time 0. The most commonly used strain-based measure of LV global systolic function is GLS. It is usually assessed by speckle-tracking echocardiography (STE) 22-24 (Table 1). On 2DE, peak GLS describes the relative length change of the LV myocardium between end-diastole and end-systole: $GLS(\%) = (MLs - MLd) / MLd;$ where ML is myocardial length at end-systole (MLs) and end-diastole (MLd). Because MLs is smaller than MLd, peak GLS is a negative number. This negative nature of GLS can lead to confusion when describing increases or decreases in strain.
ALOKA_004-002	Peak GLS of Overlay	same as above

Code Value	Code Meaning	Definition
ALOKA_004-003	Averaged Peak GLS	Averaged Peak Global Longitudinal Strain For example: $\{ \text{Peak GLS (2ch)} + \text{Peak GLS (4ch)} \} / 2,$ $\{ \text{Peak GLS (2ch)} + \text{Peak GLS (3ch)} + \text{Peak GLS (4ch)} \} / 3$ References : Lang, Roberto M., et al. "Guidelines and Standards, Recommendations for cardiac chamber quantification by echocardiography in adults: an update from the American Society of Echocardiography and the European Association of Cardiovascular Imaging." European Heart Journal-Cardiovascular Imaging 16.3 (2015): 233-271. "GLS measurements should be made in the three standard apical views and averaged. (P7)
ALOKA_004-004	GLS-Sel.	Selected value of Global Longitudinal Strain $\text{Global Longitudinal Strain} = \{ L(t) - L(\text{R-wave}) \} / L(\text{R-wave}) \times 100$ References : Definitions for a common standard for 2D speckle tracking echocardiography: consensus document of the EACVI/ASE/Industry Task Force to standardize deformation imaging, European Heart Journal Cardiovascular Imaging, 2015, 16, 1-11
ALOKA_004-005	GLS-Sel. of Overlay	Selected value of Global Longitudinal Strain by using some images. For example $[\text{GLS}_{4\text{ch}}(t) + \text{GLS}_{3\text{ch}}(t) + \text{GLS}_{2\text{ch}}(t)] / 3$
ALOKA_004-006	Time from R to Peak GLS phase	Time from R to Peak GLS phase.
ALOKA_004-007	Time from R to Peak GLS phase of Overlay	Time from R to Peak GLS phase by using some images.
ALOKA_004-008	Time from R to GLS-Sel. phase	Time from R to selected GLS phase .
ALOKA_004-009	Time from R to GLS-Sel. phase of Overlay	Time from R-wave to selected time at Global Longitudinal Strain by using some images.
ALOKA_005-001	Vessel Diameter Stenosis Ratio	%Steno Diam
ALOKA_005-002	Vessel Area Stenosis Ratio	%Steno Area

Code Value	Code Meaning	Definition
ALOKA_005-003	Gall Bladder Wall	Gall Bladder Wall Thickness
ALOKA_005-004	Common Bile Duct	CBD
ALOKA_005-005	Renal Volume	Lt/Rt.Renal Volume
ALOKA_005-006	Space Occupying Lesion1	SOL
ALOKA_005-007	Space Occupying Lesion2	SOL
ALOKA_005-008	Aorta Diameter	Aorta Diam.
ALOKA_005-009	Portal Vein Diameter	PV Diam.
ALOKA_006-001	Common Bile Duct	Common Bile Duct
ALOKA_006-002	Pancreas Head	Pancreas Head
ALOKA_006-003	Pancreas Body	Pancreas Body
ALOKA_006-004	Pancreas Tail	Pancreas Tail
ALOKA_006-005	Pancreas Duct	Pancreas Duct
ALOKA_006-006	Spleen Index	SI = Long axis × Thickness

**Private Code Definitions (Advanced Cardiac,
Coding Scheme Designator: "99IMPROOF", Coding Scheme Version: not specified)**

Code Value	Code Meaning	Definition
RWT	Relative Wall Thickness	Left Ventricle - Relative Wall Thickness
LVEI	LV Eccentricity Index	Left Ventricle - Eccentricity Index - End Diastole
LVEI	LV Eccentricity Index	Left Ventricle - Eccentricity Index - End Systole
SPWMD	Septal-to-Posterior-Wall Motion Delay	Left Ventricle - Septal-to-Posterior-Wall Motion Delay
CA_LV_VCF	Velocity of Circumferential Fiber Shortening	Left Ventricle - Velocity of Circumferential Fiber Shortening - Teichholz
RWT	Relative Wall Thickness	Left Ventricle - Relative Wall Thickness
CA_LV_DIF	Major axis length difference	Left Ventricle - Major axis length diastolic difference - End Diastole
CA_LV_DIF	Major axis length difference	Left Ventricle - Major axis length systolic difference - End Systole
THI_Myo	Myocardial Thickness	Left Ventricle - Myocardial thickness - Left Ventricle Mass by Truncated Ellipse

Code Value	Code Meaning	Definition
RADsax	Short axis radius calculated from short axis cavity area	Left Ventricle - Short axis radius calculated from short axis cavity area - Left Ventricle Mass by Truncated Ellipse
THI_Myo	Myocardial Thickness	Left Ventricle - Myocardial thickness
AOs2LAs	Aortic Root to Left Atrium Ratio	Aortic Root Diameter to Left Atrium Diameter Ratio - Systole
LAs2AOd	Left Atrium Systole to Aortic Root Diastole Diameter Ratio	Left Atrium Systole to Aortic Root Diastole Diameter Ratio
CA_DIAM_SINO	Aortic Sinotubular Junction Diameter	Aorta - Diameter - Sinotubular junction
CA_DIAM_SINUS	Aortic Sinuses of Valsalva Diameter	Aorta - Diameter - Sinus Valsava
CA_DIAM_ANN	Aortic Annulus Diameter	Aorta - Diameter - Annulus
AOs2LAs	Aortic Root to Left Atrium Ratio	Left Atrium - Aortic Root Diameter to Left Atrium Diameter Ratio
LAs2AOd	Left Atrium Systole to Aortic Root Diastole Diameter Ratio	Left Atrium Systole to Aortic Root Diastole Diameter Ratio
CA_MV_A2E	Mitral Valve A to E Ratio	Mitral Valve - A to E Ratio - Antegrade Flow
CA_MV_DUR	Mitral Valve E-Wave Duration	Mitral Valve - E-Wave Duration - Antegrade Flow
CA_MV_Dia_CC	Mitral Valve Commissural Diameter	Mitral Annulus - Diameter - Systole - commissural
CA_MV_Dia_AP	Mitral Valve Anterior-Posterior Annulus Diameter	Mitral Annulus - Diameter - Systole - anterior-posterior
CA_MV_Ea2Aa	Mitral Valve Ea to Aa Ratio	Left Ventricle Doppler Tissue Imaging - Mitral Valve Ea to Aa Ratio - Medial Mitral Annulus
CA_MV_EWaveAmpl	Amplitude E- Wave	Mitral Valve - Amplitude E- Wave
CA_MV_AWaveAmpl	Amplitude A-Wave	Mitral Valve - Amplitude A-Wave
CA_MV_DEWaveAmpl	Amplitude D-E Wave	Mitral Valve - Amplitude D-E Wave
MV_AMPL_A2E_AMPL_E	A-Wave Amplitude to E-Wave Amplitude	Mitral Valve - A-Wave Amplitude to E-Wave Amplitude
MV_AMPL_E2_AMPL_A	E-Wave Amplitude to A-Wave Amplitude	Mitral Valve - E-Wave Amplitude to A-Wave Amplitude
MAPSE	Mitral Annular Plane Systolic Excursion	Left Ventricle - Mitral Annular Plane Systolic Excursion
CA_DIAM_VC	Vena Contracta Diameter	Mitral Valve - Vena Contracta Diameter
HEIGHT_JET	Jet Height	Mitral Valve - Jet Height

Code Value	Code Meaning	Definition
CA_REG_VOL	Regurgitation Volume	Mitral Valve - Regurgitation volume - Regurgitant Flow - Doppler Mode
CA_AV_P	Aortic Valve Pressure	Aortic Valve - Pressure - End Diastole - Regurgitant Flow - Doppler Mode
dPEP	Interventricular Delay	Interventricular Dyssynchrony - Interventricular Delay - Doppler Mode
LVOTTOAoV_V P	Maximum velocity in the LVOT to Maximum velocity in AV	Aortic Valve - Maximum velocity in the LVOT to Maximum velocity in AV - Antegrade Flow - Doppler Mode
LVOTTOAoV_V TI	VTI in the LVOT to VTI in AV	Aortic Valve - VTI in the LVOT to VTI in AV - Antegrade Flow - Doppler Mode
CA_DIAM_VC	Vena Contracta Diameter	Aortic Valve - Vena Contracta Diameter
HEIGHT_JET	Jet Height	Aortic Valve - Jet Height
CA_AR_WIDT H_JET	Jet Width	Aortic Valve - Jet Width
CA_AR_WJET 2LVOTD	Jet Width to LVOT Diam	Aortic Valve - Jet Width to LVOT diam
CA_TV_A2E	Tricuspid Valve A to E Ratio	Tricuspid Valve - A to E Ratio - Antegrade Flow
TV_Dia_CC	Tricuspid Valve Commissural Diameter	Tricuspid Annulus - Diameter - Systole - commissural
TV_Dia_AP	Tricuspid Valve Anterior-Posterior Annulus Diamater	Tricuspid Annulus - Diameter - Systole - anterior-posterior
RV_Ea	Right Ventricular Peak Early Diastolic Tissue Velocity	Right Ventricle Doppler Tissue Imaging - Right Ventricular Peak Early Diastolic Tissue Velocity - Medial Tricuspid Annulus
RV_Aa	RV Peak Diastolic Tissue Velocity During Atrial Systole	Right Ventricle Doppler Tissue Imaging - RV Peak Diastolic Tissue Velocity During Atrial Systole - Medial Tricuspid Annulus
RV_Sa	Right Ventricular Peak Systolic Tissue Velocity	Right Ventricle Doppler Tissue Imaging - Right Ventricular Peak Systolic Tissue Velocity - Medial Mitral Annulus
RV_Ea	Right Ventricular Peak Early Diastolic Tissue Velocity	Right Ventricle Doppler Tissue Imaging - Right Ventricular Peak Early Diastolic Tissue Velocity - Lateral Tricuspid Annulus
RV_Aa	RV Peak Diastolic Tissue Velocity During Atrial Systole	Right Ventricle Doppler Tissue Imaging - RV Peak Diastolic Tissue Velocity During Atrial Systole - Lateral Tricuspid Annulus
RV_Sa	Right Ventricular Peak Systolic Tissue Velocity	Right Ventricle Doppler Tissue Imaging - Right Ventricular Peak Systolic Tissue Velocity - Lateral Mitral Annulus
CA_DIAM_VC	Vena Contracta Diameter	Tricuspid Valve - Vena Contracta Diameter
HEIGHT_JET	Jet Height	Tricuspid Valve - Jet Height
dPEP	Interventricular Delay	Interventricular Dyssynchrony - Interventricular Delay - Doppler Mode

Code Value	Code Meaning	Definition
CA_DIAM_VC	Vena Contracta Diameter	Pulmonary Valve - Vena Contracta Diameter - 2D mode
HEIGHT_JET	Jet Height	Pulmonary Valve - Jet Height
MPAP	Mean Pulmonary Artery Pressure	Pulmonary Artery - Mean Pulmonary Artery Pressure
CA_VTI_AWave_PVE	Pulmonary Vein A-Wave Velocity Time Integral	Pulmonary Veins - A-Wave Velocity Time Integral - Antegrade Flow
CA_PVE_SF	Systolic Fraction	Pulmonary Veins - Systolic Fraction
RVWT	Right Ventricular Wall Thickness	Right Ventricle - Wall Thickness
CA_RV_L_BaseApex	Right Ventricular base-apex Length	Right Ventricle - Right Ventricular base-apex Length - End Diastole
RVWT	Right Ventricular Wall Thickness	Right Ventricle - Wall Thickness
TAPSE	Tricuspid Annular Plane Systolic Excursion	Right Ventricle - Tricuspid Annular Plane Systolic Excursion
AREA_DIA_RV	Right Ventricular Diastolic Area	Right Ventricle - Area - End Diastole
AREA_SYS_RV	Right Ventricular Systolic Area	Right Ventricle - Area - End Systole
RVFAC	Right Ventricular Fractional Area Change	Right Ventricle - Fractional Area Change
DIAM_VC_Sup	Superior Vena Cava Diameter	Vena Cava - Diameter - Superior
DIAM_VC_Sup	Superior Vena Cava Diameter	Vena Cava - Diameter - During Inspiration - Superior
DIAM_VC_Sup	Superior Vena Cava Diameter	Vena Cava - Diameter - During Expiration - Superior
DIAM_SUPINF	Left Atrium Superior-inferior Systolic Dimension	Left Atrium - Diameter - Systole - Superior-Inferior
DIAM_MIDLAT	Left Atrium Medial-lateral Systolic Dimension	Left Atrium - Diameter - Systole - Medial-Lateral
CA_RA_MAJAXES	Right Atrium Systolic Major Axis	Right Atrium - Systolic major axis - 2D mode - Apical four chamber
CA_RA_MINAXES	Right Atrium Systolic Minor Axis	Right Atrium - Systolic minor axis - 2D mode - Apical four chamber
CA_RA_VOLes	Right Atrial End Systolic Volume	Right Atrium - End Systolic Volume - Area-Length, Single Plane
CA_RA_VOLes	Right Atrial End Systolic Volume	Right Atrium - End Systolic Volume - Area-Length, Single Plane - Body Surface Area
CA_RA_VOLes	Right Atrial End Systolic Volume	Right Atrium - End Systolic Volume - 2D mode - Method of Disks, Single Plane - Apical four chamber
CA_RA_VOLes	Right Atrial End Systolic Volume	Right Atrium - End Systolic Volume - 2D mode - Method of Disks, Single Plane - Apical four chamber - Body Surface Area

Code Value	Code Meaning	Definition
CA_MV_Ea2Aa	Mitral Valve Ea to Aa Ratio	Left Ventricle Doppler Tissue Imaging - Mitral Valve Ea to Aa Ratio - Medial Mitral Annulus
CA_MV_ATa	Mitral Valve Acceleration time of Ea	Left Ventricle Doppler Tissue Imaging - Mitral Valve Acceleration time of Ea - Medial Mitral Annulus
CA_MV_ARa	Mitral Valve Acceleration rate of Ea	Left Ventricle Doppler Tissue Imaging - Mitral Valve Acceleration rate of Ea - Medial Mitral Annulus
CA_MV_DTa	Mitral Valve Deceleration time of Ea	Left Ventricle Doppler Tissue Imaging - Mitral Valve Deceleration time of Ea - Medial Mitral Annulus
CA_MV_DRa	Mitral Valve Deceleration rate of Ea	Left Ventricle Doppler Tissue Imaging - Mitral Valve Deceleration rate of Ea - Medial Mitral Annulus
CA_MV_Ea2Aa	Mitral Valve Ea to Aa Ratio	Left Ventricle Doppler Tissue Imaging - Mitral Valve Ea to Aa Ratio - Lateral Mitral Annulus
CA_MV_ATa	Mitral Valve Acceleration time of Ea	Left Ventricle Doppler Tissue Imaging - Mitral Valve Acceleration time of Ea - Lateral Mitral Annulus
CA_MV_ARa	Mitral Valve Acceleration rate of Ea	Left Ventricle Doppler Tissue Imaging - Mitral Valve Acceleration rate of Ea - Lateral Mitral Annulus
CA_MV_DTa	Mitral Valve Deceleration time of Ea	Left Ventricle Doppler Tissue Imaging - Mitral Valve Deceleration time of Ea - Lateral Mitral Annulus
CA_MV_DRa	Mitral Valve Deceleration rate of Ea	Left Ventricle Doppler Tissue Imaging - Mitral Valve Deceleration rate of Ea - Lateral Mitral Annulus
RV_Ea	Right Ventricular Peak Early Diastolic Tissue Velocity	Right Ventricle Doppler Tissue Imaging - Right Ventricular Peak Early Diastolic Tissue Velocity - Medial Tricuspid Annulus
RV_Aa	RV Peak Diastolic Tissue Velocity During Atrial Systole	Right Ventricle Doppler Tissue Imaging - RV Peak Diastolic Tissue Velocity During Atrial Systole - Medial Tricuspid Annulus
RV_Sa	Right Ventricular Peak Systolic Tissue Velocity	Right Ventricle Doppler Tissue Imaging - Right Ventricular Peak Systolic Tissue Velocity - Medial Mitral Annulus
RV_Ea	Right Ventricular Peak Early Diastolic Tissue Velocity	Right Ventricle Doppler Tissue Imaging - Right Ventricular Peak Early Diastolic Tissue Velocity - Lateral Tricuspid Annulus
RV_Aa	RV Peak Diastolic Tissue Velocity During Atrial Systole	Right Ventricle Doppler Tissue Imaging - RV Peak Diastolic Tissue Velocity During Atrial Systole - Lateral Tricuspid Annulus
RV_Sa	Right Ventricular Peak Systolic Tissue Velocity	Right Ventricle Doppler Tissue Imaging - Right Ventricular Peak Systolic Tissue Velocity - Lateral Mitral Annulus
CA_CS_Pf-SF	Pulmonary minus Systemic Shunt Flow	Cardiac Shunt Study - Pulmonary minus Systemic Shunt Flow - Doppler Volume Flow

8.6.5 Standard Extended and Private Template Definitions (Non-Advanced Cardiac)

This section defines the Standard Extended template and Private templates.

Extended Template Definitions TID 5000 OB-GYN Ultrasound Procedure Report

	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value set Constraint
1			CONTAINER	EV(12500,DCM, "OB-GYN Ultrasound Procedure Report")	1	M		
2	>	HAS CONCEPT MOD	INCLUDE	DTID(1204)	1	U		
3	>	HAS OBS CONTEXT	INCLUDE	DTID(1001)	1	M		
4	>	CONTAINS	INCLUDE	DTID(5001)	1	U		
5	>	CONTAINS	CONTAINER	DT(111028, DCM, "Image Library")	1	U		
6	>>	CONTAINS	IMAGE	No purpose of reference	1-n	M		
7	>	CONTAINS	INCLUDE	DTID(5002)	1	U		
8	>	CONTAINS	INCLUDE	DTID(5004)	1-n	U		
9	>	CONTAINS	INCLUDE	DTID(5005)	1-n	U		
10	>	CONTAINS	INCLUDE	DTID(5006)	1-n	U		
11	>	CONTAINS	INCLUDE	DTID(5007)	1-n	U		
12	>	CONTAINS	INCLUDE	DTID(5009)	1-n	U		
13	>	CONTAINS	INCLUDE	DTID(5011)	1-n	U		
14	>	CONTAINS	INCLUDE	DTID(5010)	1	U		
15	>	CONTAINS	INCLUDE	DTID(5015)	1	U		
16	>	CONTAINS	INCLUDE	DTID(5012)	1	U		
17	>	CONTAINS	INCLUDE	DTID(5013)	1	U		\$Laterality = EV (G-A101,SRT,"Left") \$Number = EV (11879-4, LN, "Number of follicles in left ovary")

	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value set Constraint
18	>	CONTAINS	INCLUDE	DTID(5013)	1	U		\$Laterality = EV (G-A100, SRT, "Right") \$Number = EV (11880-2, LN, "Number of follicles in right ovary")
19	>	CONTAINS	CONTAINER	EV(121070, DCM, "Finding Site")	1-n	U		
20	>>	HAS CONCEPT MOD	CODE	EV(G-C0E3, SRT, "Finding Site")	1	M		EV(T-F6800, SRT, "Embryonic Vascular Structure")
21	>>	CONTAINS	INCLUDE	DTID(5025)	1	M		\$AnatomyGroup= DCID(12141)
22	>	CONTAINS	CONTAINER	EV(121070, DCM, "Findings")	1-n	U		
23	>>	HAS CONCEPT MOD	CODE	EV(G-C0E3, SRT, "Finding Site")	1	M		EV(T-D6007, SRT, "Pelvic Vascular Structure")
24	>>	CONTAINS	INCLUDE	DTID(5026)	1	M		\$AnatomyGroup= DCID(12140)
25	>	CONTAINS	INCLUDE	DTID(ALOKA-5001)	1	U		\$Anatomy = DCID(12140) \$Meas = DCID(12119)
26	>	CONTAINS	INCLUDE	DTID(ALOKA-5000)	1-n	U		\$Anatomy = DCID(12140) \$Meas = DCID(12119)
27	>	CONTAINS	INCLUDE	DTID(ALOKA-5000)	1-n	U		\$Anatomy = DCID(ALOKA_001) \$Meas = DCID(12119)
28	>	CONTAINS	INCLUDE	DTID(ALOKA-5002) OB-GYN User Defined Equation Section	1-n	U		
29	>	CONTAINS	INCLUDE	DTID(ALOKA-301) User Defined Equation Section	1-n	U		

Extended Template Definitions
TID 5100 Vascular Ultrasound Procedure Report

	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value set Constraint
1			CONTAINER	EV(12510,DCM, "Vascular Ultrasound Procedure Report")	1	M		
2	>	HAS OBS CONTEXT	CODE	EV(R-40FB8, SRT, "Temporal periods Relating to Procedure")	1	U		DCID(12102) Temporal Periods Relating To Procedure or Therapy
3	>	HAS CONCEPT MOD	INCLUDE	DTID(1204) Language of Content Item and Descendants	1	U		
4	>	HAS OBS CONTEXT	INCLUDE	DTID(1001)Observation Context	1	M		
5	>	CONTAINS	INCLUDE	DTID(5101) Vascular Patient Characteristics	1	U		
6	>	CONTAINS		EV(111028, DCM, "Image Library")	1	U		
7	>>	CONTAINS		No purpose of reference	1-n	M		
8	>	CONTAINS	INCLUDE	DTID(5102) Vascular Procedure Summary Section	1	U		
9	>	CONTAINS	INCLUDE	DTID(5103) Vascular Ultrasound Section	1	U		\$SectionScope = DT (T-40501, SRT, "Blood Vessel of Head") \$SectionLaterality = EV (G-A101, SRT, "Left") \$Anatomy = DCID (12105) Intracranial Cerebral Vessels

	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value set Constraint
10	>	CONTAINS	INCLUDE	DTID(5103) Vascular Ultrasound Section	1	U		\$SectionScope = DT (T-40501, SRT, "Blood Vessel of Head") \$SectionLaterality = EV (G-A100, SRT, "Right") \$Anatomy = DCID (12105) Intracranial Cerebral Vessels
11	>	CONTAINS	INCLUDE	DTID(5103) Vascular Ultrasound Section	1	U		\$SectionScope = DT (T-40501, SRT, "Blood Vessel of Head") \$SectionLaterality = EV (G-A103, SRT, "Unilateral") \$Anatomy = DCID (12106) Intracranial Cerebral Vessels (unilateral)
12	>	CONTAINS	INCLUDE	DTID(5103) Vascular Ultrasound Section	1	U		\$SectionScope = DT (T-45005, SRT, "Artery of neck") \$SectionLaterality = EV (G-A101, SRT, "Left") \$Anatomy = DCID (12104) Extracranial Arteries \$AnatomyRatio = DCID (12123) Carotid Ratios

	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value set Constraint
13	>	CONTAINS	INCLUDE	DTID(5103) Vascular Ultrasound Section	1	U		\$SectionScope = DT (T-45005, SRT, "Artery of neck") \$SectionLaterality = EV (G-A100, SRT, "Right") \$Anatomy = DCID (12104) xtracranial Arteries \$AnatomyRatio = DCID(12123) Carotid Ratios
14	>	CONTAINS	INCLUDE	DTID(5103) Vascular Ultrasound Section	1	U		\$SectionScope = DT (T-47040, SRT, "Artery of Lower Extremity") \$SectionLaterality = EV (G-A101, SRT, "Left") \$Anatomy = DCID (12109) Lower Extremity Arteries
15	>	CONTAINS	INCLUDE	DTID(5103) Vascular Ultrasound Section	1	U		\$SectionScope = DT (T-47040, SRT, "Artery of Lower Extremity") \$SectionLaterality = EV (G-A100, SRT, "Right") \$Anatomy = DCID (12109) Lower Extremity Arteries
16	>	CONTAINS	INCLUDE	DTID(5103) Vascular Ultrasound Section	1	U		\$SectionScope = DT (T-49403, SRT, "Vein of Lower Extremity") \$SectionLaterality = EV (G-A101, SRT, "Left") \$Anatomy = DCID (12110) Lower Extremity Veins

	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value set Constraint
17	>	CONTAINS	INCLUDE	DTID(5103) Vascular Ultrasound Section	1	U		\$SectionScope = DT (T-49403, SRT, "Vein of Lower Extremity") \$SectionLaterality = EV (G-A100, SRT, "Right") \$Anatomy = DCID (12110) Lower Extremity Veins
18	>	CONTAINS	INCLUDE	DTID(5103) Vascular Ultrasound Section	1	U		\$SectionScope = DT (T-47020, SRT, "Artery Of Upper Extremity") \$SectionLaterality = EV (G-A101, SRT, "Left") \$Anatomy = DCID (12107) Upper Extremity Arteries
19	>	CONTAINS	INCLUDE	DTID(5103) Vascular Ultrasound Section	1	U		\$SectionScope = DT (T-47020, SRT, "Artery Of Upper Extremity") \$SectionLaterality = EV (G-A100, SRT, "Right") \$Anatomy = DCID (12107) Upper Extremity Arteries
20	>	CONTAINS	INCLUDE	DTID(5103) Vascular Ultrasound Section	1	U		\$SectionScope = DT (T-49103, SRT, "Vein Of Upper Extremity") \$SectionLaterality = EV (G-A101, SRT, "Left") \$Anatomy = DCID (12108) Upper Extremity Veins

	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value set Constraint
21	>	CONTAINS	INCLUDE	DTID(5103) Vascular Ultrasound Section	1	U		\$SectionScope = DT (T-49103, SRT, "Vein Of Upper Extremity") \$SectionLaterality = EV (G-A100, SRT, "Right") \$Anatomy = DCID (12108) Upper Extremity Veins
22	>	CONTAINS	INCLUDE	DTID(5103) Vascular Ultrasound Section	1	U		\$SectionScope = DT (T-71019, SRT, "Vascular Structure Of Kidney") \$SectionLaterality = EV (G-A101, SRT, "Left") \$Anatomy = DCID (12115) Renal Vessels \$AnatomyRatio = DCID (12124) Renal Ratios
23	>	CONTAINS	INCLUDE	DTID(5103) Vascular Ultrasound Section	1	U		\$SectionScope = DT (T-71019, SRT, "Vascular Structure Of Kidney") \$SectionLaterality = EV (G-A100, SRT, "Right") \$Anatomy = DCID (12115) Renal Vessels \$AnatomyRatio = DCID (12124) Renal Ratios

	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value set Constraint
24	>	CONTAINS	INCLUDE	DTID(5103) Vascular Ultrasound Section	1	U		\$SectionScope = DT (T-46002, SRT, "Artery of Abdomen") \$SectionLaterality = EV (G-A101, SRT, "Left") \$Anatomy = DCID (12111) Abdominal Arteries (lateral)
25	>	CONTAINS	INCLUDE	DTID(5103) Vascular Ultrasound Section	1	U		\$SectionScope = DT (T-46002, SRT, "Artery of Abdomen") \$SectionLaterality = EV (G-A100, SRT, "Right") \$Anatomy = DCID (12111) Abdominal Arteries (lateral)
26	>	CONTAINS	INCLUDE	DTID(5103) Vascular Ultrasound Section	1	U		\$SectionScope = DT (T-46002, SRT, "Artery of Abdomen") \$SectionLaterality = EV (G-A103, SRT, "Unilateral") \$Anatomy = DCID (12112) Abdominal Arteries (unilateral)
27	>	CONTAINS	INCLUDE	DTID(5103) Vascular Ultrasound Section	1	U		\$SectionScope = DT (T-487A0, SRT, "Vein of Abdomen") \$SectionLaterality = EV (G-A101, SRT, "Left") \$Anatomy = DCID (12113) Abdominal Veins(lateral)

	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value set Constraint
28	>	CONTAINS	INCLUDE	DTID(5103) Vascular Ultrasound Section	1	U		\$SectionScope = DT (T-487A0, SRT, "Vein of Abdomen") \$SectionLaterality = EV (G-A100, SRT, "Right") \$Anatomy = DCID (12113) Abdominal Veins(lateral)
29	>	CONTAINS	INCLUDE	DTID(5103) Vascular Ultrasound Section	1	U		\$SectionScope = DT (T-487A0, SRT, " Vein of Abdomen") \$SectionLaterality = EV (G-A103, SRT, "Unilateral") \$Anatomy = DCID (12114) Abdominal Veins(unilateral)
30	>	CONTAINS	INCLUDE	DTID(5105) Ultrasound Graft Section	1	U		
31	>	CONTAINS	INCLUDE	DTID(ALOKA-301) User Defined Equation Section	1-n	U		
32	>	CONTAINS	INCLUDE	DTID(ALOKA-5100) Vascular Original Measurement section	1	U		\$SectionScope= DT(A-011, 99ALOKA, "Original Measurement Group of Abdomen") \$Anatomy = DCID(ALOKA_005) Abdomen Original Measurements Group \$Meas = DCID(ALOKA_006) Abdomen Original Measurement Items

Extended Template Definitions
TID 5200 Echocardiography Ultrasound Procedure Report

	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value set Constraint
1			CONTAINER	EV (125200, DCM, "Adult Echocardiography Procedure Report")	1	M		
2	>	HAS CONCEPT MOD	INCLUDE	DTID (1204) Language of Content Item and Descendants	1	U		
3	>	HAS OBS CONTEXT	INCLUDE	DTID (1001) Observation Context	1	M		
4	>	CONTAINS	CONTAINER	DT (121064, DCM, "Current Procedure Descriptions")	1	U		
5	>>	CONTAINS	CODE	DT (125203, DCM, "Acquisition Protocol")	1-n	M		BCID (12001) Ultrasound Protocol Types
6	>	CONTAINS	INCLUDE	DTID (5201) Echocardiography Patient Characteristics	1	U		
7	>	CONTAINS	CONTAINER	(111028, DCM, "Image Library")	1	U		
8	>>	CONTAINS	IMAGE	No purpose of reference	1-n	M		
9	>	CONTAINS	INCLUDE	DTID(5202) Echo Section	1	U		\$SectionSubject = EV (T-32600, SRT, "Left Ventricle") \$MeasType = DCID (12200) Echocardiography Left Ventricle
10	>	CONTAINS	INCLUDE	DTID(5202) Echo Section	1	U		\$SectionSubject = EV (T-32500, SRT, "Right Ventricle") \$MeasType = DCID (12204) Echocardiography Right Ventricle

	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value set Constraint
11	>	CONTAINS	INCLUDE	DTID(5202) Echo Section	1	U		\$SectionSubject = EV (T-32300, SRT, "Left Atrium") \$MeasType = DCID (12205) Echocardiography Left Atrium
12	>	CONTAINS	INCLUDE	DTID(5202) Echo Section	1	U		\$SectionSubject = EV (T-32200, SRT, "Right Atrium") \$MeasType = DCID (12206) Echocardiography Right Atrium
13	>	CONTAINS	INCLUDE	DTID(5202) Echo Section	1	U		\$SectionSubject = EV (T-35400, SRT, "Aortic Valve") \$MeasType = DCID (12211) Echocardiography Aortic Valve
14	>	CONTAINS	INCLUDE	DTID(5202) Echo Section	1	U		\$SectionSubject = EV (T-35300, SRT, "Mitral Valve") \$MeasType = DCID (12207) Echocardiography Mitral Valve
15	>	CONTAINS	INCLUDE	DTID(5202) Echo Section	1	U		\$SectionSubject = EV (T-35200, SRT, "Pulmonic Valve") \$MeasType = DCID (12209) Echocardiography Pulmonic Valve
16	>	CONTAINS	INCLUDE	DTID(5202) Echo Section	1	U		\$SectionSubject = EV (T-35100, SRT, "Tricuspid Valve") \$MeasType = DCID (12208) Echocardiography Tricuspid Valve

	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value set Constraint
17	>	CONTAINS	INCLUDE	DTID(5202) Echo Section	1	U		\$SectionSubject = EV (T-42000, SRT, "Aorta") \$MeasType= DCID (12212) Echocardiography Aorta
18	>	CONTAINS	INCLUDE	DTID(5202) Echo Section	1	U		\$SectionSubject = EV (T-44000, SRT, "Pulmonary artery") \$MeasType DCID (12210) = Echocardiography Pulmonary Artery
19	>	CONTAINS	INCLUDE	DTID(5202) Echo Section	1	U		\$SectionSubject = EV (T-48600, SRT, "Vena Cava" \$MeasType = DCID (12215) Echocardiography Vena Cavae
20	>	CONTAINS	INCLUDE	DTID(5202) Echo Section	1	U		\$SectionSubject = EV (T-48581, SRT, "Pulmonary Venous Structure") \$MeasType = DCID (12214) Echocardiography Pulmonary Veins
21	>	CONTAINS	INCLUDE	DTID(5202) Echo Section	1	U		\$SectionSubject = EV (T-39050,SRT, "Pericardial cavity") \$MeasType = DCID 12250 "CardiacUltrasound Common LinearMeasuremen ts"
22	>	CONTAINS	INCLUDE	DTID(5202) Echo Section	1	U		\$SectionSubject = EV (P5-30031, SRT, "Cardiac Shunt Study") \$MeasType = DCID (12217) Echocardiography Cardiac Shunt

	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value set Constraint
23	>	CONTAINS	INCLUDE	DTID(5202) Echo Section	1	U		\$SectionSubject = EV (D4-30000, SRT, "Congenital Anomaly of Cardiovascular System") \$MeasType = DCID (12218) Echocardiography Congenital
24	>	CONTAINS	INCLUDE	DTID(5204) Wall Motion Analysis	1-n	U		\$Procedure = DT (P5-B3121, SRT, "Echocardiography for Determining Ventricular Contraction")
25	>	CONTAINS	INCLUDE	DTID(ALOKA- 301) User Defined Equation Section	1-n	U		

Private Template Definitions
TID ALOKA-300 Common Original Measurement Section

	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value set Constraint
1			CONTAINER	\$AnatomyGroup	1	M		
2	>	HAS CONCEPT MOD	CODE	EV(G-C171, SRT,"Laterality")	1	U		\$SectionLaterality
3	>	CONTAINS	INCLUDE	DTID(300)	1-n	M		\$Measurement = \$ MeasType \$Derivation = DCID(3627)

Private Template Definitions
TID ALOKA-301 User Defined Equation Section

	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value set Constraint
1			CONTAINER	EV (A-004, 99ALOKA, "User Defined Equation")	1	M		
2	>	CONTAINS	TEXT	EV (A-005, 99ALOKA, "User Defined Equation Group Number")	1	U		
3	>	CONTAINS	TEXT	EV (A-006, 99ALOKA, "User Defined Equation Group Name")	1	M		
4	>>	HAS PROPERTIES	INCLUDE	DTID(ALOKA-302) User Defined Equation Items Section	1-n	M		

Private Template Definitions
TID ALOKA-302 User Defined Equation Items Section

	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value set Constraint
1			INCLUDE	DTID(300) Measurement	1	M		\$Measurement = DCID(ALOKA_002) User Defined Equation Group \$Derivation = DCID(3627) Measurement Type
2	>	INFERRED FROM	TEXT	EV(A-007, 99ALOKA, "User Defined Parameter Number")	1	U		
3	>	INFERRED FROM	TEXT	EV (A-008, 99ALOKA, "User Defined Parameter Name")	1	U		

4	>	INFERRED FROM	TEXT	EV (A-009, 99ALOKA, "User Defined Equation")	1	U		
---	---	---------------	------	--	---	---	--	--

**Private Template Definitions
TID ALOKA-5000 OB Original Fetal Measurement Section**

	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value set Constraint
1			CONTAINER	EV(A-002, 99ALOKA, "OB Original Fetal Measurement Group")	1	M		
2	>	HAS OBS CONTEXT	INCLUDE	DTID(1008)	1	MC	If this template is invoked more than once to describe more than one fetus	
3	>	CONTAINS	INCLUDE	DTID(ALOKA-300)	1-n	M		\$AnatomyGroup = \$Anatomy \$MeasType = \$Meas \$SectionLaterality = DCID(244)

**Private Template Definitions
TID ALOKA-5001 OB-GYN Original Measurement Section**

	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value set Constraint
1			CONTAINER	EV(A-003, 99ALOKA, "OB-GYN Original Fetal Measurement Group")	1	M		
2	>	CONTAINS	INCLUDE	DTID(ALOKA-300)	1-n	M		\$AnatomyGroup = \$Anatomy \$MeasType = \$Meas \$SectionLaterality = DCID(244)

Private Template Definitions
TID ALOKA-5002 OB-GYN User Defined Equation Section

	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value set Constraint
1			CONTAINER	EV(A-010, 99ALOKA, "OB-GYN User Defined Measurement Group")	1			
2	>	HAS OBS CONTEXT	INCLUDE	DTID(1008) Subject Context, Fetus	1	MC	IFF this template is invoked to describe fetus	
3	>	CONTAINS	INCLUDE	DTID(ALOKA-301) User Defined Equation Section	1-n	M		

Note: It is configurable to output TID ALOKA-5002 described in prior version of the DICOM Conformance Statement of this product for backward compatibility.

Private Template Definitions
TID ALOKA-5100 Vascular Original Measurement Section

	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value set Constraint
1			CONTAINER	\$SectionScope	1	M		
2	>	CONTAINS	INCLUDE	DTID(ALOKA-300)	1-n	M		\$AnatomyGroup = \$Anatomy \$MeasType = \$Meas \$SectionLaterality = DCID(244)