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SU-1

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- Information for Manufacturers Seeking Marketing Clearance of Diagnostic Ultrasound Systems and Transducers, FDA (1997 年)
- Medical Ultrasound Safety, 米国超音波医学会議 (AIUM) (1994 年)
- Acoustic Output Measurement Standard for Diagnostic Ultrasound Equipment (NEMA UD2-2004)
- Acoustic Output Measurement and Labeling Standard for Diagnostic Ultrasound Equipment, 米国超音波医学会議 (AIUM) (1993 年)
- Standard for Real-Time Display of Thermal and Mechanical Acoustic Output Indices on Diagnostic Ultrasound Equipment (NEMA UD3-2004)
- Guidance on the interpretation of TI and MI to be used to inform the operator, Annex HH, BS EN 60601-2-37 reprinted at P05699.

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Index Label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	
Maximum Index Value		0.72	0.21		0.39		-
Index Component Value			0.21	0.07	0.17	0.27	
Acoustic Parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	1.79					
	P (mW)		3.23	3.03	3.23	1.72	-
	P_{1x1} (mW)		3.23	3.03	3.23	1.72	
	z_s (cm)			0.65			
	z_b (cm)					0.85	
	z_{MI} (cm)	0.73					
	$z_{pii,\alpha}$ (cm)	0.73					
	f_{awf} (MHz)	6.26	6.18	6.20	6.18	6.19	-
Other Information	prr (Hz)	1536					
	srr (Hz)	60					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	75.30					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm ²)	3.17					
	I_{spta} at z_{pii} or z_{sii} (mW/cm ²)	33.28					
	p_r at z_{pii} (MPa)	2.10					
Operating Control Conditions	Control 1: 2D Frequency/THI	10MHz/On	10 MHz/On		10 MHz/On		-
	Control 2: Color Frequency	6 MHz	6 MHz		6 MHz		-
	Control 3: Color Box Size	Narrow	Narrow		Narrow		-
	Control 4: PW Frequency	6 MHz	6 MHz		6 MHz		-
	Control 5: Color PRF	5376 Hz	5376 Hz		5376 Hz		-
	Control 6: PW PRF	6579 Hz	6579 Hz		6579 Hz		-
	Control 7: Focus	10 mm	110 mm		10 mm		-
	Control 8: Acoustic Power	100%	100%		100%		-
	Control 9: Line Density	1	1		-		-

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Index Label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	
Maximum Index Value		0.72	0.19		0.37		-
Index Component Value			0.19	0.07	0.15	0.27	
Acoustic Parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	1.79					
	P (mW)		2.50	3.03	2.50	1.72	-
	$P_{1 \times 1}$ (mW)		2.50	3.03	2.50	1.72	
	z_s (cm)			0.65			
	z_b (cm)					0.85	
	z_{MI} (cm)	0.73					
	$z_{pii,\alpha}$ (cm)	0.73					
	f_{awf} (MHz)	6.26	6.22	6.20	6.22	6.19	-
Other Information	p_{rr} (Hz)	1536					
	s_{rr} (Hz)	60					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	75.30					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm ²)	3.17					
	I_{spta} at z_{pii} or z_{sii} (mW/cm ²)	33.28					
	p_r at z_{pii} (MPa)	2.10					
Operating Control Conditions	Control 1: 2D Frequency/THI	10MHz/On	10 MHz/On		10 MHz/On		-
	Control 2: Color Frequency	6 MHz	6 MHz		6 MHz		-
	Control 3: Color Box Size	Narrow	Narrow		Narrow		-
	Control 4: PW Frequency	6 MHz	6 MHz		6 MHz		-
	Control 5: Color PRF	5376 Hz	5376 Hz		5376 Hz		-
	Control 6: PW PRF	6579 Hz	6579 Hz		6579 Hz		-
	Control 7: Focus	10 mm	100 mm		10 mm		-
	Control 8: Acoustic Power	100%	100%		100%		-
	Control 9: Line Density	1	1		1		-

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Index Label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	
Maximum Index Value		0.74	0.20		0.20		-
Index Component Value			0.20	-	0.20	-	
Acoustic Parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	1.66					
	P (mW)		5.47	-	5.47	-	-
	$P_{1 \times 1}$ (mW)		5.47	-	5.47	-	
	z_s (cm)			-			
	z_b (cm)					-	
	z_{MI} (cm)	2.00					
	$z_{pii,\alpha}$ (cm)	2.00					
	f_{awf} (MHz)	4.98	6.02	-	6.02	-	-
Other Information	p_{rr} (Hz)	992					
	s_{rr} (Hz)	30					
	n_{pps} (Repeat)	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	121.95					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm ²)	3.90					
	I_{spta} at z_{pii} or z_{sii} (mW/cm ²)	49.21					
	p_r at z_{pii} (MPa)	2.34					
Operating Control Conditions	Control 1: 2D Frequency/THI	5 MHz/Off	10 MHz/On		10 MHz/On		10 MHz/On
	Control 2: Color Frequency	4 MHz	6 MHz		6 MHz		4 MHz
	Control 3: Color Box Size	Narrow	Narrow		Narrow		Narrow
	Control 4: Color PRF	11160 Hz	10929 Hz		10929 Hz		10929 Hz
	Control 5: Focus	30 mm	70 mm		70 mm		70 mm
	Control 6: Acoustic Power	100%	100%		100%		100%
	Control 7: Line Density	1	1		1		1

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Index Label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	
Maximum Index Value		0.57	0.17		0.34		-
Index Component Value			0.17	0.05	0.14	0.24	
Acoustic Parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	1.36					
	P (mW)		4.28	2.54	4.28	2.33	-
	P_{1x1} (mW)		4.28	2.54	4.28	2.33	
	z_s (cm)			1.13			
	z_b (cm)					1.38	
	z_{MI} (cm)	1.73					
	$z_{pii,\alpha}$ (cm)	1.73					
	f_{awf} (MHz)	5.71	4.03	6.11	4.03	4.00	-
Other Information	p_{rr} (Hz)	650					
	s_{rr} (Hz)	24					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	82.20					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm ²)	1.51					
	I_{spta} at z_{pii} or z_{sii} (mW/cm ²)	18.98					
	p_r at z_{pii} (MPa)	1.92					
Operating Control Conditions	Control 1: 2D Frequency/THI	10 MHz/On	10 MHz/On		10 MHz/On		-
	Control 2: Color Frequency	4 MHz	4 MHz		4 MHz		-
	Control 3: Color Box Size	Narrow	Narrow		Narrow		-
	Control 4: PW Frequency	4 MHz	6 MHz		4 MHz		-
	Control 5: Color PRF	6586 Hz	6586 Hz		6586 Hz		-
	Control 6: PW PRF	6579 Hz	6579 Hz		6579 Hz		-
	Control 7: Focus	25 mm	90mm		15 mm		-
	Control 8: Acoustic Power	100%	100%		100%		-
	Control 9: Line Density	1	1		1		-

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Index Label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	
Maximum Index Value		0.57	0.17		0.34		-
Index Component Value			0.17	0.05	0.14	0.24	
Acoustic Parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	1.36					
	P (mW)		2.27	2.54	2.27	2.33	-
	$P_{1 \times 1}$ (mW)		2.27	2.54	2.27	2.33	
	z_s (cm)			1.13			
	z_b (cm)					1.38	
	z_{MI} (cm)	1.73					
	$z_{pii,\alpha}$ (cm)	1.73					
	f_{awf} (MHz)	5.71	6.03	6.11	6.03	4.00	-
Other Information	prr (Hz)	650					
	srr (Hz)	24					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	82.20					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm ²)	1.51					
	I_{spta} at z_{pii} or z_{sii} (mW/cm ²)	18.98					
	p_r at z_{pii} (MPa)	1.92					
Operating Control Conditions	Control 1: 2D Frequency/THI	10 MHz/On	10 MHz/On		10 MHz/On		-
	Control 2: Color Frequency	6 MHz	6 MHz		6 MHz		-
	Control 3: Color Box Size	Narrow	Narrow		Narrow		-
	Control 4: PW Frequency	4 MHz	6 MHz		4 MHz		-
	Control 5: Color PRF	6586 Hz	6586 Hz		6586 Hz		-
	Control 6: PW PRF	6579 Hz	6579 Hz		6579 Hz		-
	Control 7: Focus	25 mm	80 mm		80 mm		-
	Control 8: Acoustic Power	100%	100%		100%		-
	Control 9: Line Density	1	1		-		-

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Index Label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	
Maximum Index Value		0.80	0.17		0.41		-
Index Component Value			0.17	0.04	0.13	0.31	
Acoustic Parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	1.80					
	P (mW)		2.19	2.58	2.19	2.10	-
	P_{1x1} (mW)		2.19	2.58	2.19	2.10	
	z_s (cm)			1.30			
	z_b (cm)					1.25	
	z_{MI} (cm)	1.38					
	$z_{pii,\alpha}$ (cm)	1.38					
	f_{awf} (MHz)	5.08	6.05	6.05	6.05	4.09	-
Other Information	p_{rr} (Hz)	778					
	s_{rr} (Hz)	30					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	101.28					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm ²)	2.23					
	I_{spta} at z_{pii} or z_{sii} (mW/cm ²)	26.69					
	p_r at z_{pii} (MPa)	2.29					
Operating Control Conditions	Control 1: 2D Frequency/THI	5 MHz/Off	8 MHz/On		8 MHz/On		-
	Control 2: Color Frequency	4MHz	6 MHz		6 MHz		-
	Control 3: Color Box Size	Narrow	Narrow		Narrow		-
	Control 4: PW Frequency	4 MHz	6 MHz		4 MHz		-
	Control 5: Color PRF	5589Hz	5589 Hz		5589 Hz		-
	Control 6: PW PRF	6579Hz	6579 Hz		6579 Hz		-
	Control 7: Focus	20 mm	110 mm		10 mm		-
	Control 8: Acoustic Power	100%	100%		100%		-
	Control 9: Line Density	1	1		1		-

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Index Label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	
Maximum Index Value		0.80	0.16		0.40		-
Index Component Value			0.16	0.04	0.13	0.31	
Acoustic Parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	1.80					
	P (mW)		2.00	2.58	2.00	2.10	-
	P_{1x1} (mW)		2.00	2.58	2.00	2.10	
	z_s (cm)			1.30			
	z_b (cm)					1.25	
	z_{MI} (cm)	1.38					
	$z_{pii,\alpha}$ (cm)	1.38					
	f_{awf} (MHz)	5.08	6.07	6.05	6.07	4.09	-
Other Information	pr_r (Hz)	778					
	srr (Hz)	30					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	101.28					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm ²)	2.23					
	I_{spta} at z_{pii} or z_{sii} (mW/cm ²)	26.69					
	p_r at z_{pii} (MPa)	2.29					
Operating Control Conditions	Control 1: 2D Frequency/THI	5 MHz/Off	8 MHz/On	-	8 MHz/On	-	-
	Control 2: Color Frequency	6 MHz	6 MHz	-	6 MHz	-	-
	Control 3: Color Box Size	Narrow	Narrow	-	Narrow	-	-
	Control 4: PW Frequency	4 MHz	-	6 MHz	-	4 MHz	-
	Control 5: Color PRF	5589Hz	5589 Hz	-	5589 Hz	-	-
	Control 6: PW PRF	6579Hz	-	6579 Hz	-	6579 Hz	-
	Control 7: Focus	20 mm	120 mm	110 mm	120 mm	10 mm	-
	Control 8: Acoustic Power	100%	100%	100%	100%	100%	-
	Control 9: Line Density	1	1	-	1	-	-

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Index Label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	
Maximum Index Value		0.73	0.23		0.55		-
Index Component Value			0.23	0.07	0.19	0.42	
Acoustic Parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	1.66					
	P (mW)		2.93	3.48	2.93	3.12	-
	P_{1x1} (mW)		2.93	3.48	2.93	3.12	
	z_s (cm)			1.05			
	z_b (cm)					1.40	
	z_{MI} (cm)	1.33					
	$z_{pii,\alpha}$ (cm)	1.33					
	f_{awf} (MHz)	5.14	6.13	6.13	6.13	4.07	-
Other Information	prr (Hz)	819					
	srr (Hz)	30					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	87.04					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm ²)	3.98					
	I_{spta} at z_{pii} or z_{sii} (mW/cm ²)	23.19					
	p_r at z_{pii} (MPa)	2.10					
Operating Control Conditions	Control 1: 2D Frequency/THI	5 MHz/Off	10 MHz/On		10 MHz/On		-
	Control 2: Color Frequency	4 MHz	6 MHz		6 MHz		-
	Control 3: Color Box Size	Narrow	Narrow		Narrow		-
	Control 4: PW Frequency	4 MHz	6 MHz		4 MHz		-
	Control 5: Color PRF	5530 Hz	5530 Hz		5530 Hz		-
	Control 6: PW PRF	6579 Hz	6579 Hz		6579 Hz		-
	Control 7: Focus	15 mm	35 mm		10 mm		-
	Control 8: Acoustic Power	100%	100%		100%		-
	Control 9: Line Density	1	1		-		-

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Index Label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	
Maximum Index Value		0.73	0.20		0.52		-
Index Component Value			0.20	0.07	0.16	0.42	
Acoustic Parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	1.66					
	P (mW)		1.94	3.48	1.94	3.12	-
	P_{1x1} (mW)		1.94	3.48	1.94	3.12	
	z_s (cm)			1.05			
	z_b (cm)					1.40	
	z_{MI} (cm)	1.33					
	$z_{pii,\alpha}$ (cm)	1.33					
	f_{awf} (MHz)	5.14	6.08	6.13	6.08	4.07	-
Other Information	p_{rr} (Hz)	819					
	s_{rr} (Hz)	30					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	87.04					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm ²)	3.98					
	I_{spta} at z_{pii} or z_{sii} (mW/cm ²)	23.19					
	p_r at z_{pii} (MPa)	2.10					
Operating Control Conditions	Control 1: 2D Frequency/THI	5 MHz/Off	10 MHz/On		10 MHz/On		-
	Control 2: Color Frequency	6 MHz	6 MHz		6 MHz		-
	Control 3: Color Box Size	Narrow	Narrow		Narrow		-
	Control 4: PW Frequency	4 MHz	6 MHz		4 MHz		-
	Control 5: Color PRF	5530 Hz	5530 Hz		5530 Hz		-
	Control 6: PW PRF	6579 Hz	6579 Hz		6579 Hz		-
	Control 7: Focus	15 mm	120 mm		10 mm		-
	Control 8: Acoustic Power	100%	100%		100%		-
	Control 9: Line Density	1	1		-		-

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Index Label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	
Maximum Index Value		1.12	0.79		1.19		-
Index Component Value			0.79	0.11	0.76	0.54	
Acoustic Parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	2.55					
	P (mW)		19.23	4.75	19.23	3.41	-
	$P_{1 \times 1}$ (mW)		19.23	4.75	19.23	3.41	
	z_s (cm)			0.55			
	z_b (cm)					0.85	
	z_{MI} (cm)	0.68					
	$z_{pii,\alpha}$ (cm)	0.68					
	f_{awf} (MHz)	5.19	7.27	6.12	7.27	6.20	-
Other Information	prr (Hz)	4480					
	srr (Hz)	40					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	1.62					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm ²)	1.62					
	I_{spta} at z_{pii} or z_{sii} (mW/cm ²)	263.03					
	p_r at z_{pii} (MPa)	2.88					
Operating Control Conditions	Control 1: 2D Frequency/THI	5MHz/OFF	7.5MHz/OFF		7.5MHz/OFF		-
	Control 2: Color Frequency	7MHz	7MHz		7MHz		-
	Control 3: Color Box Size	Narrow	Narrow		Narrow		-
	Control 4: PW Frequency	4 MHz	6 MHz		6 MHz		-
	Control 5: Color PRF	13158Hz	13158Hz		13158Hz		-
	Control 6: PW PRF	6579 Hz	6579 Hz		6579 Hz		-
	Control 7: Focus	5mm	60mm		10mm		-
	Control 8: Acoustic Power	100%	100%		100%		-
	Control 9: Line Density	1	1		N/A		-

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Index Label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	
Maximum Index Value		1.12	0.56		0.96		-
Index Component Value			0.56	0.11	0.53	0.54	
Acoustic Parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	2.55					
	P (mW)		12.60	4.75	12.60	3.41	-
	P_{1x1} (mW)		12.60	4.75	12.60	3.41	
	z_s (cm)			0.55			
	z_b (cm)					0.85	
	z_{MI} (cm)	0.68					
	$z_{pii,\alpha}$ (cm)	0.68					
	f_{awf} (MHz)	5.19	7.26	6.12	7.26	6.20	-
Other Information	p_{rr} (Hz)	4480					
	s_{rr} (Hz)	40					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	181.24					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm ²)	1.62					
	I_{spta} at z_{pii} or z_{sii} (mW/cm ²)	263.03					
	p_r at z_{pii} (MPa)	2.88					
Operating Control Conditions	Control 1: 2D Frequency/THI	5MHz/OFF	7.5MHz/OFF		7.5MHz/OFF		-
	Control 2: Color Frequency	7 MHz	7 MHz		7 MHz		-
	Control 3: Color Box Size	Narrow	Narrow		Narrow		-
	Control 4: PW Frequency	4 MHz	6 MHz		6 MHz		-
	Control 5: Color PRF	13158Hz	13158Hz		13158Hz		-
	Control 6: PW PRF	6579 Hz	6579 Hz		6579 Hz		-
	Control 7: Focus	5mm	100 mm		10 mm		-
	Control 8: Acoustic Power	100%	100%		100%		-
	Control 9: Line Density	1	1		N/A		-

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Index Label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	
Maximum Index Value		1.21	0.42		0.77		-
Index Component Value			0.42	0.09	0.37	0.47	
Acoustic Parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	2.61					
	P (mW)		14.68	5.95	14.68	3.58	-
	P_{1xl} (mW)		14.68	5.95	14.68	3.58	
	z_s (cm)			0.95			
	z_b (cm)					1.00	
	z_{MI} (cm)	1.18					
	$z_{pii,\alpha}$ (cm)	1.18					
	f_{awf} (MHz)	4.66	4.03	4.05	4.03	4.05	-
Other Information	prf (Hz)	1680					
	srf (Hz)	15					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	218.74					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm ²)	2.09					
	I_{spta} at z_{pii} or z_{sii} (mW/cm ²)	144.50					
	p_r at z_{pii} (MPa)	3.15					
Operating Control Conditions	Control 1: 2D Frequency/THI	8MHz/ON	10MHz/ON		10MHz/ON		-
	Control 2: Color Frequency	4 MHz	4MHz		4MHz		-
	Control 3: Color Box Size	Narrow	Narrow		Narrow		-
	Control 4: PW Frequency	4 MHz	4MHz		4 MHz		-
	Control 5: Color PRF	13158Hz	13158Hz		13158Hz		-
	Control 6: PW PRF	6579 Hz	6579 Hz		6579 Hz		-
	Control 7: Focus	10mm	110mm		10 mm		-
	Control 8: Acoustic Power	100%	100%		100%		-
	Control 9: Line Density	0.5	0.5		-		-

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Index Label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	
Maximum Index Value		1.21	0.26		0.40		-
Index Component Value			0.26	0.04	0.23	0.20	
Acoustic Parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	2.80					
	P (mW)		7.51	2.15	7.51	1.74	-
	$P_{1 \times 1}$ (mW)		7.51	2.15	7.51	1.74	
	z_s (cm)			1.00			
	z_b (cm)					1.18	
	z_{MI} (cm)	1.25					
	$z_{pii,\alpha}$ (cm)	1.25					
	f_{awf} (MHz)	5.31	5.24	6.03	5.24	4.09	-
Other Information	prr (Hz)	947					
	srr (Hz)	16.9					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	347.20					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm ²)	1.45					
	I_{spta} at z_{pii} or z_{sii} (mW/cm ²)	92.36					
	p_r at z_{pii} (MPa)	3.52					
Operating Control Conditions	Control 1: 2D Frequency/THI	8MHz/OFF	10MHz/ON		10MHz/ON		-
	Control 2: Color Frequency	4 MHz	6.15MHz		6.15MHz		-
	Control 3: Color Box Size	Narrow	Narrow		Narrow		-
	Control 4: PW Frequency	4 MHz	6.15MHz		4 MHz		-
	Control 5: Color PRF	13158Hz	13158Hz		13158Hz		-
	Control 6: PW PRF	6579 Hz	6579 Hz		6579 Hz		-
	Control 7: Focus	10mm	80mm		15mm		-
	Control 8: Acoustic Power	100%	100%		100%		-
	Control 9: Line Density	1	1		-		-

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Index Label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	
Maximum Index Value		1.21	0.30		0.65		-
Index Component Value			0.30	0.09	0.25	0.47	
Acoustic Parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	2.61					
	P (mW)		7.03	5.95	7.03	3.58	-
	$P_{1 \times 1}$ (mW)		7.03	5.95	7.03	3.58	
	z_s (cm)			0.95			
	z_b (cm)					1.00	
	z_{MI} (cm)	1.18					
	$z_{pii,\alpha}$ (cm)	1.18					
	f_{awf} (MHz)	4.66	5.16	4.05	5.16	4.05	-
Other Information	pr_r (Hz)	1680					
	sr_r (Hz)	15					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	218.74					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm ²)	2.09					
	I_{spta} at z_{pii} or z_{sii} (mW/cm ²)	144.50					
	p_r at z_{pii} (MPa)	3.15					
Operating Control Conditions	Control 1: 2D Frequency/THI	8MHz/ON	10 MHz/ON		10 MHz/ON		-
	Control 2: Color Frequency	6 MHz	6 MHz		6 MHz		-
	Control 3: Color Box Size	Narrow	Narrow		Narrow		-
	Control 4: PW Frequency	4 MHz	4MHz		4 MHz		-
	Control 5: Color PRF	13158Hz	13158Hz		13158Hz		-
	Control 6: PW PRF	6579 Hz	6579 Hz		6579 Hz		-
	Control 7: Focus	10mm	120 mm		10 mm		-
	Control 8: Acoustic Power	100%	100%		100%		-
	Control 9: Line Density	0.5	0.5		-		-

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Index Label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	
Maximum Index Value		1.21	0.22		0.36		-
Index Component Value			0.22	0.04	0.19	0.20	
Acoustic Parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	2.80					
	P (mW)		6.04	2.15	6.04	1.74	-
	$P_{1 \times 1}$ (mW)		6.04	2.15	6.04	1.74	
	z_s (cm)			1.00			
	z_b (cm)					1.18	
	z_{MI} (cm)	1.25					
	$z_{pii,\alpha}$ (cm)	1.25					
	f_{awf} (MHz)	5.31	5.24	6.03	5.24	4.09	-
Other Information	prr (Hz)	846					
	srr (Hz)	15.1					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	347.20					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm ²)	1.29					
	I_{spta} at z_{pii} or z_{sii} (mW/cm ²)	82.51					
	p_r at z_{pii} (MPa)	3.52					
Operating Control Conditions	Control 1: 2D Frequency/THI	8MHz/OFF	10 MHz/ON		10 MHz/ON		-
	Control 2: Color Frequency	6 MHz	6 MHz		6 MHz		-
	Control 3: Color Box Size	Narrow	Narrow		Narrow		-
	Control 4: PW Frequency	4 MHz	6MHz		4 MHz		-
	Control 5: Color PRF	13158Hz	13158Hz		13158Hz		-
	Control 6: PW PRF	6579 Hz	6579 Hz		6579 Hz		-
	Control 7: Focus	10mm	80mm		15mm		-
	Control 8: Acoustic Power	100%	100%		100%		-
	Control 9: Line Density	1	1		-		-

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
User Setting	
Security	Security
Page 3 / 3	
Close	

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
Security		
Security	User login Setting	Access
Page 1 / 1		
Close		

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
The 'Password Input' screen features three input fields: 'Enter', 'Re-Enter', and 'Current'. The 'Enter' field is the first and is currently active. At the bottom, there are four buttons: 'OK', 'Show', 'Page 1/1', and 'Cancel'.

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
An on-screen keyboard is displayed. The top section includes an 'Enter' field with a cursor, navigation arrows ('<<' and '>>'), and a 'Close' button. Below the input field is a grid of letters (QWERTYUIOP, ASDFGHJKL, ZXCVBNM) and special characters (aA, iü, 1@, ↑, Space). A 'Back Space' button is located to the right of the top row of letters.

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The 'Password Input' screen is shown again. The 'Enter' field now contains four black dots, indicating the password has been masked. The 'Re-Enter' and 'Current' fields are empty. The bottom buttons remain the same.

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The on-screen keyboard is shown again, identical to the one in step (6). The 'Enter' field at the top contains four black dots, indicating the user is about to re-enter the password for confirmation.

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The screenshot shows a 'Patient Info' screen with the following elements:

- Buttons at the top: Patient Info, Exam Info, Hide Info.
- Fields: ID, Accession No., Patient Name (Last [space] First), DOB (YYMMDD), AGE, Gender (U), Preset (EUS-1).
- Buttons at the bottom: Keyboard, Patient Log, Worklist, Exam List, OK.

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
The screenshot shows a 'FUJIFILM CP-1' software keyboard interface with the following elements:

- Top bar: ID, <<, >>, Close.
- Keyboard grid: Q, W, E, R, T, Y, U, I, O, P; A, S, D, F, G, H, J, K, L; Z, X, C, V, B, N, M; aA, uü, i@, ↑, Space.
- Bottom bar: Back Space, five directional arrows (up, down, left, right, center).

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
The screen displays four buttons for gender selection: M, P, O, and U. The 'U' button is highlighted with a grey background. A 'Cancel' button is located at the bottom right.

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
The screen displays two buttons for preset selection: EUS-1 and EUS-2. The 'EUS-1' button is highlighted with a grey background. A 'Close' button is located at the top right.

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The screen shows a form with two tabs: 'Patient Info' and 'Exam Info'. The 'Exam Info' tab is selected. The form contains fields for Study ID, Operator, Phys, Ref. Phys, Rep. Phys, Req. Phys, Department, Residence, Occupation, Patient History, Study Description, Patient State, and Comment. The 'Operator' through 'Req. Phys' fields are pre-filled with 'Last [space] First'. At the bottom, there are buttons for Keyboard, Patient Log, Worklist, Exam List, and OK.

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The screenshot shows the 'Patient Info' screen. At the top, there are tabs for 'Patient Info', 'Exam Info', and 'Hide Info'. Below these are input fields for 'ID:', 'Accession No.:', 'Patient Name:' (with 'Last [space] First' sub-labels), 'DOB:' (with 'YYYY', 'MM', 'DD' sub-labels), 'AGE:', 'Gender:' (with a dropdown menu showing 'U'), and 'Preset:' (with a dropdown menu showing 'EUS-1'). At the bottom, there are five buttons: 'Keyboard', 'Patient Log', 'Worklist', 'Exam List', and 'OK'.

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The screenshot shows the 'Worklist' screen. At the top, it says 'FUJIFILM CP-1' and 'Worklist'. Below this is a search bar with 'Search key:' and a dropdown menu set to 'Patient ID', and a 'String:' input field. To the right are buttons for 'Search', 'Select All', 'Find', 'Query', and 'Delete'. Below the search bar is a table with the following columns: 'Patient ID', 'Patient Name', 'Receipt Number', 'Schedule Procedure Step', 'Start Date-Time', and 'Log'. At the bottom, there are buttons for 'First/Last', 'L/R Scroll', 'Page 1/1', 'Open', and 'Close'.

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The screenshot shows the 'FUJIFILM CP-1' interface. At the top is a 'Worklist' header. Below it, there are search controls: a 'Search key:' dropdown (1) set to 'Patient ID', a 'String:' text input (2), and a 'Search' button (13). Below these are several buttons: 'Select All' (5), 'Find' (6), 'Query' (11), and 'Delete' (12). The main area is a large table (4) with columns: 'Patient ID', 'Patient Name', 'Receipt Number', 'Schedule Procedure Step Start Date-Time', and 'Sch Pro Des'. At the bottom, there are more controls: 'First/Last' (7), 'L/R Scroll' (8), 'Page 1/1' (9), an 'Open' button (10), and a 'Close' button (14).

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- Patient ID
- Patient Name
- Accession No.
- Scheduled Procedure Step Start Date-Time
- Scheduled Procedure Step Description
- Requested Procedure ID
- Requested Procedure Description

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The 'Patient Info' screen features a top navigation bar with 'Patient Info', 'Exam Info', and 'Hide Info' buttons. Below this, there are input fields for 'ID', 'Accession No.', 'Patient Name' (split into 'Last' and 'First'), 'DOB' (YYMMDD), 'AGE', and 'Gender' (a dropdown menu currently showing 'U'). A 'Preset' dropdown menu is set to 'EUS-1'. At the bottom, a row of buttons includes 'Keyboard', 'Patient Log', 'Worklist', 'Exam List', and 'OK'.

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The 'Exam List' screen has a title bar 'FUJIFILM CP-1' and a sub-header 'Exam List'. It includes a search section with a 'Search key' dropdown (set to 'Patient Name'), a 'String' input field, and 'Search', 'Load', 'Select All', and 'Delete' buttons. The main area is a table with columns 'Patient ID', 'Patient Name', and 'Receipt Number'. At the bottom, there are buttons for 'First/Last', 'L/R Scroll', 'Page 1/1', 'Open', and 'Close'.

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The screenshot shows the 'Exam List' screen on a FUJIFILM CP-1 device. The screen is titled 'Exam List' and contains a search bar at the top with a dropdown menu set to 'Patient Name' (callout 1) and a 'String:' input field (callout 2). To the right of the search bar are buttons for 'Search' (callout 9), 'Load' (callout 10), and 'Select All' (callout 11). Below the search bar is a table with three columns: 'Patient ID' (callout 3), 'Patient Name' (callout 4), and 'Receipt Number' (callout 5). The table is currently empty. To the right of the table is a 'Delete' button (callout 12). At the bottom of the screen are buttons for 'First/Last' (callout 6), 'L/R Scroll' (callout 7), 'Page 1/1' (callout 8), 'Open' (callout 9), and 'Close' (callout 10).

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- Patient ID
- Patient Name
- Accession No.

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< A >

A/B A/B Ratio
AC Alternating Current
Accel..... Acceleration
AT Delta Time

< B >

BS..... Backspace

< C >

Circ..... Circumference

< D >

dB..... Decibel
Delta T Delta Time
Delta V Delta Velocity
DICOM..... Digital Imaging and
Communications in Medicine
Dist..... Distance
D max..... Maximum Diameter
D min..... Minimum Diameter
DOB..... Date of Birth
D.RNG..... Dynamic Range

< E >

EDV End Diastolic Velocity
EMC Electro Mechanical
Compliance

< F >

FORM..... Format
FREQ..... Frequency

< H >

HR..... Heart-Rate

< I >

ID Identity, Identification
Number
IMG Image
IP Internet protocol

< J >

JPEG..... Joint Photographic Expert
Group

< L >

LAN Local area network
LED..... Light-emitting diode
L/R..... Left/Right

< M >

MDV..... Minimum Diastolic Velocity
MI Mechanical Index
MPPS..... Modality Performed
Procedure

< IN >

NTSC..... National Television Standards

< P >

PAL..... Phase Alternating Line
PI..... Pulsatility Index
PinP Picture in Picture
PSV Peak Systolic Velocity
PW..... Pulse wave Doppler

< R >

RI..... Resistance Index
ROI Region of interest

< S >

S/D..... Systolic/Diastolic
STC Sensitivity time control
SV..... Doppler Sample Volume

< T >

TAMX..... Time Average of Maximum
Velocity
TCP..... Transmission Control
Protocol
TH..... Tissue harmonic
TI..... Thermal index
TIFF..... Tagged Image File Format

< U >

US Ultrasonic

USB..... Universal serial bus

< V >

Vel Velocity

< Others >

% Stenosis Dist Percent Diameter
Stenosis

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<本製品で使用しているソフトウェアライセンス表示>

Free-Type2

This software is based in part of the work of FreeType Team.

HarfBuzz

This software is based in part of the work of FreeType Team.

libjpeg

This software is based in part on the work of the Independent JPEG Group.

< DCMTK >

Unless otherwise specified, the DCMTK software package has the following copyright:

```
/*
 * Copyright (C) 1994-2011, OFFIS e.V.
 * All rights reserved.
 *
 * This software and supporting documentation were developed by
 *
 * OFFIS e.V.
 * R&D Division Health
 * Escherweg 2
 * 26121 Oldenburg, Germany
 *
```

* Redistribution and use in source and binary forms, with or without
* modification, are permitted provided that the following conditions
* are met:
*
* - Redistributions of source code must retain the above copyright
* notice, this list of conditions and the following disclaimer.
*
* - Redistributions in binary form must reproduce the above copyright
* notice, this list of conditions and the following disclaimer in the
* documentation and/or other materials provided with the distribution.
*
* - Neither the name of OFFIS nor the names of its contributors may be
* used to endorse or promote products derived from this software
* without specific prior written permission.
*
* THIS SOFTWARE IS PROVIDED BY THE COPYRIGHT HOLDERS AND
* CONTRIBUTORS "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES,
* INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF
* MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE
* DISCLAIMED. IN NO EVENT SHALL THE COPYRIGHT HOLDER OR
* CONTRIBUTORS BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL,
* SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING, BUT
* NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES;
* LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER
* CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT,
* STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE)
* ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE, EVEN IF
* ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.
*
*/

Please note that some DCMTK modules, especially those that are not part of the free toolkit, are covered by a separate license which can be found in the COPYRIGHT file in the corresponding module directory.

Some portions of the DCMTK software package are derived from earlier versions of this software with the following copyright, and can be identified by the following copyright notice located in each source file:

/*
* Copyright (C) 1993/1994, OFFIS, Oldenburg University and CERIU
*
* This software and supporting documentation were
* developed by

*
 * Institut OFFIS
 * Bereich Kommunikationssysteme
 * Westerstr. 10-12
 * 26121 Oldenburg, Germany
 *
 * Fachbereich Informatik
 * Abteilung Prozessinformatik
 * Carl von Ossietzky Universitaet Oldenburg
 * Ammerlaender Heerstr. 114-118
 * 26111 Oldenburg, Germany
 *
 * CERIUM
 * Laboratoire SIM
 * Faculte de Medecine
 * 2 Avenue du Pr. Leon Bernard
 * 35043 Rennes Cedex, France
 *
 * for CEN/TC251/WG4 as a contribution to the Radiological
 * Society of North America (RSNA) 1993 Digital Imaging and
 * Communications in Medicine (DICOM) Demonstration.
 *
 * THIS SOFTWARE IS MADE AVAILABLE, AS IS, AND NEITHER OFFIS,
 * OLDENBURG UNIVERSITY NOR CERIUM MAKE ANY WARRANTY REGARDING
 * THE SOFTWARE, ITS PERFORMANCE, ITS MERCHANTABILITY OR
 * FITNESS FOR ANY PARTICULAR USE, FREEDOM FROM ANY COMPUTER
 * DISEASES OR ITS CONFORMITY TO ANY SPECIFICATION. THE
 * ENTIRE RISK AS TO QUALITY AND PERFORMANCE OF THE SOFTWARE
 * IS WITH THE USER.
 *
 * Copyright of the software and supporting documentation
 * is, unless otherwise stated, jointly owned by OFFIS,
 * Oldenburg University and CERIUM and free access is hereby
 * granted as a license to use this software, copy this
 * software and prepare derivative works based upon this
 * software. However, any distribution of this software
 * source code or supporting documentation or derivative
 * works (source code and supporting documentation) must
 * include the three paragraphs of this copyright notice.
 *
 */

Some other parts of this software within the dcmTk/dcmnet sub-package related to the DICOM Upper Layer Protocol are derived from software developed for the RSNA'93 DICOM demonstration and kindly made available to us by the Mallinckrodt Institute of Radiology. Such software can be identified by the following copyright notice located in each affected source file:

```

/*
 * Copyright (C) 1993, RSNA and Washington University
 *
 * The software and supporting documentation for the Radiological
 * Society of North America (RSNA) 1993 Digital Imaging and
 * Communications in Medicine (DICOM) Demonstration were developed
 * at the
 *   Electronic Radiology Laboratory
 *   Mallinckrodt Institute of Radiology
 *   Washington University School of Medicine
 *   510 S. Kingshighway Blvd.
 *   St. Louis, MO 63110
 * as part of the 1993 DICOM Central Test Node project for, and
 * under contract with, the Radiological Society of North America.
 *
 * THIS SOFTWARE IS MADE AVAILABLE, AS IS, AND NEITHER RSNA NOR
 * WASHINGTON UNIVERSITY MAKE ANY WARRANTY ABOUT THE SOFTWARE,
 * ITS PERFORMANCE, ITS MERCHANTABILITY OR FITNESS FOR ANY
 * PARTICULAR USE, FREEDOM FROM ANY COMPUTER DISEASES OR ITS
 * CONFORMITY TO ANY SPECIFICATION. THE ENTIRE RISK AS TO QUALITY
 * AND PERFORMANCE OF THE SOFTWARE IS WITH THE USER.
 *
 * Copyright of the software and supporting documentation is
 * jointly owned by RSNA and Washington University, and free access
 * is hereby granted as a license to use this software, copy this
 * software and prepare derivative works based upon this software.
 * However, any distribution of this software source code or
 * supporting documentation or derivative works (source code and
 * supporting documentation) must include the three paragraphs of
 * the copyright notice.
 */

```

The dcmjpeg sub-package includes an adapted version of the Independent JPEG Group Toolkit Version 6b, which is contained in dcmjpeg/libijg8, dcmjpeg/libijg12 and dcmjpeg/libijg16. This toolkit is covered by the following copyright. The original README file for the Independent JPEG Group Toolkit is located in dcmjpeg/docs/ijg_readme.txt.

```

/*
 * The authors make NO WARRANTY or representation, either express or
 * implied, with respect to this software, its quality, accuracy, merchantability, or
 * fitness for a particular purpose. This software is provided "AS IS", and you,
 * its user, assume the entire risk as to its quality and accuracy.
 *
 * This software is copyright (C) 1991-1998, Thomas G. Lane.
 * All Rights Reserved except as specified below.
 *
 * Permission is hereby granted to use, copy, modify, and distribute this
 * software (or portions thereof) for any purpose, without fee, subject to these
 * conditions:
 * (1) If any part of the source code for this software is distributed, then this
 * README file must be included, with this copyright and no-warranty notice
 * unaltered; and any additions, deletions, or changes to the original files
 * must be clearly indicated in accompanying documentation.
 * (2) If only executable code is distributed, then the accompanying
 * documentation must state that "this software is based in part on the work of
 * the Independent JPEG Group".
 * (3) Permission for use of this software is granted only if the user accepts
 * full responsibility for any undesirable consequences; the authors accept
 * NO LIABILITY for damages of any kind.
 *
 * These conditions apply to any software derived from or based on the IJG code,
 * not just to the unmodified library. If you use our work, you ought to
 * acknowledge us.
 *
 * Permission is NOT granted for the use of any IJG author's name or company
 * name
 * in advertising or publicity relating to this software or products derived from
 * it. This software may be referred to only as "the Independent JPEG Group's
 * software".
 *
 * We specifically permit and encourage the use of this software as the basis of
 * commercial products, provided that all warranty or liability claims are
 * assumed by the product vendor.
 */

```

The code for the `interpolatePixel()` image scaling algorithm in module `dcmimgle` has been derived from code written by Jef Poskanzer for the "Extended Portable Bitmap Toolkit" (`pbmplus10dec91`) which has the following copyright:

```

/*
 * Copyright (C) 1991 by Jef Poskanzer.
 *

```

```
* Permission to use, copy, modify, and distribute this software and its
* documentation for any purpose and without fee is hereby granted, provided
* that the above copyright notice appear in all copies and that both that
* copyright notice and this permission notice appear in supporting
* documentation. This software is provided "as is" without express or
* implied warranty.
*/
```

The color quantization code in module dcmimage (dcmquant and the related classes) has been derived from code written by Jef Poskanzer for the NetPBM toolkit which has the following copyright:

```
/*
* Copyright (C) 1989, 1991 by Jef Poskanzer.
*
* Permission to use, copy, modify, and distribute this software and its
* documentation for any purpose and without fee is hereby granted, provided
* that the above copyright notice appear in all copies and that both that
* copyright notice and this permission notice appear in supporting
* documentation. This software is provided "as is" without express or
* implied warranty.
*/
```

The code for the OFStandard::strncpy and OFStandard::strcat helper functions in ofstd/libsrc/ofstd.cc has been derived from the BSD implementation of strncpy() and strcat() and which carries the following copyright notice:

```
/*
* Copyright (c) 1998 Todd C. Miller <Todd.Miller(at)courtesan.com>
* All rights reserved.
*
* Redistribution and use in source and binary forms, with or without
* modification, are permitted provided that the following conditions
* are met:
* 1. Redistributions of source code must retain the above copyright
* notice, this list of conditions and the following disclaimer.
* 2. Redistributions in binary form must reproduce the above copyright
* notice, this list of conditions and the following disclaimer in the
* documentation and/or other materials provided with the distribution.
* 3. The name of the author may not be used to endorse or promote products
* derived from this software without specific prior written permission.
*
* THIS SOFTWARE IS PROVIDED ``AS IS" AND ANY EXPRESS OR IMPLIED
* WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED
* WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR
```

```

* PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL THE AUTHOR BE LIABLE
* FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR
* CONSEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED TO,
* PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA,
* OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY
* THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR
* TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT
* OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF
* SUCH DAMAGE.
*/

```

The code for the OFStandard::atof helper function in ofstd/libsrc/ofstd.cc has been derived from an implementation which carries the following copyright notice:

```

/*
* Copyright 1988 Regents of the University of California
* Permission to use, copy, modify, and distribute this software and
* its documentation for any purpose and without fee is hereby granted,
* provided that the above copyright notice appear in all copies. The
* University of California makes no representations about the
* suitability of this software for any purpose. It is provided "as
* is" without express or implied warranty.
*
*
* The code for OFStandard::ftoa has been derived
* from an implementation which carries the following copyright notice:
*
* Copyright (c) 1988 Regents of the University of California.
* All rights reserved.
*
* Redistribution and use in source and binary forms are permitted
* provided that the above copyright notice and this paragraph are
* duplicated in all such forms and that any documentation,
* advertising materials, and other materials related to such
* distribution and use acknowledge that the software was developed
* by the University of California, Berkeley. The name of the
* University may not be used to endorse or promote products derived
* from this software without specific prior written permission.
* THIS SOFTWARE IS PROVIDED ``AS IS'' AND WITHOUT ANY EXPRESS OR
* IMPLIED WARRANTIES, INCLUDING, WITHOUT LIMITATION, THE IMPLIED
* WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR
* PURPOSE.
*/

```

The "Base64" encoder/decoder in ofstd/libsrc/ofstd.cc has been derived from an implementation which carries the following copyright notice:

```
/*
 * Copyright (c) 1999, Bob Withers - bwit(at)pobox.com
 *
 * This code may be freely used for any purpose, either personal or commercial,
 * provided the authors copyright notice remains intact.
 */
```

The oflog sub-package is based on the log4cplus library which is covered by the following two copyright notices (for details see oflog/docs/LICENSE):

```
/*
 * Copyright (C) 1999-2009 Contributors to log4cplus project.
 * All rights reserved.
 *
 * Redistribution and use in source and binary forms, with or without modifica-
 * tion, are permitted provided that the following conditions are met:
 *
 * 1. Redistributions of source code must retain the above copyright notice,
 * this list of conditions and the following disclaimer.
 *
 * 2. Redistributions in binary form must reproduce the above copyright notice,
 * this list of conditions and the following disclaimer in the documentation
 * and/or other materials provided with the distribution.
 *
 * THIS SOFTWARE IS PROVIDED ``AS IS" AND ANY EXPRESSED OR IMPLIED
 * WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED
 * WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR
 * PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL THE APACHE SOFTWARE
 * FOUNDATION OR ITS CONTRIBUTORS BE LIABLE FOR ANY DIRECT,
 * INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL
 * DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF
 * SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR
 * BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF
 * LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT
 * (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF
 * THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF
 * SUCH DAMAGE.
 */
```

```
// Module: Log4CPLUS
// File: appender.h
// Created: 6/2001
// Author: Tad E. Smith
//
//
// Copyright 2001-2009 Tad E. Smith
//
// Licensed under the Apache License, Version 2.0 (the "License");
// you may not use this file except in compliance with the License.
// You may obtain a copy of the License at
//
// http://www.apache.org/licenses/LICENSE-2.0
//
// Unless required by applicable law or agreed to in writing, software
// distributed under the License is distributed on an "AS IS" BASIS,
// WITHOUT WARRANTIES OR CONDITIONS OF ANY KIND, either express or
// implied.
// See the License for the specific language governing permissions and
// limitations under the License.
```

The dcmjpls sub-package is based on the CharLS library, which is contained in dcmjpls/libcharls. This toolkit is covered by the following copyright:

```
/*
 * The CharLS library is available under the following license:
 * -----
 * Copyright (c) 2007-2010, Jan de Vaan
 * All rights reserved.
 *
 * Redistribution and use in source and binary forms, with or without
 * modification, are permitted provided that the following conditions are met:
 *
 * * Redistributions of source code must retain the above copyright notice, this
 *   list of conditions and the following disclaimer.
 *
 * * Redistributions in binary form must reproduce the above copyright notice,
 *   this list of conditions and the following disclaimer in the documentation
 *   and/or other materials provided with the distribution.
 *
 * * Neither the name of my employer, nor the names of its contributors may be
 *   used to endorse or promote products derived from this software without
 *   specific prior written permission.
 *
```

* THIS SOFTWARE IS PROVIDED BY THE COPYRIGHT HOLDERS AND
* CONTRIBUTORS "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES,
* INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF
* MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE
* DISCLAIMED. IN NO EVENT SHALL THE COPYRIGHT OWNER OR
* CONTRIBUTORS BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL,
* SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING, BUT
* NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES;
* LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER
* CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT,
* STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE)
* ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE, EVEN IF
* ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.
*/

< LiteSQL >

LiteSQL is distributed under the terms of BSD-license stated below:

Copyright (c) 2009, Tero Laitinen and all contributors All rights reserved.

Redistribution and use in source and binary forms, with or without modification,
are permitted provided that the following conditions are met:

Redistributions of source code must retain the above copyright notice, this list of
conditions and the following disclaimer.

Redistributions in binary form must reproduce the above copyright notice, this
list of conditions and the following disclaimer in the documentation and/or other
materials provided with the distribution.

Neither the name of the Tero Laitinen nor the names of its contributors may be
used to endorse or promote products derived from this software without specific
prior written permission.

THIS SOFTWARE IS PROVIDED BY THE COPYRIGHT HOLDERS AND
CONTRIBUTORS "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES,
INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF
MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE
ARE DISCLAIMED. IN NO EVENT SHALL THE COPYRIGHT OWNER OR
CONTRIBUTORS BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL,
EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED
TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA,
OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY
THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT
(INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE
USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH
DAMAGE.

< MkTemp >

Mktemp is distributed under the following BSD-style license:

Copyright (c) 1996, 2000, 2001 Todd C. Miller <Todd.Miller@courtesan.com>
All rights reserved.

Redistribution and use in source and binary forms, with or without modification, are permitted provided that the following conditions are met:

1. Redistributions of source code must retain the above copyright notice, this list of conditions and the following disclaimer.
2. Redistributions in binary form must reproduce the above copyright notice, this list of conditions and the following disclaimer in the documentation and/or other materials provided with the distribution.
3. The name of the author may not be used to endorse or promote products derived from this software without specific prior written permission from the author.

THIS SOFTWARE IS PROVIDED ``AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL THE AUTHOR BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.

Additionally, priv_mktemp.c bears the following UCB license:

Copyright (c) 1987, 1993

The Regents of the University of California. All rights reserved.

Redistribution and use in source and binary forms, with or without modification, are permitted provided that the following conditions are met:

1. Redistributions of source code must retain the above copyright notice, this list of conditions and the following disclaimer.
2. Redistributions in binary form must reproduce the above copyright notice, this list of conditions and the following disclaimer in the documentation and/or other materials provided with the distribution.

3. Neither the name of the University nor the names of its contributors may be used to endorse or promote products derived from this software without specific prior written permission.

THIS SOFTWARE IS PROVIDED BY THE REGENTS AND CONTRIBUTORS ``AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL THE REGENTS OR CONTRIBUTORS BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.

< Ncurses >

Copyright (c) 1998-2004,2006 Free Software Foundation, Inc.

Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the "Software"), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, distribute with modifications, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software.

THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE ABOVE COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

Except as contained in this notice, the name(s) of the above copyright holders shall not be used in advertising or otherwise to promote the sale, use or other dealings in this Software without prior written authorization.

< Popt >

Copyright (c) 1998 Red Hat Software

Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the "Software"), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software.

THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE X CONSORTIUM BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

Except as contained in this notice, the name of the X Consortium shall not be used in advertising or otherwise to promote the sale, use or other dealings in this Software without prior written authorization from the X Consortium.

< TOPPERS/JSP Kernel >

Toyohashi Open Platform for Embedded Real-Time Systems/
Just Standard Profile Kernel

Copyright (C) 2001-2003 by Embedded and Real-Time Systems Laboratory
Toyohashi Univ. of Technology, JAPAN

The above copyright holders grant permission gratis to use, duplicate, modify, or redistribute (hereafter called use) this software (including the one made by modifying this software), provided that the following four conditions (1) through (4) are satisfied.

- (1) When this software is used in the form of source code, the above copyright notice, this use conditions, and the disclaimer shown below must be retained in the source code without modification.
- (2) When this software is redistributed in the forms usable for the development of other software, such as in library form, the above copyright notice, this use conditions, and the disclaimer shown below must be shown without modification in the document provided with the redistributed software, such as the user manual.

- (3) When this software is redistributed in the forms unusable for the development of other software, such as the case when the software is embedded in a piece of equipment, either of the following two conditions must be satisfied:
 - (a) The above copyright notice, this use conditions, and the disclaimer shown below must be shown without modification in the document provided with the redistributed software, such as the user manual.
 - (b) How the software is to be redistributed must be reported to the TOPPERS Project according to the procedure described separately.
- (4) The above copyright holders and the TOPPERS Project are exempt from responsibility for any type of damage directly or indirectly caused from the use of this software and are indemnified by any users or end users of this software from any and all causes of action whatsoever.

THIS SOFTWARE IS PROVIDED "AS IS." THE ABOVE COPYRIGHT HOLDERS AND THE TOPPERS PROJECT DISCLAIM ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, ITS APPLICABILITY TO A PARTICULAR PURPOSE. IN NO EVENT SHALL THE ABOVE COPYRIGHT HOLDERS AND THE TOPPERS PROJECT BE LIABLE FOR ANY TYPE OF DAMAGE DIRECTLY OR INDIRECTLY CAUSED FROM THE USE OF THIS SOFTWARE.

This product uses the software based on GNU General Public License (hereinafter called GPL).

[Note] GPL is the free software license advocated by GNU project.

Permission to use, copy, modify and redistribute the GPL software is granted without fee. For distribution and copy of this software, the source program should be provided.

If you would like the source program of this software, contact the sales agent from whom you purchased this product or our service representatives. It will be provided at actual cost.

Note that you will be provided with only the free software based on GPL in the software of which this product consists. The software developed and created independently for this product is excluded.

Package List:

- <Qt>
- <linux-kernel>
- <alsa-lib>
- <alsa-utils>
- <bash>
- <Busybox>
- <Floppy formatting utility>
- <glibc>
- <libblkid>
- <MTD Utilities>

<GCC>
<uClibc>

GNU GENERAL PUBLIC LICENSE Version 2, June 1991

Copyright (C) 1989, 1991 Free Software Foundation, Inc.
51 Franklin St, Fifth Floor, Boston, MA 02110-1301 USA
Everyone is permitted to copy and distribute verbatim copies of this license document, but changing it is not allowed.

Preamble

The licenses for most software are designed to take away your freedom to share and change it. By contrast, the GNU General Public License is intended to guarantee your freedom to share and change free software--to make sure the software is free for all its users. This General Public License applies to most of the Free Software Foundation's software and to any other program whose authors commit to using it. (Some other Free Software Foundation software is covered by the GNU Library General Public License instead.) You can apply it to your programs, too.

When we speak of free software, we are referring to freedom, not price. Our General Public Licenses are designed to make sure that you have the freedom to distribute copies of free software (and charge for this service if you wish), that you receive source code or can get it if you want it, that you can change the software or use pieces of it in new free programs; and that you know you can do these things.

To protect your rights, we need to make restrictions that forbid anyone to deny you these rights or to ask you to surrender the rights.

These restrictions translate to certain responsibilities for you if you distribute copies of the software, or if you modify it.

For example, if you distribute copies of such a program, whether gratis or for a fee, you must give the recipients all the rights that you have. You must make sure that they, too, receive or can get the source code. And you must show them these terms so they know their rights.

We protect your rights with two steps: (1) copyright the software, and (2) offer you this license which gives you legal permission to copy, distribute and/or modify the software.

Also, for each author's protection and ours, we want to make certain that everyone understands that there is no warranty for this free software. If the software is modified by someone else and passed on, we want its recipients to know that what they have is not the original, so that any problems introduced by others will not reflect on the original authors' reputations.

Finally, any free program is threatened constantly by software patents. We wish to avoid the danger that redistributors of a free program will individually obtain patent licenses, in effect making the program proprietary. To prevent this, we have made it clear that any patent must be licensed for everyone's free use or not licensed at all.

The precise terms and conditions for copying, distribution and modification follow.

GNU GENERAL PUBLIC LICENSE TERMS AND CONDITIONS FOR COPYING, DISTRIBUTION AND MODIFICATION

0. This License applies to any program or other work which contains a notice placed by the copyright holder saying it may be distributed under the terms of this General Public License. The "Program", below, refers to any such program or work, and a "work based on the Program" means either the Program or any derivative work under copyright law: that is to say, a work containing the Program or a portion of it, either verbatim or with modifications and/or translated into another language. (Hereinafter, translation is included without limitation in the term "modification".) Each licensee is addressed as "you".

Activities other than copying, distribution and modification are not covered by this License; they are outside its scope. The act of running the Program is not restricted, and the output from the Program is covered only if its contents constitute a work based on the Program (independent of having been made by running the Program).

Whether that is true depends on what the Program does.

1. You may copy and distribute verbatim copies of the Program's source code as you receive it, in any medium, provided that you conspicuously and appropriately publish on each copy an appropriate copyright notice and disclaimer of warranty; keep intact all the notices that refer to this License and to the absence of any warranty; and give any other recipients of the Program a copy of this License along with the Program.

You may charge a fee for the physical act of transferring a copy, and you may at your option offer warranty protection in exchange for a fee.

2. You may modify your copy or copies of the Program or any portion of it, thus forming a work based on the Program, and copy and distribute such modifications or work under the terms of Section 1 above, provided that you also meet all of these conditions:
 - a) You must cause the modified files to carry prominent notices stating that you changed the files and the date of any change.
 - b) You must cause any work that you distribute or publish, that in whole or in part contains or is derived from the Program or any part thereof, to be licensed as a whole at no charge to all third parties under the terms of this License.
 - c) If the modified program normally reads commands interactively when run, you must cause it, when started running for such interactive use in the most ordinary way, to print or display an announcement including an appropriate copyright notice and a notice that there is no warranty (or else, saying that you provide a warranty) and that users may redistribute the program under these conditions, and telling the user how to view a copy of this License. (Exception: if the Program itself is interactive but does not normally print such an announcement, your work based on the Program is not required to print an announcement.)

These requirements apply to the modified work as a whole. If identifiable sections of that work are not derived from the Program, and can be reasonably considered independent and separate works in themselves, then this License, and its terms, do not apply to those sections when you distribute them as separate works. But when you distribute the same sections as part of a whole which is a work based on the Program, the distribution of the whole must be on the terms of this License, whose permissions for other licensees extend to the entire whole, and thus to each and every part regardless of who wrote it.

Thus, it is not the intent of this section to claim rights or contest your rights to work written entirely by you; rather, the intent is to exercise the right to control the distribution of derivative or collective works based on the Program.

In addition, mere aggregation of another work not based on the Program with the Program (or with a work based on the Program) on a volume of a storage or distribution medium does not bring the other work under the scope of this License.

3. You may copy and distribute the Program (or a work based on it, under Section 2) in object code or executable form under the terms of Sections 1 and 2 above provided that you also do one of the following:
 - a) Accompany it with the complete corresponding machine-readable source code, which must be distributed under the terms of Sections 1 and 2 above on a medium customarily used for software interchange; or,
 - b) Accompany it with a written offer, valid for at least three years, to give any third party, for a charge no more than your cost of physically performing source distribution, a complete machine-readable copy of the corresponding source code, to be distributed under the terms of Sections 1 and 2 above on a medium customarily used for software interchange; or,
 - c) Accompany it with the information you received as to the offer to distribute corresponding source code. (This alternative is allowed only for noncommercial distribution and only if you received the program in object code or executable form with such an offer, in accord with Subsection b above.)

The source code for a work means the preferred form of the work for making modifications to it. For an executable work, complete source code means all the source code for all modules it contains, plus any associated interface definition files, plus the scripts used to control compilation and installation of the executable. However, as a special exception, the source code distributed need not include anything that is normally distributed (in either source or binary form) with the major components (compiler, kernel, and so on) of the operating system on which the executable runs, unless that component itself accompanies the executable.

If distribution of executable or object code is made by offering access to copy from a designated place, then offering equivalent access to copy the source code from the same place counts as distribution of the source code, even though third parties are not compelled to copy the source along with the object code.

4. You may not copy, modify, sublicense, or distribute the Program except as expressly provided under this License. Any attempt otherwise to copy, modify, sublicense or distribute the Program is void, and will automatically terminate your rights under this License.
However, parties who have received copies, or rights, from you under this License will not have their licenses terminated so long as such parties remain in full compliance.

5. You are not required to accept this License, since you have not signed it. However, nothing else grants you permission to modify or distribute the Program or its derivative works. These actions are prohibited by law if you do not accept this License. Therefore, by modifying or distributing the Program (or any work based on the Program), you indicate your acceptance of this License to do so, and all its terms and conditions for copying, distributing or modifying the Program or works based on it.
6. Each time you redistribute the Program (or any work based on the Program), the recipient automatically receives a license from the original licensor to copy, distribute or modify the Program subject to these terms and conditions. You may not impose any further restrictions on the recipients' exercise of the rights granted herein.
You are not responsible for enforcing compliance by third parties to this License.
7. If, as a consequence of a court judgment or allegation of patent infringement or for any other reason (not limited to patent issues), conditions are imposed on you (whether by court order, agreement or otherwise) that contradict the conditions of this License, they do not excuse you from the conditions of this License. If you cannot distribute so as to satisfy simultaneously your obligations under this License and any other pertinent obligations, then as a consequence you may not distribute the Program at all. For example, if a patent license would not permit royalty-free redistribution of the Program by all those who receive copies directly or indirectly through you, then the only way you could satisfy both it and this License would be to refrain entirely from distribution of the Program.

If any portion of this section is held invalid or unenforceable under any particular circumstance, the balance of the section is intended to apply and the section as a whole is intended to apply in other circumstances.

It is not the purpose of this section to induce you to infringe any patents or other property right claims or to contest validity of any such claims; this section has the sole purpose of protecting the integrity of the free software distribution system, which is implemented by public license practices. Many people have made generous contributions to the wide range of software distributed through that system in reliance on consistent application of that system; it is up to the author/donor to decide if he or she is willing to distribute software through any other system and a licensee cannot impose that choice.

This section is intended to make thoroughly clear what is believed to be a consequence of the rest of this License.

8. If the distribution and/or use of the Program is restricted in certain countries either by patents or by copyrighted interfaces, the original copyright holder who places the Program under this License may add an explicit geographical distribution limitation excluding those countries, so that distribution is permitted only in or among countries not thus excluded. In such case, this License incorporates the limitation as if written in the body of this License.
9. The Free Software Foundation may publish revised and/or new versions of the General Public License from time to time. Such new versions will be similar in spirit to the present version, but may differ in detail to address new problems or concerns.

Each version is given a distinguishing version number. If the Program specifies a version number of this License which applies to it and "any later version", you have the option of following the terms and conditions either of that version or of any later version published by the Free Software Foundation. If the Program does not specify a version number of this License, you may choose any version ever published by the Free Software Foundation.

10. If you wish to incorporate parts of the Program into other free programs whose distribution conditions are different, write to the author to ask for permission. For software which is copyrighted by the Free Software Foundation, write to the Free Software Foundation; we sometimes make exceptions for this. Our decision will be guided by the two goals of preserving the free status of all derivatives of our free software and of promoting the sharing and reuse of software generally.

NO WARRANTY

11. BECAUSE THE PROGRAM IS LICENSED FREE OF CHARGE, THERE IS NO WARRANTY FOR THE PROGRAM, TO THE EXTENT PERMITTED BY APPLICABLE LAW. EXCEPT WHEN OTHERWISE STATED IN WRITING THE COPYRIGHT HOLDERS AND/OR OTHER PARTIES PROVIDE THE PROGRAM "AS IS" WITHOUT WARRANTY OF ANY KIND, EITHER EXPRESSED OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE ENTIRE RISK AS TO THE QUALITY AND PERFORMANCE OF THE PROGRAM IS WITH YOU. SHOULD THE PROGRAM PROVE DEFECTIVE, YOU ASSUME THE COST OF ALL NECESSARY SERVICING, REPAIR OR CORRECTION.
12. IN NO EVENT UNLESS REQUIRED BY APPLICABLE LAW OR AGREED TO IN WRITING WILL ANY COPYRIGHT HOLDER, OR ANY OTHER PARTY WHO MAY MODIFY AND/OR REDISTRIBUTE THE PROGRAM AS PERMITTED ABOVE, BE LIABLE TO YOU FOR DAMAGES, INCLUDING ANY GENERAL,

SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES ARISING OUT OF THE USE OR INABILITY TO USE THE PROGRAM (INCLUDING BUT NOT LIMITED TO LOSS OF DATA OR DATA BEING RENDERED INACCURATE OR LOSSES SUSTAINED BY YOU OR THIRD PARTIES OR A FAILURE OF THE PROGRAM TO OPERATE WITH ANY OTHER PROGRAMS), EVEN IF SUCH HOLDER OR OTHER PARTY HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

END OF TERMS AND CONDITIONS

How to Apply These Terms to Your New Programs

If you develop a new program, and you want it to be of the greatest possible use to the public, the best way to achieve this is to make it free software which everyone can redistribute and change under these terms.

To do so, attach the following notices to the program. It is safest to attach them to the start of each source file to most effectively convey the exclusion of warranty; and each file should have at least the "copyright" line and a pointer to where the full notice is found.

<one line to give the program's name and a brief idea of what it does.>
Copyright (C) <year> <name of author>

This program is free software; you can redistribute it and/or modify it under the terms of the GNU General Public License as published by the Free Software Foundation; either version 2 of the License, or (at your option) any later version.

This program is distributed in the hope that it will be useful, but WITHOUT ANY WARRANTY; without even the implied warranty of MERCHANTABILITY or FITNESS FOR A PARTICULAR PURPOSE. See the GNU General Public License for more details.

You should have received a copy of the GNU General Public License along with this program; if not, write to the Free Software Foundation, Inc., 51 Franklin St, Fifth Floor, Boston, MA 02110-1301 USA

Also add information on how to contact you by electronic and paper mail.

If the program is interactive, make it output a short notice like this when it starts in an interactive mode:

Gnomovision version 69, Copyright (C) year name of author
Gnomovision comes with ABSOLUTELY NO WARRANTY; for details type `show w'.

This is free software, and you are welcome to redistribute it under certain conditions; type ``show c'` for details.

The hypothetical commands ``show w'` and ``show c'` should show the appropriate parts of the General Public License. Of course, the commands you use may be called something other than ``show w'` and ``show c'`; they could even be mouse-clicks or menu items--whatever suits your program.

You should also get your employer (if you work as a programmer) or your school, if any, to sign a "copyright disclaimer" for the program, if necessary. Here is a sample; alter the names:

Yoyodyne, Inc., hereby disclaims all copyright interest in the program ``Gnomovision'` (which makes passes at compilers) written by James Hacker.

<signature of Ty Coon>, 1 April 1989
Ty Coon, President of Vice

This General Public License does not permit incorporating your program into proprietary programs. If your program is a subroutine library, you may consider it more useful to permit linking proprietary applications with the library. If this is what you want to do, use the GNU Library General Public License instead of this License.

GNU LESSER GENERAL PUBLIC LICENSE Version 2.1, February 1999

Copyright (C) 1991, 1999 Free Software Foundation, Inc.

59 Temple Place, Suite 330, Boston, MA 02111-1307 USA

Everyone is permitted to copy and distribute verbatim copies of this license document, but changing it is not allowed.

[This is the first released version of the Lesser GPL. It also counts as the successor of the GNU Library Public License, version 2, hence the version number 2.1.]

Preamble

The licenses for most software are designed to take away your freedom to share and change it. By contrast, the GNU General Public Licenses are intended to guarantee your freedom to share and change free software--to make sure the software is free for all its users.

This license, the Lesser General Public License, applies to some specially designated software packages--typically libraries--of the Free Software Foundation and other authors who decide to use it. You can use it too, but we suggest you first think carefully about whether this license or the ordinary General Public License is the better strategy to use in any particular case, based on the explanations below.

When we speak of free software, we are referring to freedom of use, not price. Our General Public Licenses are designed to make sure that you have the freedom to distribute copies of free software (and charge for this service if you wish); that you receive source code or can get it if you want it; that you can change the software and use pieces of it in new free programs; and that you are informed that you can do these things.

To protect your rights, we need to make restrictions that forbid distributors to deny you these rights or to ask you to surrender these rights. These restrictions translate to certain responsibilities for you if you distribute copies of the library or if you modify it.

For example, if you distribute copies of the library, whether gratis or for a fee, you must give the recipients all the rights that we gave you. You must make sure that they, too, receive or can get the source code. If you link other code with the library, you must provide complete object files to the recipients, so that they can relink them with the library after making changes to the library and recompiling it. And you must show them these terms so they know their rights.

We protect your rights with a two-step method: (1) we copyright the library, and (2) we offer you this license, which gives you legal permission to copy, distribute and/or modify the library.

To protect each distributor, we want to make it very clear that there is no warranty for the free library. Also, if the library is modified by someone else and passed on, the recipients should know that what they have is not the original version, so that the original author's reputation will not be affected by problems that might be introduced by others.

Finally, software patents pose a constant threat to the existence of any free program. We wish to make sure that a company cannot effectively restrict the users of a free program by obtaining a restrictive license from a patent holder. Therefore, we insist that any patent license obtained for a version of the library must be consistent with the full freedom of use specified in this license.

Most GNU software, including some libraries, is covered by the ordinary GNU General Public License. This license, the GNU Lesser General Public License, applies to certain designated libraries, and is quite different from the ordinary General Public License. We use this license for certain libraries in order to permit linking those libraries into non-free programs.

When a program is linked with a library, whether statically or using a shared library, the combination of the two is legally speaking a combined work, a derivative of the original library. The ordinary General Public License therefore permits such linking only if the entire combination fits its criteria of freedom. The Lesser General Public License permits more lax criteria for linking other code with the library.

We call this license the "Lesser" General Public License because it does Less to protect the user's freedom than the ordinary General Public License. It also provides other free software developers Less of an advantage over competing non-free programs. These disadvantages are the reason we use the ordinary General Public License for many libraries. However, the Lesser license provides advantages in certain special circumstances.

For example, on rare occasions, there may be a special need to encourage the widest possible use of a certain library, so that it becomes a de-facto standard. To achieve this, non-free programs must be allowed to use the library. A more frequent case is that a free library does the same job as widely used non-free libraries. In this case, there is little to gain by limiting the free library to free software only, so we use the Lesser General Public License.

In other cases, permission to use a particular library in non-free programs enables a greater number of people to use a large body of free software. For example, permission to use the GNU C Library in non-free programs enables many more people to use the whole GNU operating system, as well as its variant, the GNU/Linux operating system.

Although the Lesser General Public License is Less protective of the users' freedom, it does ensure that the user of a program that is linked with the Library has the freedom and the wherewithal to run that program using a modified version of the Library.

The precise terms and conditions for copying, distribution and modification follow. Pay close attention to the difference between a "work based on the library" and a "work that uses the library". The former contains code derived from the library, whereas the latter must be combined with the library in order to run.

GNU LESSER GENERAL PUBLIC LICENSE
TERMS AND CONDITIONS FOR COPYING, DISTRIBUTION AND MODIFICATION

0. This License Agreement applies to any software library or other program which contains a notice placed by the copyright holder or other authorized party saying it may be distributed under the terms of this Lesser General Public License (also called "this License").
Each licensee is addressed as "you".

A "library" means a collection of software functions and/or data prepared so as to be conveniently linked with application programs (which use some of those functions and data) to form executables.

The "Library", below, refers to any such software library or work which has been distributed under these terms. A "work based on the Library" means either the Library or any derivative work under copyright law: that is to say, a work containing the Library or a portion of it, either verbatim or with modifications and/or translated straightforwardly into another language. (Hereinafter, translation is included without limitation in the term "modification".)

"Source code" for a work means the preferred form of the work for making modifications to it. For a library, complete source code means all the source code for all modules it contains, plus any associated interface definition files, plus the scripts used to control compilation and installation of the library.

Activities other than copying, distribution and modification are not covered by this License; they are outside its scope. The act of running a program using the Library is not restricted, and output from such a program is covered only if its contents constitute a work based on the Library (independent of the use of the Library in a tool for writing it). Whether that is true depends on what the Library does and what the program that uses the Library does.

1. You may copy and distribute verbatim copies of the Library's complete source code as you receive it, in any medium, provided that you conspicuously and appropriately publish on each copy an appropriate copyright notice and disclaimer of warranty; keep intact all the notices that refer to this License and to the absence of any warranty; and distribute a copy of this License along with the Library.

You may charge a fee for the physical act of transferring a copy, and you may at your option offer warranty protection in exchange for a fee.

2. You may modify your copy or copies of the Library or any portion of it, thus forming a work based on the Library, and copy and distribute such modifications or work under the terms of Section 1 above, provided that you also meet all of these conditions:

- a) The modified work must itself be a software library.
- b) You must cause the files modified to carry prominent notices stating that you changed the files and the date of any change.
- c) You must cause the whole of the work to be licensed at no charge to all third parties under the terms of this License.
- d) If a facility in the modified Library refers to a function or a table of data to be supplied by an application program that uses the facility, other than as an argument passed when the facility is invoked, then you must make a good faith effort to ensure that, in the event an application does not supply such function or table, the facility still operates, and performs whatever part of its purpose remains meaningful.

(For example, a function in a library to compute square roots has a purpose that is entirely well-defined independent of the application. Therefore, Subsection 2d requires that any application-supplied function or table used by this function must be optional: if the application does not supply it, the square root function must still compute square roots.)

These requirements apply to the modified work as a whole. If identifiable sections of that work are not derived from the Library, and can be reasonably considered independent and separate works in themselves, then this License, and its terms, do not apply to those sections when you distribute them as separate works. But when you distribute the same sections as part of a whole which is a work based on the Library, the distribution of the whole must be on the terms of this License, whose permissions for other licensees extend to the entire whole, and thus to each and every part regardless of who wrote it.

Thus, it is not the intent of this section to claim rights or contest your rights to work written entirely by you; rather, the intent is to exercise the right to control the distribution of derivative or collective works based on the Library.

In addition, mere aggregation of another work not based on the Library with the Library (or with a work based on the Library) on a volume of a storage or distribution medium does not bring the other work under the scope of this License.

3. You may opt to apply the terms of the ordinary GNU General Public License instead of this License to a given copy of the Library. To do this, you must alter all the notices that refer to this License, so that they refer to the ordinary GNU General Public License, version 2, instead of to this License. (If a newer version than version 2 of the ordinary GNU General Public License has appeared, then you can specify that version instead if you wish.) Do not make any other change in these notices.

Once this change is made in a given copy, it is irreversible for that copy, so the ordinary GNU General Public License applies to all subsequent copies and derivative works made from that copy.

This option is useful when you wish to copy part of the code of the Library into a program that is not a library.

4. You may copy and distribute the Library (or a portion or derivative of it, under Section 2) in object code or executable form under the terms of Sections 1 and 2 above provided that you accompany it with the complete corresponding machine-readable source code, which must be distributed under the terms of Sections 1 and 2 above on a medium customarily used for software interchange.

If distribution of object code is made by offering access to copy from a designated place, then offering equivalent access to copy the source code from the same place satisfies the requirement to distribute the source code, even though third parties are not compelled to copy the source along with the object code.

5. A program that contains no derivative of any portion of the Library, but is designed to work with the Library by being compiled or linked with it, is called a "work that uses the Library". Such a work, in isolation, is not a derivative work of the Library, and therefore falls outside the scope of this License.

However, linking a "work that uses the Library" with the Library creates an executable that is a derivative of the Library (because it contains portions of the Library), rather than a "work that uses the library". The executable is therefore covered by this License. Section 6 states terms for distribution of such executables.

When a "work that uses the Library" uses material from a header file that is part of the Library, the object code for the work may be a derivative work of the Library even though the source code is not. Whether this is true is especially significant if the work can be linked without the Library, or if the work is itself a library. The threshold for this to be true is not precisely defined by law.

If such an object file uses only numerical parameters, data structure layouts and accessors, and small macros and small inline functions (ten lines or less in length), then the use of the object file is unrestricted, regardless of whether it is legally a derivative work. (Executables containing this object code plus portions of the Library will still fall under Section 6.)

Otherwise, if the work is a derivative of the Library, you may distribute the object code for the work under the terms of Section 6.

Any executables containing that work also fall under Section 6, whether or not they are linked directly with the Library itself.

6. As an exception to the Sections above, you may also combine or link a "work that uses the Library" with the Library to produce a work containing portions of the Library, and distribute that work under terms of your choice, provided that the terms permit modification of the work for the customer's own use and reverse engineering for debugging such modifications.

You must give prominent notice with each copy of the work that the Library is used in it and that the Library and its use are covered by this License. You must supply a copy of this License. If the work during execution displays copyright notices, you must include the copyright notice for the Library among them, as well as a reference directing the user to the copy of this License. Also, you must do one of these things:

- a) Accompany the work with the complete corresponding machine-readable source code for the Library including whatever changes were used in the work (which must be distributed under Sections 1 and 2 above); and, if the work is an executable linked with the Library, with the complete machine-readable "work that uses the Library", as object code and/or source code, so that the user can modify the Library and then relink to produce a modified executable containing the modified Library. (It is understood that the user who changes the contents of definitions files in the Library will not necessarily be able to recompile the application to use the modified definitions.)
- b) Use a suitable shared library mechanism for linking with the Library. A suitable mechanism is one that (1) uses at run time a copy of the library already present on the user's computer system, rather than copying library functions into the executable, and (2) will operate properly with a modified version of the library, if the user installs one, as long as the modified version is interface-compatible with the version that the work was made with.
- c) Accompany the work with a written offer, valid for at least three years, to give the same user the materials specified in Subsection 6a, above, for a charge no more than the cost of performing this distribution.

- d) If distribution of the work is made by offering access to copy from a designated place, offer equivalent access to copy the above specified materials from the same place.
- e) Verify that the user has already received a copy of these materials or that you have already sent this user a copy.

For an executable, the required form of the "work that uses the Library" must include any data and utility programs needed for reproducing the executable from it. However, as a special exception, the materials to be distributed need not include anything that is normally distributed (in either source or binary form) with the major components (compiler, kernel, and so on) of the operating system on which the executable runs, unless that component itself accompanies the executable.

It may happen that this requirement contradicts the license restrictions of other proprietary libraries that do not normally accompany the operating system. Such a contradiction means you cannot use both them and the Library together in an executable that you distribute.

7. You may place library facilities that are a work based on the Library side-by-side in a single library together with other library facilities not covered by this License, and distribute such a combined library, provided that the separate distribution of the work based on the Library and of the other library facilities is otherwise permitted, and provided that you do these two things:

- a) Accompany the combined library with a copy of the same work based on the Library, uncombined with any other library facilities. This must be distributed under the terms of the Sections above.
- b) Give prominent notice with the combined library of the fact that part of it is a work based on the Library, and explaining where to find the accompanying uncombined form of the same work.

8. You may not copy, modify, sublicense, link with, or distribute the Library except as expressly provided under this License. Any attempt otherwise to copy, modify, sublicense, link with, or distribute the Library is void, and will automatically terminate your rights under this License. However, parties who have received copies, or rights, from you under this License will not have their licenses terminated so long as such parties remain in full compliance.

9. You are not required to accept this License, since you have not signed it. However, nothing else grants you permission to modify or distribute the Library or its derivative works. These actions are prohibited by law if you do not accept this License. Therefore, by modifying or distributing the Library (or any work based on the Library), you indicate your acceptance of this License to do so, and all its terms and conditions for copying, distributing or modifying the Library or works based on it.
10. Each time you redistribute the Library (or any work based on the Library), the recipient automatically receives a license from the original licensor to copy, distribute, link with or modify the Library subject to these terms and conditions. You may not impose any further restrictions on the recipients' exercise of the rights granted herein.
You are not responsible for enforcing compliance by third parties with this License.
11. If, as a consequence of a court judgment or allegation of patent infringement or for any other reason (not limited to patent issues), conditions are imposed on you (whether by court order, agreement or otherwise) that contradict the conditions of this License, they do not excuse you from the conditions of this License. If you cannot distribute so as to satisfy simultaneously your obligations under this License and any other pertinent obligations, then as a consequence you may not distribute the Library at all. For example, if a patent license would not permit royalty-free redistribution of the Library by all those who receive copies directly or indirectly through you, then the only way you could satisfy both it and this License would be to refrain entirely from distribution of the Library.

If any portion of this section is held invalid or unenforceable under any particular circumstance, the balance of the section is intended to apply, and the section as a whole is intended to apply in other circumstances.

It is not the purpose of this section to induce you to infringe any patents or other property right claims or to contest validity of any such claims; this section has the sole purpose of protecting the integrity of the free software distribution system which is implemented by public license practices. Many people have made generous contributions to the wide range of software distributed through that system in reliance on consistent application of that system; it is up to the author/donor to decide if he or she is willing to distribute software through any other system and a licensee cannot impose that choice.

This section is intended to make thoroughly clear what is believed to be a consequence of the rest of this License.

12. If the distribution and/or use of the Library is restricted in certain countries either by patents or by copyrighted interfaces, the original copyright holder who places the Library under this License may add an explicit geographical distribution limitation excluding those countries, so that distribution is permitted only in or among countries not thus excluded. In such case, this License incorporates the limitation as if written in the body of this License.
13. The Free Software Foundation may publish revised and/or new versions of the Lesser General Public License from time to time.
Such new versions will be similar in spirit to the present version, but may differ in detail to address new problems or concerns.

Each version is given a distinguishing version number. If the Library specifies a version number of this License which applies to it and "any later version", you have the option of following the terms and conditions either of that version or of any later version published by the Free Software Foundation. If the Library does not specify a license version number, you may choose any version ever published by the Free Software Foundation.

14. If you wish to incorporate parts of the Library into other free programs whose distribution conditions are incompatible with these, write to the author to ask for permission. For software which is copyrighted by the Free Software Foundation, write to the Free Software Foundation; we sometimes make exceptions for this. Our decision will be guided by the two goals of preserving the free status of all derivatives of our free software and of promoting the sharing and reuse of software generally.

NO WARRANTY

15. BECAUSE THE LIBRARY IS LICENSED FREE OF CHARGE, THERE IS NO WARRANTY FOR THE LIBRARY, TO THE EXTENT PERMITTED BY APPLICABLE LAW.
EXCEPT WHEN OTHERWISE STATED IN WRITING THE COPYRIGHT HOLDERS AND/OR OTHER PARTIES PROVIDE THE LIBRARY "AS IS" WITHOUT WARRANTY OF ANY KIND, EITHER EXPRESSED OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE ENTIRE RISK AS TO THE QUALITY AND PERFORMANCE OF THE LIBRARY IS WITH YOU. SHOULD THE LIBRARY PROVE DEFECTIVE, YOU ASSUME THE COST OF ALL NECESSARY SERVICING, REPAIR OR CORRECTION.

16. IN NO EVENT UNLESS REQUIRED BY APPLICABLE LAW OR AGREED TO IN WRITING WILL ANY COPYRIGHT HOLDER, OR ANY OTHER PARTY WHO MAY MODIFY AND/OR REDISTRIBUTE THE LIBRARY AS PERMITTED ABOVE, BE LIABLE TO YOU FOR DAMAGES, INCLUDING ANY GENERAL, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES ARISING OUT OF THE USE OR INABILITY TO USE THE LIBRARY (INCLUDING BUT NOT LIMITED TO LOSS OF DATA OR DATA BEING RENDERED INACCURATE OR LOSSES SUSTAINED BY YOU OR THIRD PARTIES OR A FAILURE OF THE LIBRARY TO OPERATE WITH ANY OTHER SOFTWARE), EVEN IF SUCH HOLDER OR OTHER PARTY HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

END OF TERMS AND CONDITIONS

How to Apply These Terms to Your New Libraries

If you develop a new library, and you want it to be of the greatest possible use to the public, we recommend making it free software that everyone can redistribute and change. You can do so by permitting redistribution under these terms (or, alternatively, under the terms of the ordinary General Public License).

To apply these terms, attach the following notices to the library. It is safest to attach them to the start of each source file to most effectively convey the exclusion of warranty; and each file should have at least the "copyright" line and a pointer to where the full notice is found.

<one line to give the library's name and a brief idea of what it does.>
Copyright (C) <year> <name of author>

This library is free software; you can redistribute it and/or modify it under the terms of the GNU Lesser General Public License as published by the Free Software Foundation; either version 2.1 of the License, or (at your option) any later version.

This library is distributed in the hope that it will be useful, but WITHOUT ANY WARRANTY; without even the implied warranty of MERCHANTABILITY or FITNESS FOR A PARTICULAR PURPOSE. See the GNU Lesser General Public License for more details.

You should have received a copy of the GNU Lesser General Public License along with this library; if not, write to the Free Software Foundation, Inc., 59 Temple Place, Suite 330, Boston, MA 02111-1307 USA

Also add information on how to contact you by electronic and paper mail.

You should also get your employer (if you work as a programmer) or your school, if any, to sign a "copyright disclaimer" for the library, if necessary. Here is a sample; alter the names:

Yoyodyne, Inc., hereby disclaims all copyright interest in the library `Frob' (a library for tweaking knobs) written by James Random Hacker.

<signature of Ty Coon>, 1 April 1990
Ty Coon, President of Vice

That's all there is to it!

<MLA>

Copyright 2016 Microchip Technology Inc. (www.microchip.com)
Copyright 2017 Fuji Xerox Advanced Technology Co., Ltd

Licensed under the Apache License, Version 2.0 (the "License");
you may not use this file except in compliance with the License.
You may obtain a copy of the License at

https://github.com/MicrochipTech/mla_usb

Unless required by applicable law or agreed to in writing, software
distributed under the License is distributed on an "AS IS" BASIS,
WITHOUT WARRANTIES OR CONDITIONS OF ANY KIND, either express or
implied.

See the License for the specific language governing permissions and
limitations under the License.

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