FUJIFILM COVID-19 Ag Test

[Introduction]

The novel coronavirus (SARS-CoV-2) is the virus that is the cause of the 2019 coronavirus disease (COVID-19). It was first reported in Wuhan (Hubei Province, China) in December 2019. The World Health Organization (WHO) announced COVID-19 as a pandemic in March 2020. Infections rapidly increased around the world, and there is now a demand for a quick and simple way to test for SARS-CoV-2 antigens.

This product differs from nucleic acid detection tests in that it does not require any special equipment, so it allows for quick detection of SARS-CoV-2 antigens. This product is considered effective in assisting in the diagnosis of SARS-CoV-2 infection.

[Intended use]

Detection of SARS-CoV-2 antigens (target: Nucleocapsid protein) in Nasopharyngeal swabs (to assist in the diagnosis of SARS-CoV-2 infection)

For healthcare professionals (doctors, nurses and laboratory technicians, etc.) use only.

[Principle of the measurement]

The "FUJIFILM COVID-19 Ag Test" is a SARS-CoV-2 antigen detection agent which uses a nasopharyngeal swab, and is based on the principles of immunochromatographic assays.

The Mouse monoclonal anti-SARS-CoV-2 antibody exists on the membrane filter in the Test Cartridge in the solid phase. Additionally, there is a signaling particle area coated with colloidal gold conjugated to mouse monoclonal anti-SARS-CoV-2 antibody.

SARS-CoV-2 antigens in the specimen react with the colloidal gold conjugated to mouse monoclonal anti-SARS-CoV-2 antibody coating the signaling particle area to form an immune complex. This complex, which migrates as per principles of immunochromatography, is bound to the solid-phase mouse monoclonal anti-SARS-CoV-2 antibodies, forming a three-part complex sandwiching the SARS-CoV-2 antigen. Next, a reductant solution (ammonium iron sulfate) and a sensitization solution (silver nitrate) are added to sensitize the signal particles, resulting in the appearance of a line. This line can be used to indicate the presence of SARS-CoV-2 antigens in the specimen.

Additionally, the colloidal gold conjugated to mouse monoclonal anti-SARS-CoV-2 antibody is immobilized by solid-phase anti-mouse immunoglobulin antibodies (rabbit) at the control line. The line will appear after subsequent sensitization. This indicates that the reaction on the test strip has occurred successfully.

[Composition of the kit]

- 1) Test Cartridge
 - Mouse monoclonal anti-SARS-CoV-2 antibody
 - Colloidal gold conjugated to mouse monoclonal anti-SARS-CoV-2 antibody
 - Ammonium iron sulfate
 - Silver nitrate
- 2) Extraction Reagent Solution R1
Buffer solution containing surfactant

[Storage and shelf life]

- Storage: Room temperature (1-30°C)
- Expiration date: 12 months (as indicated on the package)

[Warnings and precautions]**1. Important precautions**

- (1) A negative results indicated by this product does not preclude acute SARS-CoV-2 infection.
- (2) Take the necessary precautions against biohazards when collecting and handling specimens.
- (3) This product is for *in vitro* diagnostic use only. Do not use it for anything other than its intended use.
- (4) Make an overall diagnosis which factors in other test results and clinical symptoms.
- (5) The reliability of results obtained by methods other than those described in the package insert is not guaranteed.

2. Precautions for handling (prevention for danger)

- (1) All specimens carry a risk of infection, so be careful when handling them.
- (2) Do not collect specimens with a sterile swab that has been dipped in the Extraction Reagent Solution R1.
- (3) The Test Cartridge of this kit incorporates a small amount of acidic reductant solution and sensitization solution. Do not disassemble the Test Cartridge.

Warning

- (4) In the event that Extraction Reagent Solution R1, reductant, or sensitization solution accidentally gets in one's eyes or mouth, or comes in contact with the skin, take appropriate emergency measures such as rinsing the affected area well with water, and if necessary, consult a doctor for treatment.
- (5) The Test Cartridge membrane is made out of nitrocellulose. It is an extremely flammable material. Do not use the Test Cartridge near an open flame.

3. Precautions for handling

- (1) Use all reagents as they are provided, without modification.
- (2) When keeping this product in chilled storage, make sure to let it sit in the place it will be used for a sufficient amount of time to allow all of the reagents (Test Cartridge, Extraction Reagent Solution R1, sterile swabs, and filters) to reach 15-30°C. Proceed to open the package, and use immediately after opening.
- (3) Avoid freezing the Test Cartridge and Extraction Reagent Solution R1 during storage.
- (4) When removing the Test Cartridge and Extraction Reagent Solution R1 from the aluminum packaging, be careful not to cut your hand to it.
- (5) When removing the Test Cartridge from the aluminum packaging, do not squeeze the Test Cartridge with excessive force. If button 2 or button 3 is pressed, it will discharge liquid which may render the product unable to be used for testing.
- (6) Do not use a Test Cartridge and Extraction Reagent Solution R1 tube, or filter that has been dropped.
- (7) Let the cartridge sit on a horizontal surface during use. Do not turn the Test Cartridge upside down or lie it on its side during use.
- (8) Prepare and test samples as quickly as possible after specimen collection, in accordance with the preparation methods.

- (9) When opening the aluminum seal on the Extraction Reagent Solution R1 tube, be careful of adhesion and scattering of liquid.
- (10) To collect a specimen, use your fingers to grasp the cotton swab, and move the tip of the cotton swab around the inside wall of the Extraction Reagent Solution R1 tube while applying pressure. After circling the sterile swab around the inside of the tube approximately 10 times, collect the specimen using the tip of the cotton swab.
- (11) When applying a sample to the Test Cartridge, make sure that the tip of the center filter is approximately 1 cm away from the sample well so that a drop is able to form, then drip a drop of the prescribed amount (approximately 24 µL) onto the Test Cartridge. If an amount other than the prescribed amount is used, the reaction may not occur successfully.
- (12) After dripping one drop of the prescribed amount (approximately 24 µL), do not drip any additional sample onto the Test Cartridge. Doing so may result in the reaction being unable to proceed successfully.
- (13) After dripping the sample onto the Test Cartridge, promptly press button 2 from above. Letting the test sit for an extended period of time without pressing button 2 after dripping the sample onto the Test Cartridge may cause the sample to desiccate, which may result in the reaction being unable to proceed successfully.
- (14) Pressing button 2 before dripping the sample onto the Test Cartridge may cause the reaction to be unable to proceed successfully.
- (15) Press button 3 once the "Go Next" mark completely turns from dark green to orange. The reaction will not occur normally if button 3 is pressed before the notification mark changes to orange. If the Test Cartridge is left to sit for an extended period of time after the mark turns orange, the reductant solution may desiccate, which may cause the reaction to be unable to proceed successfully.
- (16) Do not wait for an extended period of time after pressing button 3 to read the test result. Lines may become hard to see due to the effects of desiccation or light.

4. Precautions for use

- (1) Avoid freezing the reagents during storage. Do not use frozen reagents, because their quality may change and prevent correct results from being obtained.
- (2) Do not use expired reagents.
- (3) Use the dedicated Extraction Reagent Solution R1 included in the kit. Do not use an extraction solution provided with another kit.
- (4) When dripping the sample onto the sample well of the Test Cartridge, do not touch it directly with your hands.
- (5) Use the sterile swabs provided with the kit.
- (6) Be careful of the sample flying out of the tube when removing the sterile swab after preparing a sample.
- (7) If an excess amount of specimen is collected, or if the specimen is very viscous, the filter may get clogged, and be unable to dispense a suitable amount of sample. In such a case, collect a new specimen, and conduct the test using it.
- (8) Do not use the reagents, sterile swabs or filters that come with this product for any purpose other than this test.
- (9) Each Test Cartridge, Extraction Reagent Solution R1 tube, filter, and sterile swab is single use, and only to be used for a single test.
- (10) When preparing a sample, extract the specimen from the surface of the cotton swab by rubbing the tip of the sterile swab against the internal wall of the Extraction Reagent Solution R1 tube. When doing this, if you press too firmly on the external wall of the tube, the cotton swab by break, and the tip may come off. This may prevent the sample from being dispensed properly, so be careful.
- (11) When testing multiple specimens at the same time, label the surface of the cartridges with a marker with the name of the patient.

5. Cautions for disposal

- (1) Process and/or dispose of all specimens and used Test Cartridges, as well as sterile swabs and Extraction Reagent Solution R1 used for testing in accordance with your facility's regulations for infectious waste.
- (2) If autoclaving waste, autoclave it for 20 minutes or more at 121°C
- (3) Do not treat the Test Cartridge with disinfectants such as sodium hypochlorite, or glutaraldehyde, etc.

[Specimen requirements]**1. Required tools and equipment**


- Sterile swab: Use for collecting the specimen. (Included)
- Filter: Use for specimen preparation. (Included)
- Tube rack: Use for specimen preparation. (Included)

2. Specimen collection**A) Prepare for specimen collection**

- (1) Sterile swab
Use the sterile swabs which come with the kit.
- (2) Extraction Reagent Solution R1: Use as-is, at 15-30°C.


B) Specimen collection procedure**For a nasopharyngeal swab**

Insert the sterile swab deep into the nasal cavity in line with the inferior nasal concha (picture a flat plane connecting the nostril to the ear hole), and collect a mucosal epithelium specimen by rubbing the swab back and forth a few times.


3. Details of Extraction Reagent Solution R1 tube

Hold the Extraction Reagent Solution R1 tube with your hand so that it does not fall over.


4. Sample preparation


- (1) Remove the aluminum seal from the Extraction Reagent Solution R1 tube.


- (2) Insert the tip of the sterile swab which was used to collect the specimen to the bottom of the Extraction Reagent Solution R1 tube. Hold the Extraction Reagent Solution R1 tube from the outside as if to squeeze the tip of the cotton swab, then rub the cotton swab around the tube approximately 10 times.

- (3) Grip the body of the bottle, and squeeze the solution from the sterile swab while removing it from the tube.


- (4) Affix the filter onto the tube, and lightly shake the tube a few times so as to sufficiently blend the solution. The sample is now ready for testing.


[Procedure]

1. Details of Test Cartridge


2. Operation procedure


- A) Reagent preparation
Use the Test Cartridge as-is, at 15-30°C.
- B) Measurement procedure
 - (1) Remove the Test Cartridge from the aluminum package. Discard the desiccant/deoxygenating agent within the package.
 - (2) When dripping a prepared sample from the Extraction Reagent Solution R1 tube onto the sample well of the Test Cartridge, apply 1 drop (approximately 24 µL).


- (3) After dripping the sample onto the Test Cartridge, confirm that the Extraction Reagent Solution R1 has soaked into the sample well completely, then promptly press button 2 from above to start the reaction. At this time, make sure to press button 2 down completely.

- (4) Allow the Test Cartridge to sit horizontally and at room temperature (15-30°C) until the "Go Next" mark completely turns from dark green to orange.
Note) It takes approximately 10 minutes for the indicator mark to completely turn from dark green to orange.

- (5) After visually confirming that the "Go Next" mark has completely turned from dark green to orange, press button 3 down from above to start the sensitization reaction. At this time, make sure to press button 3 down completely.


- (6) Observe the presence or absence of the black line of the results window to read the test result.
Note) The sensitization reaction will take approximately 1 minute after pressing button 3 until completion.

[Interpretation of results]


1. Positive

When both a black control line (C) and black test line (T) can be seen


2. Negative

When only a black control line (C) can be seen

3. Invalid

If the control line cannot be seen, there was likely an operational error, such as insufficient sample used, so read through the operating procedures once more, then re-do the test with a new Test Cartridge.

Cautions for interpretation

- 1) This product is used to assist in the diagnosis of SARS-CoV-2 infection.
- 2) If the viral load in the specimen is below the detection sensitivity of this product, a positive specimen may have a negative result.
- 3) When collecting a specimen, if the specimen is not handled appropriately, an accurate test result might not be obtained.
- 4) Additionally, factors in the specimen may cause a non-specific reaction which causes a negative specimen to have a positive result. For the final diagnosis, perform a comprehensive judgment by considering clinical symptoms and other test results.
- 5) When reading the test result, a black vertical line may appear over the control line in the results window. This line does not have an effect on the final result. Please judge the result based on the result interpretation protocol.
In addition, if you have trouble making out the presence of either the test line or the control line because of the vertical lines, it is recommended that you re-do the test with a new Test Cartridge.

[Limitations of the examination procedure]

- 1) Interfering substances
At the concentrations listed below, the following substances will not invalidate the results given by this product.
Hemoglobin from hemolysis (450 mg/dL), blood (2%)
- 2) Cross reactivity
Cross reactivity was not observed with inactivated influenza virus (influenzavirus H1N1, influenzavirus H3N2, influenzavirus B), respiratory syncytial virus nor recombinant coronavirus antigens (MERS-CoV, HCoV-229E, HCoV-OC43, HCoV-NL63, HCoV-HKU1, SARS-CoV). No cross reactivity studies using the actual active viruses were conducted.

[Performance characteristics]

1. Performance

- (1) Accuracy
 - Testing an in-house positive control sample¹⁾ results in a positive result.
 - Testing an in-house negative control sample²⁾ results in a negative result.
 - (2) Precision
 - Testing three in-house positive control samples¹⁾ at the same time yields all positive results.
 - Testing three in-house negative control samples²⁾ at the same time yields all negative results.
- Note 1): A sample made from SARS-CoV-2 recombinant protein diluted to 25 pg/mL in Extraction Reagent Solution R1
Note 2): Extraction Reagent Solution R1
- (3) Limit of Detection
1.1×10² TCID₅₀/mL

2. Clinical performance

Clinical performance was established using retrospectively collected 225 samples (nasopharyngeal swabs).
Samples were stored in viral transport media (2-3 ml) and frozen until use. Thus these samples are diluted compared with prospective clinical samples.
Summary of the performance of the FUJIFILM COVID-19 Ag Test compared to RT-PCR for nasopharyngeal swabs is as follows.

Sensitivity (PPA): 90.9% (110/121, 95% CI: 84.3-95.4%)

Specificity (NPA): 100.0% (104/104, 95% CI: 96.5-100.0%)

		RT-PCR		
		Positive	Negative	Total
FUJIFILM COVID-19 Ag Test	Positive	110	0	110
	Negative	11	104	115
	Total	121	104	225

Sensitivity according to the RT-PCR Ct value

		RT-PCR Ct value			
		< 26	< 28	< 30	< 35
FUJIFILM COVID-19 Ag Test	Positive	100	107	110	110
	Negative	0	3	6	11
	Total	100	110	116	121
	Sensitivity (PPA)	100.0%	97.3%	94.8%	90.9%
	95% CI	96.4-100.0%	92.2-99.4%	89.1-98.1%	84.3-95.4%

EXPLANATION OF TERMS:

Sensitivity or Positive Percent Agreement (PPA) =
True Positives / (True Positives + False Negatives)

Specificity or Negative Percent Agreement (NPA) =
True Negatives / (True Negatives + False Positives)

[Contents]

FUJIFILM COVID-19 Ag Test (10 tests)

- Test Cartridge 10 tests
- Extraction Reagent Solution R1 10×0.4 mL
- Sterile Swab 10 pieces
- Tube Rack (for Extraction Reagent Solution R1) 1 piece
- Filter (for Extraction Reagent Solution R1) 10 pieces

FUJIFILM Corporation

26-30, Nishiazabu 2-Chome, Minato-ku, Tokyo 106-8620, JAPAN

EC REP FUJIFILM Europe GmbH
Heesenstrasse 31, 40549 Düsseldorf, GERMANY

<http://www.fujifilm.com/products/medical/>

[Symbols]

	In vitro diagnostic medical devices
	Catalogue number
	Consult instructions for use
	Do not reuse
	Temperature limitation
	This way up
	Contains sufficient for <n> tests
	Lot number
	Use by
	Unique Device Identifier
	Manufacturer
	Authorised representative in the European Community
	GHS pictograms (Warning)
	GHS07 Causes serious eye irritation Category 2
	GHS09 Chronic hazards to the aquatic environment Category 2

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