

List of disinfectants that can be used against COVID-19 for clinical chemical analyzer product (medium level and above)

We have confirmed that the disinfectants introduced below do not affect the performance of the equipment. If other disinfectants are used, this will be the customer's responsibility.

Area to be disinfected	Clinical Chemistry Analyzer
Product exterior	<Usable disinfectant> ■Ethanol for disinfection (excluding the transparent sampler cover)
Product interior	■Follow the instruction manual to clean or wipe with water, lukewarm water or ethanol for disinfection.