

PRO6 Conditioner Plus & Repl. For Processing Fujichrome and Process E-6 Compatible Color Reversal Films



Catalog Number: 600008663 (50 G)

DESCRIPTION

This product is a Process E-6 pre-bleach and replenisher. The product is designed for use in all equipment processing color reversal films compatible with Process E-6 and is recommended for use with Pro6® Rinse Plus & Replenisher **ONLY**.

SAFE HANDLING INSTRUCTIONS

Please refer to Material Safety Data Sheets (MSDS) for specific chemical details. For emergencies, please contact CHEMTREC at 1-800-424-9300. You can obtain an MSDS by calling the Fujifilm Technical Hotline at 1-800-526-0851 Ext. 10.

HMIS INFORMATION

Hazardous Materials Identification System. Please refer to Section 3 in the MSDS for further information.

| | Concentrate | Replenisher | Working Tank |
|---------------------|-------------|-------------|--------------|
| Health | 1 | 1 | 1 |
| Flammability | 0 | 0 | 0 |
| Reactivity | 0 | 0 | 0 |
| Personal Protection | C | C | C |

SPECS & RECOMMENDATIONS

| | | |
|-----------------------------|----------------------------------|---------|
| Temperature | 24° – 39° C 75° – 103° F | |
| Processing Time* | 2'00" Spec (1'00"- 4'00" Spread) | |
| Starting Replenishment Rate | per ft ² | 100 ml |
| | per m ² | 1076 ml |
| | per roll 135-24 | 42 ml |
| | per roll 120 | 54.7 ml |
| | per roll 220 | 109 ml |
| | per 4x5 | 13.9 ml |

*Once the process time is set the time should remain constant.

SPECIFIC GRAVITY & pH

| | pH @ 25° C | SpG @ 25° C | SpG @ 38° C |
|------------------|-------------|---------------|---------------|
| Replenisher Tank | 6.23 ± 0.10 | 1.015 ± 0.003 | 1.011 ± 0.003 |
| Fresh Tank | 6.23 ± 0.10 | 1.015 ± 0.003 | 1.011 ± 0.003 |
| Seasoned Tank | 8.70 ± 1.25 | 1.025 ± 0.006 | 1.022 ± 0.006 |

MIXING INSTRUCTIONS – REPLENISHER

| R | Water 68° – 86° F (20° – 30° C) | Concentrate |
|------------|---------------------------------------|-----------------|
| 50 Gallons | 45.0 G | 5.0 G 1 Cube |
| 1 Liter | 900 ml | 100 ml |
| 1 Gallon | 115.2 oz | 12.8 oz |

To Mix 50 Gallons of Replenisher:

1. Start with 45 gallons of warm water.
2. Add 1 cube concentrate. Rinse cube with small volume of water and add to mix tank. Mix 2 to 3 minutes or until uniform.

Do Not Over Mix

MIXING INSTRUCTIONS – WORKING TANK

| W | Water 68° – 86° F (20° – 30° C) | Concentrate |
|-------------|---------------------------------------|-------------|
| 1 Liter | 900 ml | 100 ml |
| 1 Gallon | 115.2 oz | 12.8 oz |
| Tank Volume | | |
| YOUR TANK | | |

To calculate the amount of water and conditioner concentrate that you will need, multiply your tank volume by each of the numbers listed, in either liters or gallons. For example, if you have a 10 liter tank, 10x900 ml = 9 liter of water and 10x100 ml = 1 liters of conditioner plus concentrate.

1000 ml = 1 liter 128 oz = 1 gallon

FOR BLENDER SETTINGS AND ADDITIONAL REPLENISHMENT RATES, PLEASE REFER TO BULLETIN PRO6 PROCESS SPECIFICATIONS

SHELF LIFE

| | Concentrate | Mixed Replenisher | Working Tank |
|-------|------------------------------|---------------------|---------------------|
| Time | 24 Months = U 4 Weeks = O | 2 Weeks | 2 Weeks |
| Color | Clear, Colorless | Clear, Colorless | Clear, Colorless |
| Odor | Slight Sulfur | No Odor | No Odor |

U = Unopened, O = Opened

The data contained in this photographic bulletin is believed to be true and accurate, but is offered solely for your consideration, investigation, and verification. Nothing herein shall be construed to be a warranty or guarantee by either FUJIFILM Hunt Chemicals, Inc. or any of its affiliates, and all such warranties, implied or otherwise, including any implied warranty of merchantability, are hereby expressly disclaimed. You, of course, are fully responsible for any use and / or domestic or foreign sale of the product(s) described, and nothing in this photographic bulletin shall be construed to constitute permission or a recommendation to practice any invention covered by patent or patent application or know-how owned by FUJIFILM, its affiliates, or others.