

PRO6

First Developer Starter

For Processing Fujichrome and Process E-6 Compatible Color Reversal Films



Catalog Number: 60001013 (1 G)

DESCRIPTION

This product is a Process E-6 first developer starter. The product is designed to mix a fresh working tank of Pro6® First Developer.

SAFE HANDLING INSTRUCTIONS

Please refer to Material Safety Data Sheets (MSDS) for specific chemical details. For emergencies, please contact CHEMTREC at 1-800-424-9300. You can obtain an MSDS by calling the Fujifilm Technical Hotline at 1-800-526-0851 Ext. 10.

HMIS INFORMATION

Hazardous Materials Identification System. Please refer to Section 3 in the MSDS for further information.

	Concentrate
Health	1
Flammability	0
Reactivity	0
Personal Protection	C

**WORKING TANK MIXING INSTRUCTIONS
for Pro6® First Developer**

W	First Developer Replenisher	Starter
1 Liter	1000 ml	5 ml
1 Gallon	128 oz	0.64 oz

To calculate the amount of first developer starter and developer replenisher that you will need, multiply your tank volume by each of the numbers listed, in either liters or gallons. For example, if you have a 10 liter tank, $10 \times 1000 \text{ ml} = 10 \text{ liters}$ of developer replenisher and $10 \times 5 \text{ ml} = 50 \text{ ml}$ of starter.

$$1000 \text{ ml} = 1 \text{ liter} \quad 128 \text{ oz} = 1 \text{ gallon}$$

SHELF LIFE

	Concentrate
Time	U = 36 Months, O = 6 Months
Color	Clear, Colorless
Odor	Slight Vinegar Odor

U = Unopened, O = Opened