

P

EC<sup>®</sup> -RA P1-S Developer Starter

For the Process of FUJICOLOR Fast Access FA Photographic Paper and Equivalents

Catalog Number: 600005405 (1 G)

**DESCRIPTION**

This product is a developer starter used for mixing a fresh working tank solution of EC<sup>®</sup> RA LR Developer, EC<sup>®</sup> RA 45 Developer, EC<sup>®</sup> RA 60 Developer, EC<sup>®</sup> RA 108 Developer, CP RA 95 Developer, and EC<sup>®</sup> LR 100 Developer.

**SAFE HANDLING INSTRUCTIONS**

Please refer to Material Safety Data Sheets (MSDS) for specific chemical details. For emergencies, please contact CHEMTREC at 1-800-424-9300. You can obtain an MSDS by calling the Fujifilm Technical Hotline at 1-800-526-0851 Ext. 10.

**HMIS INFORMATION**

Hazardous Materials Identification System. Please refer to Section 3 in the MSDS for further information.

	Concentrate
Health	1
Flammability	0
Reactivity	0
Personal Protection	C

**WORKING TANK MIX INSTRUCTIONS****using EC<sup>1</sup> RA 60 Developer Repl.**

W	Warm Water	Starter	EC <sup>1</sup> RA 60 Developer Replenisher
1 Liter	460 ml	140 ml	400 ml
1 Gallon	58.9 oz	17.9 oz	51.2 oz
Tank Volume			
YOUR TANK			

**WORKING TANK MIX INSTRUCTIONS****using EC<sup>1</sup> RA 45 Developer Repl.**

W	Warm Water	Starter	EC <sup>1</sup> RA VLR Developer Replenisher
1 Liter	570 ml	130 ml	300 ml
1 Gallon	73 oz	16.6 oz	38.4 oz
Tank Volume			
YOUR TANK			

**WORKING TANK MIX INSTRUCTIONS****using EC<sup>®</sup> RA 108 Developer Repl.**

W	Warm Water	Starter	EC <sup>®</sup> RA 108 Developer Replenisher
1 Liter	300 ml	85 ml	600 ml
1 Gallon	40.3 oz	10.9 oz	76.8 oz
Tank Volume			
YOUR TANK			

**WORKING TANK MIX INSTRUCTIONS****using CP RA 95 Developer Repl.**

W	Warm Water	Starter	CP RA 95 Developer Replenisher
1 Liter	350 ml	50 ml	600 ml
1 Gallon	44.8 oz	6.4 oz	76.8 oz
Tank Volume			
YOUR TANK			

**WORKING TANK MIX INSTRUCTIONS****using EC<sup>®</sup> LR 100 Developer Repl.**

W	Warm Water	Starter	EC <sup>®</sup> LR 100 Developer Replenisher
1 Liter	415 ml	85 ml	500 ml
1 Gallon	53.12 oz	10.88 oz	64 oz
Tank Volume			
YOUR TANK			

To figure out the amount of water, developer starter and developer replenisher that you will need, multiply your tank volume by each of the numbers listed, in either liters or gallons. For example, if you have a 10 liter tank, 10x415 ml = 4.15 liters of water, 10x85 ml = 850 ml of developer starter and 10x500 ml = 5.0 liters of developer replenisher.

1000 ml = 1 liter 128 oz = 1 gallon

**SHELF LIFE**

	Concentrate
Time	U = 36 Months O = 12 Months
Color	Clear, Colorless
Odor	No Odor

U = Unopened, O = Opened

The data contained in this photographic bulletin is believed to be true and accurate, but is offered solely for your consideration, investigation, and verification. Nothing herein shall be construed to be a warranty or guarantee by either FUJIFILM Hunt Chemicals, Inc. or any of its affiliates, and all such warranties, implied or otherwise, including any implied warranty of merchantability, are hereby expressly disclaimed. You, of course, are fully responsible for any use and / or domestic or foreign sale of the product(s) described, and nothing in this photographic bulletin shall be construed to constitute permission or a recommendation to practice any invention covered by patent or patent application or know-how owned by FUJIFILM, its affiliates, or others.