

DATA TAPE COURIER PRO CASE



**THERE IS NO BETTER
MEDIA CONTAINER OUT THERE.**

- Rubberized resin
- Removable media tray
- Shock resistant
- Oversized bumper guards
- Soft interior contact points
- Heavy duty compression locks
- Water resistant gasket seal
- Used extensively by the best off-site media vaults
- Holds up to 36 LTO cartridges (double width version)

Outperforms the rest!



FUJIFILM



(A) Data Tape Courier Pro Case with LTO Tray • (B) Double Width Data Tape Courier Pro Case with LTO Trays • (C) Data Tape Courier Pro Case with DLT® Tray
 Also available with: (D) DLT Tray (without cases) • (E) Soft LTO Foam Tray (without cases) • (F) Enterprise Tape Tray

DATA TAPE COURIER PRO CASE

Achieves the best balance between features, performance, compatibility and price

FEATURES:

- Removable multimedia tray system makes it easy to load and unload media.
- Individual media slots in tray allow you to safely carry a partially filled case.
- Accepted by major Disaster Recovery and Off-site vaults.
- Case and tray do not shed contaminant particles.
- Large ergonomic side grips and oversized bumpers.
- Label areas on front and sides.
- Additional trays available in various media formats.
- Holds up to 36 LTO cartridges (double width version).
- Pro Case Measurements:
 14.75" x 5.9" x 11.25" (regular width)
 27.00" x 5.9" x 11.25" (double width version)
- Cases can be stacked together.

PRODUCT LINE-UP:

Data Tape Courier Professional Case	Material Number	Data Tape Courier Professional Double Case	Material Number
Courier Pro Case - LTO W/Pcase	600006031	Courier Pro Double - LTO w/Pcase	600006572
Courier Pro Case - DLT W/Pcase	600006026	Courier Pro Double - LTO No Pcase	600006571
Courier Pro Case - HDD	600006028	Courier Pro Double - DLT w/Pcase	600006575
Data Tape Courier Pro OverPK (3 Pcs)	600006032	Courier Pro Double - DLT No Pcase	600006568
Courier Pro Case - LTO No Pcase	600006567	Courier Pro Double - Enterprise	600006574
Courier Pro Case - Enterprise	600006569	Courier Pro Double - HDD	600006573
Courier Pro Case - DLT No Pcase	600006570		

FUJIFILM

CODE: DTCP10-1 Printed on Recycled Paper

FUJIFILM Recording Media U.S.A., Inc., 200 Summit Lake Drive, Valhalla, NY 10595-1356
 For more information, please call 1-800-488-3854 or visit us at www.fujifilmusa.com/tapestorage

© 2010 FUJIFILM Recording Media U.S.A., Inc. Data Tape Courier is a trademark of FUJIFILM Recording Media U.S.A., Inc. DLT is a trademark of Quantum Corporation. Linear Tape-Open, LTO, the LTO Logo, Ultrium and the Ultrium Logo are registered trademarks of HP, IBM and Quantum in the US and other countries.

