

LTO Ultrium 4 Data Cartridge

Featuring FUJIFILM's proprietary NANOCUBIC Technology, the LTO Ultrium 4 data cartridge is FUJIFILM's latest addition to their successful LTO Ultrium family.¹ FUJIFILM's LTO Ultrium 4 tape format has been optimized for high-capacity and performance with outstanding reliability, in either stand-alone or automated environments. With an impressive 1600* GB capacity and 240* MB/sec transfer rate, the LTO Ultrium 4 data cartridge was designed to address your data storage needs by ensuring safe, accessible, and reliable data. Keep your data FUJIFILM safe.




Your logo or graphic here (optional)



FUJIFILM solely uses and is the exclusive authorized manufacturer of Tri-Optic bar code labels

Call FUJIFILM toll-free or visit www.fujibarcodes.ca for more details.



FUJIFILM's Universal Cleaning Cartridge is designed for use with all Ultrium 1, 2, 3, and 4 tape drives

Product Specifications

Product Code	Standard: 130301 Labelled: 131360
Durability	Avg. 1,000,000+ head passes
Capacity	800/1600* GB (Native/Compressed)
Transfer Rate	120/240* MB/sec (Native/Compressed)
Number of Tracks	1024
Dimensions (mm)	21.5 x 105.4 x 102 (W x H x L)

* Assumes 2:1 data compression
(1) LTO 4 data cartridge will be available Summer 2007
© 2007 FUJIFILM Canada Inc., The FUJIFILM, NANOCUBIC, and ATOMM logos are trademarks of FUJIFILM Canada Inc.
Linear Tape-Open, LTO, LTO logo, Ultrium and Ultrium logo are trademarks of Quantum, HP and IBM

