

PRODUCT LIST

2020-2021

FUJIFILM CN16 Process
MINILAB

September 2020
Version E01 / 09-20

Catn°	Product Name	Conc	Size to make
CN16S REPLENISHER CARTRIDGE			
967026	FJ CN16S NC1 Replenisher Kit x2		2 x 200 films
990077	FJ CN16S NC2 Replenisher Kit x2		2 x 1000 films
998484	FJ CN16S-ER NC1 Replenisher Kit x2		2 x 120 films
998492	FJ CN16S-ER NC2 Replenisher Kit x2		1 x 600 films
CN16S / START UP KIT			
975003	FJ CN16S N1 Developer Startup chemical		1x5,2L sz
975011	FJ CN16S N2 Bleach Startup chemical		1x3,6L sz
975029	FJ CN16S N3 Fixer Startup chemical		6x3,6L sz
990069	FJ CN16S N4 Stabilizer Startup chemical		24x1,9L sz
CN16L			
972497	<i>EnviroNeg Universal Developer Starter AC</i>	6x1L conc	
992271	EnviroNeg Developer Replenisher LR AC		2x5L sz
954743	EnviroNeg Developer Replenisher LR AC (INA Y2021)		2x10L sz
992966	FJ CN16L N2R-LQ Bleach Replenisher		4x2L sz
958694	FJ CN16L N3RS Fixer & Replenisher		4x2L sz
991596	EnviroNeg FF Superflo Stabilizer & Replenisher MB AC		2x10x10L sz
CN16Q			
954693	EnviroNeg Developer Replenisher Kit AC		2x10L sz
998518	EnviroNeg Developer Replenisher 60 AC LOW VOLUME		2x5L sz
993600	FJ CN16Q NQ2RS Bleach Replenisher		4x4L sz
931261	FJ CN16Q NQ3RS Bleach-Fix Replenisher		4x4L sz
991596	EnviroNeg FF Superflo Stabilizer & Replenisher MB AC		2x10x10L sz
CN16 SUPER RINSE			
999721	FJ Super Conditioner Tablets (FSC-100 20x 100 tablets)		20x(100x5L) sz
999632	FJ Super Conditioner Tablets (FSC-100 2x 100 tablets)		2x(100x5L) sz
999769	FHRSS-12 Cartridge		
999127	HIPURO High Purity Water	20L conc	
992156	HIPURO High Purity Water	60L conc	

Check PIS716 for FUJI HUNT replacement for **FUJI CN16L N1CR** Developer Replenisher

Check PIS715 for FUJI HUNT replacement for **FUJI CN16Q N1RS** Developer Replenisher

Check PIS717 for FUJI HUNT replacement for **FUJI CN16 N4-R** Safer Stabilizer