



PRODUCT LIST

2020-2021

RA4 Process
LARGE LAB

September 2020
Version E01 / 09-20

I. FUJI HUNT RA4 PROCESSING – REPLENISHER PRODUCTS

Catn°	Product Name	Conc	Size to make
DEVELOPER			
971657	EnviroPrint General Developer Starter AC *	6x1L conc	
999267	EnviroPrint CP48II Developer Replenisher	210L conc	
993790	EnviroPrint Developer Replenisher MP60 AC	60L conc	
994376	EnviroPrint Developer Replenisher MP60 AC	210L conc	
999557	EnviroPrint Developer Replenisher MP60 AC	640L conc	
999706	EnviroPrint Developer Replenisher MP73 AC	60L conc	
999691	EnviroPrint Developer Replenisher MP73 AC	640L conc	
990929	EnviroPrint Developer Replenisher MP108	60L conc	
991281	EnviroPrint Developer Replenisher MP160	10L conc	50L sz
991283	EnviroPrint Developer Replenisher MP160	20L conc	100L sz
999737	CPRA Digital Pro Developer Replenisher Part A	210L conc	
999740	CPRA Digital Pro Developer Replenisher Part B	200L conc	
BLEACH-FIX			
999275	EnviroPrint CP48II Bleach-Fix Replenisher Part A	210L conc	
999283	EnviroPrint CP48II Bleach-Fix Replenisher Part B	210L conc	
995845	EnviroPrint Bleach-Fix Replenisher VR AC Part A	210L conc	
995860	EnviroPrint Bleach-Fix Replenisher VR AC Part A	640L conc	
995886	EnviroPrint Bleach-Fix Replenisher VR AC Part A	1000L conc	
995852	EnviroPrint Bleach-Fix Replenisher VR AC Part B	210L conc	
995878	EnviroPrint Bleach-Fix Replenisher VR AC Part B	640L conc	
995894	EnviroPrint Bleach-Fix Replenisher VR AC Part B	1000L conc	
995951	CPRA Digital Pro Bleach-Fix Replenisher AC Part A	20L conc	
995969	CPRA Digital Pro Bleach-Fix Replenisher AC Part B	20L conc	
SUPER STABILIZER			
976050	EnviroPrint Super Stabilizer AC	60L conc	
976811	EnviroPrint Super Stabilizer AC	210L conc	
976647	EnviroPrint Super Stabilizer AC	640L conc	

* Check PIS 721 for replacement RA4 starters & TIS049 Starter substitution table for RA4 Developers

II. FUJI HUNT RA4 PROCESSING – REGENERATION PRODUCTS

Catn°	Product Name	Conc	Size to make
DEVELOPER REGENERATOR			
991802	EnviroPrint Developer. Regenerator HV 1R	60L conc	
991810	EnviroPrint Developer. Regenerator HV 1R	210L conc	
991828	EnviroPrint Developer. Regenerator HV 1R	640L conc	
972125	Multigen CD-60	60L conc	
972133	Multigen CD-60	210L conc	
972109	Multigen CD-60	640L conc	
991836	EnviroPrint Developer. Regenerator HV 3R	60L conc	
991844	EnviroPrint Developer. Regenerator HV 3R	210L conc	
991851	EnviroPrint Developer. Regenerator HV 3R	640L conc	
BLEACH-FIX REGENERATORS			
946443	EnviroPrint Electrolytic Bleach-Fix Regenerator Pt A	60L conc	
954867	EnviroPrint Electrolytic Bleach-Fix Regenerator Pt B	60L conc	
945220	EnviroPrint Electrolytic Bleach-Fix Regenerator Pt B	210L conc	
970020	EnviroPrint Electrolytic Bleach-Fix Regenerator Pt B	640L conc	
954875	EnviroPrint Electrolytic Bleach-Fix Regenerator Pt C	60L conc	
945238	EnviroPrint Electrolytic Bleach-Fix Regenerator Pt C	210L conc	
970046	EnviroPrint Electrolytic Bleach-Fix Regenerator Pt C	640L conc	
BLEACH REPLENISHER & REGENERATOR			
975763	EnviroPrint Bio Bleach Replenisher	640L conc	
FIXER REPLENISHER & REGENERATOR			
954123	Super Unilec Fixer	60L conc	
925420	Super Unilec Fixer	210L conc	
970327	Super Unilec Fixer	640L conc	
992552	Super Unilec Fixer	1000L conc	

III. MISCELLANEOUS PRODUCTS

Catn°	Product Name	Conc	Size to make
928127	Ultra Bleach-Fix Extender	10L conc	
942229	Acticide MV	5KG conc	
999668	ATS 58%	640L conc	
957571	FeAmmedta 50% w/w	60L conc	
970178	FeAmmedta 50% w/w	640L conc	
929547	Acetic Acid 60%	60L conc	
972869	ABS 70%	60L conc	
958306	ABS 70%	200L conc	
994202	Field Test Kit PE (48 x 250 ml)		