

PERSONALIZED STAINLESS STEEL

Whether you're sipping coffee at home or hydrating on the go, PPPG has got you covered with a personalized stainless steel accessory that will keep your beverage at the ideal temperature! Our stainless steel products are offered in a variety of sizes and styles, are insulated, and the appropriate drinkware items fit in most standard cup holders. Please refer to the most recent PPPG Price List for a listing of prices for the products within this sheet.



10 OZ. INSULATED STAINLESS STEEL MUG

- Photo is displayed in a single image wrapped around the mug
- Comfortable handle and splash-proof lid

PRODUCT NAME

10 oz. Insulated Stainless Steel Mug

BILL CODE

25604

PRODUCT CODE

PRGift;6656



12 OZ. INSULATED WINE TUMBLER

- Photo(s) can be displayed as single image or collage
- Features a finger slide splash-proof lid
- Great for glass-free areas

PRODUCT NAME

12 oz. Insulated Wine Tumbler

BILL CODE

25618

PRODUCT CODE

PRGift;6668



14 OZ. STAINLESS STEEL TUMBLER

- Photo(s) can be displayed as a single image or collage
- Comes with a finger-slide spill proof top

PRODUCT NAME

14 oz. Stainless Steel Tumbler

BILL CODE

4950

PRODUCT CODE

PRGift;5365



17 OZ. SLIM WHITE WATER BOTTLE

- Photo is displayed in a single image wrapped around the bottle
- Stainless steel and vacuum insulated
- Features a brushed aluminum cap with rubber seal

PRODUCT NAME

17 oz. Slim White Water Bottle

BILL CODE

24586

PRODUCT CODE

PRGift;6334



20 OZ. STAINLESS STEEL PHOTO WATER BOTTLE

- Photo is displayed in a single image wrapped around the bottle
- Comes with a spill-proof screw top and clip

PRODUCT NAME

20 oz. Stainless Steel Photo Water Bottle

BILL CODE

4927

PRODUCT CODE

PRGift;5249



PERSONALIZED DOUBLE-WALL TUMBLERS

- Photo(s) can be displayed as a single image or collage
- Available in 20 oz. & 30 oz.
- Stainless steel and comes with a plastic gasket seal top

PRODUCT NAME

20 oz. Personalized Double-Wall Tumbler
30 oz. Personalized Double-Wall Tumbler

BILL CODE

24843
24844

PRODUCT CODE

PRGift;5049
PRGift;5195



25 OZ. WHITE THERMOS

- Photo is displayed in a single image wrapped around the Thermos
- Comes with a pop up sip inner lid
- Outer lid is also stainless steel

PRODUCT NAME

25 oz. White Thermos

BILL CODE

25278

PRODUCT CODE

PRGift;6544

For more information on PPPG fulfillment, contact your FUJIFILM sales representative.