

FUJIFILM

X-T1 IR

FORENSIC IR KIT




PHOTO BY 2015 JULIO SOSA | X-T1 IR CAMERA AND XF18-55MM LENS AT 1/180 SEC F22, ISO 200
PHOTO TAKEN WITH 87C IR FILTER, #15 YELLOW FILTER, AND QUANTUM QFLASH T5DR

FUJIFILM X-T1 IR INFRARED DIGITAL CAMERA SYSTEM OFFERS LIGHT SPECTRUM RECORDED FROM APPROX. 380NM-1,000NM

ALLOWS DETAILED IMAGE CAPTURE OF BLOOD SPATTER, PLANTAR IMPRESSIONS, BONE FRAGMENTS, BRUISING, BITE MARKS, AND GSR (GUN SHOT RESIDUE)

HIGHLY EFFECTIVE TOOL FOR TATTOO IDENTIFICATION, DOCUMENT FORGERY AND NIGHTTIME SURVEILLANCE WORK

PERMITS LIVE, CONTINUOUS VIEWING OF UV/IR DETAILS BEFORE FINAL CAPTURE

WEATHER-RESISTANT STRUCTURE DESIGNED TO STAND UP TO THE RIGORS OF FIELD WORK.

IDEAL INVESTIGATIVE TOOL FOR LAW ENFORCEMENT, MEDICAL EXAMINERS, AND PUBLIC SAFETY AGENCIES

AVAILABLE IN EARLY 2016 TO APPROPRIATE BUYERS AT A LIMITED NUMBER OF AUTHORIZED FUJIFILM DEALERS ONLY. TO IDENTIFY AN APPROPRIATE AUTHORIZED DEALER, CONTACT FUJIFILM AT [DIGITALINFO@FUJIFILM.COM](mailto:DigitalInfo@Fujifilm.com)

X-T1 IR AVAILABLE IN FORENSIC KIT AND BODY ONLY OPTION.
X-T1 IR BODY ONLY OPTION REQUIRES SEPARATE X-SERIES LENS AND FILTER PURCHASE.



FUJIFILM and FUJINON are trademarks of FUJIFILM Corporation and its affiliates.
© 2015 FUJIFILM North America Corporation and its affiliates. All rights reserved.